

EUS

**Expect™ Slimline (SL) Endoscopic
Ultrasound Aspiration Needles**

➤ Expect Slimline (SL) Endoscopic Ultrasound
Aspiration Needle41

**Acquire™ Endoscopic Ultrasound
Fine Needle Biopsy (FNB) Device**

➤ Acquire Endoscopic Ultrasound Fine Needle Biopsy (FNB)
Device.....41

LumiCoil™ Platinum Fiducial Markers

➤ LumiCoil Platinum Fiducial Markers – Straight Configuration..... 42


➤ LumiCoil Platinum Fiducial Markers – Figure 8 Configuration 42

AXIOS™ System

➤ AXIOS Stent and Electrocautery Enhanced Delivery System ... 42

➤ AXIOS Stent and Delivery System 42

Expect™ Slimline (SL) Endoscopic Ultrasound Aspiration Needles



Expect Slimline (SL) Endoscopic Ultrasound Aspiration Needle*




Order Number	GTIN	Needle Size	Minimum Working Channel (mm)	Sheath Diameter (mm)	Packaging
M00555500	08714729848141	19 gauge (1.10mm)	2.8	1.83	Each
M00555510	08714729848158	22 gauge (0.72mm)	2.4	1.65	Each
M00555520	08714729848165	25 gauge (0.52mm)	2.4	1.52	Each
M00555530	08714729848172	19 gauge Flex (1.14mm)	2.8	1.73	Each
M00555501	08714729848714	19 gauge (1.10mm)	2.8	1.83	Box 5
M00555511	08714729848721	22 gauge (0.72mm)	2.4	1.65	Box 5
M00555521	08714729848738	25 gauge (0.52mm)	2.4	1.52	Box 5
M00555531	08714729848806	19 gauge Flex (1.14mm)	2.8	1.73	Box 5

*Packaging also includes: 20cc syringe and one-way stopcock

*Needle length: 0cm to 8cm, adjustable

*Working length: 137.5cm to 141.5cm, adjustable

Acquire™ Endoscopic Ultrasound Fine Needle Biopsy (FNB) Device



Acquire Endoscopic Ultrasound Fine Needle Biopsy (FNB) Device*




Order Number	GTIN	Needle Size	Minimum Working Channel (mm)	Sheath Diameter (mm)	Packaging
M00555540	08714729931782	22 gauge (0.72mm)	2.4	1.65	Please contact your Boston Scientific Sales Representative for ordering information
M00555560	08714729931805	25 gauge (0.52mm)	2.4	1.52	Please contact your Boston Scientific Sales Representative for ordering information
M00555580	08714729931829	19 gauge Flex (1.14mm)	2.8	1.73	Please contact your Boston Scientific Sales Representative for ordering information
M00555541	08714729931843	22 gauge (0.72mm)	2.4	1.65	Please contact your Boston Scientific Sales Representative for ordering information
M00555561	08714729931867	25 gauge (0.52mm)	2.4	1.52	Please contact your Boston Scientific Sales Representative for ordering information
M00555581	08714729931881	19 gauge Flex (1.14mm)	2.8	1.73	Please contact your Boston Scientific Sales Representative for ordering information

*Packaging also includes: 20cc syringe and one-way stopcock

*Needle length: 0cm to 8cm, adjustable

*Working length: 137.5cm to 141.5cm, adjustable

LumiCoil™ Platinum Fiducial Markers




LumiCoil Platinum Fiducial Markers



Order Number	GTIN	Description	Fiducial O.D. (in)	Overall Fiducial Length (mm)	Implanted Fiducial Length (mm)	Differentiator	Packaging
M00553831	08714729972877	Lumicoil Platinum Fiducial Marker	0.018	5	5	Straight	Box 3
M00553841	08714729972891	Lumicoil Platinum Fiducial Marker	0.018	10	5	Figure 8	Box 3

AXIOS™ System



AXIOS Stent and Electrocautery Enhanced Delivery System



Order Number	GTIN	Flange Diameter (mm)	Lumen Diameter (mm)	Saddle Length (mm)	Catheter O.D. (Fr)	Catheter Working Length (cm)	Catheter Total Length (cm)	Packaging
M00553640	08714729904588	21	10	10	10.8	138	144	Each
M00553650	08714729904595	24	15	10	10.8	138	144	Each
M00553660	08714729951179	29	20	10	10.8	138	144	Each

AXIOS Stent and Delivery System



Order Number	GTIN	Flange Diameter (mm)	Lumen Diameter (mm)	Saddle Length (mm)	Catheter O.D. (Fr)	Catheter Working Length (cm)	Catheter Total Length (cm)	Packaging
M00553740	08714729904625	21	10	10	10.8	138	144	Each
M00553750	08714729904632	24	15	10	10.8	138	144	Each