


Resection & Retrieval



Resection Devices

EMR for Upper GI

➤ Captivator™ EMR Device	75
--------------------------------	----

Stiff Snares

➤ Captivator COLD Single-Use Snares	75
➤ Captivator II Single-Use Snares	76
➤ Captivator Single-Use Snares	76

Specialty Snares

➤ Profile™ Single-Use Snares	77
➤ Single-Use Rotatable Snares	77

Snares

➤ Captiflex™ Single-Use Snares	78
➤ Sensation™ Single-Use Short Throw Snares	78

Active Cords

➤ Active Cords	79
----------------------	----

Submucosal Lifting

➤ ORISE™ Gel	79
--------------------	----

Tissue Retraction

➤ ORISE Tissue Retractor	79
➤ ORISE Instrument Guides	79

Retrieval Devices

Retrieval Device

➤ RescueNet™ Retrieval Device	80
-------------------------------------	----

Retrieval Grasping Forceps

➤ Rescue™ Rat Tooth Grasping Forceps	80
➤ Rescue Rat Tooth Grasping Forceps - Pediatric	80
➤ Rescue Alligator Grasping Forceps, Long	80
➤ Rescue Alligator Grasping Forceps	80
➤ Rescue Rat Tooth / Alligator Grasping Forceps	80


Polyp Traps

➤ TrapEase™ Single Chamber Polyp Trap	80
➤ TrapEase Quad Chamber Polyp Trap	80



Resection Devices

EMR for Upper GI




Captivator™ EMR Device



Order Number	GTIN	Description	Packaging
M00561600	08714729842675	Captivator EMR Device for Standard Gastrosopes	Box 1
M00561610	08714729842699	Captivator EMR Device for Large Gastrosopes	Box 1
M00561601	08714729842682	Captivator EMR Device for Standard Gastrosopes	Box 5
M00561611	08714729842705	Captivator EMR Device for Large Gastrosopes	Box 5

Stiff Snares



Captivator COLD Single-Use Snare




Order Number	GTIN	Description	Loop Width (mm)	Sheath O.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Packaging
M00561101	08714729971269	10mm Rounded – Stiff	10	2.4	2.8	240	Box 10
M00561102	08714729984979	10mm Rounded – Stiff	10	2.4	2.8	240	Box 20
M00561103	08714729984986	10mm Rounded – Stiff	10	2.4	2.8	240	Box 40



Resection Devices


Stiff Snares (Continued)



Captivator™ II Single-Use Snares



Order Number	GTIN	Description	Loop Width (mm)	Sheath O.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Packaging
M00561221	08714729855903	10mm Rounded – Stiff	10	2.4	2.8	240	Box 10
M00561222	08714729861300	10mm Rounded – Stiff	10	2.4	2.8	240	Box 20
M00561223	08714729861317	10mm Rounded – Stiff	10	2.4	2.8	240	Box 40
M00561231	08714729855910	15mm Rounded – Stiff	15	2.4	2.8	240	Box 10
M00561232	08714729861331	15mm Rounded – Stiff	15	2.4	2.8	240	Box 20
M00561233	08714729861348	15mm Rounded – Stiff	15	2.4	2.8	240	Box 40
M00561241	08714729855927	20mm Rounded – Stiff	20	2.4	2.8	240	Box 10
M00561242	08714729861362	20mm Rounded – Stiff	20	2.4	2.8	240	Box 20
M00561243	08714729861379	20mm Rounded – Stiff	20	2.4	2.8	240	Box 40
M00561191	08714729855934	25mm Rounded – Stiff	25	2.4	2.8	240	Box 10
M00561192	08714729861270	25mm Rounded – Stiff	25	2.4	2.8	240	Box 20
M00561193	08714729861287	25mm Rounded – Stiff	25	2.4	2.8	240	Box 40
M00561291	08714729645801	Extra Large Rounded – Stiff	33	2.4	2.8	240	Box 10



Captivator Single-Use Snares




Order Number	GTIN	Description	Loop Width (mm)	Sheath O.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Packaging
M00562301	08714729019312	Small Oval – Stiff	13	2.4	2.8	240	Box 10
M00562321	08714729019336	Medium Oval – Stiff	27	2.4	2.8	240	Box 10
M00561311	08714729019251	Large Oval – Medium Stiff	30	2.4	2.8	240	Box 10
M00562391	08714729019404	Large Oval (thin wire) – Flexible	30	2.4	2.8	240	Box 10
M00562451	08714729071068	Small Hexagonal – Stiff	13	2.4	2.8	240	Box 10
M00562341	08714729019350	Medium Hexagonal – Stiff	27	2.4	2.8	240	Box 10
M00562371	08714729019381	Medium Crescent – Stiff	27	2.4	2.8	240	Box 10



Resection Devices


Speciality Snares



Profile™ Single-Use Snares



Order Number	GTIN	Description	Loop Width (mm)	Sheath O.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Packaging
M00562531	08714729156628	Extra Small Oval – Flexible	11	1.9	2.0	240	Box 10
M00562551	08714729156659	Small Oval – Flexible	13	1.9	2.0	240	Box 10
M00562571	08714729156635	Medium Oval – Flexible	27	1.9	2.0	240	Box 10



Single-Use Rotatable Snares*



Order Number	GTIN	Description	Loop Width (mm)	Sheath O.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Packaging
M00561821	08714729283980	Small Oval – Medium Stiff	13	2.4	2.8	240	Box 10
M00561831	08714729283997	Medium Oval – Medium Stiff	20	2.4	2.8	240	Box 10

*With Star Catheter™ Technology.



Section Table of Contents

Prev Page ◀ 77 ▶ Next Page


Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Resection Devices


Snares



Captiflex™ Single-Use Snares



Order Number	GTIN	Description	Loop Width (mm)	Sheath O.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Packaging
M00562471	08714729019459	Extra Small Oval – Flexible	11	2.4	2.8	240	Box 10
M00562421	08714729019435	Small Oval – Flexible	13	2.4	2.8	240	Box 10
M00562422	08714729501657	Small Oval – Flexible	13	2.4	2.8	240	Box 40
M00562401	08714729019411	Medium Oval – Flexible	27	2.4	2.8	240	Box 10
M00562402	08714729501640	Medium Oval – Flexible	27	2.4	2.8	240	Box 40



Sensation™ Single-Use Short Throw Snares



Order Number	GTIN	Description	Loop Width (mm)	Sheath O.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Packaging
M00562691	08714729158158	Small Oval – Flexible	13	2.4	2.8	240	Box 5
M00562692	08714729268826	Small Oval – Flexible	13	2.4	2.8	240	Box 20
M00562693	08714729501671	Small Oval – Flexible	13	2.4	2.8	240	Box 40
M00560321	08714729283928	Small Oval – Medium Stiff	13	2.4	2.8	240	Box 5
M00562671	08714729158141	Medium Oval – Flexible	27	2.4	2.8	240	Box 5
M00562672	08714729268819	Medium Oval – Flexible	27	2.4	2.8	240	Box 20
M00562673	08714729501664	Medium Oval – Flexible	27	2.4	2.8	240	Box 40
M00560311	08714729283904	Medium Oval – Medium Stiff	27	2.4	2.8	240	Box 5
M00562651	08714729158110	Large Oval – Medium Stiff	30	2.4	2.8	240	Box 5
M00562652	08714729268802	Large Oval – Medium Stiff	30	2.4	2.8	240	Box 20
M00560331	08714729283935	Medium Crescent – Medium Stiff	27	2.4	2.8	240	Box 5



Section Table of Contents

Prev Page ◀ 78 ▶ Next Page

Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001


Resection Devices

Active Cords

Active Cords

Order Number	GTIN	Description	Packaging
M00561250	08714729019060	Active Cord, Fits Olympus™ Brand Generator	Each
M00561270	08714729053361	Active Cord, Fits Endostat™, Endostat II and Endostat III Generators, Bovie™, Valleylab™, Aspen Labs™, and ERBE™ Brand Generators	Each

Submucosal Lifting



ORISE™ Gel

Order Number	GTIN	Description	Packaging
M00519201	08714729974574	ORISE Gel Syringe Twin Pack (two 10mL syringes per tray)	Box 10
M00519210	08714729974581	ORISE Gel Syringe Kit Twin Pack with Interject™ Injection Therapy Needle Catheter 23ga Clear (two 10mL syringes per kit)	Box 1
M00519211	08714729974598	ORISE Gel Syringe Kit Twin Pack with Interject Injection Therapy Needle Catheter 23ga Clear (two 10mL syringes per kit)	Box 10

Tissue Retraction



ORISE Tissue Retractor

Order Number	GTIN	Description	Anatomy	Diameter (mm)	Cage Length (cm)	Working Length (cm)	# of channels	Scope Compatibility (mm)	Packaging
M00519510	08714729974604	ORISE Tissue Retractor 60L -OTR 60L	Left Colon	26	8	60	3	≤11.7mm outer diameter	Box 1


ORISE Instrument Guides

Order Number	GTIN	Description	Anatomy	Working Length (cm)	Tip Angles (degrees)	Tip Color	Instrument Compatibility	Packaging
M00519520	08714729974611	ORISE Instrument Guide 90B -OIG 90B	Left Colon	81	90	Blue	2.8	Box 1
M00519530	08714729974628	ORISE Instrument Guide 90G -OIG 90G	Left Colon	81	90	Gray	2.8	Box 1
M00519540	08714729974635	ORISE Instrument Guide 60B -OIG 60B	Left Colon	81	60	Blue	2.8	Box 1
M00519550	08714729974642	ORISE Instrument Guide 45B -OIG 45B	Left Colon	81	45	Blue	2.8	Box 1




Retrieval Devices

Retrieval Device



RescueNet™ Retrieval Device							
Order Number	GTIN	Description	Net Size (mm)	Sheath O.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Packaging
DGN-538-5	10840253107425	RescueNet Retrieval Device	30	2.5	2.8	230	Box 5

Retrieval Grasping Forceps



Rescue™ Retrieval Single-Use Graspers							
Order Number	GTIN	Description	Jaw Opening Width (mm)	Sheath O.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Packaging
DGR-276-5	00816849017334	Rescue Rat Tooth Grasping Forceps	7	2.4	2.8	230	Box 5
DGR-277-5	00816849017341	Rescue Rat Tooth Grasping Forceps - Pediatric	4.5	1.8	2.0	180	Box 5
DGA-275-5	00816849017310	Rescue Alligator Grasping Forceps, Long	9	2.4	2.8	230	Box 5
DGA-274-5	00816849017303	Rescue Alligator Grasping Forceps	8	2.4	2.8	230	Box 5
DGC-383-5	00816849017327	Rescue Rat Tooth / Alligator Grasping Forceps	8	2.4	2.8	230	Box 5

Polyp Traps



TrapEase™ Polyp Traps							
Order Number	GTIN	Description					
STE-297-50	00816849010182	TrapEase Single Chamber Polyp Trap					
STE-298-25	00816849010199	TrapEase Quad Chamber Polyp Trap					

