


Tissue Resection and Hemostasis Products

Boston Scientific
Advancing science for life™

Recommended Options For Endoscopic Mucosal Resection (EMR) Stiff Snares



STIFF SNARES

CAPTIVATOR™ SNARE


- **13mm** Small Hexagon – 6245
- **27mm** Medium Hexagon – 6234

CAPTIVATOR II SNARE

- **10mm** Rounded – 6122
- **15mm** Rounded – 6123
- **20mm** Rounded – 6124
- **25mm** Rounded – 6119
- **33mm** Extra Large Rounded – 6129

- Stiff braided wire design helps facilitate tactile feel and control
- Provides traction and limited tip deflection on flat lesions
- Rounded and Hexagonal Shapes
- Multiple sizes for large polyps and residual tissue

Recommended Options For COLD Polypectomy Stiff and Flexible Snares



STIFF SNARES


CAPTIVATOR II SNARE

- **10mm** Rounded – 6122
- **15mm** Rounded – 6123

CAPTIVATOR SNARE

- **13mm** Small Hexagon – 6245

- Now with smaller loop sizes for diminutive polypectomy – 10mm and 15mm



FLEXIBLE SNARES


CAPTIFLEX™ SNARE

- **11mm** Extra Small Oval – 6247
- **13mm** Small Oval – 6242

- Braided design for flexibility
- Extra small and small oval shaped snares


Tissue Resection Products for Polypectomy, Endoscopic Mucosal Resection (EMR), Injection, Closure, and Retrieval*

A combination of devices are available for use in conjunction with EMR procedures including a wide variety of snares for resection, sclerotherapy needles for injection, retrieval nets for removal of tissue and endoscopic clips for hemostasis and defect closure (less than 3cm).



Single-Use Rotatable Snare

Order Number	GTIN	Quantity	Descriptor	Loop Width Diameter (mm)	Sheath O.D. (mm)
M00561821	08714729283980	Box 10	Small Oval, Medium Stiff	13	2.8
M00561831	08714729283997	Box 10	Medium Oval , Medium Stiff	20	2.8



Captivator™ Snare

Order Number	GTIN	Quantity	Descriptor	Loop Width (mm)	Minimum Working Channel (mm)
M00562301	08714729019312	Box 10	Small Oval - Stiff	13	2.8
M00562321	08714729019312	Box 10	Medium Oval - Stiff	27	2.8
M00561311	08714729019251	Box 10	Large Oval - Medium Stiff	30	2.8
M00562391	08714729019404	Box 10	Large Oval (Thin-Wire) - Flexible	30	2.8
M00562451	08714729071068	Box 10	Small Hexagonal - Stiff	13	2.8
M00562341	08714729019350	Box 10	Medium Hexagonal - Stiff	27	2.8
M00562371	08714729019381	Box 10	Medium Crescent - Stiff	27	2.8



Captivator II Snare

Order Number	GTIN	Quantity	Descriptor	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)
M00561221	08714729855903	Box 10	10mm Rounded – Stiff	10	2.4	240	2.8
M00561222	08714729861300	Box 20	10mm Rounded – Stiff	10	2.4	240	2.8
M00561231	08714729855910	Box 10	15mm Rounded – Stiff	15	2.4	240	2.8
M00561232	08714729861331	Box 20	15mm Rounded – Stiff	15	2.4	240	2.8
M00561241	08714729855927	Box 10	20mm Rounded – Stiff	20	2.4	240	2.8
M00561242	08714729861362	Box 20	20mm Rounded – Stiff	20	2.4	240	2.8
M00561191	08714729855934	Box 10	25mm Rounded – Stiff	25	2.4	240	2.8
M00561192	08714729861270	Box 20	25mm Rounded – Stiff	25	2.4	240	2.8
M00561291	08714729645801	Box 10	Extra Large Rounded – Stiff	33	2.4	240	2.8



Captiflex™ Snare

Order Number	GTIN	Quantity	Descriptor	Loop Width (mm)	Minimum Working Channel (mm)
M00562471	08714729019459	Box 10	Extra Small Oval - Flexible	11	2.8
M00562421	08714729019435	Box 10	Small Oval - Flexible	13	2.8
M00562422	08714729501657	Box 40	Small Oval - Flexible	13	2.8
M00562401	08714729019411	Box 10	Medium Oval - Flexible	27	2.8
M00562402	08714729501640	Box 40	Medium Oval - Flexible	27	2.8




Profile™ Snare

Order Number	GTIN	Quantity	Descriptor	Loop Width (mm)	Minimum Working Channel (mm)
M00562531	08714729156628	Box 10	Mini Micro Oval	11	2.0
M00562551	08714729156659	Box 10	Micro Oval	13	2.0
M00562571	08714729156635	Box 10	Wide Oval	27	2.0




Sensation™ Short Throw Snare

Order Number	GTIN	Quantity	Descriptor	Loop Width (mm)	Minimum Working Channel (mm)
M00562691	08714729158158	Box 5	Small Oval - Flexible	13	2.8
M00562692	08714729268826	Box 20	Small Oval - Flexible	13	2.8
M00562693	08714729501671	Box 40	Small Oval - Flexible	13	2.8
M00560321	08714729283928	Box 5	Small Oval - Medium Stiff	13	2.8
M00562671	08714729158141	Box 5	Medium Oval - Flexible	27	2.8
M00562672	08714729268819	Box 20	Medium Oal - Flexible	27	2.8
M00562673	08714729501664	Box 40	Medium Oval - Flexible	27	2.8
M00560311	08714729283904	Box 5	Medium Oval - Medium Stiff	27	2.8
M00562651	08714729158110	Box 5	Large Oval - Medium Stiff	30	2.8
M00562652	08714729268802	Box 20	Large Oval - Medium Stiff	30	2.8
M00560331	08714729283935	Box 5	Medium Crescent - Medium Stiff	27	2.8



Resolution™ Clip

Order Number	GTIN	Quantity	Working Length (cm)	Minimum Working Channel (mm)	Clip Opening (mm)
M00522600	08714729504764	Each	155	2.8	11
M00522601	08714729504771	Box 10	155	2.8	11
M00522602	08714729504788	Box 20	155	2.8	11
M00522610	08714729504795	Each	235	2.8	11
M00522611	08714729504801	Box 10	235	2.8	11
M00522612	08714729504818	Box 20	235	2.8	11




Interject™ Clear Single-Use Injection Therapy Needle Catheters

Order Number	GTIN	Quantity	Sheath Design	Needle Gauge	Maximum Needle Extension Length (mm)	Sheath O.D. (mm)	Working Length (cm)
M00518101	08714729296416	Box 5	Clear	23	4	1.8	200
M00518111	08714729296423	Box 5	Clear	25	4	1.8	200
M00518301	08714729296508	Box 5	Clear	23	4	2.3	240
M00518311	08714729296515	Box 5	Clear	25	4	2.3	240



TWISTER® PLUS Rotatable Retrieval Device

Order Number	Quantity	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)
M00561411	Box 5	22	2.5	230	2.8
M00561421	Box 5	26	2.5	230	2.8
M00561412	Box 10	22	2.5	230	2.8
M00561422	Box 10	26	2.5	230	2.8



Interject Contrast Single-Use Injection Therapy Needle Catheters

Order Number	GTIN	Quantity	Sheath Design	Needle Gauge	Maximum Needle Extension Length (mm)	Sheath O.D. (mm)	Working Length (cm)
M00518151	08714729296430	Box 5	Contrast	23	4	1.8	200
M00518161	08714729296447	Box 5	Contrast	25	4	1.8	200
M00518251	08714729296478	Box 5	Contrast	23	6	1.8	200
M00518261	08714729296485	Box 5	Contrast	25	6	1.8	200
M00518351	08714729296522	Box 5	Contrast	23	4	2.3	240
M00518361	08714729296539	Box 5	Contrast	25	4	2.3	240

*For the products shown, prior to use consult the product labeling for indications, contraindications, precautions, and warnings.

All trademarks are the property of their respective owners.

TWISTER® PLUS is a trademark of Horizons International Corp. Distributed by Boston Scientific Corporation, Marlborough, Massachusetts; manufactured by Horizons International Corp., Heredia, Costa Rica.

Indications, Contraindications, Warnings and Instructions for Use can be found in the product labeling supplied with each device.

Caution: Federal (U.S.) law restricts this device to sale by or on the order of a physician.

©2015 Boston Scientific Corporation or its affiliates. All rights reserved.