



NRG™
RF Transseptal Kit

**TECHNICAL SPECIFICATIONS
& ORDERING INFORMATION**

Each Kit includes (in individual packaging)

- NRG™ Transseptal Needle
- TorFlex™ Transseptal Guiding Sheath



Choose your needle and sheath configuration

Product Code	Needle Length	Curve	Sheath Length	Sheath curve
NTK3001	71 cm	C0	63 cm	45°
NTK3003	71 cm	C0	63 cm	55°
NTK3007	71 cm	C1	63 cm	45°
NTK3009	71 cm	C1	63 cm	55°



NRG Transseptal Needle*

Needle length: 71 cm
Proximal gauge: 18 ga
Distal gauge: 21 ga



TorFlex Transseptal Guiding Sheath

French size compatibility: 8.5F (2.84 mm)
Sheath usable length: 63cm
Dilator usable length: 67cm

* NRG Transseptal Needle compatible only with RFP-100A RF Puncture Generator.

CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings, and instructions for use can be found in the product labelling supplied with each device or at www.IFU-BSCI.com. Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. This material not intended for use in France. © 2023 Boston Scientific Corporation or its affiliates. All rights reserved. EP-1631601-AA