



**NRG™**  
Transseptal Needle




**TECHNICAL SPECIFICATIONS  
& ORDERING INFORMATION**

Features	Specifications
Needle curves	C0, C1
Needle lengths	56 cm*, 71 cm, 89 cm, 98 cm
Dilator compatibility	0.032"
Needle gauge	Prox: 18 ga, Dist: 21 ga (excl. 56 cm*)

**Choose your needle configuration**

Product Code	Needle Length	Curve	Compatible Transseptal Sheaths
<b>NRG3001</b>	71 cm	C0	8F or 8.5F fixed curve – 63 cm
<b>NRG3003</b>	71 cm	C1	
<b>NRG3005</b>	89 cm	C0	8.5F fixed curve – 81 cm
<b>NRG3007</b>	89 cm	C1	
<b>NRG3009</b>	98 cm	C0	8.5F steerable curve – 71 cm
<b>NRG3011</b>	98 cm	C1	
<b>NRG3013</b>	56 cm*	C0	6F small anatomy fixed curve – 48 cm



\* Proximal gauge 19 ga, distal gauge 22 ga.