



TorFlex™
Transseptal Guiding Sheath



**TECHNICAL SPECIFICATIONS
& ORDERING INFORMATION**

Features	Specifications
French size compatibility	8.5F (2.84 mm)
Sheath usable length	63 cm, 81 cm
Sheath radiopaque marker	Platinum/iridium band
Sheath curves	45°, 55°, 90°, 135°
Dilator usable length	67 cm, 85 cm
Dilator inner diameter	0.033" (0.83 mm)
Guidewire max outer diameter	0.032" (0.81 mm)
Guidewire length	180 cm

Product Code	Sheath Usable Length	Sheath curve
TF3001	63 cm	135°
TF3005	63 cm	45°
TF3007	63 cm	55°
TF3009	63 cm	90°
TF3011	81 cm	135°
TF3015	81 cm	45°
TF3017	81 cm	55°
TF3019	81 cm	90°

