

Ranger™ Family

Paclitaxel-Coated PTA Balloon Catheter

UPN – CATALOG NUMBER				Balloon Diameter (mm)	Balloon Length (mm)	Sheath Compatibility (F)*	Rated Burst Pressure (atm) (kPa)	
Shaft Length 80 cm	Shaft Length 90 cm	Shaft Length 135 cm	Shaft Length 150 cm					
–	H74939220200890**	–	H74939220200810**	2	80	4	14	1418
–	H74939220201090**	–	H74939220201010**	2	100	4	14	1418
–	H74939220201290**	–	H74939220201210**	2	120	4	14	1418
–	H74939220201590**	–	H74939220201510**	2	150	4	14	1418
–	H74939220250890**	–	H74939220250810**	2.5	80	4	14	1418
–	H74939220251090**	–	H74939220251010**	2.5	100	4	14	1418
–	H74939220251290**	–	H74939220251210**	2.5	120	4	14	1418
–	H74939220251590**	–	H74939220251510**	2.5	150	4	14	1418
–	H74939220300890**	–	H74939220300810**	3	80	4	14	1418
–	H74939220301090**	–	H74939220301010**	3	100	4	14	1418
–	H74939220301290**	–	H74939220301210**	3	120	4	14	1418
–	H74939220301590**	–	H74939220301510**	3	150	4	14	1418
–	H74939220350890**	–	H74939220350810**	3.5	80	4	14	1418
–	H74939220351090**	–	H74939220351010**	3.5	100	4	14	1418
–	H74939220351290**	–	H74939220351210**	3.5	120	4	14	1418
–	H74939220351590**	–	H74939220351510**	3.5	150	4	14	1418
H74939219400380	–	H74939219400310	–	4	30	5	14	1418
H74939219400480	–	H74939219400410	–	4	40	5	14	1418
H74939219400680	–	H74939219400610	–	4	60	5	14	1418
–	–	–	H74939220400810**	4	80	4	14	1418
H74939219400880	–	–	–	4	80	5	14	1418
–	–	–	H74939220401010**	4	100	4	14	1418
H74939219401080	–	–	–	4	100	5	14	1418
–	H74939220401290**	–	H74939220401210**	4	120	4	14	1418
–	H74939220401590**	–	H74939220401510**	4	150	4	14	1418
H74939219500380	–	H74939219500310	–	5	30	5	14	1418
H74939219500480	–	H74939219500410	–	5	40	5	14	1418
H74939219500680	–	H74939219500610	–	5	60	5	14	1418
H74939219500880	–	H74939219500810	–	5	80	5	14	1418
H74939219501080	–	H74939219501010	–	5	100	5	14	1418
H74939219600380	–	H74939219600310	–	6	30	5	14	1418
H74939219600480	–	H74939219600410	–	6	40	5	14	1418
H74939219600680	–	H74939219600610	–	6	60	5	14	1418
H74939219600880	–	H74939219600810	–	6	80	5	14	1418
H74939219601080	–	H74939219601010	–	6	100	5	14	1418
H74939219700380	–	H74939219700310	–	7	30	5	14	1418
H74939219700480	–	H74939219700410	–	7	40	5	14	1418
H74939219700680	–	H74939219700610	–	7	60	5	14	1418
H74939219700880	–	H74939219700810	–	7	80	6	14	1418
H74939219701080	–	H74939219701010	–	7	100	6	14	1418
H74939219800380	–	H74939219800310	–	8	30	5	12	1216
H74939219800480	–	H74939219800410	–	8	40	5	12	1216
H74939219800680	–	H74939219800610	–	8	60	5	12	1216
H74939219800880	–	H74939219800810	–	8	80	6	12	1216
Ranger™	Ranger™ SL	Ranger™	Ranger™ SL					

* 4F = 1.33 mm, 5F = 1.7 mm, 6F = 2.0 mm
 ** Ranger™ SL

Ranger™ Paclitaxel-Coated PTA Balloon Catheter is owned by Hemoteg AG.

All cited trademarks are the property of their respective owners. CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device. Information for the use only in countries with applicable health authority product registrations.

Information contained herein is for distribution outside the U.S. only. Illustrations for information purposes – not indicative of actual size or clinical outcome.

PI-231101-AA OCT 2014 Printed in Germany by medicalvision.

Boston Scientific

Advancing science for life™

www.bostonscientific-international.com

© 2014 Boston Scientific Corporation or its affiliates. All rights reserved.
 DINPER4646EC