


EMBOZENE™ Microspheres

PRECISION FLEXIBILITY CONTROL




Precise Calibration

- Offered in a broad range of tightly calibrated sizes
- 95% of the microspheres are within the stated size variance, resulting in greater confidence when performing super selective and targeted embolization procedures

Procedural Flexibility

- Largest array of distinctly sized microspheres
- EMBOZENE's elastic compressibility enables to rebound to original shape and volume after passage through appropriate sized microcatheter

Procedural Control



Applied Use

EMBOZENE™ Microspheres are indicated for the embolization of:

- Hypervascular tumors (H.V.T), incl. Uterine Fibroid Embolization (UFE) and hepatocellular carcinoma
- Arteriovenous malformations (A.V.M.)
- Benign Prostate Hyperplasia (BPH)
- Tumors of head, neck, torso, and skeletal system

EMBOZENE™ Microspheres

EMBOZENE™ Microspheres are available in prefilled 1 ml and 2 ml syringes.

Size	REF 1 ml	REF 2 ml
40 µm	10410-S1	10420-S1
75 µm	10710-S1	10720-S1
100 µm	11010-S1	11020-S1
250 µm	12010-S1	12020-S1
400 µm	14010-S1	14020-S1
500 µm	15010-S1	15020-S1
700 µm	17010-S1	17020-S1
900 µm	19010-S1	19020-S1
1100 µm	111010-S1	111020-S1
1300 µm	113010-S1	113020-S1



EMBOZENE™ Microspheres

- Color-coded to enhance visualization and aid in size recognition
- Sold as single syringe units
- 3-year shelf life from date of manufacture

EMBOZENE™ Microspheres are registered or unregistered trademarks of Boston Scientific Corporation or its affiliates.

CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device. Information for the use only in countries with applicable health authority product registrations. Information not intended for use or distribution in France. Illustrations are not necessary to scale. All images are property of Boston Scientific.

PI-384807-AA APR 2016 Printed in Germany by medicalvision.

**Boston
Scientific**
Advancing science for life™

www.bostonscientific-international.com

© 2016 Boston Scientific Corporation
or its affiliates. All rights reserved.
DINONC3519EA