



ONCOZENE™ Microspheres

SIZED FOR SUCCESS




More control and confidence in your embolization procedures
with the most precisely calibrated microspheres.



ONCOZENE™ Microspheres

ONCOZENE Microspheres offer:

- Controlled embolization with precisely calibrated microspheres
- Strength in Numbers: Precisely calibrated microspheres designed to deliver more microspheres per ml
- Sized for tumor penetration with unique microspheres sizes: 40um, 75um, 100um
- High confidence when performing super selective and targeted embolization procedures as 95% of the microspheres are within the stated size variance
- Greater flexibility with single sale units in 2 ml and 3 ml sized syringes

			
2 ml Syringe	40 ± 10 µm	75 ± 15 µm	100 ± 25 µm
UPN	06-0301-00402-02	06-0301-00752-02	06-0301-01002-02
Catalog No.	10420-US1	10720-US1	11020-US1
3 ml Syringe	40 ± 10 µm	75 ± 15 µm	100 ± 25 µm
UPN	06-0301-00403-02	06-0301-00753-02	06-0301-01003-02
Catalog No.	10430-US1	10730-US1	11030-US1

ONCOZENE™ MICROSPHERES

INDICATIONS FOR USE: ONCOZENE™ Microspheres are indicated for the embolization of arteriovenous malformations and hypervascular tumors including hepatoma. **CONTRAINDICATIONS:** The contraindications of ONCOZENE™ Microspheres include the presence of vasculature where ONCOZENE™ Microspheres could pass directly into the central nervous system, central circulatory system, internal carotid artery, or other non-target territories. Procedures should not be performed if vascular anatomy precludes correct catheter placement or embolic injection. **WARNINGS AND PRECAUTIONS:** Vascular embolization is a high-risk procedure. The procedure should be performed by specialized physicians trained in vascular embolization procedures. Complications can occur at any time during or after the procedure. **CAUTION:** Federal (USA) law restricts this device to sale by or on the order of a physician.

A complete list of indications, contraindications, warnings and precautions are described in ONCOZENE™ Microspheres Instructions for Use. Please consult these before using the product.

Oncozene is a unregistered or registered trademark of Boston Scientific Corporation or its affiliates. All other trademarks are property of their respective owners.

**Boston
Scientific**
Advancing science for life™

Peripheral Interventions
300 Boston Scientific Way
Marlborough, MA 01752-1234
www.bostonscientific.com

To order product or for more information contact customer service at 1.888.272.1001.

© 2018 Boston Scientific Corporation or its affiliates. All rights reserved.

PI-363807-AC JAN2018