

Explorer ST[™]

Fixed Curve Diagnostic Catheter

Boston
Scientific

**High Performance
with Security
of a Soft Distal**



Explorer ST™ Fixed Curve Diagnostic Catheter


Clinical Benefits

Combining materials and technology to improve security and performance



- **Available in decapolar configuration** for use in the coronary sinus
- **Soft tip designed** to reduce trauma to the vasculature and cardiac tissues
- **Braided proximal shaft** provides excellent control and maneuverability
- **Quick Connector** design for intuitive connections and excellent signal quality


Technical Information

Description	Specifications	Description	Specifications
Usable length	100cm	Shaft size	5F and 6F
Catheter material	Braided Polyurethane	Distal tip	10cm soft distal tip
Electrode material	Platinum	Connector	Quick Connect



Ordering Information

Model Number	Electrode Configuration	Shaft Ø	Curve Option	Electrode Spacing	Recommended Cables
M004 5470 0	Quadripolar	6F	Josephson	5mm	M004 5446S 0
M004 5481 0	Quadripolar	6F	Courmand	5mm	M004 5446S 0
M004 5457 0	Decapolar	6F	Courmand	2/5/2mm	M004 5454S 0
M004 5459 0	Decapolar	6F	Courmand	2/8/2mm	M004 5454S 0
M004 5563 0	Quadripolar	6F	Multipurpose	5mm	M004 5446S 0
M004 5458 0	Decapolar	6F	Courmand	2mm	M004 5454S 0
M004 5477 0 	Decapolar	6F	Damato	2/8/2mm	M004 5454S 0
M004 5469 0	Quadripolar	6F	Josephson	10mm	M004 5446S 0
M004 5480 0	Quadripolar	6F	Courmand	10mm	M004 5446S 0
M004 5561 0	Quadripolar	6F	Conduction Study	5mm	M004 5446S 0
M004 5920 0 	Quadripolar	6F	K-Curve	2mm	M004 5446S 0
M004 5418 0	Quadripolar	5F	Josephson	5mm	M004 5446S 0
M004 5420 0	Quadripolar	5F	Courmand	5mm	M004 5446S 0
M004 5451 0	Quadripolar	5F	Josephson	2/5/2mm	M004 5446S 0
M004 5449 0	Quadripolar	5F	Courmand	2/5/2mm	M004 5446S 0

 Special order: minimum quantity 25 units, 8-10 weeks delivery time

Boston Scientific

Delivering what's next.™

Boston Scientific Corporation
2710 Orchard Parkway
San Jose, CA 95134
www.bostonscientific.com/electrophysiology

To order product or for more information, contact Electrophysiology customer service at 1-888-272-1001.

© 2006 by Boston Scientific Corporation or its affiliates. All rights reserved.

Printed in the US. EPT-10551_04/06

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician. Carefully read all instructions prior to use. Observe all contraindications, warnings and precautions noted in these instructions. Failure to do so may result in patient complications. Boston Scientific relies on the physician to determine, assess and communicate to each patient all foreseeable risks of the procedure.