


AFocus™ Steerable Diagnostic Catheters

AFocus I: Steerable Loop
Atrial Mapping Catheters



17 - 20mm diameter

AFocus II: Steerable Shaft
Atrial Mapping Catheters



Curve SM

Product Information

AFocus I: Steerable Loop Catheters - 5F

Catheter Model No.	Description	French Size	No. of Electrodes	Unipolar/Bipolar	Electrode Spacing	Curve Diameter
M004 81674 0	1110-5-5-AF20	5	10	n/a	5mm	17-20mm

AFocus II: Steerable Shaft Catheters - 7F and 4F (7F shaft with a 4F loop)

Catheter Model No.	Description	French Size	No. of Electrodes	Unipolar/Bipolar	Electrode Spacing	Curve Diameter
M004 81587 0	1120-7-1(2.5)-SM-AF20	7/4	20	Bipolar	1-2.5-1mm	20mm
M004 81594 0	1110-7-3.5-SM-AF15	7/4	10	Unipolar	3.5mm	15mm
M004 81595 0	1110-7-5-SM-AF20	7/4	10	Unipolar	5mm	20mm
M004 81598 0	1120-7-1(4.5)-SM-AF25	7/4	20	Bipolar	1-4.5-1mm	25mm

Note: All catheters are 110cm in length. All AFocus Catheters use Cable Model No. 85954 except M004 **81587** 0 and M004 **81598** 0, which use Cable Model No 85931.

Indications, contraindications, warnings and instructions for use can be found in the instructions for use supplied with each device. CAUTION: Federal (USA) law restricts these devices to sale by or on the order of a physician.

©2007 Boston Scientific Corporation or its affiliates. All rights reserved. EPT-0015 _ 04/07

**Boston
Scientific**

For additional product information call Customer Service at (888) 272-1001

www.bostonscientific.com/electrophysiology