

Chilli II™ Cooled Ablation Catheters



Product Information

Tip: 7F/4mm

Electrode Spacing: 2.5mm

Catheter Model No.	Type	Shaft Size	Curve Style	Distal Shaft Length
M004 9031 0	Thermocouple	7.5F	Standard	Standard
M004 9031K2 0	Thermocouple	7.5F	Large	Standard
M004 9031N4 0	Thermocouple	7.5F	Asymmetric 4	Standard

Cables and Accessories

M004 352104 0 CircuCool® Tubing Set

For use with Maestro 3000™ Controller


M004 681 0 Cable: Maestro 3000 Pod to Chilli II Catheter (10 ft.)

For use with EPT-1000XP™ Controller

M004 2071 0 Hardware: Thermocouple-Thermistor-Translator (TTT) for EPT-1000XP Controllers

M004 62062 0 Cable: Chilli II™ Catheter to TTT for EPT-1000XP Controllers (8 ft.)

Note: All Chilli II Catheters are 110cm in length.
Chilli II Catheters require the use of 1 Valleylab® Ground Pad (Model M004 354 0).



Indications, contraindications, warnings and instructions for use can be found in the instructions for use supplied with each device. CAUTION: Federal (USA) law restricts these devices to sale by or on the order of a physician.


©2007 Boston Scientific Corporation or its affiliates. All rights reserved. EPT-0015_04/07

**Boston
Scientific**

For additional product information call Customer Service at (888) 272-1001

www.bostonscientific.com/electrophysiology

Chilli II™ Cooled Ablation System



Product Information

Model Number	Description
M004 8005 0	CircuCool Fluid Circulation Pump

Cables and Accessories

M004 352104 0	CircuCool Tubing Set
M004 354 0	Valleylab® Ground Pad

For use with Maestro 3000™ Generator

M004 681 0	Cable: Maestro 3000 Pod to Chilli II™ Catheter (10ft.)
------------	--

For use with EPT-1000XP™ Generator

M004 2071 0	Hardware: Thermocouple-Thermistor-Translator (TTT) for use with EPT-1000XP Controller
M004 62062 0	Cable: Chilli II Catheter to TTT for use with EPT-1000XP Controller (8 ft.)

Technical Information: CircuCool Pump

Description	Specification
AC power requirements	100-240 VAC
Environmental conditions:	
Ambient temperature	18-35°C
Relative humidity	15-90% noncondensing relative humidity
Pump flow rate	0.6ml/second
Pressure limit	Less than 300psi
Dimensions	11in W x 11.8in D x 9in H
Weight	9.6lbs

Extended Warranty options are available, see "Service and Warranties" of the pricing catalog for details.

Indications, contraindications, warnings and instructions for use can be found in the instructions for use supplied with each device. CAUTION: Federal (USA) law restricts these devices to sale by or on the order of a physician.

©2007 Boston Scientific Corporation or its affiliates. All rights reserved. EPT-0015 _ 04/07

**Boston
Scientific**

For additional product information call Customer Service at (888) 272-1001

www.bostonscientific.com/electrophysiology