


Constellation® Mapping Catheters



Product Information

Constellation Mapping Catheters: *Catheter length 120cm, 64 electrodes, 8F shaft with 2 x 44 pin D-sub male connectors*

Catheter Model No.	Basket Diameter	Electrode Spacing
M004 US8031U 0	31mm	2mm
M004 US8038U 0	38mm	3mm
M004 US8048U 0	48mm	4mm
M004 US8060U 0	60mm	5mm
M004 US8075U 0	75mm	7mm

Constellation Cables

Model No.	Description
M004 900 0	Extension Cable (10 ft.)
M004 921B 0	Bard Adapter Cable
M004 921P 0	Prucka Adapter Cable
M004 940 0	8-pin Quick Connect (951 and 952M to recorder)
M004 626S 0	Octapolar Cable—Quick Female to Shrouded Pins (3 ft.)
M004 657 0	Quick Connect Spline Recording Cable (3 ft.)
M004 930S 0	Octapolar Adapter
M004 927 0	RPM™ System Adapter Cables (set)
M004 905S 0	SLD Ground Plug Adapter Cable

Note: Constellation Mapping Catheters are to be used with the Boston Scientific Convoy® Advanced Delivery Sheaths.

Indications, contraindications, warnings and instructions for use can be found in the instructions for use supplied with each device. CAUTION: Federal (USA) law restricts these devices to sale by or on the order of a physician.

©2007 Boston Scientific Corporation or its affiliates. All rights reserved. EPT-0015 _ 04/07

**Boston
Scientific**

For additional product information call Customer Service at (888) 272-1001

www.bostonscientific.com/electrophysiology

Constellation® Catheter Mapping System



Product Information

Model No.	Description
M004 952M 0	Constellation Catheter Mapping System
	Includes:
	• Switching/Locating Device
	• Astronomer Plus® Software
	• Model 952A Remote Display
	• RS232 computer serial cable
	• Model M004 905S 0 cable

Extended Warranty options are available, see "Service and Warranties" of the pricing catalog for details.

Indications, contraindications, warnings and instructions for use can be found in the instructions for use supplied with each device. CAUTION: Federal (USA) law restricts these devices to sale by or on the order of a physician.

©2007 Boston Scientific Corporation or its affiliates. All rights reserved. EPT-0015 _ 04/07

**Boston
Scientific**

For additional product information call Customer Service at (888) 272-1001
www.bostonscientific.com/electrophysiology