

# Explorer 360™ Fixed Curve Diagnostic Catheters



## Product Information

**Tip:** 6F


Catheter Model No.	Shaft Size	Curve Style	Electrode Configuration	Electrode Spacing	Cable Model No.
M004 5292S 0	6F	Josephson	Quadripolar	2mm	M004 5444S 0
M004 5291S 0	6F	Josephson	Quadripolar	5mm	M004 5444S 0
M004 5290S 0	6F	Josephson	Quadripolar	10mm	M004 5444S 0
M004 5461S 0	6F	Josephson	Quadripolar	2/5/2mm	M004 5444S 0
M004 5295S 0	6F	Cournand	Quadripolar	2mm	M004 5444S 0
M004 5294S 0	6F	Cournand	Quadripolar	5mm	M004 5444S 0
M004 5298S 0	6F	Cournand	Quadripolar	10mm	M004 5444S 0
M004 5229S 0	6F	Conduction Study	Quadripolar	5mm	M004 5444S 0
M004 5222S 0	6F	Conduction Study	Quadripolar	10mm	M004 5444S 0
M004 5433S 0	6F	Multipurpose	Quadripolar	5mm	M004 5444S 0

**Made to Order Item (minimum order 25 units; see Ordering Information section for details)**

M004 5429S 0	6F	Damato	Quadripolar	5mm	M004 5444S 0
--------------	----	--------	-------------	-----	--------------

Note: All Explorer 360 Diagnostic Catheters are 100cm in length.

## Explorer 360 Catheter Curve Options



Indications, contraindications, warnings and instructions for use can be found in the instructions for use supplied with each device. CAUTION: Federal (USA) law restricts these devices to sale by or on the order of a physician.

©2007 Boston Scientific Corporation or its affiliates. All rights reserved. EPT-0015 \_ 04/07

**Boston  
Scientific**

For additional product information call Customer Service at (888) 272-1001

[www.bostonscientific.com/electrophysiology](http://www.bostonscientific.com/electrophysiology)

# Explorer 360 Jr.™ Fixed Curve Diagnostic Catheters




## Product Information

**Tip:** 5F

Catheter Model No.	Shaft Size	Curve Style	Electrode Configuration	Electrode Spacing	Cable Model No.
M004 5413S 0	5F	Josephson	Quadripolar	2mm	M004 5444S 0
M004 5404S 0	5F	Josephson	Quadripolar	5mm	M004 5444S 0
M004 5414S 0	5F	Cournand	Quadripolar	5mm	M004 5444S 0
M004 5409S 0	5F	Conduction Study	Quadripolar	5mm	M004 5444S 0
M004 5415A1S 0	5F	Multipurpose	Quadripolar	10mm	M004 5446S 0

\*Note: All Explorer 360 Jr. Diagnostic Catheters are 100cm in length EXCEPT for the M004 5415A1S 0 which is 125cm in length.

## Explorer 360 Jr. Catheter Curve Options



Indications, contraindications, warnings and instructions for use can be found in the instructions for use supplied with each device. CAUTION: Federal (USA) law restricts these devices to sale by or on the order of a physician.

©2007 Boston Scientific Corporation or its affiliates. All rights reserved. EPT-0015 \_ 04/07

**Boston  
Scientific**

For additional product information call Customer Service at (888) 272-1001

[www.bostonscientific.com/electrophysiology](http://www.bostonscientific.com/electrophysiology)