


Explorer ST™ Fixed Curve Diagnostic Catheters



Product Information

Tip: 6F (10cm Soft Distal Tip)


Catheter Model No.	Shaft Size	Curve Style	Electrode Configuration	Electrode Spacing	Cable Model No.
M004 5470 0	6F	Josephson	Quadripolar	5mm	M004 5446S 0
M004 5469 0	6F	Josephson	Quadripolar	10mm	M004 5446S 0
M004 5481 0	6F	Cournand	Quadripolar	5mm	M004 5446S 0
M004 5480 0	6F	Cournand	Quadripolar	10mm	M004 5446S 0
M004 5561 0	6F	Conduction Study	Quadripolar	5mm	M004 5446S 0
M004 5563 0	6F	Multipurpose	Quadripolar	5mm	M004 5446S 0
M004 5457 0	6F	Cournand	Decapolar	2/5/2mm	M004 5454S 0
M004 5458 0	6F	Cournand	Decapolar	2mm	M004 5454S 0
M004 5459 0	6F	Cournand	Decapolar	2/8/2mm	M004 5454S 0

Made to Order Items (minimum order 25 units; see Ordering Information section for details)

M004 5477 0	6F	Damato	Decapolar	2/8/2mm	M004 5454S 0
M004 5920 0	6F	K-Curve	Quadripolar	2mm	M004 5446S 0

Note: All Explorer ST Fixed Curve Diagnostic Catheters are 100cm in length.

Explorer ST Catheter Curve Options



Indications, contraindications, warnings and instructions for use can be found in the instructions for use supplied with each device. CAUTION: Federal (USA) law restricts these devices to sale by or on the order of a physician.


©2007 Boston Scientific Corporation or its affiliates. All rights reserved. EPT-0015 _ 04/07

**Boston
Scientific**

For additional product information call Customer Service at (888) 272-1001

www.bostonscientific.com/electrophysiology

Explorer ST™ Fixed Curve Diagnostic Catheters




Product Information

Tip: 5F (10cm Soft Distal Tip)

Catheter Model No.	Shaft Size	Curve Style	Electrode Configuration	Electrode Spacing	Cable Model No.
M004 5418 0	5F	Josephson	Quadripolar	5mm	M004 5446S 0
M004 5451 0	5F	Josephson	Quadripolar	2/5/2mm	M004 5446S 0
M004 5420 0	5F	Cournand	Quadripolar	5mm	M004 5446S 0
M004 5449 0	5F	Cournand	Quadripolar	2/5/2mm	M004 5446S 0

Note: All Explorer ST Fixed Curve Diagnostic Catheters are 100cm in length.

Explorer ST (5F) Catheter Curve Options



Indications, contraindications, warnings and instructions for use can be found in the instructions for use supplied with each device. CAUTION: Federal (USA) law restricts these devices to sale by or on the order of a physician.

©2007 Boston Scientific Corporation or its affiliates. All rights reserved. EPT-0015 _ 04/07

**Boston
Scientific**

For additional product information call Customer Service at (888) 272-1001

www.bostonscientific.com/electrophysiology