

Intracardiac Ultrasound System & Catheters

iLab™ System Cart



Product Information

Model No.	Description
M004 EPiL120CART 0	iLab Ultrasound Imaging System (on cart)
M004 EPiL120INS 0	iLab Ultrasound Imaging System (installed*)



*Installation of the system is not included in the system price and will be quoted separately. Average installation is \$10,400, but will vary based on the requirements of the site.

Extended Warranty options are available, see "Service and Warranties" of the pricing catalog for details.

Catheter Information	Model No.	Description
	M004 9900 0	Ultra ICE™ Catheter Size: 9F Frequency: 9 MHz Usable Length: 110cm
	M004 9915 1	Fluid Dock® Filling Device Box of 10 individual sterile units (Filling device enables filling of the Ultra ICE Catheter tip without damage to the ultrasound sensor during catheter preparation.)

Please refer to the Sheath section for Convoy® Advanced Delivery Sheaths, recommended to deploy the Ultra ICE Catheter. For septal imaging, use Convoy Sheath Model M004 **5662ST** 0.

Indications, contraindications, warnings and instructions for use can be found in the instructions for use supplied with each device. CAUTION: Federal (USA) law restricts these devices to sale by or on the order of a physician.

©2007 Boston Scientific Corporation or its affiliates. All rights reserved. EPT-0015 _ 04/07

**Boston
Scientific**

For additional product information call Customer Service at (888) 272-1001

www.bostonscientific.com/electrophysiology