


Inquiry® Fixed Curve Diagnostic Catheters



Product Information

QUADRIPOLAR - 4F Diameter, Electrode Size (mm) Tip 1/Band 1

Catheter Model No.	Description	French Size	Electrode Spacing	Curve Type	Cable Model No.
M004 80465 0	1004-4-5-C	4F	5mm	Cournand	85953
M004 80466 0	1004-4-25-C	4F	2-5-2mm	Cournand	85953
M004 80468 0	1004-4-5-J	4F	5mm	Josephson	85953

QUADRIPOLAR - 6F Diameter, Electrode Size (mm) Tip 2/Band 2

Catheter Model No.	Description	French Size	Electrode Spacing	Curve Type	Cable Model No.
M004 80409 0	1004-6-5-C	6F	5mm	Cournand	85953

DECAPOLAR - 6F Diameter, Electrode Size (mm) Tip 4/Band 4

Catheter Model No.	Description	French Size	Electrode Spacing	Curve Type	Cable Model No.
M004 80116 0	1010-6-2-J-TE4BE4	6F	2mm	Josephson	85954

Note: All catheters are 110cm in length.

Indications, contraindications, warnings and instructions for use can be found in the instructions for use supplied with each device. CAUTION: Federal (USA) law restricts these devices to sale by or on the order of a physician.

©2007 Boston Scientific Corporation or its affiliates. All rights reserved. EPT-0015 _ 04/07

**Boston
Scientific**

For additional product information call Customer Service at (888) 272-1001

www.bostonscientific.com/electrophysiology