


Inquiry® HIS (SOFT) Diagnostic Catheters



Steerable E (HIS) Curve



Fixed C (HIS) Curve

Product Information

QUADRIPOLAR - 5F Diameter, Electrode Size (mm) Tip 1/Band 1

Catheter Model No.	Description	French Size	Electrode Spacing	Curve Type	Cable Model No.
M004 80440 0	1004-5-5-C (HIS)	5F	5mm	C (HIS) Fixed	85953
M004 81483 0	1104-5-5-E (HIS)	5F	5mm	E (HIS) Steerable	85953

Note: All catheters are 110cm in length.

Indications, contraindications, warnings and instructions for use can be found in the instructions for use supplied with each device. CAUTION: Federal (USA) law restricts these devices to sale by or on the order of a physician.

©2007 Boston Scientific Corporation or its affiliates. All rights reserved. EPT-0015 _ 04/07

**Boston
Scientific**

For additional product information call Customer Service at (888) 272-1001

www.bostonscientific.com/electrophysiology