


Inquiry® H-Curve™ Steerable Diagnostic Catheters



Product Information

Regular Curve - Electrode size (mm) Tip 1/Band 1

| Catheter Model No. | Description | French Size | Electrode Spacing | Curve Type | Tip | Cable Model No. |
|--------------------|------------------|-------------|-------------------|---|----------|-----------------|
| M004 81120 0 | 1120-7-17-H | 7F | 1-7-1mm | H-Curve; 10 pair electrodes | Inactive | 85931 |
| M004 81128 0 | 1121-7-5-H-UP | 7F | 5mm | H-Curve; 20 uni-polar electrodes plus reference at 35cm | Inactive | 85931 |
| M004 81124 0 | 1124-7-271-H | 7F | 2-7-1mm | H-Curve; 12 pair electrodes | Active | 85931 |
| M004 81150 0 | 1124-7-271-H-SCE | 7F | 2-7-1mm | H-SCE Curve; 12 pair electrodes | Active | 85931 |

Large Curve - Electrode size (mm) Tip 1/Band 1

| Catheter Model No. | Description | French Size | Electrode Spacing | Curve Type | Tip | Cable Model No. |
|--------------------|-----------------|-------------|-------------------|---|----------|-----------------|
| M004 81130 0 | 1120-7-19-HL | 7F | 1-9-1mm | HL-Curve; 10 pair electrodes | Inactive | 85931 |
| M004 81131 0 | 1121-7-19-HL-SC | 7F | 1-9-1mm | HL-SC Curve; tip plus 10 pair band electrodes | Active | 85931 |
| M004 81134 0 | 1124-7-291-HL | 7F | 2-9-1mm | HL-Curve; 12 pair electrodes | Active | 85931 |

Note: All catheters are 110cm in length.

Indications, contraindications, warnings and instructions for use can be found in the instructions for use supplied with each device. CAUTION: Federal (USA) law restricts these devices to sale by or on the order of a physician.

©2007 Boston Scientific Corporation or its affiliates. All rights reserved. EPT-0015 _ 04/07

**Boston
Scientific**

For additional product information call Customer Service at (888) 272-1001

www.bostonscientific.com/electrophysiology