


Inquiry® Steerable Curve Diagnostic Catheters



Medium- 3.2cm Ø

Large - 4.3cm Ø

Product Information

QUADRIPOLAR - 4F Diameter, Electrode Size (mm) Tip 1/Band 1

Catheter Model No.	Description	French Size	Electrode Spacing	Curve Type	Cable Model No.
M004 81540 0	1104-4-25-M	4F	2-5-2mm	Medium	85953
M004 81541 0	1104-4-2-M	4F	2mm	Medium	85953
M004 81542 0	1104-4-5-M	4F	5mm	Medium	85953
M004 81543 0	1104-4-25-L	4F	2-5-2mm	Large	85953

DECAPOLAR - 4F Diameter, Electrode Size (mm) Tip 1/Band 1

Catheter Model No.	Description	French Size	Electrode Spacing	Curve Type	Cable Model No.
M004 81530 0	1110-4-2-M	4F	2mm	Medium	85954
M004 81531 0	1110-4-25-M	4F	2-5-2mm	Medium	85954
M004 81532 0	1110-4-25-L	4F	2-5-2mm	Large	85954

OCTAPOLAR - 5F Diameter, Electrode Size (mm) Tip 1/Band 1

Catheter Model No.	Description	French Size	Electrode Spacing	Curve Type	Cable Model No.
M004 81871 0	1108-5-2-M	5F	2mm	Medium	85954

DECAPOLAR - 5F Diameter, Electrode Size (mm) Tip 1/Band 1

Catheter Model No.	Description	French Size	Electrode Spacing	Curve Type	Cable Model No.
M004 81171 0	1110-5-2-M	5F	2mm	Medium	85954
M004 81172 0	1110-5-25-M	5F	2-5-2mm	Medium	85954
M004 81174 0	1110-5-25-L	5F	2-5-2mm	Large	85954

Note: All catheters are 110cm in length.

Indications, contraindications, warnings and instructions for use can be found in the instructions for use supplied with each device. CAUTION: Federal (USA) law restricts these devices to sale by or on the order of a physician.


©2007 Boston Scientific Corporation or its affiliates. All rights reserved. EPT-0015 _ 04/07

**Boston
Scientific**

For additional product information call Customer Service at (888) 272-1001

www.bostonscientific.com/electrophysiology

Inquiry® Steerable Curve Diagnostic Catheters



Product Information

QUADRIPOlar - 5F Diameter, Electrode Size (mm) Tip 1/Band 1

Catheter Model No.	Description	French Size	Electrode Spacing	Curve Type	Cable Model No.
M004 81472 0	1104-5-25-M	5F	2-5-2mm	Medium	85953
M004 81473 0	1104-5-5-M	5F	5mm	Medium	85953
M004 81474 0	1104-5-25-L	5F	2-5-2mm	Large	85953

DECAPOlar - 6F Diameter, Electrode Size (mm) Tip 4/Band 4

Catheter Model No.	Description	French Size	Electrode Spacing	Curve Type	Cable Model No.
M004 81104 0	1110-6-25-L	6F	2-5-2mm	Large	85954
M004 81520 0	1110-6-2-XL-TE4BE4	6F	2mm	Xtra Large	85954
M004 81524 0	1110-6-2-L-TE4BE4	6F	2mm	Large	85954

DECAPOlar - 6F Diameter, Bi-directional, Electrode Size (mm) Tip 2/Band 2

Catheter Model No.	Description	French Size	Electrode Spacing	Curve Type	Cable Model No.
M004 81504 0	1110-6-5-M-TE2BE2-BD	6F	5mm	Medium (2.4cm Ø)	85954

Note: All catheters are 110cm in length.

Indications, contraindications, warnings and instructions for use can be found in the instructions for use supplied with each device. CAUTION: Federal (USA) law restricts these devices to sale by or on the order of a physician.

©2007 Boston Scientific Corporation or its affiliates. All rights reserved. EPT-0015 _ 04/07

**Boston
Scientific**

For additional product information call Customer Service at (888) 272-1001

www.bostonscientific.com/electrophysiology