


# Maestro 3000™ Cardiac Ablation System



**Controller**

**Remote Control Unit**

**Maestro 3000 Pod**

## Product Information

Model Number	Type	Description
M004 21000TC 0	Hardware	Maestro 3000 Controller
M004 21860T 0	Hardware	Maestro 3000 Pod
M004 21880 0	Hardware	Maestro 3000 Remote Control (optional)
M004 21840 0	Hardware	Maestro 3000 Controller Foot Switch

## Cables and Accessories

Model Number	Type	Description	Length
M004 613 0	Cable	Pod (or APM) to standard or high-power catheter	10 ft.
M004 651 0	Cable	Pod (or APM) to standard or high-power catheter	3 ft.
M004 681 0	Cable	Maestro 3000 Pod to Chilli II™ Catheter	10 ft.
M004 62053 0	Cable	Maestro 3000 Pod to Blazer RPM/XP RPM™ Catheter	10 ft.
M004 653S 0	Cable	Pod (or APM) to EGM	3 ft.
M004 21600 0	Cable	Maestro 3000 Pod to PAM	8 ft.
M004 SVGA25 0	Cable	Maestro 3000 Controller to Remote Control	25 ft.
M004 SVGA50 0	Cable	Maestro 3000 Controller to Remote Control	50 ft.
M004 SVGA75 0	Cable	Maestro 3000 Controller to Remote Control	75 ft.
M004 354 0	Electrode	Valleylab® Ground Pad	n/a
M004 21010 0	Manual	Maestro 3000 System Operators Manual	n/a

## Technical Information

Description	Specification
Line Power	100-120 / 220-240 V~
Output Power	100 watts maximum
Operating Frequency	460 kHz (nominal)
Operating Temperature	15°C to 40°C
Relative Humidity	30% to 75%, non-condensing
Sound Output Power Up/Sign On Tone	500 Hz - 1.5 kHz - 2.5 kHz Square Wave
<b>Maestro 3000 Controller</b>	Dimensions 7.2in H x 12.9in W x 15in D Weight 16.95 lbs
<b>Maestro 3000 Pod</b>	Dimensions 2.3in H x 6.8in W x 6.2in D Weight 2.2 lbs
<b>Maestro 3000 Remote Control</b>	Dimensions 7.2in H x 12.9in W x 7.4in D Weight 5.8 lbs

Extended Warranty options are available, see "Service and Warranties" of the pricing catalog for details.

Indications, contraindications, warnings and instructions for use can be found in the instructions for use supplied with each device. CAUTION: Federal (USA) law restricts these devices to sale by or on the order of a physician.

©2007 Boston Scientific Corporation or its affiliates. All rights reserved. EPT-0015\_04/07

**Boston  
Scientific**

For additional product information call Customer Service at (888) 272-1001

[www.bostonscientific.com/electrophysiology](http://www.bostonscientific.com/electrophysiology)