


# Polaris Dx™ Steerable Diagnostic Catheters





## Product Information

Catheter Model No.	Shaft Size	Curve Style	Electrode Configuration	Electrode Spacing	Cable Model No.
M004 5427S 0*	7F	Standard	Quadripolar	5mm	M004 5444S 0
M004 5570S 0	6F	Standard	Quadripolar	5mm	M004 5444S 0
M004 5571S 0	6F	Standard	Quadripolar	10mm	M004 5444S 0
M004 5572S 0	6F	Standard	Quadripolar	2/5/2mm	M004 5444S 0
M004 5573S 0	6F	Longest	Quadripolar	5mm	M004 5444S 0
M004 5574S 0	6F	Longest	Quadripolar	2/5/2mm	M004 5444S 0
M004 5575 0	6F	Standard	Hexapolar	2mm	M004 5440S 0
M004 5577 0	6F	Longest	Hexapolar	5mm	M004 5440S 0
M004 5579 0	6F	Standard	Octapolar	2/5/2mm	M004 5441S 0
M004 9663S 0*	6F	Standard	Quadripolar	5mm	M004 5444S 0

### Made to Order Items (minimum order 25 units; see Ordering Information section for details)

M004 5576 0	6F	Standard	Hexapolar	2/5/2mm	M004 5440S 0
M004 5578 0	6F	Longest	Hexapolar	20/5/5mm	M004 5440S 0

\*Note: All Polaris Dx Steerable Diagnostic Catheters are 100cm in length, except for M004 5427S 0 and M004 9663S 0 which are 125cm in length.



Indications, contraindications, warnings and instructions for use can be found in the instructions for use supplied with each device. CAUTION: Federal (USA) law restricts these devices to sale by or on the order of a physician.

©2007 Boston Scientific Corporation or its affiliates. All rights reserved. EPT-0015\_04/07

**Boston  
Scientific**

For additional product information call Customer Service at (888) 272-1001

[www.bostonscientific.com/electrophysiology](http://www.bostonscientific.com/electrophysiology)