


# Special Procedure SteeroCath-Dx™ Catheters




## Product Information

Catheter Model No.	Shaft Size	Distal Shaft Length	Tip	Curve Style	Electrode Configuration	Electrode Spacing	Cable Model No.
M004 2027BL 0	7F	Extended	7F/2mm	Standard	Octapolar	2.5mm	M004 626S 0
M004 2028BL 0	7F	Extended	7F/2mm	Standard	Octapolar	2.5/5/2.5mm	M004 626S 0

Note: All Special Procedure Octapolar SteeroCath-Dx Catheters are 110cm in length.

### Special Procedure Curve



Indications, contraindications, warnings and instructions for use can be found in the instructions for use supplied with each device. CAUTION: Federal (USA) law restricts these devices to sale by or on the order of a physician.

©2007 Boston Scientific Corporation or its affiliates. All rights reserved. EPT-0015\_04/07

**Boston  
Scientific**

For additional product information call Customer Service at (888) 272-1001

[www.bostonscientific.com/electrophysiology](http://www.bostonscientific.com/electrophysiology)