


Polypectomy Snares

Boston
Scientific



Captivator™ II

Captiflex™

Profile™

Captivator Crescent

Captivator Hexagonal

ROTATABLE
Boston
Scientific

Polypectomy Snares

PRODUCT DESCRIPTIONS

Captivator™ II Single-Use Snare

- ◆ Rounded shape intended to facilitate polyp capture

Captivator Single-Use Snare

- ◆ Offers a variety of snare loop configurations to accommodate a wide range of therapeutic functions

Single-Use Rotatable Snare

- ◆ Designed to provide 360° rotation to help access difficult polyps

Sensation™ Single-Use Short Throw Snare

- ◆ Short throw handle is designed for nurse and technician comfort

Captiflex™ Single-Use Snare

- ◆ Incorporates a braided wire design engineered to produce a flexible snare

Profile™ Single-Use Pediatric Snare

- ◆ Engineered to fit through most 2.0mm pediatric working channels

ORDERING INFORMATION

NEW CAPTIVATOR II SINGLE-USE SNARE

Order Number	Description	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)
M00561291	Rounded Shape	33	2.4	240

CAPTIVATOR SINGLE-USE SNARES

Order Number	Description	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)
M00561311	Jumbo Oval	30	2.4	240
M00562301	Micro Oval	13	2.4	240
M00562321	Standard Oval	27	2.4	240
M00562341	Hexagonal	27	2.4	240
M00562371	Crescent	27	2.4	240
M00562391	Thin-Wire Jumbo Oval	30	2.4	240
M00562451	Micro-Hex	13	2.4	240

SINGLE-USE ROTATABLE SNARES*

Order Number	Description	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)
M00561821	Micro Oval	13	2.4	195
M00561831	Mini-Standard Oval	20	2.4	195

*With Star Catheter Technology

SENSATION SINGLE-USE SHORT THROW SNARES

Order Number	Description	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)
M00562651	Jumbo Oval	30	2.4	240
M00562671	Standard Oval	27	2.4	240
M00562691	Micro Oval	13	2.4	240
M00562652	Jumbo Oval	30	2.4	240
M00562672	Standard Oval	27	2.4	240
M00562692	Micro Oval	13	2.4	240
M00562673	Standard Oval	27	2.4	240
M00562693	Micro Oval	13	2.4	240

SENSATION SINGLE-USE MEDIUM STIFF WIRE SHORT THROW SNARES

Order Number	Description	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)
M00560301	Jumbo Oval	30	2.4	240
M00560311	Standard Oval	27	2.4	240
M00560321	Micro Oval	13	2.4	240
M00560331	Crescent	27	2.4	240

CAPTIFLEX SINGLE-USE SNARES

Order Number	Description	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)
M00562401	Standard Oval	27	2.4	240
M00562421	Micro Oval	13	2.4	240
M00562471	Mini-Micro Oval	11	2.4	240
M00562402	Standard Oval	27	2.4	240
M00562422	Micro Oval	13	2.4	240

PROFILE SINGLE-USE PEDIATRIC SNARES

Order Number	Description	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)
M00562531	Mini Micro Oval	11	1.9	240
M00562551	Micro Oval	13	1.9	240
M00562571	Wide Oval	27	1.9	240

**Boston
Scientific**

Delivering what's next.™

Boston Scientific Corporation
One Boston Scientific Place
Natick, MA 01760-1537
www.bostonscientific.com

Ordering Information
1.800.225.3226

© 2005 Boston Scientific Corporation
or its affiliates. All rights reserved.

DVG1730 5M 9/05

CAUTION: Federal Law (USA) restricts these devices to sale by or on the order of a physician.
Refer to package insert provided with these products for complete Instructions for Use, Contraindications for Use, Potential Adverse Effects, Warnings and Precautions prior to using these products.