


Interject[®]

Sclerotherapy Needles

Boston
Scientific



Performance
you can count on



Interject®

Sclerotherapy Needles

Designed to provide **performance you can count on** for injection therapy applications.

TORTUOUS ANATOMY

Unique star-shaped inner catheter design helps promote reliable needle actuation in tortuous anatomy (even after kinking).

EFFICIENT ACCESS

Strong, pushable sheath construction facilitates advancement through difficult pathways to remote anatomy.

SMOOTH NEEDLE PENETRATION


Beveled needle designed to facilitate injection.

SECURE NEEDLE POSITIONING

Red Spacer Clip is designed to maintain needle position within the sheath during scope passage.

RANGE OF OPTIONS

- 200cm and 240cm sheath lengths
- 4mm and 6mm needle lengths
- 23 and 25 gauge needles
- Clear and contrast sheath available



Smooth beveled design facilitates injection

Packaged with new Red Spacer Clip which is designed to keep the needle inside the catheter until used



Convenient locking hub for procedural efficiency

ORDERING INFORMATION

Interject Single-Use Sclerotherapy Needles*

Order Number	O.D. (mm)	Sheath Length (cm)	Sheath Design	Needle Gauge	Maximum Needle Extension Length (mm)
M00518151	1.8	200	Contrast	23	4
M00518161	1.8	200	Contrast	25	4
M00518251	1.8	200	Contrast	23	6
M00518261	1.8	200	Contrast	25	6
M00518351	2.3	240	Contrast	23	4
M00518361	2.3	240	Contrast	25	4

Clear Interject Single-Use Sclerotherapy Needles*

Order Number	O.D. (mm)	Sheath Length (cm)	Sheath Design	Needle Gauge	Maximum Needle Extension Length (mm)
M00518111	1.8	200	Clear	25	4
M00518301	2.3	240	Clear	23	4
M00518311	2.3	240	Clear	25	4

*Packaged 5 Per Box

Indications, Contraindications, Warnings and Instructions for Use can be found in the product labeling supplied with each device.
CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Boston Scientific

Delivering what's next.™

Boston Scientific Corporation
One Boston Scientific Place
Natick, MA 01760-1537
www.bostonscientific.com

Ordering Information
1.800.225.3226

© 2007 Boston Scientific Corporation
or its affiliates. All rights reserved.

DVG1520 5M 1/07