


# *Resolution*<sup>®</sup> Clip Device

Boston  
Scientific



When positioning is critical,  
there is a ***Resolution***

# Resolution<sup>®</sup> Clip Device

Clips are being used today for a large number of indications, most of which are related to bleeding. They are designed for prompt hemostatic effect. Several studies have been published relating to the safety and effectiveness of endoscopic clip placement.

## RESOLUTION CLIP DEVICE IS INTENDED TO BE USED FOR:

- Hemostasis
- Endoscopic marking
- Closure
- Anchoring jejunal feeding tubes

*(see detailed information in the "Instructions for Use" supplied with the product)*



Resolution Clip is fully exposed and wide open, ready to grasp tissue.



Resolution Clip is closed, grasping tissue, and is now ready to be deployed.



Resolution Clip is deployed on tissue.

The image features a Resolution Clip Device, a long, thin, silver metal instrument with a blue handle and a red jaw. Two circular callouts provide close-up views of the clip's mechanism. The background is a solid teal color with faint, semi-transparent images of the device's handle and jaw.

### REOPENING CAPABILITY

The radiopaque Resolution Clip® is **engineered to enable opening and closing up to five times** prior to deployment, aiding in repositioning of the clip.\*

### PRE-LOADED

Resolution Clip Device is designed to be **ready-to-use**, which is essential for emergency bleeding situations.

### 11MM-WIDE JAW SPAN


Intended to **grasp a sizeable amount** of tissue.

### USER-FRIENDLY HANDLE

Designed to deliver a **familiar actuation** for opening and closing the jaws. Provides **tactile feedback** upon closing of jaws, indicating that the clip is about to be deployed.

\*Re-opening and closing capability may be limited by clinical circumstances and patient anatomy, among other factors.

# Resolution<sup>®</sup> Clip Device



## Resolution Clip Device

Order Number	Working Length (cm)	Minimal Working Channel (mm)	Clip Opening (mm)	Box
M00522600	155	2.8	11	1
M00522601	155	2.8	11	10
M00522602	155	2.8	11	20
M00522610	235	2.8	11	1
M00522611	235	2.8	11	10
M00522612	235	2.8	11	20

\*Endoclips versus heater probe in preventing early recurrent bleeding from peptic ulcer: a prospective and randomized trial. From Cipolletta, Naples, Italy, in *Gastrointestinal Endoscopy* 2001; 53:147-51

Dieulafoy's disease treated by endoscopic hemostatic clipping. Teo EK, Fock KM. *J Gastroenterology Hepatology* 1998; 13:320-1

<sup>†</sup>See detailed information in the "Instructions for Use" supplied with the product

**Boston  
Scientific**

*Delivering what's next.™*

Boston Scientific Corporation  
One Boston Scientific Place  
Natick, MA 01760-1537  
www.bostonscientific.com

**Ordering Information**  
**1.800.225.3226**

© 2006 Boston Scientific Corporation  
or its affiliates. All rights reserved.

BVG5030 revB 5M 3/06

**CAUTION:** Federal (USA) law and governing law outside the USA restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device.