


EndoVive®

Low Profile PEGs

Boston Scientific



Product Highlights

- Deluxe Procedure Kit
- Indicated for Adult and Pediatric Use
- Cost Effective
- Low Profile Design
- Anti-Reflux Valve
- Available in Both Pull and Push Techniques
- Universal-Sized Snare
- Multiple Sizes and Accessories

PRODUCT FEATURES


- | | |
|--|---|
| Cost Effective | ◦ A cost-effective solution to provide low profile enteral feeding in one initial procedure |
| High-Grade Silicone Button-Bolstered Low Profile PEG Tube | ◦ Designed for performance and consistent flow
◦ Radiopaque silicone designed to enhance placement confirmation |
| Low Profile Design | ◦ Low profile design helps reduce risk of accidental pull-out
◦ Product flexibility at skin level designed to enhance patient comfort
◦ Low profile devices designed to be discrete and provide greater cosmetic appeal |
| Anti-Reflux Valve | ◦ Designed to reduce gastric backflow, which may reduce tube clogging and extend the life of the implant
◦ Decompression accessory is designed to allow for periodic gastric decompression |
| Universal 1.9 mm O.D. x 240 cm Retrieval Snare | ◦ Designed to traverse scopes with a 2.0 mm or larger working channel
◦ Intended to eliminate the cost of a specialty snare |
| Multiple Sizes and Accessories | ◦ Color-coded French sizes with corresponding feeding accessories and multiple lengths for broad patient coverage
◦ Specially designed spacer discs to accommodate "in-between" sizes and short-term weight gain |

KIT CONTENTS

- | | | |
|--|---|--|
| (1) Three-Pack Iodine Swabs | (1) Optional Pediatric Bumper | (1) Press Fit Right Angle Feeding Set |
| (1) Fenestrated Drape | (1) 180 cm Insertion Wire/290 cm Guidewire | (1) Press Fit Bolus Straight Feeding Set |
| (1) 5.0 ml Syringe | (1) High-Grade Silicone Button-Bolstered Low Profile PEG Tube | (1) Press Fit Right Angle Decompression Tube |
| (1) 19 ga x 1½" Filter Needle | (4) 4" x 4" Gauze Pads | (1) 60 ml Catheter Tip Syringe |
| (1) 5.0 ml Xylocaine® 1% Ampule | (1) 5.0 g Lubricating Jelly | |
| (1) 25 ga x 1" Injection Needle | (1) Antibiotic Ointment | |
| (1) No. 11 Blade Surgical Scalpel | (4) 2" x 2" Drainage Sponges | |
| (1) 1.9 mm O.D./240 cm Retrieval Snare | (1) Spacer Disc Variety Pack | |
| (1) Percutaneous Stoma Measuring Device (PSMD) | | |

EndoVive®

Low Profile PEGs



Ordering Information

LOW PROFILE PEG KITS

Order Number	Description	Length (cm)	Color Code	Units
M00563120	18Fr Low Profile PEG Kit, Pull	1.2	Yellow	Box 1
M00563130	18Fr Low Profile PEG Kit, Pull	1.7	Yellow	Box 1
M00563140	18Fr Low Profile PEG Kit, Pull	2.4	Yellow	Box 1
M00563150	18Fr Low Profile PEG Kit, Pull	3.4	Yellow	Box 1
M00563160	18Fr Low Profile PEG Kit, Pull	4.4	Yellow	Box 1
M00563180	18Fr Low Profile PEG Kit, Push (Guidewire)	1.7	Yellow	Box 1
M00563190	18Fr Low Profile PEG Kit, Push (Guidewire)	2.4	Yellow	Box 1
M00563010	24Fr Low Profile PEG Kit, Pull	1.7	Blue	Box 1
M00563020	24Fr Low Profile PEG Kit, Pull	2.4	Blue	Box 1
M00563030	24Fr Low Profile PEG Kit, Pull	3.4	Blue	Box 1
M00563040	24Fr Low Profile PEG Kit, Pull	4.4	Blue	Box 1
M00563080	24Fr Low Profile PEG Kit, Push (Guidewire)	2.4	Blue	Box 1
M00563090	24Fr Low Profile PEG Kit, Push (Guidewire)	3.4	Blue	Box 1
M00563100	24Fr Low Profile PEG Kit, Push (Guidewire)	4.4	Blue	Box 1

= Latex-Free

Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device.
CAUTION: Federal (USA) law and governing law outside the USA restricts these devices to sale by or on the order of a physician.
Xylocaine is a trademark of APP Pharmaceuticals, LLC.

Boston Scientific

Delivering what's next.™

Boston Scientific Corporation
One Boston Scientific Place
Natick, MA 01760-1537
www.bostonscientific.com

Ordering Information
1.800.225.3226

© 2009 Boston Scientific Corporation
or its affiliates. All rights reserved.

DVG2380 500 April 2009