



Product Highlights


- Deluxe Procedure Kit
- Available in Both Pull and Push Techniques
- Universal-Sized Snare
- Multiple Sizes and Accessories

PRODUCT FEATURES

- | | |
|---|--|
| High-Grade Silicone Dome-Bolstered PEG Tube | <ul style="list-style-type: none"> ◦ Designed for performance and consistent flow ◦ Facilitates traction removal ◦ Universal 1.9mm O.D. x 240 cm retrieval snare |
| Universal 1.9 mm O.D. x 240 cm Retrieval Snare | <ul style="list-style-type: none"> ◦ Traverses scopes with a 2.0 mm or larger working channel ◦ Intended to eliminate the cost of a specialty snare |
| Double-Tier Tray Packaging | <ul style="list-style-type: none"> ◦ Provides kit organization and promotes procedural efficiency |
| External Bolsters | <ul style="list-style-type: none"> ◦ One external bolster provided to accommodate physician preference and patient needs:
Round Bolster |
| Multiple Sizes and Accessories | <ul style="list-style-type: none"> ◦ Multiple French sizes for broad patient coverage ◦ Separate medication and feeding ports ◦ C-Clamp designed to help reduce risk of tube reflux and leaking |

KIT CONTENTS

- | | | |
|--|----------------------------------|--------------------------------|
| (1) Three-Pack Iodine Swabs | (1) Trocar Catheter | (1) Curved Hemostat |
| (1) Fenestrated Drape | (1) Seldinger Needle (Push only) | (2) External Bolsters |
| (1) 5.0 ml Syringe | (1) 180 cm Insertion Wire (Pull) | (1) C-Clamp |
| (1) 5.0 ml Xylocaine® 1% Ampule | (1) 290 cm Guidewire (Push) | (1) Y-Port |
| (1) 19 ga x 1-1½" Filter Needle | (1) High-Grade Silicone PEG Tube | (1) Packet Antiseptic Ointment |
| (1) 22 ga x 1½" Injection Needle | (4) 4" x 4" Gauze Pads | (4) 2" x 2" Drainage Sponges |
| (1) No. 11 Blade Surgical Scalpel | (1) 5.0 g Lubricating Jelly | |
| (1) 1.9 mm O.D./240 cm Retrieval Snare | (1) Scissors | |



Ordering Information

STANDARD PEG KITS

Order Number	Description	Units
M00568201	20Fr Standard PEG Kit, Pull	Box 2
M00568211	20Fr Standard PEG Kit, Push	Box 2
M00568241	24Fr Standard PEG Kit, Pull	Box 2
M00568251	24Fr Standard PEG Kit, Push	Box 2

**Boston
Scientific**

Delivering what's next.™

Boston Scientific Corporation
One Boston Scientific Place
Natick, MA 01760-1537
www.bostonscientific.com

Ordering Information
1.800.225.3226

© 2009 Boston Scientific Corporation
or its affiliates. All rights reserved.

DVG2410 500 April 2009

Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device.
CAUTION: Federal (USA) law and governing law outside the USA restricts these devices to sale by or on the order of a physician.

Xylocaine is a trademark of APP Pharmaceuticals, LLP.