

Radial Jaw[®] 4 Biopsy Forceps

The next generation of our market-leading single-use biopsy forceps

RADIAL JAW 4 SINGLE-USE JUMBO BIOPSY FORCEPS

Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Color Code	Unit
M00513343	Radial Jaw 4 Jumbo	2.8	160	3.2	Yellow	Box 5
M00513353	Radial Jaw 4 Jumbo w/ Needle	2.8	160	3.2	Yellow	Box 5
M00513363	Radial Jaw 4 Jumbo	2.8	240	3.2	Orange	Box 5
M00513373	Radial Jaw 4 Jumbo w/ Needle	2.8	240	3.2	Orange	Box 5
M00513341	Radial Jaw 4 Jumbo	2.8	160	3.2	Yellow	Box 20
M00513351	Radial Jaw 4 Jumbo w/ Needle	2.8	160	3.2	Yellow	Box 20
M00513361	Radial Jaw 4 Jumbo	2.8	240	3.2	Orange	Box 20
M00513371	Radial Jaw 4 Jumbo w/ Needle	2.8	240	3.2	Orange	Box 20
M00513342	Radial Jaw 4 Jumbo	2.8	240	3.2	Yellow	Box 40
M00513352	Radial Jaw 4 Jumbo w/ Needle	2.8	240	3.2	Yellow	Box 40
M00513362	Radial Jaw 4 Jumbo	2.8	240	3.2	Orange	Box 40
M00513372	Radial Jaw 4 Jumbo w/ Needle	2.8	240	3.2	Orange	Box 40

RADIAL JAW 4 SINGLE-USE LARGE CAPACITY BIOPSY FORCEPS

Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Color Code	Unit
M00513303	Radial Jaw 4 Large Capacity	2.4	160	2.8	Yellow	Box 5
M00513313	Radial Jaw 4 Large Capacity w/ Needle	2.4	160	2.8	Yellow	Box 5
M00513323	Radial Jaw 4 Large Capacity	2.4	240	2.8	Orange	Box 5
M00513333	Radial Jaw 4 Large Capacity w/ Needle	2.4	240	2.8	Orange	Box 5
M00513301	Radial Jaw 4 Large Capacity	2.4	160	2.8	Yellow	Box 20
M00513311	Radial Jaw 4 Large Capacity w/ Needle	2.4	160	2.8	Yellow	Box 20
M00513321	Radial Jaw 4 Large Capacity	2.4	240	2.8	Orange	Box 20
M00513331	Radial Jaw 4 Large Capacity w/ Needle	2.4	240	2.8	Orange	Box 20
M00513302	Radial Jaw 4 Large Capacity	2.4	160	2.8	Yellow	Box 40
M00513312	Radial Jaw 4 Large Capacity w/ Needle	2.4	160	2.8	Yellow	Box 40
M00513322	Radial Jaw 4 Large Capacity	2.4	240	2.8	Orange	Box 40
M00513332	Radial Jaw 4 Large Capacity w/ Needle	2.4	240	2.8	Orange	Box 40

RADIAL JAW 3 SINGLE-USE STANDARD CAPACITY BIOPSY FORCEPS

Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Color Code	Unit
M00515341	Radial Jaw 3 Standard Capacity	2.2	160	2.8	Yellow	Box 5
M00515351	Radial Jaw 3 Standard Capacity w/ Needle	2.2	160	2.8	Yellow	Box 5
M00515361	Radial Jaw 3 Standard Capacity	2.2	240	2.8	Orange	Box 5
M00515371	Radial Jaw 3 Standard Capacity w/ Needle	2.2	240	2.8	Orange	Box 5
M00515342	Radial Jaw 3 Standard Capacity	2.2	160	2.8	Yellow	Box 20
M00515352	Radial Jaw 3 Standard Capacity w/ Needle	2.2	160	2.8	Yellow	Box 20
M00515362	Radial Jaw 3 Standard Capacity	2.2	240	2.8	Orange	Box 20
M00515372	Radial Jaw 3 Standard Capacity w/ Needle	2.2	240	2.8	Orange	Box 20
M00515373	Radial Jaw 3 Standard Capacity w/ Needle	2.2	240	2.8	Orange	Box 40

RADIAL JAW 3 SINGLE-USE LARGE CAPACITY BIOPSY FORCEPS

Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Color Code	Unit
M00515961	Radial Jaw 3 Large Capacity	2.2	160	2.8	Yellow	Box 5
M00515971	Radial Jaw 3 Large Capacity with Needle	2.2	160	2.8	Yellow	Box 5
M00515981	Radial Jaw 3 Large Capacity	2.2	240	2.8	Orange	Box 5
M00515991	Radial Jaw 3 Large Capacity with Needle	2.2	240	2.8	Orange	Box 5
M00515962	Radial Jaw 3 Large Capacity	2.2	160	2.8	Yellow	Box 20
M00515972	Radial Jaw 3 Large Capacity with Needle	2.2	160	2.8	Yellow	Box 20
M00515982	Radial Jaw 3 Large Capacity	2.2	240	2.8	Orange	Box 20
M00515983	Radial Jaw 3 Large Capacity	2.2	240	2.8	Orange	Box 40
M00515992	Radial Jaw 3 Large Capacity with Needle	2.2	240	2.8	Orange	Box 20
M00515993	Radial Jaw 3 Large Capacity with Needle	2.2	240	2.8	Orange	Box 40

Radial Jaw[®] 4 Biopsy Forceps

The next generation of our market-leading single-use biopsy forceps

RADIAL JAW 3 SINGLE-USE MAX CAPACITY BIOPSY FORCEPS

Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Color Code	Unit
M00515861	Radial Jaw 3 Max Capacity	3.3	160	3.7	Yellow	Box 5
M00515871	Radial Jaw 3 Max Capacity w/ Needle	3.3	160	3.7	Yellow	Box 5
M00515881	Radial Jaw 3 Max Capacity	3.3	240	3.7	Orange	Box 5
M00515891	Radial Jaw 3 Max Capacity w/ Needle	3.3	240	3.7	Orange	Box 5
M00515892	Radial Jaw 3 Max Capacity w/ Needle	3.3	240	3.7	Orange	Box 20

RADIAL JAW 3 SINGLE-USE GASTROPEDIATRIC BIOPSY FORCEPS

Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Color Code	Unit
M00515781	Radial Jaw 3 GP GastroPediatric	1.8	160	2.0	Yellow	Box 5
M00515791	Radial Jaw 3 GP GastroPediatric w/ Needle	1.8	160	2.0	Yellow	Box 5

INSTRUCTIONS FOR USE

Refer to the operator's manual for complete instructions for use.

INDICATIONS

These single-use biopsy forceps are specifically designed to collect tissue endoscopically for histologic examination. These forceps should not be used for purposes other than their intended function.

CONTRAINDICATIONS

None known.

WARNINGS

These single-use biopsy forceps should only be used to biopsy tissue where possible hemorrhage will not present a danger for the patients. Adequate plans for management of potential bleeding and appropriate airway management should be in place.

Contents supplied STERILE using a radiation process. Do not use if sterile barrier is damaged. If damage is found call your Boston Scientific representative.

For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

POTENTIAL ADVERSE EFFECTS

Complications associated with the use of the single-use biopsy forceps may include: bleeding, perforation, and infection.

Please be aware that potential adverse effects may arise even with the proper use of medical devices. Accordingly, this device should only be used by persons qualified in the procedures for which it is indicated.

CAUTIONS

Cautions can be found in the product labeling supplied with each device. CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.



Boston Scientific Endoscopy
Tel 508.650.8000
www.bostonscientific.com

Ordering Information
1.800.225.3226

©2008 Boston Scientific Corporation or its affiliates.
All rights reserved.