

Rotatable Snares

Polypectomy Snares

SINGLE-USE ROTATABLE SNARES*

Order Number	Description (Box 10)	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Unit
M00561821	Micro Oval	13	2.4	195	Box 10
M00561831	Mini-Standard Oval	20	2.4	195	Box 10

*With Star Catheter Technology

INSTRUCTIONS FOR USE

Refer to the operator's manual for complete instructions for use.

INDICATIONS

Indicated for the electrosurgical removal and cauterization of gastrointestinal tract polyps through an endoscope.

CONTRAINDICATIONS

Contraindications for these devices are for those specific to endoscopic Polypectomy.

WARNINGS

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

Ensure that the patient is properly grounded.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

POTENTIAL ADVERSE EVENTS

Potential adverse events include perforation, fulguration, hemorrhages, and transmural burns.

Please be aware that potential adverse events may arise even with the proper use of medical devices. Accordingly, this device should only be used by persons qualified in the procedures for which it is indicated.

CAUTIONS

A thorough understanding of the technical principles, clinical applications, and risks associated with Monopolar Polypectomy is necessary before using this product.

Cautions can be found in the product labeling supplied with each device. CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

**Boston
Scientific**

Delivering what's next.™

Tel 508.650.8000
www.bostonscientific.com

Ordering Information
1.800.225.3226

©2008 Boston Scientific Corporation or its affiliates. All rights reserved.