

Sensation™ Short Throw Snares

Polypectomy Snares

SENSATION SINGLE-USE SHORT THROW SNARES

Order Number	Description	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Unit
M00562651	Jumbo Oval	30	2.4	240	Box 5
M00562671	Standard Oval	27	2.4	240	Box 5
M00562691	Micro Oval	13	2.4	240	Box 5
M00562652	Jumbo Oval	30	2.4	240	Box 20
M00562672	Standard Oval	27	2.4	240	Box 20
M00562692	Micro Oval	13	2.4	240	Box 20
M00562673	Standard Oval	27	2.4	240	Box 40
M00562693	Micro Oval	13	2.4	240	Box 40

SENSATION SINGLE-USE MEDIUM STIFF WIRE SHORT THROW SNARES

Order Number	Description	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Unit
M00560301	Jumbo Oval	30	2.4	240	Box 5
M00560311	Standard Oval	27	2.4	240	Box 5
M00560321	Micro Oval	13	2.4	240	Box 5
M00560331	Crescent	27	2.4	240	Box 5

INSTRUCTIONS FOR USE

Refer to the operator's manual for complete instructions for use.

INDICATIONS

Indicated for use in the removal and cauterization of diminutive, sessile, and pedunculated polyps.

CONTRAINDICATIONS

Contraindicated for those specific to endoscopic polypectomy.

WARNINGS

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

POTENTIAL ADVERSE EVENTS

Potential adverse events include perforation, fulguration, hemorrhages, and transmural burns.

Please be aware that potential adverse events may arise even with the proper use of medical devices. Accordingly, this device should only be used by persons qualified in the procedures for which it is indicated.

CAUTIONS

A thorough understanding of the technical principles, clinical applications, and risks associated with Monopolar Polypectomy is necessary before using this product

Cautions can be found in the product labeling supplied with each device. CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Boston Scientific

Boston Scientific Endoscopy
Tel 508.650.8000
www.bostonscientific.com

Ordering Information
1.800.225.3226

©2008 Boston Scientific Corporation or its affiliates. All rights reserved.