

# Capio<sup>®</sup> Suture Capturing Devices

*Reach, Throw and Capture: One Step. One Device.*

## CAPIO<sup>®</sup> SUTURE CAPTURING DEVICES

Order Number	Description	Size
M0068312321.....	Standard - for direct vision endoscopic, open or transvaginal use.....	12mm x 32cm.....
M0068311251.....	Open Access - for direct vision sites where a trocar is not.....	25cm length.....

Packaged Four (4) Per Box

## CAPIO DEVICE SUTURES

Order Number	Description	Size	Closest Equivalent
M0068331131.....	Non-absorbable, coated braided polyester, double armed with TC tapercut needles, 48"	0.....	TI-CRON <sup>®</sup> / Ethibond <sup>®</sup>
M0068331231.....	Non-absorbable, polypropylene monofilament, double armed with TC tapercut needles, 48"	0.....	Surgilene <sup>®</sup> / Prolene <sup>®</sup>
M0068331371*.....	Absorbable, monofilament PDO, double armed with TC tapercut needle and a T 26mm 1/2 circle taper needle, 48"	0.....	PDS <sup>®</sup> II
M0068332131.....	Absorbable, coated braided PGA, double armed with TC tapercut needles, 48" ...	0.....	Dexon II <sup>®</sup> / Vicryl <sup>®</sup>
M0068331141.....	Non-absorbable, coated braided polyester, double armed with TC tapercut needle and a T 26mm 1/2 circle taper needle, 36"	0.....	TI-CRON / Ethibond
M0068331241.....	Non-absorbable, polypropylene monofilament, double armed with TC tapercut needle and a T 26mm 1/2 circle taper needle, 36"	0.....	Surgilene / Prolene

Packaged Twelve (12) Per Box

\*Please check with your Boston Scientific representative for availability.

### INSTRUCTIONS FOR USE

Refer to the operator's manual for complete instructions for use.

### INDICATIONS

The Capio Suture Capturing Device is intended for use in general suturing applications during open and endoscopic surgery to assist in the placement of suture material in tissues at the operative site and is to be used under direct visualization only.

### CONTRAINDICATIONS

The use of the Capio Suture Capturing Device is contraindicated as follows: When, in the judgement of the physician, open or endoscopic procedures would be contrary to the best interests of the patient; and for placing sutures into or through bone.

### WARNINGS

For single use only. Do not reuse, reprocess or resterilize. Use only needles and sutures specified for use with the Capio Suture Capturing Device.

### POTENTIAL ADVERSE EFFECTS

The potential adverse effects which may result from suture placement procedure include, but are not limited to, are injury to internal vessels and nerves, bleeding, hematoma, inflammatory reaction to the suture material or to the trauma of being sutured.

### CAUTIONS

Federal (USA) law and governing law outside the USA restrict these devices to sale by or on the order of a physician.

### TRADEMARKS

Capio is a trademark of Boston Scientific Corporation or its affiliates.  
Ethibond, Prolene, PDS and Vicryl are trademarks of Johnson & Johnson Corporation.  
Surgilene is a trademark of American Cyanamid Company.  
Dexon is a trademark of Sherwood Services AG.  
TI-CRON is a trademark of Tyco Group S.A.R.L.

**Boston Scientific**  
*Delivering what's next.™*

Boston Scientific  
Tel 508.650.8000  
www.bostonscientific.com

Ordering Information  
1.888.272.1001

©2005 Boston Scientific Corporation  
or its affiliates. All rights reserved.