

# *Radial Jaw<sup>®</sup> and Radial Edge<sup>®</sup>*

Single-Use Biopsy Forceps

Boston  
Scientific

Pulmonary  
Endoscopy



# Radial Jaw® and Radial Edge®

## Single-Use Biopsy Forceps

### RADIAL JAW 3 SINGLE-USE BIOPSY FORCEPS (Box 5)

Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Required Working Channel (mm)
M00515231	Radial Jaw 3	1.8	100	2.0
M00515241	Radial Jaw 3 with needle	1.8	100	2.0
M00515211	Radial Jaw 3 Large Capacity	2.2	100	2.8

### RADIAL JAW SINGLE-USE BIOPSY FORCEPS (Box 5)

Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Required Working Channel (mm)
M00512661	Radial Jaw	1.8	100	2.0
M00512691	Radial Jaw with needle	1.8	100	2.0

### RADIAL EDGE SINGLE-USE BIOPSY FORCEPS (Box 5)

Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Required Working Channel (mm)
M00515561	Radial Edge	1.8	100	2.0
M00515591	Radial Edge with needle	1.8	100	2.0

## Indications

- Designed to collect tissue endoscopically for histologic examination

## Single-Use Design

- Eliminates the need for cleaning and resterilization
- Provides first-time sharpness and helps promote patient safety

## Ergonomic Handle

- Designed to provide enhanced comfort and “feel” for end-user

## PolyGlide™ and EndoGlide™ Coatings

- Helps reduce friction to facilitate passage through a bronchoscope

## Options

- Available in jaw O.D. of 2.2mm or 1.8mm
- Available with or without needle

## RADIAL JAW BIOPSY FORCEPS

- Proprietary full-surface, micro-mesh teeth
- Designed to facilitate clean, precise resection
- Pivoting head is designed to facilitate tangential tissue acquisition

## RADIAL EDGE BIOPSY FORCEPS

- Smooth, bilaterally sharpened cup edges.
- Designed to facilitate clean resection with less risk of trauma.



Radial Jaw 3 Biopsy Forceps  
with EndoGlide Coating



Radial Jaw Biopsy Forceps  
with PolyGlide Coating



Radial Edge Biopsy Forceps  
with PolyGlide Coating

# Boston Scientific

Delivering what's next.™

Boston Scientific Corporation  
One Boston Scientific Place  
Natick, MA 01760-1537  
www.bostonscientific.com/pulmonary

**Ordering Information**  
**1.800.225.3226**

© 2005 Boston Scientific Corporation  
or its affiliates. All rights reserved.

DVG1240 5M 2/05

Indications, Contraindications, Warnings and Instructions for Use can be found in the product labeling supplied with each device.

**CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.**