

Easy Core™

Biopsy Device

Boston
Scientific



Generous notch size is designed to improve volume of sample to facilitate accurate diagnosis

Underside grip designed for comfort and security facilitates multiple grip options

Yellow caution indicator for confirmation of device arming

Needle gauge identification for proper device selection



Easy Core™

Biopsy Device

ERGONOMIC HANDLE

- Lightweight, symmetric design enhances dexterity
- Single arming thumb tab simplifies device preparation
- Thumb tab is designed to remain stationary during firing to enhance stability
- Dual firing option reduces need to alter hand position for firing
- Yellow caution indicator at handle base provides visual confirmation that device is ready for sampling




Side Firing Button




Rear Firing Button

PRECISION-GROUND, ELECTROPOLISHED STYLET AND CANNULA

- Delta Cut Needle is designed to facilitate acquisition of clean, uncrushed specimens
- Echogenic cannula tip designed to improve visibility during ultrasound guidance
- Increased notch length designed to improve sample size



Electropolished stylet and cutting cannula are designed to provide the sharpness required for clean, uncrushed specimen



ORDERING INFORMATION

ORDER NUMBER	DESCRIPTION	GAUGE	WORKING LENGTH	QUANTITY
M0065001041	Easy Core Biopsy Device	18g	21cm	5
M0065001051	Easy Core Biopsy Device	18g	15cm	5
M0065001061	Easy Core Biopsy Device	15g	15cm	5
M0065001071	Easy Core Biopsy Device	18g	25cm	5

Max-Core is a trademark of C.R. Bard Corporation.
Tru-Cut is a trademark of Allegiance Corp.

Refer to User Manual included with this product for complete instructions for use.
Caution: Federal Law (USA) restricts these devices to sale by or on the order of a physician.

*Data on file, Boston Scientific Corporation

Bench test results may not necessarily be indicative of clinical performance.

Boston Scientific

Delivering what's next.™

Boston Scientific Corporation
One Boston Scientific Place
Natick, MA 01760-1537
www.bostonscientific.com/urology

Ordering Information
1.888.272.1001

© 2005 Boston Scientific Corporation
or its affiliates. All rights reserved.

DVU1360 5M 01/05-01/07