

SlimLine EZ[™]

Single-Use Holmium Laser Fiber

Boston
Scientific

Custom Optical Fit
For Optimum
Performance



SlimLine EZ™

Single-Use Holmium Laser Fiber

New SlimLine EZ Single-Use Holmium Laser Fibers are designed and manufactured by Lumenis to offer optimum performance on Lumenis Holmium Lasers. The Single-Use configurations of the original SlimLine fibers offer new fiber performance for every case without the hassle of reprocessing. SlimLine EZ Fibers are designed to deliver a customized optical fit to Lumenis Holmium Lasers and are designed to reduce scope and console damage.

SINGLE-USE CONFIGURATION


Single use fibers ensure best performance of the laser fiber for every case. No fiber maintenance or reprocessing is required between cases with the SlimLine EZ Single-Use Holmium Laser Fibers.

OPTICAL CONNECTION

All Lumenis Laser Fibers are designed to match exactly the laser beam "spot size," the focal point distance, and the angle of focus of Lumenis Laser Systems. These matching measurements give Lumenis Laser Fibers a custom optical fit to Lumenis Lasers.

REDUCE SCOPE DAMAGE

Because of the custom optical fit and a new, single-use fiber for every case, SlimLine EZ Fibers are designed to be less likely to cause scope or laser console damage than reusable fibers.



ORDERING INFORMATION

| PRODUCT CODE | DESCRIPTION | UNITS |
|--|---|-------|
| SLIMLINE EZ SINGLE-USE HOLMIUM LASER FIBERS | | |
| M0068408921 | 200µ Single-Use Holmium Laser Fiber | Box/3 |
| M0068408931 | 365µ Single-Use Holmium Laser Fiber | Box/3 |
| M0068408941 | 550µ Single-Use Holmium Laser Fiber | Box/3 |
| SLIMLINE REUSABLE HOLMIUM LASER FIBERS | | |
| M0068408400 | 200µ Reusable Holmium Laser Fiber | Each |
| M0068408420 | 365µ Reusable Holmium Laser Fiber | Each |
| M0068408430 | 550µ Reusable Holmium Laser Fiber | Each |
| M0068408440 | 1000µ Reusable Holmium Laser Fiber | Each |
| SIDE FIRING LASER FIBERS | | |
| M0068408460 | 550µ DuoTome Single-Use Side Firing Laser Fiber | Each |

Boston Scientific

Delivering what's next.™

Boston Scientific Corporation
One Boston Scientific Place
Natick, MA 01760
www.bostonscientific.com/urology

Ordering Information
1.888.272.1001

© 2005 Boston Scientific Corporation
or its affiliates. All rights reserved.

BVU2240 7.5M 06/05-06/07



SlimLine and SlimLine EZ are trademarks of Lumenis Ltd.
SlimLine EZ is manufactured by Lumenis and distributed by Boston Scientific.

Caution: Federal Law (USA) restricts these devices to sale by or on the order of a physician.
Refer to package insert provided with product for complete instructions for use including Indications for Use, Contraindications, Potential Complications, Warning and Cautions.