


TruPath™

Biopsy Device

Boston
Scientific



Generous notch size is designed to improve volume of sample to facilitate accurate diagnosis

Yellow caution indicator for confirmation of device arming

Underside grip for comfort and security facilitates multiple grip options

TruPath
18G

Needle gauge identification for proper device selection

ERGONOMIC HANDLE

- Lightweight, symmetric design enhances dexterity
- Single arming thumb tab simplifies device preparation
- Thumb tab is designed to remain stationary during firing to enhance stability
- Dual firing option reduces need to alter hand position for firing
- Yellow caution indicator at handle base provides visual confirmation that device is ready for sampling




Side Firing Button



Rear Firing Button

PRECISION-GROUND, ELECTROPOLISHED STYLET AND CANNULA

- Designed to facilitate acquisition of clean, uncrushed specimens
- Echogenic cannula tip designed to improve visibility during ultrasound guidance
- Increased notch length designed to improve sample size



ORDERING INFORMATION

ORDER NUMBER	DESCRIPTION	GAUGE	WORKING LENGTH	QUANTITY
M0065001081	TruPath Biopsy Device	18g	25cm	10
M0065001141	TruPath Biopsy Device	18g	15cm	10
M0065001151	TruPath Biopsy Device	18g	21cm	10

Refer to User Manual included with this product for complete instructions for use.
Caution: Federal Law (USA) restricts these devices to sale by or on the order of a physician.

*Data on file, Boston Scientific Corporation

Bench test results may not necessarily be indicative of clinical performance.

Boston Scientific

Delivering what's next.™

Boston Scientific Corporation
 One Boston Scientific Place
 Natick, MA 01760-1537
www.bostonscientific.com/urology

Ordering Information
1.888.272.1001

© 2005 Boston Scientific Corporation or its affiliates. All rights reserved.

DVU1350 5M 01/05-01/07