


ImageReady™ Supported Devices

PPM and CRT-P Systems

ImageReady™ MR Conditional Pacing System Components for 1.5T¹ and 3T¹

Component	Model Number(s)	1.5T	3T
Pacemaker Pulse Generators			
ADVANTIO MRI Pulse Generator	K085, K086, K087	✓	
INGENIO MRI Pulse Generator	K185, K186, K187	✓	
VITALIO MRI Pulse Generator	K285, K286, K287	✓	
FORMIO MRI Pulse Generator	K289	✓	
ACCOLADE MRI Pulse Generator	L310, L311, L331	✓	✓
CRT-P Pulse Generators			
VISIONIST X4 Pulse Generator	U228	✓	✓
Leads and Accessories			
Right Atrial and Right Ventricular Leads and Accessories			
FINELINE II Sterox Pacing Leads	4456, 4457, 4458, 4459, 4479, 4480	✓	✓
FINELINE II Sterox EZ Pacing Leads	4469, 4470, 4471, 4472, 4473, 4474	✓	✓
Suture Sleeve for FINELINE II Leads	6220, 6221	✓	✓
INGEVITY MRI Pacing Leads (Tined Fixation)	7731, 7732, 7735, 7736	✓	✓
INGEVITY MRI Pacing Leads (Extendable / Retractable Fixation)	7740, 7741, 7742	✓	✓
INGEVITY + Pacing Leads (Extendable / Retractable Fixation)	7840, 7841, 7842	✓	✓
Suture Sleeve for INGEVITY MRI Leads	6402	✓	✓
IS-1 Lead Port Plug	7145	✓	✓
Left Ventricular Leads and Accessories			
ACUITY X4 (IS4) Pacing Leads	4671, 4672, 4674, 4675, 4677, 4678	✓	✓
Suture Sleeve for ACUITY X4 Leads	4603	✓	✓
IS4 Lead Port Plug	7148	✓	✓

1. For conditions of use and scan conditions, refer to ImageReady Conditional Pacing System Technical Guide 359489-002 EN Australia 2018-08



ImageReady™ Supported Devices

ICD, CRT-D & S-ICD Systems

ImageReady™ MR Conditional Defibrillation System Components for 1.5T¹ and 3T¹

Component	Model Number(s)	1.5T	3T
CRT-D Pulse Generators			
AUTOGEN CRT-D	G172, G173, G177, G179	✓	
DYNAGEN CRT-D	G150, G151, G156, G158	✓	
MOMENTUM CRT-D	G124, G125, G128, G138	✓	
RESONATE X4 CRT-D	G447	✓	✓
ICD Pulse Generators			
AUTOGEN EL ICD	D174, D175, D176, D177	✓	
AUTOGEN MINI ICD	D044, D045, D046, D047	✓	
DYNAGEN EL ICD	D150, D151, D152, D153	✓	
DYNAGEN MINI ICD	D020, D021, D022, D023	✓	
MOMENTUM EL ICD	D120, D121	✓	
RESONATE EL ICD	D432, D433	✓	✓
S-ICD Pulse Generators			
EMBLEM™, EMBLEM™ MRI S-ICD	A209, A219	✓	
Leads and Accesories			
Right Atrial Leads and Accessories			
FINELINE II Sterox Pacing Leads	4479, 4480	✓	✓
FINELINE II Sterox Pacing Leads	4469, 4470, 4471, 4472, 4473, 4474	✓	✓
Suture Sleeves for FINELINE II Leads	6220, 6221	✓	✓
INGEVITY MRI Pacing Leads (Tined Fixation)	7735, 7736	✓	✓
INGEVITY MRI Pacing Leads (Extendable / Retractable Fixation)	7740, 7741, 7742	✓	✓
INGEVITY + Pacing Leads (Extendable / Retractable Fixation)	7840, 7841, 7842	✓	✓
Suture Sleeve for INGEVITY MRI Leads	6402	✓	✓
IS-1 Lead Port Plug	7145	✓	✓

ImageReady™ Supported Devices

ICD, CRT-D & S-ICD Systems

ImageReady™ MR Conditional Defibrillation System Components for 1.5T¹ and 3T¹

Component	Model Number(s)	1.5T	3T
Right Ventricular Leads and Accessories			
ENDOTAK RELIANCE (IS-1) Leads	0170, 0171, 0172, 0173, 0174, 0175, 0176, 0177, 0180, 0181, 0182, 0183, 0184, 0185, 0186, 0187	✓	
ENDOTAK RELIANCE (DF4) Defibrillation Leads	0262, 0263, 0265, 0266, 0272, 0273, 0275, 0276, 0282, 0283, 0285, 0286, 0292, 0293, 0295, 0296	✓	✓
RELIANCE 4-FRONT (DF4) Defibrillation Leads	0636, 0650, 0651, 0652, 0653, 0654, 0655, 0657, 0658, 0662, 0663, 0665, 0672, 0673, 0675, 0676, 0682, 0683, 0685, 0686, 0692, 0693, 0695, 0696	✓	✓
Suture Sleeve for RELIANCE 4-FRONT Leads	6403	✓	✓
Left Ventricular Leads and Accessories			
ACUITY Spiral Leads	4591, 4592, 4593	✓	
ACUITY X4 (IS4) Pacing Leads	4671, 4672, 4674, 4675, 4677, 4678	✓	✓
Suture Sleeve for ACUITY X4 Leads	4603	✓	✓
IS4 / DF4 Lead Port Plug	7148	✓	✓
IS-1 Lead Port Plug	7145	✓	
S-ICD Electrodes and Accessories			
EMBLEM S-ICD Electrode	All Models	✓	
Cameron Health Q-TRAK S-ICD Electrode	All Models	✓	
S-ICD Electrode Suture Sleeve	All Models	✓	

1. For conditions of use and scan conditions, refer to MRI Technical Guide ImageReady MR Conditional Defibrillation System 92138524-002 EN Australia 2019-11 and ImageReady MR Conditional S-ICD System 359475-019 EN Europe 2020-01

