

# New Cardio Product Catalogue 2020





# TABLE OF CONTENTS

## 1 STENTS & DRUG TECHNOLOGIES

### Abluminal Bioabsorbable Polymer DES

- SYNERGY™ \*\* ..... 7-8
- SYNERGY MEGATRON™ \*\* ..... 8

### Drug-Eluting Stents

- Promus ELITE™ ..... 10
- Promus PREMIER™ ..... 11

### Bare Metal Stent

- REBEL™ ..... 13

### Drug Coated Balloon

- AGENT™ ..... 15

## 2 COMPLEX PCI

### Balloon Catheters

#### Pre-dilatation Balloon Catheters

- EMERGE™ Monorail™ ..... 18-19
- EMERGE™ OTW ..... 19-20
- MAVERICK2™ Monorail™ ..... 21

#### Post-dilatation Balloon Catheters

- NC QUANTUM APEX™ ..... 22-23
- NC EMERGE™ Monorail™ ..... 24-25

### Microdilatation Catheters

- THREADER™ Microdilatation Catheter ..... 27
- MAMBA™ Microdilatation Catheter ..... 27

### Guide Catheters

- Curves Overview ..... 29
- CONVEY™ 5F / 6F / 7F / 8F ..... 30-33
- MACH 1™ 6F / 7F / 8F ..... 34-38
- RUNWAY™ 6F ..... 39-41

### Guide Catheter Extension

- GUIDEZILLA™ II Guide Catheter Extension ..... 43

### Angiographic Catheters

- IMPULSE™ 5F / 6F ..... 45-48
- EXPO™ 5F / 6F / 7F ..... 49-51

### Guidewires

#### Standard Lesions / Spring Tip Guidewires

- CHOICE™ Floppy ..... 53
- LUGE™ Moderate Support Guidewire ..... 53

#### Complex Lesions / Polymer Tip Guidewires

- PT2™ LS ..... 54
- CHOICE™ PT Floppy ..... 54
- PT2™ MS ..... 54

#### Stiff Guidewires

- CHOICE™ Intermediate ..... 55
- CHOICE™ Standard ..... 55
- PT GRAPHIX™ Intermediate ..... 56
- PT GRAPHIX™ Standard ..... 56

#### Extra Support Guidewires

- CHOICE™ Extra Support ..... 57
- CHOICE™ PT Extra Support ..... 57
- MAILMAN™ ..... 58
- PT GRAPHIX™ Super Support ..... 58

#### SENTAI Workhorse Wires

- MARVEL™ Guidewire ..... 59
- SAMURAI™ Guidewire ..... 59

#### SENTAI Specialty Crossing Guidewires

- SAMURAI™ RC Guidewire ..... 60
- FIGHTER™ Guidewire ..... 60
- JUDO™ 1 Guidewire ..... 60
- JUDO™ 3 Guidewire ..... 60
- JUDO™ 6 Guidewire ..... 60
- HORNET™ 10 Guidewire ..... 61
- HORNET™ 14 Guidewire ..... 61

#### Extension Guidewires

- ADDWIRE™ Extension Wire ..... 62
- The 38 MAGNET™ Device ..... 62
- STRETCH Extension Wire ..... 62

### Embolic Protection

- FILTERWIRE EZ™ Embolic Protection System ..... 64

### Accessories

#### Inflation Devices

- ENCORE™ 26 Inflation Device ..... 66
- ENCORE™ 26 Advantage Kit ..... 66

#### Y-Adaptors

- GATEWAY™ Plus Y-Adaptor ..... 67

#### Hemostasis Valve

- WATCHDOG™ Hemostasis Valve ..... 68

#### Introducer Sheaths

- SUPER SHEATH™ 0.889 mm (0.035") compatible ..... 69
- SUPER SHEATH™ 0.965 mm (0.038") compatible ..... 69

### Cutting Balloon™

- FLEXTOME™ Cutting Balloon™ ..... 71
- WOLVERINE™ Cutting Balloon™ Dilatation Device ..... 72

### Atherectomy Devices

- ROTAPRO Rotational Atherectomy System ..... 74
- ROTAPRO Pre-Connected Burr and Advancing Device ..... 74
- ROTAWIRE™ ..... 74
- ROTALINK™ Advancer ..... 74
- ROTALINK™ Exchangeable Burr Catheter ..... 75
- ROTALINK™ Plus (Advancer & Burr) ..... 75
- ROTABLATOR™ Consoles ..... 75
- ROTAGO™ Cart ..... 76
- Accessories ..... 76

### Chronic Total Occlusions

- CROSSBOSS™ & STINGRAY™ LP ..... 78

### TRAPPER™ Exchange Device

- TRAPPER™ Exchange Device ..... 80

### ANGIOJET™ Thrombectomy System

- SPIROFLEX™ ..... 82
- SPIROFLEX™ VG ..... 82

## 3 PCI GUIDANCE

### POLARIS™ PCI Guidance System\*

- POLARIS™ Cart System ..... 85
- POLARIS™ Installed System ..... 85

### POLARIS™ PCI Guidance Accessories

- POLARIS™ Accessories ..... 87

### Imaging

- OPTICROSS™ Coronary Imaging Catheter - IVUS ..... 89
- OPTICROSS™ High Definition Imaging Catheter - HD IVUS ..... 89

### Physiology

- COMET™ ..... 91

## 4 STRUCTURAL HEART

### Transcatheter Aortic Valve Implantation (TAVI)

- ACURATE neo™ Aortic Bioprosthesis ..... 94
- ACURATE neo TF™ Transfemoral Delivery System ..... 94
- ACURATE neo™ TA Transapical Delivery System ..... 94
- LOTUS Edge™ Aortic Valve System ..... 95
- LOTUS™ Introducer Set ..... 95
- SENTINEL™ Cerebral Protection System ..... 95
- 14F iSLEEVE™ Expandable Introducer Set ..... 96
- SAFARI2™ Pre-Shaped Guidewire ..... 96

### Left Atrial Appendage Closure (LAAC)

- WATCHMAN FLX™ LAA Closure Device and Delivery System ..... 98
- WATCHMAN™ TruSeal™ Access System ..... 98
- WATCHMAN™ LAA Closure Device and Delivery System ..... 98
- WATCHMAN™ Transseptal Kits ..... 98

# ORDERING INFORMATION

The complete model number must be indicated for ordering purposes.

Example:

| Model Number                                        |                                |                                                            |
|-----------------------------------------------------|--------------------------------|------------------------------------------------------------|
| H749 1981505 0                                      |                                |                                                            |
| Code for<br>Interventional<br>Cardiology<br>Devices | Product<br>Reference<br>Number | Package Indication<br>0 = Single Unit<br>1 = Multiple Pack |

Bold characters are used to facilitate product identification

Indications, contraindications, warnings and instructions for use can be found in the instructions for use supplied with each device.

CAUTION: The law restricts these devices to sale by or on the order of a physician.

For more information, please contact your Boston Scientific sales representative.

**Illustrations** in this catalogue are not drawn to scale.

**Products** shown are subject to modification without prior notice.

1

STENTS  
AND DRUG  
TECHNOLOGIES

2

COMPLEX  
PCI

3

PCI  
GUIDANCE

4

STRUCTURAL  
HEART

# CONVERSION TABLES

| French (F) | mm    |
|------------|-------|
| 1.0        | 0.333 |
| 2.4        | 0.80  |
| 2.6        | 0.87  |
| 2.7        | 0.90  |
| 2.8        | 0.93  |
| 2.9        | 0.97  |
| 3.0        | 1.00  |
| 3.3        | 1.10  |
| 3.4        | 1.13  |
| 3.7        | 1.23  |
| 4.0        | 1.33  |
| 4.5        | 1.50  |
| 5.0        | 1.66  |
| 6.0        | 2.00  |
| 7.0        | 2.33  |
| 8.0        | 2.66  |
| 9.0        | 3.00  |
| 10.0       | 3.33  |
| 11.5       | 3.83  |
| 15.5       | 5.17  |
| 18.0       | 6.00  |
| 22.0       | 7.33  |
| 24.0       | 8.00  |
| 26.0       | 8.66  |
| 28.0       | 9.33  |
| 30.0       | 10.00 |
| 32.0       | 10.66 |
| 34.0       | 11.33 |

| inch (") | mm   |
|----------|------|
| .008     | 0.20 |
| .010     | 0.25 |
| .012     | 0.30 |
| .014     | 0.36 |
| .018     | 0.46 |
| .020     | 0.51 |
| .021     | 0.53 |
| .025     | 0.64 |
| .026     | 0.66 |
| .027     | 0.69 |
| .028     | 0.71 |
| .029     | 0.74 |
| .030     | 0.76 |
| .031     | 0.79 |
| .032     | 0.81 |
| .033     | 0.84 |
| .034     | 0.86 |
| .035     | 0.89 |
| .036     | 0.91 |
| .037     | 0.94 |
| .038     | 0.97 |
| .039     | 0.99 |
| .040     | 1.02 |
| .041     | 1.04 |
| .042     | 1.07 |

| inch (") | mm   |
|----------|------|
| .043     | 1.09 |
| .044     | 1.12 |
| .045     | 1.14 |
| .046     | 1.17 |
| .047     | 1.19 |
| .048     | 1.22 |
| .050     | 1.27 |
| .052     | 1.32 |
| .054     | 1.37 |
| .056     | 1.42 |
| .059     | 1.50 |
| .061     | 1.55 |
| .063     | 1.60 |
| .064     | 1.63 |
| .068     | 1.73 |
| .071     | 1.80 |
| .072     | 1.83 |
| .074     | 1.88 |
| .076     | 1.93 |
| .077     | 1.96 |
| .079     | 2.01 |
| .082     | 2.08 |
| .086     | 2.18 |
| .092     | 2.34 |
| .107     | 2.72 |



| French (F) | 3    | 4    | 5    | 6    | 7    | 8    | 9    | 10   | 11   | 12   | 13   | 14   | 15   | 16   | 17   |
|------------|------|------|------|------|------|------|------|------|------|------|------|------|------|------|------|
| mm ø       | 1.0  | 1.33 | 1.66 | 2.0  | 2.33 | 2.66 | 3.0  | 3.33 | 3.66 | 4.0  | 4.33 | 4.66 | 5.0  | 5.33 | 5.66 |
| inch (")   | .039 | .053 | .066 | .079 | .092 | .105 | .118 | .144 | .144 | .156 | .170 | .184 | .197 | .210 | .223 |



| French (F) | 18   | 19   | 20   | 22   | 24   | 26   | 28   | 30   | 32    | 34    |
|------------|------|------|------|------|------|------|------|------|-------|-------|
| mm ø       | 6.0  | 6.33 | 6.66 | 7.33 | 8.0  | 8.66 | 9.33 | 10.0 | 10.66 | 11.33 |
| inch (")   | .236 | .249 | .263 | .288 | .315 | .341 | .367 | .393 | .419  | .445  |



1

STENTS  
AND DRUG  
TECHNOLOGIES

2

COMPLEX  
PCI

3

PCI  
GUIDANCE

4

STRUCTURAL  
HEART



NEW CARDIO

# 1 STENTS & DRUG TECHNOLOGIES

Abluminal Bioabsorbable Polymer DES

Drug-Eluting Stents

Bare Metal Stent

Drug Coated Balloon

Abluminal  
Bioabsorbable  
Polymer DES

Drug-Eluting  
Stents

Bare Metal  
Stent

Drug Coated  
Balloon



# 1 STENTS & DRUG TECHNOLOGIES

Abluminal Bioabsorbable Polymer DES

Drug-Eluting Stents


Bare Metal Stent

Drug Coated Balloon

Drug-Eluting  
Stents

Bare Metal  
Stent

Drug Coated  
Balloon



# 1. Abluminal Bioabsorbable Polymer DES

**SYNERGY™ \***

Abluminal Bioabsorbable Polymer Drug Eluting Stent

Abluminal  
Bioabsorbable  
Polymer DES

Drug-Eluting  
Stents

Bare Metal  
Stent

Drug Coated  
Balloon

| CATALOG NUMBER /<br>ORDER NUMBER (GTIN) | LENGTH /<br>DIAMETER | PROX.SHAFT<br>DIAMETER | DIST.SHAFT<br>DIAMETER | NOMINAL<br>PRESSURE | RATED BURST<br>PRESSURE |
|-----------------------------------------|----------------------|------------------------|------------------------|---------------------|-------------------------|
| H749 392620822 0 / 08714729840923       | 8 / 2.25 mm          | 2.1F (0.70 mm)         | 2.6F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392620825 0 / 08714729840930       | 8 / 2.50 mm          | 2.1F (0.70 mm)         | 2.6F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392620827 0 / 08714729840947       | 8 / 2.75 mm          | 2.1F (0.70 mm)         | 2.6F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392620830 0 / 08714729840954       | 8 / 3.00 mm          | 2.1F (0.70 mm)         | 2.6F (0.90 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 392620835 0 / 08714729840961       | 8 / 3.50 mm          | 2.1F (0.70 mm)         | 2.6F (0.90 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 392620840 0 / 08714729840978       | 8 / 4.00 mm          | 2.1F (0.70 mm)         | 2.7F (0.95 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 392621222 0 / 08714729840985       | 12 / 2.25 mm         | 2.1F (0.70 mm)         | 2.6F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392621225 0 / 08714729840992       | 12 / 2.50 mm         | 2.1F (0.70 mm)         | 2.6F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392621227 0 / 08714729841005       | 12 / 2.75 mm         | 2.1F (0.70 mm)         | 2.6F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392621230 0 / 08714729841012       | 12 / 3.00 mm         | 2.1F (0.70 mm)         | 2.6F (0.90 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 392621235 0 / 08714729841029       | 12 / 3.50 mm         | 2.1F (0.70 mm)         | 2.6F (0.90 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 392621240 0 / 08714729841036       | 12 / 4.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.95 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 392621250 0 / 08714729974321       | 12 / 4.50 mm         | 2.1F (0.70 mm)         | 2.7F (0.95 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 392621500 0 / 08714729974406       | 12 / 5.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.95 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 392621622 0 / 08714729841043       | 16 / 2.25 mm         | 2.1F (0.70 mm)         | 2.6F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1827 kPa)       |
| H749 392621625 0 / 08714729841050       | 16 / 2.50 mm         | 2.1F (0.70 mm)         | 2.6F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1827 kPa)       |
| H749 392621627 0 / 08714729841067       | 16 / 2.75 mm         | 2.1F (0.70 mm)         | 2.6F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1827 kPa)       |
| H749 392621630 0 / 08714729841074       | 16 / 3.00 mm         | 2.1F (0.70 mm)         | 2.6F (0.90 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1620 kPa)       |
| H749 392621635 0 / 08714729841081       | 16 / 3.50 mm         | 2.1F (0.70 mm)         | 2.6F (0.90 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1620 kPa)       |
| H749 392621640 0 / 08714729841098       | 16 / 4.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.95 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1620 kPa)       |
| H749 392621645 0 / 08714729974338       | 16 / 4.50 mm         | 2.1F (0.70 mm)         | 2.6F (0.90 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1620 kPa)       |
| H749 392621650 0 / 08714729974413       | 16 / 5.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.95 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1620 kPa)       |
| H749 392622022 0 / 08714729841104       | 20 / 2.25 mm         | 2.1F (0.70 mm)         | 2.6F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1827 kPa)       |
| H749 392622025 0 / 08714729841111       | 20 / 2.50 mm         | 2.1F (0.70 mm)         | 2.6F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1827 kPa)       |
| H749 392622027 0 / 08714729841128       | 20 / 2.75 mm         | 2.1F (0.70 mm)         | 2.6F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1827 kPa)       |
| H749 392622030 0 / 08714729841135       | 20 / 3.00 mm         | 2.1F (0.70 mm)         | 2.6F (0.90 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1620 kPa)       |
| H749 392622035 0 / 08714729841142       | 20 / 3.50 mm         | 2.1F (0.70 mm)         | 2.6F (0.90 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1620 kPa)       |
| H749 392622040 0 / 08714729841159       | 20 / 4.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.95 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1620 kPa)       |
| H749 392622045 0 / 08714729974345       | 20 / 4.50 mm         | 2.1F (0.70 mm)         | 2.6F (0.90 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1620 kPa)       |
| H749 392622050 0 / 08714729974420       | 20 / 5.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.95 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1620 kPa)       |
| H749 392622422 0 / 08714729841166       | 24 / 2.25 mm         | 2.1F (0.70 mm)         | 2.6F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1827 kPa)       |
| H749 392622425 0 / 08714729841173       | 24 / 2.50 mm         | 2.1F (0.70 mm)         | 2.6F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1827 kPa)       |
| H749 392622427 0 / 08714729841180       | 24 / 2.75 mm         | 2.1F (0.70 mm)         | 2.6F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1827 kPa)       |
| H749 392622430 0 / 08714729841197       | 24 / 3.00 mm         | 2.1F (0.70 mm)         | 2.6F (0.90 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1620 kPa)       |
| H749 392622435 0 / 08714729841203       | 24 / 3.50 mm         | 2.1F (0.70 mm)         | 2.7F (0.95 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1620 kPa)       |
| H749 392622440 0 / 08714729841210       | 24 / 4.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.95 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1620 kPa)       |
| H749 392622445 0 / 08714729974352       | 24 / 4.50 mm         | 2.1F (0.70 mm)         | 2.7F (0.95 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1620 kPa)       |
| H749 392622450 0 / 08714729974437       | 24 / 5.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.95 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1620 kPa)       |
| H749 392622822 0 / 08714729841227       | 28 / 2.25 mm         | 2.1F (0.70 mm)         | 2.6F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1827 kPa)       |
| H749 392622825 0 / 08714729841234       | 28 / 2.50 mm         | 2.1F (0.70 mm)         | 2.6F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1827 kPa)       |
| H749 392622827 0 / 08714729841241       | 28 / 2.75 mm         | 2.1F (0.70 mm)         | 2.6F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1827 kPa)       |
| H749 392622830 0 / 08714729841258       | 28 / 3.00 mm         | 2.1F (0.70 mm)         | 2.6F (0.90 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1620 kPa)       |
| H749 392622835 0 / 08714729841265       | 28 / 3.50 mm         | 2.1F (0.70 mm)         | 2.7F (0.95 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1620 kPa)       |
| H749 392622840 0 / 08714729841272       | 28 / 4.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.95 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1620 kPa)       |
| H749 392622845 0 / 08714729974369       | 28 / 4.50 mm         | 2.1F (0.70 mm)         | 2.7F (0.95 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1620 kPa)       |
| H749 392622850 0 / 08714729974444       | 28 / 5.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.95 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1620 kPa)       |
| H749 392623222 0 / 08714729841289       | 32 / 2.25 mm         | 2.1F (0.70 mm)         | 2.6F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1827 kPa)       |
| H749 392623225 0 / 08714729841296       | 32 / 2.50 mm         | 2.1F (0.70 mm)         | 2.6F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1827 kPa)       |
| H749 392623227 0 / 08714729841302       | 32 / 2.75 mm         | 2.1F (0.70 mm)         | 2.6F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1827 kPa)       |
| H749 392623230 0 / 08714729841319       | 32 / 3.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.95 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1620 kPa)       |
| H749 392623235 0 / 08714729841326       | 32 / 3.50 mm         | 2.1F (0.70 mm)         | 2.7F (0.95 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1620 kPa)       |
| H749 392623240 0 / 08714729841333       | 32 / 4.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.95 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1620 kPa)       |
| H749 392623245 0 / 08714729974376       | 32 / 4.50 mm         | 2.1F (0.70 mm)         | 2.7F (0.95 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1620 kPa)       |
| H749 392623250 0 / 08714729974451       | 32 / 5.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.95 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1620 kPa)       |

Continued next page ➤



Stents & Drug Technologies - 7  
NEW CARDIO

**Boston  
Scientific**

\*Please check availability of this product with your local sales representative.



# 1. Abluminal Bioabsorbable Polymer DES

## SYNERGY™ \* Continued Abluminal Bioabsorbable Polymer Drug Eluting Stent

| CATALOG NUMBER /<br>ORDER NUMBER (GTIN) | LENGTH /<br>DIAMETER | PROX.SHAFT<br>DIAMETER | DIST.SHAFT<br>DIAMETER | NOMINAL<br>PRESSURE | RATED BURST<br>PRESSURE |
|-----------------------------------------|----------------------|------------------------|------------------------|---------------------|-------------------------|
| H749 392623822 0 / 08714729841340       | 38 / 2.25 mm         | 2.1F (0.70 mm)         | 2.6F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1827 kPa)       |
| H749 392623825 0 / 08714729841357       | 38 / 2.50 mm         | 2.1F (0.70 mm)         | 2.6F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1827 kPa)       |
| H749 392623827 0 / 08714729841364       | 38 / 2.75 mm         | 2.1F (0.70 mm)         | 2.6F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1827 kPa)       |
| H749 392623830 0 / 08714729841371       | 38 / 3.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.95 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1620 kPa)       |
| H749 392623835 0 / 08714729841388       | 38 / 3.50 mm         | 2.1F (0.70 mm)         | 2.7F (0.95 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1620 kPa)       |
| H749 392623840 0 / 08714729841395       | 38 / 4.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.95 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1620 kPa)       |
| H749 392624825 0 / 08714729935612       | 48 / 2.50 mm         | 2.1F (0.70 mm)         | 2.6F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1827 kPa)       |
| H749 392624827 0 / 08714729935629       | 48 / 2.75 mm         | 2.1F (0.70 mm)         | 2.6F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1827 kPa)       |
| H749 392624830 0 / 08714729935636       | 48 / 3.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.95 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1620 kPa)       |
| H749 392624835 0 / 08714729935643       | 48 / 3.50 mm         | 2.1F (0.70 mm)         | 2.7F (0.95 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1620 kPa)       |
| H749 392624840 0 / 08714729935650       | 48 / 4.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.95 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1620 kPa)       |

## SYNERGY MEGATRON™ \* Everolimus-Eluting Platinum Chromium Coronary Stent System

| CATALOG NUMBER /<br>ORDER NUMBER (GTIN) | LENGTH | DIAMETER | OVEREXPANSION<br>CAPABILITIES |
|-----------------------------------------|--------|----------|-------------------------------|
| H749 394270835 0 / 08714729985259       | 8 mm   | 3.50 mm  | 6.0                           |
| H749 394270840 0 / 08714729985266       | 8 mm   | 4.00 mm  | 6.0                           |
| H749 394270845 0 / 08714729985273       | 8 mm   | 4.50 mm  | 6.0                           |
| H749 394270850 0 / 08714729985280       | 8 mm   | 5.00 mm  | 6.0                           |
| H749 394271235 0 / 08714729985297       | 12 mm  | 3.50 mm  | 6.0                           |
| H749 394271240 0 / 08714729985303       | 12 mm  | 4.00 mm  | 6.0                           |
| H749 394271245 0 / 08714729985310       | 12 mm  | 4.50 mm  | 6.0                           |
| H749 394271250 0 / 08714729985327       | 12 mm  | 5.00 mm  | 6.0                           |
| H749 394271635 0 / 08714729985334       | 16 mm  | 3.50 mm  | 6.0                           |
| H749 394271640 0 / 08714729985341       | 16 mm  | 4.00 mm  | 6.0                           |
| H749 394271645 0 / 08714729985358       | 16 mm  | 4.50 mm  | 6.0                           |
| H749 394271650 0 / 08714729985365       | 16 mm  | 5.00 mm  | 6.0                           |
| H749 394272035 0 / 08714729985372       | 20 mm  | 3.50 mm  | 6.0                           |
| H749 394272040 0 / 08714729985389       | 20 mm  | 4.00 mm  | 6.0                           |
| H749 394272045 0 / 08714729985396       | 20 mm  | 4.50 mm  | 6.0                           |
| H749 394272050 0 / 08714729985402       | 20 mm  | 5.00 mm  | 6.0                           |
| H749 394272435 0 / 08714729985419       | 24 mm  | 3.50 mm  | 6.0                           |
| H749 394272440 0 / 08714729985426       | 24 mm  | 4.00 mm  | 6.0                           |
| H749 394272445 0 / 08714729985433       | 24 mm  | 4.50 mm  | 6.0                           |
| H749 394272450 0 / 08714729985440       | 24 mm  | 5.00 mm  | 6.0                           |
| H749 394272835 0 / 08714729985457       | 28 mm  | 3.50 mm  | 6.0                           |
| H749 394272840 0 / 08714729985464       | 28 mm  | 4.00 mm  | 6.0                           |
| H749 394272845 0 / 08714729985471       | 28 mm  | 4.50 mm  | 6.0                           |
| H749 394272850 0 / 08714729985488       | 28 mm  | 5.00 mm  | 6.0                           |
| H749 394273235 0 / 08714729985495       | 32 mm  | 3.50 mm  | 6.0                           |
| H749 394273240 0 / 08714729985501       | 32 mm  | 4.00 mm  | 6.0                           |
| H749 394273245 0 / 08714729985518       | 32 mm  | 4.50 mm  | 6.0                           |
| H749 394273250 0 / 08714729985525       | 32 mm  | 5.00 mm  | 6.0                           |

Abluminal  
Bioabsorbable  
Polymer DES

Drug-Eluting  
Stents

Bare Metal  
Stent

Drug Coated  
Balloon



# 1

## STENTS & DRUG TECHNOLOGIES

Abluminal Bioabsorbable Polymer DES

Drug-Eluting Stents

Bare Metal Stent

Drug Coated Balloon

Drug-Eluting  
Stents

Bare Metal  
Stent

Drug Coated  
Balloon



# 1. Drug-Eluting Stents

## Promus ELITE™ \*

Everolimus-Eluting Platinum Chromium Coronary Stent System

Abluminal  
Bioabsorbable  
Polymer DES

Drug-Eluting  
Stents

Bare Metal  
Stent

Drug Coated  
Balloon

| CATALOG NUMBER /<br>ORDER NUMBER (GTIN) | LENGTH /<br>DIAMETER | PROX. SHAFT<br>DIAMETER | DIST. SHAFT<br>DIAMETER | NOMINAL<br>PRESSURE | RATED BURST<br>PRESSURE |
|-----------------------------------------|----------------------|-------------------------|-------------------------|---------------------|-------------------------|
| H749 394130822 0 / 08714729971979       | 8 / 2.25 mm          | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 394130825 0 / 08714729971986       | 8 / 2.50 mm          | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 394130827 0 / 08714729971993       | 8 / 2.75 mm          | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 394130830 0 / 08714729972006       | 8 / 3.00 mm          | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 394130835 0 / 08714729972013       | 8 / 3.50 mm          | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 394130840 0 / 08714729972020       | 8 / 4.00 mm          | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 394131222 0 / 08714729972037       | 12 / 2.25 mm         | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 394131225 0 / 08714729972044       | 12 / 2.50 mm         | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 394131227 0 / 08714729972051       | 12 / 2.75 mm         | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 394131230 0 / 08714729972068       | 12 / 3.00 mm         | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 394131235 0 / 08714729972075       | 12 / 3.50 mm         | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 394131240 0 / 08714729972082       | 12 / 4.00 mm         | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 394131622 0 / 08714729972099       | 16 / 2.25 mm         | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 394131625 0 / 08714729972105       | 16 / 2.50 mm         | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 394131627 0 / 08714729972112       | 16 / 2.75 mm         | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 394131630 0 / 08714729972129       | 16 / 3.00 mm         | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 394131635 0 / 08714729972136       | 16 / 3.50 mm         | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 394131640 0 / 08714729972143       | 16 / 4.00 mm         | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 394132022 0 / 08714729972150       | 20 / 2.25 mm         | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 394132025 0 / 08714729972167       | 20 / 2.50 mm         | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 394132027 0 / 08714729972174       | 20 / 2.75 mm         | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 394132030 0 / 08714729972181       | 20 / 3.00 mm         | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 394132035 0 / 08714729972198       | 20 / 3.50 mm         | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 394132040 0 / 08714729972204       | 20 / 4.00 mm         | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 394132422 0 / 08714729972211       | 24 / 2.25 mm         | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 394132425 0 / 08714729972228       | 24 / 2.50 mm         | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 394132427 0 / 08714729972235       | 24 / 2.75 mm         | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 394132430 0 / 08714729972242       | 24 / 3.00 mm         | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 394132435 0 / 08714729972259       | 24 / 3.50 mm         | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 394132440 0 / 08714729972266       | 24 / 4.00 mm         | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 394132822 0 / 08714729972273       | 28 / 2.25 mm         | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 394132825 0 / 08714729972280       | 28 / 2.50 mm         | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 394132827 0 / 08714729972297       | 28 / 2.75 mm         | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 394132830 0 / 08714729972303       | 28 / 3.00 mm         | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 394132835 0 / 08714729972310       | 28 / 3.50 mm         | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 394132840 0 / 08714729972327       | 28 / 4.00 mm         | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 394133222 0 / 08714729972334       | 32 / 2.25 mm         | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 394133225 0 / 08714729972341       | 32 / 2.50 mm         | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 394133227 0 / 08714729972358       | 32 / 2.75 mm         | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 394133230 0 / 08714729972365       | 32 / 3.00 mm         | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 394133235 0 / 08714729972372       | 32 / 3.50 mm         | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 394133240 0 / 08714729972389       | 32 / 4.00 mm         | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 394133825 0 / 08714729972396       | 38 / 2.50 mm         | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 394133827 0 / 08714729972402       | 38 / 2.75 mm         | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 394133830 0 / 08714729972419       | 38 / 3.00 mm         | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 394133835 0 / 08714729972426       | 38 / 3.50 mm         | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 394133840 0 / 08714729972433       | 38 / 4.00 mm         | 2.1 F (0.70 mm)         | 2.7F (0.90 mm)          | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |





# 1. Drug-Eluting Stents

## Promus PREMIER™ \*

Everolimus-Eluting Platinum Chromium Coronary Stent System with Customized Stent Architecture

Abluminal  
Bioabsorbable  
Polymer DES

Drug-Eluting  
Stents

Bare Metal  
Stent

Drug Coated  
Balloon

| CATALOG NUMBER /<br>ORDER NUMBER (GTIN) | LENGTH /<br>DIAMETER | PROX.SHAFT<br>DIAMETER | DIST.SHAFT<br>DIAMETER | NOMINAL<br>PRESSURE | RATED BURST<br>PRESSURE |
|-----------------------------------------|----------------------|------------------------|------------------------|---------------------|-------------------------|
| H749 392510822 0 / 08714729828266       | 8 / 2.25 mm          | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392510825 0 / 08714729828273       | 8 / 2.50 mm          | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392510827 0 / 08714729828280       | 8 / 2.75 mm          | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392510830 0 / 08714729828297       | 8 / 3.00 mm          | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 392510835 0 / 08714729828303       | 8 / 3.50 mm          | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 392510840 0 / 08714729828310       | 8 / 4.00 mm          | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 392511222 0 / 08714729828327       | 12 / 2.25 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392511225 0 / 08714729828334       | 12 / 2.50 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392511227 0 / 08714729828341       | 12 / 2.75 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392511230 0 / 08714729828358       | 12 / 3.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 392511235 0 / 08714729828365       | 12 / 3.50 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 392511240 0 / 08714729828372       | 12 / 4.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 392511622 0 / 08714729828389       | 16 / 2.25 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392511625 0 / 08714729828396       | 16 / 2.50 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392511627 0 / 08714729828402       | 16 / 2.75 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392511630 0 / 08714729828419       | 16 / 3.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 392511635 0 / 08714729828426       | 16 / 3.50 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 392511640 0 / 08714729828433       | 16 / 4.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 392512022 0 / 08714729828440       | 20 / 2.25 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392512025 0 / 08714729828457       | 20 / 2.50 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392512027 0 / 08714729828464       | 20 / 2.75 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392512030 0 / 08714729828471       | 20 / 3.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 392512035 0 / 08714729828488       | 20 / 3.50 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 392512040 0 / 08714729828495       | 20 / 4.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 392512422 0 / 08714729828501       | 24 / 2.25 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392512425 0 / 08714729828518       | 24 / 2.50 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392512427 0 / 08714729828525       | 24 / 2.75 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392512430 0 / 08714729828532       | 24 / 3.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 392512435 0 / 08714729828549       | 24 / 3.50 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 392512440 0 / 08714729828556       | 24 / 4.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 392512822 0 / 08714729828563       | 28 / 2.25 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392512825 0 / 08714729828570       | 28 / 2.50 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392512827 0 / 08714729828587       | 28 / 2.75 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392512830 0 / 08714729828594       | 28 / 3.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 392512835 0 / 08714729828600       | 28 / 3.50 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 392512840 0 / 08714729828617       | 28 / 4.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 392513222 0 / 08714729828624       | 32 / 2.25 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392513225 0 / 08714729828631       | 32 / 2.50 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392513227 0 / 08714729828648       | 32 / 2.75 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392513230 0 / 08714729828655       | 32 / 3.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 392513235 0 / 08714729828662       | 32 / 3.50 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 392513240 0 / 08714729828679       | 32 / 4.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 392513825 0 / 08714729828686       | 38 / 2.50 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392513827 0 / 08714729828693       | 38 / 2.75 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392513830 0 / 08714729828709       | 38 / 3.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 392513835 0 / 08714729828716       | 38 / 3.50 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |
| H749 392513840 0 / 08714729828723       | 38 / 4.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 16 ATM (1621 kPa)       |



# 1 STENTS & DRUG TECHNOLOGIES

Abluminal Bioabsorbable Polymer DES

Drug-Eluting Stents

Bare Metal Stent

Drug Coated Balloon

Abluminal  
Bioabsorbable  
Polymer DES

Drug-Eluting  
Stents

Bare Metal  
Stent

Drug Coated  
Balloon



# 1. Bare Metal Stent

**REBEL™**

Platinum Chromium Coronary Stent System with Customized Stent Architecture

Abluminal  
Bioabsorbable  
Polymer DES

Drug-Eluting  
Stents

Bare Metal  
Stent

Drug Coated  
Balloon

| CATALOG NUMBER /<br>ORDER NUMBER (GTIN) | LENGTH /<br>DIAMETER | PROX.SHAFT<br>DIAMETER | DIST.SHAFT<br>DIAMETER | NOMINAL<br>PRESSURE | RATED BURST<br>PRESSURE |
|-----------------------------------------|----------------------|------------------------|------------------------|---------------------|-------------------------|
| H749 392730822 0 / 08714729843146       | 8 / 2.25 mm          | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392730825 0 / 08714729843153       | 8 / 2.50 mm          | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392730827 0 / 08714729843160       | 8 / 2.75 mm          | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392730830 0 / 08714729843177       | 8 / 3.00 mm          | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392730835 0 / 08714729843184       | 8 / 3.50 mm          | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392730840 0 / 08714729843191       | 8 / 4.00 mm          | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392731222 0 / 08714729843207       | 12 / 2.25 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392731225 0 / 08714729843214       | 12 / 2.50 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392731227 0 / 08714729843221       | 12 / 2.75 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392731230 0 / 08714729843238       | 12 / 3.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392731235 0 / 08714729843245       | 12 / 3.50 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392731240 0 / 08714729843252       | 12 / 4.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392731245 0 / 08714729843269       | 12 / 4.50 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392731622 0 / 08714729843276       | 16 / 2.25 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392731625 0 / 08714729843283       | 16 / 2.50 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392731627 0 / 08714729843290       | 16 / 2.75 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392731630 0 / 08714729843306       | 16 / 3.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392731635 0 / 08714729843313       | 16 / 3.50 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392731640 0 / 08714729843320       | 16 / 4.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392731645 0 / 08714729843337       | 16 / 4.50 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392732022 0 / 08714729843344       | 20 / 2.25 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392732025 0 / 08714729843351       | 20 / 2.50 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392732027 0 / 08714729843368       | 20 / 2.75 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392732030 0 / 08714729843375       | 20 / 3.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392732035 0 / 08714729843382       | 20 / 3.50 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392732040 0 / 08714729843399       | 20 / 4.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392732045 0 / 08714729843405       | 20 / 4.50 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392732422 0 / 08714729843412       | 24 / 2.25 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392732425 0 / 08714729843429       | 24 / 2.50 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392732427 0 / 08714729843436       | 24 / 2.75 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392732430 0 / 08714729843443       | 24 / 3.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392732435 0 / 08714729843450       | 24 / 3.50 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392732440 0 / 08714729843467       | 24 / 4.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392732445 0 / 08714729843474       | 24 / 4.50 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392732822 0 / 08714729843481       | 28 / 2.25 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392732825 0 / 08714729843498       | 28 / 2.50 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392732827 0 / 08714729843504       | 28 / 2.75 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392732830 0 / 08714729843511       | 28 / 3.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392732835 0 / 08714729843528       | 28 / 3.50 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392732840 0 / 08714729843535       | 28 / 4.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392732845 0 / 08714729843542       | 28 / 4.50 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392733222 0 / 08714729843559       | 32 / 2.25 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392733225 0 / 08714729843566       | 32 / 2.50 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392733227 0 / 08714729843573       | 32 / 2.75 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392733230 0 / 08714729843580       | 32 / 3.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392733235 0 / 08714729843597       | 32 / 3.50 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392733240 0 / 08714729843603       | 32 / 4.00 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |
| H749 392733245 0 / 08714729843610       | 32 / 4.50 mm         | 2.1F (0.70 mm)         | 2.7F (0.90 mm)         | 11 ATM (1117 kPa)   | 18 ATM (1824 kPa)       |




# 1 STENTS & DRUG TECHNOLOGIES

Abluminal Bioabsorbable Polymer DES

Drug-Eluting Stents

Bare Metal Stent

Drug Coated Balloon



Abluminal  
Bioabsorbable  
Polymer DES

Drug-Eluting  
Stents

Bare Metal  
Stent

Drug Coated  
Balloon



1. Drug Coated Balloon

AGENT™

Paclitaxel-Coated PTCA Balloon Catheter

| CATALOG NUMBER /<br>ORDER NUMBER (GTIN) | LENGTH /<br>DIAMETER | PROX.SHAFT<br>DIAMETER | DIST.SHAFT<br>DIAMETER | NOMINAL<br>PRESSURE | RATED BURST<br>PRESSURE |
|-----------------------------------------|----------------------|------------------------|------------------------|---------------------|-------------------------|
| H749 3922220121 0 / 08714729831402      | 12 mm                | 2.00 mm                | 0.017" / 0.43 mm       | 6 ATM (608 kPa)     | 14 ATM (1419 kPa)       |
| H749 3922220151 0 / 08714729831419      | 15 mm                | 2.00 mm                | 0.017" / 0.43 mm       | 6 ATM (608 kPa)     | 14 ATM (1419 kPa)       |
| H749 3922220201 0 / 08714729831426      | 20 mm                | 2.00 mm                | 0.017" / 0.43 mm       | 6 ATM (608 kPa)     | 14 ATM (1419 kPa)       |
| H749 3922220301 0 / 08714729831433      | 30 mm                | 2.00 mm                | 0.017" / 0.43 mm       | 6 ATM (608 kPa)     | 14 ATM (1419 kPa)       |
| H749 3922222121 0 / 08714729831457      | 12 mm                | 2.25 mm                | 0.017" / 0.43 mm       | 6 ATM (608 kPa)     | 14 ATM (1419 kPa)       |
| H749 3922222151 0 / 08714729831464      | 15 mm                | 2.25 mm                | 0.017" / 0.43 mm       | 6 ATM (608 kPa)     | 14 ATM (1419 kPa)       |
| H749 3922222201 0 / 08714729831471      | 20 mm                | 2.25 mm                | 0.017" / 0.43 mm       | 6 ATM (608 kPa)     | 14 ATM (1419 kPa)       |
| H749 3922222301 0 / 08714729831488      | 30 mm                | 2.25 mm                | 0.017" / 0.43 mm       | 6 ATM (608 kPa)     | 14 ATM (1419 kPa)       |
| H749 3922225121 0 / 08714729831501      | 12 mm                | 2.50 mm                | 0.017" / 0.43 mm       | 6 ATM (608 kPa)     | 14 ATM (1419 kPa)       |
| H749 3922225151 0 / 08714729831518      | 15 mm                | 2.50 mm                | 0.017" / 0.43 mm       | 6 ATM (608 kPa)     | 14 ATM (1419 kPa)       |
| H749 3922225201 0 / 08714729831525      | 20 mm                | 2.50 mm                | 0.017" / 0.43 mm       | 6 ATM (608 kPa)     | 14 ATM (1419 kPa)       |
| H749 3922225301 0 / 08714729831532      | 30 mm                | 2.50 mm                | 0.017" / 0.43 mm       | 6 ATM (608 kPa)     | 14 ATM (1419 kPa)       |
| H749 3922227121 0 / 08714729831556      | 12 mm                | 2.75 mm                | 0.017" / 0.43 mm       | 6 ATM (608 kPa)     | 14 ATM (1419 kPa)       |
| H749 3922227151 0 / 08714729831563      | 15 mm                | 2.75 mm                | 0.017" / 0.43 mm       | 6 ATM (608 kPa)     | 14 ATM (1419 kPa)       |
| H749 3922227201 0 / 08714729831570      | 20 mm                | 2.75 mm                | 0.017" / 0.43 mm       | 6 ATM (608 kPa)     | 14 ATM (1419 kPa)       |
| H749 3922227301 0 / 08714729831587      | 30 mm                | 2.75 mm                | 0.017" / 0.43 mm       | 6 ATM (608 kPa)     | 14 ATM (1419 kPa)       |
| H749 3922230121 0 / 08714729831600      | 12 mm                | 3.00 mm                | 0.017" / 0.43 mm       | 6 ATM (608 kPa)     | 14 ATM (1419 kPa)       |
| H749 3922230151 0 / 08714729831617      | 15 mm                | 3.00 mm                | 0.017" / 0.43 mm       | 6 ATM (608 kPa)     | 14 ATM (1419 kPa)       |
| H749 3922230201 0 / 08714729831624      | 20 mm                | 3.00 mm                | 0.017" / 0.43 mm       | 6 ATM (608 kPa)     | 14 ATM (1419 kPa)       |
| H749 3922230301 0 / 08714729831631      | 30 mm                | 3.00 mm                | 0.017" / 0.43 mm       | 6 ATM (608 kPa)     | 14 ATM (1419 kPa)       |
| H749 3922235121 0 / 08714729831709      | 12 mm                | 3.50 mm                | 0.017" / 0.43 mm       | 6 ATM (608 kPa)     | 12 ATM (1216 kPa)       |
| H749 3922235151 0 / 08714729831716      | 15 mm                | 3.50 mm                | 0.017" / 0.43 mm       | 6 ATM (608 kPa)     | 12 ATM (1216 kPa)       |
| H749 3922235201 0 / 08714729831723      | 20 mm                | 3.50 mm                | 0.017" / 0.43 mm       | 6 ATM (608 kPa)     | 12 ATM (1216 kPa)       |
| H749 3922235301 0 / 08714729831730      | 30 mm                | 3.50 mm                | 0.017" / 0.43 mm       | 6 ATM (608 kPa)     | 12 ATM (1216 kPa)       |
| H749 3922240121 0 / 08714729831808      | 12 mm                | 4.00 mm                | 0.017" / 0.43 mm       | 6 ATM (608 kPa)     | 12 ATM (1216 kPa)       |
| H749 3922240151 0 / 08714729831815      | 15 mm                | 4.00 mm                | 0.017" / 0.43 mm       | 6 ATM (608 kPa)     | 12 ATM (1216 kPa)       |
| H749 3922240201 0 / 08714729831822      | 20 mm                | 4.00 mm                | 0.017" / 0.43 mm       | 6 ATM (608 kPa)     | 12 ATM (1216 kPa)       |
| H749 3922240301 0 / 08714729831839      | 30 mm                | 4.00 mm                | 0.017" / 0.43 mm       | 6 ATM (608 kPa)     | 12 ATM (1216 kPa)       |

Abluminal  
Bioabsorbable  
Polymer DES

Drug-Eluting  
Stents

Bare Metal  
Stent

Drug Coated  
Balloon



NEW CARDIO

2

## COMPLEX PCI

Balloon Catheters

Microdilatation Catheters

Guide Catheters

Guide Catheter Extension

Angiographic Catheters

Guidewires

Embolic Protection

Accessories

Cutting Balloon™

Atherectomy Devices

Chronic Total Occlusions

TRAPPER™ Exchange Device

ANGIOJET™ Thrombectomy System



Balloon  
Catheters

Micro-  
dilatation  
Catheters

Guide  
Catheters

Guide  
Catheter  
extension

Angiographic  
Catheters

Guidewires

Embolic  
Protection

Accessories

Cutting  
Balloon™

Atherectomy  
Devices

Chronic Total  
Occlusions

Trapper™  
Exchange  
Device

ANGIOJET™  
Thrombec-  
tomy  
System

NEW CARDIO

2

## COMPLEX PCI

Balloon Catheters

Microdilatation Catheters

Guide Catheters

Guide Catheter Extension

Angiographic Catheters

Guidewires

Embolic Protection

Accessories

Cutting Balloon™

Atherectomy Devices

Chronic Total Occlusions

TRAPPER™ Exchange Device

ANGIOJET™ Thrombectomy System



Balloon  
Catheters

Micro-  
dilatation  
Catheters

Guide  
Catheters

Guide  
Catheter  
extension

Angiographic  
Catheters

Guidewires

Embolic  
Protection

Accessories

Cutting  
Balloon™

Atherectomy  
Devices

Chronic Total  
Occlusions

Trapper™  
Exchange  
Device

ANGIOJET™  
Thrombec-  
tomy  
System

## 2. Pre-dilatation Balloon Catheters

### EMERGE™ Monorail™ PTCA Dilatation Catheter

| CATALOG NUMBER / ORDER NUMBER (GTIN) | LENGTH / DIAMETER | LESION ENTRY PROFILE | NOMINAL PRESSURE | RATED BURST PRESSURE |
|--------------------------------------|-------------------|----------------------|------------------|----------------------|
| H749 391930812 0 / 08714729799764    | 8 / 1.20 mm       | 0.017" / 0.43        | 6 ATM (608 kPa)  | 18 ATM (1824 kPa)    |
| H749 391940812 0 / 08714729799689    | 8 / 1.20 mm Push  | 0.017" / 0.43        | 6 ATM (608 kPa)  | 18 ATM (1824 kPa)    |
| H749 391930815 0 / 08714729799771    | 8 / 1.50 mm       | 0.017" / 0.43        | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391940815 0 / 08714729799726    | 8 / 1.50 mm Push  | 0.017" / 0.43        | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391930820 0 / 08714729799788    | 8 / 2.00 mm       | 0.017" / 0.43        | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391930822 0 / 08714729799795    | 8 / 2.25 mm       | 0.017" / 0.43        | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391930825 0 / 08714729799801    | 8 / 2.50 mm       | 0.017" / 0.43        | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391930827 0 / 08714729799818    | 8 / 2.75 mm       | 0.017" / 0.43        | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391930830 0 / 08714729799825    | 8 / 3.00 mm       | 0.017" / 0.43        | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391930832 0 / 08714729799832    | 8 / 3.25 mm       | 0.017" / 0.43        | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391930835 0 / 08714729799849    | 8 / 3.50 mm       | 0.017" / 0.43        | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 391930837 0 / 08714729799856    | 8 / 3.75 mm       | 0.017" / 0.43        | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 391930840 0 / 08714729799863    | 8 / 4.00 mm       | 0.017" / 0.43        | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 391931212 0 / 08714729799870    | 12 / 1.20 mm      | 0.017" / 0.43        | 6 ATM (608 kPa)  | 18 ATM (1824 kPa)    |
| H749 391941212 0 / 08714729799696    | 12 / 1.20 mm Push | 0.017" / 0.43        | 6 ATM (608 kPa)  | 18 ATM (1824 kPa)    |
| H749 391931215 0 / 08714729799887    | 12 / 1.50 mm      | 0.017" / 0.43        | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391941215 0 / 08714729799733    | 12 / 1.50 mm Push | 0.017" / 0.43        | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391931220 0 / 08714729799894    | 12 / 2.00 mm      | 0.017" / 0.43        | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391931222 0 / 08714729799900    | 12 / 2.25 mm      | 0.017" / 0.43        | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391931225 0 / 08714729799917    | 12 / 2.50 mm      | 0.017" / 0.43        | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391931227 0 / 08714729799924    | 12 / 2.75 mm      | 0.017" / 0.43        | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391931230 0 / 08714729799931    | 12 / 3.00 mm      | 0.017" / 0.43        | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391931232 0 / 08714729799948    | 12 / 3.25 mm      | 0.017" / 0.43        | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391931235 0 / 08714729799955    | 12 / 3.50 mm      | 0.017" / 0.43        | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 391931237 0 / 08714729799962    | 12 / 3.75 mm      | 0.017" / 0.43        | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 391931240 0 / 08714729799979    | 12 / 4.00 mm      | 0.017" / 0.43        | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 391931512 0 / 08714729799986    | 15 / 1.20 mm      | 0.017" / 0.43        | 6 ATM (608 kPa)  | 18 ATM (1824 kPa)    |
| H749 391941512 0 / 08714729799702    | 15 / 1.20 mm Push | 0.017" / 0.43        | 6 ATM (608 kPa)  | 18 ATM (1824 kPa)    |
| H749 391931515 0 / 08714729799993    | 15 / 1.50 mm      | 0.017" / 0.43        | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391941515 0 / 08714729799740    | 15 / 1.50 mm Push | 0.017" / 0.43        | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391931520 0 / 08714729800002    | 15 / 2.00 mm      | 0.017" / 0.43        | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391931522 0 / 08714729800019    | 15 / 2.25 mm      | 0.017" / 0.43        | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391931525 0 / 08714729800026    | 15 / 2.50 mm      | 0.017" / 0.43        | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391931527 0 / 08714729800033    | 15 / 2.75 mm      | 0.017" / 0.43        | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391931530 0 / 08714729800040    | 15 / 3.00 mm      | 0.017" / 0.43        | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391931532 0 / 08714729800057    | 15 / 3.25 mm      | 0.017" / 0.43        | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391931535 0 / 08714729800064    | 15 / 3.50 mm      | 0.017" / 0.43        | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 391931537 0 / 08714729800071    | 15 / 3.75 mm      | 0.017" / 0.43        | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 391931540 0 / 08714729800088    | 15 / 4.00 mm      | 0.017" / 0.43        | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 391932012 0 / 08714729800095    | 20 / 1.20 mm      | 0.017" / 0.43        | 6 ATM (608 kPa)  | 18 ATM (1824 kPa)    |
| H749 391942012 0 / 08714729799719    | 20 / 1.20 mm Push | 0.017" / 0.43        | 6 ATM (608 kPa)  | 18 ATM (1824 kPa)    |
| H749 391932015 0 / 08714729800101    | 20 / 1.50 mm      | 0.017" / 0.43        | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391942015 0 / 08714729799757    | 20 / 1.50 mm Push | 0.017" / 0.43        | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391932020 0 / 08714729800118    | 20 / 2.00 mm      | 0.017" / 0.43        | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391932022 0 / 08714729800125    | 20 / 2.25 mm      | 0.017" / 0.43        | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391932025 0 / 08714729800132    | 20 / 2.50 mm      | 0.017" / 0.43        | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391932027 0 / 08714729800149    | 20 / 2.75 mm      | 0.017" / 0.43        | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391932030 0 / 08714729800156    | 20 / 3.00 mm      | 0.017" / 0.43        | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391932032 0 / 08714729800163    | 20 / 3.25 mm      | 0.017" / 0.43        | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391932035 0 / 08714729800170    | 20 / 3.50 mm      | 0.017" / 0.43        | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 391932037 0 / 08714729800187    | 20 / 3.75 mm      | 0.017" / 0.43        | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 391932040 0 / 08714729800194    | 20 / 4.00 mm      | 0.017" / 0.43        | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |

Continued next page >



Complex PCI - 18  
NEW CARDIO

Boston  
Scientific

Balloon  
Catheters

Micro-  
dilatation  
Catheters

Guide  
Catheters

Guide  
Catheter  
extension

Angiographic  
Catheters

Guidewires

Embolic  
Protection

Accessories

Cutting  
Balloon™

Atherectomy  
Devices

Chronic Total  
Occlusions

Trapper™  
Exchange  
Device

ANGIOJET™  
Thrombec-  
tomy  
System



## 2. Pre-dilatation Balloon Catheters

### EMERGE™ Monorail™ Continued PTCA Dilatation Catheter

| CATALOG NUMBER / ORDER NUMBER (GTIN) | LENGTH / DIAMETER | LESION ENTRY PROFILE | NOMINAL PRESSURE | RATED BURST PRESSURE |
|--------------------------------------|-------------------|----------------------|------------------|----------------------|
| H749 391933020 0 / 08714729800200    | 30 / 2.00 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391933022 0 / 08714729800217    | 30 / 2.25 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391933025 0 / 08714729800224    | 30 / 2.50 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391933027 0 / 08714729800231    | 30 / 2.75 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391933030 0 / 08714729800248    | 30 / 3.00 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391933032 0 / 08714729800255    | 30 / 3.25 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391933035 0 / 08714729800262    | 30 / 3.50 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 391933037 0 / 08714729800279    | 30 / 3.75 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 391933040 0 / 08714729800286    | 30 / 4.00 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |

### EMERGE™ OTW PTCA Dilatation Catheter

| CATALOG NUMBER / ORDER NUMBER (GTIN) | LENGTH / DIAMETER | LESION ENTRY PROFILE | NOMINAL PRESSURE | RATED BURST PRESSURE |
|--------------------------------------|-------------------|----------------------|------------------|----------------------|
| H749 391950812 0 / 08714729800378    | 8 / 1.20 mm       | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 18 ATM (1824 kPa)    |
| H749 391960812 0 / 08714729800293    | 8 / 1.20 mm Push  | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 18 ATM (1824 kPa)    |
| H749 391950815 0 / 08714729800385    | 8 / 1.50 mm       | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391960815 0 / 08714729800330    | 8 / 1.50 mm Push  | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391950820 0 / 08714729800392    | 8 / 2.00 mm       | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391950822 0 / 08714729800408    | 8 / 2.25 mm       | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391950825 0 / 08714729800415    | 8 / 2.50 mm       | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391950827 0 / 08714729800422    | 8 / 2.75 mm       | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391950830 0 / 08714729800439    | 8 / 3.00 mm       | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391950832 0 / 08714729800446    | 8 / 3.25 mm       | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391950835 0 / 08714729800453    | 8 / 3.50 mm       | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 391950837 0 / 08714729800460    | 8 / 3.75 mm       | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 391950840 0 / 08714729800477    | 8 / 4.00 mm       | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 391951212 0 / 08714729800484    | 12 / 1.20 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 18 ATM (1824 kPa)    |
| H749 391961212 0 / 08714729800309    | 12 / 1.20 mm Push | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 18 ATM (1824 kPa)    |
| H749 391951215 0 / 08714729800491    | 12 / 1.50 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391961215 0 / 08714729800347    | 12 / 1.50 mm Push | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391951220 0 / 08714729800507    | 12 / 2.00 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391951222 0 / 08714729800514    | 12 / 2.25 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391951225 0 / 08714729800521    | 12 / 2.50 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391951227 0 / 08714729800538    | 12 / 2.75 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391951230 0 / 08714729800545    | 12 / 3.00 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391951232 0 / 08714729800552    | 12 / 3.25 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391951235 0 / 08714729800569    | 12 / 3.50 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 391951237 0 / 08714729800576    | 12 / 3.75 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 391951240 0 / 08714729800583    | 12 / 4.00 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 391951512 0 / 08714729800590    | 15 / 1.20 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 18 ATM (1824 kPa)    |
| H749 391961512 0 / 08714729800316    | 15 / 1.20 mm Push | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 18 ATM (1824 kPa)    |
| H749 391951515 0 / 08714729800606    | 15 / 1.50 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391961515 0 / 08714729800354    | 15 / 1.50 mm Push | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391951520 0 / 08714729800613    | 15 / 2.00 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391951522 0 / 08714729800620    | 15 / 2.25 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391951525 0 / 08714729800637    | 15 / 2.50 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391951527 0 / 08714729800644    | 15 / 2.75 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |

Continued next page >



## 2. Pre-dilatation Balloon Catheters

### EMERGE™ OTW Continued PTCA Dilatation Catheter

| CATALOG NUMBER / ORDER NUMBER (GTIN) | LENGTH / DIAMETER | LESION ENTRY PROFILE | NOMINAL PRESSURE | RATED BURST PRESSURE |
|--------------------------------------|-------------------|----------------------|------------------|----------------------|
| H749 391951530 0 / 08714729800651    | 15 / 3.00 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391951532 0 / 08714729800668    | 15 / 3.25 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391951535 0 / 08714729800675    | 15 / 3.50 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 391951537 0 / 08714729800682    | 15 / 3.75 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 391951540 0 / 08714729800699    | 15 / 4.00 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 391952012 0 / 08714729800705    | 20 / 1.20 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 18 ATM (1824 kPa)    |
| H749 391962012 0 / 08714729800323    | 20 / 1.20 mm Push | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 18 ATM (1824 kPa)    |
| H749 391952015 0 / 08714729800712    | 20 / 1.50 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391962015 0 / 08714729800361    | 20 / 1.50 mm Push | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391952020 0 / 08714729800729    | 20 / 2.00 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391952022 0 / 08714729800736    | 20 / 2.25 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391952025 0 / 08714729800743    | 20 / 2.50 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391952027 0 / 08714729800750    | 20 / 2.75 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391952030 0 / 08714729800767    | 20 / 3.00 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391952032 0 / 08714729800774    | 20 / 3.25 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391952035 0 / 08714729800781    | 20 / 3.50 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 391952037 0 / 08714729800798    | 20 / 3.75 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 391952040 0 / 08714729800804    | 20 / 4.00 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 391953020 0 / 08714729800811    | 30 / 2.00 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391953022 0 / 08714729800828    | 30 / 2.25 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391953025 0 / 08714729800835    | 30 / 2.50 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391953027 0 / 08714729800842    | 30 / 2.75 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391953030 0 / 08714729800859    | 30 / 3.00 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391953032 0 / 08714729800866    | 30 / 3.25 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 391953035 0 / 08714729800873    | 30 / 3.50 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 391953037 0 / 08714729800880    | 30 / 3.75 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 391953040 0 / 08714729800897    | 30 / 4.00 mm      | 0.017" / 0.43 mm     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |



## 2. Pre-dilatation Balloon Catheters

### MAVERICK<sup>2</sup>™ Monorail™

| CATALOG NUMBER / ORDER NUMBER (GTIN) | LENGTH / DIAMETER | LESION ENTRY PROFILE | NOMINAL PRESSURE | RATED BURST PRESSURE |
|--------------------------------------|-------------------|----------------------|------------------|----------------------|
| H749 389280920 0 / 08714729369875    | 9 / 2.00 mm       | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 389280922 0 / 08714729369882    | 9 / 2.25 mm       | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 389280925 0 / 08714729369899    | 9 / 2.50 mm       | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 389280927 0 / 08714729369905    | 9 / 2.75 mm       | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 389280930 0 / 08714729369912    | 9 / 3.00 mm       | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 389280932 0 / 08714729369929    | 9 / 3.25 mm       | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 389280935 0 / 08714729369936    | 9 / 3.50 mm       | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 389280937 0 / 08714729369943    | 9 / 3.75 mm       | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 389280940 0 / 08714729369950    | 9 / 4.00 mm       | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 389281220 0 / 08714729369967    | 12 / 2.00 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 389281222 0 / 08714729369974    | 12 / 2.25 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 389281225 0 / 08714729369981    | 12 / 2.50 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 389281227 0 / 08714729369998    | 12 / 2.75 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 389281230 0 / 08714729370000    | 12 / 3.00 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 389281232 0 / 08714729370017    | 12 / 3.25 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 389281235 0 / 08714729370024    | 12 / 3.50 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 389281237 0 / 08714729370031    | 12 / 3.75 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 389281240 0 / 08714729370048    | 12 / 4.00 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 389281515 0 / 08714729370055    | 15 / 1.50 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 389281520 0 / 08714729370062    | 15 / 2.00 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 389281522 0 / 08714729370079    | 15 / 2.25 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 389281525 0 / 08714729370086    | 15 / 2.50 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 389281527 0 / 08714729370093    | 15 / 2.75 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 389281530 0 / 08714729370109    | 15 / 3.00 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 389281532 0 / 08714729370116    | 15 / 3.25 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 389281535 0 / 08714729370123    | 15 / 3.50 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 389281537 0 / 08714729370130    | 15 / 3.75 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 389281540 0 / 08714729370147    | 15 / 4.00 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 389282015 0 / 08714729370154    | 20 / 1.50 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 389282020 0 / 08714729370161    | 20 / 2.00 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 389282022 0 / 08714729370178    | 20 / 2.25 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 389282025 0 / 08714729370185    | 20 / 2.50 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 389282027 0 / 08714729370192    | 20 / 2.75 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 389282030 0 / 08714729370208    | 20 / 3.00 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 389282032 0 / 08714729370215    | 20 / 3.25 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 389282035 0 / 08714729370222    | 20 / 3.50 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 389282037 0 / 08714729370239    | 20 / 3.75 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 389282040 0 / 08714729370246    | 20 / 4.00 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 389282520 0 / 08714729370253    | 25 / 2.00 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 389282522 0 / 08714729370260    | 25 / 2.25 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 389282525 0 / 08714729370277    | 25 / 2.50 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 389282527 0 / 08714729370284    | 25 / 2.75 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 389282530 0 / 08714729370291    | 25 / 3.00 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 389282532 0 / 08714729370307    | 25 / 3.25 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 389282535 0 / 08714729370314    | 25 / 3.50 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 389282537 0 / 08714729370321    | 25 / 3.75 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 389282540 0 / 08714729370338    | 25 / 4.00 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 389283020 0 / 08714729370345    | 30 / 2.00 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 389283022 0 / 08714729370352    | 30 / 2.25 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 389283025 0 / 08714729370369    | 30 / 2.50 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 389283027 0 / 08714729370376    | 30 / 2.75 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 389283030 0 / 08714729370383    | 30 / 3.00 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 389283032 0 / 08714729370390    | 30 / 3.25 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 14 ATM (1419 kPa)    |
| H749 389283035 0 / 08714729370406    | 30 / 3.50 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 389283037 0 / 08714729370413    | 30 / 3.75 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |
| H749 389283040 0 / 08714729370420    | 30 / 4.00 mm      | 0.017" (0.43 mm)     | 6 ATM (608 kPa)  | 12 ATM (1216 kPa)    |



## 2. Post-dilatation Balloon Catheters

### NC QUANTUM APEX™

| CATALOG NUMBER / ORDER NUMBER (GTIN) | LENGTH / DIAMETER | LESION ENTRY PROFILE | NOMINAL PRESSURE  | RATED BURST PRESSURE |
|--------------------------------------|-------------------|----------------------|-------------------|----------------------|
| H749 391240620 0 / 08714729783008    | 6 / 2.00 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391240622 0 / 08714729783015    | 6 / 2.25 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391240625 0 / 08714729783022    | 6 / 2.50 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391240627 0 / 08714729783039    | 6 / 2.75 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391240630 0 / 08714729783046    | 6 / 3.00 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391240632 0 / 08714729783053    | 6 / 3.25 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391240635 0 / 08714729783060    | 6 / 3.50 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391240637 0 / 08714729783077    | 6 / 3.75 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391240640 0 / 08714729783084    | 6 / 4.00 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391240645 0 / 08714729783091    | 6 / 4.50 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 18 ATM (1824 kPa)    |
| H749 391240650 0 / 08714729783107    | 6 / 5.00 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 18 ATM (1824 kPa)    |
| H749 391240820 0 / 08714729783114    | 8 / 2.00 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391240822 0 / 08714729783121    | 8 / 2.25 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391240825 0 / 08714729783138    | 8 / 2.50 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391240827 0 / 08714729783145    | 8 / 2.75 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391240830 0 / 08714729783152    | 8 / 3.00 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391240832 0 / 08714729783169    | 8 / 3.25 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391240835 0 / 08714729783176    | 8 / 3.50 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391240837 0 / 08714729783183    | 8 / 3.75 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391240840 0 / 08714729783190    | 8 / 4.00 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391240845 0 / 08714729783206    | 8 / 4.50 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 18 ATM (1824 kPa)    |
| H749 391240850 0 / 08714729783213    | 8 / 5.00 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 18 ATM (1824 kPa)    |
| H749 391241220 0 / 08714729783220    | 12 / 2.00 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391241222 0 / 08714729783237    | 12 / 2.25 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391241225 0 / 08714729783244    | 12 / 2.50 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391241227 0 / 08714729783251    | 12 / 2.75 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391241230 0 / 08714729783268    | 12 / 3.00 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391241232 0 / 08714729783275    | 12 / 3.25 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391241235 0 / 08714729783282    | 12 / 3.50 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391241237 0 / 08714729783299    | 12 / 3.75 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391241240 0 / 08714729783305    | 12 / 4.00 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391241245 0 / 08714729783312    | 12 / 4.50 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 18 ATM (1824 kPa)    |
| H749 391241250 0 / 08714729783329    | 12 / 5.00 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 18 ATM (1824 kPa)    |
| H749 391241520 0 / 08714729783336    | 15 / 2.00 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391241522 0 / 08714729783343    | 15 / 2.25 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391241525 0 / 08714729783350    | 15 / 2.50 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391241527 0 / 08714729783367    | 15 / 2.75 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391241530 0 / 08714729783374    | 15 / 3.00 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391241532 0 / 08714729783381    | 15 / 3.25 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391241535 0 / 08714729783398    | 15 / 3.50 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391241537 0 / 08714729783404    | 15 / 3.75 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391241540 0 / 08714729783411    | 15 / 4.00 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391241545 0 / 08714729783428    | 15 / 4.50 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 18 ATM (1824 kPa)    |
| H749 391241550 0 / 08714729783435    | 15 / 5.00 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 18 ATM (1824 kPa)    |

Continued next page ➤





## 2. Post-dilatation Balloon Catheters

### NC QUANTUM APEX™ Continued

| CATALOG NUMBER / ORDER NUMBER (GTIN) | LENGTH / DIAMETER | LESION ENTRY PROFILE | NOMINAL PRESSURE  | RATED BURST PRESSURE |
|--------------------------------------|-------------------|----------------------|-------------------|----------------------|
| H749 391242020 0 / 08714729783336    | 20 / 2.00 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391242022 0 / 08714729783343    | 20 / 2.25 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391242025 0 / 08714729783350    | 20 / 2.50 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391242027 0 / 08714729783367    | 20 / 2.75 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391242030 0 / 08714729783374    | 20 / 3.00 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391242032 0 / 08714729783381    | 20 / 3.25 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391242035 0 / 08714729783398    | 20 / 3.50 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391242037 0 / 08714729783404    | 20 / 3.75 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391242040 0 / 08714729783411    | 20 / 4.00 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391242045 0 / 08714729783428    | 20 / 4.50 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 18 ATM (1824 kPa)    |
| H749 391242050 0 / 08714729783435    | 20 / 5.00 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 18 ATM (1824 kPa)    |
| H749 391243020 0 / 08714729783374    | 30 / 2.00 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391243022 0 / 08714729783381    | 30 / 2.25 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391243025 0 / 08714729783398    | 30 / 2.50 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391243027 0 / 08714729783404    | 30 / 2.75 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391243030 0 / 08714729783411    | 30 / 3.00 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391243032 0 / 08714729783428    | 30 / 3.25 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391243035 0 / 08714729783435    | 30 / 3.50 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391243037 0 / 08714729783428    | 30 / 3.75 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 391243040 0 / 08714729783435    | 30 / 4.00 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |



## 2. Post-dilatation Balloon Catheters

### NC EMERGE™ Monorail™

| CATALOG NUMBER / ORDER NUMBER (GTIN) | LENGTH / DIAMETER | LESION ENTRY PROFILE | NOMINAL PRESSURE  | RATED BURST PRESSURE |
|--------------------------------------|-------------------|----------------------|-------------------|----------------------|
| H749 392760620 0 / 08714729846833    | 6 / 2.00 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392760622 0 / 08714729846840    | 6 / 2.25 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392760625 0 / 08714729846857    | 6 / 2.50 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392760627 0 / 08714729846864    | 6 / 2.75 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392760630 0 / 08714729846871    | 6 / 3.00 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392760632 0 / 08714729846888    | 6 / 3.25 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392760635 0 / 08714729846895    | 6 / 3.50 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392760637 0 / 08714729846901    | 6 / 3.75 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392760640 0 / 08714729846918    | 6 / 4.00 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392760645 0 / 08714729846925    | 6 / 4.50 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 18 ATM (1824 kPa)    |
| H749 392760650 0 / 08714729846932    | 6 / 5.00 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 18 ATM (1824 kPa)    |
| H749 392760820 0 / 08714729846949    | 8 / 2.00 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392760822 0 / 08714729846956    | 8 / 2.25 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392760825 0 / 08714729846963    | 8 / 2.50 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392760827 0 / 08714729846970    | 8 / 2.75 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392760830 0 / 08714729846987    | 8 / 3.00 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392760832 0 / 08714729846994    | 8 / 3.25 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392760835 0 / 08714729847007    | 8 / 3.50 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392760837 0 / 08714729847014    | 8 / 3.75 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392760840 0 / 08714729847021    | 8 / 4.00 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392760845 0 / 08714729847038    | 8 / 4.50 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 18 ATM (1824 kPa)    |
| H749 392760850 0 / 08714729847045    | 8 / 5.00 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 18 ATM (1824 kPa)    |
| H749 392760855 0 / 08714729847052    | 8 / 5.50 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 18 ATM (1824 kPa)    |
| H749 392760860 0 / 08714729847069    | 8 / 6.00 mm       | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 18 ATM (1824 kPa)    |
| H749 392761220 0 / 08714729847076    | 12 / 2.00 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392761222 0 / 08714729847083    | 12 / 2.25 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392761225 0 / 08714729847090    | 12 / 2.50 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392761227 0 / 08714729847106    | 12 / 2.75 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392761230 0 / 08714729847113    | 12 / 3.00 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392761232 0 / 08714729847120    | 12 / 3.25 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392761235 0 / 08714729847137    | 12 / 3.50 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392761237 0 / 08714729847144    | 12 / 3.75 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392761240 0 / 08714729847151    | 12 / 4.00 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392761245 0 / 08714729847168    | 12 / 4.50 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 18 ATM (1824 kPa)    |
| H749 392761250 0 / 08714729847175    | 12 / 5.00 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 18 ATM (1824 kPa)    |
| H749 392761255 0 / 08714729847182    | 12 / 5.50 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 18 ATM (1824 kPa)    |
| H749 392761260 0 / 08714729847199    | 12 / 6.00 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 18 ATM (1824 kPa)    |
| H749 392761520 0 / 08714729847205    | 15 / 2.00 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392761522 0 / 08714729847212    | 15 / 2.25 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392761525 0 / 08714729847229    | 15 / 2.50 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392761527 0 / 08714729847236    | 15 / 2.75 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392761530 0 / 08714729847243    | 15 / 3.00 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392761532 0 / 08714729847250    | 15 / 3.25 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392761535 0 / 08714729847267    | 15 / 3.50 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392761537 0 / 08714729847274    | 15 / 3.75 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392761540 0 / 08714729847281    | 15 / 4.00 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392761545 0 / 08714729847298    | 15 / 4.50 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 18 ATM (1824 kPa)    |
| H749 392761550 0 / 08714729847304    | 15 / 5.00 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 18 ATM (1824 kPa)    |
| H749 392761555 0 / 08714729847311    | 15 / 5.50 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 18 ATM (1824 kPa)    |
| H749 392761560 0 / 08714729847328    | 15 / 6.00 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 18 ATM (1824 kPa)    |

Continued next page >



Complex PCI - 24  
NEW CARDIO

Boston  
Scientific

Balloon  
Catheters

Micro-  
dilatation  
Catheters

Guide  
Catheters

Guide  
Catheter  
extension

Angiographic  
Catheters

Guidewires

Embolic  
Protection

Accessories

Cutting  
Balloon™

Atherectomy  
Devices


Chronic Total  
Occlusions

Trapper™  
Exchange  
Device

ANGIOJET™  
Thrombec-  
tomy  
System

2. Post-dilatation Balloon Catheters

| NC EMERGE™ Monorail™ Continued       |                   |                      |                   |                      |
|--------------------------------------|-------------------|----------------------|-------------------|----------------------|
| CATALOG NUMBER / ORDER NUMBER (GTIN) | LENGTH / DIAMETER | LESION ENTRY PROFILE | NOMINAL PRESSURE  | RATED BURST PRESSURE |
| H749 392762020 0 / 08714729847335    | 20 / 2.00 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392762022 0 / 08714729847342    | 20 / 2.25 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392762025 0 / 08714729847359    | 20 / 2.50 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392762027 0 / 08714729847366    | 20 / 2.75 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392762030 0 / 08714729847373    | 20 / 3.00 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392762032 0 / 08714729847380    | 20 / 3.25 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392762035 0 / 08714729847397    | 20 / 3.50 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392762037 0 / 08714729847403    | 20 / 3.75 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392762040 0 / 08714729847410    | 20 / 4.00 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392762045 0 / 08714729847427    | 20 / 4.50 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 18 ATM (1824 kPa)    |
| H749 392762050 0 / 08714729847434    | 20 / 5.00 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 18 ATM (1824 kPa)    |
| H749 392762055 0 / 08714729847441    | 20 / 5.50 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 18 ATM (1824 kPa)    |
| H749 392762060 0 / 08714729847458    | 20 / 6.00 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 18 ATM (1824 kPa)    |
| H749 392763020 0 / 08714729847465    | 30 / 2.00 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392763022 0 / 08714729847472    | 30 / 2.25 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392763025 0 / 08714729847489    | 30 / 2.50 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392763027 0 / 08714729847496    | 30 / 2.75 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392763030 0 / 08714729847502    | 30 / 3.00 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392763032 0 / 08714729847519    | 30 / 3.25 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392763035 0 / 08714729847526    | 30 / 3.50 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392763037 0 / 08714729847533    | 30 / 3.75 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |
| H749 392763040 0 / 08714729847540    | 30 / 4.00 mm      | 0.017" (0.43 mm)     | 12 ATM (1216 kPa) | 20 ATM (2026 kPa)    |



Complex PCI - 25  
NEW CARDIO

Boston  
Scientific

Balloon  
Catheters

Micro-  
dilatation  
Catheters

Guide  
Catheters

Guide  
Catheter  
extension

Angiographic  
Catheters

Guidewires

Embolic  
Protection

Accessories

Cutting  
Balloon™

Atherectomy  
Devices

Chronic Total  
Occlusions

Trapper™  
Exchange  
Device

ANGIOJET™  
Thrombec-  
tomy  
System

NEW CARDIO

2

## COMPLEX PCI

Balloon Catheters

Microdilatation Catheters

Guide Catheters

Guide Catheter Extension

Angiographic Catheters

Guidewires

Embolic Protection

Accessories

Cutting Balloon™

Atherectomy Devices

Chronic Total Occlusions

TRAPPER™ Exchange Device

ANGIOJET™ Thrombectomy System



Balloon  
Catheters

Micro-  
dilatation  
Catheters

Guide  
Catheters

Guide  
Catheter  
extension

Angiographic  
Catheters

Guidewires

Embolic  
Protection

Accessories

Cutting  
Balloon™

Atherectomy  
Devices

Chronic Total  
Occlusions

Trapper™  
Exchange  
Device

ANGIOJET™  
Thrombec-  
tomy  
System



2. Microdilatation Catheters

THREADER™  
Microdilatation Catheter

| CATALOG NUMBER / ORDER NUMBER (GTIN) | SIZE           | TYPE          |
|--------------------------------------|----------------|---------------|
| H749 392691212 0 / 8714729848127     | 1.2 mm x 12 mm | Monorail™     |
| H749 392701212 0 / 8714729848134     | 1.2 mm x 12 mm | Over-the-Wire |

MAMBA™  
Microdilatation Catheter

| CATALOG NUMBER / ORDER NUMBER (GTIN)                                                | LESION ENTRY PROFILE | CROSSING PROFILE | PROXIMAL SHAFT PROFILE | DISTAL I.D. (IN/MM) | PROXIMAL I.D. (IN/MM) | USABLE LENGTH (CM) | COMPATIBLE GUIDEWIRE DIAMETER (IN/MM) |
|-------------------------------------------------------------------------------------|----------------------|------------------|------------------------|---------------------|-----------------------|--------------------|---------------------------------------|
| H749 392871352 0 / 8714729940265                                                    | 1.4F                 | 2.4F             | 2.9F                   | 0.018/.46           | 0.023/.58             | 135                | 0.014/.36                             |
|   |                      |                  |                        |                     |                       |                    |                                       |
| H749 392871356 0 / 08714729940302                                                   | 1.4F                 | 2.1F             | 2.9F                   | 0.018/.46           | 0.023/.58             | 135                | 0.014/.36                             |
|   |                      |                  |                        |                     |                       |                    |                                       |
| H749 392871502 0 / 08714729940340                                                   | 1.4F                 | 2.1F             | 2.9F                   | 0.018/.46           | 0.023/.58             | 150                | 0.014/.36                             |
|  |                      |                  |                        |                     |                       |                    |                                       |



NEW CARDIO

2

## COMPLEX PCI

Balloon Catheters

Microdilatation Catheters

Guide Catheters

Guide Catheter Extension

Angiographic Catheters

Guidewires

Embolic Protection

Accessories

Cutting Balloon™

Atherectomy Devices

Chronic Total Occlusions

TRAPPER™ Exchange Device

ANGIOJET™ Thrombectomy System



Balloon  
Catheters

Micro-  
dilatation  
Catheters

Guide  
Catheters

Guide  
Catheter  
extension

Angiographic  
Catheters

Guidewires

Embolic  
Protection

Accessories

Cutting  
Balloon™

Atherectomy  
Devices


Chronic Total  
Occlusions

Trapper™  
Exchange  
Device


ANGIOJET™  
Thrombec-  
tomy  
System

## Curves Overview


### Left Coronary Curves




### Right Coronary Curves



### Radial Curves



### Additional Curves



## 2. Guide Catheters

Convey Guide Catheters will be discontinued in November 2023 except in the countries where FDA certification only can be applied.

### CONVEY™ \*

#### Femoral Approach Left Curve

|                    | CURVE STYLE | 5F = 1.66 MM<br>(.058" = 1.47 MM ID) | 6F = 2.0 MM (.071" = 1.80 MM ID)<br>NO SIDE HOLES | SIDE HOLES                        | 7F = 2.33 MM (.082" = 2.08 MM ID)<br>NO SIDE HOLES | SIDE HOLES                        | 8F = 2.66 MM (.091" = 2.31 MM ID)<br>NO SIDE HOLES | SIDE HOLES                        |
|--------------------|-------------|--------------------------------------|---------------------------------------------------|-----------------------------------|----------------------------------------------------|-----------------------------------|----------------------------------------------------|-----------------------------------|
|                    |             |                                      |                                                   |                                   |                                                    |                                   |                                                    |                                   |
| Judkins Left       | JL3         | H749 39228660 0<br>08714729797876    | H749 39264660 0<br>08714729849575                 | H749 39264666 0<br>08714729849636 | -<br>-                                             | H749 39265666 0<br>08714729851233 | -<br>-                                             | -<br>-                            |
|                    | JL3.5       | H749 39228661 0<br>08714729797883    | H749 39264661 0<br>08714729849582                 | H749 39264667 0<br>08714729849643 | H749 39265661 0<br>08714729851189                  | H749 39265667 0<br>08714729851240 | H749 39266661 0<br>08714729852902                  | H749 39266667 0<br>08714729852964 |
|                    | JL4         | H749 39228662 0<br>08714729797890    | H749 39264662 0<br>08714729849599                 | H749 39264668 0<br>08714729849650 | H749 39265662 0<br>08714729851196                  | H749 39265668 0<br>08714729851257 | H749 39266662 0<br>08714729852919                  | H749 39266668 0<br>08714729852971 |
|                    | JL4.5       | H749 39228663 0<br>08714729797906    | H749 39264663 0<br>08714729849605                 | H749 39264669 0<br>08714729849667 | H749 39265663 0<br>08714729851202                  | H749 39265669 0<br>08714729851264 | H749 39266663 0<br>08714729852926                  | H749 39266669 0<br>08714729852988 |
|                    | JL5         | H749 39228664 0<br>08714729797913    | H749 39264664 0<br>08714729849612                 | H749 39264670 0<br>08714729849681 | H749 39265664 0<br>08714729851219                  | H749 39265670 0<br>08714729851288 | H749 39266664 0<br>08714729852933                  | H749 39266670 0<br>08714729853008 |
|                    | JL6         | H749 39228665 0<br>08714729797920    | H749 39264665 0<br>08714729849629                 | H749 39264671 0<br>08714729849698 | H749 39265665 0<br>08714729851226                  | H749 39265671 0<br>08714729851295 | H749 39266665 0<br>08714729852940                  | -<br>-                            |
| Femoral Left       | FL3         | H749 3922805 0<br>08714729797937     | H749 3926405 0<br>08714729890294                  | -<br>-                            | -<br>-                                             | -<br>-                            | -<br>-                                             | -<br>-                            |
|                    | FL3.5       | H749 3922806 0<br>08714729797944     | H749 3926406 0<br>08714729890300                  | -<br>-                            | -<br>-                                             | -<br>-                            | -<br>-                                             | -<br>-                            |
|                    | FL4         | H749 3922807 0<br>08714729797951     | H749 3926407 0<br>08714729890317                  | -<br>-                            | -<br>-                                             | -<br>-                            | -<br>-                                             | -<br>-                            |
|                    | FL4.5       | H749 3922808 0<br>08714729797968     | H749 3926408 0<br>08714729890324                  | -<br>-                            | -<br>-                                             | -<br>-                            | -<br>-                                             | -<br>-                            |
|                    | FL5         | H749 3922809 0<br>08714729797975     | H749 3926409 0<br>08714729890331                  | -<br>-                            | -<br>-                                             | -<br>-                            | -<br>-                                             | -<br>-                            |
| Femoral Curve Left | FCL3.5      | H749 39228272 0<br>08714729798002    | H749 39264272 0<br>08714729849148                 | -<br>-                            | -<br>-                                             | -<br>-                            | -<br>-                                             | -<br>-                            |
|                    | FCL4        | H749 39228273 0<br>08714729798019    | H749 39264273 0<br>08714729849155                 | -<br>-                            | -<br>-                                             | -<br>-                            | -<br>-                                             | -<br>-                            |
| Amplatz™ Left      | AL.75       | H749 3922822 0<br>08714729798187     | H749 3926422 0<br>08714729849063                  | H749 39264124 0<br>08714729890362 | H749 3926522 0<br>08714729850663                   | H749 39265124 0<br>08714729850564 | H749 3926622 0<br>08714729852346                   | H749 39266124 0<br>08714729852162 |
|                    | AL1         | H749 3922819 0<br>08714729798194     | H749 3926419 0<br>08714729849032                  | H749 3926468 0<br>08714729849704  | H749 3926519 0<br>08714729850632                   | H749 3926568 0<br>08714729851301  | H749 3926619 0<br>08714729852230                   | H749 3926668 0<br>08714729853022  |
|                    | AL2         | H749 3922820 0<br>08714729798217     | H749 3926420 0<br>08714729849049                  | H749 3926469 0<br>08714729849773  | H749 3926520 0<br>08714729850649                   | H749 3926569 0<br>08714729851370  | H749 3926620 0<br>08714729852285                   | H749 3926669 0<br>08714729853091  |

Continued next page >



Complex PCI - 30  
NEW CARDIO

Boston  
Scientific

\*The CONVEY™ Guiding Catheter is manufactured by PendraCare International B.V. and distributed by Boston Scientific Corporation.

Balloon  
Catheters

Micro-  
dilatation  
Catheters

Guide  
Catheters

Guide  
Catheter  
extension

Angiographic  
Catheters

Guidewires

Embolic  
Protection

Accessories

Cutting  
Balloon™

Atherectomy  
Devices

Chronic Total  
Occlusions

Trapper™  
Exchange  
Device

ANGIOJET™  
Thrombec-  
tomy  
System



## 2. Guide Catheters

Convey Guide Catheters will be discontinued in November 2023 except in the countries where FDA certification only can be applied.

### CONVEY™\* Continued Femoral Approach Left Curve

|                                       | CURVE<br>STYLE | 5F = 1.66 MM<br>(.058" = 1.47 MM ID) | 6F = 2.0 MM (.071" = 1.80 MM ID)<br>NO SIDE HOLES | 6F = 2.0 MM (.071" = 1.80 MM ID)<br>SIDE HOLES | 7F = 2.33 MM (.082" = 2.08 MM ID)<br>NO SIDE HOLES | 7F = 2.33 MM (.082" = 2.08 MM ID)<br>SIDE HOLES | 8F = 2.66 MM (.091" = 2.31 MM ID)<br>NO SIDE HOLES | 8F = 2.66 MM (.091" = 2.31 MM ID)<br>SIDE HOLES |
|---------------------------------------|----------------|--------------------------------------|---------------------------------------------------|------------------------------------------------|----------------------------------------------------|-------------------------------------------------|----------------------------------------------------|-------------------------------------------------|
| Contra-<br>lateral<br>Left<br>Support | CLS3           | H749 39228573 0<br>08714729798255    | H749 39264573 0<br>08714729849469                 | H749 39264577 0<br>08714729849506              | H749 39265573 0<br>08714729851066                  | H749 39265577 0<br>08714729851103               | H749 39266573 0<br>08714729852780                  | H749 39266577 0<br>08714729852827               |
|                                       | CLS3.5         | H749 39228574 0<br>08714729798279    | H749 39264574 0<br>08714729849476                 | H749 39264578 0<br>08714729849513              | H749 39265574 0<br>08714729851073                  | H749 39265578 0<br>08714729851110               | H749 39266574 0<br>08714729852797                  | H749 39266578 0<br>08714729852834               |
|                                       | CLS4           | H749 39228575 0<br>08714729798293    | H749 39264575 0<br>08714729849483                 | H749 39264579 0<br>08714729849520              | H749 39265575 0<br>08714729851080                  | H749 39265579 0<br>08714729851127               | H749 39266575 0<br>08714729852803                  | H749 39266579 0<br>08714729852841               |
| Left<br>Back-Up                       | LeftBU3        | H749 39228713 0<br>08714729798323    | H749 39264713 0<br>08714729849858                 | H749 39264719 0<br>08714729849919              | H749 39265713 0<br>08714729851455                  | H749 39265719 0<br>08714729851516               | H749 39266713 0<br>08714729853176                  | H749 39266719 0<br>08714729853237               |
|                                       | LeftBU3.25     | H749 39228729 0<br>08714729798330    | H749 39264729 0<br>08714729850014                 | -                                              | H749 39265729 0<br>08714729851615                  | -                                               | -                                                  | -                                               |
|                                       | LeftBU3.5      | H749 39228714 0<br>08714729798347    | H749 39264714 0<br>08714729849865                 | H749 39264720 0<br>08714729849926              | H749 39265714 0<br>08714729851462                  | H749 39265720 0<br>08714729851523               | H749 39266714 0<br>08714729853183                  | H749 39266720 0<br>08714729853244               |
|                                       | LeftBU3.75     | H749 39228715 0<br>08714729798354    | H749 39264715 0<br>08714729849872                 | -                                              | H749 39265715 0<br>08714729851479                  | -                                               | H749 39266715 0<br>08714729853190                  | -                                               |
|                                       | LeftBU4        | H749 39228716 0<br>08714729798361    | H749 39264716 0<br>08714729849889                 | H749 39264722 0<br>08714729849940              | H749 39265716 0<br>08714729851486                  | H749 39265722 0<br>08714729851547               | H749 39266716 0<br>08714729853206                  | H749 39266722 0<br>08714729853268               |
|                                       | LeftBU4.25     | -                                    | H749 39264730 0<br>08714729850021                 | -                                              | H749 39265730 0<br>08714729851622                  | -                                               | -                                                  | -                                               |
|                                       | LeftBU4.5      | H749 39228717 0<br>08714729798736    | H749 39264717 0<br>08714729849896                 | H749 39264723 0<br>08714729849957              | H749 39265717 0<br>08714729851493                  | H749 39265723 0<br>08714729851554               | H749 39266717 0<br>08714729853213                  | H749 39266723 0<br>08714729853275               |
|                                       | LeftBU4.75     | -                                    | H749 39264731 0<br>08714729850038                 | -                                              | -                                                  | -                                               | -                                                  | -                                               |
|                                       | LeftBU5        | H749 39228718 0<br>08714729798392    | H749 39264718 0<br>08714729849902                 | H749 39264724 0<br>08714729849964              | H749 39265718 0<br>08714729851509                  | H749 39265724 0<br>08714729851561               | -                                                  | -                                               |



Complex PCI - 31  
NEW CARDIO

Boston  
Scientific

\*The CONVEY™ Guiding Catheter is manufactured by PendraCare International B.V.  
and distributed by Boston Scientific Corporation.

Balloon  
Catheters

Micro-  
dilatation  
Catheters

Guide  
Catheters

Guide  
Catheter  
extension

Angiographic  
Catheters

Guidewires

Embolic  
Protection

Accessories

Cutting  
Balloon™

Atherectomy  
Devices

Chronic Total  
Occlusions

Trapper™  
Exchange  
Device

ANGIOJET™  
Thrombec-  
tomy  
System

## 2. Guide Catheters

Convey Guide Catheters will be discontinued in November 2023 except in the countries where FDA certification only can be applied.

### CONVEY™ \*

#### Femoral Approach Right Curve

| CURVE STYLE         | 5F = 1.66 MM (.058" = 1.47 MM ID) |                                   | 6F = 2.0 MM (.071" = 1.80 MM ID)  |                                   | 7F = 2.33 MM (.082" = 2.08 MM ID) |                                   | 8F = 2.66 MM (.091" = 2.31 MM ID) |                                   |
|---------------------|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|
|                     | NO SIDE HOLES                     |                                   | SIDE HOLES                        |                                   | NO SIDE HOLES                     |                                   | SIDE HOLES                        |                                   |
| JR3                 | H749 39228684 0<br>08714729798040 | H749 39264684 0<br>08714729849711 | H749 39264690 0<br>08714729849780 | -                                 | -                                 | -                                 | -                                 | -                                 |
| JR3.5               | H749 39228685 0<br>08714729798057 | H749 39264685 0<br>08714729849728 | H749 39264691 0<br>08714729849797 | H749 39265685 0<br>08714729851325 | H749 39265691 0<br>08714729851394 | H749 39266685 0<br>08714729853046 | H749 39266691 0<br>08714729853114 |                                   |
| JR4                 | H749 39228686 0<br>08714729798064 | H749 39264686 0<br>08714729849735 | H749 39264692 0<br>08714729849803 | H749 39265686 0<br>08714729851332 | H749 39265692 0<br>08714729851400 | H749 39266686 0<br>08714729853053 | H749 39266692 0<br>08714729853121 |                                   |
| JR4.5               | H749 39228687 0<br>08714729798071 | H749 39264687 0<br>08714729849742 | -                                 | H749 39265687 0<br>08714729851349 | -                                 | -                                 | -                                 |                                   |
| JR5                 | H749 39228688 0<br>08714729798088 | H749 39264688 0<br>08714729849759 | H749 39264694 0<br>08714729849827 | H749 39265688 0<br>08714729851356 | H749 39265694 0<br>08714729851424 | -                                 | H749 39266694 0<br>08714729853145 |                                   |
| JR6                 | H749 39228689 0<br>08714729798095 | H749 39264689 0<br>08714729849766 | -                                 | H749 39265689 0<br>08714729851363 | -                                 | -                                 | -                                 |                                   |
| FR3                 | -                                 | H749 39264442 0<br>08714729849193 | -                                 | -                                 | -                                 | -                                 | -                                 |                                   |
| FR3.5               | H749 3922801 0<br>08714729798118  | H749 3926401 0<br>08714729890256  | -                                 | -                                 | -                                 | -                                 | -                                 |                                   |
| FR4                 | H749 3922802 0<br>08714729798125  | H749 3926402 0<br>08714729890263  | -                                 | -                                 | -                                 | -                                 | -                                 |                                   |
| FR5                 | H749 3922803 0<br>08714729798149  | H749 3926403 0<br>08714729890270  | -                                 | -                                 | -                                 | -                                 | -                                 |                                   |
| Multi-purpose       | MPA                               | H749 39228767 0<br>08714729798163 | H749 39264767 0<br>08714729850397 | H749 39264768 0<br>08714729850403 | H749 39265767 0<br>08714729851998 | H749 39265768 0<br>08714729852001 | H749 39266767 0<br>08714729853718 | H749 39266768 0<br>08714729853725 |
| Hockey Stick        | HS                                | H749 3922841 0<br>08714729798170  | H749 3926441 0<br>08714729849186  | H749 39264141 0<br>08714729849001 | H749 3926541 0<br>08714729850786  | -                                 | H749 3926641 0<br>08714729852506  | -                                 |
| Amplatz™ Right      | AR1                               | H749 3922817 0<br>08714729798231  | H749 3926417 0<br>08714729849018  | H749 3926466 0<br>08714729849568  | H749 3926517 0<br>08714729850618  | H749 3926566 0<br>08714729851165  | H749 3926617 0<br>08714729852216  | H749 3926666 0<br>08714729852889  |
|                     | AR2                               | H749 3922818 0<br>08714729798248  | H749 3926418 0<br>08714729849025  | H749 3926467 0<br>08714729849674  | H749 3926518 0<br>08714729850625  | H749 3926567 0<br>08714729851271  | H749 3926618 0<br>08714729852223  | H749 3926667 0<br>08714729852995  |
| allRight™ Curve     | ART3                              | H749 39228496 0<br>08714729798408 | H749 39264496 0<br>08714729849292 | H749 39264497 0<br>08714729849308 | H749 39265496 0<br>08714729850892 | -                                 | -                                 | -                                 |
|                     | ART3.5                            | H749 39228498 0<br>08714729798415 | H749 39264498 0<br>08714729849315 | H749 39264499 0<br>08714729849322 | H749 39265498 0<br>08714729850915 | H749 39265499 0<br>08714729850922 | -                                 | -                                 |
|                     | ART4                              | H749 39228500 0<br>08714729798439 | H749 39264500 0<br>08714729849346 | H749 39264501 0<br>08714729849353 | H749 39265500 0<br>08714729850946 | H749 39265501 0<br>08714729850953 | H749 39266500 0<br>08714729852667 | -                                 |
| RC Shepherd's Crook | RC3.5 SC                          | H749 3922846 0<br>08714729798460  | H749 3926446 0<br>08714729849230  | -                                 | -                                 | -                                 | -                                 | -                                 |



Complex PCI - 32  
NEW CARDIO

Boston  
Scientific

\*The CONVEY™ Guiding Catheter is manufactured by PendraCare International B.V. and distributed by Boston Scientific Corporation.

Balloon  
Catheters

Micro-  
dilatation  
Catheters

Guide  
Catheters

Guide  
Catheter  
extension

Angiographic  
Catheters

Guidewires

Embolic  
Protection

Accessories

Cutting  
Balloon™

Atherectomy  
Devices

Chronic Total  
Occlusions

Trapper™  
Exchange  
Device

ANGIOJET™  
Thrombec-  
tomy  
System

2. Guide Catheters

Convey Guide Catheters will be discontinued in November 2023 except in the countries where FDA certification only can be applied.

| CONVEY™ *            |             |                                   |                                   |                                   |                                   |                                   |
|----------------------|-------------|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|
| Transradial Approach |             |                                   |                                   |                                   |                                   |                                   |
|                      | CURVE STYLE | 5F = 1.66 MM (.058" = 1.47 MM ID) | 6F = 2.0 MM (.071" = 1.80 MM ID)  |                                   | 7F = 2.33 MM (.082" = 2.08 MM ID) |                                   |
|                      |             |                                   | NO SIDE HOLES                     | SIDE HOLES                        | NO SIDE HOLES                     | SIDE HOLES                        |
| Radial               | RADIAL      | H749 39228472 0<br>08714729798491 | H749 39264472 0<br>08714729849254 | H749 39264473 0<br>08714729849261 | H749 39265472 0<br>08714729850854 | -<br>-                            |
| Radial Back-Up       | RB          | H749 39228769 0<br>08714729798507 | H749 39264769 0<br>08714729850410 | -<br>-                            | H749 39265769 0<br>08714729852018 | -<br>-                            |
| Radial Back-Up Left  | RBL3.5      | H749 39228735 0<br>08714729798514 | H749 39264735 0<br>08714729850076 | H749 39264756 0<br>08714729850281 | H749 39265735 0<br>08714729851677 | H749 39265756 0<br>08714729851882 |
|                      | RBL4        | H749 39228736 0<br>08714729798521 | H749 39264736 0<br>08714729850083 | H749 39264757 0<br>08714729850298 | H749 39265736 0<br>08714729851684 | H749 39265757 0<br>08714729851899 |
|                      | RBL4.5      | H749 39228737 0<br>08714729798538 | H749 39264737 0<br>08714729850090 | H749 39264758 0<br>08714729850304 | H749 39265737 0<br>08714729851691 | H749 39265758 0<br>08714729851905 |
| Radial Back-Up Right | RBR3        | H749 39228738 0<br>08714729798545 | H749 39264738 0<br>08714729850106 | -<br>-                            | H749 39265738 0<br>08714729851707 | -<br>-                            |
|                      | RBR3.5      | H749 39228739 0<br>08714729798552 | H749 39264739 0<br>08714729850113 | H749 39264760 0<br>08714729850328 | H749 39265739 0<br>08714729851714 | -<br>-                            |
|                      | RBR4        | H749 39228740 0<br>08714729798569 | H749 39264740 0<br>08714729850120 | H749 39264761 0<br>08714729850335 | H749 39265740 0<br>08714729851721 | H749 39265761 0<br>08714729851936 |
| TIG                  | TIG3        | H749 39228741 0<br>08714729798576 | H749 39264741 0<br>08714729850137 | H749 39264762 0<br>08714729850342 | -<br>-                            | -<br>-                            |
|                      | TIG3.5      | H749 39228742 0<br>08714729798583 | H749 39264742 0<br>08714729850144 | H749 39264763 0<br>08714729850359 | -<br>-                            | -<br>-                            |
|                      | TIG4        | H749 39228743 0<br>08714729798590 | H749 39264743 0<br>08714729850151 | H749 39264764 0<br>08714729850366 | -<br>-                            | -<br>-                            |
|                      | TIG4.5      | H749 39228744 0<br>08714729798606 | H749 39264744 0<br>08714729850168 | H749 39264765 0<br>08714729850373 | H749 39265744 0<br>08714729851769 | -<br>-                            |
|                      | TIG5        | H749 39228745 0<br>08714729798613 | H749 39264745 0<br>08714729850175 | H749 39264766 0<br>08714729850380 | -<br>-                            | -<br>-                            |
| Kimny™               | Kimny™      | H749 39228583 0<br>08714729798620 | H749 39264583 0<br>08714729849544 | -<br>-                            | -<br>-                            | -<br>-                            |

Convey Guide Catheters will be discontinued in November 2023 except in the countries where FDA certification only can be applied.

| CONVEY™ *               |             |                                   |                                  |                                   |                                   |                                   |                                   |                                   |
|-------------------------|-------------|-----------------------------------|----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|
| Mammary & Bypass Curves |             |                                   |                                  |                                   |                                   |                                   |                                   |                                   |
|                         | CURVE STYLE | 5F = 1.66 MM (.058" = 1.47 MM ID) | 6F = 2.0 MM (.070" = 1.78 MM ID) |                                   | 7F = 2.33 MM (.081" = 2.05 MM ID) |                                   | 8F = 2.66 MM (.091" = 2.31 MM ID) |                                   |
|                         |             |                                   | NO SIDE HOLES                    | SIDE HOLES                        | NO SIDE HOLES                     | SIDE HOLES                        | NO SIDE HOLES                     | SIDE HOLES                        |
| Internal Mammary        | IM          | H749 3922826 0<br>08714729798637  | H749 3926426 0<br>08714729849094 | H749 39264129 0<br>08714729848998 | H749 3926526 0<br>08714729850694  | H749 39265129 0<br>08714729850595 | H749 3926626 0<br>08714729852377  | H749 39266129 0<br>08714729852193 |
| Left Coronary Bypass    | LCB         | H749 3922825 0<br>08714729798644  | H749 3926425 0<br>08714729849087 | -<br>-                            | H749 3926525 0<br>08714729850687  | -<br>-                            | -<br>-                            | -<br>-                            |



## 2. Guide Catheters

### MACH 1™

#### Femoral Approach Left Curve

|                           | CURVE<br>STYLE | 6F = 2.0 MM (.070" = 1.78 MM ID)  |                                   | 7F = 2.33 MM (.081" = 2.05 MM ID) |                                   | 8F = 2.66 MM (.091" = 2.31 MM ID) |                                   |
|---------------------------|----------------|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|
|                           |                | NO SIDE HOLES                     | SIDE HOLES                        | NO SIDE HOLES                     | SIDE HOLES                        | NO SIDE HOLES                     | SIDE HOLES                        |
| Femoral Left              | FL3            | H749 3435605 0<br>08714729350446  | -<br>-                            |                                   |                                   |                                   |                                   |
|                           | FL3.5          | H749 3435606 0<br>08714729350453  | H749 3435652 0<br>08714729350781  | H749 3435706 0<br>08714729378846  | H749 3435752 0<br>08714729380191  | H749 3435806 0<br>08714729381501  | H749 3435852 0<br>08714729382935  |
|                           | FL4            | H749 3435607 0<br>08714729350460  | H749 3435653 0<br>-               | H749 3435707 0<br>08714729378853  | H749 3435753 0<br>08714729380306  | H749 3435807 0<br>08714729381518  | H749 3435853 0<br>08714729383048  |
|                           | FL4.5          | H749 3435608 0<br>08714729350477  | -<br>-                            | H749 3435708 0<br>08714729378860  | H749 3435754 0<br>08714729380351  | H749 3435808 0<br>08714729381525  | -<br>-                            |
|                           | FL5            | H749 3435609 0<br>08714729350484  | -<br>-                            | H749 3435709 0<br>08714729378877  | H749 3435755 0<br>08714729380368  | H749 3435809 0<br>08714729381532  | H749 3435855 0<br>08714729383109  |
|                           | FL6            | H749 3435610 0<br>08714729350491  | -<br>-                            | -<br>-                            | -<br>-                            | -<br>-                            | H749 3435856 0<br>08714729383116  |
| Femoral Left<br>Short Tip | FL3.5 ST       | H749 3435628 0<br>08714729351269  | -<br>-                            | -<br>-                            | -<br>-                            | -<br>-                            | -<br>-                            |
|                           | FL4 ST         | H749 3435629 0<br>08714729351283  | -<br>-                            | -<br>-                            | -<br>-                            | -<br>-                            | -<br>-                            |
|                           | FL4.5 ST       | H749 3435630 0<br>08714729351351  | -<br>-                            | -<br>-                            | -<br>-                            | -<br>-                            | -<br>-                            |
|                           | FL5 ST         | H749 3435631 0<br>08714729351368  | -<br>-                            | -<br>-                            | -<br>-                            | -<br>-                            | -<br>-                            |
| Femoral Curve<br>Left     | FCL3.5         | H749 34356272 0<br>08714729351238 | -<br>-                            | H749 34357272 0<br>08714729379584 | -<br>-                            | -<br>-                            | -<br>-                            |
|                           | FCL4           | H749 34356273 0<br>08714729351245 | -<br>-                            | -<br>-                            | -<br>-                            | H749 34358273 0<br>08714729382270 | -<br>-                            |
| Judkins Left              | JL3            | H749 34356660 0<br>08714729351849 | -<br>-                            | H749 34357660 0<br>08714729380818 | H749 34357666 0<br>08714729380870 | H749 34358660 0<br>08714729383536 | H749 34358666 0<br>08714729383598 |
|                           | JL3.5          | H749 34356661 0<br>08714729351856 | H749 34356667 0<br>08714729378747 | H749 34357661 0<br>08714729380825 | H749 34357667 0<br>08714729380887 | H749 34358661 0<br>08714729383543 | H749 34358667 0<br>08714729383604 |
|                           | JL4            | H749 34356662 0<br>08714729351863 | H749 34356668 0<br>08714729378754 | H749 34357662 0<br>08714729380832 | H749 34357668 0<br>08714729380894 | H749 34358662 0<br>08714729383550 | H749 34358668 0<br>08714729383611 |
|                           | JL4 90 cm      | -<br>-                            | -<br>-                            | -<br>-                            | -<br>-                            | H749 34358200 0<br>08714729842576 | H749 34358201 0<br>08714729842583 |
|                           | JL4.5          | H749 34356663 0<br>08714729351863 | -<br>-                            | H749 34357663 0<br>08714729380849 | H749 34357669 0<br>08714729380900 | H749 34358663 0<br>08714729383567 | H749 34358669 0<br>08714729383628 |
|                           | JL5            | H749 34356664 0<br>08714729351887 | -<br>-                            | H749 34357664 0<br>08714729380856 | H749 34357670 0<br>08714729380924 | H749 34358664 0<br>08714729383574 | H749 34358670 0<br>08714729383642 |
|                           | JL6            | H749 34356665 0<br>08714729351894 | -<br>-                            | H749 34357665 0<br>08714729380863 | H749 34357671 0<br>08714729380931 | H749 34358665 0<br>08714729383581 | -<br>-                            |
| Judkins Left<br>Short Tip | JL3.5 ST       | H749 34356673 0<br>08714729351917 | -<br>-                            | -<br>-                            | -<br>-                            | H749 34358673 0<br>08714729383666 | -<br>-                            |
|                           | JL4 ST         | H749 34356674 0<br>08714729351924 | -<br>-                            | H749 34357674 0<br>08714729380955 | H749 34357680 0<br>08714729381006 | H749 34358674 0<br>08714729383673 | H749 34358680 0<br>08714729383727 |

Continued next page >





## 2. Guide Catheters

### MACH 1™ Continued Femoral Approach Left Curve

|               | CURVE<br>STYLE | 6F = 2.0 MM (.070" = 1.78 MM ID)  |                                   | 7F = 2.33 MM (.081" = 2.05 MM ID) |                                   | 8F = 2.66 MM (.091" = 2.31 MM ID) |                                   |
|---------------|----------------|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|
|               |                | NO SIDE HOLES                     | SIDE HOLES                        | NO SIDE HOLES                     | SIDE HOLES                        | NO SIDE HOLES                     | SIDE HOLES                        |
| Amplatz™ Left | AL.75          | H749 3435622 0<br>08714729351894  | -<br>-                            | H749 3435722 0<br>08714729379348  | H749 34357124 0<br>08714729379041 | H749 3435822 0<br>08714729382027  | H749 34358124 0<br>08714729381716 |
|               | AL.75 90 cm    | -                                 | -                                 | -                                 | -                                 | H749 34358191 0<br>08714729842538 | H749 34358192 0<br>08714729842545 |
|               | AL1            | H749 3435619 0<br>08714729351153  | H749 3435668 0<br>08714729351955  | H749 3435719 0<br>08714729379317  | H749 3435768 0<br>08714729380993  | H749 3435819 0<br>08714729381990  | H749 3435868 0<br>08714729383710  |
|               | AL1 90 cm      | -                                 | -                                 | -                                 | -                                 | H749 34358376 0<br>08714729842613 | H749 34358398 0<br>08714729842637 |
|               | AL1.5          | H749 3435685 0<br>08714729352143  | -<br>-                            | H749 3435785 0<br>08714729381358  | -<br>-                            | H749 3435885 0<br>08714729384083  | -<br>-                            |
|               | AL2            | H749 3435620 0<br>08714729351160  | H749 3435669 0<br>08714729352020  | H749 3435720 0<br>08714729379324  | H749 3435769 0<br>08714729381099  | H749 3435820 0<br>08714729382003  | H749 3435869 0<br>08714729383819  |
|               | AL2 90 cm      | -                                 | -                                 | -                                 | -                                 | H749 34358377 0<br>08714729842620 | H749 34358399 0<br>08714729842644 |
|               | AL3            | H749 3435621 0<br>08714729378549  | -<br>-                            |                                   |                                   |                                   |                                   |
|               | AL.75 ST       | -<br>-                            | -<br>-                            | H749 34357108 0<br>08714729378938 | H749 34357119 0<br>08714729378990 | -<br>-                            | -<br>-                            |
|               | AL1 ST         | -<br>-                            | -<br>-                            | H749 34357103 0<br>08714729378891 | H749 34357114 0<br>08714729378952 | H749 34358103 0<br>08714729381556 | H749 34358114 0<br>08714729381617 |
|               | AL2 ST         | -<br>-                            | -<br>-                            | H749 34357105 0<br>08714729378914 | H749 34357116 0<br>08714729378976 | -<br>-                            | -<br>-                            |
| CLS™ Curve    | CLS3           | H749 34356573 0<br>08714729350637 | -<br>-                            | H749 34357573 0<br>08714729380382 | H749 34357577 0<br>08714729380467 | H749 34358573 0<br>08714729383130 | H749 34358577 0<br>08714729383215 |
|               | CLS3.5         | H749 34356574 0<br>08714729350644 | H749 34356578 0<br>08714729350675 | H749 34357574 0<br>08714729380405 | H749 34357578 0<br>08714729380481 | H749 34358574 0<br>08714729383154 | H749 34358578 0<br>08714729383239 |
|               | CLS4           | H749 34356575 0<br>08714729350651 | H749 34356579 0<br>08714729350682 | H749 34357575 0<br>08714729380429 | H749 34357579 0<br>08714729380504 | H749 34358575 0<br>08714729383178 | H749 34358579 0<br>08714729383253 |
|               | CLS4.5         | H749 34356576 0<br>08714729350668 | -<br>-                            | H749 34357576 0<br>08714729380443 | -<br>-                            | H749 34358576 0<br>08714729383192 | -<br>-                            |
| Kiesz™ Left   | KL3            | H749 34356611 0<br>08714729351696 | -<br>-                            | H749 34357611 0<br>08714729380597 | -<br>-                            | -<br>-                            | -<br>-                            |
|               | KL4            | H749 34356612 0<br>08714729351702 | -<br>-                            | H749 34357612 0<br>08714729380603 | -<br>-                            | -<br>-                            | -<br>-                            |
| Q-Curve™      | Q3.5           | H749 34356293 0<br>08714729351290 | H749 34356297 0<br>08714729351337 | H749 34357293 0<br>08714729350644 | H749 34357297 0<br>08714729379744 | H749 34358293 0<br>08714729382386 | -<br>-                            |
|               | Q4             | H749 34356294 0<br>08714729351306 | H749 34356298 0<br>08714729351344 | H749 34357294 0<br>08714729379713 | H749 34357298 0<br>08714729379751 | H749 34358294 0<br>08714729382393 | -<br>-                            |
|               | Q4.5           | H749 34356295 0<br>08714729351313 | -<br>-                            | H749 34357295 0<br>08714729379720 | H749 34357299 0<br>08714729379768 | H749 34358295 0<br>08714729382409 | -<br>-                            |
|               | Q5             | H749 34356296 0<br>08714729351320 | -<br>-                            | H749 34357296 0<br>08714729379737 | H749 34357300 0<br>08714729379782 | H749 34358296 0<br>08714729382416 | -<br>-                            |
| Voda Left™    | VL3            | H749 3435671 0<br>08714729352075  | -<br>-                            | H749 3435771 0<br>08714729381228  | H749 3435777 0<br>08714729381297  | H749 3435871 0<br>08714729383956  | H749 3435877 0<br>08714729384021  |
|               | VL3.5          | H749 3435672 0<br>08714729352082  | H749 3435678 0<br>08714729352129  | H749 3435772 0<br>08714729381242  | H749 3435778 0<br>08714729381303  | H749 3435872 0<br>08714729383970  | H749 3435878 0<br>08714729384038  |
|               | VL4            | H749 3435673 0<br>08714729352099  | H749 3435679 0<br>08714729352136  | H749 3435773 0<br>08714729381259  | H749 3435779 0<br>08714729381310  | H749 3435873 0<br>08714729383987  | H749 3435879 0<br>08714729384045  |
|               | VL4.5          | H749 3435674 0<br>08714729352105  | -<br>-                            | H749 3435774 0<br>08714729381266  | H749 3435780 0<br>08714729381327  | H749 3435874 0<br>08714729383994  | H749 3435880 0<br>08714729384052  |
|               | VL5            | H749 3435675 0<br>08714729352112  | -<br>-                            | H749 3435775 0<br>08714729381273  | H749 3435781 0<br>08714729381334  | H749 3435875 0<br>08714729384007  | H749 3435881 0<br>08714729384069  |
| Left Back-Up  | 90 cm/3,75     |                                   |                                   |                                   |                                   | H749 34358715 0<br>08714729842651 | H749 34358720 0<br>08714729842668 |
|               | 90 cm/3,5      |                                   |                                   |                                   |                                   | H749 34358198 0<br>08714729842552 |                                   |



## 2. Guide Catheters

### MACH 1™

#### Femoral Approach Right Curve

|                         | CURVE STYLE | 6F = 2.0 MM (.070" = 1.78 MM ID)  |                                   | 7F = 2.33 MM (.081" = 2.05 MM ID) |                                   | 8F = 2.66 MM (.091" = 2.31 MM ID) |                                   |
|-------------------------|-------------|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|
|                         |             | NO SIDE HOLES                     | SIDE HOLES                        | NO SIDE HOLES                     | SIDE HOLES                        | NO SIDE HOLES                     | SIDE HOLES                        |
| Femoral Right           | FR3         | H749 34356442 0<br>08714729351412 | -                                 | -                                 | -                                 | -                                 | -                                 |
|                         | FR3.5       | H749 3435601 0<br>08714729350408  | H749 3435647 0<br>08714729351450  | H749 3435701 0<br>08714729378730  | H749 3435747 0<br>08714729379966  | H749 3435801 0<br>08714729381457  | H749 3435847 0<br>08714729382768  |
|                         | FR4         | H749 3435602 0<br>08714729350415  | H749 3435648 0<br>08714729351498  | H749 3435702 0<br>08714729378808  | H749 3435748 0<br>08714729380030  | H749 3435802 0<br>08714729381464  | H749 3435848 0<br>08714729382775  |
|                         | FR4.5       | H749 34356514 0<br>08714729350774 | -                                 | -                                 | -                                 | -                                 | -                                 |
|                         | FR5         | H749 3435603 0<br>08714729350422  | -                                 | -                                 | -                                 | -                                 | -                                 |
|                         | FR4 90 cm   | -                                 | -                                 | -                                 | -                                 | H749 34358372 0<br>08714729382508 | H749 34358394 0<br>08714729382577 |
| Femoral Right Short Tip | FR4 ST      | H749 3435692 0<br>08714729352167  | -                                 | -                                 | -                                 | -                                 | -                                 |
| Femoral Curve Right     | FCR3.5      | H749 34356279 0<br>08714729351252 | -                                 | H749 34357279 0<br>08714729379652 | H749 34357276 0<br>08714729379621 | -                                 | -                                 |
|                         | FCR4        | -                                 | -                                 | H749 34357280 0<br>08714729379676 | H749 34357277 0<br>08714729379638 | H749 34358280 0<br>08714729382355 | -                                 |
| Judkins Right           | JR3         | H749 34356684 0<br>08714729351962 | -                                 | H749 34357684 0<br>08714729381037 | H749 34357690 0<br>08714729381105 | H749 34358684 0<br>08714729383758 | -                                 |
|                         | JR3.5       | H749 34356685 0<br>08714729351979 | H749 34356691 0<br>08714729352037 | H749 34357685 0<br>08714729381044 | H749 34357691 0<br>08714729381112 | H749 34358685 0<br>08714729383765 | H749 34358691 0<br>08714729383833 |
|                         | JR4         | H749 34356686 0<br>08714729351986 | H749 34356692 0<br>08714729352044 | H749 34357686 0<br>08714729381051 | H749 34357692 0<br>08714729381129 | H749 34358686 0<br>08714729383772 | H749 34358692 0<br>08714729383840 |
|                         | JR4.5       | H749 34356687 0<br>08714729351993 | -                                 | H749 34357687 0<br>08714729381068 | H749 34357693 0<br>08714729381136 | H749 34358687 0<br>08714729383789 | H749 34358693 0<br>08714729383857 |
|                         | JR5         | H749 34356688 0<br>08714729352006 | -                                 | H749 34357688 0<br>08714729381075 | H749 34357694 0<br>08714729381143 | H749 34358688 0<br>08714729383796 | H749 34358694 0<br>08714729383864 |
|                         | JR6         | H749 34356689 0<br>08714729352013 | -                                 | H749 34357689 0<br>08714729381082 | -                                 | -                                 | -                                 |
|                         | JR3.5 90 cm | -                                 | -                                 | -                                 | -                                 | H749 34358202 0<br>08714729842590 | -                                 |
|                         | JR4 90 cm   | -                                 | -                                 | -                                 | -                                 | H749 34358708 0<br>08714729383932 | -                                 |
| Judkins Right Short Tip | JR3.5 ST    | H749 34356697 0<br>08714729352051 | -                                 | -                                 | -                                 | -                                 | -                                 |
|                         | JR4 ST      | H749 34356698 0<br>08714729352068 | -                                 | H749 34357698 0<br>08714729381174 | H749 34357704 0<br>08714729381204 | -                                 | -                                 |
| allRight™               | ART3        | H749 34356496 0<br>08714729351504 | H749 34356497 0<br>08714729351511 | -                                 | H749 34357497 0<br>08714729380061 | H749 34358496 0<br>08714729382799 | H749 34358497 0<br>08714729382805 |
|                         | ART3.5      | H749 34356498 0<br>08714729351528 | H749 34356499 0<br>08714729351535 | H749 34357498 0<br>08714729380078 | H749 34357499 0<br>08714729380085 | H749 34358498 0<br>08714729382812 | H749 34358499 0<br>08714729382829 |
|                         | ART4        | H749 34356500 0<br>08714729350743 | H749 34356501 0<br>08714729350750 | H749 34357500 0<br>08714729380108 | H749 34357501 0<br>08714729380115 | H749 34358500 0<br>08714729382843 | H749 34358501 0<br>08714729382850 |
|                         | ART4.5      | H749 34356502 0<br>08714729350767 | -                                 | -                                 | -                                 | -                                 | -                                 |
| Amplatz™ Right          | AR1         | H749 3435617 0<br>08714729351139  | H749 3435666 0<br>08714729351832  | H749 3435717 0<br>08714729379218  | H749 3435766 0<br>08714729380801  | H749 3435817 0<br>08714729381891  | -                                 |
|                         | AR2         | H749 3435618 0<br>08714729351146  | H749 3435667 0<br>08714729351900  | H749 3435718 0<br>08714729379263  | H749 3435767 0<br>08714729380917  | H749 3435818 0<br>08714729381945  | H749 3435867 0<br>08714729383635  |

Continued next page >



## 2. Guide Catheters

### MACH 1™ Continued Femoral Approach Right Curve

|                                    | CURVE<br>STYLE | 6F = 2.0 MM (.070" = 1.78 MM ID)  |                                   | 7F = 2.33 MM (.081" = 2.05 MM ID) |                                   | 8F = 2.66 MM (.091" = 2.31 MM ID) |                                   |
|------------------------------------|----------------|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|
|                                    |                | NO SIDE HOLES                     | SIDE HOLES                        | NO SIDE HOLES                     | SIDE HOLES                        | NO SIDE HOLES                     | SIDE HOLES                        |
| Kiesz™ Right                       | KR3H           | H749 34356648 0<br>08714729351771 | H749 34356651 0<br>08714729378686 | H749 34357648 0<br>08714729380689 | -<br>-                            | H749 34358648 0<br>08714729383406 | -<br>-                            |
|                                    | KR4H           | H749 34356649 0<br>08714729351788 | H749 34356652 0<br>08714729378693 | H749 34357649 0<br>08714729380696 | -<br>-                            | H749 34358649 0<br>08714729383413 | -<br>-                            |
|                                    | KR3S           | H749 34356654 0<br>08714729351801 | -                                 | H749 34357654 0<br>08714729380740 | -<br>-                            | -                                 | -                                 |
|                                    | KR4S           | H749 34356655 0<br>08714729351818 | H749 34356658 0<br>08714729378723 | H749 34357655 0<br>08714729380757 | -<br>-                            | H749 34358655 0<br>08714729383475 | -<br>-                            |
| Right Coronary<br>Shepherd's Crook | RC3.5 SC       | H749 3435646 0<br>08714729351443  | -                                 | H749 3435746 0<br>08714729379959  | -<br>-                            | -                                 | -                                 |
|                                    | RC4 SC         | H749 34356448 0<br>08714729351436 | -                                 | H749 34357448 0<br>08714729379935 | H749 34357449 0<br>08714729379942 | H749 34358448 0<br>08714729382737 | H749 34358449 0<br>08714729382744 |
| Voda Right®                        | VR1            | H749 34356230 0<br>08714729351184 | -                                 | H749 34357230 0<br>08714729379362 | -<br>-                            | -                                 | -                                 |
|                                    | VR2            | H749 34356231 0<br>08714729351191 | -                                 | H749 34357231 0<br>08714729379379 | -<br>-                            | H749 34358231 0<br>08714729382058 | -<br>-                            |
| Multipurpose                       | MP1            | H749 3435638 0<br>08714729351375  | H749 34356138 0<br>08714729351115 | H749 3435738 0<br>08714729379812  | H749 34357138 0<br>08714729379171 | H749 3435838 0<br>08714729382515  | H749 34358138 0<br>08714729381846 |
|                                    | MP1 90 cm      | -                                 | -                                 | -                                 | -                                 | H749 34358380 0<br>08714729382522 | -                                 |
|                                    | MP2            | H749 3435639 0<br>08714729351382  | -                                 | H749 3435739 0<br>08714729379829  | H749 34357139 0<br>08714729379188 | H749 3435839 0<br>08714729382560  | H749 34358139 0<br>08714729381853 |
| Hockey Stick                       | HS             | H749 3435641 0<br>08714729351405  | H749 34356141 0<br>08714729351122 | H749 3435741 0<br>08714729379843  | H749 34357141 0<br>08714729379201 | H749 3435841 0<br>08714729382638  | H749 34358141 0<br>08714729381877 |

### MACH 1™ Mammary & Bypass Curves

|                                   | CURVE<br>STYLE | 6F = 2.0 MM (.070" = 1.78 MM ID)  |                                   | 7F = 2.33 MM (.081" = 2.05 MM ID) |                                   | 8F = 2.66 MM (.091" = 2.31 MM ID) |                                   |
|-----------------------------------|----------------|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|
|                                   |                | NO SIDE HOLES                     | SIDE HOLES                        | NO SIDE HOLES                     | SIDE HOLES                        | NO SIDE HOLES                     | SIDE HOLES                        |
| Internal Mammary                  | IM 85 cm       | H749 34356122 0<br>08714729351061 | -<br>-                            | H749 34357122 0<br>08714729379027 | H749 34357128 0<br>08714729379089 | -<br>-                            | -<br>-                            |
|                                   | IM 90 cm       | H749 34356446 0<br>08714729351429 | -                                 |                                   | H749 34357447 0<br>08714729379928 | H749 34358446 0<br>08714729382713 | -                                 |
|                                   | IM 100 cm      | H749 3435626 0<br>08714729351221  | H749 34356129 0<br>08714729351108 | H749 3435726 0<br>08714729379515  | H749 34357129 0<br>08714729379096 | H749 3435826 0<br>08714729382195  | H749 34358129 0<br>08714729381761 |
| Left and Right<br>Coronary Bypass | LCB            | H749 3435625 0<br>08714729351214  | H749 34356127 0<br>08714729351092 | H749 3435725 0<br>08714729379478  | H749 34357127 0<br>08714729379072 | H749 3435825 0<br>08714729382157  | H749 34358127 0<br>08714729381747 |
|                                   | LCB 90 cm      | -                                 | -                                 | -                                 | -                                 | H749 34358383 0<br>08714729382546 | H749 34358405 0<br>08714729382614 |
|                                   | RCB            | H749 3435624 0<br>08714729351207  | H749 34356126 0<br>08714729351085 | H749 3435724 0<br>08714729379409  | -                                 | H749 3435824 0<br>08714729382089  | H749 34358126 0<br>08714729381730 |



2. Guide Catheters

MACH 1™  
Brachial Approach

|          | CURVE<br>STYLE | 6F = 2.0 MM (.070" = 1.78 MM ID) |            | 7F = 2.33 MM (.081" = 2.05 MM ID) |            |
|----------|----------------|----------------------------------|------------|-----------------------------------|------------|
|          |                | NO SIDE HOLES                    | SIDE HOLES | NO SIDE HOLES                     | SIDE HOLES |
| Brachial | BR1 80 cm      | H749 3435640 0<br>08714729351399 | -<br>-     | -<br>-                            | -<br>-     |

MACH 1™  
Transradial Approach

|                  | CURVE<br>STYLE | 6F = 2.0 MM (.070" = 1.78 MM ID)  |                                   | 7F = 2.33 MM (.081" = 2.05 MM ID) |            |
|------------------|----------------|-----------------------------------|-----------------------------------|-----------------------------------|------------|
|                  |                | NO SIDE HOLES                     | SIDE HOLES                        | NO SIDE HOLES                     | SIDE HOLES |
| Kimny™           | Kimny™         | H749 34356583 0<br>08714729349600 | H749 34356591 0<br>08714729350712 | H749 34357583 0<br>08714729380542 | -<br>-     |
|                  | Kimny™ Mini    | H749 34356584 0<br>08714729350705 | H749 34356592 0<br>08714729350729 | -<br>-                            | -<br>-     |
| Muta™ Left       | ML3            | H749 34356536 0<br>08714729350514 | -<br>-                            | -<br>-                            | -<br>-     |
|                  | ML3.5          | H749 34356537 0<br>08714729350521 | -<br>-                            | -<br>-                            | -<br>-     |
|                  | ML4            | H749 34356538 0<br>08714729350538 | -<br>-                            | -<br>-                            | -<br>-     |
| Radial           | RADIAL         | H749 34356472 0<br>08714729351474 | H749 34356473 0<br>08714729351481 | -<br>-                            | -<br>-     |
| Radial Short Tip | RADIAL ST      | H749 34356470 0<br>08714729351467 | -<br>-                            | -<br>-                            | -<br>-     |

|                               | CURVE<br>STYLE | 6F = 2.0 MM (.066" = 1.68 MM ID)  |
|-------------------------------|----------------|-----------------------------------|
|                               |                | NO SIDE HOLES                     |
| Radial Tig Curve              | TIG3           | H749 34356741 0<br>08714729886457 |
|                               | TIG3.5         | H749 34356742 0<br>08714729886464 |
|                               | TIG4           | H749 34356743 0<br>08714729886471 |
|                               | TIG4.5         | H749 34356744 0<br>08714729886488 |
|                               | TIG5           | H749 34356745 0<br>08714729886495 |
| Radial Back-Up<br>Left Curve  | RBL3.5         | H749 34356735 0<br>08714729886396 |
|                               | RBL4           | H749 34356736 0<br>08714729886402 |
|                               | RBL4.5         | H749 34356737 0<br>08714729886419 |
| Radial Back-Up<br>Right Curve | RBR3           | H749 34356738 0<br>08714729886426 |
|                               | RBR3.5         | H749 34356739 0<br>08714729886433 |
|                               | RBR4           | H749 34356740 0<br>08714729886440 |



## 2. Guide Catheters

### RUNWAY™

#### Femoral Approach Left Curve

|                           | CURVE<br>STYLE   | 6F = 2.0 MM (.070" = 1.78 MM ID)<br>NO SIDE HOLES | SIDE HOLES                        |
|---------------------------|------------------|---------------------------------------------------|-----------------------------------|
| Femoral Left              | FL3              | H749 3896905 0<br>08714729477662                  | -<br>-                            |
|                           | FL3.5            | H749 3896906 0<br>08714729477679                  | H749 3896952 0<br>08714729478256  |
|                           | FL4              | H749 3896907 0<br>08714729477686                  | H749 3896953 0<br>08714729478263  |
|                           | FL4.5            | H749 3896908 0<br>08714729477693                  | -<br>-                            |
|                           | FL5              | H749 3896909 0<br>08714729477709                  | -<br>-                            |
|                           | FL6              | H749 3896910 0<br>08714729477716                  | -<br>-                            |
| Femoral Curve<br>Left     | FCL3.5           | H749 38969272 0<br>08714729477907                 | -<br>-                            |
|                           | FCL4             | H749 38969273 0<br>08714729477914                 | -<br>-                            |
| Judkins Left              | JL3              | H749 38969660 0<br>08714729478706                 | -<br>-                            |
|                           | JL3.5            | H749 38969661 0<br>08714729478713                 | H749 38969667 0<br>08714729478768 |
|                           | JL4              | H749 38969662 0<br>08714729478720                 | H749 38969668 0<br>08714729478775 |
|                           | JL4.5            | H749 38969663 0<br>08714729478737                 | -<br>-                            |
|                           | JL5              | H749 38969664 0<br>08714729478744                 | -<br>-                            |
|                           | JL6              | H749 38969665 0<br>08714729478751                 | -<br>-                            |
| Judkins Left<br>Short Tip | JL3.5 ST         | H749 38969673 0<br>08714729478799                 | -<br>-                            |
|                           | JL4 ST           | H749 38969674 0<br>08714729478805                 | -<br>-                            |
|                           | JL4.5 ST         | H749 38969675 0<br>08714729478812                 | -<br>-                            |
| Amplatz™ Left             | AL .75           | H749 3896922 0<br>08714729477846                  | -<br>-                            |
|                           | AL .75/<br>90 cm | H749 38969191 0<br>08714729964995                 | -<br>-                            |
|                           | AL1              | H749 3896919 0<br>08714729477815                  | H749 3896968 0<br>08714729478836  |
|                           | AL1/<br>90 cm    | H749 38969376 0<br>08714729965060                 | -<br>-                            |
|                           | AL1.5            | H749 3896985 0<br>08714729479031                  | -<br>-                            |
|                           | AL2              | H749 3896920 0<br>08714729477822                  | H749 3896969 0<br>08714729478904  |
|                           | AL3              | H749 3896921 0<br>08714729477839                  | -<br>-                            |

|                           | CURVE<br>STYLE | 6F = 2.0 MM (.070" = 1.78 MM ID)<br>NO SIDE HOLES | SIDE HOLES                        |
|---------------------------|----------------|---------------------------------------------------|-----------------------------------|
| Left Back-up              | Left BU3       | H749 38969713 0<br>08714729813118                 | H749 38969719 0<br>08714729813125 |
|                           | Left BU3.5     | H749 38969714 0<br>08714729813132                 | H749 38969720 0<br>08714729790686 |
|                           | Left BU3.75    | H749 38969715 0<br>08714729813149                 | H749 38969721 0<br>08714729820703 |
|                           | Left BU4       | H749 38969716 0<br>08714729813156                 | H749 38969722 0<br>08714729790679 |
|                           | Left BU4.5     | H749 38969717 0<br>08714729813163                 | H749 38969723 0<br>08714729820710 |
|                           | Left BU5       | H749 38969718 0<br>08714729813170                 | H749 38969724 0<br>08714729820727 |
| Femoral Left<br>Short Tip | FL3.5 ST       | H749 3896928 0<br>08714729477938                  | -<br>-                            |
|                           | FL4 ST         | H749 3896929 0<br>08714729477952                  | -<br>-                            |
|                           | FL4.5 ST       | H749 3896930 0<br>08714729478027                  | -<br>-                            |
|                           | FL5 ST         | H749 3896931 0<br>08714729478034                  | -<br>-                            |
|                           | CLS3           | H749 38969573 0<br>08714729478379                 | -<br>-                            |
| CLS™ Curve Left           | CLS3.5         | H749 38969574 0<br>08714729478386                 | H749 38969578 0<br>08714729478416 |
|                           | CLS3.75        | H749 38969577 0<br>08714729642534                 | H749 38969580 0<br>08714729642541 |
|                           | CLS4           | H749 38969575 0<br>08714729478393                 | H749 38969579 0<br>08714729478423 |
|                           | CLS4.5         | H749 38969576 0<br>08714729478409                 | -<br>-                            |
|                           | KL3H           | H749 38969611 0<br>08714729478492                 | H749 38969614 0<br>08714729478522 |
| Kiesz™ Left               | KL4H           | H749 38969612 0<br>08714729478508                 | H749 38969615 0<br>08714729478539 |
|                           | Q3.5           | H749 38969293 0<br>08714729477969                 | H749 38969297 0<br>08714729478003 |
| Q-Curve™                  | Q4             | H749 38969294 0<br>08714729477976                 | H749 38969298 0<br>08714729478010 |
|                           | Q4.5           | H749 38969295 0<br>08714729477983                 | -<br>-                            |
|                           | Q5             | H749 38969296 0<br>08714729477990                 | -<br>-                            |
| Voda Left™                | VL3            | H749 3896971 0<br>08714729478966                  | -<br>-                            |
|                           | VL3.5          | H749 3896972 0<br>08714729478973                  | H749 3896978 0<br>08714729479017  |
|                           | VL4            | H749 3896973 0<br>08714729478980                  | H749 3896979 0<br>08714729479024  |
|                           | VL4.5          | H749 3896974 0<br>08714729478997                  | -<br>-                            |
|                           | VL5            | H749 3896975 0<br>08714729479000                  | -<br>-                            |
|                           |                |                                                   |                                   |





## 2. Guide Catheters

### RUNWAY™

#### Femoral Approach Right Curve

|                                 | CURVE STYLE | 6F = 2.0 MM (.070" = 1.78 MM ID)<br>NO SIDE HOLES | SIDE HOLES                        |
|---------------------------------|-------------|---------------------------------------------------|-----------------------------------|
| Femoral Curve Right             | FR3         | H749 38969442 0<br>08714729478089                 | -<br>-                            |
|                                 | FR3.5       | H749 3896901 0<br>08714729477624                  | H749 3896947 0<br>08714729478126  |
|                                 | FR4         | H749 3896902 0<br>08714729477631                  | H749 3896948 0<br>08714729478164  |
|                                 | FR4.5       | H749 38969514 0<br>08714729478249                 | -<br>-                            |
|                                 | FR5         | H749 3896903 0<br>08714729477648                  | -<br>-                            |
|                                 | FR6         | H749 3896904 0<br>08714729477655                  | -<br>-                            |
| Femoral Right Short Tip         | FR4 ST      | H749 3896992 0<br>08714729479055                  | -<br>-                            |
| Femoral Curve Right             | FCR3.5      | H749 38969279 0<br>08714729477921                 | -<br>-                            |
|                                 | FCR4        | H749 38969280 0<br>08714729477945                 | -<br>-                            |
| Amplatz™ Right                  | AR1         | H749 3896917 0<br>08714729477792                  | H749 3896966 0<br>08714729478690  |
|                                 | AR2         | H749 3896918 0<br>08714729477808                  | H749 3896967 0<br>08714729478782  |
| Kiesz™ Right                    | KR3H        | H749 38969648 0<br>08714729478577                 | H749 38969651 0<br>08714729478607 |
|                                 | KR4H        | H749 38969649 0<br>08714729478584                 | H749 38969652 0<br>08714729478614 |
|                                 | KR5H        | H749 38969650 0<br>08714729478591                 | H749 38969653 0<br>08714729478621 |
|                                 | KR3S        | H749 38969654 0<br>08714729478638                 | H749 38969657 0<br>08714729478669 |
|                                 | KR4S        | H749 38969655 0<br>08714729478645                 | H749 38969658 0<br>08714729478676 |
|                                 | KR5S        | -<br>-                                            | H749 38969659 0<br>08714729478683 |
| Right Coronary Shepherd's Crook | RC3.5 SC    | H749 3896946 0<br>08714729478119                  | -<br>-                            |
|                                 | RC4 SC      | H749 38969448 0<br>08714729478102                 | -<br>-                            |
| Multipurpose                    | MP1         | H749 3896938 0<br>08714729478041                  | H749 38969138 0<br>08714729477778 |
|                                 | MP1/125 cm  | H749 38969953 0<br>08714729965374                 | -<br>-                            |
|                                 | MP2         | H749 3896939 0<br>08714729478058                  | -<br>-                            |
| Hockey Stick                    | HS          | H749 3896941 0<br>08714729478072                  | H749 38969141 0<br>08714729477785 |
|                                 | CURVE STYLE | 6F = 2.0 MM (.070" = 1.78 MM ID)<br>NO SIDE HOLES | SIDE HOLES                        |
| Judkins Right                   | JR3         | H749 38969684 0<br>08714729478843                 | -<br>-                            |
|                                 | JR3.5       | H749 38969685 0<br>08714729478850                 | H749 38969691 0<br>08714729478911 |
|                                 | JR4         | H749 38969686 0<br>08714729478867                 | H749 38969692 0<br>08714729478928 |
|                                 | JR4/125 cm  | H749 38969956 0<br>08714729965381                 | -<br>-                            |
|                                 | JR4.5       | H749 38969687 0<br>08714729478874                 | -<br>-                            |
|                                 | JR5         | H749 38969688 0<br>08714729478881                 | -<br>-                            |
| Judkins Right Short Tip         | JR6         | H749 38969689 0<br>08714729478898                 | -<br>-                            |
|                                 | JR4 ST      | H749 38969698 0<br>08714729478942                 | -<br>-                            |
| allRight™                       | ART3        | H749 38969496 0<br>08714729478171                 | H749 38969497 0<br>08714729478188 |
|                                 | ART3.5      | H749 38969498 0<br>08714729478195                 | H749 38969499 0<br>08714729478201 |
|                                 | ART4        | H749 38969500 0<br>08714729478218                 | H749 38969501 0<br>08714729478225 |
|                                 | ART4.5      | H749 38969502 0<br>08714729478232                 | -<br>-                            |
| Internal Mammary                | IM          | H749 3896926 0<br>08714729477891                  | H749 38969129 0<br>08714729477761 |
|                                 | IM 85 cm    | H749 38969122 0<br>08714729477723                 | -<br>-                            |
|                                 | IM 90 cm    | H749 38969446 0<br>08714729478096                 | -<br>-                            |
|                                 | IM 110 cm   | H749 38969123 0<br>08714729477730                 | -<br>-                            |
| Internal Mammary Cummings       | IMC 90 cm   | H749 38969638 0<br>08714729478560                 | -<br>-                            |
| Left and Right Coronary         | LCB         | H749 3896925 0<br>08714729477884                  | H749 38969127 0<br>08714729477754 |
|                                 | RCB         | H749 3896924 0<br>08714729477877                  | H749 38969126 0<br>08714729477747 |
| William Right Posterior         | WRP         | H749 38969609 0<br>08714729478478                 | H749 38969610 0<br>08714729478485 |



2. Guide Catheters

| RUNWAY™<br>Transradial Approach |                  |                                   |                                  |                                   |        |        |                                   |                                  |                                   |
|---------------------------------|------------------|-----------------------------------|----------------------------------|-----------------------------------|--------|--------|-----------------------------------|----------------------------------|-----------------------------------|
|                                 |                  | CURVE<br>STYLE                    | 6F = 2.0 MM (.070" = 1.78 MM ID) |                                   |        |        | CURVE<br>STYLE                    | 6F = 2.0 MM (.070" = 1.78 MM ID) |                                   |
|                                 |                  |                                   | NO SIDE HOLES                    | SIDE HOLES                        |        |        |                                   | NO SIDE HOLES                    | SIDE HOLES                        |
| Kimny™                          | Kimny™           | H749 38969583 0<br>08714729478430 |                                  | H749 38969591 0<br>08714729478454 | Radial | Radial | H749 38969472 0<br>08714729478140 |                                  | H749 38969473 0<br>08714729478157 |
|                                 | Kimny™ mini      | H749 38969584 0<br>08714729478447 |                                  | -                                 |        |        |                                   |                                  |                                   |
| Mann Internal<br>Mammary        | Mann IM<br>90 cm | H749 38969632 0<br>08714729478553 |                                  | -                                 |        |        |                                   |                                  |                                   |
|                                 | ML3              | H749 38969536 0<br>08714729478270 |                                  | -                                 |        |        |                                   |                                  |                                   |
| Muta™ Left                      | ML3.5            | H749 38969537 0<br>08714729478287 |                                  | -                                 |        |        |                                   |                                  |                                   |
|                                 | ML4              | H749 38969538 0<br>08714729478294 |                                  | -                                 |        |        |                                   |                                  |                                   |
|                                 | ML4.5            | H749 38969539 0<br>08714729478300 |                                  | -                                 |        |        |                                   |                                  |                                   |

| RUNWAY™<br>Brachial Approach |           |                                   |                                  |            |  |  |  |  |  |
|------------------------------|-----------|-----------------------------------|----------------------------------|------------|--|--|--|--|--|
|                              |           | CURVE<br>STYLE                    | 6F = 2.0 MM (.070" = 1.78 MM ID) |            |  |  |  |  |  |
|                              |           |                                   | NO SIDE HOLES                    | SIDE HOLES |  |  |  |  |  |
| Brachial                     | BR1 80 cm | H749 3896940 0<br>408714729478065 |                                  | -          |  |  |  |  |  |

Balloon  
Catheters

Micro-  
dilatation  
Catheters

Guide  
Catheters

Guide  
Catheter  
extension

Angiographic  
Catheters

Guidewires

Embolic  
Protection

Accessories

Cutting  
Balloon™

Atherectomy  
Devices

Chronic Total  
Occlusions

Trapper™  
Exchange  
Device

ANGIOJET™  
Thrombec-  
tomy  
System



NEW CARDIO

2

## COMPLEX PCI

Balloon Catheters

Microdilatation Catheters

Guide Catheters

Guide Catheter Extension

Angiographic Catheters

Guidewires

Embolic Protection

Accessories

Cutting Balloon™

Atherectomy Devices

Chronic Total Occlusions

TRAPPER™ Exchange Device

ANGIOJET™ Thrombectomy System



Balloon  
Catheters

Micro-  
dilatation  
Catheters

Guide  
Catheters

Guide  
Catheter  
extension

Angiographic  
Catheters

Guidewires

Embolic  
Protection

Accessories

Cutting  
Balloon™

Atherectomy  
Devices

Chronic Total  
Occlusions

Trapper™  
Exchange  
Device

ANGIOJET™  
Thrombec-  
tomy  
System

2. Guide Catheter Extension

GUIDEZILLA™ II  
Guide Catheter Extension

| CATALOG NUMBER / ORDER NUMBER (GTIN) | DESCRIPTION               | COMPATIBLE GUIDE CATHETER | INNER DIAMETER   | OUTER DIAMETER   |
|--------------------------------------|---------------------------|---------------------------|------------------|------------------|
| H749 393351506 0 / 08714729939450    | GUIDEZILLA™ II, 6F*       | 6F / 0.070" (1,78 mm)     | 0.057" (1,45 mm) | 0.067" (1.71 mm) |
| H749 393351507 0 / 08714729939474    | GUIDEZILLA™ II, 7F*       | 7F / 0.078" (1,98 mm)     | 0.063" (1,60 mm) | 0.073" (1.86 mm) |
| H749 393351508 0 / 08714729939481    | GUIDEZILLA™ II, 8F*       | 8F / 0.088" (2,24 mm)     | 0.072" (1,83 mm) | 0.083" (2.11 mm) |
| H749 3933515061 0 / 08714729939467   | GUIDEZILLA™ II, 6F*, Long | 6F / 0.070" (1,78 mm)     | 0.057" (1,45 mm) | 0.067" (1.71 mm) |

\* 6F = 2 mm, 7F = 2.33mm, 8F = 2.67 mm

Balloon  
Catheters

Micro-  
dilatation  
Catheters

Guide  
Catheters

Guide  
Catheter  
extension

Angiographic  
Catheters

Guidewires

Embolic  
Protection

Accessories

Cutting  
Balloon™

Atherectomy  
Devices

Chronic Total  
Occlusions

Trapper™  
Exchange  
Device

ANGIOJET™  
Thrombec-  
tomy  
System



NEW CARDIO

2

## COMPLEX PCI

Balloon Catheters

Microdilatation Catheters

Guide Catheters

Guide Catheter Extension

Angiographic Catheters

Guidewires

Embolic Protection

Accessories

Cutting Balloon™

Atherectomy Devices

Chronic Total Occlusions

TRAPPER™ Exchange Device

ANGIOJET™ Thrombectomy System



Balloon  
Catheters

Micro-  
dilatation  
Catheters

Guide  
Catheters

Guide  
Catheter  
extension

Angiographic  
Catheters

Guidewires

Embolic  
Protection

Accessories

Cutting  
Balloon™

Atherectomy  
Devices

Chronic Total  
Occlusions

Trapper™  
Exchange  
Device

ANGIOJET™  
Thrombec-  
tomy  
System



## 2. Angiographic Catheters

### IMPULSE™ Femoral Selective

|                            | CURVE<br>STYLE | LENGTH | QUANTITY | 5F = 1.66 MM<br>LEFT .047" = 1.2 MM ID<br>RIGHT .045" = 1.14 MM ID | 6F = 2.00 MM<br>(.056" = 1.42 MM ID) |
|----------------------------|----------------|--------|----------|--------------------------------------------------------------------|--------------------------------------|
| Femoral Left               | FL 3.5         | 100 cm | 1        | H749 1639121 1<br>08714729214007                                   | H749 1659921 0<br>08714729304692     |
|                            | FL 3.5         | 100 cm | 5        | H749 1639121 2<br>08714729188889                                   | H749 1659921 2<br>08714729200222     |
|                            | FL 4           | -      | 1        | H749 1639122 1<br>08714729344834                                   | -                                    |
|                            | FL 4           | 100 cm | 5        | H749 1639122 2<br>08714729188896                                   | H749 1659922 2<br>08714729200246     |
|                            | FL 4           | 125 cm | 5        | H749 16391220 2<br>08714729153047                                  | H749 16599220 2<br>08714729304739    |
|                            | FL 4.5         | 100 cm | 5        | -                                                                  | H749 1659923 2<br>08714729200260     |
|                            | FL 5           | -      | 1        | H749 1639124 1<br>08714729343431                                   | H749 1659924 0<br>08714729304784     |
|                            | FL 5           | 100 cm | 5        | H749 1639124 2<br>08714729153092                                   | H749 1659924 2<br>08714729200284     |
|                            | FL 5           | 125 cm | 5        | H749 16391235 2<br>08714729188919                                  | H749 16599235 2<br>08714729200277    |
|                            | FL 6           | 100 cm | 5        | H749 1639125 2<br>08714729153115                                   | H749 1659925 2<br>08714729200291     |
| Femoral Right              | FR 3.5         | 100 cm | 1        | H749 1639101 1<br>08714729213161                                   | H749 1659901 0<br>08714729214816     |
|                            | FR 3.5         | 100 cm | 5        | H749 1639101 2<br>08714729188803                                   | H749 1659901 2<br>08714729200093     |
|                            | FR 4           | -      | 1        | H749 1639102 1<br>08714729213178                                   | -                                    |
|                            | FR 4           | 100 cm | 1        | -                                                                  | H749 1659902 0<br>08714729345404     |
|                            | FR 4           | 100 cm | 5        | H749 1639102 2<br>08714729188810                                   | H749 1659902 2<br>08714729200109     |
|                            | FR 4           | 125 cm | 1        | H749 16391225 1<br>08714729344841                                  | -                                    |
|                            | FR 4           | 125 cm | 5        | H749 16391225 2<br>08714729188902                                  | H749 16599225 2<br>08714729200253    |
|                            | FR 4 MOD       | 100 cm | 5        | H749 1639105 2<br>08714729166627                                   | H749 1659905 2<br>08714729299431     |
|                            | FR 5           | 100 cm | 1        | H749 1639103 1<br>08714729213185                                   | -                                    |
|                            | FR 5           | -      | 1        | -                                                                  | H749 1659903 0<br>08714729214830     |
|                            | FR 5           | 100 cm | 5        | H749 1639103 2<br>08714729152873                                   | H749 1659903 2<br>08714729200116     |
|                            | FR 5           | 125 cm | 5        | H749 16391230 2<br>08714729153078                                  | H749 16599230 2<br>08714729215769    |
|                            | FR 6           | 100 cm | 5        | H749 1639104 2<br>08714729188827                                   | H749 1659904 2<br>08714729200123     |
|                            | FR 6           | 100 cm | 5        | H749 1639104 2<br>08714729188827                                   | H749 1659904 2<br>08714729200123     |
| Amplatz™ Left<br>and Right | AL 1           | 100 cm | 1        | H749 1639196 1<br>08714729214656                                   | H749 1659996 0<br>08714729305279     |
|                            | AL 1           | 100 cm | 5        | H749 1639196 2<br>08714729153214                                   | H749 1659996 2<br>08714729200383     |
|                            | AL 2           | 100 cm | 1        | H749 1639198 1<br>08714729214663                                   | -                                    |
|                            | AL 2           | 100 cm | 5        | H749 1639198 2<br>08714729190011                                   | H749 1659998 2<br>08714729200390     |
|                            | AL 3           | 100 cm | 5        | H749 16391100 2<br>08714729152897                                  | H749 16599100 2<br>08714729200130    |
|                            | AR 1           | 100 cm | 1        | H749 1639184 1<br>08714729345329                                   | -                                    |
|                            | AR 1           | 100 cm | 5        | H749 1639184 2<br>08714729153184                                   | H749 1659984 2<br>08714729200369     |
|                            | AR 2           | 100 cm | 5        | H749 1639186 2<br>08714729153191                                   | H749 1659986 2<br>08714729200376     |
|                            | AR MOD         | 100 cm | 5        | H749 1639189 2<br>08714729153207                                   | H749 1659989 2<br>08714729216278     |

Continued next page >




## 2. Angiographic Catheters




### IMPULSE™ Continued Femoral Selective

|                  | CURVE<br>STYLE         | LENGTH | QUANTITY | 5F = 1.66 MM<br>LEFT .047" = 1.2 MM ID<br>RIGHT .045" = 1.14 MM ID | 6F = 2.00 MM<br>(.056" = 1.42 MM ID) |
|------------------|------------------------|--------|----------|--------------------------------------------------------------------|--------------------------------------|
| Multipurpose     | MP A-1                 | 100 cm | 5        | H749 16391117 2<br>08714729152903                                  | H749 16599117 2<br>08714729200147    |
|                  | MP A-2                 | -      | 1        | H749 16391121 1<br>08714729342687                                  | -                                    |
|                  | MP A-2 with side holes | 100 cm | 5        | H749 16391121 2<br>08714729188834                                  | H749 16599121 2<br>08714729200154    |
|                  | MP A-2                 | 125 cm | 5        | H749 16391122 2<br>08714729152972                                  | H749 16599122 2<br>08714729200161    |
|                  | MP B-1                 | 100 cm | 5        | H749 16391142 2<br>08714729188841                                  | H749 16599142 2<br>08714729200178    |
|                  | MP B-2 with side holes | 100 cm | 5        | H749 16391146 2<br>08714729152941                                  | H749 16599146 2<br>08714729200185    |
| Coronary Bypass  | LCB                    | 100 cm | 1        | H749 16391195 1<br>08714729213888                                  | H749 16599195 0<br>08714729304593    |
|                  | LCB                    | 100 cm | 5        | H749 16391195 2<br>08714729152996                                  | H749 16599195 2<br>08714729200208    |
|                  | RCB                    | 100 cm | 1        | -                                                                  | H749 16599190 0<br>08714729304586    |
|                  | RCB                    | 100 cm | 5        | H749 16391190 2<br>08714729200086                                  | H749 16599190 2<br>08714729200192    |
| Internal Mammary | IM                     | -      | 1        | H749 16391201 1<br>08714729213918                                  | -                                    |
|                  | IM                     | 100 cm | 5        | H749 16391201 2<br>08714729188872                                  | H749 16599201 2<br>08714729200215    |
|                  | IM                     | 125 cm | 5        | H749 16391240 2<br>08714729188926                                  | H749 16599240 2<br>08714729304807    |
|                  | IM Tight Curve         | 100 cm | 5        | H749 16391202 2<br>08714729343295                                  | H749 16599202 2<br>08714729304647    |
|                  | IMC                    | 100 cm | 5        | H749 16391203 2<br>08714729343318                                  | H749 16599203 2<br>08714729304661    |



2. Angiographic Catheters

| IMPULSE™<br>Ventricular Pigtail                                                   |                         |        |          |                                                                    |                                      |
|-----------------------------------------------------------------------------------|-------------------------|--------|----------|--------------------------------------------------------------------|--------------------------------------|
|  | CURVE<br>STYLE          | LENGTH | QUANTITY | 5F = 1.66 MM<br>LEFT .047" = 1.2 MM ID<br>RIGHT .045" = 1.14 MM ID | 6F = 2.00 MM<br>(.056" = 1.42 MM ID) |
|                                                                                   | PIG with Side Holes     | 100 cm | 5        | H749 1639140 2<br>08714729153153                                   | H749 1659940 2<br>08714729200338     |
|                                                                                   | PIG                     | 110 cm | 1        | H749 1639140 1<br>08714729343639                                   | -                                    |
|                                                                                   | PIG with Side Holes     | 125 cm | 5        | H749 16391215 2<br>08714729153023                                  | H749 16599215 2<br>08714729200239    |
|                                                                                   | PIG 145 with Side Holes | 100 cm | 5        | H749 1639141 2<br>08714729153160                                   | H749 1659941 2<br>08714729200345     |
|                                                                                   | PIG 145                 | 110 cm | 1        | H749 1639141 1<br>08714729214342                                   | -                                    |
|                                                                                   | PIG 155 with Side Holes | 100 cm | 5        | H749 1639142 2<br>08714729153177                                   | H749 1659942 2<br>08714729200352     |

| IMPULSE™<br>Brachial - Radial Selective                                             |             |        |          |                                                                    |                                      |
|-------------------------------------------------------------------------------------|-------------|--------|----------|--------------------------------------------------------------------|--------------------------------------|
|    | CURVE STYLE | LENGTH | QUANTITY | 5F = 1.66 MM<br>LEFT .047" = 1.2 MM ID<br>RIGHT .045" = 1.14 MM ID | 6F = 2.00 MM<br>(.056" = 1.42 MM ID) |
|                                                                                     | TIG4        | 100 cm | 5        | H749 16391250 2<br>08714729886242                                  | H749 16599250 2<br>08714729886327    |
|  | TIG4.5      | 100 cm | 5        | H749 16391251 2<br>08714729886266                                  | H749 16599251 2<br>08714729886341    |
|                                                                                     | JIG3.5      | 100 cm | 5        | H749 16391252 2<br>08714729886280                                  | -                                    |
|  | JIG4        | 100 cm | 5        | H749 16391253 2<br>08714729886303                                  | -                                    |
|                                                                                     | CAS 1       | 100 cm | 5        | H749 16391171 2<br>08714729188858                                  | H749 16599171 2<br>08714729215325    |
|  | CAS 2       | 100 cm | 5        | H749 16391173 2<br>08714729152965                                  | H749 16599173 2<br>08714729215349    |
|                                                                                     | Radial      | 100 cm | 5        | H749 16391174 2<br>08714729213673                                  | H749 16599174 2<br>08714729215363    |
|  | Kimny™      | 100 cm | 5        | H749 16391180 2<br>08714729343073                                  | H749 16599180 2<br>08714729304456    |



2. Angiographic Catheters

| IMPULSE™<br>Multipack                                                             |                       |                  |          |                                                                    |                                      |
|-----------------------------------------------------------------------------------|-----------------------|------------------|----------|--------------------------------------------------------------------|--------------------------------------|
|  | CURVE<br>STYLE        | LENGTH           | QUANTITY | 5F = 1.66 MM<br>LEFT .047" = 1.2 MM ID<br>RIGHT .045" = 1.14 MM ID | 6F = 2.00 MM<br>(.056" = 1.42 MM ID) |
|                                                                                   | FL 4, FR 4, PIG       | 100, 100, 110 cm | 5        | H749 16391300 2<br>08714729188933                                  | H749 16599300 2<br>08714729200307    |
|                                                                                   | FL 4, FR 4, PIG 145   | 100, 100, 110 cm | 5        | H749 16391301 2<br>08714729188940                                  | H749 16599301 2<br>08714729200314    |
|                                                                                   | FL 4, AR MOD, PIG 145 | 100, 100, 110 cm | 5        | -                                                                  | H749 16599500 2<br>08714729305088    |
|                                                                                   | FL 4, FR 4 MOD, PIG   | 100, 100, 110 cm | 5        | H749 16391501 2<br>08714729214489                                  | -                                    |

Balloon  
Catheters

Micro-  
dilatation  
Catheters

Guide  
Catheters

Guide  
Catheter  
extension

Angiographic  
Catheters

Guidewires

Embolic  
Protection

Accessories

Cutting  
Balloon™

Atherectomy  
Devices

Chronic Total  
Occlusions

Trapper™  
Exchange  
Device

ANGIOJET™  
Thrombec-  
tomy  
System



## 2. Angiographic Catheters


### EXPO™ Femoral Selective


|                            | CURVE<br>STYLE        | LENGTH | QUANTITY | 5F = 1.66 MM<br>(.045" = 1.14 MM ID) | 6F = 2.00 MM<br>(.056" = 1.42 MM ID) | 7F = 2.33 MM<br>(.064" = 1.62 MM ID) |
|----------------------------|-----------------------|--------|----------|--------------------------------------|--------------------------------------|--------------------------------------|
| Femoral Left               | FL 3.5                | 100 cm | 5        | H749 0852621 2<br>08714729187813     | H749 0864121 2<br>08714729192992     | -<br>-                               |
|                            | FL 4                  | 100 cm | 5        | H749 0852622 2<br>08714729187820     | H749 0864122 2<br>08714729168607     | -<br>-                               |
|                            | FL 4                  | 125 cm | 5        | H749 08526220 2<br>08714729007531    | H749 08641220 2<br>08714729008453    | -<br>-                               |
|                            | FL 4.5                | 100 cm | 5        | H749 0852623 2<br>08714729187837     | H749 0864123 2<br>08714729168614     | -<br>-                               |
|                            | FL 5                  | 100 cm | 5        | H749 0852624 2<br>08714729187844     | H749 0864124 2<br>08714729168621     | H749 0873824 2<br>08714729193081     |
|                            | FL 5                  | 125 cm | 5        | H749 08526235 2<br>08714729007579    | H749 08641235 2<br>08714729008491    | -<br>-                               |
|                            | FL 6                  | 100 cm | 5        | H749 0852625 2<br>08714729187851     | H749 0864125 2<br>08714729193005     | H749 0873825 2<br>08714729193098     |
| Femoral Right              | FR 3.5                | 100 cm | 5        | H749 0852601 2<br>08714729187684     | H749 0864101 2<br>08714729168560     | -<br>-                               |
|                            | FR 4                  | 100 cm | 5        | H749 0852602 2<br>08714729187691     | H749 0864102 2<br>08714729168577     | -<br>-                               |
|                            | FR 4                  | 125 cm | 5        | -                                    | H749 08641225 2<br>08714729008460    | -<br>-                               |
|                            | FR 5                  | 100 cm | 5        | H749 0852603 2<br>08714729187707     | H749 0864103 2<br>08714729168584     | -<br>-                               |
|                            | FR 5                  | 125 cm | 5        | H749 08526230 2<br>08714729007562    | -                                    | -<br>-                               |
|                            | FR 6                  | 100 cm | 5        | H749 0852604 2<br>08714729187714     | H749 0864104 2<br>08714729192909     | -<br>-                               |
| Amplatz™<br>Left and Right | AL 1                  | 100 cm | 5        | -                                    | H749 0864196 2<br>08714729193029     | -<br>-                               |
|                            | AL 2                  | 100 cm | 5        | H749 0852698 2<br>08714729187943     | H749 0864198 2<br>08714729192893     | -<br>-                               |
|                            | AL 3                  | 100 cm | 5        | H749 08526100 2<br>08714729187721    | H749 08641100 2<br>08714729192916    | -<br>-                               |
|                            | AR 1                  | 100 cm | 5        | H749 0852684 2<br>08714729187929     | H749 0864184 2<br>08714729168638     | -<br>-                               |
|                            | AR 2                  | 100 cm | 5        | H749 0852686 2<br>08714729168874     | H749 0864186 2<br>08714729193012     | -<br>-                               |
|                            | AR MOD                | 100 cm | 5        | H749 0852689 2<br>08714729007692     | H749 0864189 2<br>08714729008613     | -<br>-                               |
| Williams                   | WR                    | 100 cm | 5        | H749 0852650 2<br>08714729249412     | H749 0864150 2<br>08714729250463     | -<br>-                               |
|                            | WRP                   | 100 cm | 5        | H749 0852651 2<br>08714729249504     | H749 0864151 2<br>08714729250555     | -<br>-                               |
| Multipurpose               | MPA 1                 | 100 cm | 5        | H749 08526117 2<br>08714729187738    | H749 08641117 2<br>08714729192923    | H749 08738117 2<br>08714729169239    |
|                            | MPA 2 with side holes | 100 cm | 5        | H749 08526121 2<br>08714729007418    | H749 08641121 2<br>08714729192930    | H749 08738121 2<br>08714729169246    |
|                            | MPA 2 with side holes | -      | 1        | H749 16599121 0<br>08714729304098    | -                                    | -<br>-                               |
|                            | MPA 2                 | 125 cm | 5        | H749 08526122 2<br>08714729187745    | H749 08641122 2<br>08714729008354    | -<br>-                               |
|                            | MPB 1                 | 100 cm | 5        | H749 08526142 2<br>08714729187752    | H749 08641142 2<br>08714729192947    | -<br>-                               |
|                            | MPB 2 with side holes | 100 cm | 5        | H749 08526146 2<br>08714729192800    | H749 08641146 2<br>08714729192954    | -<br>-                               |
| Coronary Bypass            | LCB                   | 100 cm | 5        | H749 08526195 2<br>08714729187790    | H749 08641195 2<br>08714729192978    | -<br>-                               |
|                            | RCB                   | 100 cm | 5        | H749 08526190 2<br>08714729168928    | H749 08641190 2<br>08714729192961    | -<br>-                               |
| Internal<br>Mammary        | IM                    | 100 cm | 5        | H749 08526201 2<br>08714729187806    | H749 08641201 2<br>08714729192985    | -<br>-                               |
|                            | IM                    | 125 cm | 5        | H749 08526240 2<br>08714729007593    | H749 08641240 2<br>08714729008514    | -<br>-                               |
|                            | IM Tight Curve        | 100 cm | 5        | H749 08526202 2<br>08714729249191    | -                                    | -<br>-                               |





2. Angiographic Catheters

| EXPO™<br>Ventricular Pigtail                                                      |                          |        |          |                                      |                                      |
|-----------------------------------------------------------------------------------|--------------------------|--------|----------|--------------------------------------|--------------------------------------|
|                                                                                   | CURVE<br>STYLE           | LENGTH | QUANTITY | 5F = 1.66 MM<br>(.045" = 1.14 MM ID) | 6F = 2.00 MM<br>(.056" = 1.42 MM ID) |
|  | Ventricular Pigtail      |        |          |                                      |                                      |
|                                                                                   | PIG with side holes      | 100 cm | 5        | H749 0852640 2<br>08714729187899     | H749 0864140 2<br>08714729169468     |
|                                                                                   | PIG with side holes      | 125 cm | 5        | H749 08526215 2<br>08714729007517    | H749 08641215 2<br>08714729195412    |
|                                                                                   | PIG 145° with side holes | 100 cm | 5        | H749 0852641 2<br>08714729187905     | H749 0864141 2<br>08714729169475     |
|                                                                                   | PIG 155° with side holes | 100 cm | 5        | H749 0852642 2<br>08714729187912     | H749 0864142 2<br>08714729169482     |

| EXPO™<br>Brachial - Radial Selective                                                |                |              |          |                                      |                                      |
|-------------------------------------------------------------------------------------|----------------|--------------|----------|--------------------------------------|--------------------------------------|
|                                                                                     | CURVE<br>STYLE | LENGTH       | QUANTITY | 5F = 1.66 MM<br>(.045" = 1.14 MM ID) | 6F = 2.00 MM<br>(.056" = 1.42 MM ID) |
|    | TIG4           | 100 cm       | 5        | H749 08526250 2<br>08714729886082    | H749 08641250 2<br>08714729886167    |
|                                                                                     | TIG4.5         | 100 cm       | 5        | H749 08526251 2<br>08714729886105    | H749 08641251 2<br>08714729886181    |
|   | JIG3.5         | 100 cm       | 5        | H749 08526252 2<br>08714729886129    | -                                    |
|                                                                                     | JIG4           | 100 cm       | 5        | H749 08526253 2<br>08714729886143    | -                                    |
|  | CAS 1          | 100 cm       | 5        | H749 08526171 2<br>08714729187776    | H749 08641171 2<br>08714729169000    |
|                                                                                     | CAS 2          | 100 cm       | 5        | H749 08526173 2<br>08714729187783    | H749 08641173 2<br>08714729169017    |
|  | Kimny™         | 102 / 103 cm | 5        | H749 08526180 2<br>08714729249030    | H749 08641180 2<br>08714729250036    |
|  | IMC            | 100 cm       | 5        | H749 08526203 2<br>08714729249214    | H749 08641203 2<br>08714729250210    |



2. Angiographic Catheters

| EXPO™<br>Multipack                                                                                     |                     |                  |          |                                     |                                      |
|--------------------------------------------------------------------------------------------------------|---------------------|------------------|----------|-------------------------------------|--------------------------------------|
| <div>Multipack</div>  | CURVE<br>STYLE      | LENGTH           | QUANTITY | 5F =1.66 MM<br>(.045" = 1.14 MM ID) | 6F = 2.00 MM<br>(.056" = 1.42 MM ID) |
|                                                                                                        | FL 4, FR 4, PIG     | 100, 100, 110 cm | 5        | H749 08526300 2<br>08714729187868   | H749 08641300 2<br>08714729169383    |
|                                                                                                        | FL 4, FR 4, PIG 145 | 100, 100, 110 cm | 5        | H749 08526301 2<br>08714729187875   | H749 08641301 2<br>08714729169390    |
|                                                                                                        | FL 4, FR 4, PIG 155 | 100, 100, 110 cm | 5        | H749 08526302 2<br>08714729187882   | H749 08641302 2<br>08714729169406    |

Balloon  
Catheters

Micro-  
dilatation  
Catheters

Guide  
Catheters

Guide  
Catheter  
extension

Angiographic  
Catheters

Guidewires

Embolic  
Protection

Accessories

Cutting  
Balloon™

Atherectomy  
Devices

Chronic Total  
Occlusions

Trapper™  
Exchange  
Device

ANGIOJET™  
Thrombec-  
tomy  
System



NEW CARDIO

2

## COMPLEX PCI

Balloon Catheters

Microdilatation Catheters

Guide Catheters

Guide Catheter Extension

Angiographic Catheters

Guidewires

Embolic Protection

Accessories

Cutting Balloon™

Atherectomy Devices

Chronic Total Occlusions

TRAPPER™ Exchange Device

ANGIOJET™ Thrombectomy System



Balloon  
Catheters

Micro-  
dilatation  
Catheters

Guide  
Catheters

Guide  
Catheter  
extension

Angiographic  
Catheters

Guidewires

Embolic  
Protection

Accessories

Cutting  
Balloon™

Atherectomy  
Devices

Chronic Total  
Occlusions

Trapper™  
Exchange  
Device

ANGIOJET™  
Thrombec-  
tomy  
System

## 2. Standard Lesions / Spring Tip Guidewires

### CHOICE™ Floppy

Diameter: .014" (0.37 mm)

| CATALOG NUMBER /<br>ORDER NUMBER (GTIN) | LENGTH | TIP<br>CONFIGURATION | TIP LOAD | TIP STYLE   | RAIL<br>SUPPORT | COATING | TIP<br>RADIOCAPACITY | COMPATIBLE  | QUANTITY |
|-----------------------------------------|--------|----------------------|----------|-------------|-----------------|---------|----------------------|-------------|----------|
| H749 1213201 M2<br>08714729877103       | 182 cm | Straight             | 0.8      | Spring Coil | Light           | Hybrid* | 3 cm                 | The MAGNET™ | 5        |
| H749 1213201 MJ2<br>108714729877110     | 182 cm | J                    | 0.8      | Spring Coil | Light           | Hybrid* | 3 cm                 | The MAGNET™ | 5        |
| H749 1213201 1<br>08714729303022        | 182 cm | Straight             | 0.8      | Spring Coil | Light           | Hybrid* | 3 cm                 | ADDWIRE™    | 1        |
| H749 1213201 2<br>08714729150541        | 182 cm | Straight             | 0.8      | Spring Coil | Light           | Hybrid* | 3 cm                 | ADDWIRE™    | 5        |
| H749 1213201 J1<br>08714729252351       | 182 cm | J                    | 0.8      | Spring Coil | Light           | Hybrid* | 3 cm                 | ADDWIRE™    | 1        |
| H749 1213201 J2<br>08714729155751       | 182 cm | J                    | 0.8      | Spring Coil | Light           | Hybrid* | 3 cm                 | ADDWIRE™    | 5        |
| H749 1211601 2<br>08714729150589        | 300 cm | Straight             | 0.8      | Spring Coil | Light           | Hybrid* | 3 cm                 | ADDWIRE™    | 5        |
| H749 1211601 J2<br>08714729193654       | 300 cm | J                    | 0.8      | Spring Coil | Light           | Hybrid* | 3 cm                 | ADDWIRE™    | 5        |



### LUGE™ Moderate Support Guidewire

Diameter: .014" (0.37 mm)

| CATALOG NUMBER /<br>ORDER NUMBER (GTIN) | LENGTH | TIP<br>CONFIGURATION | TIP LOAD | TIP STYLE                                 | RAIL<br>SUPPORT | COATING | TIP<br>RADIOCAPACITY | COMPATIBLE | QUANTITY |
|-----------------------------------------|--------|----------------------|----------|-------------------------------------------|-----------------|---------|----------------------|------------|----------|
| H749 1213001 0<br>08714729169956        | 182 cm | Straight             | 0.9      | Floppy; Spring Coil<br>w/a Polymer Sleeve | Moderate        | Hybrid* | 3 cm                 | ADDWIRE™   | 1        |
| H749 1213001 2<br>08714729252320        | 182 cm | Straight             | 0.9      | Floppy; Spring Coil<br>w/a Polymer Sleeve | Moderate        | Hybrid* | 3 cm                 | ADDWIRE™   | 5        |
| H749 1213001 J0<br>08714729169963       | 182 cm | J                    | 0.9      | Floppy; Spring Coil<br>w/a Polymer Sleeve | Moderate        | Hybrid* | 3 cm                 | ADDWIRE™   | 1        |
| H749 1213001 J2<br>08714729252337       | 182 cm | J                    | 0.9      | Floppy; Spring Coil<br>w/a Polymer Sleeve | Moderate        | Hybrid* | 3 cm                 | ADDWIRE™   | 5        |
| H749 1210001 2<br>08714729252207        | 300 cm | Straight             | 0.9      | Floppy; Spring Coil<br>w/a Polymer Sleeve | Moderate        | Hybrid* | 3 cm                 | ADDWIRE™   | 5        |
| H749 1210001 J0<br>08714729169949       | 300 cm | J                    | 0.9      | Floppy; Spring Coil<br>w/a Polymer Sleeve | Moderate        | Hybrid* | 3 cm                 | ADDWIRE™   | 1        |
| H749 1210001 J2<br>08714729252214       | 300 cm | J                    | 0.9      | Floppy; Spring Coil<br>w/a Polymer Sleeve | Moderate        | Hybrid* | 3 cm                 | ADDWIRE™   | 5        |



Complex PCI - 53  
NEW CARDIO

Boston  
Scientific

\* Hybrid, Hydrophilic coating Distal 3 cm uncoated.

Balloon  
Catheters

Micro-  
dilatation  
Catheters

Guide  
Catheters

Guide  
Catheter  
extension

Angiographic  
Catheters

Guidewires

Embolic  
Protection

Accessories

Cutting  
Balloon™

Atherectomy  
Devices

Chronic Total  
Occlusions

Trapper™  
Exchange  
Device

ANGIOJET™  
Thrombec-  
tomy  
System

## 2. Complex Lesions / Polymer Tip Guidewires

### PT2™ LS

Diameter: .014" (0.37 mm)

| CATALOG NUMBER / ORDER NUMBER (GTIN) | LENGTH | TIP CONFIGURATION | TIP LOAD | TIP STYLE   | RAIL SUPPORT | COATING     | TIP RADIOCAPACITY | COMPATIBLE | QUANTITY |
|--------------------------------------|--------|-------------------|----------|-------------|--------------|-------------|-------------------|------------|----------|
| H749 3893101 0<br>08714729471103     | 185 cm | Straight          | 2.5      | Polymer Tip | Light        | Hydrophilic | 2 cm              | ADDWIRE™   | 1        |
| H749 3893101 2<br>08714729471110     | 185 cm | Straight          | 2.5      | Polymer Tip | Light        | Hydrophilic | 2 cm              | ADDWIRE™   | 5        |
| H749 3893101 J0<br>08714729471127    | 185 cm | J                 | 2.5      | Polymer Tip | Light        | Hydrophilic | 2 cm              | ADDWIRE™   | 1        |
| H749 3893101 J2<br>08714729471134    | 185 cm | J                 | 2.5      | Polymer Tip | Light        | Hydrophilic | 2 cm              | ADDWIRE™   | 5        |
| H749 3893102 J2<br>08714729471172    | 300 cm | J                 | 2.5      | Polymer Tip | Light        | Hydrophilic | 2 cm              | ADDWIRE™   | 5        |
| H749 3893102 2<br>08714729471158     | 300 cm | Straight          | 2.5      | Polymer Tip | Light        | Hydrophilic | 2 cm              | ADDWIRE™   | 5        |



### CHOICE™ PT Floppy

Diameter: .014" (0.37 mm)

| CATALOG NUMBER / ORDER NUMBER (GTIN) | LENGTH | TIP CONFIGURATION | TIP LOAD | TIP STYLE   | RAIL SUPPORT | COATING     | TIP RADIOCAPACITY | COMPATIBLE  | QUANTITY |
|--------------------------------------|--------|-------------------|----------|-------------|--------------|-------------|-------------------|-------------|----------|
| H749 1216001 M2<br>108714729877264   | 182 cm | Straight          | 2.1      | Polymer Tip | Light        | Hydrophilic | 35 cm             | The MAGNET™ | 5        |
| H749 1216001 MJ2<br>108714729877271  | 182 cm | J                 | 2.1      | Polymer Tip | Light        | Hydrophilic | 35 cm             | The MAGNET™ | 5        |
| H749 1216001 1<br>08714729252504     | 182 cm | Straight          | 2.1      | Polymer Tip | Light        | Hydrophilic | 35 cm             | ADDWIRE™    | 1        |
| H749 1216001 2<br>08714729150626     | 182 cm | Straight          | 2.1      | Polymer Tip | Light        | Hydrophilic | 35 cm             | ADDWIRE™    | 5        |
| H749 1216001 J1<br>08714729252511    | 182 cm | J                 | 2.1      | Polymer Tip | Light        | Hydrophilic | 35 cm             | ADDWIRE™    | 1        |
| H749 1216001 J2<br>08714729193715    | 182 cm | J                 | 2.1      | Polymer Tip | Light        | Hydrophilic | 35 cm             | ADDWIRE™    | 5        |
| H749 1215401 2<br>08714729150633     | 300 cm | Straight          | 2.1      | Polymer Tip | Light        | Hydrophilic | 35 cm             | ADDWIRE™    | 1        |
| H749 1215401 J2<br>08714729193708    | 300 cm | J                 | 2.1      | Polymer Tip | Light        | Hydrophilic | 35 cm             | ADDWIRE™    | 5        |



### PT2™ MS

Diameter: .014" (0.37 mm)

| CATALOG NUMBER / ORDER NUMBER (GTIN) | LENGTH | TIP CONFIGURATION | TIP LOAD | TIP STYLE   | RAIL SUPPORT | COATING     | TIP RADIOCAPACITY | COMPATIBLE | QUANTITY |
|--------------------------------------|--------|-------------------|----------|-------------|--------------|-------------|-------------------|------------|----------|
| H749 1490201 1<br>08714729254928     | 182 cm | Straight          | 2.9      | Polymer Tip | Moderate     | Hydrophilic | 2 cm              | ADDWIRE™   | 1        |
| H749 1490201 J1<br>08714729254935    | 182 cm | J                 | 2.9      | Polymer Tip | Moderate     | Hydrophilic | 2 cm              | ADDWIRE™   | 1        |
| H749 3893103 0<br>08714729471189     | 185 cm | Straight          | 2.9      | Polymer Tip | Moderate     | Hydrophilic | 2 cm              | ADDWIRE™   | 1        |
| H749 3893103 2<br>08714729471196     | 185 cm | Straight          | 2.9      | Polymer Tip | Moderate     | Hydrophilic | 2 cm              | ADDWIRE™   | 5        |
| H749 3893103 J0<br>08714729471202    | 185 cm | J                 | 2.9      | Polymer Tip | Moderate     | Hydrophilic | 2 cm              | ADDWIRE™   | 1        |
| H749 3893103 J2<br>08714729471219    | 185 cm | J                 | 2.9      | Polymer Tip | Moderate     | Hydrophilic | 2 cm              | ADDWIRE™   | 5        |
| H749 3893104 0<br>08714729471226     | 300 cm | Straight          | 2.9      | Polymer Tip | Moderate     | Hydrophilic | 2 cm              | ADDWIRE™   | 1        |
| H749 3893104 2<br>08714729471233     | 300 cm | Straight          | 2.9      | Polymer Tip | Moderate     | Hydrophilic | 2 cm              | ADDWIRE™   | 5        |
| H749 3893104 J0<br>08714729471240    | 300 cm | J                 | 2.9      | Polymer Tip | Moderate     | Hydrophilic | 2 cm              | ADDWIRE™   | 1        |
| H749 3893104 J2<br>08714729471257    | 300 cm | J                 | 2.9      | Polymer Tip | Moderate     | Hydrophilic | 2 cm              | ADDWIRE™   | 5        |



Complex PCI - 54  
NEW CARDIO

Boston  
Scientific

Balloon  
Catheters

Micro-  
dilatation  
Catheters

Guide  
Catheters

Guide  
Catheter  
extension

Angiographic  
Catheters

Guidewires

Embolic  
Protection

Accessories

Cutting  
Balloon™

Atherectomy  
Devices

Chronic Total  
Occlusions

Trapper™  
Exchange  
Device

ANGIOJET™  
Thrombec-  
tomy  
System



## 2. Stiff Guidewires

### CHOICE™ Intermediate

Diameter: .014" (0.37 mm)

| CATALOG NUMBER /<br>ORDER NUMBER (GTIN) | LENGTH | TIP<br>CONFIGURATION | TIP LOAD | TIP STYLE       | RAIL<br>SUPPORT | COATING | TIP<br>RADIOCAPACITY | COMPATIBLE  | QUANTITY |
|-----------------------------------------|--------|----------------------|----------|-----------------|-----------------|---------|----------------------|-------------|----------|
| H749 1213301 1<br>08714729252368        | 182 cm | Straight             | 1.1      | Spring Coil Tip | Light           | Hybrid* | 3 cm                 | ADDWIRE™    | 1        |
| H749 1213301 M2<br>108714729877127      | 182 cm | Straight             | 1.1      | Spring Coil Tip | Light           | Hybrid* | 3 cm                 | The MAGNET™ | 5        |
| H749 1213301 2<br>08714729150558        | 182 cm | Straight             | 1.1      | Spring Coil Tip | Light           | Hybrid* | 3 cm                 | ADDWIRE™    | 5        |
| H749 1213301 J1<br>08714729252375       | 182 cm | J                    | 1.1      | Spring Coil Tip | Light           | Hybrid* | 3 cm                 | ADDWIRE™    | 1        |
| H749 1213301 J2<br>08714729155744       | 182 cm | J                    | 1.1      | Spring Coil Tip | Light           | Hybrid* | 3 cm                 | ADDWIRE™    | 5        |
| H749 1211701 2<br>08714729150596        | 300 cm | Straight             | 1.1      | Spring Coil Tip | Light           | Hybrid* | 3 cm                 | ADDWIRE™    | 5        |
| H749 1211701 J2<br>08714729193661       | 300 cm | J                    | 1.1      | Spring Coil Tip | Light           | Hybrid* | 3 cm                 | ADDWIRE™    | 5        |



### CHOICE™ Standard

Diameter: .014" (0.37 mm)

| CATALOG NUMBER /<br>ORDER NUMBER (GTIN) | LENGTH | TIP<br>CONFIGURATION | TIP LOAD | TIP STYLE       | RAIL<br>SUPPORT | COATING | TIP<br>RADIOCAPACITY | COMPATIBLE  | QUANTITY |
|-----------------------------------------|--------|----------------------|----------|-----------------|-----------------|---------|----------------------|-------------|----------|
| H749 1213401 J1<br>08714729252399       | 182 cm | J                    | 2.8      | Spring Coil Tip | Light           | Hybrid* | 3 cm                 | ADDWIRE™    | 1        |
| H749 1213401 MJ2<br>108714729877158     | 182 cm | J                    | 2.8      | Spring Coil Tip | Light           | Hybrid* | 3 cm                 | The MAGNET™ | 5        |
| H749 1213401 J2<br>08714729155737       | 182 cm | J                    | 2.8      | Spring Coil Tip | Light           | Hybrid* | 3 cm                 | ADDWIRE™    | 5        |
| H749 1213401 2<br>08714729150565        | 182 cm | Straight             | 2.8      | Spring Coil Tip | Light           | Hybrid* | 3 cm                 | ADDWIRE™    | 5        |
| H749 1213401 1<br>08714729252382        | 182 cm | Straight             | 2.8      | Spring Coil Tip | Light           | Hybrid* | 3 cm                 | ADDWIRE™    | 1        |
| H749 1211801 2<br>08714729150602        | 300 cm | Straight             | 2.8      | Spring Coil Tip | Light           | Hybrid* | 3 cm                 | ADDWIRE™    | 5        |



Complex PCI - 55  
NEW CARDIO

Boston  
Scientific

\* Hybrid, Hydrophilic coating Distal 3 cm uncoated.

Balloon  
Catheters

Micro-  
dilatation  
Catheters

Guide  
Catheters

Guide  
Catheter  
extension

Angiographic  
Catheters

Guidewires

Embolic  
Protection

Accessories

Cutting  
Balloon™

Atherectomy  
Devices

Chronic Total  
Occlusions

Trapper™  
Exchange  
Device

ANGIOJET™  
Thrombec-  
tomy  
System

2. Stiff Guidewires

PT GRAPHIX™ Intermediate

Diameter: .014" (0.37 mm)

| CATALOG NUMBER / ORDER NUMBER (GTIN) | LENGTH | TIP CONFIGURATION | TIP LOAD | TIP STYLE   | RAIL SUPPORT | COATING     | TIP RADIOPACITY | COMPATIBLE  | QUANTITY |
|--------------------------------------|--------|-------------------|----------|-------------|--------------|-------------|-----------------|-------------|----------|
| H749 1490201 M2<br>08714729877202    | 182 cm | Straight          | 1.7      | Polymer Tip | Moderate     | Hydrophilic | 2 cm            | The MAGNET™ | 5        |
| H749 1490201 MJ2<br>08714729877219   | 182 cm | J                 | 1.7      | Polymer Tip | Moderate     | Hydrophilic | 2 cm            | The MAGNET™ | 5        |
| H749 1490201 2<br>08714729141723     | 182 cm | Straight          | 1.7      | Polymer Tip | Moderate     | Hydrophilic | 2 cm            | ADDWIRE™    | 5        |
| H749 1490201 J2<br>08714729141730    | 182 cm | J                 | 1.7      | Polymer Tip | Moderate     | Hydrophilic | 2 cm            | ADDWIRE™    | 5        |
| H749 1491401 1<br>08714729255017     | 300 cm | Straight          | 1.7      | Polymer Tip | Moderate     | Hydrophilic | 2 cm            | ADDWIRE™    | 1        |
| H749 1491401 2<br>08714729141969     | 300 cm | Straight          | 1.7      | Polymer Tip | Moderate     | Hydrophilic | 2 cm            | ADDWIRE™    | 5        |
| H749 1491401 J2<br>08714729141976    | 300 cm | J                 | 1.7      | Polymer Tip | Moderate     | Hydrophilic | 2 cm            | ADDWIRE™    | 5        |



PT GRAPHIX™ Standard

Diameter: .014" (0.37 mm)

| CATALOG NUMBER / ORDER NUMBER (GTIN) | LENGTH | TIP CONFIGURATION | TIP LOAD | TIP STYLE   | RAIL SUPPORT | COATING     | TIP RADIOPACITY | COMPATIBLE  | QUANTITY |
|--------------------------------------|--------|-------------------|----------|-------------|--------------|-------------|-----------------|-------------|----------|
| H749 1490301 MJ2<br>08714729877219   | 182 cm | J                 | 5.3      | Polymer Tip | Moderate     | Hydrophilic | 2 cm            | The MAGNET™ | 5        |
| H749 1490301 2<br>08714729188469     | 182 cm | Straight          | 5.3      | Polymer Tip | Moderate     | Hydrophilic | 2 cm            | ADDWIRE™    | 5        |
| H749 1490301 J2<br>08714729141754    | 182 cm | J                 | 5.3      | Polymer Tip | Moderate     | Hydrophilic | 2 cm            | ADDWIRE™    | 5        |



## 2. Extra Support Guidewires

### CHOICE™ Extra Support

Diameter: .014" (0.37 mm)

| CATALOG NUMBER /<br>ORDER NUMBER (GTIN) | LENGTH | TIP<br>CONFIGURATION | TIP LOAD | TIP STYLE       | RAIL<br>SUPPORT | COATING | TIP<br>RADIOCAPACITY | COMPATIBLE  | QUANTITY |
|-----------------------------------------|--------|----------------------|----------|-----------------|-----------------|---------|----------------------|-------------|----------|
| H749 1213501 1<br>08714729252405        | 182 cm | Straight             | 0.9      | Spring Coil Tip | Extra           | Hybrid* | 3 cm                 | ADDWIRE™    | 1        |
| H749 1213501 M2<br>108714729877165      | 182 cm | Straight             | 0.9      | Spring Coil Tip | Extra           | Hybrid* | 3 cm                 | The MAGNET™ | 5        |
| H749 1213501 J1<br>08714729252412       | 182 cm | Straight             | 0.9      | Spring Coil Tip | Extra           | Hybrid* | 3 cm                 | ADDWIRE™    | 1        |
| H749 1213501 MJ2<br>108714729877172     | 182 cm | J                    | 0.9      | Spring Coil Tip | Extra           | Hybrid* | 3 cm                 | The MAGNET™ | 5        |
| H749 1213501 2<br>08714729150572        | 182 cm | Straight             | 0.9      | Spring Coil Tip | Extra           | Hybrid* | 3 cm                 | ADDWIRE™    | 5        |
| H749 1213501 J2<br>08714729193692       | 182 cm | J                    | 0.9      | Spring Coil Tip | Extra           | Hybrid* | 3 cm                 | ADDWIRE™    | 5        |
| H749 1211901 2<br>08714729150619        | 300 cm | Straight             | 0.9      | Spring Coil Tip | Extra           | Hybrid* | 3 cm                 | ADDWIRE™    | 5        |
| H749 1211901 J2<br>08714729193685       | 300 cm | J                    | 0.9      | Spring Coil Tip | Extra           | Hybrid* | 3 cm                 | ADDWIRE™    | 5        |



### CHOICE™ PT Extra Support

Diameter: .014" (0.37 mm)

| CATALOG NUMBER /<br>ORDER NUMBER (GTIN) | LENGTH | TIP<br>CONFIGURATION | TIP LOAD | TIP STYLE   | RAIL<br>SUPPORT | COATING     | TIP<br>RADIOCAPACITY | COMPATIBLE  | QUANTITY |
|-----------------------------------------|--------|----------------------|----------|-------------|-----------------|-------------|----------------------|-------------|----------|
| H749 1216101 1<br>08714729252528        | 182 cm | Straight             | 2.0      | Polymer Tip | Extra           | Hydrophilic | 35 cm                | ADDWIRE™    | 1        |
| H749 1216101 M2<br>08714729877288       | 182 cm | Straight             | 2.0      | Polymer Tip | Extra           | Hydrophilic | 35 cm                | The MAGNET™ | 5        |
| H749 1216101 J1<br>08714729252535       | 182 cm | J                    | 2.0      | Polymer Tip | Extra           | Hydrophilic | 35 cm                | ADDWIRE™    | 1        |
| H749 1216101 MJ2<br>108714729877295     | 182 cm | J                    | 2.0      | Polymer Tip | Extra           | Hydrophilic | 35 cm                | The MAGNET™ | 5        |
| H749 1216101 2<br>08714729176992        | 182 cm | Straight             | 2.0      | Polymer Tip | Extra           | Hydrophilic | 35 cm                | ADDWIRE™    | 5        |
| H749 1216101 J2<br>08714729177005       | 182 cm | J                    | 2.0      | Polymer Tip | Extra           | Hydrophilic | 35 cm                | ADDWIRE™    | 5        |
| H749 1215501 2<br>08714729177012        | 300 cm | Straight             | 2.0      | Polymer Tip | Extra           | Hydrophilic | 35 cm                | ADDWIRE™    | 5        |
| H749 1215501 J2<br>08714729176947       | 300 cm | J                    | 2.0      | Polymer Tip | Extra           | Hydrophilic | 35 cm                | ADDWIRE™    | 5        |



Complex PCI - 57  
NEW CARDIO

Boston  
Scientific

\* Hybrid, Hydrophilic coating Distal 3 cm uncoated.

Balloon  
Catheters

Micro-  
dilatation  
Catheters

Guide  
Catheters

Guide  
Catheter  
extension

Angiographic  
Catheters

Guidewires

Embolic  
Protection

Accessories

Cutting  
Balloon™

Atherectomy  
Devices

Chronic Total  
Occlusions

Trapper™  
Exchange  
Device

ANGIOJET™  
Thrombec-  
tomy  
System

2. Extra Support Guidewires

| MAILMAN™                             |        |                   |          |                 |              |         |                 |             |          |
|--------------------------------------|--------|-------------------|----------|-----------------|--------------|---------|-----------------|-------------|----------|
| Diameter: .014" (0.37 mm)            |        |                   |          |                 |              |         |                 |             |          |
| CATALOG NUMBER / ORDER NUMBER (GTIN) | LENGTH | TIP CONFIGURATION | TIP LOAD | TIP STYLE       | RAIL SUPPORT | COATING | TIP RADIOPACITY | COMPATIBLE  | QUANTITY |
| H749 1213601 M2<br>108714729877189   | 182 cm | Straight          | 0.8      | Spring Coil Tip | Super        | Hybrid* | 3 cm            | The MAGNET™ | 5        |
| H749 1213601 MJ2<br>108714729877196  | 182 cm | J                 | 0.8      | Spring Coil Tip | Super        | Hybrid* | 3 cm            | The MAGNET™ | 5        |
| H749 1213601 2<br>08714729165712     | 182 cm | Straight          | 0.8      | Spring Coil Tip | Super        | Hybrid* | 3 cm            | ADDWIRE™    | 5        |
| H749 1213601 J2<br>08714729165729    | 182 cm | J                 | 0.8      | Spring Coil Tip | Super        | Hybrid* | 3 cm            | ADDWIRE™    | 5        |
| H749 1212001 2<br>08714729188162     | 300 cm | Straight          | 0.8      | Spring Coil Tip | Super        | Hybrid* | 3 cm            | ADDWIRE™    | 5        |
| H749 1212001 J2<br>08714729165705    | 300 cm | J                 | 0.8      | Spring Coil Tip | Super        | Hybrid* | 3 cm            | ADDWIRE™    | 5        |



| PT GRAPHIX™ Super Support            |        |                   |          |             |              |             |                 |            |          |
|--------------------------------------|--------|-------------------|----------|-------------|--------------|-------------|-----------------|------------|----------|
| Diameter: .014" (0.37 mm)            |        |                   |          |             |              |             |                 |            |          |
| CATALOG NUMBER / ORDER NUMBER (GTIN) | LENGTH | TIP CONFIGURATION | TIP LOAD | TIP STYLE   | RAIL SUPPORT | COATING     | TIP RADIOPACITY | COMPATIBLE | QUANTITY |
| H749 1490501 2<br>08714729141785     | 182 cm | Straight          | 3.6      | Polymer Tip | Super        | Hydrophilic | 2 cm            | ADDWIRE™   | 5        |
| H749 1490501 J2<br>08714729188476    | 182 cm | J                 | 3.6      | Polymer Tip | Super        | Hydrophilic | 2 cm            | ADDWIRE™   | 5        |



\* Hybrid, Hydrophilic coating Distal 3 cm uncoated.

## 2. SENTAI Workhorse Wires

### MARVEL™ Wire

Diameter: .014" (0.37 mm)

| CATALOG NUMBER / ORDER NUMBER (GTIN) | LENGTH | TIP CONFIGURATION | TIP LOAD (GF) | TIP DIAMETER (INCH / MM) | TIP RADIOCAPACITY | CORE MATERIAL   | COATING     | COMPATIBLE | QUANTITY |
|--------------------------------------|--------|-------------------|---------------|--------------------------|-------------------|-----------------|-------------|------------|----------|
| H749 39301190 2<br>08714729881179    | 190 cm | Straight          | 0.9           | 0.014 / 0.37             | 3 cm              | Stainless Steel | Hydrophilic | STRETCH    | 5        |
| H749 39301300 2<br>08714729881193    | 300 cm | Straight          | 0.9           | 0.014 / 0.37             | 3 cm              | Stainless Steel | Hydrophilic | -          | 5        |
| H749 39301190 J2<br>08714729881186   | 190 cm | J                 | 0.9           | 0.014 / 0.37             | 3 cm              | Stainless Steel | Hydrophilic | STRETCH    | 5        |
| H749 39301300 J2<br>08714729881247   | 300 cm | J                 | 0.9           | 0.014 / 0.37             | 3 cm              | Stainless Steel | Hydrophilic | -          | 5        |



### SAMURAI™ Wire

Diameter: .014" (0.37 mm)

| CATALOG NUMBER / ORDER NUMBER (GTIN) | LENGTH | TIP CONFIGURATION | TIP LOAD (GF) | TIP DIAMETER (INCH / MM) | TIP RADIOCAPACITY | CORE MATERIAL   | COATING                 | COMPATIBLE | QUANTITY |
|--------------------------------------|--------|-------------------|---------------|--------------------------|-------------------|-----------------|-------------------------|------------|----------|
| H749 39302190 2<br>08714729881216    | 190 cm | Straight          | 0.5           | 0.014 / 0.37             | 4 cm              | Stainless Steel | Moderated hydrophilic** | STRETCH    | 5        |
| H749 39302300 2<br>08714729881230    | 300 cm | Straight          | 0.5           | 0.014 / 0.37             | 4 cm              | Stainless Steel | Moderated hydrophilic** | -          | 5        |
| H749 39302190 J2<br>08714729881223   | 190 cm | J                 | 0.5           | 0.014 / 0.37             | 4 cm              | Stainless Steel | Moderated hydrophilic** | STRETCH    | 5        |
| H749 39302300 J2<br>08714729881209   | 300 cm | J                 | 0.5           | 0.014 / 0.37             | 4 cm              | Stainless Steel | Moderated hydrophilic** | -          | 5        |

\*\* Reduced hydrophilic coating on distal 1 cm



Complex PCI - 59  
NEW CARDIO

Boston  
Scientific

Balloon  
Catheters

Micro-  
dilatation  
Catheters

Guide  
Catheters

Guide  
Catheter  
extension

Angiographic  
Catheters

Guidewires

Embolic  
Protection

Accessories

Cutting  
Balloon™

Atherectomy  
Devices

Chronic Total  
Occlusions

Trapper™  
Exchange  
Device

ANGIOJET™  
Thrombec-  
tomy  
System



## 2. SENTAI Specialty Crossing Wires

### SAMURAI™ RC Wire

Diameter: .014" (0.37 mm)

| CATALOG NUMBER / ORDER NUMBER (GTIN) | LENGTH | TIP CONFIGURATION | TIP LOAD (GF) | TIP DIAMETER (INCH / MM) | TIP RADIOCAPACITY | CORE MATERIAL   | COATING     | COMPATIBLE | QUANTITY |
|--------------------------------------|--------|-------------------|---------------|--------------------------|-------------------|-----------------|-------------|------------|----------|
| H749 39306190 0<br>08714729887737    | 190 cm | Straight          | 1.2           | 0.014 / 0.37             | 4 cm              | Stainless Steel | Hydrophilic | STRETCH    | 1        |
| H749 39306190 2<br>08714729886679    | 190 cm | Straight          | 1.2           | 0.014 / 0.37             | 4 cm              | Stainless Steel | Hydrophilic | STRETCH    | 5        |
| H749 39306300 2<br>08714729886686    | 300 cm | Straight          | 1.2           | 0.014 / 0.37             | 4 cm              | Stainless Steel | Hydrophilic | -          | 5        |



### FIGHTER™ Wire

Diameter: .014" (0.37 mm)

| CATALOG NUMBER / ORDER NUMBER (GTIN) | LENGTH | TIP CONFIGURATION | TIP LOAD (GF) | TIP DIAMETER (INCH / MM) | TIP RADIOCAPACITY | CORE MATERIAL   | COATING                            | COMPATIBLE | QUANTITY |
|--------------------------------------|--------|-------------------|---------------|--------------------------|-------------------|-----------------|------------------------------------|------------|----------|
| H749 39303190 2<br>08714729886617    | 190 cm | Straight          | 1.5           | 0.009 / 0.23             | 3.5 cm            | Stainless Steel | Clear polymer jacket / hydrophilic | STRETCH    | 5        |
| H749 39303300 2<br>08714729886624    | 300 cm | Straight          | 1.5           | 0.009 / 0.23             | 3.5 cm            | Stainless Steel | Clear polymer jacket / hydrophilic | -          | 5        |



### JUDO™ 1 Wire

Diameter: .014" (0.37 mm)

| CATALOG NUMBER / ORDER NUMBER (GTIN) | LENGTH | TIP CONFIGURATION | TIP LOAD (GF) | TIP DIAMETER (INCH / MM) | TIP RADIOCAPACITY | CORE MATERIAL   | COATING     | COMPATIBLE | QUANTITY |
|--------------------------------------|--------|-------------------|---------------|--------------------------|-------------------|-----------------|-------------|------------|----------|
| H749 39415190 2<br>08714729975038    | 190 cm | Straight          | 1.0           | 0.008 / 0.20             | 3.5 cm            | Stainless Steel | Hydrophilic | STRETCH    | 5        |
| H749 39415300 2<br>08714729975052    | 300 cm | Straight          | 1.0           | 0.008 / 0.20             | 3.5 cm            | Stainless Steel | Hydrophilic | -          | 5        |

### JUDO™ 3 Wire

Diameter: .014" (0.37 mm)

| CATALOG NUMBER / ORDER NUMBER (GTIN) | LENGTH | TIP CONFIGURATION | TIP LOAD (GF) | TIP DIAMETER (INCH / MM) | TIP RADIOCAPACITY | CORE MATERIAL   | COATING     | COMPATIBLE | QUANTITY |
|--------------------------------------|--------|-------------------|---------------|--------------------------|-------------------|-----------------|-------------|------------|----------|
| H749 39377190 2<br>08714729974956    | 190 cm | Straight          | 3.0           | 0.008 / 0.20             | 3.5 cm            | Stainless Steel | Hydrophilic | STRETCH    | 5        |
| H749 39377300 2<br>08714729974970    | 300 cm | Straight          | 3.0           | 0.008 / 0.20             | 3.5 cm            | Stainless Steel | Hydrophilic | -          | 5        |

### JUDO™ 6 Wire

Diameter: .014" (0.37 mm)

| CATALOG NUMBER / ORDER NUMBER (GTIN) | LENGTH | TIP CONFIGURATION | TIP LOAD (GF) | TIP DIAMETER (INCH / MM) | TIP RADIOCAPACITY | CORE MATERIAL   | COATING     | COMPATIBLE | QUANTITY |
|--------------------------------------|--------|-------------------|---------------|--------------------------|-------------------|-----------------|-------------|------------|----------|
| H749 39378190 2<br>08714729974994    | 190 cm | Straight          | 6.0           | 0.008 / 0.20             | 3.5 cm            | Stainless Steel | Hydrophilic | STRETCH    | 5        |
| H749 39378300 2<br>08714729975014    | 300 cm | Straight          | 6.0           | 0.008 / 0.20             | 3.5 cm            | Stainless Steel | Hydrophilic | -          | 5        |



Complex PCI - 60  
NEW CARDIO

Boston  
Scientific

Balloon  
Catheters

Micro-  
dilatation  
Catheters

Guide  
Catheters

Guide  
Catheter  
extension

Angiographic  
Catheters

Guidewires

Embolic  
Protection

Accessories

Cutting  
Balloon™

Atherectomy  
Devices

Chronic Total  
Occlusions

Trapper™  
Exchange  
Device

ANGIOJET™  
Thrombec-  
tomy  
System

2. SENTAI Specialty Crossing Wires

HORNET™ 10 Wire

Diameter: .014" (0.37 mm)

| CATALOG NUMBER / ORDER NUMBER (GTIN) | LENGTH | TIP CONFIGURATION | TIP LOAD (GF) | TIP DIAMETER (INCH / MM) | TIP RADIOPACITY | CORE MATERIAL   | COATING     | COMPATIBLE | QUANTITY |
|--------------------------------------|--------|-------------------|---------------|--------------------------|-----------------|-----------------|-------------|------------|----------|
| H749 39304190 2<br>08714729886631    | 190 cm | Straight          | 10            | 0.008 / 0.20             | 3.5 cm          | Stainless Steel | Hydrophilic | STRETCH    | 5        |
| H749 39304300 2<br>08714729886648    | 300 cm | Straight          | 10            | 0.008 / 0.20             | 3.5 cm          | Stainless Steel | Hydrophilic | -          | 5        |



HORNET™ 14 Wire

Diameter: .014" (0.37 mm)

| CATALOG NUMBER / ORDER NUMBER (GTIN) | LENGTH | TIP CONFIGURATION | TIP LOAD (GF) | TIP DIAMETER (INCH / MM) | TIP RADIOPACITY | CORE MATERIAL   | COATING     | COMPATIBLE | QUANTITY |
|--------------------------------------|--------|-------------------|---------------|--------------------------|-----------------|-----------------|-------------|------------|----------|
| H749 39307190 2<br>08714729886693    | 190 cm | Straight          | 14            | 0.008 / 0.20             | 3.5 cm          | Stainless Steel | Hydrophilic | STRETCH    | 5        |
| H749 39307300 2<br>08714729886709    | 300 cm | Straight          | 14            | 0.008 / 0.20             | 3.5 cm          | Stainless Steel | Hydrophilic | -          | 5        |




2. Extension Guidewires

ADDWIRE™ Extension Wire

ADDWIRE™ Extension Wire is compatible with 182-185 cm CHOICE™ Floppy, LUGE™ Moderate Support, CHOICE™ Intermediate, CHOICE™ PT Floppy, PT2™ Light Support, CHOICE™ Standard, PT Graphix™ Intermediate, PT2™ Moderate Support, CHOICE™ Extra Support, CHOICE™ PT Extra Support, PT Graphix™ Super Support, PT Graphix™ Standard and MAILMAN™ guidewires. Refer to the guidewire product label for Guidewire Extension system compatibility.

| CATALOG NUMBER / ORDER NUMBER (GTIN) | DIAMETER        | LENGTH | QUANTITY |
|--------------------------------------|-----------------|--------|----------|
| H749 2215001 2 / 08714729315940      | .014" (0.37 mm) | 145 cm | 5        |




The ADDWIRE™ is intended to be attached to a Boston Scientific Extendable Guidewire.

The 38 MAGNET™ Device

The MAGNET™ Device is compatible with 182 cm lengths of CHOICE™, CHOICE™ PT, MAILMAN™ and PT GRAPHIX™ wires (M in Order Number refers the products to be compatible with MAGNET™ Device. Short Wires without M are compatible with ADDWIRE™ Extension Wire).

| CATALOG NUMBER / ORDER NUMBER (GTIN) | DESCRIPTION | QUANTITY |
|--------------------------------------|-------------|----------|
| H749 0896601 1 / 08714729146841      | The MAGNET™ | 1        |
| H749 0896605 2 / 08714729008828      | The MAGNET™ | 5        |



STRETCH Extension Wire\*

Extension wire for 190 cm length MARVEL, SAMURAI, SAMURAI RC, FIGHTER, JUDO 1, JUDO 3, JUDO 6, HORNET, HORNET 10 and HORNET 14 wires.

Diameter: .014" (0.37 mm)

| CATALOG NUMBER / ORDER NUMBER (GTIN) | LENGTH | QUANTITY |
|--------------------------------------|--------|----------|
| H7493930802 / 08714729886716         | 150 cm | 5        |



Complex PCI - 62  
NEW CARDIO

Boston  
Scientific

\* Please check availability of this product with your local sales representative.

Balloon  
Catheters

Micro-  
dilatation  
Catheters

Guide  
Catheters

Guide  
Catheter  
extension

Angiographic  
Catheters

Guidewires

Embolic  
Protection

Accessories

Cutting  
Balloon™

Atherectomy  
Devices

Chronic Total  
Occlusions

Trapper™  
Exchange  
Device

ANGIOJET™  
Thrombec-  
tomy  
System

# 2

## COMPLEX PCI

Balloon Catheters

Microdilatation Catheters

Guide Catheters

Guide Catheter Extension

Angiographic Catheters

Guidewires

**Embolic Protection**

Accessories


Cutting Balloon™

Atherectomy Devices

Chronic Total Occlusions

TRAPPER™ Exchange Device

ANGIOJET™ Thrombectomy System



Balloon  
Catheters

Micro-  
dilatation  
Catheters

Guide  
Catheters

Guide  
Catheter  
extension

Angiographic  
Catheters

Guidewires

**Embolic  
Protection**

Accessories

Cutting  
Balloon™

Atherectomy  
Devices

Chronic Total  
Occlusions

Trapper™  
Exchange  
Device

ANGIOJET™  
Thrombec-  
tomy  
System

FILTERWIRE EZ™ Embolic Protection System

Material: Stainless steel guidewire, Platinum coated Nitinol loop, microporous 110 micron filter, Floppy wire tip 3 cm, 6F guide catheter compatible.

| CATALOG NUMBER / ORDER NUMBER (GTIN) | DESCRIPTION                  | VESSEL SIZE COVERAGE | LENGTH | CROSSING PROFILE | RETRIEVAL PROFILE |
|--------------------------------------|------------------------------|----------------------|--------|------------------|-------------------|
| H749 20100399 0<br>08714729470878    | FILTERWIRE EZ™               | 3.50 - 5.50 mm       | 190 cm | 3.2F / 1.07 mm   | 4.3F / 1.43 mm    |
| H749 50100150 0<br>08714729479109    | Bent Tip Retrieval Sheath EZ | N / A                | 150 cm | 4.3F / 1.43 mm   | 4.3F / 1.43 mm    |
| H749 20100400 0<br>08714729470885    | FILTERWIRE EZ™ 300 cm        | 3.50 - 5.50 mm       | 300 cm | 3.2F / 1.07 mm   | 4.3F / 1.43 mm    |

The 190 cm protection wire is compatible with the Boston Scientific ADDWIRE™ Extension Wire.

Balloon Catheters

Micro-dilatation Catheters

Guide Catheters

Guide Catheter extension

Angiographic Catheters

Guidewires

Embolic Protection

Accessories

Cutting Balloon™

Atherectomy Devices

Chronic Total Occlusions

Trapper™ Exchange Device

ANGIOJET™ Thrombectomy System





# 2

## COMPLEX PCI

Balloon Catheters

Microdilatation Catheters

Guide Catheters

Guide Catheter Extension

Angiographic Catheters

Guidewires

Embolic Protection

Accessories

Cutting Balloon™

Atherectomy Devices

Chronic Total Occlusions

TRAPPER™ Exchange Device

ANGIOJET™ Thrombectomy System



Balloon  
Catheters

Micro-  
dilatation  
Catheters

Guide  
Catheters

Guide  
Catheter  
extension

Angiographic  
Catheters

Guidewires

Embolic  
Protection

Accessories

Cutting  
Balloon™

Atherectomy  
Devices


Chronic Total  
Occlusions

Trapper™  
Exchange  
Device

ANGIOJET™  
Thrombec-  
tomy  
System

ENCORE™ 26 Inflation Device, 26 ATM (2634 kPa) gauge


| CATALOG NUMBER / ORDER NUMBER (GTIN) | MODEL                       | VOLUME CAPACITY | QUANTITY |
|--------------------------------------|-----------------------------|-----------------|----------|
| H749 0452601 1 / 08714729177029      | ENCORE™ 26 Inflation Device | 20 ml           | 1        |
| H749 0452605 2 / 08714729127062      | ENCORE™ 26 Inflation Device | 20 ml           | 5-pack   |



ENCORE™ 26 Advantage Kit, 26 ATM (2634 kPa) gauge

The ENCORE™ 26 Advantage Kit combines the GATEWAY™ Plus™ Y-Adaptor, the TD2™ torque device, an insertion tool, and the ENCORE™ 26 high pressure inflation device into one convenient kit.


| CATALOG NUMBER / ORDER NUMBER (GTIN) | MODEL                    | VOLUME CAPACITY | QUANTITY |
|--------------------------------------|--------------------------|-----------------|----------|
| H749 0452701 1 / 08714729180005      | ENCORE™ 26 Advantage Kit | 20 ml           | 1        |
| H749 0452705 2 / 08714729127048      | ENCORE™ 26 Advantage Kit | 20 ml           | 5-pack   |



2. Y-Adaptors

GATEWAY™ Plus Y-Adaptor

| CATALOG NUMBER / ORDER NUMBER (GTIN) | MODEL                   | QUANTITY |
|--------------------------------------|-------------------------|----------|
| H749 0448310 2 / 08714729246763      | GATEWAY™ Plus Y-Adaptor | 10-pack  |



Balloon Catheters

Micro-dilatation Catheters

Guide Catheters

Guide Catheter extension

Angiographic Catheters

Guidewires

Embolic Protection

Accessories

Cutting Balloon™

Atherectomy Devices

Chronic Total Occlusions

Trapper™ Exchange Device

ANGIOJET™ Thrombectomy System




2. Hemostasis Valve

WATCHDOG  
Hemostasis Valve

| CATALOG NUMBER / ORDER NUMBER (GTIN) | DESCRIPTION                    | GENERAL SPECIFICATIONS                         | QUANTITY             |
|--------------------------------------|--------------------------------|------------------------------------------------|----------------------|
| H749 393430 2 / 08714729996002       | WATCHDOG Hemostasis Valve      | 3.4" (8.64 cm) / 0.105" (2.67 mm) / Cross-Slit | Carton: 1, Pack: 5** |
| H749 3934302 1 / 08714729965763      | WATCHDOG Hemostasis Valve      | 3.4" (8.64 cm) / 0.105" (2.67 mm) / Cross-Slit | Carton: 1, Pack: 20  |
| H749 39343A02 2 / 08714729996019     | WATCHDOG Hemostasis Valve Kit* | 3.4" (8.64 cm) / 0.105" (2.67 mm) / Cross-Slit | Carton: 1, Pack: 5** |
| H749 39343A02 1 / 08714729965770     | WATCHDOG Hemostasis Valve Kit* | 3.4" (8.64 cm) / 0.105" (2.67 mm) / Cross-Slit | Carton: 1, Pack: 20  |

\*Kit configuration includes hemostasis valve, insertion tool and torque device    \*\*available to order starting from June, 2020



## 2. Introducer Sheaths

### SUPER SHEATH™ 0.889 mm (0.035") compatible

| CATALOG NUMBER / ORDER NUMBER (GTIN) | ID           | SHEATH LENGTH | MAX GUIDEWIRE     | MATERIALS INCLUDED                          | UNITS PER BOX |
|--------------------------------------|--------------|---------------|-------------------|---------------------------------------------|---------------|
| H749 3908601 B1<br>08714729761471    | 4F (1.33 mm) | 7 cm          | 0.035" (0.889 mm) | Introducer Sheath,<br>Dilator, Guidewire    | 10            |
| H749 3908602 B1<br>08714729761488    | 5F (1.67 mm) | 7 cm          | 0.035" (0.889 mm) |                                             | 10            |
| H749 1603404 B1<br>14543527181393    | 4F (1.33 mm) | 11 cm         | 0.035" (0.889 mm) | Introducer Sheath,<br>Dilator,<br>Guidewire | 10            |
| H749 3908613 B1<br>14543527196403    | 5F (1.67 mm) | 11 cm         | 0.035" (0.889 mm) |                                             | 10            |
| H749 3908614 B1<br>14543527196410    | 6F (2.00 mm) | 11 cm         | 0.035" (0.889 mm) |                                             | 10            |
| H749 3908615 B1<br>14543527196427    | 7F (2.33 mm) | 11 cm         | 0.035" (0.889 mm) |                                             | 10            |
| H749 3908616 B1<br>14543527196434    | 8F (2.67 mm) | 11 cm         | 0.035" (0.889 mm) |                                             | 10            |
| H749 3908617 B1<br>08714729761563    | 9F (3.00 mm) | 11 cm         | 0.035" (0.889 mm) |                                             | 10            |
| H749 3908624 B1<br>08714729761570    | 5F (1.67 mm) | 25 cm         | 0.035" (0.889 mm) | Introducer Sheath,<br>Dilator               | 10            |
| H749 3908625 B1<br>08714729761587    | 6F (2.00 mm) | 25 cm         | 0.035" (0.889 mm) |                                             | 10            |
| H749 3908626 B1<br>14543527196533    | 7F (2.33 mm) | 25 cm         | 0.035" (0.889 mm) |                                             | 10            |
| H749 3908627 B1<br>14543527196540    | 8F (2.67 mm) | 25 cm         | 0.035" (0.889 mm) |                                             | 10            |

### SUPER SHEATH™ 0.965 mm (0.038") compatible

| CATALOG NUMBER / ORDER NUMBER (GTIN) | ID           | SHEATH LENGTH | MAX GUIDEWIRE     | MATERIALS INCLUDED                          | UNITS PER BOX |
|--------------------------------------|--------------|---------------|-------------------|---------------------------------------------|---------------|
| H749 3908630 B1<br>14543527196571    | 6F (2.00 mm) | 7 cm          | 0.038" (0.965 mm) | Introducer Sheath,<br>Dilator, Guidewire    | 10            |
| H749 1603505 B1<br>14543527181409    | 5F (1.67 mm) | 11 cm         | 0.038" (0.965 mm) |                                             | 10            |
| H749 1603506 B1<br>14543527181416    | 6F (2.00 mm) | 11 cm         | 0.038" (0.965 mm) | Introducer Sheath,<br>Dilator,<br>Guidewire | 10            |
| H749 1603507 B1<br>14543527181423    | 7F (2.33 mm) | 11 cm         | 0.038" (0.965 mm) |                                             | 10            |
| H749 1603508 B1<br>14543527181430    | 8F (2.67 mm) | 11 cm         | 0.038" (0.965 mm) |                                             | 10            |
| H749 1603509 B1<br>14543527181447    | 9F (3.00 mm) | 11 cm         | 0.038" (0.965 mm) |                                             | 10            |
| H749 1603705 B1<br>08714729756712    | 5F (1.67 mm) | 25 cm         | 0.038" (0.965 mm) |                                             | 10            |
| H749 1603706 B1<br>14543527181522    | 6F (2.00 mm) | 25 cm         | 0.038" (0.965 mm) |                                             | 10            |
| H749 1603707 B1<br>14543527181539    | 7F (2.33 mm) | 25 cm         | 0.038" (0.965 mm) | Introducer Sheath,<br>Dilator               | 10            |
| H749 1603708 B1<br>14543527181546    | 8F (2.67 mm) | 25 cm         | 0.038" (0.965 mm) |                                             | 10            |
| H749 1603709 B1<br>08714729756798    | 9F (3.00 mm) | 25 cm         | 0.038" (0.965 mm) |                                             | 10            |



# 2

## COMPLEX PCI

Balloon Catheters

Microdilatation Catheters

Guide Catheters

Guide Catheter Extension

Angiographic Catheters

Guidewires

Embolic Protection

Accessories

**Cutting Balloon™**

Atherectomy Devices

Chronic Total Occlusions

TRAPPER™ Exchange Device

ANGIOJET™ Thrombectomy System



Balloon Catheters

Micro-dilatation Catheters

Guide Catheters

Guide Catheter extension

Angiographic Catheters

Guidewires

Embolic Protection

Accessories

**Cutting Balloon™**

Atherectomy Devices

Chronic Total Occlusions

Trapper™ Exchange Device

ANGIOJET™ Thrombectomy System



## 2. Cutting Balloon™

### FLEXTOME™ Cutting Balloon™

| CATALOG NUMBER / ORDER NUMBER (GTIN) | LENGTH / DIAMETER | NUMBER OF ATHEROTOMES | NOMINAL PRESSURE | RATED BURST PRESSURE | MAX. RECOMMENDED GUIDEWIRE |
|--------------------------------------|-------------------|-----------------------|------------------|----------------------|----------------------------|
| H749 RB420006 0 / 08714729719861     | 6 mm / 2.00 mm    | 3                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     | .014" (0.36 mm)            |
| H749 RB422506 0 / 08714729719878     | 6 mm / 2.25 mm    | 3                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     | .014" (0.36 mm)            |
| H749 RB425006 0 / 08714729719885     | 6 mm / 2.50 mm    | 3                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     | .014" (0.36 mm)            |
| H749 RB427506 0 / 08714729719892     | 6 mm / 2.75 mm    | 3                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     | .014" (0.36 mm)            |
| H749 RB430006 0 / 08714729719908     | 6 mm / 3.00 mm    | 3                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     | .014" (0.36 mm)            |
| H749 RB432506 0 / 08714729719915     | 6 mm / 3.25 mm    | 3                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     | .014" (0.36 mm)            |
| H749 RB435006 0 / 08714729719922     | 6 mm / 3.50 mm    | 4                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     | .014" (0.36 mm)            |
| H749 RB437506 0 / 08714729719939     | 6 mm / 3.75 mm    | 4                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     | .014" (0.36 mm)            |
| H749 RB440006 0 / 08714729719946     | 6 mm / 4.00 mm    | 4                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     | .014" (0.36 mm)            |
| H749 RB420010 0 / 08714729719953     | 10 mm / 2.00 mm   | 3                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     | .014" (0.36 mm)            |
| H749 RB422510 0 / 08714729719960     | 10 mm / 2.25 mm   | 3                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     | .014" (0.36 mm)            |
| H749 RB425010 0 / 08714729719977     | 10 mm / 2.50 mm   | 3                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     | .014" (0.36 mm)            |
| H749 RB427510 0 / 08714729719984     | 10 mm / 2.75 mm   | 3                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     | .014" (0.36 mm)            |
| H749 RB430010 0 / 08714729719991     | 10 mm / 3.00 mm   | 3                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     | .014" (0.36 mm)            |
| H749 RB432510 0 / 08714729720003     | 10 mm / 3.25 mm   | 3                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     | .014" (0.36 mm)            |
| H749 RB435010 0 / 08714729720010     | 10 mm / 3.50 mm   | 4                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     | .014" (0.36 mm)            |
| H749 RB437510 0 / 08714729720027     | 10 mm / 3.75 mm   | 4                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     | .014" (0.36 mm)            |
| H749 RB440010 0 / 08714729720034     | 10 mm / 4.00 mm   | 4                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     | .014" (0.36 mm)            |
| H749 RB420015 0 / 08714729720041     | 15 mm / 2.00 mm   | 3                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     | .014" (0.36 mm)            |
| H749 RB422515 0 / 08714729720058     | 15 mm / 2.25 mm   | 3                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     | .014" (0.36 mm)            |
| H749 RB425015 0 / 08714729720065     | 15 mm / 2.50 mm   | 3                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     | .014" (0.36 mm)            |
| H749 RB427515 0 / 08714729720072     | 15 mm / 2.75 mm   | 3                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     | .014" (0.36 mm)            |
| H749 RB430015 0 / 08714729720089     | 15 mm / 3.00 mm   | 3                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     | .014" (0.36 mm)            |
| H749 RB432515 0 / 08714729720096     | 15 mm / 3.25 mm   | 3                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     | .014" (0.36 mm)            |
| H749 RB435015 0 / 08714729720102     | 15 mm / 3.50 mm   | 4                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     | .014" (0.36 mm)            |
| H749 RB437515 0 / 08714729720119     | 15 mm / 3.75 mm   | 4                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     | .014" (0.36 mm)            |
| H749 RB440015 0 / 08714729728115     | 15 mm / 4.00 mm   | 4                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     | .014" (0.36 mm)            |

COMPLIANCE CHART

| PRESSURE ATM-(KPA)         | 2.00 MM DIAMETER | 2.25 MM DIAMETER | 2.50 MM DIAMETER | 2.75 MM DIAMETER | 3.00 MM DIAMETER | 3.25 MM DIAMETER | 3.50 MM DIAMETER | 3.75 MM DIAMETER | 4.00 MM DIAMETER |
|----------------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|
| 4.0 (405)                  | 1.99             | 2.22             | 2.44             | 2.69             | 2.94             | 3.17             | 3.40             | 3.61             | 3.89             |
| 5.0 (507)                  | 2.01             | 2.24             | 2.46             | 2.72             | 2.98             | 3.21             | 3.44             | 3.67             | 3.95             |
| 6.0 (608) NOMINAL PRESSURE | 2.03             | 2.26             | 2.49             | 2.75             | 3.01             | 3.26             | 3.48             | 3.72             | 4.01             |
| 7.0 (709)                  | 2.05             | 2.29             | 2.51             | 2.77             | 3.04             | 3.31             | 3.53             | 3.76             | 4.06             |
| 8.0 (811)                  | 2.08             | 2.32             | 2.54             | 2.80             | 3.08             | 3.35             | 3.56             | 3.80             | 4.10             |
| 9.0 (91 )                  | 2.10             | 2.34             | 2.56             | 2.82             | 3.10             | 3.38             | 3.59             | 3.83             | 4.14             |
| 10.0 (1013)                | 2.11             | 2.36             | 2.58             | 2.84             | 3.12             | 3.41             | 3.62             | 3.86             | 4.17             |
| 11.0 (1115)                | 2.13             | 2.38             | 2.60             | 2.86             | 3.15             | 3.44             | 3.65             | 3.89             | 4.20             |
| 12.0 (1216) RBP            | 2.14             | 2.40             | 2.62             | 2.88             | 3.17             | 3.47             | 3.68             | 3.91             | 4.23             |



2. Cutting Balloon™

WOLVERINE™  
Cutting Balloon™ Dilatation Device

6 mm Blade Length

| CATALOG NUMBER / ORDER NUMBER (GTIN) | LENGTH / DIAMETER | NUMBER OF ATHEROTOMES | NOMINAL PRESSURE | RATED BURST PRESSURE |
|--------------------------------------|-------------------|-----------------------|------------------|----------------------|
| H749 3940306200 0 / 08714729888406   | 6 mm / 2.00 mm    | 3                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     |
| H749 3940306225 0 / 08714729888413   | 6 mm / 2.25 mm    | 3                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     |
| H749 3940306250 0 / 08714729888420   | 6 mm / 2.50 mm    | 3                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     |
| H749 3940306275 0 / 08714729888437   | 6 mm / 2.75 mm    | 3                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     |
| H749 3940306300 0 / 08714729888444   | 6 mm / 3.00 mm    | 3                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     |
| H749 3940306325 0 / 08714729888451   | 6 mm / 3.25 mm    | 3                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     |
| H749 3940306350 0 / 08714729888468   | 6 mm / 3.50 mm    | 4                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     |
| H749 3940306375 0 / 08714729888475   | 6 mm / 3.75 mm    | 4                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     |
| H749 3940306400 0 / 08714729888482   | 6 mm / 4.00 mm    | 4                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     |

10 mm Blade Length

| CATALOG NUMBER / ORDER NUMBER (GTIN) | LENGTH / DIAMETER | NUMBER OF ATHEROTOMES | NOMINAL PRESSURE | RATED BURST PRESSURE |
|--------------------------------------|-------------------|-----------------------|------------------|----------------------|
| H749 3940310200 0 / 08714729888499   | 10 mm / 2.00 mm   | 3                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     |
| H749 3940310225 0 / 08714729888505   | 10 mm / 2.25 mm   | 3                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     |
| H749 3940310250 0 / 08714729888512   | 10 mm / 2.50 mm   | 3                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     |
| H749 3940310275 0 / 08714729888529   | 10 mm / 2.75 mm   | 3                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     |
| H749 3940310300 0 / 08714729888536   | 10 mm / 3.00 mm   | 3                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     |
| H749 3940310325 0 / 08714729888543   | 10 mm / 3.25 mm   | 3                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     |
| H749 3940310350 0 / 08714729888550   | 10 mm / 3.50 mm   | 4                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     |
| H749 3940310375 0 / 08714729888567   | 10 mm / 3.75 mm   | 4                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     |
| H749 3940310400 0 / 08714729888574   | 10 mm / 4.00 mm   | 4                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     |

15 mm Blade Length

| CATALOG NUMBER / ORDER NUMBER (GTIN) | LENGTH / DIAMETER | NUMBER OF ATHEROTOMES | NOMINAL PRESSURE | RATED BURST PRESSURE |
|--------------------------------------|-------------------|-----------------------|------------------|----------------------|
| H749 3940315200 0 / 08714729888581   | 15 mm / 2.00 mm   | 3                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     |
| H749 3940315225 0 / 08714729888598   | 15 mm / 2.25 mm   | 3                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     |
| H749 3940315250 0 / 08714729888604   | 15 mm / 2.50 mm   | 3                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     |
| H749 3940315275 0 / 08714729888611   | 15 mm / 2.75 mm   | 3                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     |
| H749 3940315300 0 / 08714729888628   | 15 mm / 3.00 mm   | 3                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     |
| H749 3940315325 0 / 08714729888635   | 15 mm / 3.25 mm   | 3                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     |
| H749 3940315350 0 / 08714729888642   | 15 mm / 3.50 mm   | 4                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     |
| H749 3940315375 0 / 08714729888659   | 15 mm / 3.75 mm   | 4                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     |
| H749 3940315400 0 / 08714729888666   | 15 mm / 4.00 mm   | 4                     | 6 ATM (608kPa)   | 12 ATM (1216kPa)     |



NEW CARDIO

2

## COMPLEX PCI

Balloon Catheters

Microdilatation Catheters

Guide Catheters

Guide Catheter Extension

Angiographic Catheters

Guidewires

Embolic Protection

Accessories

Cutting Balloon™

Atherectomy Devices

Chronic Total Occlusions

TRAPPER™ Exchange Device

ANGIOJET™ Thrombectomy System



Balloon  
Catheters

Micro-  
dilatation  
Catheters

Guide  
Catheters

Guide  
Catheter  
extension

Angiographic  
Catheters

Guidewires

Embolic  
Protection

Accessories

Cutting  
Balloon™

Atherectomy  
Devices

Chronic Total  
Occlusions

Trapper™  
Exchange  
Device

ANGIOJET™  
Thrombec-  
tomy  
System

ROTAPRO™ Rotational Atherectomy System

| CATALOG NUMBER / ORDER NUMBER (GTIN) | DESCRIPTION          |
|--------------------------------------|----------------------|
| H749 3930903 0 / 08714729975151      | ROTAPRO™ Console Kit |

ROTAPRO™ Pre-Connected Burr and Advancing Device

| CATALOG NUMBER / ORDER NUMBER (GTIN) | DIAMETER |
|--------------------------------------|----------|
| H749 39300125 0 / 08714729893356     | 1.25 mm  |
| H749 39300150 0 / 08714729893363     | 1.50 mm  |
| H749 39300175 0 / 08714729893370     | 1.75 mm  |
| H749 39300200 0 / 08714729893387     | 2.0 mm   |
| H749 39300215 0 / 08714729893394     | 2.15 mm  |
| H749 39300225 0 / 08714729893400     | 2.25 mm  |
| H749 39300238 0 / 08714729893417     | 2.38 mm  |
| H749 39300250 0 / 08714729893424     | 2.5 mm   |


ROTAWIRE™

| CATALOG NUMBER / ORDER NUMBER (GTIN) | DESCRIPTION                                   | LENGTH | TIP LENGTH | MAX. DIAMETER   | QUANTITY |
|--------------------------------------|-----------------------------------------------|--------|------------|-----------------|----------|
| H802 22824002 2 / 08714729195566     | ROTAWIRE™ Floppy w / WIRECLIP™ Torquer        | 330 cm | 2.2 cm     | .014" (0.36 mm) | 5-pack   |
| H802 23239001 2 / 08714729195573     | ROTAWIRE™ Extra Support w / WIRECLIP™ Torquer | 330 cm | 2.8 cm     | .014" (0.36 mm) | 5-pack   |



ROTALINK™ Advancer

| CATALOG NUMBER / ORDER NUMBER (GTIN) | DESCRIPTION        |
|--------------------------------------|--------------------|
| H802 22782001 A0 / 08714729317739    | ROTALINK™ Advancer |




ROTALINK™ Exchangeable Burr Catheter

| CATALOG NUMBER / ORDER NUMBER (GTIN) | LENGTH | BURR SIZE |
|--------------------------------------|--------|-----------|
| H802 22768002 0 / 08714729185864     | 135 cm | 1.25 mm   |
| H802 22768003 0 / 08714729185857     | 135 cm | 1.50 mm   |
| H802 22768004 0 / 08714729185840     | 135 cm | 1.75 mm   |
| H802 22768005 0 / 08714729126881     | 135 cm | 2.00 mm   |
| H802 22768015 0 / 08714729185826     | 135 cm | 2.15 mm   |
| H802 22768006 0 / 08714729185819     | 135 cm | 2.25 mm   |
| H802 22768016 0 / 08714729185802     | 135 cm | 2.38 mm   |
| H802 22768007 0 / 08714729185796     | 135 cm | 2.50 mm   |




ROTALINK™ Plus (Advancer & Burr)

| CATALOG NUMBER / ORDER NUMBER (GTIN) | LENGTH | BURR SIZE |
|--------------------------------------|--------|-----------|
| H749 23631002 0 / 08714729316411     | 135 cm | 1.25 mm   |
| H749 23631003 0 / 08714729228363     | 135 cm | 1.50 mm   |
| H749 23631004 0 / 08714729228370     | 135 cm | 1.75 mm   |
| H749 23631005 0 / 08714729316442     | 135 cm | 2.00 mm   |
| H749 23631015 0 / 08714729316473     | 135 cm | 2.15 mm   |
| H749 23631006 0 / 08714729228394     | 135 cm | 2.25 mm   |
| H749 23631016 0 / 08714729228424     | 135 cm | 2.38 mm   |
| H749 23631007 0 / 08714729316466     | 135 cm | 2.50 mm   |



ROTABLATOR™ Consoles

| CATALOG NUMBER / ORDER NUMBER (GTIN) | MODEL                                   |
|--------------------------------------|-----------------------------------------|
| H802 22020038 1 / 08714729353294     | ROTABLATOR™ Console 220-240V (European) |
| H802 22020036 1 / 08714729353287     | ROTABLATOR™ Console 220-240V (TUV)      |




ROTAGO CART™

| CATALOG NUMBER / ORDER NUMBER (GTIN) | DESCRIPTION  |
|--------------------------------------|--------------|
| H749 3936901 0 / 08714729974307      | ROTAGO CART™ |

Accessories

| CATALOG NUMBER / ORDER NUMBER (GTIN) | DESCRIPTION                       |
|--------------------------------------|-----------------------------------|
| H802 22196003 2 / 08714729195535     | WIRECLIP™ Torquer (5-pack)        |
| H802 21600003 1 / 08714729076964     | 15" Air Supply Hose Replacement   |
| H802 22436002 1 / 08714729353317     | DynaGlide™ Foot Pedal Replacement |





NEW CARDIO

2

## COMPLEX PCI

Balloon Catheters

Microdilatation Catheters

Guide Catheters

Guide Catheter Extension

Angiographic Catheters

Guidewires

Embolic Protection

Accessories


Cutting Balloon™

Atherectomy Devices

Chronic Total Occlusions

TRAPPER™ Exchange Device

ANGIOJET™ Thrombectomy System



Balloon  
Catheters

Micro-  
dilatation  
Catheters

Guide  
Catheters

Guide  
Catheter  
extension

Angiographic  
Catheters

Guidewires

Embolic  
Protection

Accessories

Cutting  
Balloon™

Atherectomy  
Devices

Chronic Total  
Occlusions

Trapper™  
Exchange  
Device

ANGIOJET™  
Thrombec-  
tomy  
System

2. Chronic Total Occlusions

CROSSBOSS™ & STINGRAY™ LP  
Coronary Crossing and Re-Entry Devices

| CATALOG NUMBER / ORDER NUMBER (GTIN) | DESCRIPTION                           |                   |
|--------------------------------------|---------------------------------------|-------------------|
| H749 M2000B 0 / 08714729837466       | CROSSBOSS™ Coronary Crossing Catheter |                   |
| H749 393130SR 0 / 08714729903369     | STINGRAY™ LP                          | Training Required |
| H749 M3004B 0 / 08714729837473       | STINGRAY™ Guidewire, 300 cm           | Training Required |
| H749 M3012B 0 / 08714729837497       | STINGRAY™ Guidewire, 185 cm           | Training Required |
| H749 M3010B 0 / 08714729837480       | STINGRAY™ Guidewire Extension         | Training Required |

Please contact local Boston Scientific representative for training support.

Balloon  
Catheters

Micro-  
dilatation  
Catheters

Guide  
Catheters

Guide  
Catheter  
extension

Angiographic  
Catheters

Guidewires

Embolic  
Protection

Accessories

Cutting  
Balloon™

Atherectomy  
Devices

Chronic Total  
Occlusions

Trapper™  
Exchange  
Device

ANGIOJET™  
Thrombec-  
tomy  
System



NEW CARDIO

2

## COMPLEX PCI

Balloon Catheters

Microdilatation Catheters

Guide Catheters

Guide Catheter Extension

Angiographic Catheters

Guidewires

Embolic Protection

Accessories

Cutting Balloon™

Atherectomy Devices

Chronic Total Occlusions

TRAPPER™ Exchange Device

ANGIOJET™ Thrombectomy System



Balloon  
Catheters

Micro-  
dilatation  
Catheters

Guide  
Catheters

Guide  
Catheter  
extension

Angiographic  
Catheters

Guidewires

Embolic  
Protection

Accessories

Cutting  
Balloon™

Atherectomy  
Devices

Chronic Total  
Occlusions

Trapper™  
Exchange  
Device

ANGIOJET™  
Thrombec-  
tomy  
System

2. TRAPPER™ Exchange Device

TRAPPER™ Exchange Device

| CATALOG NUMBER / ORDER NUMBER (GTIN) | DESCRIPTION              | BALLOON TYPE  | BALLOON LENGTH | COMPATIBLE GUIDE CATHETER SIZE | QUANTITY |
|--------------------------------------|--------------------------|---------------|----------------|--------------------------------|----------|
| H749 3933013 0 / 05391526210154      | TRAPPER™ Exchange Device | Non-compliant | 10 mm          | 6-8F*                          | 1        |

\* 6F = 2 mm, 7F = 2.33mm, 8F = 2.67 mm

Balloon Catheters

Micro-dilatation Catheters

Guide Catheters

Guide Catheter extension

Angiographic Catheters

Guidewires

Embolic Protection

Accessories

Cutting Balloon™

Atherectomy Devices

Chronic Total Occlusions

Trapper™ Exchange Device

ANGIOJET™ Thrombectomy System



NEW CARDIO

## 2 COMPLEX PCI

Balloon Catheters

Microdilatation Catheters

Guide Catheters

Guide Catheter Extension

Angiographic Catheters

Guidewires

Embolic Protection

Accessories

Cutting Balloon™

Atherectomy Devices

Chronic Total Occlusions

TRAPPER™ Exchange Device

ANGIOJET™ Thrombectomy System



Balloon Catheters

Micro-dilatation Catheters

Guide Catheters

Guide Catheter extension

Angiographic Catheters

Guidewires

Embolic Protection

Accessories

Cutting Balloon™

Atherectomy Devices

Chronic Total Occlusions

Trapper™ Exchange Device

ANGIOJET™ Thrombectomy System

2. ANGIOJET™ Thrombectomy System

SPIROFLEX™

| CATALOG NUMBER /<br>ORDER NUMBER (GTIN) | GUIDE<br>CATHETER | CATHETER<br>DIAMETER | MINIMUM VESSEL<br>DIAMETER | DELIVERY<br>PLATFORM |
|-----------------------------------------|-------------------|----------------------|----------------------------|----------------------|
| 106553-004 / 08714729889410             | 6F                | 4F                   | 2 mm                       | RX                   |



SPIROFLEX™ VG

| CATALOG NUMBER /<br>ORDER NUMBER (GTIN) | GUIDE<br>CATHETER | CATHETER<br>DIAMETER | MINIMUM VESSEL<br>DIAMETER | DELIVERY<br>PLATFORM |
|-----------------------------------------|-------------------|----------------------|----------------------------|----------------------|
| 106608-004 / 08714729889472             | 7F                | 5F                   | 3 mm                       | RX                   |



Balloon  
Catheters

Micro-  
dilatation  
Catheters

Guide  
Catheters

Guide  
Catheter  
extension

Angiographic  
Catheters

Guidewires

Embolic  
Protection

Accessories

Cutting  
Balloon™

Atherectomy  
Devices

Chronic Total  
Occlusions

Trapper™  
Exchange  
Device

ANGIOJET™  
Thrombec-  
tomy  
System





NEW CARDIO

## 3 PCI GUIDANCE

POLARIS™ PCI Guidance System\*

POLARIS™ PCI Guidance Accessories

Imaging

Physiology

POLARIS™  
PCI Guidance  
System\*

POLARIS™  
PCI Guidance  
Accessories

Imaging

Physiology



\* POLARIS™ 1.0 software running on iLab™ System 3.0

NEW CARDIO

3

## PCI GUIDANCE

POLARIS™ PCI Guidance System\*

POLARIS™ PCI Guidance Accessories

Imaging

Physiology

POLARIS™  
PCI Guidance  
System\*

POLARIS™  
PCI Guidance  
Accessories

Imaging

Physiology



\* POLARIS™ 1.0 software running on iLab™ System 3.0

3. POLARIS™ PCI Guidance System\*

POLARIS™ Cart System\*\*

| CATALOG NUMBER / ORDER NUMBER (GTIN)  | DESCRIPTION                        |
|---------------------------------------|------------------------------------|
| BUN30 POLARIS CART22 / 08714729909736 | POLARIS™ CART 220V 2.7 & MDU 5+*** |
| BUN30 POLARIS CART24 / 08714729909743 | POLARIS™ CART 240V 2.7 & MDU 5+*** |



POLARIS™  
PCI Guidance  
System\*

POLARIS™ Installed System\*\*

| CATALOG NUMBER / ORDER NUMBER (GTIN)  | DESCRIPTION                           |
|---------------------------------------|---------------------------------------|
| BUN30 POLARIS INST22 / 08714729909750 | POLARIS™ INSTALL 220V 2.7 & MDU 5+*** |
| BUN30 POLARIS INST24 / 08714729909767 | POLARIS™ INSTALL 240V 2.7 & MDU 5+*** |



POLARIS™  
PCI Guidance  
Accessories

Imaging

Physiology



\* POLARIS™ 1.0 software running on iLab™ System 3.0  
\*\* Please contact your Boston Scientific sales representative.  
POLARIS™ Installed will require additional special devices for installation  
requiring a visit from a Boston Scientific Field Service Engineer.  
\*\*\* Please contact your Boston Scientific representative  
to choose appropriate voltage for your country.

NEW CARDIO

3

## PCI GUIDANCE

POLARIS™ PCI Guidance System\*

POLARIS™ PCI Guidance Accessories

Imaging

Physiology

POLARIS™  
PCI Guidance  
System\*

POLARIS™  
PCI Guidance  
Accessories

Imaging

Physiology



\* POLARIS™ 1.0 software running on iLab™ System 3.0

3. POLARIS™ PCI Guidance Accessories

| POLARIS™ Accessories                 |                              |
|--------------------------------------|------------------------------|
| CATALOG NUMBER / ORDER NUMBER (GTIN) | DESCRIPTION                  |
| H749 A7020 0 / 08714729228684        | Sterile Pullback Sled        |
| H749 A7026 0 / 08714729400578        | Sterile Motordrive Bag       |
| H749 00190 0 / 08714729842262        | Mouse                        |
| H749 00191 0 / 08714729838609        | Mouse pad (cart)             |
| H749 00192 0 / 08714729838616        | Mouse pad (installed system) |
| H749 00193 0 / 08714729948803        | POLARIS™ * Cart Static Deck  |
| H749 MDU5PLUSBAGO / 08714729842255   | Extra MDU 5+ Sterile Bag     |
| H749 MDU5PLUS / 08714729842279       | Motordrive MDU5+             |

POLARIS™  
PCI Guidance  
System\*

POLARIS™  
PCI Guidance  
Accessories

Imaging

Physiology



NEW CARDIO

# 3 PCI GUIDANCE

POLARIS™ PCI Guidance System\*

POLARIS™ PCI Guidance Accessories

Imaging

Physiology

POLARIS™  
PCI Guidance  
System\*

POLARIS™  
PCI Guidance  
Accessories

Imaging

Physiology



\* POLARIS™ 1.0 software running on iLab™ System 3.0



OPTICROSS™ Coronary Imaging Catheter - IVUS

| CATALOG NUMBER / ORDER NUMBER (GTIN) | DESCRIPTION                                  | USABLE LENGTH | GUIDE CATHETER COMPATIBILITY |
|--------------------------------------|----------------------------------------------|---------------|------------------------------|
| H749 51812 0 / 08714729841531        | OPTICROSS™ 40MHz Coronary Imaging Catheter   | 135 cm        | 5F (1.6 mm)                  |
| H749 518126 0 / 08714729904779       | OPTICROSS™ 6 40MHz Coronary Imaging Catheter | 135 cm        | 6F (2 mm)                    |

OPTICROSS™ High-definition Coronary Imaging Catheter - HD IVUS

| CATALOG NUMBER / ORDER NUMBER (GTIN) | DESCRIPTION | USABLE LENGTH | GUIDE CATHETER COMPATIBILITY |
|--------------------------------------|-------------|---------------|------------------------------|
| H749 3935202 0 / 08714729960751      | HD IVUS CE  | 135 cm        | 5F (1.6 mm)                  |

POLARIS™  
PCI Guidance  
System\*

POLARIS™  
PCI Guidance  
Accessories

Imaging

Physiology



NEW CARDIO

3

## PCI GUIDANCE

POLARIS™ PCI Guidance System\*

POLARIS™ PCI Guidance Accessories

Imaging

Physiology

POLARIS™  
PCI Guidance  
System\*

POLARIS™  
PCI Guidance  
Accessories

Imaging

Physiology



\* POLARIS™ 1.0 software running on iLab™ System 3.0

COMET™ Physiological pressure guidewire -  
FFR, DFR™, Pd/Pa

| CATALOG NUMBER / ORDER NUMBER (GTIN) | DESCRIPTION               |
|--------------------------------------|---------------------------|
| H749 393243 0 / 08714729904397       | COMET™ Pressure Guidewire |
| H749 555100 0 / 08714729890010       | FFRLink™                  |



POLARIS™  
PCI Guidance  
System\*

POLARIS™  
PCI Guidance  
Accessories

Imaging

Physiology



NEW CARDIO

## 4 STRUCTURAL HEART

Transcatheter Aortic Valve Implantation (TAVI)

Left Atrial Appendage Closure (LAAC)

TAVI


LAAC



## 4 STRUCTURAL HEART

Transcatheter Aortic Valve Implantation (TAVI)

Left Atrial Appendage Closure (LAAC)



TAVI

LAAC



## ACURATE neo™ Aortic Bioprosthesis

Compatible with ACURATE neo TF™ Transfemoral Delivery System & ACURATE neo™ Transapical Delivery System


| CATALOG NUMBER / ORDER NUMBER (GTIN) | PRODUCT DESCRIPTION                 | AORTIC ANNULUS PERIMETER | AORTIC ANNULUS DIAMETER | UNITS |
|--------------------------------------|-------------------------------------|--------------------------|-------------------------|-------|
| SYM-SV23-002 / 07640168110062        | ACURATE neo™ Aortic Bioprosthesis S | S 66 - 72 mm             | S 21 - 23 mm            | 1     |
| SYM-SV25-002 / 07640168110079        | ACURATE neo™ Aortic Bioprosthesis M | M 72 - 79 mm             | M 23 - 25 mm            | 1     |
| SYM-SV27-002 / 07640168110086        | ACURATE neo™ Aortic Bioprosthesis L | L 79 - 85 mm             | L 25 - 27 mm            | 1     |



TAVI

## ACURATE neo TF™ Transfemoral Delivery System

| CATALOG NUMBER / ORDER NUMBER (GTIN) | PRODUCT DESCRIPTION                          | UNITS |
|--------------------------------------|----------------------------------------------|-------|
| SYM-DS-002 / 07640168110055          | ACURATE neo TF™ Transfemoral Delivery System | 1     |



## ACURATE neo™ TA Transapical Delivery System

| CATALOG NUMBER / ORDER NUMBER (GTIN) | PRODUCT DESCRIPTION                         | UNITS |
|--------------------------------------|---------------------------------------------|-------|
| SYM-DS-004 / 07640168110093          | ACURATE neo™ TA Transapical Delivery System | 1     |




LAAC





LOTUS Edge™ Aortic Valve System

| CATALOG NUMBER / ORDER NUMBER (GTIN) | INTRODUCER SHEATH SIZE |
|--------------------------------------|------------------------|
| H749LVS230                           | LOTUS Edge™ 23 mm      |
| H749LVS250                           | LOTUS Edge™ 25 mm      |
| H749LVS270                           | LOTUS Edge™ 27 mm      |



LOTUS™ Introducer Set\*

| CATALOG NUMBER / ORDER NUMBER (GTIN) | INTRODUCER SHEATH SIZE |
|--------------------------------------|------------------------|
| H749 NTR180 / 05391526210109         | Small (LIS-S)          |
| H749 NTR200 / 05391526210116         | Large (LIS-L)          |

\* The LOTUS™ Introducer Set is manufactured by Creganna-Tactx Medical and distributed by Boston Scientific Corporation.

SENTINEL™ Cerebral Protection System

| CATALOG NUMBER / ORDER NUMBER (GTIN) | PRODUCT DESCRIPTION                  | UNITS |
|--------------------------------------|--------------------------------------|-------|
| CMS15-10C / 00863229000028           | SENTINEL™ Cerebral Protection System | 1     |

TAVI

LAAC



14F iSLEEVE™ Expandable Introducer Set

Included with the ACURATE neo™ Aortic Valve System

| CATALOG NUMBER / ORDER NUMBER (GTIN) | PRODUCT DESCRIPTION                          |
|--------------------------------------|----------------------------------------------|
| BUN67ACURDSTFISL                     | ACURATE neo™ + 14 F iSLEEVE,™ Generic Bundle |
| BUN67ACURTFISL14FS*                  | ACURATE neo™ S + 14 F iSLEEVE™               |
| BUN67ACUTFISL14FM*                   | ACURATE neo™ M + 14 F iSLEEVE™               |
| BUN67ACUTFISL14FL*                   | ACURATE neo™ L + 14 F iSLEEVE™               |

SAFARI2™ Pre-Shaped Guidewire\*

Pre-Shaped TAVI Guidewire

| CATALOG NUMBER / ORDER NUMBER (GTIN) | OUTER            | DIAMETER LENGTH | CURVE CONFIGURATION | CORE MATERIAL   | COATING          | QUANTITY |
|--------------------------------------|------------------|-----------------|---------------------|-----------------|------------------|----------|
| H749 39406XS1 / 08714729887614       | .035" / 0.889 mm | 275 cm          | Extra Small         | Stainless Steel | LUBRIGREEN™ PTFE | 5-pack   |
| H749 39406S1 / 08714729887591        | .035" / 0.889 mm | 275 cm          | Small               | Stainless Steel | LUBRIGREEN™ PTFE | 5-pack   |
| H749 39406L1 / 08714729887577        | .035" / 0.889 mm | 275 cm          | Large               | Stainless Steel | LUBRIGREEN™ PTFE | 5-pack   |
| H749 39407XS0 / 08714729887638       | .035" / 0.889 mm | 275 cm          | Extra Small         | Stainless Steel | LUBRIGREEN™ PTFE | 1-single |
| H749 39407S0 / 08714729887652        | .035" / 0.889 mm | 275 cm          | Small               | Stainless Steel | LUBRIGREEN™ PTFE | 1-single |
| H749 39407L0 / 08714729887645        | .035" / 0.889 mm | 275 cm          | Large               | Stainless Steel | LUBRIGREEN™ PTFE | 1-single |

\* The SAFARI2™ Guidewire is manufactured by Lake Region Medical and distributed by Boston Scientific Corporation.

TAVI

LAAC




# 4 STRUCTURAL HEART

Transcatheter Aortic Valve Implantation (TAVI)

Left Atrial Appendage Closure (LAAC)

TAVI

LAAC



## WATCHMAN FLX™ LAA Closure Device and Delivery System

| CATALOG NUMBER / ORDER NUMBER (GTIN) | SIZE  | OD  | LENGTH  | QUANTITY |
|--------------------------------------|-------|-----|---------|----------|
| M635WS50200 / 08714729860433         | 20 mm | 12F | 75.5 cm | 1        |
| M635WS50240 / 08714729860440         | 24 mm | 12F | 75.5 cm | 1        |
| M635WS50270 / 08714729860457         | 27 mm | 12F | 75.5 cm | 1        |
| M635WS50310 / 08714729860464         | 31 mm | 12F | 75.5 cm | 1        |
| M635WS50350 / 08714729860471         | 35 mm | 12F | 75.5 cm | 1        |

## WATCHMAN™ TruSeal™ Access System

| CATALOG NUMBER / ORDER NUMBER (GTIN) | PRODUCT DESCRIPTION                             | I.D.         | O.D.         |
|--------------------------------------|-------------------------------------------------|--------------|--------------|
| M635TS70010 / 08714729965732         | WATCHMAN™ TruSeal™ Access System Single Curve   | 12F (4.2 mm) | 14F (4.8 mm) |
| M635TS70020 / 08714729965749         | WATCHMAN™ TruSeal™ Access System Double Curve   | 12F (4.2 mm) | 14F (4.8 mm) |
| M635TS70040 / 08714729965756         | WATCHMAN™ TruSeal™ Access System Anterior Curve | 12F (4.2 mm) | 14F (4.8 mm) |

## WATCHMAN™ LAA Closure Device and Delivery System

| CATALOG NUMBER / ORDER NUMBER (GTIN) | SIZE  | OD  | LENGTH  | QUANTITY |
|--------------------------------------|-------|-----|---------|----------|
| M635 WS21060 / 08714729797685        | 21 mm | 12F | 75.5 cm | 1        |
| M635 WS24060 / 08714729797692        | 24 mm | 12F | 75.5 cm | 1        |
| M635 WS27060 / 08714729797708        | 27 mm | 12F | 75.5 cm | 1        |
| M635 WS30060 / 08714729797715        | 30 mm | 12F | 75.5 cm | 1        |
| M635 WS33060 / 08714729797722        | 33 mm | 12F | 75.5 cm | 1        |

## WATCHMAN™ Transseptal Kits\*

WATCHMAN™ LAA Closure Device & Delivery System, WATCHMAN™ Access Sheath, Heart Span Sheath and Heart Span Needle

| CATALOG NUMBER / ORDER NUMBER (GTIN) | SIZE  | OD  | LENGTH  | QUANTITY |
|--------------------------------------|-------|-----|---------|----------|
| TNS HWS2106 0 / 08714729915300       | 21 mm | 12F | 75.5 cm | 1        |
| TNS HWS2406 0 / 08714729915317       | 24 mm | 12F | 75.5 cm | 1        |
| TNS HWS2706 0 / 08714729915324       | 27 mm | 12F | 75.5 cm | 1        |
| TNS HWS3006 0 / 08714729915331       | 30 mm | 12F | 75.5 cm | 1        |
| TNS HWS3306 0 / 08714729915348       | 33 mm | 12F | 75.5 cm | 1        |

\*Check availability in your country with your BSC sales representative.

