

JUDO™ Guidewires

Part of the SENTAI™ Guidewire Family

JUDO Guidewires are designed to tackle your toughest cases with the softest touch. Their combination of excellent penetration force, unmatched flexibility, and control make them the ultimate guidewires for intraluminal crossing in challenging lesions.

Micro EMT™ Technology

A novel guidewire technology that exactly matches the taper of the core and coil, enabling outstanding lesion **Entry**, intraluminal **Maneuverability**, and controlled **Torque**.

15 cm Hydrophilic

0.008"
(0.20 mm)

6 cm Taper

JUDO 1

Soft intraluminal crossing wire for antegrade microchannels

1 gf

JUDO 3

Intermediate intraluminal crossing wire for fibro-calcific lesions

3 gf

JUDO 6

Extra penetration with excellent steerability in tight lesions

6 gf

Maneuverability

A flexible distal shaft enables excellent intraluminal navigation in complex lesions, with minimal perforation risk

Entry

The smallest tip on the market 0.008" (0.20 mm) and the long taper (6 cm) facilitate easy lesion entry in even the most complex lesions

Torque

With exactly matched tapers, the core and coil move as one, eliminating whip and empowering precise steerability

**MASTER THE
COMPLEX™**

Product Specifications

	Tip Load (gf)	Penetration Force (gf/mm ²)	Coating	Tip Radiopacity (cm)	Clinical Application
JUDO 1 Guidewire	1.0	31	Hydrophilic	3.5	Soft intraluminal crossing wire for antegrade microchannels
JUDO 3 Guidewire	3.0	93	Hydrophilic	3.5	Intermediate intraluminal crossing wire for fibro-calcific lesions
JUDO 6 Guidewire	6.0	185	Hydrophilic	3.5	Extra penetration with excellent steerability in tight lesions

Ordering Information

	Order Number (GTIN)	Ref/Catalog Number	Length (cm)
JUDO 1 Guidewire	08714729975038 08714729975052	H749 39415190 2 H749 39415300 2	190 300
JUDO 3 Guidewire	08714729974956 08714729974970	H749 39377190 2 H749 39377300 2	190 300
JUDO 6 Guidewire	08714729974994 08714729975014	H749 39378190 2 H749 39378300 2	190 300

All JUDO Guidewires are straight-tip and come in quantities of 5.

**Boston
Scientific**
Advancing science for life™

All trademarks are the property of their respective owners. Prior to use, please review full device DFU for Indications, Contraindications, Warnings, Precautions, Adverse Events, and Operator's Instructions.

CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device. Information for use only in countries with applicable health authority registrations. Material not intended for use in France.
IC-581409-AA-MAY2019 - CreativeServices

www.bostonscientific-international.com

© 2019 Boston Scientific Corporation
or its affiliates. All rights reserved.
DINCAR2750EA