

INLIVEN™ – The *Physiological* Solution.
CRT-P

**Boston
Scientific**
Advancing science for life™

*Follows
the breath*



INLIVEN™ is the first¹ CRT-P relying on a physiological parameter - **respiration** - to:



Optimize pacing therapy thanks to RightRate™, based on Minute Ventilation Technology¹.

ACC/AHA/ESC Guidelines recommend exercise training of HF patient as part of the treatment, reducing mortality and improving Quality of Life²



Provide you with valuable screening and diagnostic information on Sleep Apnea and Respiratory Rate, for risk stratification and co-morbidity management.¹

ACC/AHA/ESC Guidelines for the Diagnosis and Treatment of Acute and Chronic Heart Failure recommend the measurement and documentation of physiological vital signs including Respiratory Rate for cardiac patients³

INLIVEN™ is designed to:



Provide your patient with optimal CRT therapy and physiological pacing thanks to RightRate™ technology.

Maximize efficiency and ease of use with auto-thresholds and user-friendly interface to provide you with **assurance, first**.

¹. To date March 2013 - 359896-020 Inliven PTM EN Europe PDF 1.0, MEDTRONIC CONSULTA CRT-P C3TR01 Clinical manual M927876A001C, BIOTRONIK EVIA HF-T ProMRI Manual Version 387513-C.

ST_JUDE MEDICAL ANTHEM Anthem_CRT_P_RF_3212_SpecSheet_Intl Ref. GMC RM459

². ACC/AHA 2005 Guideline Update for the Diagnosis and Management of Chronic Heart Failure in the Adult. A Report of the American College of Cardiology/American Heart Association Task Force on Practice Guidelines.

SA Hunt et al. Circulation. 2005;112:e154-e235 - Guidelines for the diagnosis and treatment of Chronic Heart Failure. The taskforce for the diagnosis and treatment of CHF of the ESC. K. Swedberg et al., EHJ 2005

³. ACC/AHA Heart Failure Clinical Data Standards. Circulation, Vol. 112 (12), September 20, 2005 - Dickstein et al. European Heart Journal 2008, doi: 10.1093/eurheart/ehn309 - Mont et al. Daily maximum and median respiratory rate becomes elevated several weeks before heart failure hospitalization in patients with implantable devices. Poster Europace 2009.

All cited trademarks are the property of their respective owners. CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labelling supplied with each device. Information for the use only in countries with applicable health authority product registrations.

**Boston
Scientific**

Advancing science for life™

www.bostonscientific-international.com

© 2013 Boston Scientific Corporation
or its affiliates. All rights reserved.
DINCRM0786EA