

LATITUDE™ NXT

PATIENT MANAGEMENT SYSTEM



LATITUDE Wireless Communicator
Models 6280 and 6290

Model 6290 Communicator

CRT-D

G124 MOMENTUM™ CRT-D
G125 MOMENTUM CRT-D
G126 MOMENTUM CRT-D
G128 MOMENTUM X4 CRT-D
G138 MOMENTUM X4 CRT-D
G140 INOGEN™ CRT-D
G141 INOGEN CRT-D
G146 INOGEN X4 CRT-D
G148 INOGEN X4 CRT-D
G150 DYNAGEN™ CRT-D
G151 DYNAGEN CRT-D
G154 DYNAGEN CRT-D
G156 DYNAGEN X4 CRT-D
G158 DYNAGEN X4 CRT-D
G247 VIGILANT™ X4 CRT-D
G248 VIGILANT X4 CRT-D
G424 RESONATE™ CRT-D

G437 RESONATE X4 CRT-D
G447 RESONATE X4 CRT-D
G448 RESONATE X4 CRT-D
G525 RESONATE HF CRT-D
G526 RESONATE HF CRT-D
G547 RESONATE HF CRT-D
G548 RESONATE HF CRT-D

S-ICD

A209 EMBLEM™ S-ICD A219
EMBLEM MRI S-ICD

ICD

D010 INOGEN MINI ICD VR
D011 INOGEN MINI ICD VR
D012 INOGEN MINI ICD DR
D013 INOGEN MINI ICD DR
D020 DYNAGEN MINI ICD VR
D021 DYNAGEN MINI ICD VR
D022 DYNAGEN MINI ICD DR
D023 DYNAGEN MINI ICD DR
D120 MOMENTUM EL ICD VR

D121 MOMENTUM EL ICD DR
D140 INOGEN EL ICD VR
D141 INOGEN EL ICD VR
D142 INOGEN EL ICD DR
D143 INOGEN EL ICD DR
D150 DYNAGEN EL ICD VR
D151 DYNAGEN EL ICD VR
D152 DYNAGEN EL ICD DR
D153 DYNAGEN EL ICD DR
D220 VIGILANT EL ICD VR
D221 VIGILANT EL ICD DR
D232 VIGILANT EL ICD VR
D233 VIGILANT EL ICD DR
D400 PERCIVA™ ICD VR
D401 PERCIVA ICD DR
D412 PERCIVA ICD VR
D413 PERCIVA ICD DR
D500 PERCIVA HF ICD VR
D501 PERCIVA HF ICD DR
D512 PERCIVA HF ICD DR
D513 PERCIVA HF ICD DR
D421 RESONATE EL ICD VR
D432 RESONATE EL ICD VR
D433 RESONATE EL ICD DR
D520 RESONATE HF ICD VR
D521 RESONATE HF ICD DR
D532 RESONATE HF ICD VR
D533 RESONATE HF ICD DR

Pacemakers

L100 ESSENTIO™ SR
L101 ESSENTIO DR
L110 ESSENTIO SR
L111 ESSENTIO DR
L121 ESSENTIO DR EL
L131 ESSENTIO DR EL
L300 ACCOLADE™ SR
L301 ACCOLADE DR
L310 ACCOLADE SR
L311 ACCOLADE DR
L321 ACCOLADE DR EL
L331 ACCOLADE DR EL

CRT-P

U125 VALITUDE™ CRT-P
U128 VALITUDE X4 CRT-P
U225 VISIONIST™ CRT-P
U226 VISIONIST CRT-P
U228 VISIONIST X4 CRT-P



Available for in-clinic or direct-to-patient distribution via eCommerce.

Model 6280 Communicator

CRT-D

N050 PUNCTUA™ CRT-D
N051 PUNCTUA CRT-D
N118 COGNIS™ CRT-D
N119 COGNIS CRT-D
N140 ENERGEN™ CRT-D VR
N141 ENERGEN CRT-D
N160 INCEPTA™ CRT-D
N161 INCEPTA CRT-D
N164 INCEPTA CRT-D

ICD

E050 PUNCTUA ICD VR
E051 PUNCTUA ICD VR
E052 PUNCTUA ICD DR
E053 PUNCTUA ICD DR
E102 TELIGEN™ ICD VR
E110 TELIGEN ICD DR
E140 ENERGEN ICD VR
E141 ENERGEN ICD VR
E142 ENERGEN ICD DR
E143 ENERGEN ICD DR
E160 INCEPTA ICD VR
E161 INCEPTA ICD VR
E162 INCEPTA ICD DR
E163 INCEPTA ICD DR

Pacemakers

K062 ADVANTIO™ SR
K063 ADVANTIO DR
K064 ADVANTIO DR EL
K172 INGENIOTM SR
K173 INGENIO DR
K174 INGENIO DR EL
K272 VITALIO™ SR
K273 VITALIO DR
K274 VITALIO DR EL

CRT-P

V172 INVIVE™ CRT-P
V173 INVIVE CRT-P
V272 INTUA™ CRT-P
V273 INTUA CRT-P



Available for direct-to-patient distribution via eCommerce.

Communicator Connectivity Options

USB Cellular Adapters Ethernet Adapters

Refer to the product labeling for specific instructions for use. Rx only.



LATITUDE™ NXT Patient Management
System Indications, Safety and Warnings

<https://qrco.de/latitudepatient>