

Piranha™

Ureteroscopic Biopsy Forceps



Tissue collection made simple.

The Piranha Ureteroscopic Biopsy Forceps offers a simple and reliable way to collect tissue specimens from the ureters and kidneys for histological examinations and diagnosis. The single-use Piranha Ureteroscopic Biopsy Forceps delivers.

Simplified approach

This deflectable cup biopsy device can be back-loaded through a flexible ureteroscope. There's no need to remove the scope during insertion of the forceps, helping to streamline endoscopic tissue collection procedures.

Ease of insertion

With an outer diameter of 3 Fr (1 mm), the Piranha Biopsy Forceps is designed to pass smoothly through the working channel of most flexible ureteroscopes – even in tortuous anatomies.

Reliable collection

The Piranha Biopsy Forceps has a long working length of 115 cm that enables access to the kidney and ureter. An ergonomic handle gives you the power to easily open and close the large forceps jaw, so you can quickly collect tissue specimens for analysis.



Similar to other commonly used devices, you can easily backload the Piranha Biopsy Forceps into our LithoVue™ Flexible Ureteroscope's working channel.

Piranha™ Ureteroscopic Biopsy Forceps

A large jaw capacity for tissue collection.



The jaw's 3 Fr (1 mm) outer diameter passes through the working channel of most flexible ureterscopes.



Backload through a flexible ureterscope without interrupting the procedure.



Single-use device supports your infection-control protocols.

Scissor-style handle makes it easy to open and close forceps jaws.



Working length of 115 cm enables access to the kidney and ureter.

Ordering Information

Product Code	Description	Units
M0065051600	Piranha Ureteroscopic Biopsy Forceps	Each

Please consult your sales representative for more information and ordering details.

Product available in the European Economic Area (EEA) only. Please check availability with your local sales representative or customer service. All trademarks are the property of their respective owners. All photographs taken by Boston Scientific.

CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labelling supplied with each device. Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. This material is not intended for use in France.

URO-542517-AA Printed in Germany by medicalvision.

**Boston
Scientific**
Advancing science for life™

www.bostonscientific.eu

© 2021 Boston Scientific Corporation or its affiliates. All rights reserved.