



Tegus Orca

Remote connectivity that's ready when you are

Digital access to any OR or cathlab in the world

that enables real-time communication, guidance and support between an operator in a healthcare setting and a remote viewer or viewers.



What does it require?



Which services does it support?



What are its advantages?



Scan me to visit the
ExpertLink EU page on
[bostonscientific.com](https://www.bostonscientific.com)




Tegus Orca

Remote connectivity that's ready when you are



What does it require?

- ▶ **In the OR:** Plug-and-play remote streaming technology so that when ready to use the tegus camera system, simply turn it on and connect the internet cable.
- The tegus team will work with your hospital IT department to arrange a secure network connection in the desired clinical areas.
- ▶ **Remote:** Web-based platform on desktop PC.



View of the proctor

Scan me to visit the
ExpertLink EU page on
[bostonscientific.com](https://www.bostonscientific.com)



Tegus Orca

Remote connectivity that's ready when you are



Which services does it support?

- ▶ Case support.
- ▶ Remote proctoring



Scan me to visit the
ExpertLink EU page on
[bostonscientific.com](https://www.bostonscientific.com)




Tegus Orca

Remote connectivity that's ready when you are



What are its advantages?

- ▶ Simple setup, ready in minutes - plug-and-play remote streaming technology for on-demand assistance.
- ▶ 10x optical zoom for optimized view of the procedure.
- ▶ Crystal-clear communication for two operators.
- ▶ Easily moveable within departments for enhanced flexibility.



Scan me to visit the
ExpertLink EU page on
[bostonscientific.com](https://www.bostonscientific.com)