




Endoscopy

**Boston
Scientific**
Advancing science for life™



[Future Product Updates](#)

[Product Guide](#)

[Services](#)

2017
Reference Guide

Interactive Product Guide

 Active Cords..... 14, 77	 Dreamwire™ High Performance Guidewire..... 3
 Acquire™ Endoscopic Ultrasound Fine Needle Biopsy (FNB) Device.. 44	 Dynamic™ (Y) Stent..... 79
 Advanix™ Center Bend - Single Stent..... 21	 Electrosurgical Bipolar Cable Adapter..... 55
 Advanix Center Bend Preloaded RX Stent Systems..... 20	 Encore™ 26 Inflation Device..... 36
 Advanix Double Pigtail - Single Stent..... 22	 Endostat™ Generator Accessories..... 54
 Advanix Duodenal Bend - Single Stent..... 21	 EndoVive™ Low Profile Replacement Button..... 71
 Advanix Duodenal Bend Preloaded RX Stent Systems..... 20	 EndoVive One-Step Button™ and Low Profile Replacement Button Accessories..... 68, 71
 Advanix Pancreatic Stent - Straight, No Leading Barb..... 27	 EndoVive One-Step Button Low Profile PEG Kits..... 68
 Advanix Pancreatic Stent - Straight, Leading Barb..... 28	 EndoVive Safety PEG Kit / Safety PEG Kit with ENFit™..... 67
 Advanix Pancreatic Stent - Pigtail, No Leading Barb..... 29	 EndoVive Replacement Gastrostomy Feeding Tube with ENFit, Straight and Right Angle Bolster..... 69
 Advanix Pancreatic Stent - Pigtail, Leading Barb..... 30	 EndoVive Standard Balloon Replacement Kit, Right Angle Bolster... 70
 Advanix Pancreatic Stent Kit - Pigtail, Leading Barb and No Leading Barb..... 26	 EndoVive Standard Balloon Replacement Kit, Straight Bolster..... 69
 Advanix Pancreatic Stent Kit - Straight, No Leading Barb..... 26	 EndoVive Standard PEG Kit / Standard PEG Kit with ENFit..... 67
 Alair™ Bronchial Thermoplasty System..... 79	 EndoVive Through-the-PEG (TTP) Jejunal Feeding Tube..... 70
 Alliance™ II Integrated Inflation / Lithotripsy Device..... 61, 85	 ENFit Feeding Port Transition Connector..... 70
 Alliance II Single-Use Syringe / Gauge Assembly..... 61, 85	 eXcelon™ Transbronchial Aspiration Needle..... 83
 Amplatz Super Stiff™ Guidewire..... 49	 Expect™ Endoscopic Ultrasound Aspiration Needle..... 43
 Autotome™ RX Cannulating Sphincterotome..... 8	 Expect Pulmonary Endobronchial Ultrasound Transbronchial Aspiration Needle / Adaptor..... 83
 AXIOS™ Stent and Electrocautery Enhanced Delivery System..... 44	 Expect Slimline (SL) Endoscopic Ultrasound Aspiration Needle..... 43
 AXIOS Stent and Delivery System..... 44	 Extractor™ Pro DL Retrieval Balloon Catheter..... 33
 Biliary Introducer..... 24	 Extractor Pro RX Retrieval Balloon Catheter – Distal Injection..... 32
 C-Flex™ Double Pigtail Biliary Stent..... 25	 Extractor Pro RX Retrieval Balloon Catheter – Proximal Injection..... 32
 Captiflex™ Single-Use Snares..... 76	 Extractor Pro RX-S Retrieval Balloon Catheter – Distal Injection..... 32
 Captivator™ II Single-Use Snares..... 73	 Extractor Pro RX-S Retrieval Balloon Catheter – Proximal Injection.. 32
 Captivator Single-Use Snares..... 74	 Extractor Pro XL Retrieval Balloon Catheter – Distal Injection..... 33
 Captivator EMR Device..... 73	 Extractor Pro XL Retrieval Balloon Catheter – Proximal Injection.. 33
 Celebri™ Single-Use Endoscopic Cytology Brush..... 84	 Flexima™ Duodenal Bend Biliary Stent with Delivery System..... 23
 Contour™ ERCP Cannulas..... 13	 Flexima Single Pigtail Nasobiliary Catheter..... 25
 CRE™ Fixed Wire Esophageal Balloon Dilatation Catheters..... 59	 Fluoro Tip™ ERCP Cannulas..... 13
 CRE PRO Wireguided Biliary Balloon Dilatation Catheters..... 34	 Gold Probe™ Single-Use Electrohemostasis Catheters..... 52
 CRE PRO Wireguided Esophageal / Pyloric / Biliary Balloon Dilatation Catheters..... 57	 Hurricane™ RX Biliary Balloon Dilatation Catheter..... 35
 CRE PRO Wireguided Esophageal / Pyloric / Colonic / Biliary Balloon Dilatation Catheters..... 57	 Hydra Jagwire™ High Performance Guidewire..... 3
 CRE SteriFlate Disposable Inflation Device..... 61	 Hydratome™ RX Cannulating Sphincterotome..... 6
 CRE Wireguided Esophageal / Pyloric / Balloon Dilatation Catheters..... 58	 Inflation / Deflation Devices and Accessories..... 61
 CRE Wireguided Esophageal / Pyloric / Colonic / Balloon Dilatation Catheters..... 58	 Initial Placement PEG Accessories..... 68
 CRE Pulmonary Balloon Dilatation Catheters..... 85	 Injection Gold Probe™ Single-Use Electrohemostasis / Injection Catheters..... 52
 Distal End ENFit Transition Connector..... 70	 Interject™ Clear Single-Use Injection Therapy Needle Catheters.. 53
 Dreamtome™ RX Cannulating Sphincterotome..... 6	 Interject Contrast Single-Use Injection Therapy Needle Catheters.. 53
	 Jagtail™ High Performance Guidewire Extension..... 4



*Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

Note: RX devices are short wire devices.

 [Terms & Conditions](#)



[Section Table of Contents](#)

Prev Page

1

Next Page

[Main Table of Contents](#)




Interactive Product Guide

 Jagtome™ RX Cannulating Sphincterotome..... 7	 SpyGlass DS Direct Visualization System – Capital Equipment 39
 Jagwire™ High Performance Guidewire..... 4	 Stonetome™ Stone Removal Device..... 11
 MaxForce TTS™ Esophageal High Performance Balloon Dilatation Catheter..... 60	 Tandem™ XL Triple-Lumen ERCP Cannula 13
 MaxForce™ Biliary Balloon Dilatation Catheter 36	 Trapezoid™ RX Wireguided Retrieval Baskets..... 33
 MicroKnife™ XL Triple-Lumen Needle Knife..... 10	 TRUEtome™ Jag Cannulating Sphincterotome 9
 Multibite™ Single-Use Multiple Sample Biopsy Forceps..... 65	 TRUEtome Hydra Cannulating Sphincterotome..... 9
 NaviFlex™ RX Delivery System 22	 TRUEtome Dreamwire™ Cannulating Sphincterotome 9
 NaviFlex RX Pancreatic Delivery System & NaviFlex Pushers 31	 TRUEtome Cannulating Sphincterotome 10
 NaviPro™ Hydrophilic Guidewire 5	 Ultraflex™ Single-Use Covered Large Esophageal NG Stent System – Distal Release 49
 NovaGold™ High Performance Guidewire..... 5	 Ultraflex Single-Use Covered Large Esophageal NG Stent System – Proximal Release..... 49
 Optiflo™ Single-Use Hemostasis Catheters 53	 Ultraflex Single-Use Uncovered Esophageal NG Stent System – Distal Release 49
 Pathfinder™ Exchange Guidewire 5	 Ultraflex Single-Use Uncovered Esophageal NG Stent System – Proximal Release..... 49
 Percuflex™ Duodenal Bend Biliary Stent With Introducer Kit 24	 Ultraflex Precision Single-Use Colonic Stent System..... 46
 Percuflex Duodenal Bend Biliary Stent 24	 Ultraflex Single-Use Covered Esophageal NG Stent System – Distal Release 48
 Polyflex™ Single-Use Esophageal Stent System..... 48	 Ultraflex Single-Use Covered Esophageal NG Stent System – Proximal Release..... 48
 Polyflex Single-Use Self-Expanding Silicone Airway Stent System 82	 Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Distal Release 80
 Profile™ Single-Use Snares..... 75	 Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release 81
 Pulmonary Jagwire™ Pulmonary Guidewire..... 85	 Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Proximal Release..... 80
 Radial Jaw™ 4 Standard Capacity Single-Use Biopsy Forceps 63	 Ultratome™ Double-Lumen Sphincterotome 12
 Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps 63	 Ultratome RX Cannulating Sphincterotome..... 8
 Radial Jaw 4 Single-Use Hot Biopsy Forceps..... 64	 Ultratome XL Triple-Lumen Sphincterotome..... 11
 Radial Jaw 4 Single-Use Jumbo Biopsy Forceps 63	 WallFlex™ Biliary RX Fully Covered Stent System RMV 17
 Radial Jaw 4 Single-Use Pediatric Biopsy Forceps 64	 WallFlex Biliary RX Partially Covered Stent System 17
 Radial Jaw 4 Single-Use Pulmonary Biopsy Forceps 83	 WallFlex Biliary RX Uncovered Stent System..... 17
 Rapid Exchange XL Cannula 12	 WallFlex Fully Covered Esophageal Stent 47
 Rapid Refill Continuous Injection System..... 14	 WallFlex Partially Covered Esophageal Stent 47
 Resolution 360™ Clip..... 51	 WallFlex Colonic Stent System with Anchor Lock Delivery System .. 45
 Resolution™ Clip..... 51	 WallFlex Duodenal Stent System with Anchor Lock Delivery System . 45
 Rigiflex™ II Single-Use Achalasia Balloon Dilatation Catheter..... 60	 WALLSTENT™ Endoscopic Biliary Endoprosthesis Stents 19
 RX Cholangiogram Kits..... 8	 WALLSTENT RX Biliary Endoprosthesis Stent System..... 18
 RX Cytology Brush Wireguided Cytology Brush..... 36	 WALLSTENT Enteral Endoprosthesis Colonic / Duodenal Stent with UNISTEP™ Plus Delivery System 46
 RX Duodenal Bend Biliary Stent with RX Delivery System 23	 WALLSTENT Super Stiff Guidewire..... 46
 RX Locking Devices and Biopsy Cap Systems 14	 Zero Tip™ Airway Retrieval Basket..... 84
 RX Needle Knife XL Triple-Lumen Needle Knife..... 8	
 RX ERCP Cannulas 12	
 Sensation™ Single-Use Short Throw Snares 76	
 Single-Use Rotatable Snares 75	
 Speedband Superview Super 7™ Multiple Band Ligators..... 54	
 SpyBite™ Biopsy Forceps 39, 40	
 SpyGlass™ Direct Visualization System – Disposable and Reusable Devices..... 40	
 SpyGlass DS Direct Visualization System – Disposable Devices.... 39	

Note: RX devices are short wire devices.

 [Terms & Conditions](#)

Biliary Access



Biliary Access

Guidewires

➤ Dreamwire™ High Performance Guidewire.....	3
➤ Hydra Jagwire™ High Performance Guidewire.....	3
➤ Jagtail™ High Performance Guidewire Extension.....	4
➤ Jagwire™ High Performance Guidewire	4
➤ NaviPro™ Hydrophilic Guidewire.....	5
➤ NovaGold™ High Performance Guidewire.....	5
➤ Pathfinder™ Exchange Guidewire.....	5

RX Preloaded Sphincterotomes

➤ Dreamtome™ RX Cannulating Sphincterotome (Preloaded with Dreamwire High Performance Guidewire).....	6
➤ Hydratome™ RX Cannulating Sphincterotome (Preloaded with Hydra Jagwire High Performance Guidewire)..	6
➤ Jagtome™ RX Cannulating Sphincterotome (Preloaded with Jagwire High Performance Guidewire).....	7

RX Sphincterotomes

➤ Autotome™ RX Cannulating Sphincterotome.....	8
➤ RX Cholangiogram Kits	8
➤ RX Needle Knife XL Triple-Lumen Needle Knife	8
➤ Ultratome™ RX Cannulating Sphincterotome	8

Preloaded Sphincterotomes

➤ TRUEtome™ Jag Cannulating Sphincterotome (Preloaded with Jagwire High Performance Guidewire).....	9
➤ TRUEtome Hydra Cannulating Sphincterotome (Preloaded with Hydra Jagwire High Performance Guidewire)..	9
➤ TRUEtome Dreamwire Cannulating Sphincterotome (Preloaded with Dreamwire High Performance Guidewire).....	9

Sphincterotomes

➤ MicroKnife™ XL Triple-Lumen Needle Knife	10
➤ Stonetome™ Stone Removal Device	11
➤ TRUEtome Cannulating Sphincterotome	10
➤ Ultratome Double-Lumen Sphincterotome.....	12
➤ Ultratome XL Triple-Lumen Sphincterotome	11

RX Cannulas

➤ Rapid Exchange XL Cannula	12
➤ RX ERCP Cannulas.....	12

Cannulas

➤ Contour™ ERCP Cannulas	13
➤ Fluoro Tip™ ERCP Cannulas	13
➤ Tandem™ XL Triple-Lumen ERCP Cannula	13

Accessories

➤ Active Cords (Fits Boston Scientific Endoscopy Sphincterotome)	14
➤ RX Locking Devices and Biopsy Cap Systems.....	14


Rapid Refill Continuous Injection System

➤ Rapid Refill Continuous Injection System.....	14
---	----



Biliary

Guidewires




Dreamwire™ High Performance Guidewire



Order Number	GTIN	Description	Distal Tip O.D. (in)	Working Length (cm)	Shaft Style	Dream Tip Style	Tip	Packaging
M00556101	08714729777816	Dreamwire Guidewire ST	.035	260	Standard	Straight	Round	Box 2
M00556111	08714729777847	Dreamwire Guidewire ANG	.035	260	Standard	Angled	Round	Box 2
M00556121	08714729777830	Dreamwire Guidewire ST Stiff	.035	260	Stiff	Straight	Round	Box 2
M00556131	08714729777823	Dreamwire Guidewire ANG Stiff	.035	260	Stiff	Angled	Round	Box 2
M00556141	08714729777885	Dreamwire Guidewire ST	.035	450	Standard	Straight	Round	Box 2
M00556151	08714729777878	Dreamwire Guidewire ANG	.035	450	Standard	Angled	Round	Box 2
M00556161	08714729777861	Dreamwire Guidewire ST Stiff	.035	450	Stiff	Straight	Round	Box 2
M00556171	08714729777854	Dreamwire Guidewire ANG Stiff	.035	450	Stiff	Angled	Round	Box 2

ST = Straight, ANG = Angled



Hydra Jagwire™ High Performance Guidewire



Order Number	GTIN	Description	Distal Tip O.D. (in)	Working Length (cm)	Shaft Style	Dream Tip Style	Tip	Packaging
M00556001	08714729648437	Hydra Jagwire Guidewire ST	.035	260	Standard	Straight	Round	Box 2
M00556011	08714729648444	Hydra Jagwire Guidewire ANG	.035	260	Standard	Angled	Round	Box 2
M00556021	08714729648451	Hydra Jagwire Guidewire ST Stiff	.035	260	Stiff	Straight	Round	Box 2
M00556031	08714729648468	Hydra Jagwire Guidewire ANG Stiff	.035	260	Stiff	Angled	Round	Box 2
M00556041	08714729648475	Hydra Jagwire Guidewire ST	.035	450	Standard	Straight	Round	Box 2
M00556051	08714729648482	Hydra Jagwire Guidewire ANG	.035	450	Standard	Angled	Round	Box 2
M00556061	08714729648499	Hydra Jagwire Guidewire ST Stiff	.035	450	Stiff	Straight	Round	Box 2
M00556071	08714729648505	Hydra Jagwire Guidewire ANG Stiff	.035	450	Stiff	Angled	Round	Box 2

Hydra Jagwire Guidewire may remain in place during sphincterotomy. See Instructions for Use for details.

ST = Straight, ANG = Angled



Section Table of Contents

Prev Page

3


Next Page

Main Table of Contents



Biliary

Guidewires (Continued)



Jagwire™ High Performance Guidewire



Order Number	GTIN	Description	Distal Tip O.D. (in)	Working Length (cm)	Shaft Style	Tip Style	Tip	Packaging
M00556511	08714729874492	Jagwire Guidewire ANG	.025	450	Standard	Angled	Round	Box 2
M00556521	08714729874508	Jagwire Guidewire ST Stiff	.025	450	Stiff	Straight	Round	Box 2
M00556531	08714729874515	Jagwire Guidewire ANG Stiff	.025	450	Stiff	Angled	Round	Box 2
M00556561	08714729501374	Jagwire Guidewire ST	.025	450	Standard	Straight	Trim	Box 2
M00556581	08714729339151	Jagwire Guidewire ST	.035	450	Standard	Straight	Trim	Box 2
M00556591	08714729501398	Jagwire Guidewire ANG	.035	450	Standard	Angled	Trim	Box 2
M00556601	08714729172154	Jagwire Guidewire ST Stiff	.035	450	Stiff	Straight	Trim	Box 2
M00556611	08714729172161	Jagwire Guidewire ANG Stiff	.035	450	Stiff	Angled	Trim	Box 2
M00556621	08714729172178	Jagwire Guidewire ST	.038	260	Standard	Straight	Trim	Box 2
M00556631	08714729172185	Jagwire Guidewire ST Stiff	.035	260	Stiff	Straight	Trim	Box 2
M00556641	08714729172192	Jagwire Guidewire ST	.035	260	Standard	Straight	Trim	Box 2
M00556651	08714729172208	Jagwire Guidewire ANG	.035	260	Standard	Angled	Trim	Box 2
M00556661	08714729172215	Jagwire Guidewire ANG Stiff	.035	260	Stiff	Angled	Trim	Box 2
M00556461	08714729415183	Jagwire Guidewire ST	.025	260	Standard	Straight	Round	Box 2
M00556481	08714729415190	Jagwire Guidewire ST Stiff	.025	260	Stiff	Straight	Round	Box 2
M00556491	08714729415206	Jagwire Guidewire ANG Stiff	.025	260	Stiff	Angled	Round	Box 2
M00556841	08714729268680	Extendable Jagwire Guidewire ST	.035	260	Standard	Straight	Round	Box 2
M00556851	08714729268697	Extendable Jagwire Guidewire ANG	.035	260	Standard	Angled	Round	Box 2
M00556861	08714729268703	Extendable Jagwire Guidewire ST Stiff	.035	260	Stiff	Straight	Round	Box 2
M00556871	08714729268710	Extendable Jagwire Guidewire ANG Stiff	.035	260	Stiff	Angled	Round	Box 2
M00556902	08714729268734	Jagtail™ Guidewire Extension	.035	200	—	Guidewire extension	—	Box 5

Jagwire Guidewire may remain in place during sphincterotomy. See Instructions for Use for details.

ST = Straight, ANG = Angled



Section Table of Contents

Prev Page

4


Next Page

Main Table of Contents



Biliary


Guidewires (Continued)



NaviPro™ Hydrophilic Guidewire




Order Number	GTIN	Description	Guidewire Diameter (in)	Length (cm)	Shaft Style	Tip Style	Packaging
M00556191	08714729839873	NaviPro Hydrophilic Guidewire	0.018	260	Standard	Angled	Box 5
M00556181	08714729839866	NaviPro Hydrophilic Guidewire	0.018	260	Standard	Straight	Box 5
M00556211	08714729839897	NaviPro Hydrophilic Guidewire	0.025	260	Standard	Angled	Box 5
M00556201	08714729839880	NaviPro Hydrophilic Guidewire	0.025	260	Standard	Straight	Box 5
M00556231	08714729839910	NaviPro Hydrophilic Guidewire	0.035	260	Standard	Angled	Box 5
M00556221	08714729839903	NaviPro Hydrophilic Guidewire	0.035	260	Standard	Straight	Box 5
M00556251	08714729839934	NaviPro Hydrophilic Guidewire	0.035	260	Stiff	Angled	Box 5
M00556241	08714729839927	NaviPro Hydrophilic Guidewire	0.035	260	Stiff	Straight	Box 5



NovaGold™ High Performance Guidewire



Order Number	GTIN	Description	Distal Tip O.D. (in)	Working Length (cm)	Tip Style	Tip Length (cm)	Packaging
M00552001	08714729854418	NovaGold High Performance Guidewire	.018	260	Straight (shapeable)	6	Box 2
M00552011	08714729854425	NovaGold High Performance Guidewire	.018	480	Straight (shapeable)	6	Box 2



Pathfinder™ Exchange Guidewire



Order Number	GTIN	Description	Distal Tip O.D. (in)	Working Length (cm)	Tip Style	Tip Length (cm)	Packaging
M00551591	08714729156116	Pathfinder Exchange Guidewire Stiff Shaft	.018	450	Shapeable	6.5	Box 2



Section Table of Contents

Prev Page

5



Next Page

Main Table of Contents



Biliary

RX Preloaded Sphincterotome





Dreamtome™ RX Cannulating Sphincterotome (Preloaded with Dreamwire™ High Performance Guidewire)



Order Number	GTIN	Description	Tip Length (mm)	Cut Wire Length (mm)	Preloaded Guidewire (in / cm)	Guidewire Tip Style	Packaging
M00584000	08714729778226	Dreamtome RX 49 Sphincterotome	5	20	.035 / 260 Dreamwire Guidewire	Straight	Each
M00584010	08714729778233	Dreamtome RX 49 Sphincterotome	5	30	.035 / 260 Dreamwire Guidewire	Straight	Each
M00584020	08714729778240	Dreamtome RX 49 Sphincterotome	5	20	.035 / 450 Dreamwire Guidewire	Straight	Each
M00584030	08714729778257	Dreamtome RX 49 Sphincterotome	5	30	.035 / 450 Dreamwire Guidewire	Straight	Each
M00584040	08714729778264	Dreamtome RX 44 Sphincterotome	5	20	.035 / 260 Dreamwire Guidewire	Straight	Each
M00584050	08714729778271	Dreamtome RX 44 Sphincterotome	5	30	.035 / 260 Dreamwire Guidewire	Straight	Each
M00584060	08714729778288	Dreamtome RX 44 Sphincterotome	5	20	.035 / 450 Dreamwire Guidewire	Straight	Each
M00584070	08714729778295	Dreamtome RX 44 Sphincterotome	5	30	.035 / 450 Dreamwire Guidewire	Straight	Each

*Preloaded with straight-tip guidewires only.
Active Cords sold separately [▶](#) (Click here to view page 14).



Hydratome™ RX Cannulating Sphincterotome (Preloaded with Hydra Jagwire™ High Performance Guidewire)



Order Number	GTIN	Description	Cut Wire Length (mm)	Tip Length (mm)	Distal Tip O.D. (Fr)	Preloaded Guidewire (cm)	Guidewire Tip Style	Packaging
M00583000	08714729755159	Hydratome RX 49 Sphincterotome	20	5	4.9	260 Hydra Jagwire Guidewire	Straight	Each
M00583010	08714729755166	Hydratome RX 49 Sphincterotome	30	5	4.9	260 Hydra Jagwire Guidewire	Straight	Each
M00583020	08714729755173	Hydratome RX 49 Sphincterotome	20	5	4.9	450 Hydra Jagwire Guidewire	Straight	Each
M00583030	08714729755180	Hydratome RX 49 Sphincterotome	30	5	4.9	450 Hydra Jagwire Guidewire	Straight	Each
M00583040	08714729755197	Hydratome RX 44 Sphincterotome	20	5	4.4	260 Hydra Jagwire Guidewire	Straight	Each
M00583050	08714729755203	Hydratome RX 44 Sphincterotome	30	5	4.4	260 Hydra Jagwire Guidewire	Straight	Each
M00583060	08714729755210	Hydratome RX 44 Sphincterotome	20	5	4.4	450 Hydra Jagwire Guidewire	Straight	Each
M00583070	08714729755227	Hydratome RX 44 Sphincterotome	30	5	4.4	450 Hydra Jagwire Guidewire	Straight	Each

Active Cords sold separately [▶](#) (Click here to view page 14).



Section Table of Contents

Prev Page

6

Next Page

Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Biliary

RX Preloaded Sphincterotome (Continued)



Jagtome™ RX Cannulating Sphincterotome (Preloaded with Jagwire™ High Performance Guidewire)



Order Number	GTIN	Description	Cut Wire (mm)	Tip Length (mm)	Distal Tip O.D. (Fr)	Preloaded Guidewire (in / cm)	Guidewire Tip Style	Packaging
M00573000	08714729776710	Jagtome RX 49 Sphincterotome	20	5	4.9	.035 / 260 Jagwire Guidewire	Straight	Each
M00573010	08714729776727	Jagtome RX 49 Sphincterotome	20	5	4.9	.035 / 450 Jagwire Guidewire	Straight	Each
M00573020	08714729776734	Jagtome RX 49 Sphincterotome	30	5	4.9	.035 / 260 Jagwire Guidewire	Straight	Each
M00573030	08714729776741	Jagtome RX 49 Sphincterotome	30	5	4.9	.035 / 450 Jagwire Guidewire	Straight	Each
M00573040	08714729776659	Jagtome RX 44 Sphincterotome	20	5	4.4	.035 / 260 Jagwire Guidewire	Straight	Each
M00573050	08714729776666	Jagtome RX 44 Sphincterotome	20	5	4.4	.035 / 450 Jagwire Guidewire	Straight	Each
M00573060	08714729776673	Jagtome RX 44 Sphincterotome	30	5	4.4	.035 / 260 Jagwire Guidewire	Straight	Each
M00573070	08714729776543	Jagtome RX 44 Sphincterotome	30	5	4.4	.035 / 450 Jagwire Guidewire	Straight	Each
M00573080	08714729776642	Jagtome RX 39 Sphincterotome	20	5	3.9	.025 / 260 Jagwire Guidewire	Straight	Each
M00573090	08714729776550	Jagtome RX 39 Sphincterotome	20	5	3.9	.025 / 450 Jagwire Guidewire	Straight	Each
M00573100	08714729776536	Jagtome RX 39 Sphincterotome	30	5	3.9	.025 / 260 Jagwire Guidewire	Straight	Each
M00573110	08714729776529	Jagtome RX 39 Sphincterotome	30	5	3.9	.025 / 450 Jagwire Guidewire	Straight	Each

Active Cords sold separately [▶ \(Click here to view page 14\).](#)



Section Table of Contents

Prev Page

7


Next Page

Main Table of Contents



Biliary

RX Sphincterotome



Autotome™ RX Cannulating Sphincterotome



Order Number	GTIN	Description	Cut Wire (mm)	Tip Length (mm)	Distal Tip O.D. (Fr)	Recommended Guidewire (in)	Packaging
M00545150	08714729415138	Autotome RX 49 Cannulating Sphincterotome	20	5	4.9	.035	Each
M00545160	08714729415152	Autotome RX 49 Cannulating Sphincterotome	30	5	4.9	.035	Each
M00545170	08714729444749	Autotome RX 44 Cannulating Sphincterotome	20	5	4.4	.035	Each
M00545180	08714729444756	Autotome RX 44 Cannulating Sphincterotome	30	5	4.4	.035	Each
M00545190	08714729444763	Autotome RX 39 Cannulating Sphincterotome	20	5	3.9	.025	Each
M00545200	08714729444770	Autotome RX 39 Cannulating Sphincterotome	30	5	3.9	.025	Each

Active Cords sold separately [▶](#) (Click here to view page 14).

Ultratome™ RX Cannulating Sphincterotome



Order Number	GTIN	Description	Cut Wire (mm)	Tip Length (mm)	Distal Tip O.D. (Fr)	Recommended Guidewire (in)	Packaging
M00545050	08714729251903	Ultratome RX Cannulating Sphincterotome, Short Nose	20	5	5.5	.035	Each
M00545060	08714729251927	Ultratome RX Cannulating Sphincterotome, Short Nose	30	5	5.5	.035	Each

Active Cords sold separately [▶](#) (Click here to view page 14).

RX Cholangiogram Kits*



Order Number	Description	Sphincterotome	Guidewire	Accessory
M005XXXX1	Cholangiogram Kit	Autotome or Ultratome XL Sphincterotome	Jagwire™ Guidewire	N / A
M005XXXX1	Cholangiogram Kit-Preloaded	Autotome Sphincterotome	Jagwire Guidewire	RX Locking Device

*To order RX Cholangiogram Kits please contact your local Boston Scientific Endoscopy Sales Representative.

RX Needle Knife XL Triple-Lumen Needle Knife



Order Number	GTIN	Description	Length (cm)	Cut Wire (mm)	Distal Tip O.D. (Fr)	Recommended Guidewire (in)	Packaging
M00545840	08714729283799	RX Needle Knife	200	5	5.5	.035	Each



Section Table of Contents

Prev Page

8

Next Page

Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Biliary


Sphincterotome

TRUEtome™ Jag Cannulating Sphincterotome (Preloaded with Jagwire™ High Performance Guidewire)								
Order Number	GTIN	Description	Tip Length (mm)	Cut Wire Length (mm)	Preloaded Guidewire (in / cm)	Tip O.D. (Fr)	Packaging	
M00583080	08714729858621	TRUEtome 39 Sphincterotome	5	20	.025 / 450 Jagwire Guidewire	3.9	Each	
M00583090	08714729858638	TRUEtome 44 Sphincterotome	5	20	.025 / 450 Jagwire Guidewire	4.4	Each	
M00583100	08714729858645	TRUEtome 44 Sphincterotome	5	20	.035 / 450 Jagwire Guidewire	4.4	Each	
M00583110	08714729858652	TRUEtome 44 Sphincterotome	5	20	.035 / Stiff 450 Jagwire Guidewire	4.4	Each	
M00583120	08714729858669	TRUEtome 39 Sphincterotome	5	30	.025 / 450 Jagwire Guidewire	3.9	Each	
M00583130	08714729858676	TRUEtome 44 Sphincterotome	5	30	.025 / 450 Jagwire Guidewire	4.4	Each	
M00583140	08714729858683	TRUEtome 44 Sphincterotome	5	30	.035 / 450 Jagwire Guidewire	4.4	Each	
M00583150	08714729858690	TRUEtome 44 Sphincterotome	5	30	.035 / Stiff 450 Jagwire Guidewire	4.4	Each	
TRUEtome Hydra Cannulating Sphincterotome (Preloaded with Hydra Jagwire™ High Performance Guidewire)								
M00583160	08714729858706	TRUEtome 44 Sphincterotome	5	20	.035 / 450 Hydra Jagwire Guidewire	4.4	Each	
M00583170	08714729858713	TRUEtome 44 Sphincterotome	5	30	.035 / 450 Hydra Jagwire Guidewire	4.4	Each	
TRUEtome Dreamwire™ Cannulating Sphincterotome (Preloaded with Dreamwire High Performance Guidewire)								
M00583180	08714729858720	TRUEtome 44 Sphincterotome	5	20	.035 / 260 Dreamwire Guidewire	4.4	Each	
M00583190	08714729858737	TRUEtome 44 Sphincterotome	5	20	.035 / 450 Dreamwire Guidewire	4.4	Each	
M00583200	08714729858744	TRUEtome 44 Sphincterotome	5	30	.035 / 260 Dreamwire Guidewire	4.4	Each	
M00583210	08714729858751	TRUEtome 44 Sphincterotome	5	30	.035 / 450 Dreamwire Guidewire	4.4	Each	

*Preloaded with straight-tip guidewires only.
Active Cords sold separately  (Click here to view page 14).

Biliary

Sphincterotome (Continued)




TRUEtome™ Cannulating Sphincterotome



Order Number	GTIN	Description	Cut Wire Length (mm)	Tip Length (mm)	Distal Tip O.D. (Fr)	Recommended Guidewire (in)	Packaging
M00584190	08714729827658	TRUEtome 39 Sphincterotome	20	5	3.9	.025	Each
M00584120	08714729827771	TRUEtome 39 Sphincterotome	30	5	3.9	.025	Each
M00584170	08714729827634	TRUEtome 44 Sphincterotome	20	5	4.4	.035	Each
M00584180	08714729827641	TRUEtome 44 Sphincterotome	30	5	4.4	.035	Each
M00584150	08714729827610	TRUEtome 49 Sphincterotome	20	5	4.9	.035	Each
M00584160	08714729827627	TRUEtome 49 Sphincterotome	30	5	4.9	.035	Each

Active Cords sold separately [▶](#) (Click here to view page 14).



MicroKnife™ XL Triple-Lumen Needle Knife



Order Number	GTIN	Description	Shaft O.D. (Fr)	Catheter Length (cm)	Recommended Guidewire (in)	Packaging
M00532810	08714729172017	MicroKnife XL Triple-Lumen Needle Knife	7.0 - 5.5	200	.035	Each

Active Cords sold separately [▶](#) (Click here to view page 14).



Section Table of Contents


Prev Page **10** Next Page

Main Table of Contents



Biliary

Sphincterotome (Continued)




Stonetome™ Stone Removal Device



Order Number	GTIN	Description	Balloon Diameter (mm)	Catheter Length (cm)	Tip Length (mm)	Cut Wire (mm)	Distal Tip O.D. (Fr)	Recommended Guidewire (in)	Packaging
M00535150	08714729146667	Stonetome Device (below cut wire)	11.5	200	5	20	5.5	.035	Each
M00535170	08714729146629	Stonetome Device (below cut wire)	11.5	200	20	20	5.5	.035	Each
M00535190	08714729146612	Stonetome Device (below cut wire)	11.5	200	5	30	5.5	.035	Each
M00535210	08714729146605	Stonetome Device (below cut wire)	11.5	200	20	30	5.5	.035	Each
M00535110	08714729146599	Stonetome Device (above cut wire)	11.5	200	20	20	5.5	.035	Each
M00535130	08714729146575	Stonetome Device (above cut wire)	11.5	200	20	30	5.5	.035	Each

Active Cords sold separately [▶](#) (Click here to view page 14).



Ultratome™ XL Triple-Lumen Sphincterotome



Order Number	GTIN	Description	Wire Length (mm)	Shaft O.D. (Fr)	Catheter Length (cm)	Tip Length (mm)	Recommended Guidewire (in)	Packaging
M00535900	08714729103257	Ultratome XL Short Nose Sphincterotome	20	5.5	200	5	.035	Each
M00535910	08714729103301	Ultratome XL Long Nose Sphincterotome	20	5.5	200	20	.035	Each
M00535920	08714729103318	Ultratome XL Short Nose Sphincterotome	30	5.5	200	5	.035	Each
M00535930	08714729103349	Ultratome XL Long Nose Sphincterotome	30	5.5	200	20	.035	Each

Active Cords sold separately [▶](#) (Click here to view page 14).



Section Table of Contents

Prev Page

11


Next Page

Main Table of Contents



Biliary

Sphincterotome (Continued)




Ultratome™ Double-Lumen Sphincterotome



Order Number	GTIN	Description	Wire Length (mm)	Distal Tip O.D.		Catheter Length (cm)	Tip Length (mm)	Recommended Guidewire (in)	Packaging
				(mm)	(Fr)				
M00530800	08714729163664	Ultratome Short Nose Sphincterotome	20	1.8	5.5	200	5	.035	Each
M00530810	08714729183143	Ultratome Long Nose Sphincterotome	20	1.8	5.5	200	20	.035	Each
M00530820	08714729193418	Ultratome Short Nose Sphincterotome	30	1.8	5.5	200	5	.035	Each
M00530830	08714729193425	Ultratome Long Nose Sphincterotome	30	1.8	5.5	200	20	.035	Each

Active Cords sold separately [▶ \(Click here to view page 14\).](#)


RX Cannulas



Rapid Exchange XL Cannula



Order Number	GTIN	Description	Working Length (cm)	Tip Style	Recommended Guidewire (in)	Packaging
M00545820	08714729301905	Rapid Exchange XL Cannula	195	Tapered	.035	Each



RX ERCP Cannulas



Order Number	GTIN	Description	Working Length (cm)	Recommended Guidewire (in)	Packaging
M00545760	08714729268642	RX ERCP Cannula Standard Tip	210	.035	Each
M00545770	08714729268659	RX ERCP Cannula Tapered Tip	210	.035	Each
M00545780	08714729268673	RX ERCP Cannula Ball Tip	210	.035	Each



Section Table of Contents

Prev Page **12** Next Page


Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Biliary


Cannulas



Tandem™ XL Triple-Lumen ERCP Cannula




Order Number	GTIN	Distal Tip O.D.		Catheter Size (Fr)	Working Length (cm)	Tip Style	Recommended Guidewire (in)	Packaging
		(Fr)	(mm)					
M00535700	08714729146636	5.5	1.8	7.0 – 5.5	200	Tapered	.035	Each



Contour™ ERCP Cannulas



Order Number	GTIN	Description	Distal Tip O.D.		Working Length (cm)	Tip Style	Recommended Guidewire (in)	Packaging
			(Fr)	(mm)				
M00530850	08714729099789	Contour ERCP Cannula	5	1.7	210	Standard	.035	Each
M00530860	08714729099819	Contour ERCP Cannula	5	1.7	210	Tapered	.035	Each
M00530870	08714729099826	Contour ERCP Cannula	5	1.7	210	Ultra Tapered	.025	Each
M00530880	08714729099833	Contour ERCP Cannula	5	1.7	210	5-4-3	.018	Each
M00530890	08714729099840	Contour ERCP Cannula	5	1.7	210	Ball Tip	.035	Each



Fluoro Tip™ ERCP Cannulas



Order Number	GTIN	Description	Distal Tip O.D.		Working Length (cm)	Recommended Guidewire (in)	Packaging
			(mm)	(Fr)			
M00530960	08714729172277	Fluoro Tip ERCP Cannula Standard Tip	1.7	5	210	.035	Each
M00530970	08714729172284	Fluoro Tip ERCP Cannula Tapered Tip	1.7	5	210	.035	Each
M00530980	08714729193432	Fluoro Tip ERCP Cannula Ultra Tapered Tip	1.7	5	210	.025	Each



Section Table of Contents

Prev Page **13** Next Page

Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Biliary


Accessories

RX Locking Devices and Biopsy Cap Systems

Order Number	GTIN	Description	Packaging
M00545471	08714729337768	RX Locking Device (Biopsy Cap not included)	Box 10
M00545211	08714729440604	RX Locking Device – Olympus™ / Fujinon™ Brand Endoscope Compatible	Box 10
M00545261	08714729296614	RX Locking Device and Biopsy Cap System – Olympus Brand Endoscope Compatible	Box 10
M00545271	08714729296621	RX Locking Device and Biopsy Cap System – Pentax™ Brand Endoscope Compatible	Box 10
M00545281	08714729296638	RX Locking Device and Biopsy Cap System – Fujinon Brand Endoscope Compatible	Box 10

Active Cords (Fit Boston Scientific Endoscopy Sphincterotome)

Order Number	GTIN	Description	Packaging
M00561250	08714729019060	Active Cord, Fits Olympus Brand Generator	Each
M00561270	08714729053361	Active Cord, Fits Endostat™, Endostat II and Endostat III Generators, Bovie™, Valleylab™, Aspen Labs™, and ERBE™ Brand Generators	Each



Rapid Refill Continuous Injection System



Order Number	GTIN	Description	Packaging
M00566001	08714729285151	Rapid Refill Continuous Injection System	Box 10



Notes and Future Product Updates



Notes and Future Product Updates



Biliary Stricture and Stone Management

Biliary and Pancreatic Stents

Biliary Metal Stent Systems

- WallFlex™ Biliary RX Fully Covered Stent System RMV..... 17
- WallFlex Biliary RX Partially Covered Stent System..... 17
- WallFlex Biliary RX Uncovered Stent System 17
- WALLSTENT™ RX Biliary Endoprosthesis Stent System 18
- WALLSTENT Endoscopic Biliary Endoprosthesis Stents..... 19

Biliary Plastic Stents

- Advanix™ Center Bend Preloaded RX Stent Systems..... 20
- Advanix Center Bend - Single Stent..... 21
- Advanix Double Pigtail - Single Stent 22
- Advanix Duodenal Bend Preloaded RX Stent Systems..... 20
- Advanix Duodenal Bend - Single Stent 21
- Biliary Introducer..... 24
- C-Flex™ Double Pigtail Biliary Stent..... 25
- Flexima™ Duodenal Bend Biliary Stent with Delivery System .. 23
- Flexima Single Pigtail Nasobiliary Catheter..... 25
- NaviFlex™ RX Delivery System..... 22
- Percuflex™ Duodenal Bend Biliary Stent..... 24
- Percuflex Duodenal Bend Biliary Stent with Introducer Kit 24
- RX Duodenal Bend Biliary Stent with RX Delivery System..... 23

Pancreatic Plastic Stents

- Advanix Pancreatic Stent Kit - Straight, No Leading Barb..... 26
- Advanix Pancreatic Stent Kit - Pigtail, Leading Barb and No Leading Barb..... 26
- Advanix Pancreatic Stent - Straight, No Leading Barb 27
- Advanix Pancreatic Stent - Straight, Leading Barb..... 28
- Advanix Pancreatic Stent - Pigtail, No Leading Barb..... 29
- Advanix Pancreatic Stent - Pigtail, Leading Barb..... 30
- NaviFlex RX Pancreatic Delivery System & NaviFlex Pushers.. 31

Biliary Stone Management

RX Stone Management

- Extractor™ Pro RX Retrieval Balloon Catheters – Distal Injection .. 32
- Extractor Pro RX Retrieval Balloon Catheters – Proximal Injection . 32
- Extractor Pro RX-S Retrieval Balloon Catheters – Distal Injection . 32
- Extractor Pro RX-S Retrieval Balloon Catheters – Proximal Injection . 32

Stone Management

- Extractor Pro DL Retrieval Balloon Catheters..... 33
- Extractor Pro XL Retrieval Balloon Catheters – Distal Injection..... 33
- Extractor Pro XL Retrieval Balloon Catheters – Proximal Injection 33
- Trapezoid™ RX Wireguided Retrieval Baskets 33

CRE™ PRO Wireguided Balloon Dilatation Catheters

- CRE PRO Wireguided Balloon Dilatation Catheters..... 34

Biliary Stricture Management

RX Dilatation

- Hurricane™ RX Biliary Balloon Dilatation Catheter..... 35

Biliary Dilatation

- MaxForce™ Biliary Balloon Dilatation Catheter..... 36

RX Accessories

- Encore™ 26 Inflation Device 36

RX Cytology

- RX Cytology Brush Wireguided Cytology Brush 36



Biliary Metal Stents

WallFlex™ Biliary RX Stents



WallFlex Biliary RX Fully Covered Stent System RMV*

Order Number	GTIN	Diameter (mm)	Length (mm)	Covered Length (mm) (PC Only)	Catheter Diameter (Fr / mm)	Guidewire Diameter (in / mm)	Working Length (cm)	Packaging
M00570340	08714729953159	8	60	—	8.5 / 2.83	.035 / .89	194	Each
M00570350	08714729953166	8	80	—	8.5 / 2.83	.035 / .89	194	Each
M00570360	08714729953173	10	40	—	8.5 / 2.83	.035 / .89	194	Each
M00570370	08714729953180	10	60	—	8.5 / 2.83	.035 / .89	194	Each
M00570380	08714729953197	10	80	—	8.5 / 2.83	.035 / .89	194	Each

*Compatible with Boston Scientific .035" guidewires, including Dreamwire™ Guidewire (Click here to view page 3).

Above (WallFlex Biliary RX Partially Covered Stent System) with Permalume™ Covering.

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



WallFlex Biliary RX Partially Covered Stent System*

Order Number	GTIN	Diameter (mm)	Length (mm)	Covered Length (mm) (PC Only)	Catheter Diameter (Fr / mm)	Guidewire Diameter (in / mm)	Working Length (cm)	Packaging
M00570700	08714729759294	8	60	48	8.5 / 2.83	.035 / .89	194	Each
M00570710	08714729759300	8	80	68	8.5 / 2.83	.035 / .89	194	Each
M00570720	08714729759317	10	40	28	8.5 / 2.83	.035 / .89	194	Each
M00570730	08714729759324	10	60	48	8.5 / 2.83	.035 / .89	194	Each
M00570740	08714729759331	10	80	68	8.5 / 2.83	.035 / .89	194	Each

*Compatible with Boston Scientific .035" guidewires, including Dreamwire Guidewire (Click here to view page 3).

Above (WallFlex Biliary RX Partially Covered Stent System) with Permalume Covering.

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



WallFlex Biliary RX Uncovered Stent System*

Order Number	GTIN	Diameter (mm)	Length (mm)	Covered Length (mm) (PC Only)	Catheter Diameter (Fr / mm)	Guidewire Diameter (in / mm)	Working Length (cm)	Packaging
M00570600	08714729759249	8	40	—	8.0 / 2.67	.035 / .89	194	Each
M00570610	08714729759256	8	60	—	8.0 / 2.67	.035 / .89	194	Each
M00570620	08714729759263	8	80	—	8.0 / 2.67	.035 / .89	194	Each
M00570630	08714729759270	8	100	—	8.0 / 2.67	.035 / .89	194	Each
M00570890	08714729765332	10	40	—	8.0 / 2.67	.035 / .89	194	Each
M00570640	08714729756231	10	60	—	8.0 / 2.67	.035 / .89	194	Each
M00570650	08714729758150	10	80	—	8.0 / 2.67	.035 / .89	194	Each
M00570660	08714729759287	10	100	—	8.0 / 2.67	.035 / .89	194	Each

*Compatible with Boston Scientific .035" guidewires, including Dreamwire Guidewire (Click here to view page 3).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

Biliary Metal Stents

WALLSTENT™ Biliary RX Stents



 Magnetic Resonance Conditional¹

WALLSTENT RX Biliary Endoprosthesis Stent System*



Order Number	GTIN	Description	Unconstrained Diameter (mm)	Nominal Unconstrained Length (mm)	Delivery System (Fr)	Cover Length (mm)	Working Length (cm)	Packaging
M00569670	08714729296874	WALLSTENT RX Biliary Partially Covered Stent System	8	40	8	30	194	Each
M00569680	08714729296881	WALLSTENT RX Biliary Partially Covered Stent System	8	60	8	50	194	Each
M00569690	08714729296898	WALLSTENT RX Biliary Partially Covered Stent System	8	80	8	70	194	Each
M00569700	08714729296904	WALLSTENT RX Biliary Partially Covered Stent System	10	40	8	30	194	Each
M00569710	08714729296911	WALLSTENT RX Biliary Partially Covered Stent System	10	60	8	50	194	Each
M00569720	08714729296928	WALLSTENT RX Biliary Partially Covered Stent System	10	80	8	70	194	Each

*Recommended Guidewire: .035", 260cm, including Dreamwire™ Guidewire [\(Click here to view page 3\)](#).

Above (Wallstent RX Biliary Endoprosthesis Stent System) with Permalume™ Covering.

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



 Magnetic Resonance Conditional¹

WALLSTENT RX Biliary Endoprosthesis Stent System*



Order Number	GTIN	Description	Unconstrained Diameter (mm)	Nominal Unconstrained Length (mm)	Delivery System (Fr)	Working Length (cm)	Packaging
M00569610	08714729296812	WALLSTENT RX Biliary Uncovered Stent System	8	40	8	194	Each
M00569620	08714729296829	WALLSTENT RX Biliary Uncovered Stent System	8	60	8	194	Each
M00569630	08714729296836	WALLSTENT RX Biliary Uncovered Stent System	8	80	8	194	Each
M00569800	08714729456711	WALLSTENT RX Biliary Uncovered Stent System	8	100	8	194	Each
M00569640	08714729296843	WALLSTENT RX Biliary Uncovered Stent System	10	40	8	194	Each
M00569650	08714729296850	WALLSTENT RX Biliary Uncovered Stent System	10	60	8	194	Each
M00569660	08714729296867	WALLSTENT RX Biliary Uncovered Stent System	10	80	8	194	Each
M00569810	08714729456728	WALLSTENT RX Biliary Uncovered Stent System	10	100	8	194	Each

*Recommended Guidewire: .035", 260cm, including Dreamwire Guidewire [\(Click here to view page 3\)](#).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



Section Table of Contents

Prev Page

18

Next Page


Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Biliary Metal Stents

WALLSTENT™ Stents



 Magnetic Resonance Conditional¹

WALLSTENT Endoscopic Biliary Endoprosthesis Stents*




Order Number	GTIN	Description	Unconstrained Diameter (mm)	Nominal Unconstrained Length (mm)	Delivery System (Fr)	Cover Length (mm)	Working Length (cm)	Packaging
H965430700	08714729229858	WALLSTENT Biliary Partially Covered Stent	8	40	8	30	194	Each
H965430800	08714729317906	WALLSTENT Biliary Partially Covered Stent	8	60	8	50	194	Each
H965430900	08714729229872	WALLSTENT Biliary Partially Covered Stent	8	80	8	70	194	Each
H965431000	08714729317920	WALLSTENT Biliary Partially Covered Stent	10	40	8	30	194	Each
H965431100	08714729317937	WALLSTENT Biliary Partially Covered Stent	10	60	8	50	194	Each
H965431200	08714729229902	WALLSTENT Biliary Partially Covered Stent	10	80	8	70	194	Each

*Recommended Guidewire: .035", 450cm, including Dreamwire™ Guidewire [▶](#) (Click here to view page 3).

Above (Wallstent Biliary Endoprosthesis) with Permalume™ Covering.

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



 Magnetic Resonance Conditional¹

WALLSTENT Endoscopic Biliary Endoprosthesis Stents*



Order Number	GTIN	Description	Unconstrained Diameter (mm)	Nominal Unconstrained Length (mm)	Delivery System Diameter (Fr)	Working Length (cm)	Packaging
H965430100	08714729204633	WALLSTENT Biliary Uncovered Stent	8	40	7.5	194	Each
H965430200	08714729204640	WALLSTENT Biliary Uncovered Stent	8	60	7.5	194	Each
H965430300	08714729204657	WALLSTENT Biliary Uncovered Stent	8	80	7.5	194	Each
H965430400	08714729204664	WALLSTENT Biliary Uncovered Stent	10	40	7.5	194	Each
H965430500	08714729204671	WALLSTENT Biliary Uncovered Stent	10	60	7.5	194	Each
H965430600	08714729204688	WALLSTENT Biliary Uncovered Stent	10	80	7.5	194	Each

*Recommended Guidewire: .035", 450cm, including Dreamwire Guidewire [▶](#) (Click here to view page 3).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



Section Table of Contents

Prev Page

19

Next Page

Main Table of Contents



Biliary Plastic Stents

Advanix™ Biliary Stents



MR Magnetic Resonance Safe

Advanix Duodenal Bend Preloaded RX Stent Systems



Order Number	GTIN	Description	Stent Diameter (Fr)	Length (cm) Barb to Barb	Packaging
M00534200	08714729787297	Duodenal Bend	7	5	Each
M00534210	08714729787303	Duodenal Bend	7	7	Each
M00534220	08714729787310	Duodenal Bend	7	9	Each
M00534230	08714729787327	Duodenal Bend	7	12	Each
M00534240	08714729787334	Duodenal Bend	7	15	Each
M00534250	08714729787341	Duodenal Bend	7	18	Each
M00534260	08714729787358	Duodenal Bend	8.5	5	Each
M00534270	08714729787365	Duodenal Bend	8.5	7	Each
M00534280	08714729787372	Duodenal Bend	8.5	9	Each
M00534290	08714729787389	Duodenal Bend	8.5	12	Each
M00534300	08714729787396	Duodenal Bend	8.5	15	Each
M00534310	08714729787402	Duodenal Bend	8.5	18	Each
M00534320	08714729787419	Duodenal Bend	10	5	Each
M00534330	08714729787426	Duodenal Bend	10	7	Each
M00534340	08714729787433	Duodenal Bend	10	9	Each
M00534350	08714729787440	Duodenal Bend	10	12	Each
M00534360	08714729787457	Duodenal Bend	10	15	Each
M00534370	08714729787464	Duodenal Bend	10	18	Each

Scope Compatibility: 7Fr and 8.5Fr stents require a 3.2mm working channel; 10Fr stents require a 4.2mm working channel.



MR Magnetic Resonance Safe

Advanix Center Bend Preloaded RX Stent Systems



Order Number	GTIN	Description	Stent Diameter (Fr)	Length (cm) Barb to Barb	Packaging
M00534560	08714729787112	Center Bend	7	5	Each
M00534570	08714729787129	Center Bend	7	7	Each
M00534580	08714729787136	Center Bend	7	9	Each
M00534590	08714729787143	Center Bend	7	12	Each
M00534600	08714729787150	Center Bend	7	15	Each
M00534610	08714729787167	Center Bend	7	18	Each
M00534620	08714729787174	Center Bend	8.5	5	Each
M00534630	08714729787181	Center Bend	8.5	7	Each
M00534640	08714729787198	Center Bend	8.5	9	Each
M00534650	08714729787204	Center Bend	8.5	12	Each
M00534660	08714729787211	Center Bend	8.5	15	Each
M00534670	08714729787228	Center Bend	8.5	18	Each
M00534680	08714729787235	Center Bend	10	5	Each
M00534690	08714729787242	Center Bend	10	7	Each
M00534700	08714729787259	Center Bend	10	9	Each
M00534710	08714729787266	Center Bend	10	12	Each
M00534720	08714729787273	Center Bend	10	15	Each
M00534730	08714729787280	Center Bend	10	18	Each

Scope Compatibility: 7Fr and 8.5Fr stents require a 3.2mm working channel; 10Fr stents require a 4.2mm working channel.



Section Table of Contents

Prev Page

20

Next Page

Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Biliary Plastic Stents

Advanix™ Biliary Stents (Continued)



MR Magnetic Resonance Safe

Advanix Duodenal Bend - Single Stent



Order Number	GTIN	Description	Stent Diameter (Fr)	Length (cm) Barb to Barb	Packaging
M00532870	08714729786788	Duodenal Bend	7	5	Each
M00532880	08714729786795	Duodenal Bend	7	7	Each
M00532890	08714729786801	Duodenal Bend	7	9	Each
M00532900	08714729786818	Duodenal Bend	7	12	Each
M00532910	08714729786825	Duodenal Bend	7	15	Each
M00532920	08714729786832	Duodenal Bend	7	18	Each
M00532930	08714729786849	Duodenal Bend	8.5	5	Each
M00532940	08714729786856	Duodenal Bend	8.5	7	Each
M00532950	08714729786863	Duodenal Bend	8.5	9	Each
M00532960	08714729786870	Duodenal Bend	8.5	12	Each
M00532970	08714729786887	Duodenal Bend	8.5	15	Each
M00532980	08714729786894	Duodenal Bend	8.5	18	Each
M00532990	08714729786900	Duodenal Bend	10	5	Each
M00533000	08714729786917	Duodenal Bend	10	7	Each
M00533010	08714729786924	Duodenal Bend	10	9	Each
M00533020	08714729786931	Duodenal Bend	10	12	Each
M00533030	08714729786948	Duodenal Bend	10	15	Each
M00533040	08714729786955	Duodenal Bend	10	18	Each

Scope Compatibility: 7Fr and 8.5Fr stents require a 3.2mm working channel; 10Fr stents require a 4.2mm working channel.



MR Magnetic Resonance Safe

Advanix Center Bend - Single Stent



Order Number	GTIN	Description	Stent Diameter (Fr)	Length (cm) Barb to Barb	Packaging
M00533230	08714729786603	Center Bend	7	5	Each
M00533240	08714729786610	Center Bend	7	7	Each
M00533250	08714729786627	Center Bend	7	9	Each
M00533260	08714729786634	Center Bend	7	12	Each
M00533270	08714729786641	Center Bend	7	15	Each
M00533280	08714729786658	Center Bend	7	18	Each
M00533290	08714729786665	Center Bend	8.5	5	Each
M00533300	08714729786672	Center Bend	8.5	7	Each
M00533310	08714729786689	Center Bend	8.5	9	Each
M00533320	08714729786696	Center Bend	8.5	12	Each
M00533330	08714729786702	Center Bend	8.5	15	Each
M00533340	08714729786719	Center Bend	8.5	18	Each
M00533350	08714729786726	Center Bend	10	5	Each
M00533360	08714729786733	Center Bend	10	7	Each
M00533370	08714729786740	Center Bend	10	9	Each
M00533380	08714729786757	Center Bend	10	12	Each
M00533390	08714729786764	Center Bend	10	15	Each
M00533400	08714729786771	Center Bend	10	18	Each

Scope Compatibility: 7Fr and 8.5Fr stents require a 3.2mm working channel; 10Fr stents require a 4.2mm working channel.



Section Table of Contents

Prev Page ◀ 21 ▶ Next Page

Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Biliary Plastic Stents

Advanix™ Biliary Stents (Continued)



MR Magnetic Resonance Safe

Advanix Double Pigtail - Single Stent



Order Number	GTIN	Description	Stent Diameter (Fr)	Length (cm) Between Pigtails	Packaging
M00532160	08714729786962	Double Pigtail	7	3	Each
M00532180	08714729786979	Double Pigtail	7	5	Each
M00532190	08714729786986	Double Pigtail	7	7	Each
M00532210	08714729786993	Double Pigtail	7	10	Each
M00532220	08714729787006	Double Pigtail	7	12	Each
M00532230	08714729787013	Double Pigtail	7	15	Each
M00532240	08714729787020	Double Pigtail	10	3	Each
M00532260	08714729787037	Double Pigtail	10	5	Each
M00532270	08714729787044	Double Pigtail	10	7	Each
M00532290	08714729787051	Double Pigtail	10	10	Each
M00532300	08714729787068	Double Pigtail	10	12	Each
M00532310	08714729787075	Double Pigtail	10	15	Each

Scope Compatibility: 7Fr and 8.5Fr stents require a 3.2mm working channel; 10Fr stents require a 4.2mm working channel.

Naviflex™ RX Delivery System



Naviflex RX Delivery System*



Order Number	GTIN	Description	Catheter Diameter (Fr)	Length (cm)	Packaging
M00533560	08714729787082	RX Delivery System	7	202.5	Each
M00533570	08714729787099	RX Delivery System	8.5	202.5	Each
M00533580	08714729787105	RX Delivery System	10	202.5	Each

*Compatible with Advanix Biliary Stents.



Section Table of Contents

Prev Page

22

Next Page

Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Biliary Plastic Stents

Duodenal Bend Stents

MR Magnetic Resonance Safe



RX Duodenal Bend Biliary Stent with RX Delivery System

Order Number	GTIN	Description	O.D. (Fr)	Length (cm) Barb to Barb	Recommended Guidewire (in)	Packaging
M00545550	08714729337850	Duodenal Bend Biliary Stent with RX Delivery System	7	5	.035	Each
M00545560	08714729337867	Duodenal Bend Biliary Stent with RX Delivery System	7	7	.035	Each
M00545570	08714729337874	Duodenal Bend Biliary Stent with RX Delivery System	7	10	.035	Each
M00545580	08714729337881	Duodenal Bend Biliary Stent with RX Delivery System	7	12	.035	Each
M00545590	08714729337898	Duodenal Bend Biliary Stent with RX Delivery System	7	15	.035	Each
M00545650	08714729252221	Duodenal Bend Biliary Stent with RX Delivery System	8.5	5	.035	Each
M00545660	08714729301837	Duodenal Bend Biliary Stent with RX Delivery System	8.5	7	.035	Each
M00545670	08714729301844	Duodenal Bend Biliary Stent with RX Delivery System	8.5	10	.035	Each
M00545680	08714729301851	Duodenal Bend Biliary Stent with RX Delivery System	8.5	12	.035	Each
M00545690	08714729301868	Duodenal Bend Biliary Stent with RX Delivery System	8.5	15	.035	Each
M00545600	08714729337904	Duodenal Bend Biliary Stent with RX Delivery System	10	5	.035	Each
M00545610	08714729337911	Duodenal Bend Biliary Stent with RX Delivery System	10	7	.035	Each
M00545620	08714729337928	Duodenal Bend Biliary Stent with RX Delivery System	10	10	.035	Each
M00545630	08714729301813	Duodenal Bend Biliary Stent with RX Delivery System	10	12	.035	Each
M00545640	08714729301820	Duodenal Bend Biliary Stent with RX Delivery System	10	15	.035	Each
M00545700	08714729301875	Duodenal Bend Biliary Stent with RX Delivery System	11.5	5	.035	Each
M00545710	08714729301882	Duodenal Bend Biliary Stent with RX Delivery System	11.5	7	.035	Each
M00545720	08714729301899	Duodenal Bend Biliary Stent with RX Delivery System	11.5	10	.035	Each
M00545730	08714729337942	Duodenal Bend Biliary Stent with RX Delivery System	11.5	12	.035	Each
M00545740	08714729337959	Duodenal Bend Biliary Stent with RX Delivery System	11.5	15	.035	Each



MR Magnetic Resonance Safe



Flexima™ Duodenal Bend Biliary Stent with Delivery System

Order Number	GTIN	O.D. (Fr)	Length (cm) Barb to Barb	Recommended Guidewire (in)	Packaging
M00539200	08714729337362	7	5	.035	Each
M00539210	08714729162520	7	7	.035	Each
M00539220	08714729162537	7	10	.035	Each
M00539230	08714729162544	7	12	.035	Each
M00539240	08714729162551	7	15	.035	Each
M00539250	08714729162568	8.5	5	.035	Each
M00539260	08714729162575	8.5	7	.035	Each
M00539270	08714729162582	8.5	10	.035	Each
M00539280	08714729162599	8.5	12	.035	Each
M00539290	08714729162605	8.5	15	.035	Each
M00539300	08714729162612	10	5	.035	Each
M00539310	08714729162629	10	7	.035	Each
M00539320	08714729162636	10	10	.035	Each
M00539330	08714729162643	10	12	.035	Each
M00539340	08714729162650	10	15	.035	Each
M00539350	08714729162667	11.5	5	.035	Each
M00539360	08714729162674	11.5	7	.035	Each
M00539370	08714729162681	11.5	10	.035	Each
M00539380	08714729162698	11.5	12	.035	Each
M00539390	08714729162704	11.5	15	.035	Each



Section Table of Contents

Prev Page 23 Next Page

Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Biliary Plastic Stents

Duodenal Bend Stents *(Continued)*



MR Magnetic Resonance Safe

Percuflex™ Duodenal Bend Biliary Stent



Order Number	GTIN	O.D. (Fr)	Length (cm) Barb to Barb	Recommended Guidewire (in)	Packaging
M00533610	08714729181217	7	5	.035	Each
M00533620	08714729181224	7	7	.035	Each
M00533630	08714729181231	7	10	.035	Each
M00533640	08714729181248	7	12	.035	Each
M00533650	08714729181255	7	15	.035	Each
M00533660	08714729181262	10	5	.035	Each
M00533670	08714729181279	10	7	.035	Each
M00533680	08714729181286	10	10	.035	Each
M00533690	08714729181293	10	12	.035	Each
M00533700	08714729181309	10	15	.035	Each

MR Magnetic Resonance Safe

Percuflex Duodenal Bend Biliary Stent with Introducer Kit



Order Number	GTIN	Description	O.D. (Fr)	Length (cm) Barb to Barb	Packaging
M00533760	08714729181088	Percuflex Stent w/ push catheter	7	5	Each
M00533770	08714729181095	Percuflex Stent w/ push catheter	7	7	Each
M00533780	08714729181101	Percuflex Stent w/ push catheter	7	10	Each
M00533790	08714729181118	Percuflex Stent w/ push catheter	7	12	Each
M00533800	08714729181125	Percuflex Stent w/ push catheter	7	15	Each
M00533810	08714729181132	Percuflex Stent w/ push and guide catheter	10	5	Each
M00533820	08714729181149	Percuflex Stent w/ push and guide catheter	10	7	Each
M00533830	08714729181156	Percuflex Stent w/ push and guide catheter	10	10	Each
M00533840	08714729193449	Percuflex Stent w/ push and guide catheter	10	12	Each
M00533850	08714729181163	Percuflex Stent w/ push and guide catheter	10	15	Each

Biliary Introducer*



Order Number	GTIN	Description	O.D. (Fr)	Recommended Guidewire (in)	Packaging
M00533910	08714729181439	Push catheter	7	.035	Each
M00533920	08714729193463	Push and guide catheter	10	.035	Each

*Stents sold separately.



Section Table of Contents

Prev Page ◀ 24 ▶ Next Page

Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Biliary Plastic Stents

Pigtail Stents




MR Magnetic Resonance Safe



C-Flex™ Double Pigtail Biliary Stent

Order Number	GTIN	O.D. (Fr)	Length (cm) Between Pigtails	Color Code	Packaging
M00532010	08714729019077	7	5	Blue	Each
M00532020	08714729019084	7	10	Blue	Each
M00532030	08714729019091	7	15	Blue	Each
M00532040	08714729019107	10	5	Grey	Each
M00532050	08714729019114	10	10	Grey	Each



MR Magnetic Resonance Safe



Flexima™ Single Pigtail Nasobiliary Catheter*

Order Number	GTIN	O.D. (Fr)	Working Length (cm)	Packaging
M00540110	08714729337379	5	250	Each
M00540120	08714729337386	6	250	Each
M00540130	08714729337393	7.5	250	Each
M00540140	08714729337409	8.5	250	Each

*Recommended Guidewire .035" Dreamwire™ Guidewire, Order #5614, [\(Click here to view page 3\)](#).



Section Table of Contents

Prev Page ◀ 25 ▶ Next Page

Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Pancreatic Plastic Stents

Advanix™ Pancreatic Stents



 Magnetic Resonance Conditional¹

Advanix Pancreatic Stent Kit - Straight, No Leading Barb



Order Number	GTIN	Description	Diameter (Fr)	Length (cm)	Packaging
M00537460	08714729856245	Straight Stent – No Leading Barb	4	3	Each
M00537470	08714729856252	Straight Stent – No Leading Barb	4	4	Each
M00537480	08714729856269	Straight Stent – No Leading Barb	4	5	Each
M00537490	08714729862956	Straight Stent – No Leading Barb	4	7	Each
M00537500	08714729856276	Straight Stent – No Leading Barb	5	3	Each
M00537510	08714729856283	Straight Stent – No Leading Barb	5	4	Each
M00537520	08714729856290	Straight Stent – No Leading Barb	5	5	Each
M00537530	08714729856306	Straight Stent – No Leading Barb	5	7	Each
M00537540	08714729856313	Straight Stent – Leading Barb	5	3	Each
M00537550	08714729856320	Straight Stent – Leading Barb	5	4	Each
M00537560	08714729856337	Straight Stent – Leading Barb	5	5	Each

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



 Magnetic Resonance Conditional¹

Advanix Pancreatic Stent Kit - Pigtail, Leading Barb and No Leading Barb



Order Number	GTIN	Description	Diameter (Fr)	Length (cm)	Packaging
M00537350	08714729856344	Single Pigtail Stent – No Leading Barb	4	3	Each
M00537360	08714729856351	Single Pigtail Stent – No Leading Barb	4	4	Each
M00537370	08714729856368	Single Pigtail Stent – No Leading Barb	4	5	Each
M00537380	08714729856375	Single Pigtail Stent – No Leading Barb	4	7	Each
M00537390	08714729856382	Single Pigtail Stent – No Leading Barb	5	3	Each
M00537400	08714729856399	Single Pigtail Stent – No Leading Barb	5	4	Each
M00537410	08714729856405	Single Pigtail Stent – No Leading Barb	5	5	Each
M00537420	08714729856412	Single Pigtail Stent – No Leading Barb	5	7	Each
M00537430	08714729856429	Single Pigtail Stent – Leading Barb	5	3	Each
M00537440	08714729856436	Single Pigtail Stent – Leading Barb	5	4	Each
M00537450	08714729856443	Single Pigtail Stent – Leading Barb	5	5	Each

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



Section Table of Contents

Prev Page ◀ 26 ▶ Next Page

Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Pancreatic Plastic Stents

Advanix™ Pancreatic Stents (Continued)



 Magnetic Resonance Conditional¹

Advanix Pancreatic Stent - Straight, No Leading Barb 					
Order Number	GTIN	Description	Diameter (Fr)	Length (cm)	Packaging
M00536000	08714729857686	Straight – No Leading Barb	3	2	Each
M00536010	08714729857693	Straight – No Leading Barb	3	3	Each
M00536020	08714729857709	Straight – No Leading Barb	3	4	Each
M00536030	08714729857716	Straight – No Leading Barb	3	5	Each
M00536040	08714729857723	Straight – No Leading Barb	3	6	Each
M00536050	08714729857730	Straight – No Leading Barb	3	7	Each
M00536060	08714729857747	Straight – No Leading Barb	3	8	Each
M00536070	08714729857754	Straight – No Leading Barb	3	9	Each
M00536080	08714729857761	Straight – No Leading Barb	3	10	Each
M00536090	08714729857778	Straight – No Leading Barb	3	11	Each
M00536100	08714729857785	Straight – No Leading Barb	3	12	Each
M00536110	08714729857792	Straight – No Leading Barb	3	13	Each
M00536120	08714729857068	Straight – No Leading Barb	4	2	Each
M00536130	08714729857075	Straight – No Leading Barb	4	3	Each
M00536140	08714729857082	Straight – No Leading Barb	4	4	Each
M00536150	08714729857099	Straight – No Leading Barb	4	5	Each
M00536160	08714729857105	Straight – No Leading Barb	4	6	Each
M00536170	08714729857112	Straight – No Leading Barb	4	7	Each
M00536180	08714729857129	Straight – No Leading Barb	4	8	Each
M00536190	08714729857136	Straight – No Leading Barb	4	9	Each
M00536200	08714729857143	Straight – No Leading Barb	4	10	Each
M00536210	08714729857150	Straight – No Leading Barb	4	11	Each
M00536220	08714729857167	Straight – No Leading Barb	4	12	Each
M00536230	08714729857174	Straight – No Leading Barb	4	13	Each
M00536240	08714729857181	Straight – No Leading Barb	5	2	Each
M00536250	08714729857198	Straight – No Leading Barb	5	3	Each
M00536260	08714729857204	Straight – No Leading Barb	5	4	Each
M00536270	08714729857211	Straight – No Leading Barb	5	5	Each
M00536280	08714729857228	Straight – No Leading Barb	5	6	Each
M00536290	08714729857235	Straight – No Leading Barb	5	7	Each
M00536300	08714729857242	Straight – No Leading Barb	5	8	Each
M00536310	08714729857259	Straight – No Leading Barb	5	9	Each
M00536320	08714729857266	Straight – No Leading Barb	5	10	Each
M00536330	08714729857273	Straight – No Leading Barb	5	11	Each
M00536340	08714729857280	Straight – No Leading Barb	5	12	Each
M00536350	08714729857297	Straight – No Leading Barb	5	13	Each

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



Section Table of Contents

Prev Page

27

Next Page

Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Pancreatic Plastic Stents

Advanix™ Pancreatic Stents (Continued)



 Magnetic Resonance
Conditional¹

Advanix Pancreatic Stent - Straight, Leading Barb 					
Order Number	GTIN	Description	Diameter (Fr)	Length (cm)	Packaging
M00536360	08714729857303	Straight – Leading Barb	5	2	Each
M00536370	08714729857310	Straight – Leading Barb	5	3	Each
M00536380	08714729857327	Straight – Leading Barb	5	4	Each
M00536390	08714729857334	Straight – Leading Barb	5	5	Each
M00536400	08714729857341	Straight – Leading Barb	5	6	Each
M00536410	08714729857358	Straight – Leading Barb	5	7	Each
M00536420	08714729857365	Straight – Leading Barb	5	8	Each
M00536430	08714729857372	Straight – Leading Barb	5	9	Each
M00536440	08714729857389	Straight – Leading Barb	5	10	Each
M00536450	08714729857396	Straight – Leading Barb	5	11	Each
M00536460	08714729857402	Straight – Leading Barb	5	12	Each
M00536470	08714729857419	Straight – Leading Barb	5	13	Each
M00536480	08714729857426	Straight – Leading Barb	7	3	Each
M00536490	08714729857433	Straight – Leading Barb	7	4	Each
M00536500	08714729857440	Straight – Leading Barb	7	5	Each
M00536510	08714729857457	Straight – Leading Barb	7	6	Each
M00536520	08714729857464	Straight – Leading Barb	7	7	Each
M00536530	08714729857471	Straight – Leading Barb	7	8	Each
M00536540	08714729857488	Straight – Leading Barb	7	9	Each
M00536550	08714729857495	Straight – Leading Barb	7	10	Each
M00536560	08714729857501	Straight – Leading Barb	7	11	Each
M00536570	08714729857518	Straight – Leading Barb	7	12	Each
M00536580	08714729857525	Straight – Leading Barb	7	13	Each
M00536590	08714729857532	Straight – Leading Barb	7	15	Each
M00536600	08714729857549	Straight – Leading Barb	7	18	Each
M00536610	08714729857556	Straight – Leading Barb	10	3	Each
M00536620	08714729857563	Straight – Leading Barb	10	4	Each
M00536630	08714729857570	Straight – Leading Barb	10	5	Each
M00536640	08714729857587	Straight – Leading Barb	10	6	Each
M00536650	08714729857594	Straight – Leading Barb	10	7	Each
M00536660	08714729857600	Straight – Leading Barb	10	8	Each
M00536670	08714729857617	Straight – Leading Barb	10	9	Each
M00536680	08714729857624	Straight – Leading Barb	10	10	Each
M00536690	08714729857631	Straight – Leading Barb	10	11	Each
M00536700	08714729857648	Straight – Leading Barb	10	12	Each
M00536710	08714729857655	Straight – Leading Barb	10	13	Each
M00536720	08714729857662	Straight – Leading Barb	10	15	Each
M00536730	08714729857679	Straight – Leading Barb	10	18	Each

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



Section Table of Contents

Prev Page ◀ 28 ▶ Next Page

Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Pancreatic Plastic Stents

Advanix™ Pancreatic Stents (Continued)



 Magnetic Resonance Conditional¹

Advanix Pancreatic Stent - Pigtail, No Leading Barb

Order Number	GTIN	Description	Diameter (Fr)	Length (cm)	Packaging
M00536740	08714729856948	Single Pigtail – No Leading Barb	3	2	Each
M00536750	08714729856955	Single Pigtail – No Leading Barb	3	3	Each
M00536760	08714729856962	Single Pigtail – No Leading Barb	3	4	Each
M00536770	08714729856979	Single Pigtail – No Leading Barb	3	5	Each
M00536780	08714729856986	Single Pigtail – No Leading Barb	3	6	Each
M00536790	08714729856993	Single Pigtail – No Leading Barb	3	7	Each
M00536800	08714729857006	Single Pigtail – No Leading Barb	3	8	Each
M00536810	08714729857013	Single Pigtail – No Leading Barb	3	9	Each
M00536820	08714729857020	Single Pigtail – No Leading Barb	3	10	Each
M00536830	08714729857037	Single Pigtail – No Leading Barb	3	11	Each
M00536840	08714729857044	Single Pigtail – No Leading Barb	3	12	Each
M00536850	08714729857051	Single Pigtail – No Leading Barb	3	13	Each
M00536860	08714729856450	Single Pigtail – No Leading Barb	4	2	Each
M00536870	08714729856467	Single Pigtail – No Leading Barb	4	3	Each
M00536880	08714729856474	Single Pigtail – No Leading Barb	4	4	Each
M00536890	08714729856481	Single Pigtail – No Leading Barb	4	5	Each
M00536900	08714729856498	Single Pigtail – No Leading Barb	4	6	Each
M00536910	08714729856504	Single Pigtail – No Leading Barb	4	7	Each
M00536920	08714729856511	Single Pigtail – No Leading Barb	4	8	Each
M00536930	08714729856528	Single Pigtail – No Leading Barb	4	9	Each
M00536940	08714729856535	Single Pigtail – No Leading Barb	4	10	Each
M00536950	08714729856542	Single Pigtail – No Leading Barb	4	11	Each
M00536960	08714729856559	Single Pigtail – No Leading Barb	4	12	Each
M00536970	08714729856566	Single Pigtail – No Leading Barb	4	13	Each
M00536980	08714729856696	Single Pigtail – No Leading Barb	5	2	Each
M00536990	08714729856702	Single Pigtail – No Leading Barb	5	3	Each
M00537000	08714729856719	Single Pigtail – No Leading Barb	5	4	Each
M00537010	08714729856726	Single Pigtail – No Leading Barb	5	5	Each
M00537020	08714729856733	Single Pigtail – No Leading Barb	5	6	Each
M00537030	08714729856740	Single Pigtail – No Leading Barb	5	7	Each
M00537040	08714729856757	Single Pigtail – No Leading Barb	5	8	Each
M00537050	08714729856764	Single Pigtail – No Leading Barb	5	9	Each
M00537060	08714729856771	Single Pigtail – No Leading Barb	5	10	Each
M00537070	08714729856788	Single Pigtail – No Leading Barb	5	11	Each
M00537080	08714729856795	Single Pigtail – No Leading Barb	5	12	Each
M00537090	08714729856801	Single Pigtail – No Leading Barb	5	13	Each

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



Section Table of Contents

Prev Page

29

Next Page

Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Pancreatic Plastic Stents

Advanix™ Pancreatic Stents (Continued)



 Magnetic Resonance Conditional¹

Advanix Pancreatic Stent - Pigtail, Leading Barb

Order Number	GTIN	Description	Diameter (Fr)	Length (cm)	Packaging
M00537100	08714729856573	Single Pigtail – Leading Barb	5	2	Each
M00537110	08714729856580	Single Pigtail – Leading Barb	5	3	Each
M00537120	08714729856597	Single Pigtail – Leading Barb	5	4	Each
M00537130	08714729856603	Single Pigtail – Leading Barb	5	5	Each
M00537140	08714729856610	Single Pigtail – Leading Barb	5	6	Each
M00537150	08714729856627	Single Pigtail – Leading Barb	5	7	Each
M00537160	08714729856634	Single Pigtail – Leading Barb	5	8	Each
M00537170	08714729856641	Single Pigtail – Leading Barb	5	9	Each
M00537180	08714729856658	Single Pigtail – Leading Barb	5	10	Each
M00537190	08714729856665	Single Pigtail – Leading Barb	5	11	Each
M00537200	08714729856672	Single Pigtail – Leading Barb	5	12	Each
M00537210	08714729856689	Single Pigtail – Leading Barb	5	13	Each
M00537220	08714729856818	Single Pigtail – Leading Barb	7	3	Each
M00537230	08714729856825	Single Pigtail – Leading Barb	7	4	Each
M00537240	08714729856832	Single Pigtail – Leading Barb	7	5	Each
M00537250	08714729856849	Single Pigtail – Leading Barb	7	6	Each
M00537260	08714729856856	Single Pigtail – Leading Barb	7	7	Each
M00537270	08714729856863	Single Pigtail – Leading Barb	7	8	Each
M00537280	08714729856870	Single Pigtail – Leading Barb	7	9	Each
M00537290	08714729856887	Single Pigtail – Leading Barb	7	10	Each
M00537300	08714729856894	Single Pigtail – Leading Barb	7	11	Each
M00537310	08714729856900	Single Pigtail – Leading Barb	7	12	Each
M00537320	08714729856917	Single Pigtail – Leading Barb	7	13	Each
M00537330	08714729856924	Single Pigtail – Leading Barb	7	15	Each
M00537340	08714729856931	Single Pigtail – Leading Barb	7	18	Each

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



Section Table of Contents

Prev Page

30

Next Page


Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Pancreatic Plastic Stents

Naviflex™ RX Pancreatic Delivery System & Naviflex Pushers



Naviflex RX Pancreatic Delivery System & Naviflex Pushers					
Order Number	GTIN	Description	Diameter (Fr)	Length	Packaging
M00535060	08714729856221	Naviflex RX Pancreatic Stent Delivery System	7	N/A	Each
M00535080	08714729856238	Naviflex RX Pancreatic Stent Delivery System	10	N/A	Each
M00535000	08714729856207	Pancreatic Stent RX Pusher	3	N/A	Each
M00535010	08714729856214	Pancreatic Stent RX Pusher	4 / 5	N/A	Each
M00535020	08714729856184	Pancreatic Stent Pusher	3	N/A	Each
M00535030	08714729856191	Pancreatic Stent Pusher	4 / 5	N/A	Each



Section Table of Contents

Prev Page

31

Next Page


Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Biliary Stone Management

RX Stone Management



Extractor™ Pro RX-S Retrieval Balloon Catheter – Distal Injection



Order Number	GTIN	Description	Balloon Inflated O.D. (mm)	Working Length (cm)	Catheter O.D. (Fr)	Recommended Guidewire (in)	Packaging
M00547300	08714729907138	Extractor Pro RX-S Retrieval Balloon	9.0 – 12.0	200	7 - 6	.035	Each
M00547310	08714729907152	Extractor Pro RX-S Retrieval Balloon	12.0 – 15.0	200	7 - 6	.035	Each
M00547320	08714729907176	Extractor Pro RX-S Retrieval Balloon	15.0 – 18.0	200	7 - 6	.035	Each



Extractor™ Pro RX-S Retrieval Balloon Catheter – Proximal Injection



Order Number	GTIN	Description	Balloon Inflated O.D. (mm)	Working Length (cm)	Catheter O.D. (Fr)	Recommended Guidewire (in)	Packaging
M00547330	08714729907145	Extractor Pro RX-S Retrieval Balloon	9.0 – 12.0	200	7 - 6	.035	Each
M00547340	08714729907169	Extractor Pro RX-S Retrieval Balloon	12.0 – 15.0	200	7 - 6	.035	Each
M00547350	08714729907183	Extractor Pro RX-S Retrieval Balloon	15.0 – 18.0	200	7 - 6	.035	Each



Extractor Pro RX Retrieval Balloon Catheter – Distal Injection



Order Number	GTIN	Description	Balloon Inflated O.D. (mm)	Working Length (cm)	Catheter O.D. (Fr)	Recommended Guidewire (in)	Packaging
M00547000	08714729790266	Extractor Pro RX Retrieval Balloon	9.0 – 12.0	200	7 - 6	.035	Each
M00547010	08714729790297	Extractor Pro RX Retrieval Balloon	12.0 – 15.0	200	7 - 6	.035	Each
M00547020	08714729790310	Extractor Pro RX Retrieval Balloon	15.0 – 18.0	200	7 - 6	.035	Each



Extractor Pro RX Retrieval Balloon Catheter – Proximal Injection




Order Number	GTIN	Description	Balloon Inflated O.D. (mm)	Working Length (cm)	Catheter O.D. (Fr)	Recommended Guidewire (in)	Packaging
M00547030	08714729790280	Extractor Pro RX Retrieval Balloon	9.0 – 12.0	200	7 - 6	.035	Each
M00547040	08714729790273	Extractor Pro RX Retrieval Balloon	12.0 – 15.0	200	7 - 6	.035	Each
M00547050	08714729790303	Extractor Pro RX Retrieval Balloon	15.0 – 18.0	200	7 - 6	.035	Each



Biliary Stone Management


Stone Management



Extractor™ Pro XL Retrieval Balloon Catheter – Distal Injection




Order Number	GTIN	Description	Balloon Inflated O.D. (mm)	Working Length (cm)	Catheter O.D. (Fr)	Recommended Guidewire (in)	Packaging
M00547100	08714729790389	Extractor Pro XL Retrieval Balloon	9.0 – 12.0	200	7 - 6	.035	Each
M00547110	08714729790402	Extractor Pro XL Retrieval Balloon	12.0 – 15.0	200	7 - 6	.035	Each
M00547120	08714729790648	Extractor Pro XL Retrieval Balloon	15.0 – 18.0	200	7 - 6	.035	Each



Extractor Pro XL Retrieval Balloon Catheter – Proximal Injection




Order Number	GTIN	Description	Balloon Inflated O.D. (mm)	Working Length (cm)	Catheter O.D. (Fr)	Recommended Guidewire (in)	Packaging
M00547130	08714729790396	Extractor Pro XL Retrieval Balloon	9.0 – 12.0	200	7 - 6	.035	Each
M00547140	08714729790365	Extractor Pro XL Retrieval Balloon	12.0 – 15.0	200	7 - 6	.035	Each
M00547150	08714729790372	Extractor Pro XL Retrieval Balloon	15.0 – 18.0	200	7 - 6	.035	Each



Extractor Pro DL Retrieval Balloon Catheter



Order Number	GTIN	Description	Balloon Inflated O.D. (mm)	Working Length (cm)	Catheter O.D. (Fr)	Recommended Guidewire (in)	Packaging
M00547200	08714729790631	Extractor Pro DL Retrieval Balloon	9.0 – 12.0	200	7 - 6	.035	Each
M00547210	08714729790655	Extractor Pro DL Retrieval Balloon	12.0 – 15.0	200	7 - 6	.035	Each
M00547220	08714729790662	Extractor Pro DL Retrieval Balloon	15.0 – 18.0	200	7 - 6	.035	Each



Trapezoid™ RX Wireguided Retrieval Baskets



Order Number	GTIN	Description	Open Basket Diameter (cm)	Recommended Guidewire (in)	Required Working Channel (mm)	Packaging
M00510860	08714729296379	Trapezoid RX Wireguided Retrieval Basket	1.5	.035	3.2	Each
M00510870	08714729296386	Trapezoid RX Wireguided Retrieval Basket	2.0	.035	3.2	Each
M00510880	08714729296393	Trapezoid RX Wireguided Retrieval Basket	2.5	.035	3.2	Each
M00510890*	08714729296409	Trapezoid RX Wireguided Retrieval Basket	3.0	.035	3.2	Each

Recommend Alliance™ II Inflation System for use as Mechanical Lithotripter; Order #5062, [\(Click here to view page 61\)](#).

*Refer to Instructions for Use when using this device; recommended for use with Biliary Calculi of >1.5cm.



Section Table of Contents

Prev Page

33

Next Page


Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Biliary Stone Management

CRE™ PRO Wireguided Balloon Dilatation Catheters



CRE PRO Wireguided Balloon Dilatation Catheters

Order Number	GTIN	Balloon Inflated O.D.		Balloon Length (cm)	Inflation Pressures		Minimum Working Channel (mm)	Catheter Size (Fr)	Working Length (cm)	Packaging
		(mm)	(Fr)		(atm)	(kPa)				
M00558600	08714729797487	6-7-8	18-21-24	5.5	3-6-10	304-608-1013	2.8 / 3.2*	7.5	180	Each
M00558610	08714729797494	8-9-10	24-27-30	5.5	3-5.5-9	304-557-912	2.8 / 3.2*	7.5	180	Each
M00558620	08714729797500	10-11-12	30-33-36	5.5	3-5-8	304-507-811	2.8 / 3.2*	7.5	180	Each
M00558630	08714729797517	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	2.8 / 3.2*	7.5	180	Each
M00558640	08714729797524	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	2.8 / 3.7*	7.5	180	Each
M00558650	08714729797531	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	2.8 / 3.7*	7.5	180	Each
M00558660	08714729797548	6-7-8	18-21-24	5.5	3-6-10	304-608-1013	2.8 / 3.2*	7.5	240	Each
M00558670	08714729797555	8-9-10	24-27-30	5.5	3-5.5-9	304-557-912	2.8 / 3.2*	7.5	240	Each
M00558680	08714729797562	10-11-12	30-33-36	5.5	3-5-8	304-507-811	2.8 / 3.2*	7.5	240	Each
M00558690	08714729797579	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	2.8 / 3.2*	7.5	240	Each
M00558700	08714729797586	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	2.8 / 3.7*	7.5	240	Each
M00558710	08714729797593	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	2.8 / 3.7*	7.5	240	Each

Recommended Inflation Device: Alliance™ II Integrated Inflation / Lithotripsy Device Order #5062, Alliance II Single-Use Syringe / Gauge Assembly Order #5060,

[\(Click here to view page 61\).](#)


Guidewire included and preloaded in balloon dilatation catheter.

*Duodenoscope channel size.



Biliary Stricture Management

RX Dilatation



Hurricane™ RX Biliary Balloon Dilatation Catheter



Order Number	GTIN	Balloon Inflated O.D. (mm)	Balloon Length (cm)	Catheter O.D. (Fr)	Catheter Length (cm)	Maximum Inflation Pressure (atm)	Minimum Channel Size (mm)	Recommended Guidewire (in)	Packaging
M00545890	08714729283805	4	2	5.8	180	12	3.2	.035	Each
M00545900	08714729283812	4	4	5.8	180	12	3.2	.035	Each
M00545910	08714729283829	6	2	5.8	180	11	3.2	.035	Each
M00545920	08714729283836	6	4	5.8	180	11	3.2	.035	Each
M00545930	08714729283843	8	2	5.8	180	11	3.2	.035	Each
M00545940	08714729283850	8	4	5.8	180	11	3.2	.035	Each
M00545950	08714729283867	10	2	5.8	180	8	3.2	.035	Each
M00545960	08714729283874	10	4	5.8	180	8	3.2	.035	Each

Recommend inflation device: ENCORE™ 26 Inflation Device Order #6667, [▶ \(Click here to view page 36\).](#)



Section Table of Contents

Prev Page

35

Next Page

Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Biliary Stricture Management

RX Cytology




RX Cytology Brush Wireguided Cytology Brush



Order Number	GTIN	Description	Brush Width (mm)	Catheter O.D. (Fr)	Working Length (cm)	Recommended Guidewire (in)	Packaging
M00545000	08714729268628	RX Cytology Brush	2.1	8	200	.035	Each

Biliary Dilatation




MaxForce™ Biliary Balloon Dilatation Catheter



Order Number	GTIN	Balloon Inflated O.D.		Balloon Length (cm)	Catheter Size		Working Length (cm)	Required Working Channel (mm)	Recommended Guidewire (in)	Packaging
		(mm)	(Fr)		(mm)	(Fr)				
M00567330	08714729097037	4	12	2	1.67	5	180	2.8	.035	Each
M00567340	08714729097051	4	12	4	1.67	5	180	2.8	.035	Each
M00567370	08714729097068	6	18	2	1.67	5	180	2.8	.035	Each
M00567380	08714729097075	6	18	4	1.67	5	180	2.8	.035	Each
M00567410	08714729097082	8	24	2	1.67	5	180	2.8	.035	Each
M00567420	08714729097099	8	24	3	1.67	5	180	2.8	.035	Each

RX Accessories



Encore™ 26 Inflation Device

Order Number	GTIN	Description	Packaging
M00566670	08714729755241	Encore 26 Inflation Device	Each



Notes and Future Product Updates



Notes and Future Product Updates



SpyGlass™ DS Direct Visualization System


➤ Capital Equipment.....	39
➤ Disposable Devices.....	39

SpyGlass Direct Visualization System

➤ Disposable and Reusable Devices.....	40
--	----

SpyGlass™ DS Direct Visualization System

Capital Equipment & Devices



Disposable Devices

Order Number	GTIN	Description	Packaging
M00546600	08714729863236	 SpyScope™ DS Access & Delivery Catheter	Each
M00546270	08714729764458	 SpyBite™ Biopsy Forceps	Each
M00546451	10677964060905	 SpyGlass Irrigation Tube Set	Box 10

Capital Equipment



Order Number	GTIN	Description	Packaging
M00546650	08714729874348	SpyGlass DS Digital Controller	Please contact your Boston Scientific Sales Representative for ordering information.

For more information on ordering a compatible EHL or Laser Probe for the SpyGlass DS Direct Visualization System, please contact your Boston Scientific Representative.



SpyGlass™ Direct Visualization System

Capital Equipment & Devices



Disposable and Reusable Devices¹

Order Number	GTIN	Description	Packaging
M00546030	08714729764441	SpyGlass Direct Visualization Probe	Each
M00546230	08714729764434	 SpyScope™ Access & Delivery Catheter	Each
M00546270	08714729764458	 SpyBite™ Biopsy Forceps	Each
M00546451	10677964060905	 SpyGlass Irrigation Tube Set	Box 10

¹ The SpyScope Access and Delivery Catheter, SpyBite Forceps and Irrigation Tubing are single-use devices. The SpyGlass Direct Visualization Probe is a reusable device that has been validated up to 20 reprocessing cycles of high level disinfection. The true number of procedural uses, however, will depend on how the SpyGlass Probe is handled. Fiberoptics are very sensitive and need delicate handling.

[Section Table of Contents](#)[Prev Page](#) ◀ 40 ▶ [Next Page](#)[Main Table of Contents](#)

Notes and Future Product Updates



Notes and Future Product Updates



EUS

Expect™ Endoscopic Ultrasound Aspiration Needles

- Expect Endoscopic Ultrasound Aspiration Needle 43
- Expect Slimline (SL) Endoscopic Ultrasound Aspiration Needle 43



Acquire™ Endoscopic Ultrasound Fine Needle Biopsy (FNB) Device

- Acquire Endoscopic Ultrasound Fine Needle Biopsy (FNB) Device..... 44

AXIOS™ System

- AXIOS Stent and Electrocautery Enhanced Delivery System ... 44
- AXIOS Stent and Delivery System..... 44

Expect™ Endoscopic Ultrasound Aspiration Needles



Expect Endoscopic Ultrasound Aspiration Needle*





Order Number	GTIN	Needle Size	Minimum Working Channel (mm)	Sheath Diameter (mm)	Packaging
M00550000	08714729790051	19 gauge	2.8	1.83	Each
M00550010	08714729790044	22 gauge	2.4	1.65	Each
M00550020	08714729790020	25 gauge	2.4	1.52	Each
M00550040	08714729797463	19 gauge Flex	2.8	1.73	Each
M00550001	08714729790068	19 gauge	2.8	1.83	Box 5
M00550011	08714729790013	22 gauge	2.4	1.65	Box 5
M00550021	08714729790037	25 gauge	2.4	1.52	Box 5
M00550041	08714729797470	19 gauge Flex	2.8	1.73	Box 5

*Packaging also includes: 20cc syringe and one-way stopcock

*Needle length: 0cm to 8cm, adjustable

*Working length: 137.5cm to 141.5cm, adjustable



Expect Slimline (SL) Endoscopic Ultrasound Aspiration Needle*



Order Number	GTIN	Needle Size	Minimum Working Channel (mm)	Sheath Diameter (mm)	Packaging
M00555500	08714729848141	19 gauge (1.10mm)	2.8	1.83	Each
M00555510	08714729848158	22 gauge (0.72mm)	2.4	1.65	Each
M00555520	08714729848165	25 gauge (0.52mm)	2.4	1.52	Each
M00555530	08714729848172	19 gauge Flex (1.14mm)	2.8	1.73	Each
M00555501	08714729848714	19 gauge (1.10mm)	2.8	1.83	Box 5
M00555511	08714729848721	22 gauge (0.72mm)	2.4	1.65	Box 5
M00555521	08714729848738	25 gauge (0.52mm)	2.4	1.52	Box 5
M00555531	08714729848806	19 gauge Flex (1.14mm)	2.8	1.73	Box 5

*Packaging also includes: 20cc syringe and one-way stopcock

*Needle length: 0cm to 8cm, adjustable

*Working length: 137.5cm to 141.5cm, adjustable

Acquire™ Endoscopic Ultrasound Fine Needle Biopsy (FNB) Device



Acquire Endoscopic Ultrasound Fine Needle Biopsy (FNB) Device*




Order Number	GTIN	Needle Size	Minimum Working Channel (mm)	Sheath Diameter (mm)	Packaging
M00555540	08714729931782	22 gauge (0.72mm)	2.4	1.65	Please contact your Boston Scientific Sales Representative for ordering information
M00555560	08714729931805	25 gauge (0.52mm)	2.4	1.52	Please contact your Boston Scientific Sales Representative for ordering information
M00555541	08714729931843	22 gauge (0.72mm)	2.4	1.65	Please contact your Boston Scientific Sales Representative for ordering information
M00555561	08714729931867	25 gauge (0.52mm)	2.4	1.52	Please contact your Boston Scientific Sales Representative for ordering information

*Packaging also includes: 20cc syringe and one-way stopcock

*Needle length: 0cm to 8cm, adjustable

*Working length: 137.5cm to 141.5cm, adjustable

AXIOS™ System



AXIOS Stent and Electrocautery Enhanced Delivery System



Order Number	Flange Diameter (mm)	Lumen Diameter (mm)	Saddle Length (mm)	Catheter O.D. (Fr)	Catheter Working Length (cm)	Catheter Total Length (cm)	Packaging
M00553640	21	10	10	10.8	138	146	Each
M00553650	24	15	10	10.8	138	146	Each

AXIOS Stent and Delivery System



Order Number	Flange Diameter (mm)	Lumen Diameter (mm)	Saddle Length (mm)	Catheter O.D. (Fr)	Catheter Working Length (cm)	Catheter Total Length (cm)	Packaging
M00553740	21	10	10	10.8	138	146	Each
M00553750	24	15	10	10.8	138	146	Each

Enteral & Esophageal Stents



Stents (Enteral and Esophageal)

Colonic and Duodenal Stents

➤ Ultraflex™ Precision Single-Use Colonic Stent System	46
➤ WallFlex™ Colonic Stent System with Anchor Lock Delivery System	45
➤ WallFlex Duodenal Stent System with Anchor Lock Delivery System	45
➤ WALLSTENT™ Enteral Endoprosthesis Colonic / Duodenal Stent with UNISTEP™ Plus Delivery System.....	46
➤ WALLSTENT Super Stiff Guidewire.....	46

Esophageal Stents

➤ Amplatz Super Stiff™ Guidewire.....	49
➤ Polyflex™ Single-Use Esophageal Stent System.....	48
➤ Ultraflex Single-Use Covered Esophageal NG Stent System – Distal Release	48
➤ Ultraflex Single-Use Covered Esophageal NG Stent System – Proximal Release.....	48
➤ Ultraflex Single-Use Covered Large Esophageal NG Stent System – Distal Release.....	49
➤ Ultraflex Single-Use Covered Large Esophageal NG Stent System – Proximal Release	49
➤ Ultraflex Single-Use Uncovered Esophageal NG Stent System – Distal Release.....	49
➤ Ultraflex Single-Use Uncovered Esophageal NG Stent System – Proximal Release	49
➤ WallFlex Fully Covered Esophageal Stent.....	47
➤ WallFlex Partially Covered Esophageal Stent.....	47



Enteral Stents

WallFlex™ Colonic and Duodenal Stents



 Magnetic Resonance Conditional¹

WallFlex Colonic Stent System with Anchor Lock Delivery System*



Order Number	GTIN	Stent Body Diameter (mm)	Stent Flare Diameter (mm)	Expanded Stent Length (mm)	Working Length (cm)	Delivery System Length (cm)	Delivery System Diameter (Fr)	Packaging
M00565040	08714729456513	25	30	60	230	270	10	Each
M00565050	08714729456520	25	30	90	230	270	10	Each
M00565060	08714729456537	25	30	120	230	270	10	Each
M00565070	08714729456544	25	30	60	135	175	10	Each
M00565080	08714729456551	25	30	90	135	175	10	Each
M00565090	08714729456568	25	30	120	135	175	10	Each
M00565100	08714729456575	22	27	60	230	270	10	Each
M00565110	08714729456582	22	27	90	230	270	10	Each
M00565120	08714729456599	22	27	120	230	270	10	Each
M00565130	08714729456605	22	27	60	135	175	10	Each
M00565140	08714729456612	22	27	90	135	175	10	Each
M00565150	08714729456629	22	27	120	135	175	10	Each

*Recommended Guidewire .035" WALLSTENT™ Super Stiff Guidewire, Order #18001, [\(Click here to view page 46\)](#).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



 Magnetic Resonance Conditional¹

WallFlex Duodenal Stent System with Anchor Lock Delivery System*



Order Number	GTIN	Stent Body Diameter (mm)	Flare Diameter (mm)	Expanded Stent Length (mm)	Working Length (cm)	Delivery System Length (cm)	Delivery System Diameter (Fr)	Packaging
M00565010	08714729456483	22	27	60	230	270	10	Each
M00565020	08714729456490	22	27	90	230	270	10	Each
M00565030	08714729456506	22	27	120	230	270	10	Each

*Recommended Guidewire .035" Hydra Jagwire™ Guidewire Stiff Shaft; or 450cm Order #5606, [\(Click here to view page 3\)](#).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



Section Table of Contents

Prev Page

45


Next Page

Main Table of Contents



Enteral Stents

WALLSTENT™ Colonic and Duodenal Stents (Continued)



 Magnetic Resonance Conditional¹

WALLSTENT Enteral Endoprosthesis Colonic / Duodenal Stent with UNISTEP™ Plus Delivery System*



Order Number	GTIN	Expanded Stent Diameter (mm)	Expanded Stent Length (mm)	Working Length (cm)	Delivery System Length (cm)	Delivery System Diameter (Fr)	Packaging
M00565560	08714729285052	20	60	230	255	10	Each
M00565570	08714729285069	20	90	230	255	10	Each
M00565580	08714729285076	22	60	230	255	10	Each
M00565590	08714729285083	22	90	230	255	10	Each
M00565620	08714729285113	20	60	135	160	10	Each
M00565630	08714729285120	20	90	135	160	10	Each
M00565640	08714729285137	22	60	135	160	10	Each
M00565650	08714729285144	22	90	135	160	10	Each



*Recommended Guidewire .035" WALLSTENT Super Stiff Guidewire, Order #18001, (see below).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

WALLSTENT Super Stiff Guidewire



Order Number	GTIN	Description	Distal Tip O.D. (in)	Length (cm)	Tip Style	Packaging
H965180011	08714729203070	WALLSTENT Super Stiff Guidewire	.035	500	Tapered	Box 5



 Magnetic Resonance Conditional¹

Ultraflex™ Precision Single-Use Colonic Stent System*



Order Number	GTIN	Stent Body Diameter (mm)	Stent Flare Diameter (mm)	Expanded Stent Length (mm)	Working Length (cm)	Delivery System Length (cm)	Delivery Tip Diameter (Fr)	Packaging
M00557330	08714729339175	25	30	57	100	105	22	Each
M00557360	08714729339199	25	30	87	100	105	22	Each
M00557380	08714729339205	25	30	117	100	105	22	Each

*Recommended Guidewire .038" Jagwire™ Guidewire, Order #5662, [\(Click here to view page 4\)](#).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



Section Table of Contents

Prev Page

46


Next Page

Main Table of Contents



Esophageal Stents

WallFlex™ Esophageal Stents



 Magnetic Resonance Conditional¹

WallFlex Fully Covered Esophageal Stent*



Order Number	GTIN	Stent O.D. (mm)	Proximal / Distal Flares O.D. (mm)	Stent Length (cm)	Catheter Diameter (Fr / mm)	Working Length (cm)	Packaging
M00516700	08714729778035	18	25 / 23	10.3	18.5 / 6.2	78	Each
M00516710	08714729778042	18	25 / 23	12.3	18.5 / 6.2	78	Each
M00516720	08714729778059	18	25 / 23	15.3	18.5 / 6.2	78	Each
M00516730	08714729778066	23	28 / 28	10.5	18.5 / 6.2	78	Each
M00516740	08714729778073	23	28 / 28	12.5	18.5 / 6.2	78	Each
M00516750	08714729778080	23	28 / 28	15.5	18.5 / 6.2	78	Each

*Recommended Guidewire .038" Jagwire™ Guidewire, Order #5662, [\(Click here to view page 4\)](#).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



 Magnetic Resonance Conditional¹

WallFlex Partially Covered Esophageal Stent*



Order Number	GTIN	Stent O.D. (mm)	Proximal / Distal Flares O.D. (mm)	Stent Length (cm)	Covered Length (cm)	Catheter Diameter (Fr / mm)	Working Length (cm)	Packaging
M00516900	08714729765240	18	23 / 23	10.3	7.1	18.5 / 6.2	78	Each
M00516910	08714729765257	18	23 / 23	12.3	9.1	18.5 / 6.2	78	Each
M00516920	08714729765264	18	23 / 23	15.3	12.1	18.5 / 6.2	78	Each
M00516930	08714729765271	23	28 / 28	10.5	7.3	18.5 / 6.2	78	Each
M00516940	08714729765288	23	28 / 28	12.5	9.3	18.5 / 6.2	78	Each
M00516950	08714729765295	23	28 / 28	15.5	12.3	18.5 / 6.2	78	Each

*Recommended Guidewire .038" Jagwire Guidewire, Order #5662, [\(Click here to view page 4\)](#).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



Section Table of Contents


Prev Page  47  Next Page

Main Table of Contents 

For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Esophageal Stents

Polyflex™ Stents




Polyflex™ Single-Use Esophageal Stent System*



Order Number	GTIN	Stent I.D. (mm)	Proximal Flare I.D. (mm)	Stent Length (mm)	Delivery System Diameter (mm)	Packaging
M00514270	08714729457015	16	20	90	12	Each
M00514280	08714729457022	16	20	120	12	Each
M00514290	08714729457039	16	20	150	12	Each
M00514300	08714729457046	18	23	90	13	Each
M00514310	08714729457053	18	23	120	13	Each
M00514320	08714729457060	18	23	150	13	Each
M00514330	08714729457077	21	25	90	14	Each
M00514340	08714729457084	21	25	120	14	Each
M00514350	08714729457091	21	25	150	14	Each

*Recommended Guidewire .038" Jagwire™ Guidewire, Order #5662, [\(Click here to view page 4\)](#).

Ultraflex™ Stents



Ultraflex Single-Use Covered Esophageal NG Stent System – Distal Release*



Order Number	GTIN	Stent O.D. (mm)	Proximal Flare O.D. (mm)	Stent Length (cm)	Covered Length (cm)	Packaging
M00513730	08714729716150	18	23	10	7	Each
M00513740	08714729716167	18	23	12	9	Each
M00513750	08714729649113	18	23	15	12	Each

*Recommended Guidewire .038" Jagwire Guidewire, Order #5662, [\(Click here to view page 4\)](#).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



Ultraflex Single-Use Covered Esophageal NG Stent System – Proximal Release*



Order Number	GTIN	Stent O.D. (mm)	Proximal Flare O.D. (mm)	Stent Length (cm)	Covered Length (cm)	Packaging
M00513840	08714729716099	18	23	10	7	Each
M00513850	08714729716075	18	23	12	9	Each
M00513860	08714729716082	18	23	15	12	Each


*Recommended Guidewire .038" Jagwire Guidewire, Order #5662, [\(Click here to view page 4\)](#).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



Esophageal Stents

Ultraflex™ Stents (Continued)



 Magnetic Resonance Conditional¹

Ultraflex Single-Use Covered Large Esophageal NG Stent System – Distal Release*



Order Number	GTIN	Stent O.D. (mm)	Proximal Flare O.D. (mm)	Stent Length (cm)	Covered Length (cm)	Packaging
M00514200	08714729121190	23	28	10	7	Each
M00514210	08714729121206	23	28	12	9	Each

*Recommended Guidewire .038" Jagwire™ Guidewire, Order #5662, [\(Click here to view page 4\)](#).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

 Magnetic Resonance Conditional¹


Ultraflex Single-Use Covered Large Esophageal NG Stent System – Proximal Release*



Order Number	GTIN	Stent O.D. (mm)	Proximal Flare O.D. (mm)	Stent Length (cm)	Covered Length (cm)	Packaging
M00514240	08714729716105	23	28	10	7	Each
M00514250	08714729716112	23	28	12	9	Each

*Recommended Guidewire .038" Jagwire Guidewire, Order #5662, [\(Click here to view page 4\)](#).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



 Magnetic Resonance Conditional¹

Ultraflex Single-Use Uncovered Esophageal NG Stent System – Distal Release*



Order Number	GTIN	Stent O.D. (mm)	Proximal Flare O.D. (mm)	Stent Length (cm)	Packaging
M00513700	08714729716129	18	23	7	Each
M00513710	08714729716136	18	23	10	Each
M00513720	08714729716143	18	23	15	Each

*Recommended Guidewire .038" Jagwire Guidewire, Order #5662, [\(Click here to view page 4\)](#).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

 Magnetic Resonance Conditional¹

Ultraflex Single-Use Uncovered Esophageal NG Stent System – Proximal Release*



Order Number	GTIN	Stent O.D. (mm)	Proximal Flare O.D. (mm)	Stent Length (cm)	Packaging
M00513800	08714729716174	18	23	7	Each
M00513810	08714729716181	18	23	10	Each
M00513830	08714729716204	18	23	15	Each

*Recommended Guidewire .038" Jagwire Guidewire, Order #5662, [\(Click here to view page 4\)](#).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

Amplatz Super Stiff™ Guidewire



Order Number	GTIN	Description	O.D. (in)	Length (cm)	Tip Style	Packaging
M00550090	08714729301929	Amplatz Super Stiff Guidewire	.038	260	Straight	Each

Notes and Future Product Updates



Hemostasis

Clipping Devices

- ▶ Resolution 360™ Clip 51
- ▶ Resolution™ Clip 51

Bipolar Electrohemostasis Catheters

- ▶ Gold Probe™ Single-Use Electrohemostasis Catheters 52
- ▶ Injection Gold Probe™ Single-Use
Electrohemostasis / Injection Catheters 52

Needles

- ▶ Interject™ Clear Single-Use Injection Therapy
Needle Catheters 53
- ▶ Interject Contrast Single-Use Injection Therapy
Needle Catheters 53
- ▶ Optiflo™ Single-Use Hemostasis Catheters 53

Multiple Band Ligators

- ▶ Speedband Superview Super 7™ Multiple Band Ligators 54

Generator Accessories

- ▶ Endostat™ Generator Accessories 54

Electrosurgical Bipolar Cable Adapter

- ▶ Electrosurgical Bipolar Cable Adapter 55

Hemostasis

Clipping Devices



 Magnetic Resonance Conditional¹

Resolution 360™ Clip



Order Number	GTIN	Working Length (cm)	Minimum Working Channel (mm)	Clip Opening (mm)	Packaging
M00521240	08714729875666	155*	2.8	11	Each
M00521241	08714729875673	155*	2.8	11	Box 10
M00521242	08714729875680	155*	2.8	11	Box 20
M00521230	08714729875628	235	2.8	11	Each
M00521231	08714729875635	235	2.8	11	Box 10
M00521232	08714729875642	235	2.8	11	Box 20

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

*Product not yet for sale.



 Magnetic Resonance Conditional¹

Resolution™ Clip



Order Number	GTIN	Working Length (cm)	Minimum Working Channel (mm)	Clip Opening (mm)	Packaging
M00522600	08714729504764	155	2.8	11	Each
M00522601	08714729504771	155	2.8	11	Box 10
M00522602	08714729504788	155	2.8	11	Box 20
M00522610	08714729504795	235	2.8	11	Each
M00522611	08714729504801	235	2.8	11	Box 10
M00522612	08714729504818	235	2.8	11	Box 20

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



Section Table of Contents

Prev Page

51


Next Page

Main Table of Contents



Hemostasis

Bipolar Electrohemostasis Catheters




Injection Gold Probe™ Single-Use Electrohemostasis / Injection Catheters



Order Number	GTIN	Description	O.D. (Fr)	Needle Size Gauge	Working Length (cm)	Type Generator / Connector	Packaging
M00560150	08714729056492	Injection Gold Probe Catheter	7	25	210	Standard Plug*	Each
M00560160	08714729056508	Injection Gold Probe Catheter	10	25	210	Standard Plug*	Each

*Compatible with Endostat™ Bipolar-Monopolar Electrosurgical Generators.

For connections to other generators, please refer to the Electrosurgical Bipolar Cable Adapters, [▶](#) (Click here to view page 55).



Gold Probe™ Single-Use Electrohemostasis Catheters



Order Number	GTIN	Description	O.D. (Fr)	Working Length (cm)	Type Generator / Connector	Packaging
M00560070	08714729174691	Gold Probe Catheter	7	300	Standard Plug*	Each
M00560100	08714729174684	Gold Probe Catheter	10	300	Standard Plug*	Each
M00560220	08714729099574	Gold Probe Catheter	7	350	Standard Plug*	Each
M00560071	08714729071013	Gold Probe Catheter	7	300	Standard Plug*	Box 5
M00560101	08714729071020	Gold Probe Catheter	10	300	Standard Plug*	Box 5


*Compatible with Endostat Bipolar-Monopolar Electrosurgical Generators.

For connections to other generators, please refer to the Electrosurgical Bipolar Cable Adapters, [▶](#) (Click here to view page 55).



Hemostasis

Needles




Interject™ Contrast Single-Use Injection Therapy Needle Catheters*



Order Number	GTIN	Sheath Design	Needle Gauge	Maximum Needle Extension Length (mm)	Sheath O.D. (mm)	Working Length (cm)	Packaging
M00518151	08714729296430	Contrast	23	4	1.8	200	Box 5
M00518161	08714729296447	Contrast	25	4	1.8	200	Box 5
M00518251	08714729296478	Contrast	23	6	1.8	200	Box 5
M00518261	08714729296485	Contrast	25	6	1.8	200	Box 5
M00518351	08714729296522	Contrast	23	4	2.3	240	Box 5
M00518361	08714729296539	Contrast	25	4	2.3	240	Box 5

*With Star Catheter™ Technology.




Interject Clear Single-Use Injection Therapy Needle Catheters*



Order Number	GTIN	Sheath Design	Needle Gauge	Maximum Needle Extension Length (mm)	Sheath O.D. (mm)	Working Length (cm)	Packaging
M00518101	08714729296416	Clear	23	4	1.8	200	Box 5
M00518111	08714729296423	Clear	25	4	1.8	200	Box 5
M00518301	08714729296508	Clear	23	4	2.3	240	Box 5
M00518311	08714729296515	Clear	25	4	2.3	240	Box 5

*With Star Catheter Technology.



Optiflo™ Single-Use Hemostasis Catheters



Order Number	GTIN	Description	Needle Gauge	Sheath O.D. (mm)	Working Length (cm)	Working Channel (mm)	Packaging
M00511731	08714729144182	Optiflo Catheter	23	2.3	200	2.8	Box 5
M00511751	08714729145264	Optiflo Catheter	25	2.3	200	2.8	Box 5



Section Table of Contents

Prev Page

53

Next Page



Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Hemostasis

Multiple Band Ligators



Speedband Superview Super 7™ Multiple Band Ligators



Order Number	GTIN	Description	Required Working Channel (mm)	Scope O.D. (mm)	Bands / Ligating Unit	Packaging
M00542251	08714729201960	Speedband Superview Super 7 Band Ligator	2.8	8.6 – 11.5	7	Box 2
M00542253	08714729504832	Speedband Superview Super 7 Band Ligator	2.8	8.6 – 11.5	7	Box 4

Generator Accessories


Endostat™ Generator Accessories

Order Number	GTIN	Description	Packaging
M00540720	08714729337485	Water Bottle	Each
M00540240	08714729804925	Endostat III Dual Footswitch	Each
M00540760	08714729337522	Power Cord w/ I.E.C. Molded Plug	Each
M00540771	08714729337539	E-Z Load Pump Tubing – Endostat III Generator (for generators w/ E-Z Load Pump)	Box 5



Hemostasis

Electrosurgical Bipolar Cable Adapter



Electrosurgical Bipolar Cable Adapter*

Order Number	GTIN	Description	Pin Spacing (mm)	Generator Compatibility**	Packaging
M00561440	08714729797609	Bipolar Cable Adapter	18	EMC (Coagulasem Dolley)	Each
M00561450	08714729797616	Bipolar Cable Adapter	22	ERBE, EMC Dolley, Olympus™ (ESG100 & PSD60)	Each
M00561460	08714729797623	Bipolar Cable Adapter	25	Lamidey (Optima 4)	Each
M00561470	08714729797630	Bipolar Cable Adapter	28	Valleylab™, Lamidey (Duo), Bovie™ (Aaron 3250™), Olympus (UES10, UES30 & Surgemaster), Bowa, Soring	Each
M00561480	08714729797647	Bipolar Cable Adapter	28 / Third Pin	ConMed™ (System 5000), Bovie (Aaron 2250™), Ethicon™	Each
M00561490	08714729797654	Bipolar Cable Adapter	32 / Third Pin	ConMed (Sable 180, Excalibur, System 7550)	Each
M00561500	08714729797661	Bipolar Cable Adapter	39	Codman, Malis™	Each

*Bipolar Cable Adapters are used to connect Gold Probe™ and Injection Gold Probe™ Catheters with the above mentioned generators.

Bipolar Cable Adapters are not required for the connection of Gold Probe and Injection Gold Probe to Endostat™ II or III Generators.

** For generators and their associated cable adaptors not listed above, please contact your Boston Scientific Representative.

Note: Monopolar active cords are in the Polypectomy Section, [Click here to view page 77](#).

Note: Some variations in generator faceplates may exist within a generator model. We recommend that you confirm pin spacing prior to ordering.

Notes and Future Product Updates



Dilation

CRE™ Balloon Dilatation Catheter

- CRE PRO Wireguided Esophageal / Pyloric / Biliary Balloon Dilatation Catheters..... 57
- CRE PRO Wireguided Esophageal / Pyloric / Colonic / Biliary Balloon Dilatation Catheters..... 57
- CRE Wireguided Esophageal / Pyloric Balloon Dilatation Catheters..... 58
- CRE Wireguided Esophageal / Pyloric / Colonic Balloon Dilatation Catheters..... 58
- CRE Fixed Wire Esophageal Balloon Dilatation Catheters 59

MaxForce TTS™ Esophageal High Performance Balloon Dilatation Catheter

- MaxForce TTS Esophageal High Performance Balloon Dilatation Catheters..... 60

Rigiflex™ II Single-Use Achalasia Balloon

- Rigiflex II Single-Use Achalasia Balloon 60

Inflation System


- Alliance™ II Integrated Inflation / Lithotripsy Device..... 61
- Alliance II Single-Use Syringe / Gauge Assembly 61
- CRE SteriFlate Disposable Inflation Device..... 61

Dilation Accessories

- Inflation / Deflation Devices and Accessories 61

Dilation

CRE™ PRO Balloon Dilatation Catheters



CRE PRO Wireguided Esophageal / Pyloric / Biliary Balloon Dilatation Catheters



Order Number	GTIN	Balloon Inflated O.D.		Balloon Length (cm)	Inflation Pressures		Minimum Working Channel (mm)	Catheter Size (Fr)	Working Length (cm)	Packaging
		(mm)	(Fr)		(atm)	(kPa)				
M00558600	08714729797487	6-7-8	18-21-24	5.5	3-6-10	304-608-1013	2.8 / 3.2*	7.5	180	Each
M00558610	08714729797494	8-9-10	24-27-30	5.5	3-5-9	304-557-912	2.8 / 3.2*	7.5	180	Each
M00558620	08714729797500	10-11-12	30-33-36	5.5	3-5-8	304-507-811	2.8 / 3.2*	7.5	180	Each
M00558630	08714729797517	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	2.8 / 3.2*	7.5	180	Each
M00558640	08714729797524	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	2.8 / 3.7*	7.5	180	Each
M00558650**	08714729797531	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	2.8 / 3.7*	7.5	180	Each

Recommended Inflation Device: Alliance™ II Integrated Inflation / Lithotripsy Device Order #5062, Alliance II Syringe / Gauge Assembly Order #5060,

[\(Click here to view page 61\).](#)

Guidewire included and preloaded in balloon dilatation catheter.

*Duodenoscope channel size.

**Esophageal and Biliary application only.



CRE PRO Wireguided Esophageal / Pyloric / Colonic / Biliary Balloon Dilatation Catheters



Order Number	GTIN	Balloon Inflated O.D.		Balloon Length (cm)	Inflation Pressures		Minimum Working Channel (mm)	Catheter Size (Fr)	Working Length (cm)	Packaging
		(mm)	(Fr)		(atm)	(kPa)				
M00558660	08714729797548	6-7-8	18-21-24	5.5	3-6-10	304-608-1013	2.8 / 3.2*	7.5	240	Each
M00558670	08714729797555	8-9-10	24-27-30	5.5	3-5-9	304-557-912	2.8 / 3.2*	7.5	240	Each
M00558680	08714729797562	10-11-12	30-33-36	5.5	3-5-8	304-507-811	2.8 / 3.2*	7.5	240	Each
M00558690	08714729797579	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	2.8 / 3.2*	7.5	240	Each
M00558700	08714729797586	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	2.8 / 3.7*	7.5	240	Each
M00558710**	08714729797593	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	2.8 / 3.7*	7.5	240	Each

Recommended Inflation Device: Alliance II Integrated Inflation / Lithotripsy Device Order #5062, Alliance II Syringe / Gauge Assembly Order #5060,

[\(Click here to view page 61\).](#)

Guidewire included and preloaded in balloon dilatation catheter.

*Duodenoscope channel size.

**Esophageal, Colonic and Biliary application only.



Section Table of Contents

Prev Page

57


Next Page

Main Table of Contents



Dilation

CRE™ Balloon Dilatation Catheters



CRE Wireguided Esophageal / Pyloric Balloon Dilatation Catheters



Order Number	GTIN	Balloon Inflated O.D.		Balloon Length (cm)	Inflation Pressures		Minimum Working Channel (mm)	Catheter Size (Fr)	Working Length (cm)	Packaging
		(mm)	(Fr)		(atm)	(kPa)				
M00558390	08714729339328	6-7-8	18-21-24	5.5	3-6-10	304-608-1013	2.8	7.5	180	Each
M00558400	08714729339335	8-9-10	24-27-30	5.5	3-5.5-9	304-557-912	2.8	7.5	180	Each
M00558410	08714729339342	10-11-12	30-33-36	5.5	3-5-8	304-507-811	2.8	7.5	180	Each
M00558420	08714729339359	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	2.8	7.5	180	Each
M00558430	08714729339366	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	2.8	7.5	180	Each
M00558440*	08714729339373	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	2.8	7.5	180	Each

Recommended Inflation Device: Alliance™ II Integrated Inflation / Lithotripsy Device Order #5062, Alliance II Single-Use Syringe / Gauge Assembly Order #5060, [\(Click here to view page 61\).](#)

Guidewire included and preloaded in balloon dilatation catheter.

*Esophageal application only.

CRE Wireguided Esophageal / Pyloric / Colonic Balloon Dilatation Catheters



Order Number	GTIN	Balloon Inflated O.D.		Balloon Length (cm)	Inflation Pressures		Minimum Working Channel (mm)	Catheter Size (Fr)	Working Length (cm)	Packaging
		(mm)	(Fr)		(atm)	(kPa)				
M00558450	08714729339380	6-7-8	18-21-24	5.5	3-6-10	304-608-1013	2.8	7.5	240	Each
M00558460	08714729339397	8-9-10	24-27-30	5.5	3-5.5-9	304-557-912	2.8	7.5	240	Each
M00558470	08714729339403	10-11-12	30-33-36	5.5	3-5-8	304-507-811	2.8	7.5	240	Each
M00558480	08714729339410	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	2.8	7.5	240	Each
M00558490	08714729339427	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	2.8	7.5	240	Each
M00558500*	08714729339434	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	2.8	7.5	240	Each

Recommended Inflation Device: Alliance II Integrated Inflation / Lithotripsy Device Order #5062, Alliance II Single-Use Syringe / Gauge Assembly Order #5060, [\(Click here to view page 61\).](#)

Guidewire included and preloaded in balloon dilatation catheter.

*Esophageal and Colonic application only.



Section Table of Contents

Prev Page **58** Next Page


Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Dilation

CRE™ Balloon Dilatation Catheters (Continued)



CRE Fixed Wire Esophageal Balloon Dilatation Catheters



Order Number	GTIN	Balloon Inflated O.D.		Balloon Length (cm)	Inflation Pressures		Minimum Working Channel (mm)	Catheter Size (Fr)	Working Length (cm)	Packaging
		(mm)	(Fr)		(atm)	(kPa)				
M00558330	08714729201984	6-7-8	18-21-24	8	3-6-10	304-608-1013	2.8	6	180	Each
M00558340	08714729201991	8-9-10	24-27-30	8	3-5.5-9	304-557-912	2.8	6	180	Each
M00558350	08714729202004	10-11-12	30-33-36	8	3-5-8	304-507-811	2.8	6	180	Each
M00558360	08714729195986	12-13.5-15	36-40.5-45	8	3-4.5-8	304-456-811	2.8	6	180	Each
M00558370	08714729195993	15-16.5-18	45-49.5-54	8	3-4.5-7	304-456-709	2.8	6	180	Each
M00558380	08714729202011	18-19-20	54-57-60	8	3-4.5-6	304-456-608	2.8	6	180	Each
M00558361	08714729339298	12-13.5-15	36-40.5-45	8	3-4.5-8	304-456-811	2.8	6	180	Box 5
M00558371	08714729339304	15-16.5-18	45-49.5-54	8	3-4.5-7	304-456-709	2.8	6	180	Box 5
M00558381	08714729339311	18-19-20	54-57-60	8	3-4.5-6	304-456-608	2.8	6	180	Box 5

Recommended Inflation Device: Alliance™ II Integrated Inflation / Lithotripsy Device Order #5062, Alliance II Single-Use Syringe / Gauge Assembly Order #5060,

[\(Click here to view page 61\).](#)



Section Table of Contents

Prev Page

59

Next Page


Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Dilation

MaxForce TTS™ Esophageal High Performance Balloon Dilatation Catheter



MaxForce TTS™ Esophageal High Performance Balloon Dilatation Catheter




Order Number	GTIN	Balloon Inflated O.D.		Balloon Length (cm)	Inflation Pressures		Minimum Working Channel (mm)	Catheter Size		Working Length (cm)	Packaging
		(mm)	(Fr)		(atm)	(kPa)		(mm)	(Fr)		
M00589060	08714729099581	6	18	6	10	1013	2.8	2.0	6	180	Each
M00589080	08714729099598	8	24	6	10	1013	2.8	2.0	6	180	Each
M00589100	08714729099895	10	30	6	7	709	2.8	2.0	6	180	Each
M00589120	08714729099901	12	36	6	7	709	2.8	2.0	6	180	Each
M00589140	08714729080602	14	42	6	7	709	2.8	2.0	6	180	Each
M00589150	08714729099918	15	45	6	5	507	2.8	2.0	6	180	Each
M00589160	08714729080619	16	48	6	5	507	2.8	2.0	6	180	Each
M00589180	08714729099925	18	54	6	5	507	2.8	2.0	6	180	Each

Recommended Inflation Device: Alliance™ II Integrated Inflation / Lithotripsy Device Order #5062, Sterile Syringe Order #5060,

[\(Click here to view page 61\).](#)

Rigiflex™ II Single-Use Achalasia Balloon Dilatation Catheter



Rigiflex II Single-Use Achalasia Balloon Dilators*



Order Number	GTIN	Balloon Inflated O.D.		Balloon Length (cm)	Inflation Pressures		Working Length (cm)	Catheter Size		Recommended Guidewire** (in)	Packaging
		(mm)	(Fr)		(psi)	(kPa)		(mm)	(Fr)		
M00554500	08714729719106	30	90	10	20	138	90	4.6	14	.038	Each
M00554510	08714729719113	35	105	10	20	138	90	4.6	14	.038	Each
M00554520	08714729719120	40	120	10	20	138	90	4.6	14	.038	Each

*Product Shipped Sterile. Recommended use with Achalasia Pneumatic Hand Pump and Monitor Order #5320,

[\(Click here to view page 61\).](#)

**Guidewire included with Rigiflex II Balloon Dilatation Catheter.



Section Table of Contents

Prev Page **60** Next Page


Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Dilation

Inflation System



Alliance™ II Single-Use Syringe / Gauge Assembly / Integrated Inflation / Lithotripsy Device

Order Number	GTIN	Description	Packaging
M00550620	08714729283881	Alliance II Integrated Inflation / Lithotripsy Device	Each
M00550601	08714729129332	 Alliance II Single-Use Syringe / Gauge Assembly	Box 5

Dilation Accessories

Inflation / Deflation Devices and Accessories

Order Number	GTIN	Description	Packaging
M00553200	08714729195450	Pneumatic Inflator	Each
M00550630	00856877002506	CRE™ SteriFlate Disposable Inflation Device	Each



Notes and Future Product Updates




Biopsy

Forceps

➤ Multibite™ Single-Use Multiple Sample Biopsy Forceps	65	➤ Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps.....	63
➤ Radial Jaw™ 4 Single-Use Hot Biopsy Forceps	64	➤ Radial Jaw 4 Single-Use Pediatric Biopsy Forceps.....	64
➤ Radial Jaw 4 Single-Use Jumbo Biopsy Forceps.....	63	➤ Radial Jaw 4 Standard Capacity Single-Use Biopsy Forceps...	63

Biopsy

Forceps



Radial Jaw™ 4 Single-Use Jumbo Biopsy Forceps



Order Number	GTIN	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Color Code	Packaging
M00513343	08714729761891	Radial Jaw 4 Jumbo Forceps	2.8	160	3.2	Yellow	Box 5
M00513353	08714729761907	Radial Jaw 4 Jumbo Forceps w/ Needle	2.8	160	3.2	Yellow	Box 5
M00513363	08714729761914	Radial Jaw 4 Jumbo Forceps	2.8	240	3.2	Orange	Box 5
M00513373	08714729761938	Radial Jaw 4 Jumbo Forceps w/ Needle	2.8	240	3.2	Orange	Box 5
M00513361	08714729648871	Radial Jaw 4 Jumbo Forceps	2.8	240	3.2	Orange	Box 20
M00513371	08714729648895	Radial Jaw 4 Jumbo Forceps w/ Needle	2.8	240	3.2	Orange	Box 20
M00513362	08714729648888	Radial Jaw 4 Jumbo Forceps	2.8	240	3.2	Orange	Box 40
M00513372	08714729648901	Radial Jaw 4 Jumbo Forceps w/ Needle	2.8	240	3.2	Orange	Box 40

Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps



Order Number	GTIN	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Color Code	Packaging
M00513303	08714729761839	Radial Jaw 4 Large Capacity Forceps	2.4	160	2.8	Yellow	Box 5
M00513313	08714729761853	Radial Jaw 4 Large Capacity Forceps w/ Needle	2.4	160	2.8	Yellow	Box 5
M00513323	08714729761877	Radial Jaw 4 Large Capacity Forceps	2.4	240	2.8	Orange	Box 5
M00513333	08714729755494	Radial Jaw 4 Large Capacity Forceps w/ Needle	2.4	240	2.8	Orange	Box 5
M00513301	08714729648758	Radial Jaw 4 Large Capacity Forceps	2.4	160	2.8	Yellow	Box 20
M00513311	08714729648772	Radial Jaw 4 Large Capacity Forceps w/ Needle	2.4	160	2.8	Yellow	Box 20
M00513321	08714729648796	Radial Jaw 4 Large Capacity Forceps	2.4	240	2.8	Orange	Box 20
M00513331	08714729648819	Radial Jaw 4 Large Capacity Forceps w/ Needle	2.4	240	2.8	Orange	Box 20
M00513302	08714729648765	Radial Jaw 4 Large Capacity Forceps	2.4	160	2.8	Yellow	Box 40
M00513312	08714729648789	Radial Jaw 4 Large Capacity Forceps w/ Needle	2.4	160	2.8	Yellow	Box 40
M00513322	08714729648802	Radial Jaw 4 Large Capacity Forceps	2.4	240	2.8	Orange	Box 40
M00513332	08714729648826	Radial Jaw 4 Large Capacity Forceps w/ Needle	2.4	240	2.8	Orange	Box 40

Radial Jaw 4 Standard Capacity Single-Use Biopsy Forceps



Order Number	GTIN	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Color Code	Packaging
M00513383	08714729777915	Radial Jaw 4 Standard Capacity Forceps	2.2	160	2.8	Yellow	Box 5
M00513393	08714729777939	Radial Jaw 4 Standard Capacity Forceps w/ Needle	2.2	160	2.8	Yellow	Box 5
M00513403	08714729777922	Radial Jaw 4 Standard Capacity Forceps	2.2	240	2.8	Orange	Box 5
M00513413	08714729777946	Radial Jaw 4 Standard Capacity Forceps w/ Needle	2.2	240	2.8	Orange	Box 5
M00513381	08714729777953	Radial Jaw 4 Standard Capacity Forceps	2.2	160	2.8	Yellow	Box 20
M00513391	08714729777977	Radial Jaw 4 Standard Capacity Forceps w/ Needle	2.2	160	2.8	Yellow	Box 20
M00513401	08714729777960	Radial Jaw 4 Standard Capacity Forceps	2.2	240	2.8	Orange	Box 20
M00513411	08714729777984	Radial Jaw 4 Standard Capacity Forceps w/ Needle	2.2	240	2.8	Orange	Box 20
M00513382	08714729778011	Radial Jaw 4 Standard Capacity Forceps	2.2	160	2.8	Yellow	Box 40
M00513392	08714729777991	Radial Jaw 4 Standard Capacity Forceps w/ Needle	2.2	160	2.8	Yellow	Box 40
M00513402	08714729778028	Radial Jaw 4 Standard Capacity Forceps	2.2	240	2.8	Orange	Box 40
M00513412	08714729778004	Radial Jaw 4 Standard Capacity Forceps w/ Needle	2.2	240	2.8	Orange	Box 40



Biopsy



Forceps (Continued)

Radial Jaw™ 4 Single-Use Pediatric Biopsy Forceps



Order Number	GTIN	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Packaging
M00513443	08714729790211	Radial Jaw 4 Pediatric Biopsy Forceps	1.8	160	≥ 2.0	Box 5
M00513453	08714729790228	Radial Jaw 4 Pediatric Biopsy Forceps w/ Needle	1.8	160	≥ 2.0	Box 5

Hot Biopsy Forceps



Radial Jaw 4 Single-Use Hot Biopsy Forceps




Order Number	GTIN	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Packaging
M00515031	08714729790136	Radial Jaw 4 Hot Biopsy Forceps	2.2	240	≥ 2.8	Box 5
M00515032	08714729790105	Radial Jaw 4 Hot Biopsy Forceps	2.2	240	≥ 2.8	Box 20
M00515033	08714729790143	Radial Jaw 4 Hot Biopsy Forceps	2.2	240	≥ 2.8	Box 40



Biopsy

Forceps (Continued)



Multibite™ Single-Use Multiple Sample Biopsy Forceps

Order Number	GTIN	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Color Code	Packaging
M00510101	08714729164975	Multibite Multiple Sample Biopsy Forceps	2.4	160	2.8	Yellow	Box 5
M00510121	08714729164982	Multibite Multiple Sample Biopsy Forceps	2.4	240	2.8	Orange	Box 5

Notes and Future Product Updates



Enteral Access



EndoVive™ Enteral Access Devices

Percutaneous Endoscopic Gastrostomy (PEG) Devices – Initial Placement Gastrostomy Devices

➤ Safety PEG Kit with ENFit™	67
➤ Standard PEG Kit with ENFit.....	67
➤ Safety PEG Kit	67
➤ Standard PEG Kit.....	67
➤ Initial Placement PEG Accessories	68

One-Step Button™ Low Profile PEG Kits

➤ One-Step Button Low Profile PEG Kits.....	68
➤ One-Step Button and Low Profile Replacement Button Accessories	68

Gastrostomy Tubes

➤ Replacement Gastrostomy Feeding Tube with ENFit.....	69
➤ Standard Balloon Replacement Kit, Straight Bolster	69
➤ Standard Balloon Replacement Kit, Right Angle Bolster.....	70

Jejunal Feeding Tube Devices

➤ Through-the-PEG (TTP) Jejunal Feeding Tube	70
--	----

Transition Connectors

➤ Distal and Feeding Transition Connectors.....	70
---	----

Accessories & Replacements


➤ Low Profile Replacement Button	71
➤ One-Step Button and Low Profile Replacement Button Accessories	71



EndoVive™ Enteral Access Devices

Percutaneous Endoscopic Gastrostomy (PEG) Devices

Initial Placement Gastrostomy Devices




Safety PEG* Kit with ENFit™



Order Number	GTIN	Description	Packaging
M00509001	08714729880394	20Fr Safety PEG Kit, Pull with ENFit	Box 2
M00509011	08714729880417	20Fr Safety PEG Kit, Push with ENFit	Box 2
M00509021	08714729880431	24Fr Safety PEG Kit, Pull with ENFit	Box 2
M00509031	08714729880455	24Fr Safety PEG Kit, Push with ENFit	Box 2

*Includes local anesthetic.




Standard PEG* Kit with ENFit



Order Number	GTIN	Description	Packaging
M00509041	08714729880479	20Fr Standard PEG Kit, Pull with ENFit	Box 2
M00509051	08714729880493	20Fr Standard PEG Kit, Push with ENFit	Box 2
M00509061	08714729880516	24Fr Standard PEG Kit, Pull with ENFit	Box 2
M00509071	08714729880530	24Fr Standard PEG Kit, Push with ENFit	Box 2

*Includes local anesthetic.




Safety PEG* Kit



Order Number	GTIN	Description	Packaging
M00566461	08714729802624	20Fr Safety PEG Kit, Pull	Box 2
M00566471	08714729802631	20Fr Safety PEG Kit, Push	Box 2
M00566481	08714729802648	24Fr Safety PEG Kit, Pull	Box 2
M00566491	08714729802655	24Fr Safety PEG Kit, Push	Box 2

*Includes local anesthetic.



Standard PEG* Kit



Order Number	GTIN	Description	Packaging
M00568201	08714729285199	20Fr Standard PEG Kit, Pull	Box 2
M00568211	08714729285205	20Fr Standard PEG Kit, Push	Box 2
M00568241	08714729296799	24Fr Standard PEG Kit, Pull	Box 2
M00568251	08714729296805	24Fr Standard PEG Kit, Push	Box 2

*Includes local anesthetic.



Section Table of Contents

Prev Page ◀ 67 ▶ Next Page

Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001


EndoVive™ Enteral Access Devices

Percutaneous Endoscopic Gastrostomy (PEG) Devices

Initial Placement Gastrostomy Devices

Initial Placement PEG Accessories				
Order Number	GTIN	Description		Packaging
M00501261	08714729863564	20Fr Y-Ports with ENFit		Box 10
M00501271	08714729863588	24Fr Y-Ports with ENFit		Box 10
M00580641	08714729169116	20Fr Y-Ports		Box 10
M00580651	08714729169123	24Fr Y-Ports		Box 10
M00580731	08714729415343	20Fr Round Bolsters		Box 10
M00580781	08714729415398	20Fr / 24Fr C Clamps		Box 10

One-Step Button™ Low Profile PEG Kits



One-Step Button Low Profile PEG Kits					
Order Number	GTIN	Description	Length (cm)		Packaging
M00563120	08714729284246	18Fr Low Profile PEG Kit, Pull	1.2		Each
M00563130	08714729284253	18Fr Low Profile PEG Kit, Pull	1.7		Each
M00563140	08714729284260	18Fr Low Profile PEG Kit, Pull	2.4		Each
M00563150	08714729284277	18Fr Low Profile PEG Kit, Pull	3.4		Each
M00563160	08714729284284	18Fr Low Profile PEG Kit, Pull	4.4		Each
M00563180	08714729284307	18Fr Low Profile PEG Kit, Push	1.7		Each
M00563190	08714729284314	18Fr Low Profile PEG Kit, Push	2.4		Each
M00563010	08714729284147	24Fr Low Profile PEG Kit, Pull	1.7		Each
M00563020	08714729284154	24Fr Low Profile PEG Kit, Pull	2.4		Each
M00563030	08714729284161	24Fr Low Profile PEG Kit, Pull	3.4		Each
M00563040	08714729268833	24Fr Low Profile PEG Kit, Pull	4.4		Each
M00563080	08714729284208	24Fr Low Profile PEG Kit, Push	2.4		Each
M00563090	08714729284215	24Fr Low Profile PEG Kit, Push	3.4		Each
M00563100	08714729284222	24Fr Low Profile PEG Kit, Push	4.4		Each

One-Step Button and Low Profile Replacement Button Accessories					
Order Number	GTIN	Description	Length (cm)		Packaging
M00580401	08714729166443	18Fr Decompression Tube	1.2		Box 10
M00580411	08714729166450	18Fr Decompression Tube	1.7		Box 10
M00580421	08714729166467	18Fr Decompression Tube	2.4		Box 10
M00580431	08714729166474	18Fr Decompression Tube	3.4		Box 10
M00580441	08714729166481	18Fr Decompression Tube	4.4		Box 10
M00580461	08714729166504	24Fr Decompression Tube	1.7		Box 10
M00580471	08714729166511	24Fr Decompression Tube	2.4		Box 10
M00580481	08714729339649	24Fr Decompression Tube	3.4		Box 10
M00580491	08714729166535	24Fr Decompression Tube	4.4		Box 10
M00580171	08714729142232	18Fr Bolus Feeding Set Straight Adapter	—		Box 10
M00580181	08714729142249	18Fr Right Angle Feeding Set Continuous Adapter	—		Box 10
M00580231	08714729166405	24Fr Bolus Feeding Set Straight Adapter	—		Box 10
M00580241	08714729166412	24Fr Right Angle Feeding Set Continuous Adapter	—		Box 10



EndoVive™ Enteral Access Devices

Gastrostomy Tubes




Replacement Gastrostomy Feeding Tube with ENFit™, Straight Angle Bolster

Order Number	GTIN	Description	Packaging
M00509881	00813939023862	12Fr Replacement G-Tube, Straight with ENFit	Box 5
M00509891	00813939023879	14Fr Replacement G-Tube, Straight with ENFit	Box 5
M00509901	00813939023886	16Fr Replacement G-Tube, Straight with ENFit	Box 5
M00509911	00813939023893	18Fr Replacement G-Tube, Straight with ENFit	Box 5
M00509921	00813939023909	20Fr Replacement G-Tube, Straight with ENFit	Box 5
M00509931	00813939023916	22Fr Replacement G-Tube, Straight with ENFit	Box 5
M00509941	00813939023923	24Fr Replacement G-Tube, Straight with ENFit	Box 5



Replacement Gastrostomy Feeding Tube with ENFit™, Right Angle Bolster

Order Number	GTIN	Description	Packaging
M00509961	00813939024005	12Fr Replacement G-Tube, Right Angle with ENFit	Box 5
M00509971	00813939024012	14Fr Replacement G-Tube, Right Angle with ENFit	Box 5
M00509981	00813939024029	16Fr Replacement G-Tube, Right Angle with ENFit	Box 5
M00509991	00813939024036	18Fr Replacement G-Tube, Right Angle with ENFit	Box 5
M00510001	00813939024043	20Fr Replacement G-Tube, Right Angle with ENFit	Box 5
M00510011	00813939024050	22Fr Replacement G-Tube, Right Angle with ENFit	Box 5
M00510021	00813939024067	24Fr Replacement G-Tube, Right Angle with ENFit	Box 5



Standard Balloon Replacement Kit, Straight Bolster

Order Number	GTIN	Description	Packaging
M00582040	00813939020243	12Fr Replacement G-Tube, Straight	Each
M00582050	00813939020250	14Fr Replacement G-Tube, Straight	Each
M00582060	00813939020267	16Fr Replacement G-Tube, Straight	Each
M00582070	00813939020274	18Fr Replacement G-Tube, Straight	Each
M00582080	00813939020281	20Fr Replacement G-Tube, Straight	Each
M00582860	00813939020298	22Fr Replacement G-Tube, Straight	Each
M00582090	00813939020304	24Fr Replacement G-Tube, Straight	Each
M00582100	00813939020311	28Fr Replacement G-Tube, Straight	Each



Section Table of Contents

Prev Page

69


Next Page

Main Table of Contents



EndoVive™ Enteral Access Devices

Gastrostomy Tubes (Continued)




Standard Balloon Replacement Kit, Right Angle Bolster



Order Number	GTIN	Description	Packaging
M00582110	00813939020007	12Fr Replacement G-Tube, Right Angle	Each
M00582120	00813939020014	14Fr Replacement G-Tube, Right Angle	Each
M00582130	00813939020021	16Fr Replacement G-Tube, Right Angle	Each
M00582140	00813939020038	18Fr Replacement G-Tube, Right Angle	Each
M00582150	00813939020045	20Fr Replacement G-Tube, Right Angle	Each
M00582870	00813939020052	22Fr Replacement G-Tube, Right Angle	Each
M00582160	00813939020069	24Fr Replacement G-Tube, Right Angle	Each
M00582170	00813939020076	28Fr Replacement G-Tube, Right Angle	Each

Jejunal Feeding Tube Devices



Through-The-PEG (TTP) Jejunal Feeding Tube



Order Number	GTIN	Description	Tip Style	Packaging
M00566320	08714729786313	8.5Fr TTP Jejunal Feeding Tube	Bent Tip	Each
M00566330	08714729789727	8.5Fr TTP Jejunal Feeding Tube	Pigtail Tip	Each
M00566340	08714729789666	12Fr TTP Jejunal Feeding Tube	Bent Tip	Each
M00566350	08714729786252	12Fr TTP Jejunal Feeding Tube	Pigtail Tip	Each

Transition Connectors



Distal and Feeding Transition Connectors




Order Number	GTIN	Description	Packaging
M00509221	N/A	Distal End ENFit™ Transition Connector	Box 10
M00509241	N/A	ENFit Feeding Port Transition Connector	Box 10



EndoVive™ Enteral Access Devices

Accessories & Replacements



Low Profile Replacement Button



Order Number	GTIN	Description	Length (cm)	Packaging
M00562760	08714729284017	18Fr Low Profile Replacement Button	1.2	Each
M00562770	08714729142201	18Fr Low Profile Replacement Button	1.7	Each
M00562780	08714729284024	18Fr Low Profile Replacement Button	2.4	Each
M00562790	08714729284031	18Fr Low Profile Replacement Button	3.4	Each
M00562800	08714729284048	18Fr Low Profile Replacement Button	4.4	Each
M00562820	08714729284062	24Fr Low Profile Replacement Button	1.7	Each
M00562830	08714729142218	24Fr Low Profile Replacement Button	2.4	Each
M00562840	08714729284079	24Fr Low Profile Replacement Button	3.4	Each
M00562850	08714729284086	24Fr Low Profile Replacement Button	4.4	Each

One-Step Button™ and Low Profile Replacement Button Accessories



Order Number	GTIN	Description	Length (cm)	Packaging
M00580401	08714729166443	18Fr Decompression Tube	1.2	Box 10
M00580411	08714729166450	18Fr Decompression Tube	1.7	Box 10
M00580421	08714729166467	18Fr Decompression Tube	2.4	Box 10
M00580431	08714729166474	18Fr Decompression Tube	3.4	Box 10
M00580441	08714729166481	18Fr Decompression Tube	4.4	Box 10
M00580461	08714729166504	24Fr Decompression Tube	1.7	Box 10
M00580471	08714729166511	24Fr Decompression Tube	2.4	Box 10
M00580481	08714729339649	24Fr Decompression Tube	3.4	Box 10
M00580491	08714729166535	24Fr Decompression Tube	4.4	Box 10
M00580171	08714729142232	18Fr Bolus Feeding Set Straight Adapter	—	Box 10
M00580181	08714729142249	18Fr Right Angle Feeding Set Continuous Adapter	—	Box 10
M00580231	08714729166405	24Fr Bolus Feeding Set Straight Adapter	—	Box 10
M00580241	08714729166412	24Fr Right Angle Feeding Set Continuous Adapter	—	Box 10



Section Table of Contents

Prev Page

71

Next Page

Main Table of Contents



Notes and Future Product Updates



Tissue Resection

Tissue Resection

EMR for Upper GI

➤ Captivator™ EMR Devices	73
---------------------------------	----

Stiff Snares

➤ Captivator II Single-Use Snare	73
➤ Captivator Single-Use Snares.....	74

Speciality Snares

➤ Profile™ Single-Use Snares	75
➤ Single-Use Rotatable Snares.....	75

Snares

➤ Captiflex™ Single-Use Snares	76
➤ Sensation™ Single-Use Short Throw Snares.....	76


Active Cords

➤ Active Cords	77
----------------------	----



Tissue Resection

EMR for Upper GI




Captivator™ EMR Device



Order Number	GTIN	Description	Packaging
M00561600	08714729842675	Captivator EMR Device for Standard Gastrosopes	Box 1
M00561610	08714729842699	Captivator EMR Device for Large Gastrosopes	Box 1
M00561601	08714729842682	Captivator EMR Device for Standard Gastrosopes	Box 5
M00561611	08714729842705	Captivator EMR Device for Large Gastrosopes	Box 5

Stiff Snares



Captivator II Single-Use Snare



Order Number	GTIN	Description	Loop Width (mm)	Sheath O.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Packaging
M00561221	08714729855903	10mm Rounded – Stiff	10	2.4	2.8	240	Box 10
M00561222	08714729861300	10mm Rounded – Stiff	10	2.4	2.8	240	Box 20
M00561223	08714729861317	10mm Rounded – Stiff	10	2.4	2.8	240	Box 40
M00561231	08714729855910	15mm Rounded – Stiff	15	2.4	2.8	240	Box 10
M00561232	08714729861331	15mm Rounded – Stiff	15	2.4	2.8	240	Box 20
M00561233	08714729861348	15mm Rounded – Stiff	15	2.4	2.8	240	Box 40
M00561241	08714729855927	20mm Rounded – Stiff	20	2.4	2.8	240	Box 10
M00561242	08714729861362	20mm Rounded – Stiff	20	2.4	2.8	240	Box 20
M00561243	08714729861379	20mm Rounded – Stiff	20	2.4	2.8	240	Box 40
M00561191	08714729855934	25mm Rounded – Stiff	25	2.4	2.8	240	Box 10
M00561192	08714729861270	25mm Rounded – Stiff	25	2.4	2.8	240	Box 20
M00561193	08714729861287	25mm Rounded – Stiff	25	2.4	2.8	240	Box 40
M00561291	08714729645801	Extra Large Rounded – Stiff	33	2.4	2.8	240	Box 10



Section Table of Contents

Prev Page

73


Next Page

Main Table of Contents



Tissue Resection

Stiff Snares (Continued)



Captivator™ Single-Use Snares



Order Number	GTIN	Description	Loop Width (mm)	Sheath O.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Packaging
M00562301	08714729019312	Small Oval – Stiff	13	2.4	2.8	240	Box 10
M00562321	08714729019312	Medium Oval – Stiff	27	2.4	2.8	240	Box 10
M00561311	08714729019251	Large Oval – Medium Stiff	30	2.4	2.8	240	Box 10
M00562391	08714729019404	Large Oval (thin wire) – Flexible	30	2.4	2.8	240	Box 10
M00562451	08714729071068	Small Hexagonal – Stiff	13	2.4	2.8	240	Box 10
M00562341	08714729019350	Medium Hexagonal – Stiff	27	2.4	2.8	240	Box 10
M00562371	08714729019381	Medium Crescent – Stiff	27	2.4	2.8	240	Box 10



Section Table of Contents

Prev Page ◀ 74 ▶ Next Page


Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Tissue Resection


Speciality Snares



Profile™ Single-Use Snares



Order Number	GTIN	Description	Loop Width (mm)	Sheath O.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Packaging
M00562531	08714729156628	Extra Small Oval – Flexible	11	1.9	2.0	240	Box 10
M00562551	08714729156659	Small Oval – Flexible	13	1.9	2.0	240	Box 10
M00562571	08714729156635	Medium Oval – Flexible	27	1.9	2.0	240	Box 10



Single-Use Rotatable Snares*




Order Number	GTIN	Description	Loop Width (mm)	Sheath O.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Packaging
M00561821	08714729283980	Small Oval – Medium Stiff	13	2.4	2.8	240	Box 10
M00561831	08714729283997	Medium Oval – Medium Stiff	20	2.4	2.8	240	Box 10

*With Star Catheter™ Technology.



Tissue Resection


Snares



Captiflex™ Single-Use Snares



Order Number	GTIN	Description	Loop Width (mm)	Sheath O.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Packaging
M00562471	08714729019459	Extra Small Oval – Flexible	11	2.4	2.8	240	Box 10
M00562421	08714729019435	Small Oval – Flexible	13	2.4	2.8	240	Box 10
M00562422	08714729501657	Small Oval – Flexible	13	2.4	2.8	240	Box 40
M00562401	08714729019411	Medium Oval – Flexible	27	2.4	2.8	240	Box 10
M00562402	08714729501640	Medium Oval – Flexible	27	2.4	2.8	240	Box 40



Sensation™ Single-Use Short Throw Snares



Order Number	GTIN	Description	Loop Width (mm)	Sheath O.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Packaging
M00562691	08714729158158	Small Oval – Flexible	13	2.4	2.8	240	Box 5
M00562692	08714729268826	Small Oval – Flexible	13	2.4	2.8	240	Box 20
M00562693	08714729501671	Small Oval – Flexible	13	2.4	2.8	240	Box 40
M00560321	08714729283928	Small Oval – Medium Stiff	13	2.4	2.8	240	Box 5
M00562671	08714729158141	Medium Oval – Flexible	27	2.4	2.8	240	Box 5
M00562672	08714729268819	Medium Oval – Flexible	27	2.4	2.8	240	Box 20
M00562673	08714729501664	Medium Oval – Flexible	27	2.4	2.8	240	Box 40
M00560311	08714729283904	Medium Oval – Medium Stiff	27	2.4	2.8	240	Box 5
M00562651	08714729158110	Large Oval – Medium Stiff	30	2.4	2.8	240	Box 5
M00562652	08714729268802	Large Oval – Medium Stiff	30	2.4	2.8	240	Box 20
M00560331	08714729283935	Medium Crescent – Medium Stiff	27	2.4	2.8	240	Box 5



Section Table of Contents

Prev Page **76** Next Page

Main Table of Contents



Tissue Resection

Active Cords

Active Cords			
Order Number	GTIN	Description	Packaging
M00561250	08714729019060	Active Cord, Fits Olympus™ Brand Generator	Each
M00561270	08714729053361	Active Cord, Fits Endostat™, Endostat II and Endostat III Generators, Bovie™, Valleylab™, Aspen Labs™, and ERBE™ Brand Generators	Each



Notes and Future Product Updates



Pulmonary Endoscopy

Alair™ Bronchial Thermoplasty System

- Alair Bronchial Thermoplasty System 79

Airway Stents

- Dynamic™ (Y) Stent..... 79
- Polyflex™ Single-Use Self-Expanding Silicone Airway Stent.... 82
- Ultraflex™ Single-Use Partially Covered Tracheobronchial Stent System – Distal Release..... 81
- Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Distal Release 80
- Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Proximal Release..... 80

Biopsy Forceps

- Radial Jaw™ 4 Single-Use Pulmonary Biopsy Forceps..... 83

Endobronchial Ultrasound Transbronchial Aspiration Needle (EBUS)

- Expect™ Pulmonary Endobronchial Ultrasound Transbronchial Aspiration Needle / Adaptor 83

Transbronchial Aspiration Needles

- eXcelon™ Transbronchial Aspiration Needle 83

Cytology Brushes

- Cellebri™ Single-Use Endoscopic Cytology Brush..... 84

Airway Retrieval Baskets

- Zero Tip™ Airway Retrieval Basket 84

Controlled Radial Expansion Pulmonary Balloon Dilatation Catheter

- CRE™ Pulmonary Balloon Dilatation Catheter..... 85

Pulmonary Guidewires


- Pulmonary Jagwire™ Pulmonary Guidewire 85

Inflation System

- Alliance™ II Integrated Inflation / Lithotripsy Device..... 85
- Alliance II Single-Use Syringe / Gauge Assembly 85

Pulmonary

Alair™ Bronchial Thermoplasty System



Alair Bronchial Thermoplasty System*


Order Number	GTIN	Model Number	Description	Packaging
M005 ATS20000	08714729796695	Alair ATS 200	Alair Bronchial Thermoplasty Controller	Each
M005 ATS25010	08714729802792	Alair ATS 2-5	⊗ Alair Bronchial Thermoplasty Catheter**	Each
M005 ATS20110	08714729802778	Alair ATS 201	Alair RF Controller Accessory Kit	Each

*Service Agreement available with 1, 2, 3, and 4 year terms.

****Note:** Initial stocking order requires a minimum order of 6 catheters (covering the complete treatment of 2 patients).

UPN	Description
SRVM005 ATS1001	ENDO POS 1 Year Service Alair
SRVM005 ATS1002	ENDO POS 2 Year Service Alair
SRVM005 ATS1003	ENDO POS 3 Year Service Alair
SRVM005 ATS1004	ENDO POS 4 Year Service Alair

Airway Stents



Dynamic™ (Y) Stent




Order Number	GTIN	Tracheal Width (ID) (mm)	Bronchial Width (ID) (mm)	Tracheal Length (mm)	Bronchial Lengths (R / L) (mm)	Packaging
M005 70670	08714729764809	11	8	110	25 / 40	Each
M005 70680	08714729764816	13	10	110	25 / 40	Each
M005 70690	08714729764823	15	12	110	25 / 40	Each



Pulmonary

Airway Stents (Continued)



 Magnetic Resonance Conditional¹

Ultraflex™ Single-Use Uncovered Tracheobronchial Stent System – Distal Release*



Order Number	GTIN	Description	Expanded Stent O.D. (mm)	Expanded Stent Length (mm)	Compressed Stent O.D. (Fr)**	Tip Max O.D. (mm)	Packaging
M00576300	08714729837688	Uncovered Distal Release	10	30	15	4.1	Each
M00576310	08714729837695	Uncovered Distal Release	12	30	17	4.1	Each
M00576320	08714729837701	Uncovered Distal Release	14	30	18	4.1	Each
M00576330	08714729837718	Uncovered Distal Release	16	40	20	5.3	Each
M00576340	08714729837725	Uncovered Distal Release	16	60	20	5.3	Each
M00576350	08714729837732	Uncovered Distal Release	18	40	20	5.3	Each
M00576360	08714729837749	Uncovered Distal Release	18	60	20	5.3	Each
M00576370	08714729837756	Uncovered Distal Release	20	40	20	5.3	Each
M00576380	08714729837763	Uncovered Distal Release	20	60	20	5.3	Each

*Recommended Guidewire .035" Pulmonary Jagwire™ Pulmonary Guidewire, Order #1517, [\(Click here to view page 85\)](#).
Delivery System working length is 95cm.

**Data on File at Boston Scientific.



 Magnetic Resonance Conditional¹

Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Proximal Release*



Order Number	GTIN	Description	Expanded Stent O.D. (mm)	Expanded Stent Length (mm)	Compressed Stent O.D. (Fr)**	Tip Max O.D. (mm)	Packaging
M00576410	08714729837770	Uncovered Proximal Release	8	20	15	4.1	Each
M00576420	08714729837787	Uncovered Proximal Release	8	40	15	4.1	Each
M00576430	08714729837794	Uncovered Proximal Release	10	20	15	4.1	Each
M00576440	08714729837800	Uncovered Proximal Release	10	40	15	4.1	Each
M00576450	08714729837817	Uncovered Proximal Release	12	20	17	4.1	Each
M00576460	08714729837824	Uncovered Proximal Release	12	40	17	4.1	Each
M00576470	08714729837831	Uncovered Proximal Release	14	20	18	4.1	Each
M00576480	08714729837848	Uncovered Proximal Release	14	40	18	4.1	Each
M00576490	08714729837855	Uncovered Proximal Release	14	60	18	4.1	Each

*Recommended Guidewire .035" Pulmonary Jagwire Pulmonary Guidewire, Order #1517, [\(Click here to view page 85\)](#).
Delivery System working length is 95cm.

**Data on File at Boston Scientific.



Section Table of Contents


Prev Page  80  Next Page

Main Table of Contents 

For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Pulmonary

Airway Stents (Continued)



 Magnetic Resonance
Conditional¹

Ultraflex™ Single-Use Partially Covered Tracheobronchial Stent System – Distal Release*



Order Number	GTIN	Description	Expanded Stent O.D. (mm)	Expanded Stent Length (mm)	Covered Compressed Stent O.D. (Fr)**	Cover Length (mm)	Tip Max O.D. (mm)	Packaging
M00576520	08714729842323	Partially Covered Distal Release	8	40	16	25	4.1	Each
M00576530	08714729842330	Partially Covered Distal Release	10	30	18	15	4.1	Each
M00576540	08714729842347	Partially Covered Distal Release	10	40	18	25	4.1	Each
M00576550	08714729842354	Partially Covered Distal Release	12	30	20	15	4.1	Each
M00576560	08714729842361	Partially Covered Distal Release	12	40	20	25	4.1	Each
M00576570	08714729842378	Partially Covered Distal Release	14	30	24	15	4.1	Each
M00576580	08714729842385	Partially Covered Distal Release	14	40	24	25	4.1	Each
M00576590	08714729842392	Partially Covered Distal Release	14	60	24	45	4.1	Each
M00576600	08714729842408	Partially Covered Distal Release	14	80	24	65	4.1	Each
M00576610	08714729842415	Partially Covered Distal Release	16	40	25	25	5.3	Each
M00576620	08714729842422	Partially Covered Distal Release	16	60	25	45	5.3	Each
M00576630	08714729842439	Partially Covered Distal Release	16	80	25	65	5.3	Each
M00576640	08714729842446	Partially Covered Distal Release	18	40	27	25	5.3	Each
M00576650	08714729842453	Partially Covered Distal Release	18	60	27	45	5.3	Each
M00576660	08714729842460	Partially Covered Distal Release	18	80	27	65	5.3	Each
M00576670	08714729842477	Partially Covered Distal Release	20	40	27	25	5.3	Each
M00576680	08714729842484	Partially Covered Distal Release	20	60	27	45	5.3	Each
M00576690	08714729842491	Partially Covered Distal Release	20	80	27	65	5.3	Each

*Recommended Guidewire .035" Pulmonary Jagwire™ Pulmonary Guidewire, Order #1517, [\(Click here to view page 85\)](#).
Delivery System working length is 95cm.

**Data on File at Boston Scientific.



Section Table of Contents

Prev Page **81** Next Page

Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Pulmonary

Airway Stents (Continued)



Polyflex™ Single-Use Self-Expanding Silicone Airway Stent




Order Number	GTIN	Stent ID (mm)	Length (mm)	Delivery System Diameter (mm)	Packaging
M00570000	08714729457138	8	20	7	Each
M00570010	08714729457145	8	30	7	Each
M00570020	08714729457152	10	20	8	Each
M00570030	08714729457169	10	30	8	Each
M00570040	08714729457176	10	40	8	Each
M00570050	08714729457183	10	50	8	Each
M00570060	08714729457190	12	20	9	Each
M00570070	08714729457206	12	30	9	Each
M00570080	08714729457213	12	40	9	Each
M00570090	08714729457220	12	50	9	Each
M00570100	08714729457237	14	20	9	Each
M00570110	08714729457244	14	30	9	Each
M00570120	08714729457251	14	40	9	Each
M00570130	08714729457268	14	50	9	Each
M00570140	08714729457275	14	60	9	Each
M00570150	08714729457282	16	30	10	Each
M00570160	08714729457299	16	40	10	Each
M00570170	08714729457305	16	50	10	Each
M00570180	08714729457312	16	60	10	Each
M00570190	08714729457329	16	70	10	Each
M00570200	08714729457336	18	30	11	Each
M00570210	08714729457343	18	40	11	Each
M00570220	08714729457350	18	50	11	Each
M00570230	08714729457367	18	60	11	Each
M00570240	08714729457374	18	70	11	Each
M00570250	08714729457381	18	80	11	Each
M00570260	08714729457398	20	40	12	Each
M00570270	08714729457404	20	50	12	Each
M00570280	08714729457411	20	60	12	Each
M00570290	08714729457428	20	70	12	Each
M00570300	08714729457435	20	80	12	Each
M00570310	08714729457442	22	50	13	Each
M00570320	08714729457459	22	60	13	Each
M00570330	08714729457466	22	80	13	Each



Pulmonary

Biopsy Forceps




Radial Jaw™ 4 Single-Use Pulmonary Biopsy Forceps



Order Number	GTIN	Description	Differentiator	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Packaging
M00515181	08714729792864	Single-Use Biopsy Forceps	Pulmonary Standard Capacity	1.8	100	2.0	Box 5
M00515182	08714729799337	Single-Use Biopsy Forceps	Pulmonary Standard Capacity	1.8	100	2.0	Box 20
M00515191	08714729792871	Single-Use Biopsy Forceps	Pulmonary Standard Capacity w/ Needle	1.8	100	2.0	Box 5
M00515192	08714729818083	Single-Use Biopsy Forceps	Pulmonary Standard Capacity w/ Needle	1.8	100	2.0	Box 20
M00515201	08714729793281	Single-Use Biopsy Forceps	Pulmonary Large Capacity	2.2	100	2.8	Box 5
M00515202	08714729818076	Single-Use Biopsy Forceps	Pulmonary Large Capacity	2.2	100	2.8	Box 20

Expect™ Pulmonary Endobronchial Ultrasound Transbronchial Aspiration Needle




Expect Pulmonary Endobronchial Ultrasound Transbronchial Aspiration Needle



Order Number	GTIN	Needle Size	Minimum Working Channel (mm)	Sheath O.D. (mm)	Packaging
M00558220	08714729861409	Expect Pulmonary 22Ga – Olympus scope compatible	2.0	1.6	Please contact your Boston Scientific Sales Representative for ordering information
M00558250	08714729861416	Expect Pulmonary 25Ga – Olympus scope compatible	2.0	1.4	Please contact your Boston Scientific Sales Representative for ordering information
M00558731	08714729893189	Expect Pulmonary Needle Adaptor – Olympus scope compatible	N/A	N/A	Please contact your Boston Scientific Sales Representative for ordering information

Transbronchial Aspiration Needles



eXcelon™ Transbronchial Aspiration Needle*



Order Number	GTIN	Description	Needle Gauge	Needle Length (mm)	Catheter Length (cm)	Sheath O.D. (mm)	Packaging
M00564101	08714729456452	eXcelon Transbronchial Aspiration Needle	19	15	130	1.8	Box 5
M00564111	08714729456469	eXcelon Transbronchial Aspiration Needle	20	15	130	1.8	Box 5
M00564121	08714729456476	eXcelon Transbronchial Aspiration Needle	21	15	130	1.8	Box 5

*Packaging also includes a 20ml syringe.



Section Table of Contents

Prev Page ◀ 83 ▶ Next Page

Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Pulmonary


Cytology Brushes

Cellebri™ Single-Use Endoscopic Cytology Brush



Order Number	GTIN	Required Working Channel (mm)	Bristle O.D. (mm)	Sheath Length (cm)	Packaging
M00516001	08714729138983	2.0	1.0	140	Box 10
M00516011	08714729138990	2.0	1.5	140	Box 10
M00516071	08714729139089	2.0	1.9	100	Box 10
M00516151	08714729014003	2.0	1.9	150	Box 10

Airway Retrieval Baskets



Zero Tip™ Airway Retrieval Basket



Order Number	GTIN	Description	Sheath O.D. (mm)	Working Length (cm)	Basket Opening (mm)	Packaging
M00513200	08714729414995	Zero Tip Airway Retrieval Basket	0.8	120	12	Each
M00513210	08714729415008	Zero Tip Airway Retrieval Basket	1.0	120	16	Each



Section Table of Contents

Prev Page

84

Next Page

Main Table of Contents



Pulmonary

Controlled Radial Expansion Pulmonary Balloon Dilatation Catheter



CRE™ Pulmonary Balloon Dilatation Catheter*



Order Number	GTIN	Balloon Inflated O.D.		Balloon Length (cm)	Inflation Pressures		Catheter Length (cm)	Packaging
		(mm)	(Fr)		(atm)	(kPa)		
M00550300	08714729456186	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	75	Each
M00550310	08714729456193	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	75	Each
M00550320	08714729456209	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	75	Each
M00550330	08714729456216	8-9-10	24-27-30	3.0	3-5.5-9	304-557-912	75	Each
M00550340	08714729456223	10-11-12	30-33-36	3.0	3-5-8	304-507-811	75	Each
M00550350	08714729456230	12-13.5-15	36-40.5-45	3.0	3-4.5-8	304-456-811	75	Each

*Recommended Inflation Device: Alliance™ II Integrated Inflation / Lithotripsy Device Order #5062, Alliance II Single-Use Gauge / Syringe Assembly Order #5060 (see below).


Pulmonary Guidewires

Pulmonary Jagwire™ Pulmonary Guidewire



Order Number	GTIN	Description	O.D. (in)	Length (cm)	Packaging
M00515171	08714729455813	Pulmonary Jagwire Pulmonary Guidewire	.035	180	Box 2

Inflation System



Alliance™ II Integrated Inflation / Lithotripsy Device and Alliance II Single-Use Syringe / Gauge Assembly

Order Number	GTIN	Description	Packaging
M00550620	08714729283881	Alliance II Integrated Inflation / Lithotripsy Device	Each
M00550601	08714729129332	 Alliance II Single-Use Syringe / Gauge Assembly	Box 5



Section Table of Contents

Prev Page **85** Next Page


Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Notes and Future Product Updates





Services

- **GuidePoint Reimbursement Solutions**..87
- **Performance Optimization Services**88

Services

GuidePoint Reimbursement Resources and Solutions

The Endoscopy Health Economics and Market Access (HEMA) Team offers a suite of reimbursement services and support for healthcare professionals seeking a better understanding of reimbursement for Boston Scientific's endoscopy and pulmonary products and procedures. We can help you develop a base level of knowledge in coding, coverage and payment for these procedures and offer educational presentations and webinars on changes in healthcare policies and Medicare reimbursement and changes to Medicare reimbursement impacting providers and facilities.

About Our Team

Boston Scientific Endoscopy's HEMA Team is comprised of highly skilled professionals specializing in reimbursement, health economics and health policy issues. The team has a broad background with experience working for public and private payers, medical associations, healthcare consulting agencies and academic research organizations. We use clinical and economic outcomes data and collaborate with outside experts to promote appropriate coding, coverage and reimbursement for Boston Scientific endoscopy and pulmonary products and procedures.

Contact Information

For more information about Boston Scientific Endoscopy's GuidePoint Reimbursement Resources and Solutions, please contact our **Health Economics and Reimbursement department at 800-876-9960 x54510 or ENDOREIMBURSEMENT@bsci.com**.

www.bostonscientific.com/reimbursement

GUIDEPOINT
Reimbursement Resources

Health economic and reimbursement information provided by Boston Scientific Corporation is gathered from third-party sources and is subject to change without notice as a result of complex and frequently changing laws, regulations, rules and policies. This information is presented for illustrative purposes only and does not constitute reimbursement or legal advice. Boston Scientific encourages providers to submit accurate and appropriate claims for services. It is always the provider's responsibility to determine medical necessity, the proper site for delivery of any services and to submit appropriate codes, charges, and modifiers for services that are rendered. Boston Scientific recommends that you consult with your payers, reimbursement specialists and/or legal counsel regarding coding, coverage and reimbursement matters. Boston Scientific does not promote the use of its products outside their FDA-approved label.



Services

Boston Scientific's Endoscopy business develops less invasive medical devices for treating gastrointestinal diseases. We are committed to providing clinically proven technologies, services and customer support to help healthcare providers improve patient outcomes while reducing the total cost of patient care.

Performance Optimization Services

Boston Scientific's Supply Chain Optimization services are designed to **improve operational efficiencies** and cost savings. Working with a hospital's internal supply chain team is a key component of developing solutions for inventory management, end-to-end order management and more.

ADVANTICS™ Innovative Healthcare Solutions is a new way to collaborate to meet your specific challenges in the changing healthcare landscape. ADVANTICS recognizes that every provider is unique and that only by working together can we find truly effective solutions.

ADVANTICS™
Innovative Healthcare Solutions

To learn more visit, www.bostonscientific.com/gastroservices or email endoscopysolutions@bsci.com.



Future Product Updates




Notes and Future Product Updates



Notes and Future Product Updates





Product Code Index	
➤ Begins on Page.....	91

Product Code Index

Product Code	Description	Page #	Product Code	Description	Page #
H965180011	WALLSTENT™ Super Stiff Guidewire.....	46	M00510011	EndoVive Replacement Gastrostomy Feeding Tube, with ENFit, Right Angle	69
H965430100	WALLSTENT Endoscopic Biliary Endoprosthesis Stents	19	M00510021	EndoVive Replacement Gastrostomy Feeding Tube, with ENFit, Right Angle	69
H965430200	WALLSTENT Endoscopic Biliary Endoprosthesis Stents	19	M00510101	Multibite™ Single-Use Multiple Sample Biopsy Forceps.....	65
H965430300	WALLSTENT Endoscopic Biliary Endoprosthesis Stents	19	M00510121	Multibite Single-Use Multiple Sample Biopsy Forceps.....	65
H965430400	WALLSTENT Endoscopic Biliary Endoprosthesis Stents	19	M00510860	Trapezoid™ RX Wireguided Retrieval Basket	33
H965430500	WALLSTENT Endoscopic Biliary Endoprosthesis Stents	19	M00510870	Trapezoid RX Wireguided Retrieval Basket.....	33
H965430600	WALLSTENT Endoscopic Biliary Endoprosthesis Stents	19	M00510880	Trapezoid RX Wireguided Retrieval Basket.....	33
H965430700	WALLSTENT Endoscopic Biliary Endoprosthesis Stents	19	M00510890	Trapezoid RX Wireguided Retrieval Basket.....	33
H965430800	WALLSTENT Endoscopic Biliary Endoprosthesis Stents	19	M00511731	Optiflo™ Single-Use Hemostasis Catheters	53
H965430900	WALLSTENT Endoscopic Biliary Endoprosthesis Stents	19	M00511751	Optiflo Single-Use Hemostasis Catheters	53
H965431000	WALLSTENT Endoscopic Biliary Endoprosthesis Stents	19	M00513200	Zero Tip™ Airway Retrieval Basket	84
H965431100	WALLSTENT Endoscopic Biliary Endoprosthesis Stents	19	M00513210	Zero Tip Airway Retrieval Basket	84
H965431200	WALLSTENT Endoscopic Biliary Endoprosthesis Stents	19	M00513301	Radial Jaw™ 4 Single-Use Large Capacity Biopsy Forceps	63
M00501261	EndoVive™ Initial Placement PEG Accessories – Y-Ports with ENFit™	68	M00513302	Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps	63
M00501271	EndoVive Initial Placement PEG Accessories – Y-Ports with ENFit.....	68	M00513303	Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps	63
M00509001	EndoVive Safety PEG Kit, Pull with ENFit.....	67	M00513311	Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps w/ Needle.....	63
M00509011	EndoVive Safety PEG Kit, Push with ENFit	67	M00513312	Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps w/ Needle.....	63
M00509021	EndoVive Safety PEG Kit, Pull with ENFit.....	67	M00513313	Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps w/ Needle.....	63
M00509031	EndoVive Safety PEG Kit, Push with ENFit	67	M00513321	Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps	63
M00509041	EndoVive Standard PEG Kit, Pull with ENFit.....	67	M00513322	Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps	63
M00509051	EndoVive Standard PEG Kit, Push with ENFit.....	67	M00513323	Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps	63
M00509061	EndoVive Standard PEG Kit, Pull with ENFit.....	67	M00513331	Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps w/ Needle.....	63
M00509071	EndoVive Standard PEG Kit, Push with ENFit.....	67	M00513332	Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps w/ Needle.....	63
M00509221	Distal End ENFit Transition Connector.....	70	M00513333	Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps w/ Needle.....	63
M00509241	ENFit Feeding Port Transition Connector.....	70	M00513343	Radial Jaw 4 Single-Use Jumbo Biopsy Forceps	63
M00509881	EndoVive Replacement Gastrostomy Feeding Tube, with ENFit, Straight.....	69	M00513353	Radial Jaw 4 Single-Use Jumbo Biopsy Forceps w/ Needle..	63
M00509891	EndoVive Replacement Gastrostomy Feeding Tube, with ENFit, Straight.....	69	M00513361	Radial Jaw 4 Single-Use Jumbo Biopsy Forceps	63
M00509901	EndoVive Replacement Gastrostomy Feeding Tube, with ENFit, Straight.....	69	M00513362	Radial Jaw 4 Single-Use Jumbo Biopsy Forceps	63
M00509911	EndoVive Replacement Gastrostomy Feeding Tube, with ENFit, Straight.....	69	M00513363	Radial Jaw 4 Single-Use Jumbo Biopsy Forceps	63
M00509921	EndoVive Replacement Gastrostomy Feeding Tube, with ENFit, Straight.....	69	M00513371	Radial Jaw 4 Single-Use Jumbo Biopsy Forceps w/ Needle..	63
M00509931	EndoVive Replacement Gastrostomy Feeding Tube, with ENFit, Straight.....	69	M00513372	Radial Jaw 4 Single-Use Jumbo Biopsy Forceps w/ Needle..	63
M00509941	EndoVive Replacement Gastrostomy Feeding Tube, with ENFit, Straight.....	69	M00513373	Radial Jaw 4 Single-Use Jumbo Biopsy Forceps w/ Needle..	63
M00509961	EndoVive Replacement Gastrostomy Feeding Tube, with ENFit, Right Angle	69	M00513381	Radial Jaw 4 Standard Capacity Single-Use Biopsy Forceps..	63
M00509971	EndoVive Replacement Gastrostomy Feeding Tube, with ENFit, Right Angle	69	M00513382	Radial Jaw 4 Standard Capacity Single-Use Biopsy Forceps..	63
M00509981	EndoVive Replacement Gastrostomy Feeding Tube, with ENFit, Right Angle	69	M00513383	Radial Jaw 4 Standard Capacity Single-Use Biopsy Forceps..	63
M00509991	EndoVive Replacement Gastrostomy Feeding Tube, with ENFit, Right Angle	69	M00513391	Radial Jaw 4 Standard Capacity Single-Use Biopsy Forceps w/ Needle.....	63
M00510001	EndoVive Replacement Gastrostomy Feeding Tube, with ENFit, Right Angle	69	M00513392	Radial Jaw 4 Standard Capacity Single-Use Biopsy Forceps w/ Needle.....	63
			M00513393	Radial Jaw 4 Standard Capacity Single-Use Biopsy Forceps w/ Needle.....	63
			M00513401	Radial Jaw 4 Standard Capacity Single-Use Biopsy Forceps..	63



Product Code Index

Product Code	Description	Page #	Product Code	Description	Page #
M00513402	Radial Jaw™ 4 Standard Capacity Single-Use Biopsy Forceps	63	M00515031	Radial Jaw 4 Single-Use Hot Biopsy Forceps	64
M00513403	Radial Jaw 4 Standard Capacity Single-Use Biopsy Forceps..	63	M00515032	Radial Jaw 4 Single-Use Hot Biopsy Forceps	64
M00513411	Radial Jaw 4 Standard Capacity Single-Use Biopsy Forceps w/ Needle.....	63	M00515033	Radial Jaw 4 Single-Use Hot Biopsy Forceps	64
M00513412	Radial Jaw 4 Standard Capacity Single-Use Biopsy Forceps w/ Needle.....	63	M00515171	Pulmonary Jagwire™ Pulmonary Guidewire.....	85
M00513413	Radial Jaw 4 Standard Capacity Single-Use Biopsy Forceps w/ Needle.....	63	M00515181	Radial Jaw 4 Single-Use Pulmonary Forceps – Standard Capacity	83
M00513443	Radial Jaw 4 Single-Use Pediatric Biopsy Forceps	64	M00515182	Radial Jaw 4 Single-Use Pulmonary Forceps – Standard Capacity	83
M00513453	Radial Jaw 4 Single-Use Pediatric Biopsy Forceps w/ Needle....	64	M00515191	Radial Jaw 4 Single-Use Pulmonary Forceps – Standard Capacity with Needle.....	83
M00513700	Ultraflex™ Single-Use Uncovered Esophageal NG Stent System – Distal Release	49	M00515192	Radial Jaw 4 Single-Use Pulmonary Forceps – Standard Capacity with Needle.....	83
M00513710	Ultraflex Single-Use Uncovered Esophageal NG Stent System – Distal Release	49	M00515201	Radial Jaw 4 Single-Use Pulmonary Forceps – Large Capacity	83
M00513720	Ultraflex Single-Use Uncovered Esophageal NG Stent System – Distal Release	49	M00515202	Radial Jaw 4 Single-Use Pulmonary Forceps – Large Capacity	83
M00513730	Ultraflex Single-Use Covered Esophageal NG Stent System – Distal Release	48	M00516001	Celebrity™ Single-Use Endoscopic Cytology Brush	84
M00513740	Ultraflex Single-Use Covered Esophageal NG Stent System – Distal Release	48	M00516011	Celebrity Single-Use Endoscopic Cytology Brush	84
M00513750	Ultraflex Single-Use Covered Esophageal NG Stent System – Distal Release	48	M00516071	Celebrity Single-Use Endoscopic Cytology Brush	84
M00513800	Ultraflex Single-Use Uncovered Esophageal NG Stent System – Proximal Release	49	M00516151	Celebrity Single-Use Endoscopic Cytology Brush	84
M00513810	Ultraflex Single-Use Uncovered Esophageal NG Stent System – Proximal Release	49	M00516700	WallFlex™ Fully Covered Esophageal Stent.....	47
M00513830	Ultraflex Single-Use Uncovered Esophageal NG Stent System – Proximal Release	49	M00516710	WallFlex Fully Covered Esophageal Stent.....	47
M00513840	Ultraflex Single-Use Covered Esophageal NG Stent System – Proximal Release	48	M00516720	WallFlex Fully Covered Esophageal Stent.....	47
M00513850	Ultraflex Single-Use Covered Esophageal NG Stent System – Proximal Release	48	M00516730	WallFlex Fully Covered Esophageal Stent.....	47
M00513860	Ultraflex Single-Use Covered Esophageal NG Stent System – Proximal Release	48	M00516740	WallFlex Fully Covered Esophageal Stent.....	47
M00514200	Ultraflex Single-Use Covered Large Esophageal NG Stent System – Distal Release	49	M00516750	WallFlex Fully Covered Esophageal Stent.....	47
M00514210	Ultraflex Single-Use Covered Large Esophageal NG Stent System – Distal Release	49	M00516900	WallFlex Partially Covered Esophageal Stent	47
M00514240	Ultraflex Single-Use Covered Large Esophageal NG Stent System – Proximal Release	49	M00516910	WallFlex Partially Covered Esophageal Stent	47
M00514250	Ultraflex Single-Use Covered Large Esophageal NG Stent System – Proximal Release	49	M00516920	WallFlex Partially Covered Esophageal Stent	47
M00514270	Polyflex™ Single-Use Esophageal Stent System.....	48	M00516930	WallFlex Partially Covered Esophageal Stent	47
M00514280	Polyflex Single-Use Esophageal Stent System.....	48	M00516940	WallFlex Partially Covered Esophageal Stent	47
M00514290	Polyflex Single-Use Esophageal Stent System.....	48	M00516950	WallFlex Partially Covered Esophageal Stent	47
M00514300	Polyflex Single-Use Esophageal Stent System.....	48	M00518101	Interject™ Clear Single-Use Injection Therapy Needle Catheters.....	53
M00514310	Polyflex Single-Use Esophageal Stent System.....	48	M00518111	Interject Clear Single-Use Injection Therapy Needle Catheters.....	53
M00514320	Polyflex Single-Use Esophageal Stent System.....	48	M00518151	Interject Contrast Single-Use Injection Therapy Needle Catheters.....	53
M00514330	Polyflex Single-Use Esophageal Stent System.....	48	M00518161	Interject Contrast Single-Use Injection Therapy Needle Catheters.....	53
M00514340	Polyflex Single-Use Esophageal Stent System.....	48	M00518251	Interject Contrast Single-Use Injection Therapy Needle Catheters.....	53
M00514350	Polyflex Single-Use Esophageal Stent System.....	48	M00518261	Interject Contrast Single-Use Injection Therapy Needle Catheters.....	53
			M00518301	Interject Clear Single-Use Injection Therapy Needle Catheters.....	53
			M00518311	Interject Clear Single-Use Injection Therapy Needle Catheters.....	53



Product Code Index

Product Code	Description	Page #	Product Code	Description	Page #
M00518351	Interject™ Contrast Single-Use Injection Therapy Needle Catheters.....	53	M00532870	Advanix Duodenal Bend - Single Stent	21
M00518361	Interject Contrast Single-Use Injection Therapy Needle Catheters.....	53	M00532880	Advanix Duodenal Bend - Single Stent	21
M00521230	Resolution 360™ Clip.....	51	M00532890	Advanix Duodenal Bend - Single Stent	21
M00521231	Resolution 360 Clip.....	51	M00532900	Advanix Duodenal Bend - Single Stent	21
M00521232	Resolution 360 Clip.....	51	M00532910	Advanix Duodenal Bend - Single Stent	21
M00521240	Resolution 360 Clip.....	51	M00532920	Advanix Duodenal Bend - Single Stent	21
M00521241	Resolution 360 Clip.....	51	M00532930	Advanix Duodenal Bend - Single Stent	21
M00521242	Resolution 360 Clip.....	51	M00532940	Advanix Duodenal Bend - Single Stent	21
M00522600	Resolution™ Clip.....	51	M00532950	Advanix Duodenal Bend - Single Stent	21
M00522601	Resolution Clip.....	51	M00532960	Advanix Duodenal Bend - Single Stent	21
M00522602	Resolution Clip.....	51	M00532970	Advanix Duodenal Bend - Single Stent	21
M00522610	Resolution Clip.....	51	M00532980	Advanix Duodenal Bend - Single Stent	21
M00522611	Resolution Clip.....	51	M00532990	Advanix Duodenal Bend - Single Stent	21
M00522612	Resolution Clip.....	51	M00533000	Advanix Duodenal Bend - Single Stent	21
M00530800	Ultratome™ Double-Lumen Sphincterotome – Short Nose.....	12	M00533010	Advanix Duodenal Bend - Single Stent	21
M00530810	Ultratome Double-Lumen Sphincterotome – Long Nose.....	12	M00533020	Advanix Duodenal Bend - Single Stent	21
M00530820	Ultratome Double-Lumen Sphincterotome – Short Nose.....	12	M00533030	Advanix Duodenal Bend - Single Stent	21
M00530830	Ultratome Double-Lumen Sphincterotome – Long Nose.....	12	M00533040	Advanix Duodenal Bend - Single Stent	21
M00530850	Contour™ ERCP Cannulas – Standard	13	M00533230	Advanix Center Bend - Single Stent	21
M00530860	Contour ERCP Cannulas – Tapered.....	13	M00533240	Advanix Center Bend - Single Stent	21
M00530870	Contour ERCP Cannulas – Ultra Tapered.....	13	M00533250	Advanix Center Bend - Single Stent	21
M00530880	Contour ERCP Cannulas – 5-4-3	13	M00533260	Advanix Center Bend - Single Stent	21
M00530890	Contour ERCP Cannulas – Ball Tip.....	13	M00533270	Advanix Center Bend - Single Stent	21
M00530960	Fluoro Tip™ ERCP Cannulas – Standard Tip	13	M00533280	Advanix Center Bend - Single Stent	21
M00530970	Fluoro Tip ERCP Cannulas – Tapered Tip	13	M00533290	Advanix Center Bend - Single Stent	21
M00530980	Fluoro Tip ERCP Cannulas – Ultra Tapered Tip	13	M00533300	Advanix Center Bend - Single Stent	21
M00532010	C-Flex™ Double Pigtail Biliary Stent.....	25	M00533310	Advanix Center Bend - Single Stent	21
M00532020	C-Flex Double Pigtail Biliary Stent	25	M00533320	Advanix Center Bend - Single Stent	21
M00532030	C-Flex Double Pigtail Biliary Stent	25	M00533330	Advanix Center Bend - Single Stent	21
M00532040	C-Flex Double Pigtail Biliary Stent	25	M00533340	Advanix Center Bend - Single Stent	21
M00532050	C-Flex Double Pigtail Biliary Stent	25	M00533350	Advanix Center Bend - Single Stent	21
M00532160	Advanix™ Double Pigtail - Single Stent	22	M00533360	Advanix Center Bend - Single Stent	21
M00532180	Advanix Double Pigtail - Single Stent	22	M00533370	Advanix Center Bend - Single Stent	21
M00532190	Advanix Double Pigtail - Single Stent	22	M00533380	Advanix Center Bend - Single Stent	21
M00532210	Advanix Double Pigtail - Single Stent	22	M00533390	Advanix Center Bend - Single Stent	21
M00532220	Advanix Double Pigtail - Single Stent	22	M00533400	Advanix Center Bend - Single Stent	21
M00532230	Advanix Double Pigtail - Single Stent	22	M00533560	NaviFlex™ RX Delivery System.....	22
M00532240	Advanix Double Pigtail - Single Stent	22	M00533570	NaviFlex RX Delivery System	22
M00532260	Advanix Double Pigtail - Single Stent	22	M00533580	NaviFlex RX Delivery System	22
M00532270	Advanix Double Pigtail - Single Stent	22	M00533610	Percuflex™ Duodenal Bend Biliary Stent.....	24
M00532290	Advanix Double Pigtail - Single Stent	22	M00533620	Percuflex Duodenal Bend Biliary Stent.....	24
M00532300	Advanix Double Pigtail - Single Stent	22	M00533630	Percuflex Duodenal Bend Biliary Stent.....	24
M00532310	Advanix Double Pigtail - Single Stent	22	M00533640	Percuflex Duodenal Bend Biliary Stent.....	24
M00532810	MicroKnife™ XL Triple-Lumen Needle Knife	10	M00533650	Percuflex Duodenal Bend Biliary Stent.....	24
			M00533660	Percuflex Duodenal Bend Biliary Stent.....	24
			M00533670	Percuflex Duodenal Bend Biliary Stent.....	24



Product Code Index

Product Code	Description	Page #	Product Code	Description	Page #
M00533680	Percuflex™ Duodenal Bend Biliary Stent.....	24	M00534650	Advanix Center Bend Preloaded RX Stent Systems	20
M00533690	Percuflex Duodenal Bend Biliary Stent.....	24	M00534660	Advanix Center Bend Preloaded RX Stent Systems	20
M00533700	Percuflex Duodenal Bend Biliary Stent.....	24	M00534670	Advanix Center Bend Preloaded RX Stent Systems	20
M00533760	Percuflex Duodenal Bend Biliary Stent w/ Push Catheter.....	24	M00534680	Advanix Center Bend Preloaded RX Stent Systems	20
M00533770	Percuflex Duodenal Bend Biliary Stent w/ Push Catheter	24	M00534690	Advanix Center Bend Preloaded RX Stent Systems	20
M00533780	Percuflex Duodenal Bend Biliary Stent w/ Push Catheter.....	24	M00534700	Advanix Center Bend Preloaded RX Stent Systems	20
M00533790	Percuflex Duodenal Bend Biliary Stent w/ Push Catheter.....	24	M00534710	Advanix Center Bend Preloaded RX Stent Systems	20
M00533800	Percuflex Duodenal Bend Biliary Stent w/ Push Catheter.....	24	M00534720	Advanix Center Bend Preloaded RX Stent Systems	20
M00533810	Percuflex Duodenal Bend Biliary Stent w/ Push and Guide Catheter.....	24	M00534730	Advanix Center Bend Preloaded RX Stent Systems	20
M00533820	Percuflex Duodenal Bend Biliary Stent w/ Push and Guide Catheter.....	24	M00535000	NaviFlex™ Pancreatic Stent RX Pusher	31
M00533830	Percuflex Duodenal Bend Biliary Stent w/ Push and Guide Catheter.....	24	M00535010	NaviFlex Pancreatic Stent RX Pusher.....	31
M00533840	Percuflex Duodenal Bend Biliary Stent w/ Push and Guide Catheter.....	24	M00535020	NaviFlex Pancreatic Stent Pusher.....	31
M00533850	Percuflex Duodenal Bend Biliary Stent w/ Push and Guide Catheter.....	24	M00535030	NaviFlex Pancreatic Stent Pusher.....	31
M00533910	Biliary Introducer – Push Catheter.....	24	M00535060	NaviFlex RX Pancreatic Stent Delivery System.....	31
M00533920	Biliary Introducer – Push and Guide Catheter	24	M00535080	NaviFlex RX Pancreatic Stent Delivery System.....	31
M00534200	Advanix™ Duodenal Bend Preloaded RX Stent Systems.....	20	M00535110	Stonetome™ Stone Removal Device	11
M00534210	Advanix Duodenal Bend Preloaded RX Stent Systems.....	20	M00535130	Stonetome Stone Removal Device.....	11
M00534220	Advanix Duodenal Bend Preloaded RX Stent Systems.....	20	M00535150	Stonetome Stone Removal Device.....	11
M00534230	Advanix Duodenal Bend Preloaded RX Stent Systems.....	20	M00535170	Stonetome Stone Removal Device	11
M00534240	Advanix Duodenal Bend Preloaded RX Stent Systems.....	20	M00535190	Stonetome Stone Removal Device	11
M00534250	Advanix Duodenal Bend Preloaded RX Stent Systems.....	20	M00535210	Stonetome Stone Removal Device.....	11
M00534260	Advanix Duodenal Bend Preloaded RX Stent Systems.....	20	M00535700	Tandem™ XL Triple-Lumen ERCP Cannula – Tapered	13
M00534270	Advanix Duodenal Bend Preloaded RX Stent Systems.....	20	M00535900	Ultratome™ XL Triple-Lumen Sphincterotome – Short Nose.....	11
M00534280	Advanix Duodenal Bend Preloaded RX Stent Systems.....	20	M00535910	Ultratome XL Triple-Lumen Sphincterotome – Long Nose	11
M00534290	Advanix Duodenal Bend Preloaded RX Stent Systems.....	20	M00535920	Ultratome XL Triple-Lumen Sphincterotome – Short Nose	11
M00534300	Advanix Duodenal Bend Preloaded RX Stent Systems.....	20	M00535930	Ultratome XL Triple-Lumen Sphincterotome – Long Nose	11
M00534310	Advanix Duodenal Bend Preloaded RX Stent Systems.....	20	M00536000	Advanix Pancreatic Stent – Straight, No Leading Barb.....	27
M00534320	Advanix Duodenal Bend Preloaded RX Stent Systems.....	20	M00536010	Advanix Pancreatic Stent – Straight, No Leading Barb.....	27
M00534330	Advanix Duodenal Bend Preloaded RX Stent Systems.....	20	M00536020	Advanix Pancreatic Stent – Straight, No Leading Barb.....	27
M00534340	Advanix Duodenal Bend Preloaded RX Stent Systems.....	20	M00536030	Advanix Pancreatic Stent – Straight, No Leading Barb.....	27
M00534350	Advanix Duodenal Bend Preloaded RX Stent Systems.....	20	M00536040	Advanix Pancreatic Stent – Straight, No Leading Barb.....	27
M00534360	Advanix Duodenal Bend Preloaded RX Stent Systems.....	20	M00536050	Advanix Pancreatic Stent – Straight, No Leading Barb.....	27
M00534370	Advanix Duodenal Bend Preloaded RX Stent Systems.....	20	M00536060	Advanix Pancreatic Stent – Straight, No Leading Barb.....	27
M00534560	Advanix Center Bend Preloaded RX Stent Systems	20	M00536070	Advanix Pancreatic Stent – Straight, No Leading Barb.....	27
M00534570	Advanix Center Bend Preloaded RX Stent Systems	20	M00536080	Advanix Pancreatic Stent – Straight, No Leading Barb.....	27
M00534580	Advanix Center Bend Preloaded RX Stent Systems	20	M00536090	Advanix Pancreatic Stent – Straight, No Leading Barb.....	27
M00534590	Advanix Center Bend Preloaded RX Stent Systems	20	M00536100	Advanix Pancreatic Stent – Straight, No Leading Barb.....	27
M00534600	Advanix Center Bend Preloaded RX Stent Systems	20	M00536110	Advanix Pancreatic Stent – Straight, No Leading Barb.....	27
M00534610	Advanix Center Bend Preloaded RX Stent Systems	20	M00536120	Advanix Pancreatic Stent – Straight, No Leading Barb.....	27
M00534620	Advanix Center Bend Preloaded RX Stent Systems	20	M00536130	Advanix Pancreatic Stent – Straight, No Leading Barb.....	27
M00534630	Advanix Center Bend Preloaded RX Stent Systems	20	M00536140	Advanix Pancreatic Stent – Straight, No Leading Barb.....	27
M00534640	Advanix Center Bend Preloaded RX Stent Systems	20	M00536150	Advanix Pancreatic Stent – Straight, No Leading Barb.....	27
			M00536160	Advanix Pancreatic Stent – Straight, No Leading Barb.....	27
			M00536170	Advanix Pancreatic Stent – Straight, No Leading Barb.....	27
			M00536180	Advanix Pancreatic Stent – Straight, No Leading Barb.....	27



Product Code Index

[illegible]

Product Code Index

Product Code	Description	Page #	Product Code	Description	Page #
M00537110	Advanix™ Pancreatic Stent – Pigtail, Leading Barb	30	M00539210	Flexima Duodenal Bend Biliary Stent with Delivery System	23
M00537120	Advanix Pancreatic Stent – Pigtail, Leading Barb	30	M00539220	Flexima Duodenal Bend Biliary Stent with Delivery System	23
M00537130	Advanix Pancreatic Stent – Pigtail, Leading Barb	30	M00539230	Flexima Duodenal Bend Biliary Stent with Delivery System	23
M00537140	Advanix Pancreatic Stent – Pigtail, Leading Barb	30	M00539240	Flexima Duodenal Bend Biliary Stent with Delivery System	23
M00537150	Advanix Pancreatic Stent – Pigtail, Leading Barb	30	M00539250	Flexima Duodenal Bend Biliary Stent with Delivery System	23
M00537160	Advanix Pancreatic Stent – Pigtail, Leading Barb	30	M00539260	Flexima Duodenal Bend Biliary Stent with Delivery System	23
M00537170	Advanix Pancreatic Stent – Pigtail, Leading Barb	30	M00539270	Flexima Duodenal Bend Biliary Stent with Delivery System	23
M00537180	Advanix Pancreatic Stent – Pigtail, Leading Barb	30	M00539280	Flexima Duodenal Bend Biliary Stent with Delivery System	23
M00537190	Advanix Pancreatic Stent – Pigtail, Leading Barb	30	M00539290	Flexima Duodenal Bend Biliary Stent with Delivery System	23
M00537200	Advanix Pancreatic Stent – Pigtail, Leading Barb	30	M00539300	Flexima Duodenal Bend Biliary Stent with Delivery System	23
M00537210	Advanix Pancreatic Stent – Pigtail, Leading Barb	30	M00539310	Flexima Duodenal Bend Biliary Stent with Delivery System	23
M00537220	Advanix Pancreatic Stent – Pigtail, Leading Barb	30	M00539320	Flexima Duodenal Bend Biliary Stent with Delivery System	23
M00537230	Advanix Pancreatic Stent – Pigtail, Leading Barb	30	M00539330	Flexima Duodenal Bend Biliary Stent with Delivery System	23
M00537240	Advanix Pancreatic Stent – Pigtail, Leading Barb	30	M00539340	Flexima Duodenal Bend Biliary Stent with Delivery System	23
M00537250	Advanix Pancreatic Stent – Pigtail, Leading Barb	30	M00539350	Flexima Duodenal Bend Biliary Stent with Delivery System	23
M00537260	Advanix Pancreatic Stent – Pigtail, Leading Barb	30	M00539360	Flexima Duodenal Bend Biliary Stent with Delivery System	23
M00537270	Advanix Pancreatic Stent – Pigtail, Leading Barb	30	M00539370	Flexima Duodenal Bend Biliary Stent with Delivery System	23
M00537280	Advanix Pancreatic Stent – Pigtail, Leading Barb	30	M00539380	Flexima Duodenal Bend Biliary Stent with Delivery System	23
M00537290	Advanix Pancreatic Stent – Pigtail, Leading Barb	30	M00539390	Flexima Duodenal Bend Biliary Stent with Delivery System	23
M00537300	Advanix Pancreatic Stent – Pigtail, Leading Barb	30	M00540110	Flexima Single Pigtail Nasobiliary Catheter	25
M00537310	Advanix Pancreatic Stent – Pigtail, Leading Barb	30	M00540120	Flexima Single Pigtail Nasobiliary Catheter	25
M00537320	Advanix Pancreatic Stent – Pigtail, Leading Barb	30	M00540130	Flexima Single Pigtail Nasobiliary Catheter	25
M00537330	Advanix Pancreatic Stent – Pigtail, Leading Barb	30	M00540140	Flexima Single Pigtail Nasobiliary Catheter	25
M00537340	Advanix Pancreatic Stent – Pigtail, Leading Barb	30	M00540240	Endostat™ Generator Accessories – Endostat III Dual Footswitch	54
M00537350	Advanix Pancreatic Stent Kit – Pigtail, No Leading Barb	26	M00540720	Endostat Generator Accessories – Water Bottle	54
M00537360	Advanix Pancreatic Stent Kit – Pigtail, No Leading Barb	26	M00540760	Endostat Generator Accessories – Power Cord w/ I.E.C. Molded Plug	54
M00537370	Advanix Pancreatic Stent Kit – Pigtail, No Leading Barb	26	M00540771	Endostat Generator Accessories – E-Z Load Pump Tubing – Endostat III	54
M00537380	Advanix Pancreatic Stent Kit – Pigtail, No Leading Barb	26	M00542251	Speedband Superview Super 7™ Multiple Band Ligators	54
M00537390	Advanix Pancreatic Stent Kit – Pigtail, No Leading Barb	26	M00542253	Speedband Superview Super 7 Multiple Band Ligators	54
M00537400	Advanix Pancreatic Stent Kit – Pigtail, No Leading Barb	26	M00545000	RX Cytology Brush Wireguided Cytology Brush	36
M00537410	Advanix Pancreatic Stent Kit – Pigtail, No Leading Barb	26	M00545050	Ultratome™ RX Cannulating Sphincterotome – Short Nose	8
M00537420	Advanix Pancreatic Stent Kit – Pigtail, No Leading Barb	26	M00545060	Ultratome RX Cannulating Sphincterotome – Short Nose	8
M00537430	Advanix Pancreatic Stent Kit – Pigtail, Leading Barb	26	M00545150	Autotome™ RX 49 Cannulating Sphincterotome	8
M00537440	Advanix Pancreatic Stent Kit – Pigtail, Leading Barb	26	M00545160	Autotome RX 49 Cannulating Sphincterotome	8
M00537450	Advanix Pancreatic Stent Kit – Pigtail, Leading Barb	26	M00545170	Autotome RX 44 Cannulating Sphincterotome	8
M00537460	Advanix Pancreatic Stent Kit – Straight, No Leading Barb	26	M00545180	Autotome RX 44 Cannulating Sphincterotome	8
M00537470	Advanix Pancreatic Stent Kit – Straight, No Leading Barb	26	M00545190	Autotome RX 39 Cannulating Sphincterotome	8
M00537480	Advanix Pancreatic Stent Kit – Straight, No Leading Barb	26	M00545200	Autotome RX 39 Cannulating Sphincterotome	8
M00537490	Advanix Pancreatic Stent Kit – Straight, No Leading Barb	26	M00545211	RX Locking Device	14
M00537500	Advanix Pancreatic Stent Kit – Straight, No Leading Barb	26	M00545261	RX Locking Device and Biopsy Cap System – Olympus™ Brand	14
M00537510	Advanix Pancreatic Stent Kit – Straight, No Leading Barb	26	M00545271	RX Locking Device and Biopsy Cap System – Pentax™ Brand	14
M00537520	Advanix Pancreatic Stent Kit – Straight, No Leading Barb	26	M00545281	RX Locking Device and Biopsy Cap System – Fujinon™ Brand	14
M00537530	Advanix Pancreatic Stent Kit – Straight, No Leading Barb	26			
M00537540	Advanix Pancreatic Stent Kit – Straight, Leading Barb	26			
M00537550	Advanix Pancreatic Stent Kit – Straight, Leading Barb	26			
M00537560	Advanix Pancreatic Stent Kit – Straight, Leading Barb	26			
M00539200	Flexima™ Duodenal Bend Biliary Stent with Delivery System	23			



Product Code Index

Product Code	Description	Page #	Product Code	Description	Page #
M00545471	RX Locking Device	14	M00547110	Extractor Pro XL Retrieval Balloon Catheter – Distal Injection	33
M00545550	RX Duodenal Bend Biliary Stent with RX Delivery System	23	M00547120	Extractor Pro XL Retrieval Balloon Catheter – Distal Injection	33
M00545560	RX Duodenal Bend Biliary Stent with RX Delivery System	23	M00547130	Extractor Pro XL Retrieval Balloon Catheter – Proximal Injection	33
M00545570	RX Duodenal Bend Biliary Stent with RX Delivery System	23	M00547140	Extractor Pro XL Retrieval Balloon Catheter – Proximal Injection	33
M00545580	RX Duodenal Bend Biliary Stent with RX Delivery System	23	M00547150	Extractor Pro XL Retrieval Balloon Catheter – Proximal Injection	33
M00545590	RX Duodenal Bend Biliary Stent with RX Delivery System	23	M00547200	Extractor Pro DL Retrieval Balloon Catheter	33
M00545600	RX Duodenal Bend Biliary Stent with RX Delivery System	23	M00547210	Extractor Pro DL Retrieval Balloon Catheter	33
M00545610	RX Duodenal Bend Biliary Stent with RX Delivery System	23	M00547220	Extractor Pro DL Retrieval Balloon Catheter	33
M00545620	RX Duodenal Bend Biliary Stent with RX Delivery System	23	M00547300	Extractor Pro RX-S Retrieval Balloon Catheter – Distal Injection	32
M00545630	RX Duodenal Bend Biliary Stent with RX Delivery System	23	M00547310	Extractor Pro RX-S Retrieval Balloon Catheter – Distal Injection	32
M00545640	RX Duodenal Bend Biliary Stent with RX Delivery System	23	M00547320	Extractor Pro RX-S Retrieval Balloon Catheter – Distal Injection	32
M00545650	RX Duodenal Bend Biliary Stent with RX Delivery System	23	M00547330	Extractor Pro RX-S Retrieval Balloon Catheter – Proximal Injection	32
M00545660	RX Duodenal Bend Biliary Stent with RX Delivery System	23	M00547340	Extractor Pro RX-S Retrieval Balloon Catheter – Proximal Injection	32
M00545670	RX Duodenal Bend Biliary Stent with RX Delivery System	23	M00547350	Extractor Pro RX-S Retrieval Balloon Catheter – Proximal Injection	32
M00545680	RX Duodenal Bend Biliary Stent with RX Delivery System	23	M00550000	Expect™ Endoscopic Ultrasound Aspiration Needle	43
M00545690	RX Duodenal Bend Biliary Stent with RX Delivery System	23	M00550001	Expect Endoscopic Ultrasound Aspiration Needle	43
M00545700	RX Duodenal Bend Biliary Stent with RX Delivery System	23	M00550010	Expect Endoscopic Ultrasound Aspiration Needle	43
M00545710	RX Duodenal Bend Biliary Stent with RX Delivery System	23	M00550011	Expect Endoscopic Ultrasound Aspiration Needle	43
M00545720	RX Duodenal Bend Biliary Stent with RX Delivery System	23	M00550020	Expect Endoscopic Ultrasound Aspiration Needle	43
M00545730	RX Duodenal Bend Biliary Stent with RX Delivery System	23	M00550021	Expect Endoscopic Ultrasound Aspiration Needle	43
M00545740	RX Duodenal Bend Biliary Stent with RX Delivery System	23	M00550040	Expect Endoscopic Ultrasound Aspiration Needle	43
M00545760	RX ERCP Cannula – Standard Tip	12	M00550041	Expect Endoscopic Ultrasound Aspiration Needle	43
M00545770	RX ERCP Cannula – Tapered Tip	12	M00550090	Amplatz Super Stiff™ Guidewires – Straight Tip	49
M00545780	RX ERCP Cannula – Ball Tip	12	M00550300	CRE™ Pulmonary Balloon Dilatation Catheter	85
M00545820	Rapid Exchange XL Cannula	12	M00550310	CRE Pulmonary Balloon Dilatation Catheter	85
M00545840	RX Needle Knife XL Triple-Lumen Needle Knife	8	M00550320	CRE Pulmonary Balloon Dilatation Catheter	85
M00545890	Hurricane™ RX Biliary Balloon Dilatation Catheter	35	M00550330	CRE Pulmonary Balloon Dilatation Catheter	85
M00545900	Hurricane RX Biliary Balloon Dilatation Catheter	35	M00550340	CRE Pulmonary Balloon Dilatation Catheter	85
M00545910	Hurricane RX Biliary Balloon Dilatation Catheter	35	M00550350	CRE Pulmonary Balloon Dilatation Catheter	85
M00545920	Hurricane RX Biliary Balloon Dilatation Catheter	35	M00550601	Alliance™ II Integrated Inflation / Lithotripsy Device and Alliance II Syringe / Gauge Assembly	61, 85
M00545930	Hurricane RX Biliary Balloon Dilatation Catheter	35	M00550620	Alliance II Integrated Inflation / Lithotripsy Device – Inflation Handle	61, 85
M00545940	Hurricane RX Biliary Balloon Dilatation Catheter	35	M00550630	CRE SteriFlate Disposable Inflation Device	61
M00545950	Hurricane RX Biliary Balloon Dilatation Catheter	35	M00551591	Pathfinder™ Exchange Guidewire – Stiff Shaft	5
M00545960	Hurricane RX Biliary Balloon Dilatation Catheter	35	M00552001	NovaGold™ High Performance Guidewire	5
M00546030	SpyGlass™ Direct Visualization System – Probe	40	M00552011	NovaGold High Performance Guidewire	5
M00546230	SpyScope™ Access & Delivery Catheter	40	M00553200	Inflation / Deflation Devices and Accessories – Pneumatic Inflator	61
M00546270	SpyBite™ Biopsy Forceps	39, 40	M00553640	AXIOS™ Stent and Electrocautery Enhanced Delivery System	44
M00546451	SpyGlass Direct Visualization System – Irrigation Tube Set	39, 40	M00553650	AXIOS Stent and Electrocautery Enhanced Delivery System	44
M00546600	SpyScope DS Access & Delivery Catheter	39	M00553740	AXIOS Stent and Delivery System	44
M00546650	SpyGlass DS Digital Controller	39	M00553750	AXIOS Stent and Delivery System	44
M00547000	Extractor™ Pro RX Retrieval Balloon Catheter – Distal Injection	32			
M00547010	Extractor Pro RX Retrieval Balloon Catheter – Distal Injection	32			
M00547020	Extractor Pro RX Retrieval Balloon Catheter – Distal Injection	32			
M00547030	Extractor Pro RX Retrieval Balloon Catheter – Proximal Injection	32			
M00547040	Extractor Pro RX Retrieval Balloon Catheter – Proximal Injection	32			
M00547050	Extractor Pro RX Retrieval Balloon Catheter – Proximal Injection	32			
M00547100	Extractor Pro XL Retrieval Balloon Catheter – Distal Injection	33			



Product Code Index

Product Code	Description	Page #	Product Code	Description	Page #
M00554500	Rigiflex™ II Single-Use Achalasia Balloon Dilators.....	60	M00556531	Jagwire High Performance Guidewire – Stiff Shaft – Angled ...	4
M00554510	Rigiflex II Single-Use Achalasia Balloon Dilators.....	60	M00556561	Jagwire High Performance Guidewire – Stiff Shaft – Straight..	4
M00554520	Rigiflex II Single-Use Achalasia Balloon Dilators.....	60	M00556581	Jagwire High Performance Guidewire – Stiff Shaft – Straight..	4
M00555500	Expect™ Slimline (SL) Endoscopic Ultrasound Aspiration Needle..	43	M00556591	Jagwire High Performance Guidewire – Stiff Shaft – Angled ...	4
M00555510	Expect Slimline (SL) Endoscopic Ultrasound Aspiration Needle ..	43	M00556601	Jagwire High Performance Guidewire – Stiff Shaft – Straight.....	4
M00555520	Expect Slimline (SL) Endoscopic Ultrasound Aspiration Needle ..	43	M00556611	Jagwire High Performance Guidewire – Stiff Shaft – Angled	4
M00555530	Expect Slimline (SL) Endoscopic Ultrasound Aspiration Needle ..	43	M00556621	Jagwire High Performance Guidewire – Straight	4
M00555501	Expect Slimline (SL) Endoscopic Ultrasound Aspiration Needle ..	43	M00556631	Jagwire High Performance Guidewire – Stiff Shaft – Straight.....	4
M00555511	Expect Slimline (SL) Endoscopic Ultrasound Aspiration Needle ..	43	M00556641	Jagwire High Performance Guidewire – Straight	4
M00555521	Expect Slimline (SL) Endoscopic Ultrasound Aspiration Needle ..	43	M00556651	Jagwire High Performance Guidewire – Angled.....	4
M00555531	Expect Slimline (SL) Endoscopic Ultrasound Aspiration Needle ..	43	M00556661	Jagwire High Performance Guidewire – Stiff Shaft – Angled	4
M00555540	Acquire™ Endoscopic Ultrasound Fine Needle Biopsy (FNB) Device..	44	M00556841	Jagwire Extendable High Performance Guidewire – Straight....	4
M00555541	Acquire Endoscopic Ultrasound Fine Needle Biopsy (FNB) Device....	44	M00556851	Jagwire Extendable High Performance Guidewire – Angled	4
M00555560	Acquire Endoscopic Ultrasound Fine Needle Biopsy (FNB) Device....	44	M00556861	Jagwire Extendable High Performance Guidewire – Stiff Shaft – Straight.....	4
M00555561	Acquire Endoscopic Ultrasound Fine Needle Biopsy (FNB) Device....	44	M00556871	Jagwire Extendable High Performance Guidewire – Stiff Shaft – Angled.....	4
M00556001	Hydra Jagwire™ High Performance Guidewire – Straight	3	M00556902	Jagtail™ High Performance Guidewire Extension.....	4
M00556011	Hydra Jagwire High Performance Guidewire – Angled	3	M00557330	Ultraflex™ Precision Single-Use Colonic Stent System.....	46
M00556021	Hydra Jagwire High Performance Guidewire – Stiff Shaft – Straight.....	3	M00557360	Ultraflex Precision Single-Use Colonic Stent System.....	46
M00556031	Hydra Jagwire High Performance Guidewire – Stiff Shaft – Angled	3	M00557380	Ultraflex Precision Single-Use Colonic Stent System.....	46
M00556041	Hydra Jagwire High Performance Guidewire – Straight.....	3	M00558220	Expect Pulmonary Endobronchial Ultrasound Transbronchial Aspiration Needle.....	83
M00556051	Hydra Jagwire High Performance Guidewire – Angled	3	M00558250	Expect Pulmonary Endobronchial Ultrasound Transbronchial Aspiration Needle.....	83
M00556061	Hydra Jagwire High Performance Guidewire – Stiff Shaft – Straight.....	3	M00558330	CRE™ Fixed Wire Esophageal Balloon Dilatation Catheters ...	59
M00556071	Hydra Jagwire High Performance Guidewire – Stiff Shaft – Angled	3	M00558340	CRE Fixed Wire Esophageal Balloon Dilatation Catheters.....	59
M00556101	Dreamwire™ High Performance Guidewire – Standard Shaft ..	3	M00558350	CRE Fixed Wire Esophageal Balloon Dilatation Catheters.....	59
M00556111	Dreamwire High Performance Guidewire – Standard Shaft	3	M00558360	CRE Fixed Wire Esophageal Balloon Dilatation Catheters.....	59
M00556121	Dreamwire High Performance Guidewire – Stiff Shaft.....	3	M00558361	CRE Fixed Wire Esophageal Balloon Dilatation Catheters.....	59
M00556131	Dreamwire High Performance Guidewire – Stiff Shaft.....	3	M00558370	CRE Fixed Wire Esophageal Balloon Dilatation Catheters.....	59
M00556141	Dreamwire High Performance Guidewire – Standard Shaft	3	M00558371	CRE Fixed Wire Esophageal Balloon Dilatation Catheters.....	59
M00556151	Dreamwire High Performance Guidewire – Standard Shaft	3	M00558380	CRE Fixed Wire Esophageal Balloon Dilatation Catheters.....	59
M00556161	Dreamwire High Performance Guidewire – Stiff Shaft.....	3	M00558381	CRE Fixed Wire Esophageal Balloon Dilatation Catheters.....	59
M00556171	Dreamwire High Performance Guidewire – Stiff Shaft.....	3	M00558390	CRE Wireguided Esophageal / Pyloric Balloon Dilatation Catheters	58
M00556181	NaviPro™ Hydrophilic Guidewire.....	5	M00558400	CRE Wireguided Esophageal / Pyloric Balloon Dilatation Catheters	58
M00556191	NaviPro Hydrophilic Guidewire.....	5	M00558410	CRE Wireguided Esophageal / Pyloric Balloon Dilatation Catheters	58
M00556201	NaviPro Hydrophilic Guidewire.....	5	M00558420	CRE Wireguided Esophageal / Pyloric Balloon Dilatation Catheters	58
M00556211	NaviPro Hydrophilic Guidewire.....	5	M00558430	CRE Wireguided Esophageal / Pyloric Balloon Dilatation Catheters	58
M00556221	NaviPro Hydrophilic Guidewire.....	5	M00558440	CRE Wireguided Esophageal / Pyloric Balloon Dilatation Catheters	58
M00556231	NaviPro Hydrophilic Guidewire.....	5			
M00556241	NaviPro Hydrophilic Guidewire.....	5			
M00556251	NaviPro Hydrophilic Guidewire.....	5			
M00556461	Jagwire™ High Performance Guidewire – Straight	4			
M00556481	Stiff Shaft Jagwire™ Guidewire – Straight	4			
M00556491	Stiff Shaft Jagwire Guidewire – Angled.....	4			
M00556511	Jagwire High Performance Guidewire – Stiff Shaft – Angled ...	4			
M00556521	Jagwire High Performance Guidewire – Stiff Shaft – Straight..	4			



Product Code Index

Product Code	Description	Page #	Product Code	Description	Page #
M00558450	CRE™ Wireguided Esophageal / Pyloric / Colonic Balloon Dilatation Catheters	58	M00560101	Gold Probe Single-Use Electrohemostasis Catheters – Standard Plug	52
M00558460	CRE Wireguided Esophageal / Pyloric / Colonic Balloon Dilatation Catheters	58	M00560150	Injection Gold Probe™ Single-Use Electrohemostasis / Injection Catheters – Standard Plug	52
M00558470	CRE Wireguided Esophageal / Pyloric / Colonic Balloon Dilatation Catheters	58	M00560160	Injection Gold Probe Single-Use Electrohemostasis / Injection Catheters – Standard Plug	52
M00558480	CRE Wireguided Esophageal / Pyloric / Colonic Balloon Dilatation Catheters	58	M00560220	Gold Probe Single-Use Electrohemostasis Catheters – Standard Plug	52
M00558490	CRE Wireguided Esophageal / Pyloric / Colonic Balloon Dilatation Catheters	58	M00560311	Sensation™ Single-Use Medium Stiff Wire Short Throw Snares – Medium Oval	76
M00558500	CRE Wireguided Esophageal / Pyloric / Colonic Balloon Dilatation Catheters	58	M00560321	Sensation Single-Use Medium Stiff Wire Short Throw Snares – Small Oval	76
M00558600	CRE PRO Wireguided Biliary Balloon Dilators CRE PRO Wireguided Esophageal / Pyloric / Biliary Balloon Dilatation Catheters	34 57	M00560331	Sensation Single-Use Medium Stiff Wire Short Throw Snares – Medium Crescent	76
M00558610	CRE PRO Wireguided Biliary Balloon Dilatation Catheters CRE PRO Wireguided Esophageal / Pyloric / Biliary Balloon Dilatation Catheters	34 57	M00561191	Captivator™ II Single-Use Snares – Rounded Stiff	73
M00558620	CRE PRO Wireguided Biliary Balloon Dilatation Catheters CRE PRO Wireguided Esophageal / Pyloric / Biliary Balloon Dilatation Catheters	34 57	M00561192	Captivator II Single-Use Snares – Rounded Stiff	73
M00558630	CRE PRO Wireguided Biliary Balloon Dilatation Catheters CRE PRO Wireguided Esophageal / Pyloric / Biliary Balloon Dilatation Catheters	34 57	M00561193	Captivator II Single-Use Snares – Rounded Stiff	73
M00558640	CRE PRO Wireguided Biliary Balloon Dilatation Catheters CRE PRO Wireguided Esophageal / Pyloric / Biliary Balloon Dilatation Catheters	34 57	M00561221	Captivator II Single-Use Snares – Rounded Stiff	73
M00558650	CRE PRO Wireguided Biliary Balloon Dilatation Catheters CRE PRO Wireguided Esophageal / Pyloric / Biliary Balloon Dilatation Catheters	34 57	M00561222	Captivator II Single-Use Snares – Rounded Stiff	73
M00558660	CRE PRO Wireguided Biliary Balloon Dilatation Catheters CRE PRO Wireguided Esophageal / Pyloric / Colonic / Biliary Balloon Dilatation Catheters	34 57	M00561223	Captivator II Single-Use Snares – Rounded Stiff	73
M00558670	CRE PRO Wireguided Biliary Balloon Dilatation Catheters CRE PRO Wireguided Esophageal / Pyloric / Colonic / Biliary Balloon Dilatation Catheters	34 57	M00561231	Captivator II Single-Use Snares – Rounded Stiff	73
M00558680	CRE PRO Wireguided Biliary Balloon Dilatation Catheters CRE PRO Wireguided Esophageal / Pyloric / Colonic / Biliary Balloon Dilatation Catheters	34 57	M00561232	Captivator II Single-Use Snares – Rounded Stiff	73
M00558690	CRE PRO Wireguided Biliary Balloon Dilatation Catheters CRE PRO Wireguided Esophageal / Pyloric / Colonic / Biliary Balloon Dilatation Catheters	34 57	M00561233	Captivator II Single-Use Snares – Rounded Stiff	73
M00558700	CRE PRO Wireguided Biliary Balloon Dilatation Catheters CRE PRO Wireguided Esophageal / Pyloric / Colonic / Biliary Balloon Dilatation Catheters	34 57	M00561241	Captivator II Single-Use Snares – Rounded Stiff	73
M00558710	CRE PRO Wireguided Biliary Balloon Dilatation Catheters CRE PRO Wireguided Esophageal / Pyloric / Colonic / Biliary Balloon Dilatation Catheters	34 57	M00561242	Captivator II Single-Use Snares – Rounded Stiff	73
M00558731	Expect™ Pulmonary Endobronchial Ultrasound Transbronchial Aspiration Needle Adaptor	83	M00561243	Captivator II Single-Use Snares – Rounded Stiff	73
M00560070	Gold Probe™ Single-Use Electrohemostasis Catheters – Standard Plug	52	M00561250	Active Cords – Fits Olympus™ Brand Generator	14, 77
M00560071	Gold Probe Single-Use Electrohemostasis Catheters – Standard Plug	52	M00561270	Active Cord, Fits Endostat™, Endostat II and Endostat III Generators, Bovie™, Valleylab™, Aspen Labs™, and ERBE™ Brand Generators	14, 77
M00560100	Gold Probe Single-Use Electrohemostasis Catheters – Standard Plug	52	M00561291	Captivator II Single-Use Snare – Extra Large Rounded	73
			M00561311	Captivator Single-Use Snares – Large Oval Medium	74
			M00561440	Electrosurgical Bipolar Cable Adapter	55
			M00561450	Electrosurgical Bipolar Cable Adapter	55
			M00561460	Electrosurgical Bipolar Cable Adapter	55
			M00561470	Electrosurgical Bipolar Cable Adapter	55
			M00561480	Electrosurgical Bipolar Cable Adapter	55
			M00561490	Electrosurgical Bipolar Cable Adapter	55
			M00561500	Electrosurgical Bipolar Cable Adapter	55
			M00561600	Captivator EMR Device for Standard Gastrosopes	73
			M00561601	Captivator EMR Device for Standard Gastrosopes	73
			M00561610	Captivator EMR Device for Large Gastrosopes	73
			M00561611	Captivator EMR Device for Large Gastrosopes	73
			M00561821	Single-Use Rotatable Snares – Small Oval	75
			M00561831	Single-Use Rotatable Snares – Medium Oval	75
			M00562301	Captivator Single-Use Snares – Small Oval	74
			M00562321	Captivator Single-Use Snares – Medium Oval	74
			M00562341	Captivator Single-Use Snares – Medium Hexagonal	74
			M00562371	Captivator Single-Use Snares – Medium Crescent	74
			M00562391	Captivator Single-Use Snares – Large Oval	74



Product Code Index

Product Code	Description	Page #	Product Code	Description	Page #
M00562401	Captiflex™ Single-Use Snares – Medium Oval	76	M00564101	eXcelon™ Transbronchial Aspiration Needle	83
M00562402	Captiflex Single-Use Snares – Medium Oval	76	M00564111	eXcelon Transbronchial Aspiration Needle	83
M00562421	Captiflex Single-Use Snares – Small Oval	76	M00564121	eXcelon Transbronchial Aspiration Needle	83
M00562422	Captiflex Single-Use Snares – Small Oval	76	M00565010	WallFlex™ Duodenal Stent System with Anchor Lock Delivery System	45
M00562451	Captivator™ Single-Use Snares – Small Hexagonal	74	M00565020	WallFlex Duodenal Stent System with Anchor Lock Delivery System	45
M00562471	Captiflex Single-Use Snares – Extra Small Oval	76	M00565030	WallFlex Duodenal Stent System with Anchor Lock Delivery System	45
M00562531	Profile™ Single-Use Snares – Extra Small Oval	75	M00565040	WallFlex Colonic Stent System with Anchor Lock Delivery System	45
M00562551	Profile Single-Use Snares – Small Oval	75	M00565050	WallFlex Colonic Stent System with Anchor Lock Delivery System	45
M00562571	Profile Single-Use Snares – Medium Oval	75	M00565060	WallFlex Colonic Stent System with Anchor Lock Delivery System	45
M00562651	Sensation™ Single-Use Short Throw Snares – Large Oval / Medium Stiff	76	M00565070	WallFlex Colonic Stent System with Anchor Lock Delivery System	45
M00562652	Sensation Single-Use Short Throw Snares – Large Oval / Medium Stiff	76	M00565080	WallFlex Colonic Stent System with Anchor Lock Delivery System	45
M00562671	Sensation Single-Use Short Throw Snares – Medium Oval / Flexible	76	M00565090	WallFlex Colonic Stent System with Anchor Lock Delivery System	45
M00562672	Sensation Single-Use Short Throw Snares – Medium Oval / Flexible	76	M00565100	WallFlex Colonic Stent System with Anchor Lock Delivery System	45
M00562673	Sensation Single-Use Short Throw Snares – Medium Oval / Flexible	76	M00565110	WallFlex Colonic Stent System with Anchor Lock Delivery System	45
M00562691	Sensation Single-Use Short Throw Snares – Small Oval / Flexible	76	M00565120	WallFlex Colonic Stent System with Anchor Lock Delivery System	45
M00562692	Sensation Single-Use Short Throw Snares – Small Oval / Flexible	76	M00565130	WallFlex Colonic Stent System with Anchor Lock Delivery System	45
M00562693	Sensation Single-Use Short Throw Snares – Small Oval / Flexible	76	M00565140	WallFlex Colonic Stent System with Anchor Lock Delivery System	45
M00562760	EndoVive™ Low Profile Replacement Button	71	M00565150	WallFlex Colonic Stent System with Anchor Lock Delivery System	45
M00562770	EndoVive Low Profile Replacement Button	71	M00565560	WALLSTENT™ Enteral Endoprosthesis Colonic / Duodenal Stent with UNISTEP™ Plus Delivery System	46
M00562780	EndoVive Low Profile Replacement Button	71	M00565570	WALLSTENT Enteral Endoprosthesis Colonic / Duodenal Stent with UNISTEP Plus Delivery System	46
M00562790	EndoVive Low Profile Replacement Button	71	M00565580	WALLSTENT Enteral Endoprosthesis Colonic / Duodenal Stent with UNISTEP Plus Delivery System	46
M00562800	EndoVive Low Profile Replacement Button	71	M00565590	WALLSTENT Enteral Endoprosthesis Colonic / Duodenal Stent with UNISTEP Plus Delivery System	46
M00562820	EndoVive Low Profile Replacement Button	71	M00565620	WALLSTENT Enteral Endoprosthesis Colonic / Duodenal Stent with UNISTEP Plus Delivery System	46
M00562830	EndoVive Low Profile Replacement Button	71	M00565630	WALLSTENT Enteral Endoprosthesis Colonic / Duodenal Stent with UNISTEP Plus Delivery System	46
M00562840	EndoVive Low Profile Replacement Button	71	M00565640	WALLSTENT Enteral Endoprosthesis Colonic / Duodenal Stent with UNISTEP Plus Delivery System	46
M00562850	EndoVive Low Profile Replacement Button	71	M00565650	WALLSTENT Enteral Endoprosthesis Colonic / Duodenal Stent with UNISTEP Plus Delivery System	46
M00563010	EndoVive One Step Button™ Low Profile PEG Kit, Pull	68	M00566001	RX Accessories – Rapid Refill Continuous Injection System..	14
M00563020	EndoVive One Step Button Low Profile PEG Kit, Pull	68	M00566320	EndoVive Through-the-PEG (TTP) Jejunal Feeding Tube – Bent Tip	70
M00563030	EndoVive One Step Button Low Profile PEG Kit, Pull	68			
M00563040	EndoVive One Step Button Low Profile PEG Kit, Pull	68			
M00563080	EndoVive One Step Button Low Profile PEG Kit, Push	68			
M00563090	EndoVive One Step Button Low Profile PEG Kit, Push	68			
M00563100	EndoVive One Step Button Low Profile PEG Kit, Push	68			
M00563120	EndoVive One Step Button Low Profile PEG Kit, Pull	68			
M00563130	EndoVive One Step Button Low Profile PEG Kit, Pull	68			
M00563140	EndoVive One Step Button Low Profile PEG Kit, Pull	68			
M00563150	EndoVive One Step Button Low Profile PEG Kit, Pull	68			
M00563160	EndoVive One Step Button Low Profile PEG Kit, Pull	68			
M00563180	EndoVive One Step Button Low Profile PEG Kit, Push	68			
M00563190	EndoVive One Step Button Low Profile PEG Kit, Push	68			



Product Code Index

Product Code	Description	Page #	Product Code	Description	Page #
M00566330	EndoVive™ Through-the-PEG (TTP) Jejunal Feeding Tube – Pigtail Tip.....	70	M00570120	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	82
M00566340	EndoVive Through-the-PEG (TTP) Jejunal Feeding Tube – Bent Tip	70	M00570130	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	82
M00566350	EndoVive Through-the-PEG (TTP) Jejunal Feeding Tube – Pigtail Tip.....	70	M00570140	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	82
M00566461	EndoVive Safety PEG Kit – Pull.....	67	M00570150	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	82
M00566471	EndoVive Safety PEG Kit – Push.....	67	M00570160	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	82
M00566481	EndoVive Safety PEG Kit – Pull.....	67	M00570170	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	82
M00566491	EndoVive Safety PEG Kit – Push.....	67	M00570180	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	82
M00566670	Encore™ 26 Inflation Device	36	M00570190	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	82
M00567330	MaxForce™ Biliary Balloon Dilatation Catheter	36	M00570200	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	82
M00567340	MaxForce Biliary Balloon Dilatation Catheter.....	36	M00570210	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	82
M00567370	MaxForce Biliary Balloon Dilatation Catheter.....	36	M00570220	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	82
M00567380	MaxForce Biliary Balloon Dilatation Catheter.....	36	M00570230	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	82
M00567410	MaxForce Biliary Balloon Dilatation Catheter.....	36	M00570240	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	82
M00567420	MaxForce Biliary Balloon Dilatation Catheter.....	36	M00570250	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	82
M00568201	EndoVive Standard PEG Kit – Pull	67	M00570260	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	82
M00568211	EndoVive Standard PEG Kit – Push.....	67	M00570270	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	82
M00568241	EndoVive Standard PEG Kit – Pull	67	M00570280	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	82
M00568251	EndoVive Standard PEG Kit – Push.....	67	M00570290	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	82
M00569610	WALLSTENT™ RX Biliary Endoprosthesis Stent System.....	18	M00570300	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	82
M00569620	WALLSTENT RX Biliary Endoprosthesis Stent System.....	18	M00570310	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	82
M00569630	WALLSTENT RX Biliary Endoprosthesis Stent System.....	18	M00570320	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	82
M00569640	WALLSTENT RX Biliary Endoprosthesis Stent System.....	18	M00570330	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	82
M00569650	WALLSTENT RX Biliary Endoprosthesis Stent System.....	18	M00570340	WallFlex™ Biliary RX Fully Covered Stent System RMV	17
M00569660	WALLSTENT RX Biliary Endoprosthesis Stent System.....	18	M00570350	WallFlex Biliary RX Fully Covered Stent System RMV	17
M00569670	WALLSTENT RX Biliary Endoprosthesis Stent System.....	18	M00570360	WallFlex Biliary RX Fully Covered Stent System RMV	17
M00569680	WALLSTENT RX Biliary Endoprosthesis Stent System.....	18	M00570370	WallFlex Biliary RX Fully Covered Stent System RMV	17
M00569690	WALLSTENT RX Biliary Endoprosthesis Stent System.....	18	M00570380	WallFlex Biliary RX Fully Covered Stent System RMV	17
M00569700	WALLSTENT RX Biliary Endoprosthesis Stent System.....	18	M00570600	WallFlex Biliary RX Uncovered Stent System.....	17
M00569710	WALLSTENT RX Biliary Endoprosthesis Stent System.....	18	M00570610	WallFlex Biliary RX Uncovered Stent System.....	17
M00569720	WALLSTENT RX Biliary Endoprosthesis Stent System.....	18	M00570620	WallFlex Biliary RX Uncovered Stent System.....	17
M00569800	WALLSTENT RX Biliary Endoprosthesis Stent System.....	18	M00570630	WallFlex Biliary RX Uncovered Stent System.....	17
M00569810	WALLSTENT RX Biliary Endoprosthesis Stent System.....	18	M00570640	WallFlex Biliary RX Uncovered Stent System.....	17
M00570000	Polyflex™ Single-Use Self-Expanding Silicone Airway Stent..	82	M00570650	WallFlex Biliary RX Uncovered Stent System.....	17
M00570010	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	82	M00570660	WallFlex Biliary RX Uncovered Stent System.....	17
M00570020	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	82	M00570670	Dynamic™ (Y) Stent	79
M00570030	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	82	M00570680	Dynamic (Y) Stent	79
M00570040	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	82	M00570690	Dynamic (Y) Stent	79
M00570050	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	82	M00570700	WallFlex Biliary RX Partially Covered Stent System	17
M00570060	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	82	M00570710	WallFlex Biliary RX Partially Covered Stent System	17
M00570070	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	82	M00570720	WallFlex Biliary RX Partially Covered Stent System	17
M00570080	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	82	M00570730	WallFlex Biliary RX Partially Covered Stent System	17
M00570090	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	82	M00570740	WallFlex Biliary RX Partially Covered Stent System	17
M00570100	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	82	M00570890	WallFlex Biliary RX Uncovered Stent System.....	17
M00570110	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	82	M00573000	Jagtome™ RX 49 Sphincterotome (Preloaded with Jagwire™ Guidewire)	7
			M00573010	Jagtome RX 49 Sphincterotome (Preloaded with Jagwire Guidewire).....	7



Product Code Index

Product Code	Description	Page #	Product Code	Description	Page #
M00573020	Jagtome™ RX 49 Sphincterotome (Preloaded with Jagwire Guidewire).....	7	M00576490	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Proximal Release.....	80
M00573030	Jagtome RX 49 Sphincterotome (Preloaded with Jagwire Guidewire).....	7	M00576520	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release.....	81
M00573040	Jagtome RX 44 Sphincterotome (Preloaded with Jagwire Guidewire).....	7	M00576530	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release.....	81
M00573050	Jagtome RX 44 Sphincterotome (Preloaded with Jagwire Guidewire).....	7	M00576540	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release.....	81
M00573060	Jagtome RX 44 Sphincterotome (Preloaded with Jagwire Guidewire).....	7	M00576550	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release.....	81
M00573070	Jagtome RX 44 Sphincterotome (Preloaded with Jagwire Guidewire).....	7	M00576560	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release.....	81
M00573080	Jagtome RX 39 Sphincterotome (Preloaded with Jagwire Guidewire).....	7	M00576570	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release.....	81
M00573090	Jagtome RX 39 Sphincterotome (Preloaded with Jagwire Guidewire).....	7	M00576580	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release.....	81
M00573100	Jagtome RX 39 Sphincterotome (Preloaded with Jagwire Guidewire).....	7	M00576590	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release.....	81
M00573110	Jagtome RX 39 Sphincterotome (Preloaded with Jagwire Guidewire).....	7	M00576600	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release.....	81
M00576300	Ultraflex™ Single-Use Uncovered Tracheobronchial Stent System – Distal Release.....	80	M00576610	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release.....	81
M00576310	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Distal Release.....	80	M00576620	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release.....	81
M00576320	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Distal Release.....	80	M00576630	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release.....	81
M00576330	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Distal Release.....	80	M00576640	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release.....	81
M00576340	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Distal Release.....	80	M00576650	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release.....	81
M00576350	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Distal Release.....	80	M00576660	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release.....	81
M00576360	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Distal Release.....	80	M00576670	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release.....	81
M00576370	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Distal Release.....	80	M00576680	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release.....	81
M00576380	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Distal Release.....	80	M00576690	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release.....	81
M00576410	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Proximal Release.....	80	M00580171	EndoVive™ One-Step Button™ and Low Profile Replacement Button Accessories – 18Fr Bolus Feeding Set, Straight Adapter.....	68, 71
M00576420	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Proximal Release.....	80	M00580181	EndoVive One-Step Button and Low Profile Replacement Button Accessories – 18Fr Right Angle Feeding Set, Continuous Adapter.....	68, 71
M00576430	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Proximal Release.....	80	M00580231	EndoVive One-Step Button and Low Profile Replacement Button Accessories – 24Fr Bolus Feeding Set, Straight Adapter..	68, 71
M00576440	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Proximal Release.....	80	M00580241	EndoVive One-Step Button and Low Profile Replacement Button Accessories – 24Fr Right Angle Feeding Set, Continuous Adapter.....	68, 71
M00576450	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Proximal Release.....	80	M00580401	EndoVive One-Step Button and Low Profile Replacement Button Accessories – 18Fr Decompression Tube.....	68, 71
M00576460	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Proximal Release.....	80	M00580411	EndoVive One-Step Button and Low Profile Replacement Button Accessories – 18Fr Decompression Tube.....	68, 71
M00576470	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Proximal Release.....	80	M00580421	EndoVive One-Step Button and Low Profile Replacement Button Accessories – 18Fr Decompression Tube.....	68, 71
M00576480	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Proximal Release.....	80			



Product Code Index

Product Code	Description	Page #	Product Code	Description	Page #
M00580431	EndoVive™ One-Step Button™ and Low Profile Replacement Button Accessories – 18Fr Decompression Tube	68, 71	M00583060	Hydratome RX 44 Sphincterotome (Preloaded with Hydra Jagwire Guidewire)	6
M00580441	EndoVive One-Step Button and Low Profile Replacement Button Accessories – 18Fr Decompression Tube	68, 71	M00583070	Hydratome RX 44 Sphincterotome (Preloaded with Hydra Jagwire Guidewire)	6
M00580461	EndoVive One-Step Button and Low Profile Replacement Button Accessories – 24Fr Decompression Tube	68, 71	M00583080	TRUEtome™ Jag 39 Cannulating Sphincterotome (Preloaded with High Performance Jagwire™ Guidewire)	9
M00580471	EndoVive One-Step Button and Low Profile Replacement Button Accessories – 24Fr Decompression Tube	68, 71	M00583090	TRUEtome Jag 44 Cannulating Sphincterotome (Preloaded with High Performance Jagwire Guidewire)	9
M00580481	EndoVive One-Step Button and Low Profile Replacement Button Accessories – 24Fr Decompression Tube	68, 71	M00583100	TRUEtome Jag 44 Cannulating Sphincterotome (Preloaded with High Performance Jagwire Guidewire)	9
M00580491	EndoVive One-Step Button and Low Profile Replacement Button Accessories – 24Fr Decompression Tube	68, 71	M00583110	TRUEtome Jag 44 Cannulating Sphincterotome (Preloaded with High Performance Jagwire™ Guidewire)	9
M00580641	EndoVive Initial Placement PEG Accessories – Y-Ports	68	M00583120	TRUEtome Jag 39 Cannulating Sphincterotome (Preloaded with High Performance Jagwire Guidewire)	9
M00580651	EndoVive Initial Placement PEG Accessories – Y-Ports	68	M00583130	TRUEtome Jag 44 Cannulating Sphincterotome (Preloaded with High Performance Jagwire Guidewire)	9
M00580731	EndoVive Initial Placement PEG Accessories – Round Bolsters	68	M00583140	TRUEtome Jag 44 Cannulating Sphincterotome (Preloaded with High Performance Jagwire Guidewire)	9
M00580781	EndoVive Initial Placement PEG Accessories – C Clamps	68	M00583150	TRUEtome Jag 44 Cannulating Sphincterotome (Preloaded with High Performance Jagwire Guidewire)	9
M00582040	EndoVive Standard Balloon Replacement Kit, Straight Bolster	69	M00583160	TRUEtome Hydra 44 Cannulating Sphincterotome (Preloaded with High Performance Hydra Guidewire)	9
M00582050	EndoVive Standard Balloon Replacement Kit, Straight Bolster	69	M00583170	TRUEtome Hydra 44 Cannulating Sphincterotome (Preloaded with High Performance Hydra Guidewire)	9
M00582060	EndoVive Standard Balloon Replacement Kit, Straight Bolster	69	M00583180	TRUEtome Dreamwire 44 Cannulating Sphincterotome (Preloaded with High Performance Dreamwire™ Guidewire)	9
M00582070	EndoVive Standard Balloon Replacement Kit, Straight Bolster	69	M00583190	TRUEtome Dreamwire 44 Cannulating Sphincterotome (Preloaded with High Performance Dreamwire Guidewire)	9
M00582080	EndoVive Standard Balloon Replacement Kit, Straight Bolster	69	M00583200	TRUEtome Dreamwire 44 Cannulating Sphincterotome (Preloaded with High Performance Dreamwire Guidewire)	9
M00582090	EndoVive Standard Balloon Replacement Kit, Straight Bolster	69	M00583210	TRUEtome Dreamwire 44 Cannulating Sphincterotome (Preloaded with High Performance Dreamwire Guidewire)	9
M00582100	EndoVive Standard Balloon Replacement Kit, Straight Bolster	69	M00584000	Dreamtome™ RX 49 Sphincterotome (Preloaded with Dreamwire High Performance Guidewire)	6
M00582110	EndoVive Standard Balloon Replacement Kit, Right Angle Bolster	70	M00584010	Dreamtome RX 49 Sphincterotome (Preloaded with Dreamwire High Performance Guidewire)	6
M00582120	EndoVive Standard Balloon Replacement Kit, Right Angle Bolster	70	M00584020	Dreamtome RX 49 Sphincterotome (Preloaded with Dreamwire High Performance Guidewire)	6
M00582130	EndoVive Standard Balloon Replacement Kit, Right Angle Bolster	70	M00584030	Dreamtome RX 49 Sphincterotome (Preloaded with Dreamwire High Performance Guidewire)	6
M00582140	EndoVive Standard Balloon Replacement Kit, Right Angle Bolster	70	M00584040	Dreamtome RX 44 Sphincterotome (Preloaded with Dreamwire High Performance Guidewire)	6
M00582150	EndoVive Standard Balloon Replacement Kit, Right Angle Bolster	70	M00584050	Dreamtome RX 44 Sphincterotome (Preloaded with Dreamwire High Performance Guidewire)	6
M00582160	EndoVive Standard Balloon Replacement Kit, Right Angle Bolster	70	M00584060	Dreamtome RX 44 Sphincterotome (Preloaded with Dreamwire High Performance Guidewire)	6
M00582170	EndoVive Standard Balloon Replacement Kit, Right Angle Bolster	70	M00584070	Dreamtome RX 44 Sphincterotome (Preloaded with Dreamwire High Performance Guidewire)	6
M00582860	EndoVive Standard Balloon Replacement Kit, Straight Bolster	69	M00584120	TRUEtome 39 Sphincterotome	10
M00582870	EndoVive Standard Balloon Replacement Kit, Right Angle Bolster	70	M00584150	TRUEtome 49 Sphincterotome	10
M00583000	Hydratome™ RX 49 Sphincterotome (Preloaded with Hydra Jagwire™ Guidewire)	6	M00584160	TRUEtome 49 Sphincterotome	10
M00583010	Hydratome RX 49 Sphincterotome (Preloaded with Hydra Jagwire Guidewire)	6	M00584170	TRUEtome 44 Sphincterotome	10
M00583020	Hydratome RX 49 Sphincterotome (Preloaded with Hydra Jagwire Guidewire)	6	M00584180	TRUEtome 44 Sphincterotome	10
M00583030	Hydratome RX 49 Sphincterotome (Preloaded with Hydra Jagwire Guidewire)	6	M00584190	TRUEtome 39 Sphincterotome	10
M00583040	Hydratome RX 44 Sphincterotome (Preloaded with Hydra Jagwire Guidewire)	6			
M00583050	Hydratome RX 44 Sphincterotome (Preloaded with Hydra Jagwire Guidewire)	6			



Product Code Index

Product Code	Description	Page #	Product Code	Description	Page #
M00589060	MaxForce TTS™ Esophageal High Performance Balloon Dilatation Catheter	60	M005ATS20000	Alair™ Bronchial Thermoplasty System Controller.....	79
M00589080	MaxForce TTS Esophageal High Performance Balloon Dilatation Catheter	60	M005ATS25010	Alair Bronchial Thermoplasty System Catheter	79
M00589100	MaxForce TTS Esophageal High Performance Balloon Dilatation Catheter	60	M005ATS20110	Alair RF Controller Accessory Kit	79
M00589120	MaxForce TTS Esophageal High Performance Balloon Dilatation Catheter	60	M005XXX1	RX Cholangiogram Kit-Preloaded – Autotome™ Sphincterotome – RX Locking Device.....	8
M00589140	MaxForce TTS Esophageal High Performance Balloon Dilatation Catheter	60	M005XXX1	RX Cholangiogram Kits – Autotome or Ultratome™ XL Sphincterotome	8
M00589150	MaxForce TTS Esophageal High Performance Balloon Dilatation Catheter	60			
M00589160	MaxForce TTS Esophageal High Performance Balloon Dilatation Catheter	60			
M00589180	MaxForce TTS Esophageal High Performance Balloon Dilatation Catheter	60			



Terms & Conditions

▶ Begins on Page..... 105

Certain Terms and Conditions

Ordering Information

All orders should be placed directly with Boston Scientific Endoscopy Customer Service
by mail, telephone, fax or Electronic Data Interchange (EDI).

Boston Scientific Endoscopy
500 Commander Shea Blvd.
Quincy, MA 02171
(888) 272-1001 Phone
(888) 272-9444 Fax
(877) 272-3340 EDI

Please contact EDI if interested in ordering your endoscopy products via the internet.

Please note that the above address is our mailing address only.

Customer Service is available between the hours of 8:30 am and 9:30 pm (Eastern Standard Time) weekdays excluding national holidays. Our toll-free fax number, (888) 272-9444, is available at all times. In addition, technical and emergency service is available 24 hours a day, 7 days a week, call (888) 272-1001 for further instructions.

Standing Orders

Boston Scientific Endoscopy accepts standing orders on all products. This convenience provides you with scheduled shipments using a single purchase order. Contact the Customer Service Database Department at (888) 325-1060 for details.

Terms

Terms are Net 30 Days. All shipments are F.O.B. Quincy, Massachusetts. Prices are subject to change without notice.

Returned Goods

Permission to return merchandise is required for all returns and, except as provided herein, may be withheld in the sole discretion of Boston Scientific Endoscopy. This authorization must be requested by calling or writing:

Customer Service

Boston Scientific Endoscopy

500 Commander Shea Blvd.

Quincy, MA 02171

(800) 862-1284

The request for authorization must include:

- Reason for the return – “Ordered in Error,” “Shipped in Error,” etc.
- Catalog number, size, quantity, lot number
- Purchase order number, Boston Scientific Endoscopy invoice number and date of invoice

All merchandise returned to Boston Scientific Endoscopy must be accompanied by a copy of the packing slip or invoice enclosed with the original shipment. The Returned Goods Authorization (RGA) number must be clearly marked on the shipping container. Merchandise returned without an RGA number will not be accepted or considered as returned product.

Merchandise that will not be accepted for return:

- Expired product or products with less than 6 months shelf life
- Merchandise manufactured as SPECIALS or CUSTOM
- Merchandise in unsaleable condition (opened, damaged, tampered, etc.)
- Merchandise not listed in this Price List or any Supplement
- Merchandise sterilized or used by customer
- Merchandise that is not packaged properly for shipping when returned
- Incomplete kits or boxes (unless defective)
- Unauthorized returns



Certain Terms and Conditions

Returned Goods (continued)

Allowance Schedule for Authorized Returns:

Allowance Schedule for Capital Equipment:

- Capital Equipment processed in error by Boston Scientific Endoscopy will receive full credit if returned in saleable condition
- Capital Equipment cannot be returned once the customer has taken title of the equipment

Allowance Schedule for Metal and Polyflex™ Stents and Alair™ Bronchial Thermoplasty Products:

- Merchandise processed in error by Boston Scientific Endoscopy will receive full credit if returned in saleable condition
- Merchandise returned within 30 days of invoice date in saleable condition will receive 100% credit
- Merchandise returned from 31 days to 90 days in saleable condition will receive full credit less a 30% restocking charge
- Merchandise returned from 91 days to 150 days in saleable condition will receive full credit less a 50% restocking fee
- Merchandise returned past 150 days will not receive any credit

Allowance Schedule for All Other Products:

- Merchandise processed in error by Boston Scientific Endoscopy will receive full credit if returned in saleable condition
- Other merchandise returned within 30 days of invoice date in saleable condition will receive full credit less a 15% restocking charge
- Merchandise returned past 31 days in saleable condition will receive full credit less a 40% restocking charge provided there is at least 6 months left on the expiration date

Pricing and Specification Changes

As we are continually introducing products and improving our existing products, this Price List does not purport to be complete, and the equipment supplied may not agree in all details with the description and specifications shown in the catalog and / or price list. All product offerings, prices and specifications are subject to change without notice.

Notice

Indications, contraindications and instructions for use, as well as product warranty information, if applicable, can be found in the product labeling supplied with each device. Caution: Federal law restricts this device to sale by or on the order of a physician.

Warning Re: Single-Use Devices

Products marked as “single-use” are for single use only. Do not reuse, reprocess or resterilize single-use products. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and / or lead to device failure which in turn may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and / or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. **BSC assumes no liability with respect to single-use devices reused, reprocessed or resterilized and makes no warranties, express or implied, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose, with respect to any such instrument.**

©2017 Boston Scientific Corporation or its Affiliates. All rights reserved.



Customer Ordering Information

Universal Product Number		
Radial Jaw™ 3 Single-Use Max Capacity Biopsy Forceps		
Order Number	GTIN	Description
M005 1586 1	08714729096092	Radial Jaw 3 Max Capacity Biopsy Forceps
M005 1587 1	08714729085799	Radial Jaw 3 Max Capacity Biopsy Forceps
M005 1588 1	08714729085805	Radial Jaw 3 Max Capacity Biopsy Forceps
M005 1589 1	08714729096108	Radial Jaw 3 Max Capacity Biopsy Forceps
M005 1589 2	08714729718802	Radial Jaw 3 Max Capacity Biopsy Forceps

Bold type
is old order
number

For customer convenience, the Boston Scientific Endoscopy ordering system accepts the UPN (Universal Product Number) e.g., **M00515892** found on product labeling as well as Boston Scientific's old 4 digit order numbers e.g., **1589**.

Please use either one of these numbers when placing an order via phone, fax, internet, or Electronic Data Interchange (EDI).

For those customers who scan product today, Boston Scientific is moving to the GS1 GTIN barcode and labeling standards as part of our ongoing commitment to improve supply chain efficiency. For more information go to the website: www.bostonscientific.com/GS1.

Boston Scientific Endoscopy Customer Service Representatives will be happy to help you with any questions. We are eager to make the ordering procedure as simple and efficient as possible.

All trademarks are the property of their respective owners.



Boston Scientific Corporation
300 Boston Scientific Way
Marlborough, MA 01752-1234
<http://www.bostonscientific.com/gastro>

Ordering Information
1.888.272.1001

©2017 Boston Scientific Corporation
or its affiliates. All rights reserved.

ENDO-443010-AA E January 2017



**Magnetic Resonance
Conditional¹**

¹Product with this symbol may be used in Magnetic Resonance environment
according to the conditions described in Directions for Use.

Indications, Contraindications, Warnings and Instructions for Use can be found in the product labeling supplied with each device.
CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.