







Endoscopy

**Boston
Scientific**
Advancing science for life™



Future Product Updates

Product Guide

2018
Reference Guide



TAP OR CLICK ON PRODUCT NAME FOR INFORMATION

Interactive Product Guide

12
79 Active Cords

41 Acquire™ Endoscopic Ultrasound Fine Needle Biopsy (FNB) Device

20 Advanix™ Center Bend - Single Stent

19 Advanix Center Bend Preloaded RX Stent Systems

21 Advanix Double Pigtail - Single Stent

20 Advanix Duodenal Bend - Single Stent

19 Advanix Duodenal Bend Preloaded RX Stent Systems

26 Advanix Pancreatic Stent - Straight, No Leading Barb

27 Advanix Pancreatic Stent - Straight, Leading Barb

28 Advanix Pancreatic Stent - Pigtail, No Leading Barb

29 Advanix Pancreatic Stent - Pigtail, Leading Barb

25 Advanix Pancreatic Stent Kit - Pigtail, Leading Barb and No Leading Barb

25 Advanix Pancreatic Stent Kit - Straight, No Leading Barb and Leading Barb

85 Alair™ Bronchial Thermoplasty System

60
91 Alliance™ II Integrated Inflation / Lithotripsy Device

60
91 Alliance II Single-Use Syringe / Gauge Assembly

47
91 Amplatz Super Stiff™ Guidewire

38 Autolith™ Touch Electrohydraulic Lithotripter

38 Autolith Touch Extender Cable

38 Autolith Touch Foot Switch

6 Autotome™ RX Cannulating Sphincterotome

42 AXIOS™ Stent and Electrocautery Enhanced Delivery System

42 AXIOS Stent and Delivery System

38 Biliary EHL Probe

23 Biliary Introducer

81 Blox™ Endoscopic Bite Block with Elastic Cloth Strap (38 FR)

81 Blox Endoscopic Bite Block with Elastic Cloth Strap (54 FR)

81 Blox Endoscopic Bite Block with Elastic Cloth Strap (54 FR) (w/o Dental Retention)

81 Blox™ Endoscopic Bite Block with Elastic Cloth Strap (60 FR)

24 C-Flex™ Double Pigtail Biliary Stent

78 Captiflex™ Single-Use Snares

75 Captivator™ II Single-Use Snares

76 Captivator Single-Use Snares

75 Captivator EMR Device

90 Cellebrity™ Single-Use Endoscopic Cytology Brush

83 CinchPad™ Endoscope Transport System

83 Compliance EndoKit™ - Colon Kit

83 Compliance EndoKit - EGD Kit

11 Contour™ ERCP Cannulas

58 CRE™ Fixed Wire Esophageal Balloon Dilatation Catheters

33 CRE PRO Wireguided Balloon Dilatation Catheters

56 CRE PRO Wireguided Esophageal / Pyloric / Biliary Balloon Dilatation Catheters

56 CRE PRO Wireguided Esophageal / Pyloric / Colonic / Biliary Balloon Dilatation Catheters

60 CRE SteriFlate Disposable Inflation Device

57 CRE Wireguided Esophageal / Pyloric / Balloon Dilatation Catheters

Note: RX devices are short wire devices.

Prev Page



Next Page

Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at **1.888.272.1001**



TAP OR CLICK ON PRODUCT NAME FOR INFORMATION

Interactive Product Guide

57 CRE™ Wireguided Esophageal / Pyloric / Colonic / Balloon Dilatation Catheters

91 CRE Pulmonary Balloon Dilatation Catheters

55 CRE RX Biliary Balloon Dilatation Catheters

68 Distal End ENFit™ Transition Connector

4 Dreamtome™ RX Cannulating Sphincterotome

1 Dreamwire™ High Performance Guidewire

85 Dynamic™ (Y) Stent

53 Electrosurgical Bipolar Cable Adapter

35 Encore™ 26 Inflation Device

82 EndoArmor™ Procedure Gown

82 EndoArmor™+ Procedure Gown

82 Endoglide™+ Lubricant

69 EndoVive Initial Placement PEG Accessories

65 EndoVive Initial Placement PEG Accessories with ENFit

72 EndoVive Low Profile Replacement Button

72 EndoVive Low Profile Replacement Button Accessories

67 EndoVive Low Profile Replacement Button with ENFit

68 EndoVive Low Profile Replacement Button Accessories with ENFit

70 EndoVive One-Step Button Low Profile PEG Kits

66 EndoVive One-Step Button™ Low Profile PEG with ENFit

66 EndoVive One-Step Button and Low Profile Replacement Accessories with ENFit

70 EndoVive One-Step Button™ and Low Profile Replacement Button Accessories

71 EndoVive™ Replacement Gastrostomy Feeding Tubes, Straight and Right Angle Bolster

67 EndoVive Replacement Gastrostomy Feeding Tubes with ENFit, Straight and Right Angle Bolster

69 EndoVive Safety PEG Kit

65 EndoVive Safety PEG Kit with ENFit™

65 EndoVive Standard PEG Kit with ENFit

69 EndoVive Standard PEG Kit

73 EndoVive Through-the-PEG (TTP) Jejunal Feeding Tube

68 ENFit Transition Connectors

16 Epic™ Endoscopic Biliary Stent System

89 eXcelon™ Transbronchial Aspiration Needle

89 Expect Pulmonary Endobronchial Ultrasound Transbronchial Aspiration Needle / Adaptor

41 Expect Slimline (SL) Endoscopic Ultrasound Aspiration Needle

32 Extractor™ Pro DL Retrieval Balloon Catheter

31 Extractor Pro RX Retrieval Balloon Catheter – Distal Injection

31 Extractor Pro RX Retrieval Balloon Catheter – Proximal Injection

31 Extractor Pro RX-S Retrieval Balloon Catheter – Distal Injection

31 Extractor Pro RX-S Retrieval Balloon Catheter – Proximal Injection

32 Extractor Pro XL Retrieval Balloon Catheter – Distal Injection

32 Extractor Pro XL Retrieval Balloon Catheter – Proximal Injection

22 Flexima™ Duodenal Bend Biliary Stent with Delivery System

24 Flexima Single Pigtail Nasobiliary Catheter

50 Gold Probe™ Single-Use Electrohemostasis Catheters

81 HedgeHog™ DedgeHog Double Bundle (Double-End Brush + Valve Brush)

Note: RX devices are short wire devices.

Prev Page



Next Page

Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at **1.888.272.1001**

Interactive Product Guide

81 HedgeHog™ Double-End Channel Brush	2 Jagtail™ High Performance Guidewire Extension	3 Percuflex™ Duodenal Bend Biliary Stent With Introducer Kit
81 HedgeHog Valve Brush	5 Jagtome™ RX Cannulating Sphincterotome	51 Percuflex Duodenal Bend Biliary Stent
81 HedgeHog Dual-End Channel / Valve Brush	2 Jagwire™ High Performance Guidewire	83 Polyflex™ Single-Use Esophageal Stent System
81 HedgeHog Single-End Channel Brush	83 Koala™ Cleaning Sponge Plus Pure Enzymatic	23 Polyflex Single-Use Self-Expanding Silicone Airway Stent System
81 HedgeHog Sweeper Channel Brush	83 Koala Cleaning Sponge	23 Profile™ Single-Use Snares
34 Hurricane™ RX Biliary Balloon Dilatation Catheter	59 MaxForce TTS™ Esophageal High Performance Balloon Dilatation Catheter	46 Pulmonary Jagwire™ Pulmonary Guidewire
82 Hydra™ Disposable Auxiliary Water Jet Connector	35 MaxForce™ Biliary Balloon Dilatation Catheter	88 Pure™ Enzymatic Detergent - Fragrance Free
82 Hydra Water Bottle Cap	8 MicroKnife™ XL Triple-Lumen Needle Knife	77 Pure Enzymatic Detergent - Spearmint Fresh Scent
82 Hydra Irrigation Tubing	63 Multibite™ Single-Use Multiple Sample Biopsy Forceps	91 Radial Jaw™ 4 Standard Capacity Single-Use Biopsy Forceps
1 Hydra Jagwire™ High Performance Guidewire	21 NaviFlex™ RX Delivery System	81 Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps
4 Hydratome™ RX Cannulating Sphincterotome	30 NaviFlex RX Pancreatic Delivery System & NaviFlex Pushers	81 Radial Jaw 4 Single-Use Hot Biopsy Forceps
60 Inflation / Deflation Devices and Accessories	3 NaviPro™ Hydrophilic Guidewire	61 Radial Jaw 4 Single-Use Jumbo Biopsy Forceps
50 Injection Gold Probe™ Single-Use Electrohemostasis / Injection Catheters	3 NovaGold™ High Performance Guidewire	61 Radial Jaw 4 Single-Use Pediatric Biopsy Forceps
51 Interject™ Clear Single-Use Injection Therapy Needle Catheters	51 Optiflo™ Single-Use Hemostasis Catheters	62 Radial Jaw 4 Single-Use Pulmonary Biopsy Forceps
51 Interject Contrast Single-Use Injection Therapy Needle Catheters	83 Orca™ Sterile Single-Use Air/Water & Suction Valves	61 Rapid Exchange XL Cannula

Note: RX devices are short wire devices.





Interactive Product Guide

12 Rapid Refill Continuous Injection System

49 Resolution 360™ Clip

80 Rescue Alligator Grasping Forceps, Long

80 Rescue Alligator Grasping Forceps

80 Rescue Rat Tooth / Alligator Grasping Forceps

80 RescueNet™ Retrieval Net

80 Rescue™ Rat Tooth Grasping Forceps

80 Rescue Rat Tooth Grasping Forceps - Pediatric

49 Resolution™ Clip

59 RigiFlex™ II Single-Use Achalasia Balloon Dilatation Catheter

6 RX Cholangiogram Kits

35 RX Cytology Brush Wireguided Cytology Brush

22 RX Duodenal Bend Biliary Stent with RX Delivery System

12 RX Locking Devices and Biopsy Cap Systems

6 RX Needle Knife XL Triple-Lumen Needle Knife

10 RX ERCP Cannulas

83 SafeStart™ Cleaning Kit

83 SafeStart™ CP Cleaning Kit

78 Sensation™ Single-Use Short Throw Snares

77 Single-Use Rotatable Snares

52 Speedband Superview Super 7™ Multiple Band Ligators

37 SpyBite™ Biopsy Forceps

37 SpyGlass DS Direct Visualization System – Capital Equipment

37 SpyScope™ DS Access & Delivery Catheter

9 Stonetome™ Stone Removal Device

11 Tandem™ XL Triple-Lumen ERCP Cannula

80 TrapEase™ Polyp Traps

32 Trapezoid™ RX Wireguided Retrieval Baskets

7 TRUEtome™ Jag Cannulating Sphincterotome

7 TRUEtome Hydra Cannulating Sphincterotome

7 TRUEtome Dreamwire™ Cannulating Sphincterotome

8 TRUEtome Cannulating Sphincterotome

47 Ultraflex™ Single-Use Covered Large Esophageal NG Stent System – Distal Release

47 Ultraflex Single-Use Covered Large Esophageal NG Stent System – Proximal Release

47 Ultraflex Single-Use Uncovered Esophageal NG Stent System – Distal Release

47 Ultraflex Single-Use Uncovered Esophageal NG Stent System – Proximal Release

44 Ultraflex Precision Single-Use Colonic Stent System

46 Ultraflex Single-Use Covered Esophageal NG Stent System – Distal Release

46 Ultraflex Single-Use Covered Esophageal NG Stent System – Proximal Release

86 Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Distal Release

87 Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release

86 Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Proximal Release

10 Ultratome™ Double-Lumen Sphincterotome

9 Ultratome XL Triple-Lumen Sphincterotome

15 WallFlex™ Biliary RX Fully Covered Stent System RMV

Note: RX devices are short wire devices.



Interactive Product Guide

15

WallFlex™ Biliary RX Partially Covered Stent System

15

WallFlex Biliary RX Uncovered Stent System

45

WallFlex Fully Covered Esophageal Stent

45

WallFlex Partially Covered Esophageal Stent

43

WallFlex Colonic Stent System with Anchor Lock Delivery System

43

WallFlex Duodenal Stent System with Anchor Lock Delivery System

18

WALLSTENT™ Endoscopic Biliary Endoprosthesis Stents

17

WALLSTENT RX Biliary Endoprosthesis Stent System

44

WALLSTENT Enteral Endoprosthesis Colonic / Duodenal Stent with UNISTEP™ Plus Delivery System

44


WALLSTENT Super Stiff Guidewire

90

Zero Tip™ Airway Retrieval Basket



Biliary Access



Biliary Access

Guidewires

- Dreamwire™ High Performance Guidewire..... 1
- Hydra Jagwire™ High Performance Guidewire..... 1
- Jagtail™ High Performance Guidewire Extension..... 2
- Jagwire™ High Performance Guidewire 2
- NaviPro™ Hydrophilic Guidewire..... 3
- NovaGold™ High Performance Guidewire..... 3

RX Preloaded Sphincterotomes

- Dreamtome™ RX Cannulating Sphincterotome
(Preloaded with Dreamwire High Performance Guidewire)..... 4
- Hydratome™ RX Cannulating Sphincterotome
(Preloaded with Hydra Jagwire High Performance Guidewire).. 4
- Jagtome™ RX Cannulating Sphincterotome
(Preloaded with Jagwire High Performance Guidewire)..... 5

RX Sphincterotomes

- Autotome™ RX Cannulating Sphincterotome..... 6
- RX Cholangiogram Kits 6
- RX Needle Knife XL Triple-Lumen Needle Knife 6

Preloaded Sphincterotomes

- TRUEtome™ Jag Cannulating Sphincterotome
(Preloaded with Jagwire High Performance Guidewire)..... 7
- TRUEtome Hydra Cannulating Sphincterotome
(Preloaded with Hydra Jagwire High Performance Guidewire)... 7
- TRUEtome Dreamwire Cannulating Sphincterotome
(Preloaded with Dreamwire High Performance Guidewire)..... 7

Sphincterotomes

- MicroKnife™ XL Triple-Lumen Needle Knife 8
- Stonetome™ Stone Removal Device 9
- TRUEtome Cannulating Sphincterotome 8
- Ultratome Double-Lumen Sphincterotome..... 10
- Ultratome XL Triple-Lumen Sphincterotome 9

RX Cannulas

- Rapid Exchange XL Cannula 10
- RX ERCP Cannulas..... 10

Cannulas

- Contour™ ERCP Cannulas 11
- Tandem™ XL Triple-Lumen ERCP Cannula 11

Accessories

- Active Cords
(Fits Boston Scientific Endoscopy Sphincterotome) 12
- RX Locking Devices and Biopsy Cap Systems..... 12


Rapid Refill Continuous Injection System

- Rapid Refill Continuous Injection System..... 12



Biliary

Guidewires




Dreamwire™ High Performance Guidewire



Order Number	GTIN	Description	Distal Tip O.D. (in)	Working Length (cm)	Shaft Style	Dream Tip Style	Tip	Packaging
M00556101	08714729777816	Dreamwire Guidewire ST	.035	260	Standard	Straight	Round	Box 2
M00556111	08714729777847	Dreamwire Guidewire ANG	.035	260	Standard	Angled	Round	Box 2
M00556121	08714729777830	Dreamwire Guidewire ST Stiff	.035	260	Stiff	Straight	Round	Box 2
M00556131	08714729777823	Dreamwire Guidewire ANG Stiff	.035	260	Stiff	Angled	Round	Box 2
M00556141	08714729777885	Dreamwire Guidewire ST	.035	450	Standard	Straight	Round	Box 2
M00556151	08714729777878	Dreamwire Guidewire ANG	.035	450	Standard	Angled	Round	Box 2
M00556161	08714729777861	Dreamwire Guidewire ST Stiff	.035	450	Stiff	Straight	Round	Box 2
M00556171	08714729777854	Dreamwire Guidewire ANG Stiff	.035	450	Stiff	Angled	Round	Box 2

ST = Straight, ANG = Angled



Hydra Jagwire™ High Performance Guidewire



Order Number	GTIN	Description	Distal Tip O.D. (in)	Working Length (cm)	Shaft Style	Dream Tip Style	Tip	Packaging
M00556001	08714729648437	Hydra Jagwire Guidewire ST	.035	260	Standard	Straight	Round	Box 2
M00556011	08714729648444	Hydra Jagwire Guidewire ANG	.035	260	Standard	Angled	Round	Box 2
M00556021	08714729648451	Hydra Jagwire Guidewire ST Stiff	.035	260	Stiff	Straight	Round	Box 2
M00556031	08714729648468	Hydra Jagwire Guidewire ANG Stiff	.035	260	Stiff	Angled	Round	Box 2
M00556041	08714729648475	Hydra Jagwire Guidewire ST	.035	450	Standard	Straight	Round	Box 2
M00556051	08714729648482	Hydra Jagwire Guidewire ANG	.035	450	Standard	Angled	Round	Box 2
M00556061	08714729648499	Hydra Jagwire Guidewire ST Stiff	.035	450	Stiff	Straight	Round	Box 2
M00556071	08714729648505	Hydra Jagwire Guidewire ANG Stiff	.035	450	Stiff	Angled	Round	Box 2

Hydra Jagwire Guidewire may remain in place during sphincterotomy. See Instructions for Use for details.

ST = Straight, ANG = Angled



Section Table of Contents

Prev Page

1

Next Page


Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Biliary

Guidewires (Continued)



Jagwire™ High Performance Guidewire



Order Number	GTIN	Description	Distal Tip O.D. (in)	Working Length (cm)	Shaft Style	Tip Style	Tip	Packaging
M00556511	08714729874492	Jagwire Guidewire ANG	.025	450	Standard	Angled	Round	Box 2
M00556521	08714729874508	Jagwire Guidewire ST Stiff	.025	450	Stiff	Straight	Round	Box 2
M00556531	08714729874515	Jagwire Guidewire ANG Stiff	.025	450	Stiff	Angled	Round	Box 2
M00556561	08714729501374	Jagwire Guidewire ST	.025	450	Standard	Straight	Trim	Box 2
M00556581	08714729339151	Jagwire Guidewire ST	.035	450	Standard	Straight	Trim	Box 2
M00556591	08714729501398	Jagwire Guidewire ANG	.035	450	Standard	Angled	Trim	Box 2
M00556601	08714729172154	Jagwire Guidewire ST Stiff	.035	450	Stiff	Straight	Trim	Box 2
M00556611	08714729172161	Jagwire Guidewire ANG Stiff	.035	450	Stiff	Angled	Trim	Box 2
M00556621	08714729172178	Jagwire Guidewire ST	.038	260	Standard	Straight	Trim	Box 2
M00556631	08714729172185	Jagwire Guidewire ST Stiff	.035	260	Stiff	Straight	Trim	Box 2
M00556641	08714729172192	Jagwire Guidewire ST	.035	260	Standard	Straight	Trim	Box 2
M00556651	08714729172208	Jagwire Guidewire ANG	.035	260	Standard	Angled	Trim	Box 2
M00556661	08714729172215	Jagwire Guidewire ANG Stiff	.035	260	Stiff	Angled	Trim	Box 2
M00556461	08714729415183	Jagwire Guidewire ST	.025	260	Standard	Straight	Round	Box 2
M00556481	08714729415190	Jagwire Guidewire ST Stiff	.025	260	Stiff	Straight	Round	Box 2
M00556491	08714729415206	Jagwire Guidewire ANG Stiff	.025	260	Stiff	Angled	Round	Box 2
M00556841	08714729268680	Extendable Jagwire Guidewire ST	.035	260	Standard	Straight	Round	Box 2
M00556851	08714729268697	Extendable Jagwire Guidewire ANG	.035	260	Standard	Angled	Round	Box 2
M00556861	08714729268703	Extendable Jagwire Guidewire ST Stiff	.035	260	Stiff	Straight	Round	Box 2
M00556871	08714729268710	Extendable Jagwire Guidewire ANG Stiff	.035	260	Stiff	Angled	Round	Box 2
M00556902	08714729268734	Jagtail™ Guidewire Extension	.035	200	—	Guidewire extension	—	Box 5

Jagwire Guidewire may remain in place during sphincterotomy. See Instructions for Use for details.

ST = Straight, ANG = Angled



Section Table of Contents

Prev Page

2


Next Page

Main Table of Contents



Biliary


Guidewires (Continued)



NaviPro™ Hydrophilic Guidewire



Order Number	GTIN	Description	Guidewire Diameter (in)	Length (cm)	Shaft Style	Tip Style	Packaging
M00556191	08714729839873	NaviPro Hydrophilic Guidewire	0.018	260	Standard	Angled	Box 5
M00556181	08714729839866	NaviPro Hydrophilic Guidewire	0.018	260	Standard	Straight	Box 5
M00556211	08714729839897	NaviPro Hydrophilic Guidewire	0.025	260	Standard	Angled	Box 5
M00556201	08714729839880	NaviPro Hydrophilic Guidewire	0.025	260	Standard	Straight	Box 5
M00556231	08714729839910	NaviPro Hydrophilic Guidewire	0.035	260	Standard	Angled	Box 5
M00556221	08714729839903	NaviPro Hydrophilic Guidewire	0.035	260	Standard	Straight	Box 5
M00556251	08714729839934	NaviPro Hydrophilic Guidewire	0.035	260	Stiff	Angled	Box 5
M00556241	08714729839927	NaviPro Hydrophilic Guidewire	0.035	260	Stiff	Straight	Box 5



NovaGold™ High Performance Guidewire





Order Number	GTIN	Description	Distal Tip O.D. (in)	Working Length (cm)	Tip Style	Tip Length (cm)	Packaging
M00552001	08714729854418	NovaGold High Performance Guidewire	.018	260	Straight (shapeable)	6	Box 2
M00552011	08714729854425	NovaGold High Performance Guidewire	.018	480	Straight (shapeable)	6	Box 2



Biliary

RX Preloaded Sphincterotome





Dreamtome™ RX Cannulating Sphincterotome (Preloaded with Dreamwire™ High Performance Guidewire)



Order Number	GTIN	Description	Tip Length (mm)	Cut Wire Length (mm)	Preloaded Guidewire (in / cm)	Guidewire Tip Style	Packaging
M00584000	08714729778226	Dreamtome RX 49 Sphincterotome	5	20	.035 / 260 Dreamwire Guidewire	Straight	Each
M00584010	08714729778233	Dreamtome RX 49 Sphincterotome	5	30	.035 / 260 Dreamwire Guidewire	Straight	Each
M00584020	08714729778240	Dreamtome RX 49 Sphincterotome	5	20	.035 / 450 Dreamwire Guidewire	Straight	Each
M00584030	08714729778257	Dreamtome RX 49 Sphincterotome	5	30	.035 / 450 Dreamwire Guidewire	Straight	Each
M00584040	08714729778264	Dreamtome RX 44 Sphincterotome	5	20	.035 / 260 Dreamwire Guidewire	Straight	Each
M00584050	08714729778271	Dreamtome RX 44 Sphincterotome	5	30	.035 / 260 Dreamwire Guidewire	Straight	Each
M00584060	08714729778288	Dreamtome RX 44 Sphincterotome	5	20	.035 / 450 Dreamwire Guidewire	Straight	Each
M00584070	08714729778295	Dreamtome RX 44 Sphincterotome	5	30	.035 / 450 Dreamwire Guidewire	Straight	Each

*Preloaded with straight-tip guidewires only.
Active Cords sold separately (See page 12).



Hydratome™ RX Cannulating Sphincterotome (Preloaded with Hydra Jagwire™ High Performance Guidewire)



Order Number	GTIN	Description	Cut Wire (mm)	Tip Length (mm)	Distal Tip O.D. (Fr)	Preloaded Guidewire (cm)	Guidewire Tip Style	Packaging
M00583000	08714729755159	Hydratome RX 49 Sphincterotome	20	5	4.9	260 Hydra Jagwire Guidewire	Straight	Each
M00583010	08714729755166	Hydratome RX 49 Sphincterotome	30	5	4.9	260 Hydra Jagwire Guidewire	Straight	Each
M00583020	08714729755173	Hydratome RX 49 Sphincterotome	20	5	4.9	450 Hydra Jagwire Guidewire	Straight	Each
M00583030	08714729755180	Hydratome RX 49 Sphincterotome	30	5	4.9	450 Hydra Jagwire Guidewire	Straight	Each
M00583040	08714729755197	Hydratome RX 44 Sphincterotome	20	5	4.4	260 Hydra Jagwire Guidewire	Straight	Each
M00583050	08714729755203	Hydratome RX 44 Sphincterotome	30	5	4.4	260 Hydra Jagwire Guidewire	Straight	Each
M00583060	08714729755210	Hydratome RX 44 Sphincterotome	20	5	4.4	450 Hydra Jagwire Guidewire	Straight	Each
M00583070	08714729755227	Hydratome RX 44 Sphincterotome	30	5	4.4	450 Hydra Jagwire Guidewire	Straight	Each

Active Cords sold separately (See page 12).



Section Table of Contents

Prev Page

4

Next Page



Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Biliary

RX Preloaded Sphincterotome (Continued)



Jagtome™ RX Cannulating Sphincterotome (Preloaded with Jagwire™ High Performance Guidewire)



Order Number	GTIN	Description	Cut Wire (mm)	Tip Length (mm)	Distal Tip O.D. (Fr)	Preloaded Guidewire (in / cm)	Guidewire Tip Style	Packaging
M00573000	08714729776710	Jagtome RX 49 Sphincterotome	20	5	4.9	.035 / 260 Jagwire Guidewire	Straight	Each
M00573010	08714729776727	Jagtome RX 49 Sphincterotome	20	5	4.9	.035 / 450 Jagwire Guidewire	Straight	Each
M00573020	08714729776734	Jagtome RX 49 Sphincterotome	30	5	4.9	.035 / 260 Jagwire Guidewire	Straight	Each
M00573030	08714729776741	Jagtome RX 49 Sphincterotome	30	5	4.9	.035 / 450 Jagwire Guidewire	Straight	Each
M00573040	08714729776659	Jagtome RX 44 Sphincterotome	20	5	4.4	.035 / 260 Jagwire Guidewire	Straight	Each
M00573050	08714729776666	Jagtome RX 44 Sphincterotome	20	5	4.4	.035 / 450 Jagwire Guidewire	Straight	Each
M00573060	08714729776673	Jagtome RX 44 Sphincterotome	30	5	4.4	.035 / 260 Jagwire Guidewire	Straight	Each
M00573070	08714729776543	Jagtome RX 44 Sphincterotome	30	5	4.4	.035 / 450 Jagwire Guidewire	Straight	Each
M00573080	08714729776642	Jagtome RX 39 Sphincterotome	20	5	3.9	.025 / 260 Jagwire Guidewire	Straight	Each
M00573090	08714729776550	Jagtome RX 39 Sphincterotome	20	5	3.9	.025 / 450 Jagwire Guidewire	Straight	Each
M00573100	08714729776536	Jagtome RX 39 Sphincterotome	30	5	3.9	.025 / 260 Jagwire Guidewire	Straight	Each
M00573110	08714729776529	Jagtome RX 39 Sphincterotome	30	5	3.9	.025 / 450 Jagwire Guidewire	Straight	Each

Active Cords sold separately (See page 12).



Section Table of Contents

Prev Page

5

Next Page


Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Biliary

RX Sphincterotome



Autotome™ RX Cannulating Sphincterotome



Order Number	GTIN	Description	Cut Wire (mm)	Tip Length (mm)	Distal Tip O.D. (Fr)	Recommended Guidewire (in)	Packaging
M00545150	08714729415138	Autotome RX 49 Cannulating Sphincterotome	20	5	4.9	.035	Each
M00545160	08714729415152	Autotome RX 49 Cannulating Sphincterotome	30	5	4.9	.035	Each
M00545170	08714729444749	Autotome RX 44 Cannulating Sphincterotome	20	5	4.4	.035	Each
M00545180	08714729444756	Autotome RX 44 Cannulating Sphincterotome	30	5	4.4	.035	Each
M00545190	08714729444763	Autotome RX 39 Cannulating Sphincterotome	20	5	3.9	.025	Each
M00545200	08714729444770	Autotome RX 39 Cannulating Sphincterotome	30	5	3.9	.025	Each

Active Cords sold separately ➤ (See page 12).

RX Cholangiogram Kits*



Order Number	Description	Sphincterotome	Guidewire	Accessory
M005XXXX1	Cholangiogram Kit	Autotome or Ultratome XL Sphincterotome	Jagwire™ Guidewire	N / A
M005XXXX1	Cholangiogram Kit-Preloaded	Autotome Sphincterotome	Jagwire Guidewire	RX Locking Device

*To order RX Cholangiogram Kits please contact your local Boston Scientific Endoscopy Sales Representative.

RX Needle Knife XL Triple-Lumen Needle Knife



Order Number	GTIN	Description	Length (cm)	Cut Wire (mm)	Distal Tip O.D. (Fr)	Recommended Guidewire (in)	Packaging
M00545840	08714729283799	RX Needle Knife	200	5	5.5	.035	Each



Section Table of Contents

Prev Page

6

Next Page

Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Biliary

Sphincterotome

TRUEtome™ Jag Cannulating Sphincterotome (Preloaded with Jagwire™ High Performance Guidewire)									
Order Number	GTIN	Description	Tip Length (mm)	Cut Wire Length (mm)	Preloaded Guidewire (in / cm)	Tip O.D. (Fr)	Packaging		
M00583080	08714729858621	TRUEtome 39 Sphincterotome	5	20	.025 / 450 Jagwire Guidewire	3.9	Each		
M00583090	08714729858638	TRUEtome 44 Sphincterotome	5	20	.025 / 450 Jagwire Guidewire	4.4	Each		
M00583100	08714729858645	TRUEtome 44 Sphincterotome	5	20	.035 / 450 Jagwire Guidewire	4.4	Each		
M00583110	08714729858652	TRUEtome 44 Sphincterotome	5	20	.035 / Stiff 450 Jagwire Guidewire	4.4	Each		
M00583120	08714729858669	TRUEtome 39 Sphincterotome	5	30	.025 / 450 Jagwire Guidewire	3.9	Each		
M00583130	08714729858676	TRUEtome 44 Sphincterotome	5	30	.025 / 450 Jagwire Guidewire	4.4	Each		
M00583140	08714729858683	TRUEtome 44 Sphincterotome	5	30	.035 / 450 Jagwire Guidewire	4.4	Each		
M00583150	08714729858690	TRUEtome 44 Sphincterotome	5	30	.035 / Stiff 450 Jagwire Guidewire	4.4	Each		
TRUEtome Hydra Cannulating Sphincterotome (Preloaded with Hydra Jagwire™ High Performance Guidewire)									
M00583160	08714729858706	TRUEtome 44 Sphincterotome	5	20	.035 / 450 Hydra Jagwire Guidewire	4.4	Each		
M00583170	08714729858713	TRUEtome 44 Sphincterotome	5	30	.035 / 450 Hydra Jagwire Guidewire	4.4	Each		
TRUEtome Dreamwire™ Cannulating Sphincterotome (Preloaded with Dreamwire High Performance Guidewire)									
M00583180	08714729858720	TRUEtome 44 Sphincterotome	5	20	.035 / 260 Dreamwire Guidewire	4.4	Each		
M00583190	08714729858737	TRUEtome 44 Sphincterotome	5	20	.035 / 450 Dreamwire Guidewire	4.4	Each		
M00583200	08714729858744	TRUEtome 44 Sphincterotome	5	30	.035 / 260 Dreamwire Guidewire	4.4	Each		
M00583210	08714729858751	TRUEtome 44 Sphincterotome	5	30	.035 / 450 Dreamwire Guidewire	4.4	Each		

*Preloaded with straight-tip guidewires only.
Active Cords sold separately  (See page 12).



Section Table of Contents

Prev Page

7


Next Page

Main Table of Contents



Biliary

Sphincterotome (Continued)




TRUEtome™ Cannulating Sphincterotome



Order Number	GTIN	Description	Cut Wire Length (mm)	Tip Length (mm)	Distal Tip O.D. (Fr)	Recommended Guidewire (in)	Packaging
M00584190	08714729827658	TRUEtome 39 Sphincterotome	20	5	3.9	.025	Each
M00584120	08714729827771	TRUEtome 39 Sphincterotome	30	5	3.9	.025	Each
M00584170	08714729827634	TRUEtome 44 Sphincterotome	20	5	4.4	.035	Each
M00584180	08714729827641	TRUEtome 44 Sphincterotome	30	5	4.4	.035	Each
M00584150	08714729827610	TRUEtome 49 Sphincterotome	20	5	4.9	.035	Each
M00584160	08714729827627	TRUEtome 49 Sphincterotome	30	5	4.9	.035	Each

Active Cords sold separately [▶](#) (See page 12).



MicroKnife™ XL Triple-Lumen Needle Knife



Order Number	GTIN	Description	Shaft O.D. (Fr)	Catheter Length (cm)	Recommended Guidewire (in)	Packaging
M00532810	08714729172017	MicroKnife XL Triple-Lumen Needle Knife	7.0 - 5.5	200	.035	Each

Active Cords sold separately [▶](#) (See page 12).



Section Table of Contents

Prev Page [◀](#)

8


[▶](#) Next Page

Main Table of Contents



Biliary

Sphincterotome (Continued)




Stonetome™ Stone Removal Device



Order Number	GTIN	Description	Balloon Diameter (mm)	Catheter Length (cm)	Tip Length (mm)	Cut Wire (mm)	Distal Tip O.D. (Fr)	Recommended Guidewire (in)	Packaging
M00535150	08714729146667	Stonetome Device (below cut wire)	11.5	200	5	20	5.5	.035	Each
M00535170	08714729146629	Stonetome Device (below cut wire)	11.5	200	20	20	5.5	.035	Each
M00535190	08714729146612	Stonetome Device (below cut wire)	11.5	200	5	30	5.5	.035	Each
M00535210	08714729146605	Stonetome Device (below cut wire)	11.5	200	20	30	5.5	.035	Each
M00535110	08714729146599	Stonetome Device (above cut wire)	11.5	200	20	20	5.5	.035	Each
M00535130	08714729146575	Stonetome Device (above cut wire)	11.5	200	20	30	5.5	.035	Each

Active Cords sold separately (See page 12).



Ultratome™ XL Triple-Lumen Sphincterotome



Order Number	GTIN	Description	Wire Length (mm)	Shaft O.D. (Fr)	Catheter Length (cm)	Tip Length (mm)	Recommended Guidewire (in)	Packaging
M00535900	08714729103257	Ultratome XL Short Nose Sphincterotome	20	5.5	200	5	.035	Each
M00535910	08714729103301	Ultratome XL Long Nose Sphincterotome	20	5.5	200	20	.035	Each
M00535920	08714729103318	Ultratome XL Short Nose Sphincterotome	30	5.5	200	5	.035	Each
M00535930	08714729103349	Ultratome XL Long Nose Sphincterotome	30	5.5	200	20	.035	Each

Active Cords sold separately (See page 12).



Section Table of Contents

Prev Page

9

Next Page


Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Biliary

Sphincterotome (Continued)




Ultratome™ Double-Lumen Sphincterotome



Order Number	GTIN	Description	Wire Length (mm)	Distal Tip O.D.		Catheter Length (cm)	Tip Length (mm)	Recommended Guidewire (in)	Packaging
				(mm)	(Fr)				
M00530800	08714729163664	Ultratome Short Nose Sphincterotome	20	1.8	5.5	200	5	.035	Each
M00530810	08714729183143	Ultratome Long Nose Sphincterotome	20	1.8	5.5	200	20	.035	Each
M00530820	08714729193418	Ultratome Short Nose Sphincterotome	30	1.8	5.5	200	5	.035	Each
M00530830	08714729193425	Ultratome Long Nose Sphincterotome	30	1.8	5.5	200	20	.035	Each

Active Cords sold separately (See page 12).


RX Cannulas



Rapid Exchange XL Cannula



Order Number	GTIN	Description	Working Length (cm)	Tip Style	Recommended Guidewire (in)	Packaging
M00545820	08714729301905	Rapid Exchange XL Cannula	195	Tapered	.035	Each



RX ERCP Cannulas



Order Number	GTIN	Description	Working Length (cm)	Recommended Guidewire (in)	Packaging
M00545760	08714729268642	RX ERCP Cannula Standard Tip	210	.035	Each
M00545770	08714729268659	RX ERCP Cannula Tapered Tip	210	.035	Each
M00545780	08714729268673	RX ERCP Cannula Ball Tip	210	.035	Each



Section Table of Contents

Prev Page **10** Next Page


Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Biliary


Cannulas



Tandem™ XL Triple-Lumen ERCP Cannula



Order Number	GTIN	Distal Tip O.D.		Catheter Size (Fr)	Working Length (cm)	Tip Style	Recommended Guidewire (in)	Packaging
		(Fr)	(mm)					
M00535700	08714729146636	5.5	1.8	7.0 – 5.5	200	Tapered	.035	Each



Contour™ ERCP Cannulas



Order Number	GTIN	Description	Distal Tip O.D.		Working Length (cm)	Tip Style	Recommended Guidewire (in)	Packaging
			(Fr)	(mm)				
M00530850	08714729099789	Contour ERCP Cannula	5	1.7	210	Standard	.035	Each
M00530860	08714729099819	Contour ERCP Cannula	5	1.7	210	Tapered	.035	Each
M00530870	08714729099826	Contour ERCP Cannula	5	1.7	210	Ultra Tapered	.025	Each
M00530880	08714729099833	Contour ERCP Cannula	5	1.7	210	5-4-3	.018	Each
M00530890	08714729099840	Contour ERCP Cannula	5	1.7	210	Ball Tip	.035	Each



Section Table of Contents

Prev Page **11** Next Page

Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at **1.888.272.1001**

Biliary

Accessories


RX Locking Devices and Biopsy Cap Systems



Order Number	GTIN	Description	Packaging
M00545471	08714729337768	RX Locking Device (Biopsy Cap not included)	Box 10
M00545211	08714729440604	RX Locking Device Olympus™ / FujiFilm™ Brand Endoscope Compatible	Box 10
M00545261	08714729296614	RX Locking Device and Biopsy Cap System – Olympus / FujiFilm G5 Series Brand Endoscope Compatible	Box 10
M00545271	08714729296621	RX Locking Device and Biopsy Cap System – Pentax™ Brand Endoscope Compatible	Box 10
M00545281	08714729296638	RX Locking Device and Biopsy Cap System – FujiFilm Type 3 Series Brand Endoscope Compatible	Box 10

Active Cords (Fit Boston Scientific Endoscopy Sphincterotome)

Order Number	GTIN	Description	Packaging
M00561250	08714729019060	Active Cord, Fits Olympus Brand Generator	Each
M00561270	08714729053361	Active Cord, Fits Endostat™, Endostat II and Endostat III Generators, Bovie™, Valleylab™, Aspen Labs™, and ERBE™ Brand Generators	Each



Rapid Refill Continuous Injection System



Order Number	GTIN	Description	Packaging
M00566001	08714729285151	Rapid Refill Continuous Injection System	Box 10



Notes and Future Product Updates

Notes and Future Product Updates



Biliary Stricture and Stone Management

Biliary and Pancreatic Stents

Biliary Metal Stent Systems

➤ WallFlex™ Biliary RX Fully Covered Stent System RMV	15
➤ WallFlex Biliary RX Partially Covered Stent System	15
➤ WallFlex Biliary RX Uncovered Stent System	15
➤ Epic™ Endoscopic Biliary Stent System	16
➤ WALLSTENT™ RX Biliary Endoprosthesis Stent System	17
➤ WALLSTENT Endoscopic Biliary Endoprosthesis Stents	18

Biliary Plastic Stents

➤ Advanix™ Center Bend Preloaded RX Stent Systems	19
➤ Advanix Center Bend - Single Stent	20
➤ Advanix Double Pigtail - Single Stent	21
➤ Advanix Duodenal Bend Preloaded RX Stent Systems	19
➤ Advanix Duodenal Bend - Single Stent	20
➤ Biliary Introducer	23
➤ C-Flex™ Double Pigtail Biliary Stent	24
➤ Flexima™ Duodenal Bend Biliary Stent with Delivery System ..	22
➤ Flexima Single Pigtail Nasobiliary Catheter	24
➤ NaviFlex™ RX Delivery System	21
➤ Percuflex™ Duodenal Bend Biliary Stent	23
➤ Percuflex Duodenal Bend Biliary Stent with Introducer Kit	23
➤ RX Duodenal Bend Biliary Stent with RX Delivery System	22

Pancreatic Plastic Stents

➤ Advanix Pancreatic Stent Kit - Straight, No Leading Barb and Leading Barb	25
➤ Advanix Pancreatic Stent Kit - Pigtail, Leading Barb and No Leading Barb	25
➤ Advanix Pancreatic Stent - Straight, No Leading Barb	26
➤ Advanix Pancreatic Stent - Straight, Leading Barb	27
➤ Advanix Pancreatic Stent - Pigtail, No Leading Barb	28
➤ Advanix Pancreatic Stent - Pigtail, Leading Barb	29
➤ NaviFlex RX Pancreatic Delivery System & NaviFlex Pushers ..	30

Biliary Stone Management

RX Stone Management

➤ Extractor™ Pro RX Retrieval Balloon Catheters – Distal Injection	31
➤ Extractor Pro RX Retrieval Balloon Catheters – Proximal Injection	31
➤ Extractor Pro RX-S Retrieval Balloon Catheters – Distal Injection	31
➤ Extractor Pro RX-S Retrieval Balloon Catheters – Proximal Injection	31

Stone Management

➤ Extractor Pro DL Retrieval Balloon Catheters	32
➤ Extractor Pro XL Retrieval Balloon Catheters – Distal Injection	32
➤ Extractor Pro XL Retrieval Balloon Catheters – Proximal Injection	32
➤ Trapezoid™ RX Wireguided Retrieval Baskets	32

CRE™ PRO Wireguided Balloon Dilatation Catheters

➤ CRE PRO Wireguided Balloon Dilatation Catheters	33
---	----

Biliary Stricture Management

RX Dilatation

➤ Hurricane™ RX Biliary Balloon Dilatation Catheter	34
---	----

Biliary Dilatation

➤ MaxForce™ Biliary Balloon Dilatation Catheter	35
---	----

RX Accessories

➤ Encore™ 26 Inflation Device	35
-------------------------------------	----

RX Cytology

➤ RX Cytology Brush Wireguided Cytology Brush	35
---	----



Biliary Metal Stents

WallFlex™ Biliary RX Stents



WallFlex Biliary RX Fully Covered Stent System RMV*



Order Number	GTIN	Diameter (mm)	Length (mm)	Covered Length (mm) (PC Only)	Catheter Diameter (Fr / mm)	Guidewire Diameter (in / mm)	Working Length (cm)	Packaging
M00570340	08714729953159	8	60	—	8.5 / 2.83	.035 / .89	194	Each
M00570350	08714729953166	8	80	—	8.5 / 2.83	.035 / .89	194	Each
M00570360	08714729953173	10	40	—	8.5 / 2.83	.035 / .89	194	Each
M00570370	08714729953180	10	60	—	8.5 / 2.83	.035 / .89	194	Each
M00570380	08714729953197	10	80	—	8.5 / 2.83	.035 / .89	194	Each

*Compatible with Boston Scientific .035" guidewires, including Dreamwire™ High Performance Guidewire (See page 1).

Above (WallFlex Biliary RX Fully Covered Stent System) with Permalume™ Covering.

†Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



WallFlex Biliary RX Partially Covered Stent System*



Order Number	GTIN	Diameter (mm)	Length (mm)	Covered Length (mm) (PC Only)	Catheter Diameter (Fr / mm)	Guidewire Diameter (in / mm)	Working Length (cm)	Packaging
M00570700	08714729759294	8	60	48	8.5 / 2.83	.035 / .89	194	Each
M00570710	08714729759300	8	80	68	8.5 / 2.83	.035 / .89	194	Each
M00570720	08714729759317	10	40	28	8.5 / 2.83	.035 / .89	194	Each
M00570730	08714729759324	10	60	48	8.5 / 2.83	.035 / .89	194	Each
M00570740	08714729759331	10	80	68	8.5 / 2.83	.035 / .89	194	Each

*Compatible with Boston Scientific .035" guidewires, including Dreamwire High Performance Guidewire (See page 1).

Above (WallFlex Biliary RX Partially Covered Stent System) with Permalume Covering.

†Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



WallFlex Biliary RX Uncovered Stent System*



Order Number	GTIN	Diameter (mm)	Length (mm)	Covered Length (mm) (PC Only)	Catheter Diameter (Fr / mm)	Guidewire Diameter (in / mm)	Working Length (cm)	Packaging
M00570600	08714729759249	8	40	—	8.0 / 2.67	.035 / .89	194	Each
M00570610	08714729759256	8	60	—	8.0 / 2.67	.035 / .89	194	Each
M00570620	08714729759263	8	80	—	8.0 / 2.67	.035 / .89	194	Each
M00570630	08714729759270	8	100	—	8.0 / 2.67	.035 / .89	194	Each
M00570890	08714729765332	10	40	—	8.0 / 2.67	.035 / .89	194	Each
M00570640	08714729756231	10	60	—	8.0 / 2.67	.035 / .89	194	Each
M00570650	08714729758150	10	80	—	8.0 / 2.67	.035 / .89	194	Each
M00570660	08714729759287	10	100	—	8.0 / 2.67	.035 / .89	194	Each

*Compatible with Boston Scientific .035" guidewires, including Dreamwire High Performance Guidewire (See page 1).

†Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



Biliary Metal Stents

Epic™ Endoscopic Biliary Stent System



 Magnetic Resonance Conditional

Epic Endoscopic Biliary Stent

Order Number	GTIN	Description	Stent Diameter (mm)	Stent Length (mm)	Catheter Diameter (Fr)	Guidewire Diameter (in)	Catheter Working Length (cm)	Packaging
M00572230	08714729896852	EPIC Biliary Stent	6	40	6	0.035	220	Each
M00572240	08714729896869	EPIC Biliary Stent	6	60	6	0.035	220	Each
M00572250	08714729896876	EPIC Biliary Stent	6	80	6	0.035	220	Each
M00572260	08714729896883	EPIC Biliary Stent	6	100	6	0.035	220	Each
M00572270	08714729896890	EPIC Biliary Stent	8	40	6	0.035	220	Each
M00572280	08714729896906	EPIC Biliary Stent	8	60	6	0.035	220	Each
M00572290	08714729896913	EPIC Biliary Stent	8	80	6	0.035	220	Each
M00572300	08714729896920	EPIC Biliary Stent	8	100	6	0.035	220	Each
M00572310	08714729896937	EPIC Biliary Stent	10	40	6	0.035	220	Each
M00572320	08714729896944	EPIC Biliary Stent	10	60	6	0.035	220	Each
M00572330	08714729896951	EPIC Biliary Stent	10	80	6	0.035	220	Each
M00572340	08714729896968	EPIC Biliary Stent	10	100	6	0.035	220	Each



Section Table of Contents

Prev Page

16

Next Page

Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Biliary Metal Stents

WALLSTENT™ Biliary RX Stents



 Magnetic Resonance Conditional¹

WALLSTENT RX Biliary Endoprosthesis Stent System*



Order Number	GTIN	Description	Unconstrained Diameter (mm)	Nominal Unconstrained Length (mm)	Delivery System (Fr)	Cover Length (mm)	Working Length (cm)	Packaging
M00569670	08714729296874	WALLSTENT RX Biliary Partially Covered Stent System	8	40	8	30	194	Each
M00569680	08714729296881	WALLSTENT RX Biliary Partially Covered Stent System	8	60	8	50	194	Each
M00569690	08714729296898	WALLSTENT RX Biliary Partially Covered Stent System	8	80	8	70	194	Each
M00569700	08714729296904	WALLSTENT RX Biliary Partially Covered Stent System	10	40	8	30	194	Each
M00569710	08714729296911	WALLSTENT RX Biliary Partially Covered Stent System	10	60	8	50	194	Each
M00569720	08714729296928	WALLSTENT RX Biliary Partially Covered Stent System	10	80	8	70	194	Each

*Recommended Guidewire: .035", 260cm, including Dreamwire™ High Performance Guidewire  (See page 1).

Above (WALLSTENT RX Biliary Endoprosthesis Stent System) with Permalume™ Covering.

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



 Magnetic Resonance Conditional¹

WALLSTENT RX Biliary Endoprosthesis Stent System*



Order Number	GTIN	Description	Unconstrained Diameter (mm)	Nominal Unconstrained Length (mm)	Delivery System (Fr)	Working Length (cm)	Packaging
M00569610	08714729296812	WALLSTENT RX Biliary Uncovered Stent System	8	40	8	194	Each
M00569620	08714729296829	WALLSTENT RX Biliary Uncovered Stent System	8	60	8	194	Each
M00569630	08714729296836	WALLSTENT RX Biliary Uncovered Stent System	8	80	8	194	Each
M00569800	08714729456711	WALLSTENT RX Biliary Uncovered Stent System	8	100	8	194	Each
M00569640	08714729296843	WALLSTENT RX Biliary Uncovered Stent System	10	40	8	194	Each
M00569650	08714729296850	WALLSTENT RX Biliary Uncovered Stent System	10	60	8	194	Each
M00569660	08714729296867	WALLSTENT RX Biliary Uncovered Stent System	10	80	8	194	Each
M00569810	08714729456728	WALLSTENT RX Biliary Uncovered Stent System	10	100	8	194	Each

*Recommended Guidewire: .035", 260cm, including Dreamwire High Performance Guidewire  (See page 1).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



Section Table of Contents

Prev Page  17  Next Page


Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Biliary Metal Stents

WALLSTENT™ Stents



 Magnetic Resonance Conditional¹

WALLSTENT Endoscopic Biliary Endoprosthesis Stents*




Order Number	GTIN	Description	Unconstrained Diameter (mm)	Nominal Unconstrained Length (mm)	Delivery System (Fr)	Cover Length (mm)	Working Length (cm)	Packaging
H965430700	08714729229858	WALLSTENT Biliary Partially Covered Stent	8	40	8	30	194	Each
H965430800	08714729317906	WALLSTENT Biliary Partially Covered Stent	8	60	8	50	194	Each
H965430900	08714729229872	WALLSTENT Biliary Partially Covered Stent	8	80	8	70	194	Each
H965431000	08714729317920	WALLSTENT Biliary Partially Covered Stent	10	40	8	30	194	Each
H965431100	08714729317937	WALLSTENT Biliary Partially Covered Stent	10	60	8	50	194	Each
H965431200	08714729229902	WALLSTENT Biliary Partially Covered Stent	10	80	8	70	194	Each

*Recommended Guidewire: .035", 450cm, including Dreamwire™ High Performance Guidewire  (See page 1).

Above (WALLSTENT Biliary Endoprosthesis) with Permalume™ Covering.

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



 Magnetic Resonance Conditional¹

WALLSTENT Endoscopic Biliary Endoprosthesis Stents*



Order Number	GTIN	Description	Unconstrained Diameter (mm)	Nominal Unconstrained Length (mm)	Delivery System Diameter (Fr)	Working Length (cm)	Packaging
H965430100	08714729204633	WALLSTENT Biliary Uncovered Stent	8	40	7.5	194	Each
H965430200	08714729204640	WALLSTENT Biliary Uncovered Stent	8	60	7.5	194	Each
H965430300	08714729204657	WALLSTENT Biliary Uncovered Stent	8	80	7.5	194	Each
H965430400	08714729204664	WALLSTENT Biliary Uncovered Stent	10	40	7.5	194	Each
H965430500	08714729204671	WALLSTENT Biliary Uncovered Stent	10	60	7.5	194	Each
H965430600	08714729204688	WALLSTENT Biliary Uncovered Stent	10	80	7.5	194	Each

*Recommended Guidewire: .035", 450cm, including Dreamwire High Performance Guidewire  (See page 1).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



Section Table of Contents

Prev Page

18

Next Page

Main Table of Contents



Biliary Plastic Stents

Advanix™ Biliary Stents



MR Magnetic Resonance Safe

Advanix Duodenal Bend Preloaded RX Stent Systems



Order Number	GTIN	Description	Stent Diameter (Fr)	Length (cm) Barb to Barb	Packaging
M00534200	08714729787297	Duodenal Bend	7	5	Each
M00534210	08714729787303	Duodenal Bend	7	7	Each
M00534220	08714729787310	Duodenal Bend	7	9	Each
M00534230	08714729787327	Duodenal Bend	7	12	Each
M00534240	08714729787334	Duodenal Bend	7	15	Each
M00534250	08714729787341	Duodenal Bend	7	18	Each
M00534260	08714729787358	Duodenal Bend	8.5	5	Each
M00534270	08714729787365	Duodenal Bend	8.5	7	Each
M00534280	08714729787372	Duodenal Bend	8.5	9	Each
M00534290	08714729787389	Duodenal Bend	8.5	12	Each
M00534300	08714729787396	Duodenal Bend	8.5	15	Each
M00534310	08714729787402	Duodenal Bend	8.5	18	Each
M00534320	08714729787419	Duodenal Bend	10	5	Each
M00534330	08714729787426	Duodenal Bend	10	7	Each
M00534340	08714729787433	Duodenal Bend	10	9	Each
M00534350	08714729787440	Duodenal Bend	10	12	Each
M00534360	08714729787457	Duodenal Bend	10	15	Each
M00534370	08714729787464	Duodenal Bend	10	18	Each

Scope Compatibility: 7Fr and 8.5Fr stents require a 3.2mm working channel; 10Fr stents require a 4.2mm working channel.



MR Magnetic Resonance Safe

Advanix Center Bend Preloaded RX Stent Systems



Order Number	GTIN	Description	Stent Diameter (Fr)	Length (cm) Barb to Barb	Packaging
M00534560	08714729787112	Center Bend	7	5	Each
M00534570	08714729787129	Center Bend	7	7	Each
M00534580	08714729787136	Center Bend	7	9	Each
M00534590	08714729787143	Center Bend	7	12	Each
M00534600	08714729787150	Center Bend	7	15	Each
M00534610	08714729787167	Center Bend	7	18	Each
M00534620	08714729787174	Center Bend	8.5	5	Each
M00534630	08714729787181	Center Bend	8.5	7	Each
M00534640	08714729787198	Center Bend	8.5	9	Each
M00534650	08714729787204	Center Bend	8.5	12	Each
M00534660	08714729787211	Center Bend	8.5	15	Each
M00534670	08714729787228	Center Bend	8.5	18	Each
M00534680	08714729787235	Center Bend	10	5	Each
M00534690	08714729787242	Center Bend	10	7	Each
M00534700	08714729787259	Center Bend	10	9	Each
M00534710	08714729787266	Center Bend	10	12	Each
M00534720	08714729787273	Center Bend	10	15	Each
M00534730	08714729787280	Center Bend	10	18	Each

Scope Compatibility: 7Fr and 8.5Fr stents require a 3.2mm working channel; 10Fr stents require a 4.2mm working channel.



Section Table of Contents

Prev Page

19

Next Page

Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Biliary Plastic Stents

Advanix™ Biliary Stents (Continued)



MR Magnetic Resonance Safe



Advanix Duodenal Bend - Single Stent

Order Number	GTIN	Description	Stent Diameter (Fr)	Length (cm) Barb to Barb	Packaging
M00532870	08714729786788	Duodenal Bend	7	5	Each
M00532880	08714729786795	Duodenal Bend	7	7	Each
M00532890	08714729786801	Duodenal Bend	7	9	Each
M00532900	08714729786818	Duodenal Bend	7	12	Each
M00532910	08714729786825	Duodenal Bend	7	15	Each
M00532920	08714729786832	Duodenal Bend	7	18	Each
M00532930	08714729786849	Duodenal Bend	8.5	5	Each
M00532940	08714729786856	Duodenal Bend	8.5	7	Each
M00532950	08714729786863	Duodenal Bend	8.5	9	Each
M00532960	08714729786870	Duodenal Bend	8.5	12	Each
M00532970	08714729786887	Duodenal Bend	8.5	15	Each
M00532980	08714729786894	Duodenal Bend	8.5	18	Each
M00532990	08714729786900	Duodenal Bend	10	5	Each
M00533000	08714729786917	Duodenal Bend	10	7	Each
M00533010	08714729786924	Duodenal Bend	10	9	Each
M00533020	08714729786931	Duodenal Bend	10	12	Each
M00533030	08714729786948	Duodenal Bend	10	15	Each
M00533040	08714729786955	Duodenal Bend	10	18	Each

Scope Compatibility: 7Fr and 8.5Fr stents require a 3.2mm working channel; 10Fr stents require a 4.2mm working channel.



MR Magnetic Resonance Safe



Advanix Center Bend - Single Stent

Order Number	GTIN	Description	Stent Diameter (Fr)	Length (cm) Barb to Barb	Packaging
M00533230	08714729786603	Center Bend	7	5	Each
M00533240	08714729786610	Center Bend	7	7	Each
M00533250	08714729786627	Center Bend	7	9	Each
M00533260	08714729786634	Center Bend	7	12	Each
M00533270	08714729786641	Center Bend	7	15	Each
M00533280	08714729786658	Center Bend	7	18	Each
M00533290	08714729786665	Center Bend	8.5	5	Each
M00533300	08714729786672	Center Bend	8.5	7	Each
M00533310	08714729786689	Center Bend	8.5	9	Each
M00533320	08714729786696	Center Bend	8.5	12	Each
M00533330	08714729786702	Center Bend	8.5	15	Each
M00533340	08714729786719	Center Bend	8.5	18	Each
M00533350	08714729786726	Center Bend	10	5	Each
M00533360	08714729786733	Center Bend	10	7	Each
M00533370	08714729786740	Center Bend	10	9	Each
M00533380	08714729786757	Center Bend	10	12	Each
M00533390	08714729786764	Center Bend	10	15	Each
M00533400	08714729786771	Center Bend	10	18	Each

Scope Compatibility: 7Fr and 8.5Fr stents require a 3.2mm working channel; 10Fr stents require a 4.2mm working channel.



Section Table of Contents

Prev Page

20

Next Page

Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Biliary Plastic Stents

Advanix™ Biliary Stents (Continued)



MR Magnetic Resonance Safe

Advanix Double Pigtail - Single Stent



Order Number	GTIN	Description	Stent Diameter (Fr)	Length (cm) Between Pigtails	Packaging
M00532160	08714729786962	Double Pigtail	7	3	Each
M00532180	08714729786979	Double Pigtail	7	5	Each
M00532190	08714729786986	Double Pigtail	7	7	Each
M00532210	08714729786993	Double Pigtail	7	10	Each
M00532220	08714729787006	Double Pigtail	7	12	Each
M00532230	08714729787013	Double Pigtail	7	15	Each
M00532240	08714729787020	Double Pigtail	10	3	Each
M00532260	08714729787037	Double Pigtail	10	5	Each
M00532270	08714729787044	Double Pigtail	10	7	Each
M00532290	08714729787051	Double Pigtail	10	10	Each
M00532300	08714729787068	Double Pigtail	10	12	Each
M00532310	08714729787075	Double Pigtail	10	15	Each

Scope Compatibility: 7Fr and 8.5Fr stents require a 3.2mm working channel; 10Fr stents require a 4.2mm working channel.

Naviflex™ RX Delivery System



NaviFlex RX Delivery System*



Order Number	GTIN	Description	Catheter Diameter (Fr)	Length (cm)	Packaging
M00533560	08714729787082	RX Delivery System	7	202.5	Each
M00533570	08714729787099	RX Delivery System	8.5	202.5	Each
M00533580	08714729787105	RX Delivery System	10	202.5	Each

*Compatible with Advanix Biliary Stents.



Section Table of Contents

Prev Page ◀ 21 ▶ Next Page

Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Biliary Plastic Stents

Duodenal Bend Stents

MR Magnetic Resonance Safe



RX Duodenal Bend Biliary Stent with RX Delivery System

Order Number	GTIN	Description	O.D. (Fr)	Length (cm) Barb to Barb	Recommended Guidewire (in)	Packaging
M00545550	08714729337850	Duodenal Bend Biliary Stent with RX Delivery System	7	5	.035	Each
M00545560	08714729337867	Duodenal Bend Biliary Stent with RX Delivery System	7	7	.035	Each
M00545570	08714729337874	Duodenal Bend Biliary Stent with RX Delivery System	7	10	.035	Each
M00545580	08714729337881	Duodenal Bend Biliary Stent with RX Delivery System	7	12	.035	Each
M00545590	08714729337898	Duodenal Bend Biliary Stent with RX Delivery System	7	15	.035	Each
M00545650	08714729252221	Duodenal Bend Biliary Stent with RX Delivery System	8.5	5	.035	Each
M00545660	08714729301837	Duodenal Bend Biliary Stent with RX Delivery System	8.5	7	.035	Each
M00545670	08714729301844	Duodenal Bend Biliary Stent with RX Delivery System	8.5	10	.035	Each
M00545680	08714729301851	Duodenal Bend Biliary Stent with RX Delivery System	8.5	12	.035	Each
M00545690	08714729301868	Duodenal Bend Biliary Stent with RX Delivery System	8.5	15	.035	Each
M00545600	08714729337904	Duodenal Bend Biliary Stent with RX Delivery System	10	5	.035	Each
M00545610	08714729337911	Duodenal Bend Biliary Stent with RX Delivery System	10	7	.035	Each
M00545620	08714729337928	Duodenal Bend Biliary Stent with RX Delivery System	10	10	.035	Each
M00545630	08714729301813	Duodenal Bend Biliary Stent with RX Delivery System	10	12	.035	Each
M00545640	08714729301820	Duodenal Bend Biliary Stent with RX Delivery System	10	15	.035	Each
M00545700	08714729301875	Duodenal Bend Biliary Stent with RX Delivery System	11.5	5	.035	Each
M00545710	08714729301882	Duodenal Bend Biliary Stent with RX Delivery System	11.5	7	.035	Each
M00545720	08714729301899	Duodenal Bend Biliary Stent with RX Delivery System	11.5	10	.035	Each
M00545730	08714729337942	Duodenal Bend Biliary Stent with RX Delivery System	11.5	12	.035	Each
M00545740	08714729337959	Duodenal Bend Biliary Stent with RX Delivery System	11.5	15	.035	Each



MR Magnetic Resonance Safe



Flexima™ Duodenal Bend Biliary Stent with Delivery System

Order Number	GTIN	O.D. (Fr)	Length (cm) Barb to Barb	Recommended Guidewire (in)	Packaging
M00539200	08714729337362	7	5	.035	Each
M00539210	08714729162520	7	7	.035	Each
M00539220	08714729162537	7	10	.035	Each
M00539230	08714729162544	7	12	.035	Each
M00539240	08714729162551	7	15	.035	Each
M00539250	08714729162568	8.5	5	.035	Each
M00539260	08714729162575	8.5	7	.035	Each
M00539270	08714729162582	8.5	10	.035	Each
M00539280	08714729162599	8.5	12	.035	Each
M00539290	08714729162605	8.5	15	.035	Each
M00539300	08714729162612	10	5	.035	Each
M00539310	08714729162629	10	7	.035	Each
M00539320	08714729162636	10	10	.035	Each
M00539330	08714729162643	10	12	.035	Each
M00539340	08714729162650	10	15	.035	Each
M00539350	08714729162667	11.5	5	.035	Each
M00539360	08714729162674	11.5	7	.035	Each
M00539370	08714729162681	11.5	10	.035	Each
M00539380	08714729162698	11.5	12	.035	Each
M00539390	08714729162704	11.5	15	.035	Each



Section Table of Contents

Prev Page ◀ 22 ▶ Next Page

Main Table of Contents



Biliary Plastic Stents

Duodenal Bend Stents *(Continued)*



MR Magnetic Resonance Safe

Percuflex™ Duodenal Bend Biliary Stent



Order Number	GTIN	O.D. (Fr)	Length (cm) Barb to Barb	Recommended Guidewire (in)	Packaging
M00533610	08714729181217	7	5	.035	Each
M00533620	08714729181224	7	7	.035	Each
M00533630	08714729181231	7	10	.035	Each
M00533640	08714729181248	7	12	.035	Each
M00533650	08714729181255	7	15	.035	Each
M00533660	08714729181262	10	5	.035	Each
M00533670	08714729181279	10	7	.035	Each
M00533680	08714729181286	10	10	.035	Each
M00533690	08714729181293	10	12	.035	Each
M00533700	08714729181309	10	15	.035	Each

MR Magnetic Resonance Safe

Percuflex Duodenal Bend Biliary Stent with Introducer Kit



Order Number	GTIN	Description	O.D. (Fr)	Length (cm) Barb to Barb	Packaging
M00533760	08714729181088	Percuflex Stent w/ push catheter	7	5	Each
M00533770	08714729181095	Percuflex Stent w/ push catheter	7	7	Each
M00533780	08714729181101	Percuflex Stent w/ push catheter	7	10	Each
M00533790	08714729181118	Percuflex Stent w/ push catheter	7	12	Each
M00533800	08714729181125	Percuflex Stent w/ push catheter	7	15	Each
M00533810	08714729181132	Percuflex Stent w/ push and guide catheter	10	5	Each
M00533820	08714729181149	Percuflex Stent w/ push and guide catheter	10	7	Each
M00533830	08714729181156	Percuflex Stent w/ push and guide catheter	10	10	Each
M00533840	08714729193449	Percuflex Stent w/ push and guide catheter	10	12	Each
M00533850	08714729181163	Percuflex Stent w/ push and guide catheter	10	15	Each

Biliary Introducer*



Order Number	GTIN	Description	O.D. (Fr)	Recommended Guidewire (in)	Packaging
M00533910	08714729181439	Push catheter	7	.035	Each
M00533920	08714729193463	Push and guide catheter	10	.035	Each

*Stents sold separately.



Section Table of Contents

Prev Page **23** Next Page

Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Biliary Plastic Stents

Pigtail Stents




MR Magnetic Resonance Safe



C-Flex™ Double Pigtail Biliary Stent

Order Number	GTIN	O.D. (Fr)	Length (cm) Between Pigtails	Color Code	Packaging
M00532010	08714729019077	7	5	Blue	Each
M00532020	08714729019084	7	10	Blue	Each
M00532030	08714729019091	7	15	Blue	Each
M00532040	08714729019107	10	5	Grey	Each
M00532050	08714729019114	10	10	Grey	Each



MR Magnetic Resonance Safe



Flexima™ Single Pigtail Nasobiliary Catheter*

Order Number	GTIN	O.D. (Fr)	Working Length (cm)	Packaging
M00540110	08714729337379	5	250	Each
M00540120	08714729337386	6	250	Each
M00540130	08714729337393	7.5	250	Each
M00540140	08714729337409	8.5	250	Each

*Recommended Guidewire .035" Dreamwire™ High Performance Guidewire, Order #5614, [See page 1](#).



Section Table of Contents

Prev Page

24

Next Page

Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Pancreatic Plastic Stents

Advanix™ Pancreatic Stents



 Magnetic Resonance Conditional¹

Advanix Pancreatic Stent Kit - Straight, No Leading Barb and Leading Barb



Order Number	GTIN	Description	Diameter (Fr)	Length (cm)	Packaging
M00537460	08714729856245	Straight Stent – No Leading Barb	4	3	Each
M00537470	08714729856252	Straight Stent – No Leading Barb	4	4	Each
M00537480	08714729856269	Straight Stent – No Leading Barb	4	5	Each
M00537490	08714729862956	Straight Stent – No Leading Barb	4	7	Each
M00537500	08714729856276	Straight Stent – No Leading Barb	5	3	Each
M00537510	08714729856283	Straight Stent – No Leading Barb	5	4	Each
M00537520	08714729856290	Straight Stent – No Leading Barb	5	5	Each
M00537530	08714729856306	Straight Stent – No Leading Barb	5	7	Each
M00537540	08714729856313	Straight Stent – Leading Barb	5	3	Each
M00537550	08714729856320	Straight Stent – Leading Barb	5	4	Each
M00537560	08714729856337	Straight Stent – Leading Barb	5	5	Each

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



 Magnetic Resonance Conditional¹

Advanix Pancreatic Stent Kit - Pigtail, Leading Barb and No Leading Barb



Order Number	GTIN	Description	Diameter (Fr)	Length (cm)	Packaging
M00537350	08714729856344	Single Pigtail Stent – No Leading Barb	4	3	Each
M00537360	08714729856351	Single Pigtail Stent – No Leading Barb	4	4	Each
M00537370	08714729856368	Single Pigtail Stent – No Leading Barb	4	5	Each
M00537380	08714729856375	Single Pigtail Stent – No Leading Barb	4	7	Each
M00537390	08714729856382	Single Pigtail Stent – No Leading Barb	5	3	Each
M00537400	08714729856399	Single Pigtail Stent – No Leading Barb	5	4	Each
M00537410	08714729856405	Single Pigtail Stent – No Leading Barb	5	5	Each
M00537420	08714729856412	Single Pigtail Stent – No Leading Barb	5	7	Each
M00537430	08714729856429	Single Pigtail Stent – Leading Barb	5	3	Each
M00537440	08714729856436	Single Pigtail Stent – Leading Barb	5	4	Each
M00537450	08714729856443	Single Pigtail Stent – Leading Barb	5	5	Each

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



Section Table of Contents

Prev Page  25  Next Page

Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Pancreatic Plastic Stents

Advanix™ Pancreatic Stents (Continued)



 Magnetic Resonance Conditional¹

Advanix Pancreatic Stent - Straight, No Leading Barb 					
Order Number	GTIN	Description	Diameter (Fr)	Length (cm)	Packaging
M00536000	08714729857686	Straight – No Leading Barb	3	2	Each
M00536010	08714729857693	Straight – No Leading Barb	3	3	Each
M00536020	08714729857709	Straight – No Leading Barb	3	4	Each
M00536030	08714729857716	Straight – No Leading Barb	3	5	Each
M00536040	08714729857723	Straight – No Leading Barb	3	6	Each
M00536050	08714729857730	Straight – No Leading Barb	3	7	Each
M00536060	08714729857747	Straight – No Leading Barb	3	8	Each
M00536070	08714729857754	Straight – No Leading Barb	3	9	Each
M00536080	08714729857761	Straight – No Leading Barb	3	10	Each
M00536090	08714729857778	Straight – No Leading Barb	3	11	Each
M00536100	08714729857785	Straight – No Leading Barb	3	12	Each
M00536110	08714729857792	Straight – No Leading Barb	3	13	Each
M00536120	08714729857068	Straight – No Leading Barb	4	2	Each
M00536130	08714729857075	Straight – No Leading Barb	4	3	Each
M00536140	08714729857082	Straight – No Leading Barb	4	4	Each
M00536150	08714729857099	Straight – No Leading Barb	4	5	Each
M00536160	08714729857105	Straight – No Leading Barb	4	6	Each
M00536170	08714729857112	Straight – No Leading Barb	4	7	Each
M00536180	08714729857129	Straight – No Leading Barb	4	8	Each
M00536190	08714729857136	Straight – No Leading Barb	4	9	Each
M00536200	08714729857143	Straight – No Leading Barb	4	10	Each
M00536210	08714729857150	Straight – No Leading Barb	4	11	Each
M00536220	08714729857167	Straight – No Leading Barb	4	12	Each
M00536230	08714729857174	Straight – No Leading Barb	4	13	Each
M00536240	08714729857181	Straight – No Leading Barb	5	2	Each
M00536250	08714729857198	Straight – No Leading Barb	5	3	Each
M00536260	08714729857204	Straight – No Leading Barb	5	4	Each
M00536270	08714729857211	Straight – No Leading Barb	5	5	Each
M00536280	08714729857228	Straight – No Leading Barb	5	6	Each
M00536290	08714729857235	Straight – No Leading Barb	5	7	Each
M00536300	08714729857242	Straight – No Leading Barb	5	8	Each
M00536310	08714729857259	Straight – No Leading Barb	5	9	Each
M00536320	08714729857266	Straight – No Leading Barb	5	10	Each
M00536330	08714729857273	Straight – No Leading Barb	5	11	Each
M00536340	08714729857280	Straight – No Leading Barb	5	12	Each
M00536350	08714729857297	Straight – No Leading Barb	5	13	Each

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

Pancreatic Plastic Stents

Advanix™ Pancreatic Stents (Continued)



 Magnetic Resonance
Conditional¹

Advanix Pancreatic Stent - Straight, Leading Barb



Order Number	GTIN	Description	Diameter (Fr)	Length (cm)	Packaging
M00536360	08714729857303	Straight – Leading Barb	5	2	Each
M00536370	08714729857310	Straight – Leading Barb	5	3	Each
M00536380	08714729857327	Straight – Leading Barb	5	4	Each
M00536390	08714729857334	Straight – Leading Barb	5	5	Each
M00536400	08714729857341	Straight – Leading Barb	5	6	Each
M00536410	08714729857358	Straight – Leading Barb	5	7	Each
M00536420	08714729857365	Straight – Leading Barb	5	8	Each
M00536430	08714729857372	Straight – Leading Barb	5	9	Each
M00536440	08714729857389	Straight – Leading Barb	5	10	Each
M00536450	08714729857396	Straight – Leading Barb	5	11	Each
M00536460	08714729857402	Straight – Leading Barb	5	12	Each
M00536470	08714729857419	Straight – Leading Barb	5	13	Each
M00536480	08714729857426	Straight – Leading Barb	7	3	Each
M00536490	08714729857433	Straight – Leading Barb	7	4	Each
M00536500	08714729857440	Straight – Leading Barb	7	5	Each
M00536510	08714729857457	Straight – Leading Barb	7	6	Each
M00536520	08714729857464	Straight – Leading Barb	7	7	Each
M00536530	08714729857471	Straight – Leading Barb	7	8	Each
M00536540	08714729857488	Straight – Leading Barb	7	9	Each
M00536550	08714729857495	Straight – Leading Barb	7	10	Each
M00536560	08714729857501	Straight – Leading Barb	7	11	Each
M00536570	08714729857518	Straight – Leading Barb	7	12	Each
M00536580	08714729857525	Straight – Leading Barb	7	13	Each
M00536590	08714729857532	Straight – Leading Barb	7	15	Each
M00536600	08714729857549	Straight – Leading Barb	7	18	Each
M00536610	08714729857556	Straight – Leading Barb	10	3	Each
M00536620	08714729857563	Straight – Leading Barb	10	4	Each
M00536630	08714729857570	Straight – Leading Barb	10	5	Each
M00536640	08714729857587	Straight – Leading Barb	10	6	Each
M00536650	08714729857594	Straight – Leading Barb	10	7	Each
M00536660	08714729857600	Straight – Leading Barb	10	8	Each
M00536670	08714729857617	Straight – Leading Barb	10	9	Each
M00536680	08714729857624	Straight – Leading Barb	10	10	Each
M00536690	08714729857631	Straight – Leading Barb	10	11	Each
M00536700	08714729857648	Straight – Leading Barb	10	12	Each
M00536710	08714729857655	Straight – Leading Barb	10	13	Each
M00536720	08714729857662	Straight – Leading Barb	10	15	Each
M00536730	08714729857679	Straight – Leading Barb	10	18	Each

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



Section Table of Contents

Prev Page ◀ 27 ▶ Next Page

Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Pancreatic Plastic Stents

Advanix™ Pancreatic Stents (Continued)



 Magnetic Resonance Conditional¹

Advanix Pancreatic Stent - Pigtail, No Leading Barb

Order Number	GTIN	Description	Diameter (Fr)	Length (cm)	Packaging
M00536740	08714729856948	Single Pigtail – No Leading Barb	3	2	Each
M00536750	08714729856955	Single Pigtail – No Leading Barb	3	3	Each
M00536760	08714729856962	Single Pigtail – No Leading Barb	3	4	Each
M00536770	08714729856979	Single Pigtail – No Leading Barb	3	5	Each
M00536780	08714729856986	Single Pigtail – No Leading Barb	3	6	Each
M00536790	08714729856993	Single Pigtail – No Leading Barb	3	7	Each
M00536800	08714729857006	Single Pigtail – No Leading Barb	3	8	Each
M00536810	08714729857013	Single Pigtail – No Leading Barb	3	9	Each
M00536820	08714729857020	Single Pigtail – No Leading Barb	3	10	Each
M00536830	08714729857037	Single Pigtail – No Leading Barb	3	11	Each
M00536840	08714729857044	Single Pigtail – No Leading Barb	3	12	Each
M00536850	08714729857051	Single Pigtail – No Leading Barb	3	13	Each
M00536860	08714729856450	Single Pigtail – No Leading Barb	4	2	Each
M00536870	08714729856467	Single Pigtail – No Leading Barb	4	3	Each
M00536880	08714729856474	Single Pigtail – No Leading Barb	4	4	Each
M00536890	08714729856481	Single Pigtail – No Leading Barb	4	5	Each
M00536900	08714729856498	Single Pigtail – No Leading Barb	4	6	Each
M00536910	08714729856504	Single Pigtail – No Leading Barb	4	7	Each
M00536920	08714729856511	Single Pigtail – No Leading Barb	4	8	Each
M00536930	08714729856528	Single Pigtail – No Leading Barb	4	9	Each
M00536940	08714729856535	Single Pigtail – No Leading Barb	4	10	Each
M00536950	08714729856542	Single Pigtail – No Leading Barb	4	11	Each
M00536960	08714729856559	Single Pigtail – No Leading Barb	4	12	Each
M00536970	08714729856566	Single Pigtail – No Leading Barb	4	13	Each
M00536980	08714729856696	Single Pigtail – No Leading Barb	5	2	Each
M00536990	08714729856702	Single Pigtail – No Leading Barb	5	3	Each
M00537000	08714729856719	Single Pigtail – No Leading Barb	5	4	Each
M00537010	08714729856726	Single Pigtail – No Leading Barb	5	5	Each
M00537020	08714729856733	Single Pigtail – No Leading Barb	5	6	Each
M00537030	08714729856740	Single Pigtail – No Leading Barb	5	7	Each
M00537040	08714729856757	Single Pigtail – No Leading Barb	5	8	Each
M00537050	08714729856764	Single Pigtail – No Leading Barb	5	9	Each
M00537060	08714729856771	Single Pigtail – No Leading Barb	5	10	Each
M00537070	08714729856788	Single Pigtail – No Leading Barb	5	11	Each
M00537080	08714729856795	Single Pigtail – No Leading Barb	5	12	Each
M00537090	08714729856801	Single Pigtail – No Leading Barb	5	13	Each

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



Section Table of Contents

Prev Page ◀ 28 ▶ Next Page

Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Pancreatic Plastic Stents

Advanix™ Pancreatic Stents (Continued)



 Magnetic Resonance Conditional¹

Advanix Pancreatic Stent - Pigtail, Leading Barb

Order Number	GTIN	Description	Diameter (Fr)	Length (cm)	Packaging
M00537100	08714729856573	Single Pigtail – Leading Barb	5	2	Each
M00537110	08714729856580	Single Pigtail – Leading Barb	5	3	Each
M00537120	08714729856597	Single Pigtail – Leading Barb	5	4	Each
M00537130	08714729856603	Single Pigtail – Leading Barb	5	5	Each
M00537140	08714729856610	Single Pigtail – Leading Barb	5	6	Each
M00537150	08714729856627	Single Pigtail – Leading Barb	5	7	Each
M00537160	08714729856634	Single Pigtail – Leading Barb	5	8	Each
M00537170	08714729856641	Single Pigtail – Leading Barb	5	9	Each
M00537180	08714729856658	Single Pigtail – Leading Barb	5	10	Each
M00537190	08714729856665	Single Pigtail – Leading Barb	5	11	Each
M00537200	08714729856672	Single Pigtail – Leading Barb	5	12	Each
M00537210	08714729856689	Single Pigtail – Leading Barb	5	13	Each
M00537220	08714729856818	Single Pigtail – Leading Barb	7	3	Each
M00537230	08714729856825	Single Pigtail – Leading Barb	7	4	Each
M00537240	08714729856832	Single Pigtail – Leading Barb	7	5	Each
M00537250	08714729856849	Single Pigtail – Leading Barb	7	6	Each
M00537260	08714729856856	Single Pigtail – Leading Barb	7	7	Each
M00537270	08714729856863	Single Pigtail – Leading Barb	7	8	Each
M00537280	08714729856870	Single Pigtail – Leading Barb	7	9	Each
M00537290	08714729856887	Single Pigtail – Leading Barb	7	10	Each
M00537300	08714729856894	Single Pigtail – Leading Barb	7	11	Each
M00537310	08714729856900	Single Pigtail – Leading Barb	7	12	Each
M00537320	08714729856917	Single Pigtail – Leading Barb	7	13	Each
M00537330	08714729856924	Single Pigtail – Leading Barb	7	15	Each
M00537340	08714729856931	Single Pigtail – Leading Barb	7	18	Each

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



Section Table of Contents

Prev Page

29

Next Page


Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Pancreatic Plastic Stents

Naviflex™ RX Pancreatic Delivery System & Naviflex Pushers



Naviflex RX Pancreatic Delivery System & Naviflex Pushers



Order Number	GTIN	Description	Diameter (Fr)	Length	Packaging
M00535060	08714729856221	Naviflex RX Pancreatic Stent Delivery System	7	N/A	Each
M00535080	08714729856238	Naviflex RX Pancreatic Stent Delivery System	10	N/A	Each
M00535000	08714729856207	Pancreatic Stent RX Pusher	3	N/A	Each
M00535010	08714729856214	Pancreatic Stent RX Pusher	4 / 5	N/A	Each
M00535020	08714729856184	Pancreatic Stent Pusher	3	N/A	Each
M00535030	08714729856191	Pancreatic Stent Pusher	4 / 5	N/A	Each



Section Table of Contents

Prev Page ◀ 30 ▶ Next Page


Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Biliary Stone Management

RX Stone Management



Extractor™ Pro RX-S Retrieval Balloon Catheter – Distal Injection




Order Number	GTIN	Description	Balloon Inflated O.D. (mm)	Working Length (cm)	Catheter O.D. (Fr)	Recommended Guidewire (in)	Packaging
M00547300	08714729907138	Extractor Pro RX-S Retrieval Balloon	9.0 – 12.0	200	7 - 6	.035	Each
M00547310	08714729907152	Extractor Pro RX-S Retrieval Balloon	12.0 – 15.0	200	7 - 6	.035	Each
M00547320	08714729907176	Extractor Pro RX-S Retrieval Balloon	15.0 – 18.0	200	7 - 6	.035	Each



Extractor™ Pro RX-S Retrieval Balloon Catheter – Proximal Injection




Order Number	GTIN	Description	Balloon Inflated O.D. (mm)	Working Length (cm)	Catheter O.D. (Fr)	Recommended Guidewire (in)	Packaging
M00547330	08714729907145	Extractor Pro RX-S Retrieval Balloon	9.0 – 12.0	200	7 - 6	.035	Each
M00547340	08714729907169	Extractor Pro RX-S Retrieval Balloon	12.0 – 15.0	200	7 - 6	.035	Each
M00547350	08714729907183	Extractor Pro RX-S Retrieval Balloon	15.0 – 18.0	200	7 - 6	.035	Each



Extractor Pro RX Retrieval Balloon Catheter – Distal Injection



Order Number	GTIN	Description	Balloon Inflated O.D. (mm)	Working Length (cm)	Catheter O.D. (Fr)	Recommended Guidewire (in)	Packaging
M00547000	08714729790266	Extractor Pro RX Retrieval Balloon	9.0 – 12.0	200	7 - 6	.035	Each
M00547010	08714729790297	Extractor Pro RX Retrieval Balloon	12.0 – 15.0	200	7 - 6	.035	Each
M00547020	08714729790310	Extractor Pro RX Retrieval Balloon	15.0 – 18.0	200	7 - 6	.035	Each



Extractor Pro RX Retrieval Balloon Catheter – Proximal Injection




Order Number	GTIN	Description	Balloon Inflated O.D. (mm)	Working Length (cm)	Catheter O.D. (Fr)	Recommended Guidewire (in)	Packaging
M00547030	08714729790280	Extractor Pro RX Retrieval Balloon	9.0 – 12.0	200	7 - 6	.035	Each
M00547040	08714729790273	Extractor Pro RX Retrieval Balloon	12.0 – 15.0	200	7 - 6	.035	Each
M00547050	08714729790303	Extractor Pro RX Retrieval Balloon	15.0 – 18.0	200	7 - 6	.035	Each



Biliary Stone Management


Stone Management



Extractor™ Pro XL Retrieval Balloon Catheter – Distal Injection




Order Number	GTIN	Description	Balloon Inflated O.D. (mm)	Working Length (cm)	Catheter O.D. (Fr)	Recommended Guidewire (in)	Packaging
M00547100	08714729790389	Extractor Pro XL Retrieval Balloon	9.0 – 12.0	200	7 - 6	.035	Each
M00547110	08714729790402	Extractor Pro XL Retrieval Balloon	12.0 – 15.0	200	7 - 6	.035	Each
M00547120	08714729790648	Extractor Pro XL Retrieval Balloon	15.0 – 18.0	200	7 - 6	.035	Each



Extractor Pro XL Retrieval Balloon Catheter – Proximal Injection




Order Number	GTIN	Description	Balloon Inflated O.D. (mm)	Working Length (cm)	Catheter O.D. (Fr)	Recommended Guidewire (in)	Packaging
M00547130	08714729790396	Extractor Pro XL Retrieval Balloon	9.0 – 12.0	200	7 - 6	.035	Each
M00547140	08714729790365	Extractor Pro XL Retrieval Balloon	12.0 – 15.0	200	7 - 6	.035	Each
M00547150	08714729790372	Extractor Pro XL Retrieval Balloon	15.0 – 18.0	200	7 - 6	.035	Each



Extractor Pro DL Retrieval Balloon Catheter



Order Number	GTIN	Description	Balloon Inflated O.D. (mm)	Working Length (cm)	Catheter O.D. (Fr)	Recommended Guidewire (in)	Packaging
M00547200	08714729790631	Extractor Pro DL Retrieval Balloon	9.0 – 12.0	200	7 - 6	.035	Each
M00547210	08714729790655	Extractor Pro DL Retrieval Balloon	12.0 – 15.0	200	7 - 6	.035	Each
M00547220	08714729790662	Extractor Pro DL Retrieval Balloon	15.0 – 18.0	200	7 - 6	.035	Each



Trapezoid™ RX Wireguided Retrieval Baskets



Order Number	GTIN	Description	Open Basket Diameter (cm)	Recommended Guidewire (in)	Required Working Channel (mm)	Packaging
M00510860	08714729296379	Trapezoid RX Wireguided Retrieval Basket	1.5	.035	3.2	Each
M00510870	08714729296386	Trapezoid RX Wireguided Retrieval Basket	2.0	.035	3.2	Each
M00510880	08714729296393	Trapezoid RX Wireguided Retrieval Basket	2.5	.035	3.2	Each
M00510890*	08714729296409	Trapezoid RX Wireguided Retrieval Basket	3.0	.035	3.2	Each

Recommend Alliance™ II Inflation System for use as Mechanical Lithotripter; Order #5062, (See page 60).

*Refer to Instructions for Use when using this device; recommended for use with Biliary Calculi of >1.5cm.



Section Table of Contents

Prev Page

32

Next Page


Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Biliary Stone Management

CRE™ PRO Wireguided Balloon Dilatation Catheters



CRE PRO Wireguided Balloon Dilatation Catheters

Order Number	GTIN	Balloon Inflated O.D.		Balloon Length (cm)	Inflation Pressures		Minimum Working Channel (mm)	Catheter Size (Fr)	Working Length (cm)	Packaging
		(mm)	(Fr)		(atm)	(kPa)				
M00558600	08714729797487	6-7-8	18-21-24	5.5	3-6-10	304-608-1013	2.8 / 3.2*	7.5	180	Each
M00558610	08714729797494	8-9-10	24-27-30	5.5	3-5.5-9	304-557-912	2.8 / 3.2*	7.5	180	Each
M00558620	08714729797500	10-11-12	30-33-36	5.5	3-5-8	304-507-811	2.8 / 3.2*	7.5	180	Each
M00558630	08714729797517	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	2.8 / 3.2*	7.5	180	Each
M00558640	08714729797524	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	2.8 / 3.7*	7.5	180	Each
M00558650	08714729797531	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	2.8 / 3.7*	7.5	180	Each
M00558660	08714729797548	6-7-8	18-21-24	5.5	3-6-10	304-608-1013	2.8 / 3.2*	7.5	240	Each
M00558670	08714729797555	8-9-10	24-27-30	5.5	3-5.5-9	304-557-912	2.8 / 3.2*	7.5	240	Each
M00558680	08714729797562	10-11-12	30-33-36	5.5	3-5-8	304-507-811	2.8 / 3.2*	7.5	240	Each
M00558690	08714729797579	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	2.8 / 3.2*	7.5	240	Each
M00558700	08714729797586	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	2.8 / 3.7*	7.5	240	Each
M00558710	08714729797593	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	2.8 / 3.7*	7.5	240	Each

Recommended Inflation Device: Alliance™ II Integrated Inflation / Lithotripsy Device Order #5062, Alliance II Single-Use Syringe / Gauge Assembly Order #5060,

 (See page 60).

Guidewire included and preloaded in balloon dilatation catheter.

*Duodenoscope channel size.



Section Table of Contents

Prev Page

33

Next Page


Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Biliary Stricture Management

RX Dilatation



Hurricane™ RX Biliary Balloon Dilatation Catheter



Order Number	GTIN	Balloon Inflated O.D. (mm)	Balloon Length (cm)	Catheter O.D. (Fr)	Catheter Length (cm)	Maximum Inflation Pressure (atm)	Minimum Channel Size (mm)	Recommended Guidewire (in)	Packaging
M00545890	08714729283805	4	2	5.8	180	12	3.2	.035	Each
M00545900	08714729283812	4	4	5.8	180	12	3.2	.035	Each
M00545910	08714729283829	6	2	5.8	180	11	3.2	.035	Each
M00545920	08714729283836	6	4	5.8	180	11	3.2	.035	Each
M00545930	08714729283843	8	2	5.8	180	11	3.2	.035	Each
M00545940	08714729283850	8	4	5.8	180	11	3.2	.035	Each
M00545950	08714729283867	10	2	5.8	180	8	3.2	.035	Each
M00545960	08714729283874	10	4	5.8	180	8	3.2	.035	Each

Recommend inflation device: Encore™ 26 Inflation Device Order #6667, [▶](#) (See page 35).



Section Table of Contents

Prev Page

34

Next Page


Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Biliary Stricture Management

RX Cytology




RX Cytology Brush Wireguided Cytology Brush



Order Number	GTIN	Description	Brush Width (mm)	Catheter O.D. (Fr)	Working Length (cm)	Recommended Guidewire (in)	Packaging
M00545000	08714729268628	RX Cytology Brush	2.1	8	200	.035	Each

Biliary Dilatation




MaxForce™ Biliary Balloon Dilatation Catheter



Order Number	GTIN	Balloon Inflated O.D.		Balloon Length (cm)	Catheter Size		Working Length (cm)	Required Working Channel (mm)	Recommended Guidewire (in)	Packaging
		(mm)	(Fr)		(mm)	(Fr)				
M00567330	08714729097037	4	12	2	1.67	5	180	2.8	.035	Each
M00567340	08714729097051	4	12	4	1.67	5	180	2.8	.035	Each
M00567370	08714729097068	6	18	2	1.67	5	180	2.8	.035	Each
M00567380	08714729097075	6	18	4	1.67	5	180	2.8	.035	Each
M00567410	08714729097082	8	24	2	1.67	5	180	2.8	.035	Each
M00567420	08714729097099	8	24	3	1.67	5	180	2.8	.035	Each

RX Accessories



Encore™ 26 Inflation Device

Order Number	GTIN	Description	Packaging
M00566670	08714729755241	Encore 26 Inflation Device	Each



Section Table of Contents

Prev Page ◀ 35 ▶ Next Page

Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Notes and Future Product Updates



Cholangioscopy



SpyGlass™ DS Direct Visualization System

➤ Capital Equipment.....	37	➤ Disposable Devices.....	37
--------------------------	----	---------------------------	----


Autolith® Touch Biliary EHL System

➤ Autolith Touch System.....	38	➤ Biliary EHL Probe.....	38
➤ Autolith Touch System Accessories	38		



SpyGlass™ DS Direct Visualization System

Capital Equipment & Devices



Disposable Devices



Order Number	GTIN	Description	Packaging
M00546600	08714729863236	SpyScope™ DS Access & Delivery Catheter	Each
M00546270	08714729764458	SpyBite™ Biopsy Forceps	Each

Capital Equipment

Order Number	GTIN	Description	Packaging
M00546650	08714729874348	SpyGlass DS Digital Controller	Please contact your Boston Scientific Sales Representative for ordering information.

For more information on ordering a compatible EHL or Laser Probe for the SpyGlass DS Direct Visualization System, please contact your Boston Scientific Representative.



Autolith™ Touch Biliary EHL System

Autolith Touch Biliary EHL System

Biliary EHL Probe

Order Number	Description	Packaging
M00546620	1.9, 375cm Biliary Electrohydraulic Lithotripter Probe	Each

Autolith Touch System

M005466800	Autolith Touch Bipolar Electrohydraulic Lithotripter	Each
------------	--	------

Autolith Touch System Accessories

M00546750	Autolith Touch Extender Cable	Each
M00546760	Autolith Touch Foot Switch	Each



Notes and Future Product Updates



Notes and Future Product Updates



EUS

**Expect™ Slimline (SL) Endoscopic
Ultrasound Aspiration Needles**

➤ Expect Slimline (SL) Endoscopic Ultrasound
Aspiration Needle 41


**Acquire™ Endoscopic Ultrasound
Fine Needle Biopsy (FNB) Device**

➤ Acquire Endoscopic Ultrasound Fine Needle Biopsy (FNB)
Device 41

AXIOS™ System

➤ AXIOS Stent and Electrocautery Enhanced Delivery System ... 42
➤ AXIOS Stent and Delivery System 42

Expect™ Slimline (SL) Endoscopic Ultrasound Aspiration Needles



Expect Slimline (SL) Endoscopic Ultrasound Aspiration Needle*




Order Number	GTIN	Needle Size	Minimum Working Channel (mm)	Sheath Diameter (mm)	Packaging
M00555500	08714729848141	19 gauge (1.10mm)	2.8	1.83	Each
M00555510	08714729848158	22 gauge (0.72mm)	2.4	1.65	Each
M00555520	08714729848165	25 gauge (0.52mm)	2.4	1.52	Each
M00555530	08714729848172	19 gauge Flex (1.14mm)	2.8	1.73	Each
M00555501	08714729848714	19 gauge (1.10mm)	2.8	1.83	Box 5
M00555511	08714729848721	22 gauge (0.72mm)	2.4	1.65	Box 5
M00555521	08714729848738	25 gauge (0.52mm)	2.4	1.52	Box 5
M00555531	08714729848806	19 gauge Flex (1.14mm)	2.8	1.73	Box 5

*Packaging also includes: 20cc syringe and one-way stopcock

*Needle length: 0cm to 8cm, adjustable

*Working length: 137.5cm to 141.5cm, adjustable

Acquire™ Endoscopic Ultrasound Fine Needle Biopsy (FNB) Device



Acquire Endoscopic Ultrasound Fine Needle Biopsy (FNB) Device*




Order Number	GTIN	Needle Size	Minimum Working Channel (mm)	Sheath Diameter (mm)	Packaging
M00555540	08714729931782	22 gauge (0.72mm)	2.4	1.65	Please contact your Boston Scientific Sales Representative for ordering information
M00555560	08714729931805	25 gauge (0.52mm)	2.4	1.52	Please contact your Boston Scientific Sales Representative for ordering information
M00555580	08714729931829	19 gauge Flex (1.14mm)	2.8	1.73	Please contact your Boston Scientific Sales Representative for ordering information
M00555541	08714729931843	22 gauge (0.72mm)	2.4	1.65	Please contact your Boston Scientific Sales Representative for ordering information
M00555561	08714729931867	25 gauge (0.52mm)	2.4	1.52	Please contact your Boston Scientific Sales Representative for ordering information
M00555581	08714729931881	19 gauge Flex (1.14mm)	2.8	1.73	Please contact your Boston Scientific Sales Representative for ordering information

*Packaging also includes: 20cc syringe and one-way stopcock

*Needle length: 0cm to 8cm, adjustable

*Working length: 137.5cm to 141.5cm, adjustable

AXIOS™ System



AXIOS Stent and Electrocautery Enhanced Delivery System



Order Number	Flange Diameter (mm)	Lumen Diameter (mm)	Saddle Length (mm)	Catheter O.D. (Fr)	Catheter Working Length (cm)	Catheter Total Length (cm)	Packaging
M00553640	21	10	10	10.8	138	146	Each
M00553650	24	15	10	10.8	138	146	Each

AXIOS Stent and Delivery System



Order Number	Flange Diameter (mm)	Lumen Diameter (mm)	Saddle Length (mm)	Catheter O.D. (Fr)	Catheter Working Length (cm)	Catheter Total Length (cm)	Packaging
M00553740	21	10	10	10.8	138	146	Each
M00553750	24	15	10	10.8	138	146	Each

Enteral & Esophageal Stents



Stents (Enteral and Esophageal)

Colonic and Duodenal Stents

- Ultraflex™ Precision Single-Use Colonic Stent System 44
- WallFlex™ Colonic Stent System with Anchor Lock
Delivery System 43
- WallFlex Duodenal Stent System with Anchor Lock
Delivery System 43
- WALLSTENT™ Enteral Endoprosthesis Colonic / Duodenal
Stent with UNISTEP™ Plus Delivery System 44
- WALLSTENT Super Stiff Guidewire..... 44

Esophageal Stents

- Amplatz Super Stiff™ Guidewire..... 47
- Polyflex™ Single-Use Esophageal Stent System..... 46
- Ultraflex Single-Use Covered Esophageal NG Stent System –
Distal Release 46
- Ultraflex Single-Use Covered Esophageal NG Stent System –
Proximal Release..... 46
- Ultraflex Single-Use Covered Large Esophageal NG Stent
System – Distal Release..... 47
- Ultraflex Single-Use Covered Large Esophageal NG Stent
System – Proximal Release 47
- Ultraflex Single-Use Uncovered Esophageal NG Stent
System – Distal Release..... 47
- Ultraflex Single-Use Uncovered Esophageal NG Stent
System – Proximal Release 47
- WallFlex Fully Covered Esophageal Stent..... 45
- WallFlex Partially Covered Esophageal Stent 45



Enteral Stents

WallFlex™ Colonic and Duodenal Stents



 Magnetic Resonance Conditional¹

WallFlex Colonic Stent System with Anchor Lock Delivery System*



Order Number	GTIN	Stent Body Diameter (mm)	Stent Flare Diameter (mm)	Expanded Stent Length (mm)	Working Length (cm)	Delivery System Length (cm)	Delivery System Diameter (Fr)	Packaging
M00565040	08714729456513	25	30	60	230	270	10	Each
M00565050	08714729456520	25	30	90	230	270	10	Each
M00565060	08714729456537	25	30	120	230	270	10	Each
M00565070	08714729456544	25	30	60	135	175	10	Each
M00565080	08714729456551	25	30	90	135	175	10	Each
M00565090	08714729456568	25	30	120	135	175	10	Each
M00565100	08714729456575	22	27	60	230	270	10	Each
M00565110	08714729456582	22	27	90	230	270	10	Each
M00565120	08714729456599	22	27	120	230	270	10	Each
M00565130	08714729456605	22	27	60	135	175	10	Each
M00565140	08714729456612	22	27	90	135	175	10	Each
M00565150	08714729456629	22	27	120	135	175	10	Each

*Recommended Guidewire .035" WALLSTENT™ Super Stiff Guidewire, Order #18001,  (See page 44).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



 Magnetic Resonance Conditional¹

WallFlex Duodenal Stent System with Anchor Lock Delivery System*



Order Number	GTIN	Stent Body Diameter (mm)	Flare Diameter (mm)	Expanded Stent Length (mm)	Working Length (cm)	Delivery System Length (cm)	Delivery System Diameter (Fr)	Packaging
M00565010	08714729456483	22	27	60	230	270	10	Each
M00565020	08714729456490	22	27	90	230	270	10	Each
M00565030	08714729456506	22	27	120	230	270	10	Each

*Recommended Guidewire .035" Hydra Jagwire™ Guidewire Stiff Shaft 450cm Order #5606,  (See page 1).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



Section Table of Contents

Prev Page

43



Next Page

Main Table of Contents



Enteral Stents

WALLSTENT™ Colonic and Duodenal Stents (Continued)



 Magnetic Resonance Conditional¹

WALLSTENT Enteral Endoprosthesis Colonic / Duodenal Stent with UNISTEP™ Plus Delivery System*



Order Number	GTIN	Expanded Stent Diameter (mm)	Expanded Stent Length (mm)	Working Length (cm)	Delivery System Length (cm)	Delivery System Diameter (Fr)	Packaging
M00565560	08714729285052	20	60	230	255	10	Each
M00565570	08714729285069	20	90	230	255	10	Each
M00565580	08714729285076	22	60	230	255	10	Each
M00565590	08714729285083	22	90	230	255	10	Each
M00565620	08714729285113	20	60	135	160	10	Each
M00565630	08714729285120	20	90	135	160	10	Each
M00565640	08714729285137	22	60	135	160	10	Each
M00565650	08714729285144	22	90	135	160	10	Each


*Recommended Guidewire .035" WALLSTENT Super Stiff Guidewire, Order #18001, (see below).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

WALLSTENT Super Stiff Guidewire



Order Number	GTIN	Description	Distal Tip O.D. (in)	Length (cm)	Tip Style	Packaging
H965180011	08714729203070	WALLSTENT Super Stiff Guidewire	.035	500	Tapered	Box 5



 Magnetic Resonance Conditional¹

Ultraflex™ Precision Single-Use Colonic Stent System*



Order Number	GTIN	Stent Body Diameter (mm)	Stent Flare Diameter (mm)	Expanded Stent Length (mm)	Working Length (cm)	Delivery System Length (cm)	Delivery Tip Diameter (Fr)	Packaging
M00557330	08714729339175	25	30	57	100	105	22	Each
M00557360	08714729339199	25	30	87	100	105	22	Each
M00557380	08714729339205	25	30	117	100	105	22	Each

*Recommended Guidewire .038" Jagwire™ Guidewire, Order #5662,  (See page 2).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



Section Table of Contents


Prev Page  44  Next Page

Main Table of Contents



Esophageal Stents

WallFlex™ Esophageal Stents



 Magnetic Resonance Conditional¹


WallFlex Fully Covered Esophageal Stent*



Order Number	GTIN	Stent O.D. (mm)	Proximal / Distal Flares O.D. (mm)	Stent Length (cm)	Catheter Diameter (Fr / mm)	Working Length (cm)	Packaging
M00516700	08714729778035	18	25 / 23	10.3	18.5 / 6.2	78	Each
M00516710	08714729778042	18	25 / 23	12.3	18.5 / 6.2	78	Each
M00516720	08714729778059	18	25 / 23	15.3	18.5 / 6.2	78	Each
M00516730	08714729778066	23	28 / 28	10.5	18.5 / 6.2	78	Each
M00516740	08714729778073	23	28 / 28	12.5	18.5 / 6.2	78	Each
M00516750	08714729778080	23	28 / 28	15.5	18.5 / 6.2	78	Each

*Recommended Guidewire .038" Jagwire™ Guidewire, Order #5662,  (See page 2).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



 Magnetic Resonance Conditional¹

WallFlex Partially Covered Esophageal Stent*



Order Number	GTIN	Stent O.D. (mm)	Proximal / Distal Flares O.D. (mm)	Stent Length (cm)	Covered Length (cm)	Catheter Diameter (Fr / mm)	Working Length (cm)	Packaging
M00516900	08714729765240	18	23 / 23	10.3	7.1	18.5 / 6.2	78	Each
M00516910	08714729765257	18	23 / 23	12.3	9.1	18.5 / 6.2	78	Each
M00516920	08714729765264	18	23 / 23	15.3	12.1	18.5 / 6.2	78	Each
M00516930	08714729765271	23	28 / 28	10.5	7.3	18.5 / 6.2	78	Each
M00516940	08714729765288	23	28 / 28	12.5	9.3	18.5 / 6.2	78	Each
M00516950	08714729765295	23	28 / 28	15.5	12.3	18.5 / 6.2	78	Each

*Recommended Guidewire .038" Jagwire Guidewire, Order #5662,  (See page 2).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



Section Table of Contents

Prev Page

45


Next Page

Main Table of Contents



Esophageal Stents

Polyflex™ Stents




Polyflex Single-Use Esophageal Stent System*



Order Number	GTIN	Stent I.D. (mm)	Proximal Flare I.D. (mm)	Stent Length (mm)	Delivery System Diameter (mm)	Packaging
M00514270	08714729457015	16	20	90	12	Each
M00514280	08714729457022	16	20	120	12	Each
M00514290	08714729457039	16	20	150	12	Each
M00514300	08714729457046	18	23	90	13	Each
M00514310	08714729457053	18	23	120	13	Each
M00514320	08714729457060	18	23	150	13	Each
M00514330	08714729457077	21	25	90	14	Each
M00514340	08714729457084	21	25	120	14	Each
M00514350	08714729457091	21	25	150	14	Each

*Recommended Guidewire .038" Jagwire™ Guidewire, Order #5662, (See page 2).

Ultraflex™ Stents



Magnetic Resonance Conditional¹

Ultraflex Single-Use Covered Esophageal NG Stent System – Distal Release*



Order Number	GTIN	Stent O.D. (mm)	Proximal Flare O.D. (mm)	Stent Length (cm)	Covered Length (cm)	Packaging
M00513730	08714729716150	18	23	10	7	Each
M00513740	08714729716167	18	23	12	9	Each
M00513750	08714729649113	18	23	15	12	Each

*Recommended Guidewire .038" Jagwire Guidewire, Order #5662, (See page 2).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

Magnetic Resonance Conditional¹

Ultraflex Single-Use Covered Esophageal NG Stent System – Proximal Release*



Order Number	GTIN	Stent O.D. (mm)	Proximal Flare O.D. (mm)	Stent Length (cm)	Covered Length (cm)	Packaging
M00513840	08714729716099	18	23	10	7	Each
M00513850	08714729716075	18	23	12	9	Each
M00513860	08714729716082	18	23	15	12	Each

*Recommended Guidewire .038" Jagwire Guidewire, Order #5662, (See page 2).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



Section Table of Contents

Prev Page

46

Next Page


Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Esophageal Stents

Ultraflex™ Stents (Continued)



 Magnetic Resonance Conditional¹

Ultraflex Single-Use Covered Large Esophageal NG Stent System – Distal Release*



Order Number	GTIN	Stent O.D. (mm)	Proximal Flare O.D. (mm)	Stent Length (cm)	Covered Length (cm)	Packaging
M00514200	08714729121190	23	28	10	7	Each
M00514210	08714729121206	23	28	12	9	Each

*Recommended Guidewire .038" Jagwire™ Guidewire, Order #5662,  (See page 2).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

 Magnetic Resonance Conditional¹


Ultraflex Single-Use Covered Large Esophageal NG Stent System – Proximal Release*



Order Number	GTIN	Stent O.D. (mm)	Proximal Flare O.D. (mm)	Stent Length (cm)	Covered Length (cm)	Packaging
M00514240	08714729716105	23	28	10	7	Each
M00514250	08714729716112	23	28	12	9	Each

*Recommended Guidewire .038" Jagwire Guidewire, Order #5662,  (See page 2).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



 Magnetic Resonance Conditional¹

Ultraflex Single-Use Uncovered Esophageal NG Stent System – Distal Release*



Order Number	GTIN	Stent O.D. (mm)	Proximal Flare O.D. (mm)	Stent Length (cm)	Packaging
M00513700	08714729716129	18	23	7	Each
M00513710	08714729716136	18	23	10	Each
M00513720	08714729716143	18	23	15	Each

*Recommended Guidewire .038" Jagwire Guidewire, Order #5662,  (See page 2).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

 Magnetic Resonance Conditional¹

Ultraflex Single-Use Uncovered Esophageal NG Stent System – Proximal Release*



Order Number	GTIN	Stent O.D. (mm)	Proximal Flare O.D. (mm)	Stent Length (cm)	Packaging
M00513800	08714729716174	18	23	7	Each
M00513810	08714729716181	18	23	10	Each
M00513830	08714729716204	18	23	15	Each

*Recommended Guidewire .038" Jagwire Guidewire, Order #5662,  (See page 2).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

Amplatz Super Stiff™ Guidewire



Order Number	GTIN	Description	O.D. (in)	Length (cm)	Tip Style	Packaging
M00550090	08714729301929	Amplatz Super Stiff Guidewire	.038	260	Straight	Each

Notes and Future Product Updates





Hemostasis

Clipping Devices

- Resolution 360™ Clip 49
- Resolution™ Clip 49

Bipolar Electrohemostasis Catheters

- Gold Probe™ Single-Use Electrohemostasis Catheters 50
- Injection Gold Probe™ Single-Use 50
- Electrohemostasis / Injection Catheters 50

Needles

- Interject™ Clear Single-Use Injection Therapy..... 51
- Needle Catheters 51
- Interject Contrast Single-Use Injection Therapy
Needle Catheters 51
- Optiflo™ Single-Use Hemostasis Catheters 51

Multiple Band Ligators

- Speedband Superview Super 7™ Multiple Band Ligators 52

Electrosurgical Bipolar Cable Adapter

- Electrosurgical Bipolar Cable Adapter..... 53



Hemostasis

Clipping Devices



 Magnetic Resonance Conditional¹

Resolution 360™ Clip



Order Number	GTIN	Working Length (cm)	Minimum Working Channel (mm)	Clip Opening (mm)	Packaging
M00521240	08714729875666	155	2.8	11	Each
M00521241	08714729875673	155	2.8	11	Box 10
M00521242	08714729875680	155	2.8	11	Box 20
M00521230	08714729875628	235	2.8	11	Each
M00521231	08714729875635	235	2.8	11	Box 10
M00521232	08714729875642	235	2.8	11	Box 20
M00521233	08714729875659	235	2.8	11	Box 40

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



 Magnetic Resonance Conditional¹

Resolution™ Clip



Order Number	GTIN	Working Length (cm)	Minimum Working Channel (mm)	Clip Opening (mm)	Packaging
M00522600	08714729504764	155	2.8	11	Each
M00522601	08714729504771	155	2.8	11	Box 10
M00522602	08714729504788	155	2.8	11	Box 20
M00522610	08714729504795	235	2.8	11	Each
M00522611	08714729504801	235	2.8	11	Box 10
M00522612	08714729504818	235	2.8	11	Box 20

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



Section Table of Contents

Prev Page

49


Next Page

Main Table of Contents



Hemostasis

Bipolar Electrohemostasis Catheters




Injection Gold Probe™ Single-Use Electrohemostasis / Injection Catheters



Order Number	GTIN	Description	O.D. (Fr)	Needle Size Gauge	Working Length (cm)	Type Generator / Connector	Packaging
M00560150	08714729056492	Injection Gold Probe Catheter	7	25	210	Standard Plug*	Each
M00560160	08714729056508	Injection Gold Probe Catheter	10	25	210	Standard Plug*	Each

*Compatible with Endostat™ Bipolar-Monopolar Electrosurgical Generators.

For connections to other generators, please refer to the Electrosurgical Bipolar Cable Adapters, (See page 53).



Gold Probe™ Single-Use Electrohemostasis Catheters



Order Number	GTIN	Description	O.D. (Fr)	Working Length (cm)	Type Generator / Connector	Packaging
M00560070	08714729174691	Gold Probe Catheter	7	300	Standard Plug*	Each
M00560100	08714729174684	Gold Probe Catheter	10	300	Standard Plug*	Each
M00560220	08714729099574	Gold Probe Catheter	7	350	Standard Plug*	Each
M00560071	08714729071013	Gold Probe Catheter	7	300	Standard Plug*	Box 5
M00560101	08714729071020	Gold Probe Catheter	10	300	Standard Plug*	Box 5


*Compatible with Endostat Bipolar-Monopolar Electrosurgical Generators.

For connections to other generators, please refer to the Electrosurgical Bipolar Cable Adapters, (See page 53).



Hemostasis

Needles




Interject™ Contrast Single-Use Injection Therapy Needle Catheters*



Order Number	GTIN	Sheath Design	Needle Gauge	Maximum Needle Extension Length (mm)	Sheath O.D. (mm)	Working Length (cm)	Packaging
M00518151	08714729296430	Contrast	23	4	1.8	200	Box 5
M00518161	08714729296447	Contrast	25	4	1.8	200	Box 5
M00518251	08714729296478	Contrast	23	6	1.8	200	Box 5
M00518261	08714729296485	Contrast	25	6	1.8	200	Box 5
M00518351	08714729296522	Contrast	23	4	2.3	240	Box 5
M00518361	08714729296539	Contrast	25	4	2.3	240	Box 5

*With Star Catheter™ Technology.




Interject Clear Single-Use Injection Therapy Needle Catheters*



Order Number	GTIN	Sheath Design	Needle Gauge	Maximum Needle Extension Length (mm)	Sheath O.D. (mm)	Working Length (cm)	Packaging
M00518101	08714729296416	Clear	23	4	1.8	200	Box 5
M00518111	08714729296423	Clear	25	4	1.8	200	Box 5
M00518301	08714729296508	Clear	23	4	2.3	240	Box 5
M00518311	08714729296515	Clear	25	4	2.3	240	Box 5

*With Star Catheter Technology.



Optiflo™ Single-Use Hemostasis Catheters





Order Number	GTIN	Description	Needle Gauge	Sheath O.D. (mm)	Working Length (cm)	Working Channel (mm)	Packaging
M00511731	08714729144182	Optiflo Catheter	23	2.3	200	2.8	Box 5
M00511751	08714729145264	Optiflo Catheter	25	2.3	200	2.8	Box 5



Hemostasis

Multiple Band Ligators



Speedband Superview Super 7™ Multiple Band Ligators




Order Number	GTIN	Description	Required Working Channel (mm)	Scope O.D. (mm)	Bands / Ligating Unit	Packaging
M00542250	08714729201953	Speedband Superview Super 7 Band Ligator	2.8	8.6 – 11.5	7	Each
M00542251	08714729201960	Speedband Superview Super 7 Band Ligator	2.8	8.6 – 11.5	7	Box 2
M00542253	08714729504832	Speedband Superview Super 7 Band Ligator	2.8	8.6 – 11.5	7	Box 4



Hemostasis

Electrosurgical Bipolar Cable Adapter



Electrosurgical Bipolar Cable Adapter*

Order Number	GTIN	Description	Pin Spacing (mm)	Generator Compatibility**	Packaging
M00561440	08714729797609	Bipolar Cable Adapter	18	EMC (Coagulasem Dolley)	Each
M00561450	08714729797616	Bipolar Cable Adapter	22	ERBE, EMC Dolley, Olympus™ (ESG100 & PSD60)	Each
M00561460	08714729797623	Bipolar Cable Adapter	25	Lamidey (Optima 4)	Each
M00561470	08714729797630	Bipolar Cable Adapter	28	Valleylab™, Lamidey (Duo), Bovie™ (Aaron 3250™), Olympus (UES10, UES30 & Surgemaster), Bowa, Soring	Each
M00561480	08714729797647	Bipolar Cable Adapter	28 / Third Pin	ConMed™ (System 5000), Bovie (Aaron 2250™), Ethicon™	Each
M00561490	08714729797654	Bipolar Cable Adapter	32 / Third Pin	ConMed (Sable 180, Excalibur, System 7550)	Each
M00561500	08714729797661	Bipolar Cable Adapter	39	Codman, Malis™	Each

*Bipolar Cable Adapters are used to connect Gold Probe™ and Injection Gold Probe™ Catheters with the above mentioned generators.

Bipolar Cable Adapters are not required for the connection of Gold Probe and Injection Gold Probe to Endostat™ II or III Generators.

** For generators and their associated cable adaptors not listed above, please contact your Boston Scientific Representative.

Note: Monopolar active cords are in the Polypectomy Section, [▶](#) (See page 79).

Note: Some variations in generator faceplates may exist within a generator model. We recommend that you confirm pin spacing prior to ordering.

Notes and Future Product Updates





Dilation

CRE™ Balloon Dilatation Catheter

- CRE RX Biliary Balloon Dilatation Catheters 55
- CRE PRO Wireguided Esophageal / Pyloric / Biliary Balloon Dilatation Catheters..... 56
- CRE PRO Wireguided Esophageal / Pyloric / Colonic / Biliary Balloon Dilatation Catheters..... 56
- CRE Wireguided Esophageal / Pyloric Balloon Dilatation Catheters 57
- CRE Wireguided Esophageal / Pyloric / Colonic Balloon Dilatation Catheters 57
- CRE Fixed Wire Esophageal Balloon Dilatation Catheters 58

MaxForce TTS™ Esophageal High Performance Balloon Dilatation Catheter

- MaxForce TTS Esophageal High Performance Balloon Dilatation Catheters 59

Rigiflex™ II Single-Use Achalasia Balloon

- Rigiflex II Single-Use Achalasia Balloon 59

Inflation System



- Alliance™ II Integrated Inflation / Lithotripsy Device..... 60
- Alliance II Single-Use Syringe / Gauge Assembly 60
- CRE SteriFlate Disposable Inflation Device..... 60

Dilation Accessories

- Inflation / Deflation Devices and Accessories 60

Dilation

CRE™ RX Biliary Balloon Dilatation Catheters



CRE RX Biliary Balloon Dilatation Catheters

Order Number	GTIN	Balloon Inflated O.D.		Balloon Length (cm)	Inflation Pressures		Minimum Working Channel (mm)	Catheter Size (Fr)	Working Length (cm)	Packaging
		(mm)	(Fr)		(atm)	(kPa)				
M00558860	08714729965411	6-7-8	18-21-24	5.5	3-6-10	304-608-1013	3.7	7.5	180	Each
M00558870	08714729965428	8-9-10	24-27-30	5.5	3-5.5-9	304-557-912	3.7	7.5	180	Each
M00558880	08714729965435	10-11-12	30-33-36	5.5	3-5-8	304-507-811	3.7	7.5	180	Each
M00558890	08714729965442	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	3.7	7.5	180	Each
M00558900	08714729965459	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	3.7	7.5	180	Each
M00558910	08714729965466	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	3.7	7.5	180	Each
M00558920	08714729965473	8-9-10	24-27-30	3.0	3-5.5-9	304-608-1013	3.7	7.5	180	Each
M00558930	08714729965480	10-11-12	30-33-36	3.0	3-5-8	304-557-912	3.7	7.5	180	Each
M00558940	08714729965497	12-13.5-15	36-40.5-45	3.0	3-4.5-8	304-507-811	3.7	7.5	180	Each

Recommended Inflation Device: Alliance™ II Integrated Inflation / Lithotripsy Device Order #5062, Alliance II Single-Use Syringe / Gauge Assembly Order #5060.
No Pre-loaded Guidewire in the Balloon Dilatation Catheter.



Section Table of Contents

Prev Page

55

Next Page


Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Dilation

CRE™ PRO Balloon Dilatation Catheters



CRE PRO Wireguided Esophageal / Pyloric / Biliary Balloon Dilatation Catheters



Order Number	GTIN	Balloon Inflated O.D.		Balloon Length (cm)	Inflation Pressures		Minimum Working Channel (mm)	Catheter Size (Fr)	Working Length (cm)	Packaging
		(mm)	(Fr)		(atm)	(kPa)				
M00558600	08714729797487	6-7-8	18-21-24	5.5	3-6-10	304-608-1013	2.8 / 3.2*	7.5	180	Each
M00558610	08714729797494	8-9-10	24-27-30	5.5	3-5-9	304-557-912	2.8 / 3.2*	7.5	180	Each
M00558620	08714729797500	10-11-12	30-33-36	5.5	3-5-8	304-507-811	2.8 / 3.2*	7.5	180	Each
M00558630	08714729797517	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	2.8 / 3.2*	7.5	180	Each
M00558640	08714729797524	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	2.8 / 3.7*	7.5	180	Each
M00558650**	08714729797531	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	2.8 / 3.7*	7.5	180	Each

Recommended Inflation Device: Alliance™ II Integrated Inflation / Lithotripsy Device Order #5062, Alliance II Syringe / Gauge Assembly Order #5060,

(See page 60).

Guidewire included and preloaded in balloon dilatation catheter.

*Duodenoscope channel size.

**Esophageal and Biliary application only.



CRE PRO Wireguided Esophageal / Pyloric / Colonic / Biliary Balloon Dilatation Catheters



Order Number	GTIN	Balloon Inflated O.D.		Balloon Length (cm)	Inflation Pressures		Minimum Working Channel (mm)	Catheter Size (Fr)	Working Length (cm)	Packaging
		(mm)	(Fr)		(atm)	(kPa)				
M00558660	08714729797548	6-7-8	18-21-24	5.5	3-6-10	304-608-1013	2.8 / 3.2*	7.5	240	Each
M00558670	08714729797555	8-9-10	24-27-30	5.5	3-5-9	304-557-912	2.8 / 3.2*	7.5	240	Each
M00558680	08714729797562	10-11-12	30-33-36	5.5	3-5-8	304-507-811	2.8 / 3.2*	7.5	240	Each
M00558690	08714729797579	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	2.8 / 3.2*	7.5	240	Each
M00558700	08714729797586	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	2.8 / 3.7*	7.5	240	Each
M00558710**	08714729797593	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	2.8 / 3.7*	7.5	240	Each

Recommended Inflation Device: Alliance II Integrated Inflation / Lithotripsy Device Order #5062, Alliance II Syringe / Gauge Assembly Order #5060,

(See page 60).

Guidewire included and preloaded in balloon dilatation catheter.

*Duodenoscope channel size.

**Esophageal, Colonic and Biliary application only.



Section Table of Contents

Prev Page

56


Next Page

Main Table of Contents



Dilation

CRE™ Balloon Dilatation Catheters



CRE Wireguided Esophageal / Pyloric Balloon Dilatation Catheters



Order Number	GTIN	Balloon Inflated O.D.		Balloon Length (cm)	Inflation Pressures		Minimum Working Channel (mm)	Catheter Size (Fr)	Working Length (cm)	Packaging
		(mm)	(Fr)		(atm)	(kPa)				
M00558390	08714729339328	6-7-8	18-21-24	5.5	3-6-10	304-608-1013	2.8	7.5	180	Each
M00558400	08714729339335	8-9-10	24-27-30	5.5	3-5.5-9	304-557-912	2.8	7.5	180	Each
M00558410	08714729339342	10-11-12	30-33-36	5.5	3-5-8	304-507-811	2.8	7.5	180	Each
M00558420	08714729339359	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	2.8	7.5	180	Each
M00558430	08714729339366	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	2.8	7.5	180	Each
M00558440*	08714729339373	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	2.8	7.5	180	Each

Recommended Inflation Device: Alliance™ II Integrated Inflation / Lithotripsy Device Order #5062, Alliance II Single-Use Syringe / Gauge Assembly Order #5060, (See page 60).

Guidewire included and preloaded in balloon dilatation catheter.

*Esophageal application only.

CRE Wireguided Esophageal / Pyloric / Colonic Balloon Dilatation Catheters



Order Number	GTIN	Balloon Inflated O.D.		Balloon Length (cm)	Inflation Pressures		Minimum Working Channel (mm)	Catheter Size (Fr)	Working Length (cm)	Packaging
		(mm)	(Fr)		(atm)	(kPa)				
M00558450	08714729339380	6-7-8	18-21-24	5.5	3-6-10	304-608-1013	2.8	7.5	240	Each
M00558460	08714729339397	8-9-10	24-27-30	5.5	3-5.5-9	304-557-912	2.8	7.5	240	Each
M00558470	08714729339403	10-11-12	30-33-36	5.5	3-5-8	304-507-811	2.8	7.5	240	Each
M00558480	08714729339410	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	2.8	7.5	240	Each
M00558490	08714729339427	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	2.8	7.5	240	Each
M00558500*	08714729339434	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	2.8	7.5	240	Each

Recommended Inflation Device: Alliance II Integrated Inflation / Lithotripsy Device Order #5062, Alliance II Single-Use Syringe / Gauge Assembly Order #5060, (See page 60).

Guidewire included and preloaded in balloon dilatation catheter.

*Esophageal and Colonic application only.



Section Table of Contents


Prev Page **57** Next Page

Main Table of Contents



Dilation

CRE™ Balloon Dilatation Catheters (Continued)



CRE Fixed Wire Esophageal Balloon Dilatation Catheters



Order Number	GTIN	Balloon Inflated O.D.		Balloon Length (cm)	Inflation Pressures		Minimum Working Channel (mm)	Catheter Size (Fr)	Working Length (cm)	Packaging
		(mm)	(Fr)		(atm)	(kPa)				
M00558330	08714729201984	6-7-8	18-21-24	8	3-6-10	304-608-1013	2.8	6	180	Each
M00558340	08714729201991	8-9-10	24-27-30	8	3-5.5-9	304-557-912	2.8	6	180	Each
M00558350	08714729202004	10-11-12	30-33-36	8	3-5-8	304-507-811	2.8	6	180	Each
M00558360	08714729195986	12-13.5-15	36-40.5-45	8	3-4.5-8	304-456-811	2.8	6	180	Each
M00558370	08714729195993	15-16.5-18	45-49.5-54	8	3-4.5-7	304-456-709	2.8	6	180	Each
M00558380	08714729202011	18-19-20	54-57-60	8	3-4.5-6	304-456-608	2.8	6	180	Each
M00558361	08714729339298	12-13.5-15	36-40.5-45	8	3-4.5-8	304-456-811	2.8	6	180	Box 5
M00558371	08714729339304	15-16.5-18	45-49.5-54	8	3-4.5-7	304-456-709	2.8	6	180	Box 5
M00558381	08714729339311	18-19-20	54-57-60	8	3-4.5-6	304-456-608	2.8	6	180	Box 5

Recommended Inflation Device: Alliance™ II Integrated Inflation / Lithotripsy Device Order #5062, Alliance II Single-Use Syringe / Gauge Assembly Order #5060,

 (See page 60).



Section Table of Contents

Prev Page

58

Next Page


Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Dilation

MaxForce TTS™ Esophageal High Performance Balloon Dilatation Catheter



MaxForce TTS™ Esophageal High Performance Balloon Dilatation Catheter




Order Number	GTIN	Balloon Inflated O.D.		Balloon Length (cm)	Inflation Pressures		Minimum Working Channel (mm)	Catheter Size		Working Length (cm)	Packaging
		(mm)	(Fr)		(atm)	(kPa)		(mm)	(Fr)		
M00589060	08714729099581	6	18	6	10	1013	2.8	2.0	6	180	Each
M00589080	08714729099598	8	24	6	10	1013	2.8	2.0	6	180	Each
M00589100	08714729099895	10	30	6	7	709	2.8	2.0	6	180	Each
M00589120	08714729099901	12	36	6	7	709	2.8	2.0	6	180	Each
M00589140	08714729080602	14	42	6	7	709	2.8	2.0	6	180	Each
M00589150	08714729099918	15	45	6	5	507	2.8	2.0	6	180	Each
M00589160	08714729080619	16	48	6	5	507	2.8	2.0	6	180	Each
M00589180	08714729099925	18	54	6	5	507	2.8	2.0	6	180	Each

Recommended Inflation Device: Alliance™ II Integrated Inflation / Lithotripsy Device Order #5062, Sterile Syringe Order #5060,



(See page 60).

Rigiflex™ II Single-Use Achalasia Balloon Dilatation Catheter



Rigiflex II Single-Use Achalasia Balloon Dilators*



Order Number	GTIN	Balloon Inflated O.D.		Balloon Length (cm)	Inflation Pressures		Working Length (cm)	Catheter Size		Recommended Guidewire** (in)	Packaging
		(mm)	(Fr)		(psi)	(kPa)		(mm)	(Fr)		
M00554500	08714729719106	30	90	10	20	138	90	4.6	14	.038	Each
M00554510	08714729719113	35	105	10	20	138	90	4.6	14	.038	Each
M00554520	08714729719120	40	120	10	20	138	90	4.6	14	.038	Each

*Product Shipped Sterile. Recommended use with Achalasia Pneumatic Hand Pump and Monitor Order #5320,



(See page 60).

**Guidewire included with Rigiflex II Balloon Dilatation Catheter.



Section Table of Contents

Prev Page

59

Next Page


Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Dilation

Inflation System



Alliance™ II Single-Use Syringe / Gauge Assembly / Integrated Inflation / Lithotripsy Device

Order Number	GTIN	Description	Packaging
M00550620	08714729283881	Alliance II Integrated Inflation / Lithotripsy Device	Each
M00550601	08714729129332	 Alliance II Single-Use Syringe / Gauge Assembly	Box 5

Dilation Accessories

Inflation / Deflation Devices and Accessories

Order Number	GTIN	Description	Packaging
M00553200	08714729195450	Pneumatic Inflator	Each
M00550630	00856877002506	CRE™ SteriFlate Disposable Inflation Device	Each

[Section Table of Contents](#)[Prev Page](#) ◀ 60 ▶ [Next Page](#)[Main Table of Contents](#)




Biopsy

Forceps

➤ Multibite™ Single-Use Multiple Sample Biopsy Forceps	63
➤ Radial Jaw™ 4 Single-Use Hot Biopsy Forceps	62
➤ Radial Jaw 4 Single-Use Jumbo Biopsy Forceps	61
➤ Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps	61
➤ Radial Jaw 4 Single-Use Pediatric Biopsy Forceps	62
➤ Radial Jaw 4 Single-Use Standard Capacity Biopsy Forceps ...	61

Biopsy

Forceps



Radial Jaw™ 4 Single-Use Jumbo Biopsy Forceps



Order Number	GTIN	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Color Code	Packaging
M00513343	08714729761891	Radial Jaw 4 Jumbo Forceps	2.8	160	3.2	Yellow	Box 5
M00513353	08714729761907	Radial Jaw 4 Jumbo Forceps w/ Needle	2.8	160	3.2	Yellow	Box 5
M00513363	08714729761914	Radial Jaw 4 Jumbo Forceps	2.8	240	3.2	Orange	Box 5
M00513373	08714729761938	Radial Jaw 4 Jumbo Forceps w/ Needle	2.8	240	3.2	Orange	Box 5
M00513361	08714729648871	Radial Jaw 4 Jumbo Forceps	2.8	240	3.2	Orange	Box 20
M00513371	08714729648895	Radial Jaw 4 Jumbo Forceps w/ Needle	2.8	240	3.2	Orange	Box 20
M00513362	08714729648888	Radial Jaw 4 Jumbo Forceps	2.8	240	3.2	Orange	Box 40
M00513372	08714729648901	Radial Jaw 4 Jumbo Forceps w/ Needle	2.8	240	3.2	Orange	Box 40

Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps



Order Number	GTIN	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Color Code	Packaging
M00513303	08714729761839	Radial Jaw 4 Large Capacity Forceps	2.4	160	2.8	Yellow	Box 5
M00513313	08714729761853	Radial Jaw 4 Large Capacity Forceps w/ Needle	2.4	160	2.8	Yellow	Box 5
M00513323	08714729761877	Radial Jaw 4 Large Capacity Forceps	2.4	240	2.8	Orange	Box 5
M00513333	08714729755494	Radial Jaw 4 Large Capacity Forceps w/ Needle	2.4	240	2.8	Orange	Box 5
M00513301	08714729648758	Radial Jaw 4 Large Capacity Forceps	2.4	160	2.8	Yellow	Box 20
M00513311	08714729648772	Radial Jaw 4 Large Capacity Forceps w/ Needle	2.4	160	2.8	Yellow	Box 20
M00513321	08714729648796	Radial Jaw 4 Large Capacity Forceps	2.4	240	2.8	Orange	Box 20
M00513331	08714729648819	Radial Jaw 4 Large Capacity Forceps w/ Needle	2.4	240	2.8	Orange	Box 20
M00513302	08714729648765	Radial Jaw 4 Large Capacity Forceps	2.4	160	2.8	Yellow	Box 40
M00513312	08714729648789	Radial Jaw 4 Large Capacity Forceps w/ Needle	2.4	160	2.8	Yellow	Box 40
M00513322	08714729648802	Radial Jaw 4 Large Capacity Forceps	2.4	240	2.8	Orange	Box 40
M00513332	08714729648826	Radial Jaw 4 Large Capacity Forceps w/ Needle	2.4	240	2.8	Orange	Box 40

Radial Jaw 4 Single-Use Standard Capacity Biopsy Forceps



Order Number	GTIN	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Color Code	Packaging
M00513383	08714729777915	Radial Jaw 4 Standard Capacity Forceps	2.2	160	2.8	Yellow	Box 5
M00513393	08714729777939	Radial Jaw 4 Standard Capacity Forceps w/ Needle	2.2	160	2.8	Yellow	Box 5
M00513403	08714729777922	Radial Jaw 4 Standard Capacity Forceps	2.2	240	2.8	Orange	Box 5
M00513413	08714729777946	Radial Jaw 4 Standard Capacity Forceps w/ Needle	2.2	240	2.8	Orange	Box 5
M00513381	08714729777953	Radial Jaw 4 Standard Capacity Forceps	2.2	160	2.8	Yellow	Box 20
M00513391	08714729777977	Radial Jaw 4 Standard Capacity Forceps w/ Needle	2.2	160	2.8	Yellow	Box 20
M00513401	08714729777960	Radial Jaw 4 Standard Capacity Forceps	2.2	240	2.8	Orange	Box 20
M00513411	08714729777984	Radial Jaw 4 Standard Capacity Forceps w/ Needle	2.2	240	2.8	Orange	Box 20
M00513382	08714729778011	Radial Jaw 4 Standard Capacity Forceps	2.2	160	2.8	Yellow	Box 40
M00513392	08714729777991	Radial Jaw 4 Standard Capacity Forceps w/ Needle	2.2	160	2.8	Yellow	Box 40
M00513402	08714729778028	Radial Jaw 4 Standard Capacity Forceps	2.2	240	2.8	Orange	Box 40
M00513412	08714729778004	Radial Jaw 4 Standard Capacity Forceps w/ Needle	2.2	240	2.8	Orange	Box 40



Biopsy



Forceps (Continued)

Radial Jaw™ 4 Single-Use Pediatric Biopsy Forceps



Order Number	GTIN	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Packaging
M00513443	08714729790211	Radial Jaw 4 Pediatric Biopsy Forceps	1.8	160	≥ 2.0	Box 5
M00513453	08714729790228	Radial Jaw 4 Pediatric Biopsy Forceps w/ Needle	1.8	160	≥ 2.0	Box 5

Hot Biopsy Forceps



Radial Jaw 4 Single-Use Hot Biopsy Forceps



Order Number	GTIN	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Packaging
M00515031	08714729790136	Radial Jaw 4 Hot Biopsy Forceps	2.2	240	≥ 2.8	Box 5
M00515032	08714729790105	Radial Jaw 4 Hot Biopsy Forceps	2.2	240	≥ 2.8	Box 20
M00515033	08714729790143	Radial Jaw 4 Hot Biopsy Forceps	2.2	240	≥ 2.8	Box 40



Section Table of Contents

Prev Page

62


Next Page

Main Table of Contents



Biopsy

Forceps (Continued)



Multibite™ Single-Use Multiple Sample Biopsy Forceps ⓧ

Order Number	GTIN	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Color Code	Packaging
M00510101	08714729164975	Multibite Multiple Sample Biopsy Forceps	2.4	160	2.8	Yellow	Box 5
M00510121	08714729164982	Multibite Multiple Sample Biopsy Forceps	2.4	240	2.8	Orange	Box 5



Notes and Future Product Updates



Enteral Access



EndoVive™ Enteral Access Devices

Percutaneous Endoscopic Gastrostomy (PEG) Devices with ENFit™

Initial Placement PEG Kits with ENFit

- Safety PEG Kits with ENFit 65
- Standard PEG Kits with ENFit 65
- Initial Placement PEG Accessories with ENFit 65

One-Step Button Low Profile PEG Kits with ENFit

- One-Step Button™ Low Profile PEG Kits with ENFit 66
- One-Step Button Accessories with ENFit 66

Replacement Gastrostomy Tubes (G-Tube) with ENFit

- Replacement G-Tubes, Straight Bolster, with ENFit 67
- Replacement G-Tubes, Angle Bolster, with ENFit 67

Replacement Low Profile Buttons and Accessories with ENFit

- Low Profile Replacement Buttons with ENFit 67
- Low Profile Replacement Button Accessories with ENFit 68

Transition Connectors

- ENFit Transition Connectors 68

Percutaneous Endoscopic Gastrostomy (PEG) Devices

Initial Placement PEG Kits

- Safety PEG Kits 69
- Standard PEG Kits 69
- Initial Placement PEG Accessories 69

Initial Placement Low Profile PEG Kits

- One-Step Button Low Profile PEG Kits 70
- One-Step Button Accessories 70

Replacement Gastrostomy Tubes (G-Tube)

- Replacement G-Tubes, Straight Bolster 71
- Replacement G-Tubes, Angle Bolster 71

Replacement Low Profile Buttons and Accessories

- Low Profile Replacement Buttons with ENFit 72
- Low Profile Replacement Button Accessories with ENFit 72


Jejunal Feeding Tubes

- Through-the-PEG (TTP) Jejunal Feeding Tube 73



EndoVive™ Enteral Access Devices

Initial Placement PEG Kits with ENFit†




Safety PEG* Kit with ENFit



Order Number	GTIN	Description	Packaging
M00509001	08714729880394	20Fr Safety PEG Kit, Pull with ENFit	Box 2
M00509011	08714729880417	20Fr Safety PEG Kit, Push with ENFit	Box 2
M00509021	08714729880431	24Fr Safety PEG Kit, Pull with ENFit	Box 2
M00509031	08714729880455	24Fr Safety PEG Kit, Push with ENFit	Box 2

*Includes local anesthetic.



Standard PEG* Kit with ENFit



Order Number	GTIN	Description	Packaging
M00509041	08714729880479	20Fr Standard PEG Kit, Pull with ENFit	Box 2
M00509051	08714729880493	20Fr Standard PEG Kit, Push with ENFit	Box 2
M00509061	08714729880516	24Fr Standard PEG Kit, Pull with ENFit	Box 2
M00509071	08714729880530	24Fr Standard PEG Kit, Push with ENFit	Box 2

*Includes local anesthetic.

Initial Placement PEG Accessories with ENFit



Order Number	GTIN	Description	Packaging
M00501261	08714729863564	20Fr Y-Ports with ENFit	Box 10
M00501271	08714729863588	24Fr Y-Ports with ENFit	Box 10
M00580731	08714729415343	20Fr Round Bolsters	Box 10
M00580781	08714729415398	20Fr / 24Fr C Clamps	Box 10

†NOTE: Make sure all appropriate service lines in your healthcare system are ready to use ENFit enteral connectors. Please contact your Boston Scientific Sales Representative for support with ENFit enteral connectors.

**Recommended ENFit Transition Connectors on [Page 68](#).



Section Table of Contents

Prev Page **65** Next Page



Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at **1.888.272.1001**

EndoVive™ Enteral Access Devices

Initial Placement Low Profile PEG Kits with ENFit†



One-Step Button™ Low Profile PEG Kits with ENFit



Order Number	GTIN	Description	Length (cm)	Packaging
M00509690	08714729897415	18Fr Low Profile PEG Kit, Pull with ENFit	1.2	Each
M00509700	08714729897422	18Fr Low Profile PEG Kit, Pull with ENFit	1.7	Each
M00509710	08714729897439	18Fr Low Profile PEG Kit, Pull with ENFit	2.4	Each
M00509720	08714729897446	18Fr Low Profile PEG Kit, Pull with ENFit	3.4	Each
M00509730	08714729897453	18Fr Low Profile PEG Kit, Pull with ENFit	4.4	Each
M00509740	08714729897491	18Fr Low Profile PEG Kit, Push with ENFit	1.7	Each
M00509750	08714729897507	18Fr Low Profile PEG Kit, Push with ENFit	2.4	Each
M00509580	08714729897378	24Fr Low Profile PEG Kit, Pull with ENFit	1.7	Each
M00509590	08714729897385	24Fr Low Profile PEG Kit, Pull with ENFit	2.4	Each
M00509600	08714729897392	24Fr Low Profile PEG Kit, Pull with ENFit	3.4	Each
M00509610	08714729897408	24Fr Low Profile PEG Kit, Pull with ENFit	4.4	Each
M00509630	08714729897460	24Fr Low Profile PEG Kit, Push with ENFit	2.4	Each
M00509620	08714729897477	24Fr Low Profile PEG Kit, Push with ENFit	3.4	Each
M00509680	08714729897484	24Fr Low Profile PEG Kit, Push with ENFit	4.4	Each

One-Step Button and Low Profile Replacement Button Accessories with ENFit



Order Number	GTIN	Description	Length (cm)	Packaging
M00501121	08714729863380	18Fr Decompression Tube with ENFit	1.2	Box 10
M00501131	08714729863403	18Fr Decompression Tube with ENFit	1.7	Box 10
M00501141	08714729863427	18Fr Decompression Tube with ENFit	2.4	Box 10
M00501151	08714729863441	18Fr Decompression Tube with ENFit	3.4	Box 10
M00501161	08714729863465	18Fr Decompression Tube with ENFit	4.4	Box 10
M00501181	08714729863489	24Fr Decompression Tube with ENFit	1.7	Box 10
M00501191	08714729863502	24Fr Decompression Tube with ENFit	2.4	Box 10
M00501201	08714729863526	24Fr Decompression Tube with ENFit	3.4	Box 10
M00501211	08714729863540	24Fr Decompression Tube with ENFit	4.4	Box 10
M00501091	08714729863335	18Fr Bolus Feeding Set Straight Adapter with ENFit	—	Box 10
M00501061	08714729863281	18Fr Right Angle Feeding Set Continuous Adapter with ENFit	—	Box 10
M00501101	08714729863359	24Fr Bolus Feeding Set Straight Adapter with ENFit	—	Box 10
M00501071	08714729863304	24Fr Right Angle Feeding Set Continuous Adapter with ENFit	—	Box 10

†NOTE: Make sure all appropriate service lines in your healthcare system are ready to use ENFit enteral connectors. Please contact your Boston Scientific Sales Representative for support with ENFit enteral connectors.

**Recommended ENFit Transition Connectors on [Page 68](#).



Section Table of Contents

Prev Page

66

Next Page

Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

EndoVive™ Enteral Access Devices

Replacement Gastrostomy Tubes with ENFit™†



Replacement Gastrostomy Feeding Tube with ENFit, Straight Bolster



Order Number	GTIN	Description	Packaging
M00509881	00813939023862	12Fr Replacement G-Tube, Straight with ENFit	Box 5
M00509891	00813939023879	14Fr Replacement G-Tube, Straight with ENFit	Box 5
M00509901	00813939023886	16Fr Replacement G-Tube, Straight with ENFit	Box 5
M00509911	00813939023893	18Fr Replacement G-Tube, Straight with ENFit	Box 5
M00509921	00813939023909	20Fr Replacement G-Tube, Straight with ENFit	Box 5
M00509931	00813939023916	22Fr Replacement G-Tube, Straight with ENFit	Box 5
M00509941	00813939023923	24Fr Replacement G-Tube, Straight with ENFit	Box 5




Replacement Gastrostomy Feeding Tube with ENFit, Right Angle Bolster



Order Number	GTIN	Description	Packaging
M00509961	00813939024005	12Fr Replacement G-Tube, Right Angle with ENFit	Box 5
M00509971	00813939024012	14Fr Replacement G-Tube, Right Angle with ENFit	Box 5
M00509981	00813939024029	16Fr Replacement G-Tube, Right Angle with ENFit	Box 5
M00509991	00813939024036	18Fr Replacement G-Tube, Right Angle with ENFit	Box 5
M00510001	00813939024043	20Fr Replacement G-Tube, Right Angle with ENFit	Box 5
M00510011	00813939024050	22Fr Replacement G-Tube, Right Angle with ENFit	Box 5
M00510021	00813939024067	24Fr Replacement G-Tube, Right Angle with ENFit	Box 5

Low Profile Accessories & Replacements with ENFit



Low Profile Replacement Button with ENFit



Order Number	GTIN	Description	Length (cm)	Packaging
M00509830	08714729902614	18Fr Low Profile Replacement Button with ENFit	1.2	Each
M00509840	08714729902621	18Fr Low Profile Replacement Button with ENFit	1.7	Each
M00509850	08714729902638	18Fr Low Profile Replacement Button with ENFit	2.4	Each
M00509860	08714729902645	18Fr Low Profile Replacement Button with ENFit	3.4	Each
M00509870	08714729902652	18Fr Low Profile Replacement Button with ENFit	4.4	Each
M00509410	08714729902669	24Fr Low Profile Replacement Button with ENFit	1.7	Each
M00509420	08714729902676	24Fr Low Profile Replacement Button with ENFit	2.4	Each
M00509430	08714729902683	24Fr Low Profile Replacement Button with ENFit	3.4	Each
M00509440	08714729902690	24Fr Low Profile Replacement Button with ENFit	4.4	Each

†NOTE: Make sure all appropriate service lines in your healthcare system are ready to use ENFit enteral connectors. Please contact your Boston Scientific Sales Representative for support with ENFit enteral connectors.

**Recommended ENFit Transition Connectors on [Page 68](#).



Section Table of Contents

Prev Page

67

Next Page


Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

EndoVive™ Enteral Access Devices

Low Profile Replacement Button Accessories with ENFit™†




One-Step Button and Low Profile Replacement Button Accessories with ENFit



Order Number	GTIN	Description	Length (cm)	Packaging
M00501121	08714729863380	18Fr Decompression Tube with ENFit	1.2	Box 10
M00501131	08714729863403	18Fr Decompression Tube with ENFit	1.7	Box 10
M00501141	08714729863427	18Fr Decompression Tube with ENFit	2.4	Box 10
M00501151	08714729863441	18Fr Decompression Tube with ENFit	3.4	Box 10
M00501161	08714729863465	18Fr Decompression Tube with ENFit	4.4	Box 10
M00501181	08714729863489	24Fr Decompression Tube with ENFit	1.7	Box 10
M00501191	08714729863502	24Fr Decompression Tube with ENFit	2.4	Box 10
M00501201	08714729863526	24Fr Decompression Tube with ENFit	3.4	Box 10
M00501211	08714729863540	24Fr Decompression Tube with ENFit	4.4	Box 10
M00501091	08714729863335	18Fr Bolus Feeding Set Straight Adapter with ENFit	—	Box 10
M00501061	08714729863281	18Fr Right Angle Feeding Set Continuous Adapter with ENFit	—	Box 10
M00501101	08714729863359	24Fr Bolus Feeding Set Straight Adapter with ENFit	—	Box 10
M00501071	08714729863304	24Fr Right Angle Feeding Set Continuous Adapter with ENFit	—	Box 10

ENFit Transition Connectors



ENFit Transition Connectors



Order Number	GTIN	Description	Packaging
M00509221	N/A	Distal End ENFit Transition Connector	Box 10
M00509241	N/A	ENFit Feeding Port Transition Connector	Box 10

†NOTE: Make sure all appropriate service lines in your healthcare system are ready to use ENFit enteral connectors. Please contact your Boston Scientific Sales Representative for support with ENFit enteral connectors.



Section Table of Contents

Prev Page

68

Next Page


Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001


EndoVive™ Enteral Access Devices

Initial Placement PEG Kits



Safety PEG* Kit				
Order Number	GTIN	Description		Packaging
M00566461	08714729802624	20Fr Safety PEG Kit, Pull		Box 2
M00566471	08714729802631	20Fr Safety PEG Kit, Push		Box 2
M00566481	08714729802648	24Fr Safety PEG Kit, Pull		Box 2
M00566491	08714729802655	24Fr Safety PEG Kit, Push		Box 2

*Includes local anesthetic.



Standard PEG* Kit				
Order Number	GTIN	Description		Packaging
M00568201	08714729285199	20Fr Standard PEG Kit, Pull		Box 2
M00568211	08714729285205	20Fr Standard PEG Kit, Push		Box 2
M00568241	08714729296799	24Fr Standard PEG Kit, Pull		Box 2
M00568251	08714729296805	24Fr Standard PEG Kit, Push		Box 2

*Includes local anesthetic.

Initial Placement PEG Accessories				
Order Number	GTIN	Description		Packaging
M00580641	08714729169116	20Fr Y-Ports		Box 10
M00580651	08714729169123	24Fr Y-Ports		Box 10
M00580731	08714729415343	20Fr Round Bolsters		Box 10
M00580781	08714729415398	20Fr / 24Fr C Clamps		Box 10



Section Table of Contents

Prev Page

69

Next Page



Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

EndoVive™ Enteral Access Devices

Initial Placement Low Profile PEG Kits




One-Step Button™ Low Profile PEG Kits				
Order Number	GTIN	Description	Length (cm)	Packaging
M00563120	08714729284246	18Fr Low Profile PEG Kit, Pull	1.2	Each
M00563130	08714729284253	18Fr Low Profile PEG Kit, Pull	1.7	Each
M00563140	08714729284260	18Fr Low Profile PEG Kit, Pull	2.4	Each
M00563150	08714729284277	18Fr Low Profile PEG Kit, Pull	3.4	Each
M00563160	08714729284284	18Fr Low Profile PEG Kit, Pull	4.4	Each
M00563180	08714729284307	18Fr Low Profile PEG Kit, Push	1.7	Each
M00563190	08714729284314	18Fr Low Profile PEG Kit, Push	2.4	Each
M00563010	08714729284147	24Fr Low Profile PEG Kit, Pull	1.7	Each
M00563020	08714729284154	24Fr Low Profile PEG Kit, Pull	2.4	Each
M00563030	08714729284161	24Fr Low Profile PEG Kit, Pull	3.4	Each
M00563040	08714729268833	24Fr Low Profile PEG Kit, Pull	4.4	Each
M00563080	08714729284208	24Fr Low Profile PEG Kit, Push	2.4	Each
M00563090	08714729284215	24Fr Low Profile PEG Kit, Push	3.4	Each
M00563100	08714729284222	24Fr Low Profile PEG Kit, Push	4.4	Each

One-Step Button and Low Profile Replacement Button Accessories				
Order Number	GTIN	Description	Length (cm)	Packaging
M00580401	08714729166443	18Fr Decompression Tube	1.2	Box 10
M00580411	08714729166450	18Fr Decompression Tube	1.7	Box 10
M00580421	08714729166467	18Fr Decompression Tube	2.4	Box 10
M00580431	08714729166474	18Fr Decompression Tube	3.4	Box 10
M00580441	08714729166481	18Fr Decompression Tube	4.4	Box 10
M00580461	08714729166504	24Fr Decompression Tube	1.7	Box 10
M00580471	08714729166511	24Fr Decompression Tube	2.4	Box 10
M00580481	08714729339649	24Fr Decompression Tube	3.4	Box 10
M00580491	08714729166535	24Fr Decompression Tube	4.4	Box 10
M00580171	08714729142232	18Fr Bolus Feeding Set Straight Adapter	—	Box 10
M00580181	08714729142249	18Fr Right Angle Feeding Set Continuous Adapter	—	Box 10
M00580231	08714729166405	24Fr Bolus Feeding Set Straight Adapter	—	Box 10
M00580241	08714729166412	24Fr Right Angle Feeding Set Continuous Adapter	—	Box 10



EndoVive™ Enteral Access Devices


Replacement Gastrostomy Tubes



Replacement Gastrostomy Feeding Tubes, Straight Bolster



Order Number	GTIN	Description	Packaging
M00582040	00813939020243	12Fr Replacement G-Tube, Straight	Each
M00582050	00813939020250	14Fr Replacement G-Tube, Straight	Each
M00582060	00813939020267	16Fr Replacement G-Tube, Straight	Each
M00582070	00813939020274	18Fr Replacement G-Tube, Straight	Each
M00582080	00813939020281	20Fr Replacement G-Tube, Straight	Each
M00582860	00813939020298	22Fr Replacement G-Tube, Straight	Each
M00582090	00813939020304	24Fr Replacement G-Tube, Straight	Each
M00582100	00813939020311	28Fr Replacement G-Tube, Straight	Each



Replacement Gastrostomy Feeding Tubes, Right Angle Bolster




Order Number	GTIN	Description	Packaging
M00582110	00813939020007	12Fr Replacement G-Tube, Right Angle	Each
M00582120	00813939020014	14Fr Replacement G-Tube, Right Angle	Each
M00582130	00813939020021	16Fr Replacement G-Tube, Right Angle	Each
M00582140	00813939020038	18Fr Replacement G-Tube, Right Angle	Each
M00582150	00813939020045	20Fr Replacement G-Tube, Right Angle	Each
M00582870	00813939020052	22Fr Replacement G-Tube, Right Angle	Each
M00582160	00813939020069	24Fr Replacement G-Tube, Right Angle	Each
M00582170	00813939020076	28Fr Replacement G-Tube, Right Angle	Each



EndoVive™ Enteral Access Devices

Low Profile Replacement Buttons & Accessories



Low Profile Replacement Button



Order Number	GTIN	Description	Length (cm)	Packaging
M00562760	08714729284017	18Fr Low Profile Replacement Button	1.2	Each
M00562770	08714729142201	18Fr Low Profile Replacement Button	1.7	Each
M00562780	08714729284024	18Fr Low Profile Replacement Button	2.4	Each
M00562790	08714729284031	18Fr Low Profile Replacement Button	3.4	Each
M00562800	08714729284048	18Fr Low Profile Replacement Button	4.4	Each
M00562820	08714729284062	24Fr Low Profile Replacement Button	1.7	Each
M00562830	08714729142218	24Fr Low Profile Replacement Button	2.4	Each
M00562840	08714729284079	24Fr Low Profile Replacement Button	3.4	Each
M00562850	08714729284086	24Fr Low Profile Replacement Button	4.4	Each

One-Step Button and Low Profile Replacement Button Accessories




Order Number	GTIN	Description	Length (cm)	Packaging
M00580401	08714729166443	18Fr Decompression Tube	1.2	Box 10
M00580411	08714729166450	18Fr Decompression Tube	1.7	Box 10
M00580421	08714729166467	18Fr Decompression Tube	2.4	Box 10
M00580431	08714729166474	18Fr Decompression Tube	3.4	Box 10
M00580441	08714729166481	18Fr Decompression Tube	4.4	Box 10
M00580461	08714729166504	24Fr Decompression Tube	1.7	Box 10
M00580471	08714729166511	24Fr Decompression Tube	2.4	Box 10
M00580481	08714729339649	24Fr Decompression Tube	3.4	Box 10
M00580491	08714729166535	24Fr Decompression Tube	4.4	Box 10
M00580171	08714729142232	18Fr Bolus Feeding Set Straight Adapter	—	Box 10
M00580181	08714729142249	18Fr Right Angle Feeding Set Continuous Adapter	—	Box 10
M00580231	08714729166405	24Fr Bolus Feeding Set Straight Adapter	—	Box 10
M00580241	08714729166412	24Fr Right Angle Feeding Set Continuous Adapter	—	Box 10



EndoVive™ Enteral Access Devices

Jejunal Feeding Tube Devices




Through-The-PEG (TTP) Jejunal Feeding Tube

Order Number	GTIN	Description	Tip Style	Packaging
M00566320	08714729786313	8.5Fr TTP Jejunal Feeding Tube	Bent Tip	Each
M00566330	08714729789727	8.5Fr TTP Jejunal Feeding Tube	Pigtail Tip	Each
M00566340	08714729789666	12Fr TTP Jejunal Feeding Tube	Bent Tip	Each
M00566350	08714729786252	12Fr TTP Jejunal Feeding Tube	Pigtail Tip	Each

Notes and Future Product Updates



Resection & Retrieval



Resection Devices

EMR for Upper GI

➤ Captivator™ EMR Devices 75

Stiff Snares

➤ Captivator II Single-Use Snares..... 75

➤ Captivator Single-Use Snares..... 76

Speciality Snares

➤ Profile™ Single-Use Snares 77

➤ Single-Use Rotatable Snares..... 77

Snares

➤ Captiflex™ Single-Use Snares 78

➤ Sensation™ Single-Use Short Throw Snares..... 78

Active Cords

➤ Active Cords 79

Retrieval Devices

Retrieval Net

➤ RescueNet™ Retrieval Net..... 80

Retrieval Grasping Forceps

➤ Rescue Rat Tooth Grasping Forceps..... 80

➤ Rescue Rat Tooth Grasping Forceps - Pediatric..... 80

➤ Rescue Alligator Grasping Forceps, Long 80

➤ Rescue Alligator Grasping Forceps..... 80

➤ Rescue Rat Tooth / Alligator Grasping Forceps..... 80

Polyp Traps


➤ TrapEase Single Chamber Poly Trap..... 80

➤ TrapEase Quad Chamber Poly Trap 80



Resection Devices

EMR for Upper GI




Captivator™ EMR Device



Order Number	GTIN	Description	Packaging
M00561600	08714729842675	Captivator EMR Device for Standard Gastrosopes	Box 1
M00561610	08714729842699	Captivator EMR Device for Large Gastrosopes	Box 1
M00561601	08714729842682	Captivator EMR Device for Standard Gastrosopes	Box 5
M00561611	08714729842705	Captivator EMR Device for Large Gastrosopes	Box 5

Stiff Snares



Captivator II Single-Use Snares




Order Number	GTIN	Description	Loop Width (mm)	Sheath O.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Packaging
M00561221	08714729855903	10mm Rounded – Stiff	10	2.4	2.8	240	Box 10
M00561222	08714729861300	10mm Rounded – Stiff	10	2.4	2.8	240	Box 20
M00561223	08714729861317	10mm Rounded – Stiff	10	2.4	2.8	240	Box 40
M00561231	08714729855910	15mm Rounded – Stiff	15	2.4	2.8	240	Box 10
M00561232	08714729861331	15mm Rounded – Stiff	15	2.4	2.8	240	Box 20
M00561233	08714729861348	15mm Rounded – Stiff	15	2.4	2.8	240	Box 40
M00561241	08714729855927	20mm Rounded – Stiff	20	2.4	2.8	240	Box 10
M00561242	08714729861362	20mm Rounded – Stiff	20	2.4	2.8	240	Box 20
M00561243	08714729861379	20mm Rounded – Stiff	20	2.4	2.8	240	Box 40
M00561191	08714729855934	25mm Rounded – Stiff	25	2.4	2.8	240	Box 10
M00561192	08714729861270	25mm Rounded – Stiff	25	2.4	2.8	240	Box 20
M00561193	08714729861287	25mm Rounded – Stiff	25	2.4	2.8	240	Box 40
M00561291	08714729645801	Extra Large Rounded – Stiff	33	2.4	2.8	240	Box 10



Resection Devices


Stiff Snares (Continued)



Captivator™ Single-Use Snares									
Order Number	GTIN	Description	Loop Width (mm)	Sheath O.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Packaging		
M00562301	08714729019312	Small Oval – Stiff	13	2.4	2.8	240	Box 10		
M00562321	08714729019312	Medium Oval – Stiff	27	2.4	2.8	240	Box 10		
M00561311	08714729019251	Large Oval – Medium Stiff	30	2.4	2.8	240	Box 10		
M00562391	08714729019404	Large Oval (thin wire) – Flexible	30	2.4	2.8	240	Box 10		
M00562451	08714729071068	Small Hexagonal – Stiff	13	2.4	2.8	240	Box 10		
M00562341	08714729019350	Medium Hexagonal – Stiff	27	2.4	2.8	240	Box 10		
M00562371	08714729019381	Medium Crescent – Stiff	27	2.4	2.8	240	Box 10		

Resection Devices


Speciality Snares



Profile™ Single-Use Snares



Order Number	GTIN	Description	Loop Width (mm)	Sheath O.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Packaging
M00562531	08714729156628	Extra Small Oval – Flexible	11	1.9	2.0	240	Box 10
M00562551	08714729156659	Small Oval – Flexible	13	1.9	2.0	240	Box 10
M00562571	08714729156635	Medium Oval – Flexible	27	1.9	2.0	240	Box 10



Single-Use Rotatable Snares*




Order Number	GTIN	Description	Loop Width (mm)	Sheath O.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Packaging
M00561821	08714729283980	Small Oval – Medium Stiff	13	2.4	2.8	240	Box 10
M00561831	08714729283997	Medium Oval – Medium Stiff	20	2.4	2.8	240	Box 10

*With Star Catheter™ Technology.



Resection Devices


Snares



Captiflex™ Single-Use Snares



Order Number	GTIN	Description	Loop Width (mm)	Sheath O.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Packaging
M00562471	08714729019459	Extra Small Oval – Flexible	11	2.4	2.8	240	Box 10
M00562421	08714729019435	Small Oval – Flexible	13	2.4	2.8	240	Box 10
M00562422	08714729501657	Small Oval – Flexible	13	2.4	2.8	240	Box 40
M00562401	08714729019411	Medium Oval – Flexible	27	2.4	2.8	240	Box 10
M00562402	08714729501640	Medium Oval – Flexible	27	2.4	2.8	240	Box 40



Sensation™ Single-Use Short Throw Snares



Order Number	GTIN	Description	Loop Width (mm)	Sheath O.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Packaging
M00562691	08714729158158	Small Oval – Flexible	13	2.4	2.8	240	Box 5
M00562692	08714729268826	Small Oval – Flexible	13	2.4	2.8	240	Box 20
M00562693	08714729501671	Small Oval – Flexible	13	2.4	2.8	240	Box 40
M00560321	08714729283928	Small Oval – Medium Stiff	13	2.4	2.8	240	Box 5
M00562671	08714729158141	Medium Oval – Flexible	27	2.4	2.8	240	Box 5
M00562672	08714729268819	Medium Oval – Flexible	27	2.4	2.8	240	Box 20
M00562673	08714729501664	Medium Oval – Flexible	27	2.4	2.8	240	Box 40
M00560311	08714729283904	Medium Oval – Medium Stiff	27	2.4	2.8	240	Box 5
M00562651	08714729158110	Large Oval – Medium Stiff	30	2.4	2.8	240	Box 5
M00562652	08714729268802	Large Oval – Medium Stiff	30	2.4	2.8	240	Box 20
M00560331	08714729283935	Medium Crescent – Medium Stiff	27	2.4	2.8	240	Box 5



Resection Devices


Active Cords

Active Cords			
Order Number	GTIN	Description	Packaging
M00561250	08714729019060	Active Cord, Fits Olympus™ Brand Generator	Each
M00561270	08714729053361	Active Cord, Fits Endostat™, Endostat II and Endostat III Generators, Bovie™, Valleylab™, Aspen Labs™, and ERBE™ Brand Generators	Each




Retrieval Devices

Retrieval Net



RescueNet™ Retrieval Net							
Order Number	GTIN	Description	Net Size (mm)	Sheath O.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Packaging
DGN-538-5	10840253107425	RescueNet Retrieval Net	30	2.5	2.8	230	Box 5

Retrieval Grasping Forceps



Rescue™ Retrieval Single-Use Graspers							
Order Number	GTIN	Description	Jaw Opening Width (mm)	Sheath O.D. (mm)	Minimum Working Channel (mm)	Working Length (cm)	Packaging
DGR-276-5	00816849017334	Rescue Rat Tooth Grasping Forceps	7	2.4	2.8	230	Box 5
DGR-277-5	00816849017341	Rescue Rat Tooth Grasping Forceps - Pediatric	4.5	1.8	2.0	180	Box 5
DGA-275-5	00816849017310	Rescue Alligator Grasping Forceps, Long	9	2.4	2.8	230	Box 5
DGA-274-5	00816849017303	Rescue Alligator Grasping Forceps	8	2.4	2.8	230	Box 5
DGC-383-5	00816849017327	Rescue Rat Tooth / Alligator Grasping Forceps	8	2.4	2.8	230	Box 5


Polyp Traps



TrapEase™ Polyp Traps							
Order Number	GTIN	Description					
STE-297-50	00816849010182	TrapEase Single Chamber Polyp Trap					
STE-298-25	00816849010199	TrapEase Quad Chamber Polyp Trap					



Infection Prevention



Infection Prevention

Cleaning Brushes

- HedgeHog™ Valve Brush..... 81
- HedgeHog Double Bundle (Double-End Brush + Valve Brush) .. 81
- HedgeHog Double-End Channel Brush 81
- HedgeHog Dual-End Brush Channel / Valve..... 81
- HedgeHog Sweeper Channel Brush 81
- HedgeHog Single-End Channel Brush 81

Detergent

- Pure™ Enzymatic Detergent - Spearmint Fresh Scent 81
- Pure Enzymatic Detergent - Fragrance Free 81

Bite Blocks

- Blox™ Endoscopic Bite Block with Elastic Cloth Strap (38 FR) .. 81
- Blox Endoscopic Bite Block with Elastic Cloth Strap (54 FR) 81
- Blox Endoscopic Bite Block with Elastic Cloth Strap (54 FR)
(w/o Dental Retention)..... 81
- Blox Endoscopic Bite Block with Elastic Cloth Strap (60 FR) 81

Lubricant

- Endoglide™+ Sterile Bacteriostatic Lubricating Gel..... 82

Procedure Gowns

- EndoArmor™ Procedure Gown..... 82
- EndoArmor+ Procedure Gown..... 82

Irrigation & Insufflation

- Hydra™ Disposable Auxiliary Water Jet Connector..... 82
- Hydra Water Bottle Cap..... 82
- Hydra Irrigation Tubing..... 82

Scope Transport

- CinchPad™ Endoscope Transport System..... 83

Sponges

- Koala™ Cleaning Sponge..... 83
- Koala Cleaning Sponge plus Pure enzymatic 83
- SafeStart™ Cleaning Kit 83
- SafeStart CP Cleaning Kit..... 83

Valves

- Orca™ Sterile Single-Use Air / Water & Suction Valves 83

Compliance EndoKit™

- Compliance EndoKit - Colon Kit..... 83
- Compliance EndoKit - EGD Kit 83



Infection Prevention

Cleaning Brushes

HedgeHog™ Cleaning Brushes



Order Number	GTIN	Description	Packaging
SBB-119-50	00816849010274	HedgeHog Valve Brush	Box 50
SBB-382-50	00816849016900	HedgeHog Double Bundle (Double-End Brush + Valve Brush)	Box 50
SBD-228-50	00816849010083	HedgeHog Double-End Channel Brush	Box 50
SBD-289-50	00816849010090	HedgeHog Dual-End Channel / Valve Brush	Box 50
SBD-291-50	00816849016443	HedgeHog Sweeper Channel Brush	Box 50
SBS-227-50	00816849010076	HedgeHog Single-End Channel Brush	Box 50

Detergent

Pure™ Enzymatic Detergents

Order Number	GTIN	Description	Packaging
SEE-220-4	10840253105308	Pure Enzymatic Detergent - Spearmint Fresh Scent	Box 4
SEE-573-4	10840253105292	Pure Enzymatic Detergent - Fragrance Free	Box 4

Bite Blocks



Blox™ Bite Blocks with Adjustable Straps



Order Number	GTIN	Description	Packaging
SBP-236-40	00816849010175	Blox Endoscopic Bite Block with Elastic Cloth Strap (38 FR)	Box 40
SBT-114-100	00816849010168	Blox Endoscopic Bite Block with Elastic Cloth Strap (54 FR)	Box 100
SBT-321-100	00816849016764	Blox Endoscopic Bite Block with Elastic Cloth Strap (54 FR) (w/o Dental Retention).	Box 100
SBT-546-100	00840253104809	Blox Endoscopic Bite Block with Elastic Cloth Strap (60 FR)	Box 100



Infection Prevention

Lubricant

Endoglide™ + Lubricant



Order Number	GTIN	Description	Packaging
SLT-394-100	00816849017167	Endoglide+ Sterile Bacteriostatic Lubricating Gel (32 g / 1.1 oz)	Box 100
SLT-395-50	00840253103055	Endoglide+ Sterile Bacteriostatic Lubricating Gel (56.7 g / 2 oz)	Box 50


Procedure Gowns

EndoArmor™ Procedure Gown



Order Number	GTIN	Description	Packaging
SGE-368-100	00816849010281	EndoArmor Procedure Gown	Box 100
SGE-389-100	00816849016481	EndoArmor+ Procedure Gown	Box 100

Irrigation & Insufflation



Hydra™ Irrigation Portfolio

Order Number	GTIN	Description	Scope Compatibility	CO2	Air/Water Tubing Length (in)	CO2 Tubing Length (in)	Packaging
SIT-361-100	00816849016955	Hydra Disposable Auxiliary Water Jet Connector	Olympus	N/A	N/A	N/A	Box 100
SIT-470-100	00816849019895	Hydra Disposable Auxiliary Water Jet Connector	Fuse	N/A	N/A	N/A	Box 100
SCT-350-10	00816849016887	Hydra Water Bottle Cap	Olympus	No	16	N/A	Box 10
SCT-351-10	00816849016894	Hydra Water Bottle Cap	Olympus	Yes	16	18	Box 10
SCT-466-10	00816849019864	Hydra Water Bottle Cap	Fuse	No	16	N/A	Box 10
SCT-467-10	00816849019833	Hydra Water Bottle Cap	Fuse	Yes	16	18	Box 10
SIT-355-15	00840253112200	Hydra Irrigation Tubing (compatible with ERBE, Medivators pumps)	N/A	N/A	N/A	N/A	Box 15
SIT-576-15	00840253108401	Hydra Irrigation Tubing (compatible with Olympus pumps)	N/A	N/A	N/A	N/A	Box 15



Section Table of Contents

Prev Page

82

Next Page

Main Table of Contents



Infection Prevention

Scope Transport

CinchPad™ Endoscope Transport System



Order Number	GTIN	Description	Packaging
SPX-017-25	00816849010045	CinchPad Endoscope Transport System	Box 25

Sponges

Koala Sponge



SafeStart Cleaning Kit



Koala™ Sponges



Order Number	GTIN	Description	Packaging
SSB-210-50	00816849015613	Koala Cleaning Sponge	Box 50
SES-240-50	00816849015569	Koala Cleaning Sponge plus Pure enzymatic	Box 50
SIK-242-24	00816849015637	SafeStart™ Cleaning Kit	Box 24
SIK-243-24	00816849015620	SafeStart CP Cleaning Kit	Box 24

Valves



Orca™ Sterile Single-Use Air / Water & Suction Valves

Orca™ Sterile Single-Use Air / Water & Suction Valves



Order Number	GTIN	Description	Scope Compatibility	Packaging
SUV-617-50	00840253111012	Orca Sterile Single-Use Air / Water & Suction Valves	Olympus	Box 50

Compliance EndoKit™

Compliance EndoKit



Order Number	GTIN	Description	Packaging
CEK-721-20	00840253113337	Compliance EndoKit - Colon Kit	Box 20
CEK-733-20	00840253113450	Compliance EndoKit - EGD Kit	Box 20

Please consult your Boston Scientific Territory Manager on options to customize your Compliance EndoKit™ to meet your individual endoscope reprocessing needs.



Section Table of Contents

Prev Page

83

Next Page

Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Notes and Future Product Updates



Pulmonary Endoscopy

Alair™ Bronchial Thermoplasty System

- Alair Bronchial Thermoplasty System 85

Airway Stents

- Dynamic™ (Y) Stent..... 85
- Polyflex™ Single-Use Self-Expanding Silicone Airway Stent 88
- Ultraflex™ Single-Use Partially Covered Tracheobronchial Stent System – Distal Release 87
- Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Distal Release..... 86
- Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Proximal Release 86

Biopsy Forceps

- Radial Jaw™ 4 Single-Use Pulmonary Biopsy Forceps..... 89

Endobronchial Ultrasound Transbronchial Aspiration Needle (EBUS)

- Expect™ Pulmonary Endobronchial Ultrasound Transbronchial Aspiration Needle / Adaptor 89

Transbronchial Aspiration Needles

- eXcelon™ Transbronchial Aspiration Needle 89

Cytology Brushes

- Cellebri™ Single-Use Endoscopic Cytology Brush..... 90

Airway Retrieval Baskets

- Zero Tip™ Airway Retrieval Basket 90

Controlled Radial Expansion Pulmonary Balloon Dilatation Catheter

- CRE™ Pulmonary Balloon Dilatation Catheter..... 91

Pulmonary Guidewires


- Amplatz Super Stiff Guidewire 91
- Pulmonary Jagwire™ Pulmonary Guidewire 91

Inflation System

- Alliance™ II Integrated Inflation / Lithotripsy Device..... 91
- Alliance II Single-Use Syringe / Gauge Assembly 91

Pulmonary

Alair™ Bronchial Thermoplasty System



Alair Bronchial Thermoplasty System*


Order Number	GTIN	Model Number	Description	Packaging
M005ATS20000	08714729796695	Alair ATS 200	Alair Bronchial Thermoplasty Controller	Each
M005ATS25010	08714729802792	Alair ATS 2-5	⊗ Alair Bronchial Thermoplasty Catheter**	Each
M005ATS20110	08714729802778	Alair ATS 201	Alair RF Controller Accessory Kit	Each

*Service Agreement available with 1, 2, 3, and 4 year terms.

**Note: Initial stocking order requires a minimum order of 6 catheters (covering the complete treatment of 2 patients).

UPN	Description
SRVM005ATS1001	ENDO POS 1 Year Service Alair
SRVM005ATS1002	ENDO POS 2 Year Service Alair
SRVM005ATS1003	ENDO POS 3 Year Service Alair
SRVM005ATS1004	ENDO POS 4 Year Service Alair

Airway Stents



Dynamic™ (Y) Stent



Order Number	GTIN	Tracheal Width (ID) (mm)	Bronchial Width (ID) (mm)	Tracheal Length (mm)	Bronchial Lengths (R / L) (mm)	Packaging
M00570670	08714729764809	11	8	110	25 / 40	Each
M00570680	08714729764816	13	10	110	25 / 40	Each
M00570690	08714729764823	15	12	110	25 / 40	Each



Section Table of Contents

Prev Page

85


Next Page

Main Table of Contents



Pulmonary

Airway Stents (Continued)



 Magnetic Resonance Conditional¹

Ultraflex™ Single-Use Uncovered Tracheobronchial Stent System – Distal Release*



Order Number	GTIN	Description	Expanded Stent O.D. (mm)	Expanded Stent Length (mm)	Compressed Stent O.D. (Fr)**	Tip Max O.D. (mm)	Packaging
M00576300	08714729837688	Uncovered Distal Release	10	30	15	4.1	Each
M00576310	08714729837695	Uncovered Distal Release	12	30	17	4.1	Each
M00576320	08714729837701	Uncovered Distal Release	14	30	18	4.1	Each
M00576330	08714729837718	Uncovered Distal Release	16	40	20	5.3	Each
M00576340	08714729837725	Uncovered Distal Release	16	60	20	5.3	Each
M00576350	08714729837732	Uncovered Distal Release	18	40	20	5.3	Each
M00576360	08714729837749	Uncovered Distal Release	18	60	20	5.3	Each
M00576370	08714729837756	Uncovered Distal Release	20	40	20	5.3	Each
M00576380	08714729837763	Uncovered Distal Release	20	60	20	5.3	Each

*Recommended Guidewire .035" Pulmonary Jagwire™ Pulmonary Guidewire and Amplatz Super Stiff Guidewire, Order #1517,  (See page 91).
Delivery System working length is 95cm.

**Data on File at Boston Scientific.



 Magnetic Resonance Conditional¹

Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Proximal Release*



Order Number	GTIN	Description	Expanded Stent O.D. (mm)	Expanded Stent Length (mm)	Compressed Stent O.D. (Fr)**	Tip Max O.D. (mm)	Packaging
M00576410	08714729837770	Uncovered Proximal Release	8	20	15	4.1	Each
M00576420	08714729837787	Uncovered Proximal Release	8	40	15	4.1	Each
M00576430	08714729837794	Uncovered Proximal Release	10	20	15	4.1	Each
M00576440	08714729837800	Uncovered Proximal Release	10	40	15	4.1	Each
M00576450	08714729837817	Uncovered Proximal Release	12	20	17	4.1	Each
M00576460	08714729837824	Uncovered Proximal Release	12	40	17	4.1	Each
M00576470	08714729837831	Uncovered Proximal Release	14	20	18	4.1	Each
M00576480	08714729837848	Uncovered Proximal Release	14	40	18	4.1	Each
M00576490	08714729837855	Uncovered Proximal Release	14	60	18	4.1	Each

*Recommended Guidewire .035" Pulmonary Jagwire Pulmonary Guidewire and Amplatz Super Stiff Guidewire, Order #1517,  (See page 91).
Delivery System working length is 95cm.

**Data on File at Boston Scientific.



Section Table of Contents


Prev Page  86  Next Page

Main Table of Contents 

For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Pulmonary

Airway Stents (Continued)



 Magnetic Resonance
Conditional[†]

Ultraflex™ Single-Use Partially Covered Tracheobronchial Stent System – Distal Release*



Order Number	GTIN	Description	Expanded Stent O.D. (mm)	Expanded Stent Length (mm)	Covered Compressed Stent O.D. (Fr)**	Cover Length (mm)	Tip Max O.D. (mm)	Packaging
M00576520	08714729842323	Partially Covered Distal Release	8	40	16	25	4.1	Each
M00576530	08714729842330	Partially Covered Distal Release	10	30	18	15	4.1	Each
M00576540	08714729842347	Partially Covered Distal Release	10	40	18	25	4.1	Each
M00576550	08714729842354	Partially Covered Distal Release	12	30	20	15	4.1	Each
M00576560	08714729842361	Partially Covered Distal Release	12	40	20	25	4.1	Each
M00576570	08714729842378	Partially Covered Distal Release	14	30	24	15	4.1	Each
M00576580	08714729842385	Partially Covered Distal Release	14	40	24	25	4.1	Each
M00576590	08714729842392	Partially Covered Distal Release	14	60	24	45	4.1	Each
M00576600	08714729842408	Partially Covered Distal Release	14	80	24	65	4.1	Each
M00576610	08714729842415	Partially Covered Distal Release	16	40	25	25	5.3	Each
M00576620	08714729842422	Partially Covered Distal Release	16	60	25	45	5.3	Each
M00576630	08714729842439	Partially Covered Distal Release	16	80	25	65	5.3	Each
M00576640	08714729842446	Partially Covered Distal Release	18	40	27	25	5.3	Each
M00576650	08714729842453	Partially Covered Distal Release	18	60	27	45	5.3	Each
M00576660	08714729842460	Partially Covered Distal Release	18	80	27	65	5.3	Each
M00576670	08714729842477	Partially Covered Distal Release	20	40	27	25	5.3	Each
M00576680	08714729842484	Partially Covered Distal Release	20	60	27	45	5.3	Each
M00576690	08714729842491	Partially Covered Distal Release	20	80	27	65	5.3	Each

*Recommended Guidewire .035" Pulmonary Jagwire™ Pulmonary Guidewire and Amplatz Super Stiff Guidewire, Order #1517,  (See page 91).
Delivery System working length is 95cm.

**Data on File at Boston Scientific.



Section Table of Contents

Prev Page  87  Next Page


Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Pulmonary

Airway Stents (Continued)




Polyflex™ Single-Use Self-Expanding Silicone Airway Stent					
Order Number	GTIN	Stent ID (mm)	Length (mm)	Delivery System Diameter (mm)	Packaging
M00570000	08714729457138	8	20	7	Each
M00570010	08714729457145	8	30	7	Each
M00570020	08714729457152	10	20	8	Each
M00570030	08714729457169	10	30	8	Each
M00570040	08714729457176	10	40	8	Each
M00570050	08714729457183	10	50	8	Each
M00570060	08714729457190	12	20	9	Each
M00570070	08714729457206	12	30	9	Each
M00570080	08714729457213	12	40	9	Each
M00570090	08714729457220	12	50	9	Each
M00570100	08714729457237	14	20	9	Each
M00570110	08714729457244	14	30	9	Each
M00570120	08714729457251	14	40	9	Each
M00570130	08714729457268	14	50	9	Each
M00570140	08714729457275	14	60	9	Each
M00570150	08714729457282	16	30	10	Each
M00570160	08714729457299	16	40	10	Each
M00570170	08714729457305	16	50	10	Each
M00570180	08714729457312	16	60	10	Each
M00570190	08714729457329	16	70	10	Each
M00570200	08714729457336	18	30	11	Each
M00570210	08714729457343	18	40	11	Each
M00570220	08714729457350	18	50	11	Each
M00570230	08714729457367	18	60	11	Each
M00570240	08714729457374	18	70	11	Each
M00570250	08714729457381	18	80	11	Each
M00570260	08714729457398	20	40	12	Each
M00570270	08714729457404	20	50	12	Each
M00570280	08714729457411	20	60	12	Each
M00570290	08714729457428	20	70	12	Each
M00570300	08714729457435	20	80	12	Each
M00570310	08714729457442	22	50	13	Each
M00570320	08714729457459	22	60	13	Each
M00570330	08714729457466	22	80	13	Each



Pulmonary

Biopsy Forceps



Radial Jaw™ 4 Single-Use Pulmonary Biopsy Forceps



Order Number	GTIN	Description	Differentiator	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Packaging
M00515181	08714729792864	Single-Use Biopsy Forceps	Pulmonary Standard Capacity	1.8	100	2.0	Box 5
M00515182	08714729799337	Single-Use Biopsy Forceps	Pulmonary Standard Capacity	1.8	100	2.0	Box 20
M00515191	08714729792871	Single-Use Biopsy Forceps	Pulmonary Standard Capacity w/ Needle	1.8	100	2.0	Box 5
M00515192	08714729818083	Single-Use Biopsy Forceps	Pulmonary Standard Capacity w/ Needle	1.8	100	2.0	Box 20
M00515201	08714729793281	Single-Use Biopsy Forceps	Pulmonary Large Capacity	2.2	100	2.8	Box 5
M00515202	08714729818076	Single-Use Biopsy Forceps	Pulmonary Large Capacity	2.2	100	2.8	Box 20

Expect™ Pulmonary Endobronchial Ultrasound Transbronchial Aspiration Needle




Expect Pulmonary Endobronchial Ultrasound Transbronchial Aspiration Needle



Order Number	GTIN	Needle Size	Minimum Working Channel (mm)	Sheath O.D. (mm)	Packaging
M00558220	08714729861409	Expect Pulmonary 22Ga – Olympus scope compatible	2.0	1.6	Please contact your Boston Scientific Sales Representative for ordering information
M00558250	08714729861416	Expect Pulmonary 25Ga – Olympus scope compatible	2.0	1.4	Please contact your Boston Scientific Sales Representative for ordering information
M00558731	08714729893189	Expect Pulmonary Needle Adaptor – Olympus scope compatible	N/A	N/A	Please contact your Boston Scientific Sales Representative for ordering information

Transbronchial Aspiration Needles



eXcelon™ Transbronchial Aspiration Needle*



Order Number	GTIN	Description	Needle Gauge	Needle Length (mm)	Catheter Length (cm)	Sheath O.D. (mm)	Packaging
M00564101	08714729456452	eXcelon Transbronchial Aspiration Needle	19	15	130	1.8	Box 5
M00564111	08714729456469	eXcelon Transbronchial Aspiration Needle	20	15	130	1.8	Box 5
M00564121	08714729456476	eXcelon Transbronchial Aspiration Needle	21	15	130	1.8	Box 5

*Packaging also includes a 20ml syringe.



Section Table of Contents

Prev Page ◀ 89 ▶ Next Page

Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Pulmonary


Cytology Brushes

Cellebrity™ Single-Use Endoscopic Cytology Brush



Order Number	GTIN	Required Working Channel (mm)	Bristle O.D. (mm)	Sheath Length (cm)	Packaging
M00516001	08714729138983	2.0	1.0	140	Box 10
M00516011	08714729138990	2.0	1.5	140	Box 10
M00516071	08714729139089	2.0	1.9	100	Box 10
M00516151	08714729014003	2.0	1.9	150	Box 10

Airway Retrieval Baskets



Zero Tip™ Airway Retrieval Basket



Order Number	GTIN	Description	Sheath O.D. (mm)	Working Length (cm)	Basket Opening (mm)	Packaging
M00513200	08714729414995	Zero Tip Airway Retrieval Basket	0.8	120	12	Each
M00513210	08714729415008	Zero Tip Airway Retrieval Basket	1.0	120	16	Each



Pulmonary

Controlled Radial Expansion Pulmonary Balloon Dilatation Catheter



CRE™ Pulmonary Balloon Dilatation Catheter*



Order Number	GTIN	Balloon Inflated O.D.		Balloon Length (cm)	Inflation Pressures		Catheter Length (cm)	Packaging
		(mm)	(Fr)		(atm)	(kPa)		
M00550300	08714729456186	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	110	Each
M00550310	08714729456193	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	110	Each
M00550320	08714729456209	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	110	Each
M00550330	08714729456216	8-9-10	24-27-30	3.0	3-5.5-9	304-557-912	110	Each
M00550340	08714729456223	10-11-12	30-33-36	3.0	3-5-8	304-507-811	110	Each
M00550350	08714729456230	12-13.5-15	36-40.5-45	3.0	3-4.5-8	304-456-811	110	Each

*Recommended Inflation Device: Alliance™ II Integrated Inflation / Lithotripsy Device Order #5062, Alliance II Single-Use Gauge / Syringe Assembly Order #5060 (see below).


Pulmonary Guidewires

Pulmonary Jagwire™ Pulmonary Guidewire



Order Number	GTIN	Description	O.D. (in)	Length (cm)	Packaging
M00515171	08714729455813	Pulmonary Jagwire Pulmonary Guidewire	.035	180	Box 2
M00550090	08714729301929	Amplatz Super Stiff Guidewire	.038	260	Each

Inflation System



Alliance™ II Integrated Inflation / Lithotripsy Device and Alliance II Single-Use Syringe / Gauge Assembly

Order Number	GTIN	Description	Packaging
M00550620	08714729283881	Alliance II Integrated Inflation / Lithotripsy Device	Each
M00550601	08714729129332	Alliance II Single-Use Syringe / Gauge Assembly	Box 5



Section Table of Contents

Prev Page **91** Next Page

Main Table of Contents




For pricing information contact your Boston Scientific Sales Representative or Customer Service at **1.888.272.1001**

Notes and Future Product Updates



Solutions



Solutions

- **GuidePoint Reimbursement Solutions**.. 93
- **ADVANTICS™ Innovative Healthcare Solutions**..... 94
- **GI Specialized Pathology**..... 94
- **Scope Repair**..... 94
- **GI Lab Optimization**..... 94

Solutions

GuidePoint Reimbursement Resources and Solutions

The Endoscopy Health Economics and Market Access (HEMA) Team offers a suite of reimbursement services and support for healthcare professionals seeking a better understanding of reimbursement for Boston Scientific's endoscopy and pulmonary products and procedures. We can help you develop a base level of knowledge in coding, coverage and payment for these procedures and offer educational presentations and webinars on changes in healthcare policies and Medicare reimbursement and changes to Medicare reimbursement impacting providers and facilities.

About Our Team

Boston Scientific Endoscopy's HEMA Team is comprised of highly skilled professionals specializing in reimbursement, health economics and health policy issues. The team has a broad background with experience working for public and private payers, medical associations, healthcare consulting agencies and academic research organizations. We use clinical and economic outcomes data and collaborate with outside experts to promote appropriate coding, coverage and reimbursement for Boston Scientific endoscopy and pulmonary products and procedures.

Contact Information

For more information about Boston Scientific Endoscopy's GuidePoint Reimbursement Resources and Solutions, please contact our **Health Economics and Reimbursement department at 800-876-9960 x54510 or ENDOREIMBURSEMENT@bsci.com**.

www.bostonscientific.com/reimbursement

GUIDEPOINT
Reimbursement Resources

Health economic and reimbursement information provided by Boston Scientific Corporation is gathered from third-party sources and is subject to change without notice as a result of complex and frequently changing laws, regulations, rules and policies. This information is presented for illustrative purposes only and does not constitute reimbursement or legal advice. Boston Scientific encourages providers to submit accurate and appropriate claims for services. It is always the provider's responsibility to determine medical necessity, the proper site for delivery of any services and to submit appropriate codes, charges, and modifiers for services that are rendered. Boston Scientific recommends that you consult with your payers, reimbursement specialists and/or legal counsel regarding coding, coverage and reimbursement matters. Boston Scientific does not promote the use of its products outside their FDA-approved label.



Solutions

Boston Scientific's Endoscopy business develops less invasive medical devices for treating gastrointestinal diseases. We are committed to providing clinically proven technologies, services and customer support to help healthcare providers improve patient outcomes while reducing the total cost of patient care.

ADVANTICS™ Innovative Healthcare Solutions

ADVANTICS™ Innovative Healthcare Solutions is a new way to collaborate to meet your specific challenges in the changing healthcare landscape. ADVANTICS recognizes that every provider is unique and that only by working together can we find truly effective solutions.



GI Specialized Pathology

Diagnostics for GI Diseases. Getting an accurate and timely diagnosis is critical to patient care. GI SPECIALIZED PATHOLOGY services include GI subspecialty trained pathologists, GI trained histotechnicians to help ensure proper slide oriented, and customer service teams dedicated to improving the diagnosis process for physicians and their patients. Today, of the top diagnostic labs in the U.S., ours is the only pathology lab focused solely on GI diseases.



Scope Repair

Capital Equipment Management. As an alternative to original equipment manufacturer (OEM) and costly service contracts, SCOPE REPAIR services are offered through a partnership with Total Scope, Inc. Established in 1993, Total Scope is focused on providing excellent customer service, including scope maintenance education and training for its customers. The company is certified to the quality standard ISO 13485, for the service, repair and maintenance of medical instrumentation. Total Scope is a woman-owned business.

GI Lab Optimization

Endoscopy Lab and Continuous Improvement. GI Lab Optimization services are designed to address specific customer issues related to quality, cost, safety, and staff and patient satisfaction. A mix of Lean Six Sigma, project management and Human Centered Design expertise can be used to tailor solutions to your needs for endoscopy lab optimization.

To learn more visit, [bostonscientific.com/endoscopysolutions](https://www.bostonscientific.com/endoscopysolutions)




Future Product Updates



Notes and Future Product Updates



Product Code Index



Product Code Index

➤ Begins on Page..... 97

Product Code Index

Product Code	Description	Page #	Product Code	Description	Page #
CEK-721-20.....	Compliance EndoKit™ - Colon Kit	83	M00501271	EndoVive Initial Placement PEG Accessories – Y-Ports with ENFit.....	65
CEK-733-20.....	Compliance EndoKit - EGD Kit.....	83	M00509001	EndoVive Safety PEG Kit, Pull with ENFit.....	65
DGA-274-5.....	Rescue Alligator Short Grasping Forceps.....	80	M00509011	EndoVive Safety PEG Kit, Push with ENFit.....	65
DGA-275-5.....	Rescue™ Alligator Long Grasping Forceps.....	80	M00509021	EndoVive Safety PEG Kit, Pull with ENFit.....	65
DGC-383-5.....	Rescue Combo Rat / Alligator Long Grasping Forceps.....	80	M00509031	EndoVive Safety PEG Kit, Push with ENFit.....	65
DGN-538-5	RescueNet™ Retrieval Net.....	80	M00509041	EndoVive Standard PEG Kit, Pull with ENFit.....	65
DGR-276-5.....	Rescue Rat Tooth Grasping Forceps	80	M00509051	EndoVive Standard PEG Kit, Push with ENFit.....	65
DGR-277-5.....	Rescue Rat Tooth PED Grasping Forceps.....	80	M00509061	EndoVive Standard PEG Kit, Pull with ENFit.....	65
H965180011	WALLSTENT™ Super Stiff Guidewire.....	44	M00509071	EndoVive Standard PEG Kit, Push with ENFit.....	65
H965430100	WALLSTENT Endoscopic Biliary Endoprosthesis Stents	18	M00509221	Distal End ENFit Transition Connector.....	68
H965430200	WALLSTENT Endoscopic Biliary Endoprosthesis Stents	18	M00509231	ENFit Medication Port Transition Connector.....	68
H965430300	WALLSTENT Endoscopic Biliary Endoprosthesis Stents	18	M00509241	ENFit Feeding Port Transition Connector.....	68
H965430400	WALLSTENT Endoscopic Biliary Endoprosthesis Stents	18	M00509580	EndoVive One-Step Button Low Profile PEG Kit, Pull with ENFit.....	66
H965430500	WALLSTENT Endoscopic Biliary Endoprosthesis Stents	18	M00509590	EndoVive One-Step Button Low Profile PEG Kit, Pull with ENFit.....	66
H965430600	WALLSTENT Endoscopic Biliary Endoprosthesis Stents	18	M00509610	EndoVive One-Step Button Low Profile PEG Kit, Pull with ENFit.....	66
H965430700	WALLSTENT Endoscopic Biliary Endoprosthesis Stents	18	M00509620	EndoVive One-Step Button Low Profile PEG Kit, Push with ENFit.....	66
H965430800	WALLSTENT Endoscopic Biliary Endoprosthesis Stents	18	M00509630	EndoVive One-Step Button Low Profile PEG Kit, Push with ENFit.....	66
H965430900	WALLSTENT Endoscopic Biliary Endoprosthesis Stents	18	M00509680	EndoVive One-Step Button Low Profile PEG Kit, Push with ENFit.....	66
H965431000	WALLSTENT Endoscopic Biliary Endoprosthesis Stents	18	M00509690	EndoVive One-Step Button™ Low Profile PEG Kit, Pull with ENFit.....	66
H965431100	WALLSTENT Endoscopic Biliary Endoprosthesis Stents	18	M00509700	EndoVive One-Step Button Low Profile PEG Kit, Pull with ENFit.....	66
H965431200	WALLSTENT Endoscopic Biliary Endoprosthesis Stents	18	M00509710	EndoVive One-Step Button Low Profile PEG Kit, Pull with ENFit.....	66
M00501061	EndoVive™ Low Profile Replacement Button Accessories with ENFit – 18Fr Right Angle Feeding Set, Continuous Adapter	66, 68	M00509720	EndoVive One-Step Button Low Profile PEG Kit, Pull with ENFit.....	66
M00501071	EndoVive Low Profile Replacement Button Accessories with ENFit – 24Fr Right Angle Feeding Set, Continuous Adapter.....	66, 68	M00509730	EndoVive One-Step Button Low Profile PEG Kit, Pull with ENFit.....	66
M00501091	EndoVive Low Profile Replacement Button Accessories with ENFit™ – 18Fr Bolus Feeding Set, Straight Adapter..	66, 68	M00509740	EndoVive One-Step Button Low Profile PEG Kit, Push with ENFit.....	66
M00501101	EndoVive Low Profile Replacement Button Accessories with ENFit – 24Fr Bolus Feeding Set, Straight Adapter.....	66, 68	M00509750	EndoVive One-Step Button Low Profile PEG Kit, Push with ENFit.....	66
M00501121	EndoVive Low Profile Replacement Button Accessories with ENFit – 18Fr Decompression Tube.....	66, 68	M00509780	EndoVive One-Step Button Low Profile PEG Kit, Pull with ENFit.....	66
M00501131	EndoVive Low Profile Replacement Button Accessories with ENFit – 18Fr Decompression Tube.....	66, 68	M00509830	EndoVive Low Profile Replacement Button with ENFit.....	67
M00501141	EndoVive Low Profile Replacement Button Accessories with ENFit – 18Fr Decompression Tube.....	66, 68	M00509840	EndoVive Low Profile Replacement Button with ENFit.....	67
M00501151	EndoVive Low Profile Replacement Button Accessories with ENFit – 18Fr Decompression Tube.....	66, 68	M00509850	EndoVive Low Profile Replacement Button with ENFit.....	67
M00501161	EndoVive Low Profile Replacement Button Accessories with ENFit – 18Fr Decompression Tube.....	66, 68	M00509860	EndoVive Low Profile Replacement Button with ENFit.....	67
M00501181	EndoVive Low Profile Replacement Button Accessories with ENFit – 24Fr Decompression Tube.....	66, 68	M00509870	EndoVive Low Profile Replacement Button with ENFit.....	67
M00501191	EndoVive Low Profile Replacement Button Accessories with ENFit – 24Fr Decompression Tube.....	66, 68	M00509410	EndoVive Low Profile Replacement Button with ENFit.....	67
M00501201	EndoVive Low Profile Replacement Button Accessories with ENFit – 24Fr Decompression Tube.....	66, 68	M00509420	EndoVive Low Profile Replacement Button with ENFit.....	67
M00501211	EndoVive Low Profile Replacement Button Accessories with ENFit – 24Fr Decompression Tube.....	66, 68	M00509430	EndoVive Low Profile Replacement Button with ENFit.....	67
M00501261	EndoVive Initial Placement PEG Accessories – Y-Ports with ENFit.....	65	M00509440	EndoVive Low Profile Replacement Button with ENFit.....	67



Product Code Index

Product Code	Description	Page #	Product Code	Description	Page #
M00509881	EndoVive™ Replacement Gastrostomy Feeding Tube, with ENFit, Straight	67	M00513332	Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps w/ Needle	61
M00509891	EndoVive Replacement Gastrostomy Feeding Tube, with ENFit, Straight	67	M00513333	Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps w/ Needle	61
M00509901	EndoVive Replacement Gastrostomy Feeding Tube, with ENFit, Straight	67	M00513343	Radial Jaw 4 Single-Use Jumbo Biopsy Forceps	61
M00509911	EndoVive Replacement Gastrostomy Feeding Tube, with ENFit, Straight	67	M00513353	Radial Jaw 4 Single-Use Jumbo Biopsy Forceps w/ Needle	61
M00509921	EndoVive Replacement Gastrostomy Feeding Tube, with ENFit, Straight	67	M00513361	Radial Jaw 4 Single-Use Jumbo Biopsy Forceps	61
M00509931	EndoVive Replacement Gastrostomy Feeding Tube, with ENFit, Straight	67	M00513362	Radial Jaw 4 Single-Use Jumbo Biopsy Forceps	61
M00509941	EndoVive Replacement Gastrostomy Feeding Tube, with ENFit, Straight	67	M00513363	Radial Jaw 4 Single-Use Jumbo Biopsy Forceps	61
M00509961	EndoVive Replacement Gastrostomy Feeding Tube, with ENFit, Right Angle	67	M00513371	Radial Jaw 4 Single-Use Jumbo Biopsy Forceps w/ Needle	61
M00509971	EndoVive Replacement Gastrostomy Feeding Tube, with ENFit, Right Angle	67	M00513372	Radial Jaw 4 Single-Use Jumbo Biopsy Forceps w/ Needle	61
M00509981	EndoVive Replacement Gastrostomy Feeding Tube, with ENFit, Right Angle	67	M00513373	Radial Jaw 4 Single-Use Jumbo Biopsy Forceps w/ Needle	61
M00509991	EndoVive Replacement Gastrostomy Feeding Tube, with ENFit, Right Angle	67	M00513381	Radial Jaw 4 Standard Capacity Single-Use Biopsy Forceps	61
M00510001	EndoVive Replacement Gastrostomy Feeding Tube, with ENFit, Right Angle	67	M00513382	Radial Jaw 4 Standard Capacity Single-Use Biopsy Forceps	61
M00510011	EndoVive Replacement Gastrostomy Feeding Tube, with ENFit, Right Angle	67	M00513383	Radial Jaw 4 Standard Capacity Single-Use Biopsy Forceps	61
M00510021	EndoVive Replacement Gastrostomy Feeding Tube, with ENFit, Right Angle	67	M00513391	Radial Jaw 4 Standard Capacity Single-Use Biopsy Forceps w/ Needle	61
M00510101	Multibite™ Single-Use Multiple Sample Biopsy Forceps	63	M00513392	Radial Jaw 4 Standard Capacity Single-Use Biopsy Forceps w/ Needle	61
M00510121	Multibite Single-Use Multiple Sample Biopsy Forceps	63	M00513393	Radial Jaw 4 Standard Capacity Single-Use Biopsy Forceps w/ Needle	61
M00510860	Trapezoid™ RX Wireguided Retrieval Basket	32	M00513401	Radial Jaw 4 Standard Capacity Single-Use Biopsy Forceps	61
M00510870	Trapezoid RX Wireguided Retrieval Basket	32	M00513402	Radial Jaw 4 Standard Capacity Single-Use Biopsy Forceps	61
M00510880	Trapezoid RX Wireguided Retrieval Basket	32	M00513403	Radial Jaw 4 Standard Capacity Single-Use Biopsy Forceps	61
M00510890	Trapezoid RX Wireguided Retrieval Basket	32	M00513411	Radial Jaw 4 Standard Capacity Single-Use Biopsy Forceps w/ Needle	61
M00511731	Optiflo™ Single-Use Hemostasis Catheters	51	M00513412	Radial Jaw 4 Standard Capacity Single-Use Biopsy Forceps w/ Needle	61
M00511751	Optiflo Single-Use Hemostasis Catheters	51	M00513413	Radial Jaw 4 Standard Capacity Single-Use Biopsy Forceps w/ Needle	61
M00513200	Zero Tip™ Airway Retrieval Basket	90	M00513443	Radial Jaw 4 Single-Use Pediatric Biopsy Forceps	62
M00513210	Zero Tip Airway Retrieval Basket	90	M00513453	Radial Jaw 4 Single-Use Pediatric Biopsy Forceps w/ Needle	62
M00513301	Radial Jaw™ 4 Single-Use Large Capacity Biopsy Forceps	61	M00513700	Ultraflex™ Single-Use Uncovered Esophageal NG Stent System – Distal Release	47
M00513302	Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps	61	M00513710	Ultraflex Single-Use Uncovered Esophageal NG Stent System – Distal Release	47
M00513303	Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps	61	M00513720	Ultraflex Single-Use Uncovered Esophageal NG Stent System – Distal Release	47
M00513311	Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps w/ Needle	61	M00513730	Ultraflex Single-Use Covered Esophageal NG Stent System – Distal Release	46
M00513312	Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps w/ Needle	61	M00513740	Ultraflex Single-Use Covered Esophageal NG Stent System – Distal Release	46
M00513313	Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps w/ Needle	61	M00513750	Ultraflex Single-Use Covered Esophageal NG Stent System – Distal Release	46
M00513321	Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps	61	M00513800	Ultraflex Single-Use Uncovered Esophageal NG Stent System – Proximal Release	47
M00513322	Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps	61	M00513810	Ultraflex Single-Use Uncovered Esophageal NG Stent System – Proximal Release	47
M00513323	Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps	61	M00513830	Ultraflex Single-Use Uncovered Esophageal NG Stent System – Proximal Release	47
M00513331	Radial Jaw 4 Single-Use Large Capacity Biopsy Forceps w/ Needle	61	M00513840	Ultraflex Single-Use Covered Esophageal NG Stent System – Proximal Release	46



Product Code Index

Product Code	Description	Page #	Product Code	Description	Page #
M00513850	Ultraflex™ Single-Use Covered Esophageal NG Stent System – Proximal Release	46	M00516930	WallFlex Partially Covered Esophageal Stent	45
M00513860	Ultraflex Single-Use Covered Esophageal NG Stent System – Proximal Release	46	M00516940	WallFlex Partially Covered Esophageal Stent	45
M00514200	Ultraflex Single-Use Covered Large Esophageal NG Stent System – Distal Release	47	M00516950	WallFlex Partially Covered Esophageal Stent	45
M00514210	Ultraflex Single-Use Covered Large Esophageal NG Stent System – Distal Release	47	M00518101	Interject™ Clear Single-Use Injection Therapy Needle Catheters	51
M00514240	Ultraflex Single-Use Covered Large Esophageal NG Stent System – Proximal Release	47	M00518111	Interject Clear Single-Use Injection Therapy Needle Catheters	51
M00514250	Ultraflex Single-Use Covered Large Esophageal NG Stent System – Proximal Release	47	M00518151	Interject Contrast Single-Use Injection Therapy Needle Catheters	51
M00514270	Polyflex™ Single-Use Esophageal Stent System	46	M00518161	Interject Contrast Single-Use Injection Therapy Needle Catheters	51
M00514280	Polyflex Single-Use Esophageal Stent System	46	M00518251	Interject Contrast Single-Use Injection Therapy Needle Catheters	51
M00514290	Polyflex Single-Use Esophageal Stent System	46	M00518261	Interject Contrast Single-Use Injection Therapy Needle Catheters	51
M00514300	Polyflex Single-Use Esophageal Stent System	46	M00518301	Interject Clear Single-Use Injection Therapy Needle Catheters	51
M00514310	Polyflex Single-Use Esophageal Stent System	46	M00518311	Interject Clear Single-Use Injection Therapy Needle Catheters	51
M00514320	Polyflex Single-Use Esophageal Stent System	46	M00518351	Interject™ Contrast Single-Use Injection Therapy Needle Catheters	51
M00514330	Polyflex Single-Use Esophageal Stent System	46	M00518361	Interject Contrast Single-Use Injection Therapy Needle Catheters	51
M00514340	Polyflex Single-Use Esophageal Stent System	46	M00521230	Resolution 360™ Clip	49
M00514350	Polyflex Single-Use Esophageal Stent System	46	M00521231	Resolution 360 Clip	49
M00515031	Radial Jaw™ 4 Single-Use Hot Biopsy Forceps	62	M00521232	Resolution 360 Clip	49
M00515032	Radial Jaw 4 Single-Use Hot Biopsy Forceps	62	M00521233	Resolution 360 Clip	49
M00515033	Radial Jaw 4 Single-Use Hot Biopsy Forceps	62	M00521240	Resolution 360 Clip	49
M00515171	Pulmonary Jagwire™ Pulmonary Guidewire	91	M00521241	Resolution 360 Clip	49
M00515181	Radial Jaw 4 Single-Use Pulmonary Forceps – Standard Capacity	89	M00521242	Resolution 360 Clip	49
M00515182	Radial Jaw 4 Single-Use Pulmonary Forceps – Standard Capacity	89	M00522600	Resolution™ Clip	49
M00515191	Radial Jaw 4 Single-Use Pulmonary Forceps – Standard Capacity with Needle	89	M00522601	Resolution Clip	49
M00515192	Radial Jaw 4 Single-Use Pulmonary Forceps – Standard Capacity with Needle	89	M00522602	Resolution Clip	49
M00515201	Radial Jaw 4 Single-Use Pulmonary Forceps – Large Capacity	89	M00522610	Resolution Clip	49
M00515202	Radial Jaw 4 Single-Use Pulmonary Forceps – Large Capacity	89	M00522611	Resolution Clip	49
M00516001	Celebrity™ Single-Use Endoscopic Cytology Brush	90	M00522612	Resolution Clip	49
M00516011	Celebrity Single-Use Endoscopic Cytology Brush	90	M00530800	Ultratome™ Double-Lumen Sphincterotome – Short Nose	10
M00516071	Celebrity Single-Use Endoscopic Cytology Brush	90	M00530810	Ultratome Double-Lumen Sphincterotome – Long Nose	10
M00516151	Celebrity Single-Use Endoscopic Cytology Brush	90	M00530820	Ultratome Double-Lumen Sphincterotome – Short Nose	10
M00516700	WallFlex™ Fully Covered Esophageal Stent	45	M00530830	Ultratome Double-Lumen Sphincterotome – Long Nose	10
M00516710	WallFlex Fully Covered Esophageal Stent	45	M00530850	Contour™ ERCP Cannulas – Standard	11
M00516720	WallFlex Fully Covered Esophageal Stent	45	M00530860	Contour ERCP Cannulas – Tapered	11
M00516730	WallFlex Fully Covered Esophageal Stent	45	M00530870	Contour ERCP Cannulas – Ultra Tapered	11
M00516740	WallFlex Fully Covered Esophageal Stent	45	M00530880	Contour ERCP Cannulas – 5-4-3	11
M00516750	WallFlex Fully Covered Esophageal Stent	45	M00530890	Contour ERCP Cannulas – Ball Tip	11
M00516900	WallFlex Partially Covered Esophageal Stent	45	M00532010	C-Flex™ Double Pigtail Biliary Stent	24
M00516910	WallFlex Partially Covered Esophageal Stent	45	M00532020	C-Flex Double Pigtail Biliary Stent	24
M00516920	WallFlex Partially Covered Esophageal Stent	45	M00532030	C-Flex Double Pigtail Biliary Stent	24
			M00532040	C-Flex Double Pigtail Biliary Stent	24
			M00532050	C-Flex Double Pigtail Biliary Stent	24



Product Code Index

Product Code	Description	Page #	Product Code	Description	Page #
M00532160	Advanix™ Double Pigtail - Single Stent	21	M00533390	Advanix Center Bend - Single Stent	20
M00532180	Advanix Double Pigtail - Single Stent	21	M00533400	Advanix Center Bend - Single Stent	20
M00532190	Advanix Double Pigtail - Single Stent	21	M00533560	NaviFlex™ RX Delivery System	21
M00532210	Advanix Double Pigtail - Single Stent	21	M00533570	NaviFlex RX Delivery System	21
M00532220	Advanix Double Pigtail - Single Stent	21	M00533580	NaviFlex RX Delivery System	21
M00532230	Advanix Double Pigtail - Single Stent	21	M00533610	Percuflex™ Duodenal Bend Biliary Stent	23
M00532240	Advanix Double Pigtail - Single Stent	21	M00533620	Percuflex Duodenal Bend Biliary Stent	23
M00532260	Advanix Double Pigtail - Single Stent	21	M00533630	Percuflex Duodenal Bend Biliary Stent	23
M00532270	Advanix Double Pigtail - Single Stent	21	M00533640	Percuflex Duodenal Bend Biliary Stent	23
M00532290	Advanix Double Pigtail - Single Stent	21	M00533650	Percuflex Duodenal Bend Biliary Stent	23
M00532300	Advanix Double Pigtail - Single Stent	21	M00533660	Percuflex Duodenal Bend Biliary Stent	23
M00532310	Advanix Double Pigtail - Single Stent	21	M00533670	Percuflex Duodenal Bend Biliary Stent	23
M00532810	MicroKnife™ XL Triple-Lumen Needle Knife	8	M00533680	Percuflex™ Duodenal Bend Biliary Stent	23
M00532870	Advanix Duodenal Bend - Single Stent	20	M00533690	Percuflex Duodenal Bend Biliary Stent	23
M00532880	Advanix Duodenal Bend - Single Stent	20	M00533700	Percuflex Duodenal Bend Biliary Stent	23
M00532890	Advanix Duodenal Bend - Single Stent	20	M00533760	Percuflex Duodenal Bend Biliary Stent w/ Push Catheter	23
M00532900	Advanix Duodenal Bend - Single Stent	20	M00533770	Percuflex Duodenal Bend Biliary Stent w/ Push Catheter	23
M00532910	Advanix Duodenal Bend - Single Stent	20	M00533780	Percuflex Duodenal Bend Biliary Stent w/ Push Catheter	23
M00532920	Advanix Duodenal Bend - Single Stent	20	M00533790	Percuflex Duodenal Bend Biliary Stent w/ Push Catheter	23
M00532930	Advanix Duodenal Bend - Single Stent	20	M00533800	Percuflex Duodenal Bend Biliary Stent w/ Push Catheter	23
M00532940	Advanix Duodenal Bend - Single Stent	20	M00533810	Percuflex Duodenal Bend Biliary Stent w/ Push and Guide Catheter	23
M00532950	Advanix Duodenal Bend - Single Stent	20	M00533820	Percuflex Duodenal Bend Biliary Stent w/ Push and Guide Catheter	23
M00532960	Advanix Duodenal Bend - Single Stent	20	M00533830	Percuflex Duodenal Bend Biliary Stent w/ Push and Guide Catheter	23
M00532970	Advanix Duodenal Bend - Single Stent	20	M00533840	Percuflex Duodenal Bend Biliary Stent w/ Push and Guide Catheter	23
M00532980	Advanix Duodenal Bend - Single Stent	20	M00533850	Percuflex Duodenal Bend Biliary Stent w/ Push and Guide Catheter	23
M00532990	Advanix Duodenal Bend - Single Stent	20	M00533910	Biliary Introducer – Push Catheter	23
M00533000	Advanix Duodenal Bend - Single Stent	20	M00533920	Biliary Introducer – Push and Guide Catheter	23
M00533010	Advanix Duodenal Bend - Single Stent	20	M00534200	Advanix Duodenal Bend Preloaded RX Stent Systems	19
M00533020	Advanix Duodenal Bend - Single Stent	20	M00534210	Advanix Duodenal Bend Preloaded RX Stent Systems	19
M00533030	Advanix Duodenal Bend - Single Stent	20	M00534220	Advanix Duodenal Bend Preloaded RX Stent Systems	19
M00533040	Advanix Duodenal Bend - Single Stent	20	M00534230	Advanix Duodenal Bend Preloaded RX Stent Systems	19
M00533230	Advanix Center Bend - Single Stent	20	M00534240	Advanix Duodenal Bend Preloaded RX Stent Systems	19
M00533240	Advanix Center Bend - Single Stent	20	M00534250	Advanix Duodenal Bend Preloaded RX Stent Systems	19
M00533250	Advanix Center Bend - Single Stent	20	M00534260	Advanix Duodenal Bend Preloaded RX Stent Systems	19
M00533260	Advanix Center Bend - Single Stent	20	M00534270	Advanix Duodenal Bend Preloaded RX Stent Systems	19
M00533270	Advanix Center Bend - Single Stent	20	M00534280	Advanix Duodenal Bend Preloaded RX Stent Systems	19
M00533280	Advanix Center Bend - Single Stent	20	M00534290	Advanix Duodenal Bend Preloaded RX Stent Systems	19
M00533290	Advanix Center Bend - Single Stent	20	M00534300	Advanix Duodenal Bend Preloaded RX Stent Systems	19
M00533300	Advanix Center Bend - Single Stent	20	M00534310	Advanix Duodenal Bend Preloaded RX Stent Systems	19
M00533310	Advanix Center Bend - Single Stent	20	M00534320	Advanix Duodenal Bend Preloaded RX Stent Systems	19
M00533320	Advanix Center Bend - Single Stent	20	M00534330	Advanix Duodenal Bend Preloaded RX Stent Systems	19
M00533330	Advanix Center Bend - Single Stent	20	M00534340	Advanix Duodenal Bend Preloaded RX Stent Systems	19
M00533340	Advanix Center Bend - Single Stent	20	M00534350	Advanix Duodenal Bend Preloaded RX Stent Systems	19
M00533350	Advanix Center Bend - Single Stent	20			
M00533360	Advanix Center Bend - Single Stent	20			
M00533370	Advanix Center Bend - Single Stent	20			
M00533380	Advanix Center Bend - Single Stent	20			



Product Code Index

Product Code	Description	Page #	Product Code	Description	Page #
M00534360	Advanix™ Duodenal Bend Preloaded RX Stent Systems	19	M00536090	Advanix Pancreatic Stent – Straight, No Leading Barb	26
M00534370	Advanix Duodenal Bend Preloaded RX Stent Systems	19	M00536100	Advanix Pancreatic Stent – Straight, No Leading Barb	26
M00534560	Advanix Center Bend Preloaded RX Stent Systems	19	M00536110	Advanix Pancreatic Stent – Straight, No Leading Barb	26
M00534570	Advanix Center Bend Preloaded RX Stent Systems	19	M00536120	Advanix Pancreatic Stent – Straight, No Leading Barb	26
M00534580	Advanix Center Bend Preloaded RX Stent Systems	19	M00536130	Advanix Pancreatic Stent – Straight, No Leading Barb	26
M00534590	Advanix Center Bend Preloaded RX Stent Systems	19	M00536140	Advanix Pancreatic Stent – Straight, No Leading Barb	26
M00534600	Advanix Center Bend Preloaded RX Stent Systems	19	M00536150	Advanix Pancreatic Stent – Straight, No Leading Barb	26
M00534610	Advanix Center Bend Preloaded RX Stent Systems	19	M00536160	Advanix Pancreatic Stent – Straight, No Leading Barb	26
M00534620	Advanix Center Bend Preloaded RX Stent Systems	19	M00536170	Advanix Pancreatic Stent – Straight, No Leading Barb	26
M00534630	Advanix Center Bend Preloaded RX Stent Systems	19	M00536180	Advanix Pancreatic Stent – Straight, No Leading Barb	26
M00534640	Advanix Center Bend Preloaded RX Stent Systems	19	M00536190	Advanix Pancreatic Stent – Straight, No Leading Barb	26
M00534650	Advanix Center Bend Preloaded RX Stent Systems	19	M00536200	Advanix Pancreatic Stent – Straight, No Leading Barb	26
M00534660	Advanix Center Bend Preloaded RX Stent Systems	19	M00536210	Advanix Pancreatic Stent – Straight, No Leading Barb	26
M00534670	Advanix Center Bend Preloaded RX Stent Systems	19	M00536220	Advanix Pancreatic Stent – Straight, No Leading Barb	26
M00534680	Advanix Center Bend Preloaded RX Stent Systems	19	M00536230	Advanix Pancreatic Stent – Straight, No Leading Barb	26
M00534690	Advanix Center Bend Preloaded RX Stent Systems	19	M00536240	Advanix Pancreatic Stent – Straight, No Leading Barb	26
M00534700	Advanix Center Bend Preloaded RX Stent Systems	19	M00536250	Advanix Pancreatic Stent – Straight, No Leading Barb	26
M00534710	Advanix Center Bend Preloaded RX Stent Systems	19	M00536260	Advanix Pancreatic Stent – Straight, No Leading Barb	26
M00534720	Advanix Center Bend Preloaded RX Stent Systems	19	M00536270	Advanix Pancreatic Stent – Straight, No Leading Barb	26
M00534730	Advanix Center Bend Preloaded RX Stent Systems	19	M00536280	Advanix Pancreatic Stent – Straight, No Leading Barb	26
M00535000	NaviFlex™ Pancreatic Stent RX Pusher	30	M00536290	Advanix Pancreatic Stent – Straight, No Leading Barb	26
M00535010	NaviFlex Pancreatic Stent RX Pusher	30	M00536300	Advanix Pancreatic Stent – Straight, No Leading Barb	26
M00535020	NaviFlex Pancreatic Stent Pusher	30	M00536310	Advanix Pancreatic Stent – Straight, No Leading Barb	26
M00535030	NaviFlex Pancreatic Stent Pusher	30	M00536320	Advanix Pancreatic Stent – Straight, No Leading Barb	26
M00535060	NaviFlex RX Pancreatic Stent Delivery System	30	M00536330	Advanix Pancreatic Stent – Straight, No Leading Barb	26
M00535080	NaviFlex RX Pancreatic Stent Delivery System	30	M00536340	Advanix Pancreatic Stent – Straight, No Leading Barb	26
M00535110	Stonetome™ Stone Removal Device	9	M00536350	Advanix Pancreatic Stent – Straight, No Leading Barb	26
M00535130	Stonetome Stone Removal Device	9	M00536360	Advanix Pancreatic Stent – Straight, Leading Barb	27
M00535150	Stonetome Stone Removal Device	9	M00536370	Advanix Pancreatic Stent – Straight, Leading Barb	27
M00535170	Stonetome Stone Removal Device	9	M00536380	Advanix Pancreatic Stent – Straight, Leading Barb	27
M00535190	Stonetome Stone Removal Device	9	M00536390	Advanix Pancreatic Stent – Straight, Leading Barb	27
M00535210	Stonetome Stone Removal Device	9	M00536400	Advanix Pancreatic Stent – Straight, Leading Barb	27
M00535700	Tandem™ XL Triple-Lumen ERCP Cannula – Tapered	11	M00536410	Advanix Pancreatic Stent – Straight, Leading Barb	27
M00535900	Ultratome™ XL Triple-Lumen Sphincterotome – Short Nose	9	M00536420	Advanix Pancreatic Stent – Straight, Leading Barb	27
M00535910	Ultratome XL Triple-Lumen Sphincterotome – Long Nose	9	M00536430	Advanix Pancreatic Stent – Straight, Leading Barb	27
M00535920	Ultratome XL Triple-Lumen Sphincterotome – Short Nose	9	M00536440	Advanix Pancreatic Stent – Straight, Leading Barb	27
M00535930	Ultratome XL Triple-Lumen Sphincterotome – Long Nose	9	M00536450	Advanix Pancreatic Stent – Straight, Leading Barb	27
M00536000	Advanix Pancreatic Stent – Straight, No Leading Barb	26	M00536460	Advanix Pancreatic Stent – Straight, Leading Barb	27
M00536010	Advanix Pancreatic Stent – Straight, No Leading Barb	26	M00536470	Advanix Pancreatic Stent – Straight, Leading Barb	27
M00536020	Advanix Pancreatic Stent – Straight, No Leading Barb	26	M00536480	Advanix Pancreatic Stent – Straight, Leading Barb	27
M00536030	Advanix Pancreatic Stent – Straight, No Leading Barb	26	M00536490	Advanix Pancreatic Stent – Straight, Leading Barb	27
M00536040	Advanix Pancreatic Stent – Straight, No Leading Barb	26	M00536500	Advanix Pancreatic Stent – Straight, Leading Barb	27
M00536050	Advanix Pancreatic Stent – Straight, No Leading Barb	26	M00536510	Advanix Pancreatic Stent – Straight, Leading Barb	27
M00536060	Advanix Pancreatic Stent – Straight, No Leading Barb	26	M00536520	Advanix Pancreatic Stent – Straight, Leading Barb	27
M00536070	Advanix Pancreatic Stent – Straight, No Leading Barb	26	M00536530	Advanix Pancreatic Stent – Straight, Leading Barb	27
M00536080	Advanix Pancreatic Stent – Straight, No Leading Barb	26	M00536540	Advanix Pancreatic Stent – Straight, Leading Barb	27
			M00536550	Advanix Pancreatic Stent – Straight, Leading Barb	27



Product Code Index

[illegible]

Product Code Index

Product Code	Description	Page #	Product Code	Description	Page #
M00537500	Advanix™ Pancreatic Stent Kit – Straight, No Leading Barb ..	25	M00545281	RX Locking Device and Biopsy Cap System – FujiFilm Type 3 Brand Endoscope Compatible	12
M00537510	Advanix Pancreatic Stent Kit – Straight, No Leading Barb	25	M00545471	RX Locking Device (Biopsy Cap not included).....	12
M00537520	Advanix Pancreatic Stent Kit – Straight, No Leading Barb	25	M00545550	RX Duodenal Bend Biliary Stent with RX Delivery System	22
M00537530	Advanix Pancreatic Stent Kit – Straight, No Leading Barb	25	M00545560	RX Duodenal Bend Biliary Stent with RX Delivery System	22
M00537540	Advanix Pancreatic Stent Kit – Straight, Leading Barb.....	25	M00545570	RX Duodenal Bend Biliary Stent with RX Delivery System	22
M00537550	Advanix Pancreatic Stent Kit – Straight, Leading Barb.....	25	M00545580	RX Duodenal Bend Biliary Stent with RX Delivery System	22
M00537560	Advanix Pancreatic Stent Kit – Straight, Leading Barb.....	25	M00545590	RX Duodenal Bend Biliary Stent with RX Delivery System	22
M00539200	Flexima™ Duodenal Bend Biliary Stent with Delivery System.....	22	M00545600	RX Duodenal Bend Biliary Stent with RX Delivery System	22
M00539210	Flexima Duodenal Bend Biliary Stent with Delivery System..	22	M00545610	RX Duodenal Bend Biliary Stent with RX Delivery System	22
M00539220	Flexima Duodenal Bend Biliary Stent with Delivery System..	22	M00545620	RX Duodenal Bend Biliary Stent with RX Delivery System	22
M00539230	Flexima Duodenal Bend Biliary Stent with Delivery System..	22	M00545630	RX Duodenal Bend Biliary Stent with RX Delivery System	22
M00539240	Flexima Duodenal Bend Biliary Stent with Delivery System..	22	M00545640	RX Duodenal Bend Biliary Stent with RX Delivery System	22
M00539250	Flexima Duodenal Bend Biliary Stent with Delivery System..	22	M00545650	RX Duodenal Bend Biliary Stent with RX Delivery System	22
M00539260	Flexima Duodenal Bend Biliary Stent with Delivery System..	22	M00545660	RX Duodenal Bend Biliary Stent with RX Delivery System	22
M00539270	Flexima Duodenal Bend Biliary Stent with Delivery System..	22	M00545670	RX Duodenal Bend Biliary Stent with RX Delivery System	22
M00539280	Flexima Duodenal Bend Biliary Stent with Delivery System..	22	M00545680	RX Duodenal Bend Biliary Stent with RX Delivery System	22
M00539290	Flexima Duodenal Bend Biliary Stent with Delivery System..	22	M00545690	RX Duodenal Bend Biliary Stent with RX Delivery System	22
M00539300	Flexima Duodenal Bend Biliary Stent with Delivery System..	22	M00545700	RX Duodenal Bend Biliary Stent with RX Delivery System	22
M00539310	Flexima Duodenal Bend Biliary Stent with Delivery System..	22	M00545710	RX Duodenal Bend Biliary Stent with RX Delivery System	22
M00539320	Flexima Duodenal Bend Biliary Stent with Delivery System....	22	M00545720	RX Duodenal Bend Biliary Stent with RX Delivery System	22
M00539330	Flexima Duodenal Bend Biliary Stent with Delivery System..	22	M00545730	RX Duodenal Bend Biliary Stent with RX Delivery System	22
M00539340	Flexima Duodenal Bend Biliary Stent with Delivery System..	22	M00545740	RX Duodenal Bend Biliary Stent with RX Delivery System	22
M00539350	Flexima Duodenal Bend Biliary Stent with Delivery System..	22	M00545760	RX ERCP Cannula – Standard Tip	10
M00539360	Flexima Duodenal Bend Biliary Stent with Delivery System..	22	M00545770	RX ERCP Cannula – Tapered Tip	10
M00539370	Flexima Duodenal Bend Biliary Stent with Delivery System..	22	M00545780	RX ERCP Cannula – Ball Tip.....	10
M00539380	Flexima Duodenal Bend Biliary Stent with Delivery System..	22	M00545820	Rapid Exchange XL Cannula	10
M00539390	Flexima Duodenal Bend Biliary Stent with Delivery System..	22	M00545840	RX Needle Knife XL Triple-Lumen Needle Knife.....	6
M00540110	Flexima Single Pigtail Nasobiliary Catheter.....	24	M00545890	Hurricane™ RX Biliary Balloon Dilatation Catheter.....	34
M00540120	Flexima Single Pigtail Nasobiliary Catheter.....	24	M00545900	Hurricane RX Biliary Balloon Dilatation Catheter.....	34
M00540130	Flexima Single Pigtail Nasobiliary Catheter.....	24	M00545910	Hurricane RX Biliary Balloon Dilatation Catheter.....	34
M00540140	Flexima Single Pigtail Nasobiliary Catheter.....	24	M00545920	Hurricane RX Biliary Balloon Dilatation Catheter.....	34
M00542250	Speedband Superview Super 7™ Multiple Band Ligators.....	52	M00545930	Hurricane RX Biliary Balloon Dilatation Catheter.....	34
M00542251	Speedband Superview Super 7™ Multiple Band Ligators.....	52	M00545940	Hurricane RX Biliary Balloon Dilatation Catheter.....	34
M00542253	Speedband Superview Super 7 Multiple Band Ligators.....	52	M00545950	Hurricane RX Biliary Balloon Dilatation Catheter.....	34
M00545000	RX Cytology Brush Wireguided Cytology Brush	35	M00545960	Hurricane RX Biliary Balloon Dilatation Catheter.....	34
M00545150	Autotome™ RX 49 Cannulating Sphincterotome.....	6	M00546270	SpyBite™ Biopsy Forceps.....	37
M00545160	Autotome RX 49 Cannulating Sphincterotome.....	6	M00546600	SpyScope™ DS Access & Delivery Catheter.....	37
M00545170	Autotome RX 44 Cannulating Sphincterotome.....	6	M00546620	Biliary Electrohydraulic Lithotripter Probe	38
M00545180	Autotome RX 44 Cannulating Sphincterotome.....	6	M00546650	SpyGlass DS Digital Controller	37
M00545190	Autotome RX 39 Cannulating Sphincterotome.....	6	M005466800	Autolith™ Touch Bipolar Electrohydraulic Lithotripter	38
M00545200	Autotome RX 39 Cannulating Sphincterotome.....	6	M00546750	Autolith Touch Extender Cable	38
M00545211	RX Locking Device Olympus™ / FujiFilm™ Brand Endoscope Compatible.....	12	M00546760	Autolith Touch Foot Switch.....	38
M00545261	RX Locking Device and Biopsy Cap System – Olympus FujiFilm G5 Series Brand Endoscope Compatible.....	12	M00547000	Extractor™ Pro RX Retrieval Balloon Catheter – Distal Injection.....	31
M00545271	RX Locking Device and Biopsy Cap System – Pentax™ Brand Endoscope Compatible.....	12	M00547010	Extractor Pro RX Retrieval Balloon Catheter – Distal Injection	31
			M00547020	Extractor Pro RX Retrieval Balloon Catheter – Distal Injection	31



Product Code Index

Product Code	Description	Page #	Product Code	Description	Page #
M00547030	Extractor™ Pro RX Retrieval Balloon Catheter – Proximal Injection	31	M00555520	Expect Slimline (SL) Endoscopic Ultrasound Aspiration Needle	41
M00547040	Extractor Pro RX Retrieval Balloon Catheter – Proximal Injection	31	M00555530	Expect Slimline (SL) Endoscopic Ultrasound Aspiration Needle	41
M00547050	Extractor Pro RX Retrieval Balloon Catheter – Proximal Injection	31	M00555501	Expect Slimline (SL) Endoscopic Ultrasound Aspiration Needle	41
M00547100	Extractor Pro XL Retrieval Balloon Catheter – Distal Injection	32	M00555511	Expect Slimline (SL) Endoscopic Ultrasound Aspiration Needle	41
M00547110	Extractor Pro XL Retrieval Balloon Catheter – Distal Injection	32	M00555521	Expect Slimline (SL) Endoscopic Ultrasound Aspiration Needle	41
M00547120	Extractor Pro XL Retrieval Balloon Catheter – Distal Injection	32	M00555531	Expect Slimline (SL) Endoscopic Ultrasound Aspiration Needle	41
M00547130	Extractor Pro XL Retrieval Balloon Catheter – Proximal Injection	32	M00555540	Acquire™ Endoscopic Ultrasound Fine Needle Biopsy (FNB) Device	41
M00547140	Extractor Pro XL Retrieval Balloon Catheter – Proximal Injection	32	M00555541	Acquire Endoscopic Ultrasound Fine Needle Biopsy (FNB) Device	41
M00547150	Extractor Pro XL Retrieval Balloon Catheter – Proximal Injection	32	M00555560	Acquire Endoscopic Ultrasound Fine Needle Biopsy (FNB) Device	41
M00547200	Extractor Pro DL Retrieval Balloon Catheter	32	M00555561	Acquire Endoscopic Ultrasound Fine Needle Biopsy (FNB) Device	41
M00547210	Extractor Pro DL Retrieval Balloon Catheter	32	M00555580	Acquire Endoscopic Ultrasound Fine Needle Biopsy (FNB) Device	41
M00547220	Extractor Pro DL Retrieval Balloon Catheter	32	M00555581	Acquire Endoscopic Ultrasound Fine Needle Biopsy (FNB) Device	41
M00547300	Extractor Pro RX-S Retrieval Balloon Catheter – Distal Injection	31	M00556001	Hydra Jagwire™ High Performance Guidewire – Straight	1
M00547310	Extractor Pro RX-S Retrieval Balloon Catheter – Distal Injection	31	M00556011	Hydra Jagwire High Performance Guidewire – Angled	1
M00547320	Extractor Pro RX-S Retrieval Balloon Catheter – Distal Injection	31	M00556021	Hydra Jagwire High Performance Guidewire – Stiff Shaft – Straight	1
M00547330	Extractor Pro RX-S Retrieval Balloon Catheter – Proximal Injection	31	M00556031	Hydra Jagwire High Performance Guidewire – Stiff Shaft – Angled	1
M00547340	Extractor Pro RX-S Retrieval Balloon Catheter – Proximal Injection	31	M00556041	Hydra Jagwire High Performance Guidewire – Straight	1
M00547350	Extractor Pro RX-S Retrieval Balloon Catheter – Proximal Injection	31	M00556051	Hydra Jagwire High Performance Guidewire – Angled	1
M00550090	Amplatz Super Stiff™ Guidewires – Straight Tip	47, 91	M00556061	Hydra Jagwire High Performance Guidewire – Stiff Shaft – Straight	1
M00550300	CRE™ Pulmonary Balloon Dilatation Catheter	91	M00556071	Hydra Jagwire High Performance Guidewire – Stiff Shaft – Angled	1
M00550310	CRE Pulmonary Balloon Dilatation Catheter	91	M00556101	Dreamwire™ High Performance Guidewire – Standard Shaft	1
M00550320	CRE Pulmonary Balloon Dilatation Catheter	91	M00556111	Dreamwire High Performance Guidewire – Standard Shaft	1
M00550330	CRE Pulmonary Balloon Dilatation Catheter	91	M00556121	Dreamwire High Performance Guidewire – Stiff Shaft	1
M00550340	CRE Pulmonary Balloon Dilatation Catheter	91	M00556131	Dreamwire High Performance Guidewire – Stiff Shaft	1
M00550350	CRE Pulmonary Balloon Dilatation Catheter	91	M00556141	Dreamwire High Performance Guidewire – Standard Shaft	1
M00550601	Alliance™ II Integrated Inflation / Lithotripsy Device and Alliance II Syringe / Gauge Assembly	60, 91	M00556151	Dreamwire High Performance Guidewire – Standard Shaft	1
M00550620	Alliance II Integrated Inflation / Lithotripsy Device – Inflation Handle	60, 91	M00556161	Dreamwire High Performance Guidewire – Stiff Shaft	1
M00550630	CRE SteriFlate Disposable Inflation Device	60	M00556171	Dreamwire High Performance Guidewire – Stiff Shaft	1
M00552001	NovaGold™ High Performance Guidewire	3	M00556181	NaviPro™ Hydrophilic Guidewire	3
M00552011	NovaGold High Performance Guidewire	3	M00556191	NaviPro Hydrophilic Guidewire	3
M00553200	Inflation / Deflation Devices and Accessories – Pneumatic Inflator	60	M00556201	NaviPro Hydrophilic Guidewire	3
M00553640	AXIOS™ Stent and Electrocautery Enhanced Delivery System	42	M00556211	NaviPro Hydrophilic Guidewire	3
M00553650	AXIOS Stent and Electrocautery Enhanced Delivery System	42	M00556221	NaviPro Hydrophilic Guidewire	3
M00553740	AXIOS Stent and Delivery System	42	M00556231	NaviPro Hydrophilic Guidewire	3
M00553750	AXIOS Stent and Delivery System	42	M00556241	NaviPro Hydrophilic Guidewire	3
M00554500	Rigiflex™ II Single-Use Achalasia Balloon Dilators	59	M00556251	NaviPro Hydrophilic Guidewire	3
M00554510	Rigiflex II Single-Use Achalasia Balloon Dilators	59	M00556461	Jagwire High Performance Guidewire – Straight	2
M00554520	Rigiflex II Single-Use Achalasia Balloon Dilators	59	M00556481	Stiff Shaft Jagwire Guidewire – Straight	2
M00555500	Expect™ Slimline (SL) Endoscopic Ultrasound Aspiration Needle	41	M00556491	Stiff Shaft Jagwire Guidewire – Angled	2
M00555510	Expect Slimline (SL) Endoscopic Ultrasound Aspiration Needle	41	M00556511	Jagwire High Performance Guidewire – Stiff Shaft – Angled	2
			M00556521	Jagwire High Performance Guidewire – Stiff Shaft – Straight	2
			M00556531	Jagwire High Performance Guidewire – Stiff Shaft – Angled	2
			M00556561	Jagwire High Performance Guidewire – Stiff Shaft – Straight	2
			M00556581	Jagwire High Performance Guidewire – Stiff Shaft – Straight	2



Product Code Index

Product Code	Description	Page #	Product Code	Description	Page #
M00556591	Jagwire™ High Performance Guidewire – Stiff Shaft – Angled..	2	M00558470	CRE Wireguided Esophageal / Pyloric / Colonic Balloon Dilatation Catheters	57
M00556601	Jagwire High Performance Guidewire – Stiff Shaft – Straight.....	2	M00558480	CRE Wireguided Esophageal / Pyloric / Colonic Balloon Dilatation Catheters	57
M00556611	Jagwire High Performance Guidewire – Stiff Shaft – Angled	2	M00558490	CRE Wireguided Esophageal / Pyloric / Colonic Balloon Dilatation Catheters	57
M00556621	Jagwire High Performance Guidewire – Straight	2	M00558500	CRE Wireguided Esophageal / Pyloric / Colonic Balloon Dilatation Catheters	57
M00556631	Jagwire High Performance Guidewire – Stiff Shaft – Straight.....	2	M00558600	CRE PRO Wireguided Biliary Balloon Dilators	33
M00556641	Jagwire High Performance Guidewire – Straight	2		CRE PRO Wireguided Esophageal / Pyloric / Biliary Balloon Dilatation Catheters	56
M00556651	Jagwire High Performance Guidewire – Angled.....	2	M00558610	CRE PRO Wireguided Biliary Balloon Dilatation Catheters	33
M00556661	Jagwire High Performance Guidewire – Stiff Shaft – Angled	2		CRE PRO Wireguided Esophageal / Pyloric / Biliary Balloon Dilatation Catheters	56
M00556841	Jagwire Extendable High Performance Guidewire – Straight....	2	M00558620	CRE PRO Wireguided Biliary Balloon Dilatation Catheters	33
M00556851	Jagwire Extendable High Performance Guidewire – Angled	2		CRE PRO Wireguided Esophageal / Pyloric / Biliary Balloon Dilatation Catheters	56
M00556861	Jagwire Extendable High Performance Guidewire – Stiff Shaft – Straight	2	M00558630	CRE PRO Wireguided Biliary Balloon Dilatation Catheters	33
M00556871	Jagwire Extendable High Performance Guidewire – Stiff Shaft – Angled.....	2		CRE PRO Wireguided Esophageal / Pyloric / Biliary Balloon Dilatation Catheters	56
M00556902	Jagtail™ High Performance Guidewire Extension	2	M00558640	CRE PRO Wireguided Biliary Balloon Dilatation Catheters	33
M00557330	Ultraflex™ Precision Single-Use Colonic Stent System.....	44		CRE PRO Wireguided Esophageal / Pyloric / Biliary Balloon Dilatation Catheters	56
M00557360	Ultraflex Precision Single-Use Colonic Stent System.....	44	M00558650	CRE PRO Wireguided Biliary Balloon Dilatation Catheters	33
M00557380	Ultraflex Precision Single-Use Colonic Stent System.....	44		CRE PRO Wireguided Esophageal / Pyloric / Biliary Balloon Dilatation Catheters	56
M00558220	Expect™ Pulmonary Endobronchial Ultrasound Transbronchial Aspiration Needle.....	89	M00558660	CRE PRO Wireguided Biliary Balloon Dilatation Catheters	33
M00558250	Expect Pulmonary Endobronchial Ultrasound Transbronchial Aspiration Needle.....	89		CRE PRO Wireguided Esophageal / Pyloric / Colonic / Biliary Balloon Dilatation Catheters	56
M00558330	CRE™ Fixed Wire Esophageal Balloon Dilatation Catheters....	58	M00558670	CRE PRO Wireguided Biliary Balloon Dilatation Catheters	33
M00558340	CRE Fixed Wire Esophageal Balloon Dilatation Catheters....	58		CRE PRO Wireguided Esophageal / Pyloric / Colonic / Biliary Balloon Dilatation Catheters	56
M00558350	CRE Fixed Wire Esophageal Balloon Dilatation Catheters....	58	M00558680	CRE PRO Wireguided Biliary Balloon Dilatation Catheters	33
M00558360	CRE Fixed Wire Esophageal Balloon Dilatation Catheters....	58		CRE PRO Wireguided Esophageal / Pyloric / Colonic / Biliary Balloon Dilatation Catheters	56
M00558361	CRE Fixed Wire Esophageal Balloon Dilatation Catheters....	58	M00558690	CRE PRO Wireguided Biliary Balloon Dilatation Catheters	33
M00558370	CRE Fixed Wire Esophageal Balloon Dilatation Catheters....	58		CRE PRO Wireguided Esophageal / Pyloric / Colonic / Biliary Balloon Dilatation Catheters	56
M00558371	CRE Fixed Wire Esophageal Balloon Dilatation Catheters....	58	M00558700	CRE PRO Wireguided Biliary Balloon Dilatation Catheters	33
M00558380	CRE Fixed Wire Esophageal Balloon Dilatation Catheters....	58		CRE PRO Wireguided Esophageal / Pyloric / Colonic / Biliary Balloon Dilatation Catheters	56
M00558381	CRE Fixed Wire Esophageal Balloon Dilatation Catheters....	58	M00558710	CRE PRO Wireguided Biliary Balloon Dilatation Catheters	33
M00558390	CRE Wireguided Esophageal / Pyloric Balloon Dilatation Catheters	57		CRE PRO Wireguided Esophageal / Pyloric / Colonic / Biliary Balloon Dilatation Catheters	56
M00558400	CRE Wireguided Esophageal / Pyloric Balloon Dilatation Catheters	57	M00558731	Expect Pulmonary Endobronchial Ultrasound Transbronchial Aspiration Needle Adaptor.....	89
M00558410	CRE Wireguided Esophageal / Pyloric Balloon Dilatation Catheters	57	M00558860	CRE RX Biliary Balloon Dilatation Catheters	55
M00558420	CRE Wireguided Esophageal / Pyloric Balloon Dilatation Catheters	57	M00558870	CRE RX Biliary Balloon Dilatation Catheters	55
M00558430	CRE Wireguided Esophageal / Pyloric Balloon Dilatation Catheters	57	M00558880	CRE RX Biliary Balloon Dilatation Catheters	55
M00558440	CRE Wireguided Esophageal / Pyloric Balloon Dilatation Catheters	57	M00558890	CRE RX Biliary Balloon Dilatation Catheters	55
M00558450	CRE Wireguided Esophageal / Pyloric / Colonic Balloon Dilatation Catheters	57	M00558900	CRE RX Biliary Balloon Dilatation Catheters	55
M00558460	CRE Wireguided Esophageal / Pyloric / Colonic Balloon Dilatation Catheters	57	M00558910	CRE RX Biliary Balloon Dilatation Catheters	55
			M00558920	CRE RX Biliary Balloon Dilatation Catheters	55
			M00558930	CRE RX Biliary Balloon Dilatation Catheters	55
			M00558940	CRE RX Biliary Balloon Dilatation Catheters	55



Product Code Index

Product Code	Description	Page #	Product Code	Description	Page #
M00560070	Gold Probe™ Single-Use Electrohemostasis Catheters – Standard Plug	50	M00561831	Single-Use Rotatable Snares – Medium Oval	77
M00560071	Gold Probe Single-Use Electrohemostasis Catheters – Standard Plug	50	M00562301	Captivator Single-Use Snares – Small Oval	76
M00560100	Gold Probe Single-Use Electrohemostasis Catheters – Standard Plug	50	M00562321	Captivator Single-Use Snares – Medium Oval	76
M00560101	Gold Probe Single-Use Electrohemostasis Catheters – Standard Plug	50	M00562341	Captivator Single-Use Snares – Medium Hexagonal	76
M00560150	Injection Gold Probe™ Single-Use Electrohemostasis / Injection Catheters – Standard Plug	50	M00562371	Captivator Single-Use Snares – Medium Crescent	76
M00560160	Injection Gold Probe Single-Use Electrohemostasis / Injection Catheters – Standard Plug	50	M00562391	Captivator Single-Use Snares – Large Oval	76
M00560220	Gold Probe Single-Use Electrohemostasis Catheters – Standard Plug	50	M00562401	Captiflex™ Single-Use Snares – Medium Oval	78
M00560311	Sensation™ Single-Use Medium Stiff Wire Short Throw Snares – Medium Oval	78	M00562402	Captiflex Single-Use Snares – Medium Oval	78
M00560321	Sensation Single-Use Medium Stiff Wire Short Throw Snares – Small Oval	78	M00562421	Captiflex Single-Use Snares – Small Oval	78
M00560331	Sensation Single-Use Medium Stiff Wire Short Throw Snares – Medium Crescent	78	M00562422	Captiflex Single-Use Snares – Small Oval	78
M00561191	Captivator™ II Single-Use Snares – Rounded Stiff	75	M00562451	Captivator™ Single-Use Snares – Small Hexagonal	76
M00561192	Captivator II Single-Use Snares – Rounded Stiff	75	M00562471	Captiflex Single-Use Snares – Extra Small Oval	78
M00561193	Captivator II Single-Use Snares – Rounded Stiff	75	M00562531	Profile™ Single-Use Snares – Extra Small Oval	77
M00561221	Captivator II Single-Use Snares – Rounded Stiff	75	M00562551	Profile Single-Use Snares – Small Oval	77
M00561222	Captivator II Single-Use Snares – Rounded Stiff	75	M00562571	Profile Single-Use Snares – Medium Oval	77
M00561223	Captivator II Single-Use Snares – Rounded Stiff	75	M00562651	Sensation Single-Use Short Throw Snares – Large Oval / Medium Stiff	78
M00561231	Captivator II Single-Use Snares – Rounded Stiff	75	M00562652	Sensation Single-Use Short Throw Snares – Large Oval / Medium Stiff	78
M00561232	Captivator II Single-Use Snares – Rounded Stiff	75	M00562671	Sensation Single-Use Short Throw Snares – Medium Oval / Flexible	78
M00561233	Captivator II Single-Use Snares – Rounded Stiff	75	M00562672	Sensation Single-Use Short Throw Snares – Medium Oval / Flexible	78
M00561241	Captivator II Single-Use Snares – Rounded Stiff	75	M00562673	Sensation Single-Use Short Throw Snares – Medium Oval / Flexible	78
M00561242	Captivator II Single-Use Snares – Rounded Stiff	75	M00562691	Sensation Single-Use Short Throw Snares – Small Oval / Flexible	78
M00561243	Captivator II Single-Use Snares – Rounded Stiff	75	M00562692	Sensation Single-Use Short Throw Snares – Small Oval / Flexible	78
M00561250	Active Cords – Fits Olympus™ Brand Generator	12, 79	M00562693	Sensation Single-Use Short Throw Snares – Small Oval / Flexible	78
M00561270	Active Cord, Fits Endostat™, Endostat II and Endostat III Generators, Bovie™, Valleylab™, Aspen Labs™, and ERBE™ Brand Generators	12, 79	M00562760	EndoVive™ Low Profile Replacement Button	72
M00561291	Captivator II Single-Use Snare – Extra Large Rounded	75	M00562770	EndoVive Low Profile Replacement Button	72
M00561311	Captivator Single-Use Snares – Large Oval Medium	76	M00562780	EndoVive Low Profile Replacement Button	72
M00561440	Electrosurgical Bipolar Cable Adapter	53	M00562790	EndoVive Low Profile Replacement Button	72
M00561450	Electrosurgical Bipolar Cable Adapter	53	M00562800	EndoVive Low Profile Replacement Button	72
M00561460	Electrosurgical Bipolar Cable Adapter	53	M00562820	EndoVive Low Profile Replacement Button	72
M00561470	Electrosurgical Bipolar Cable Adapter	53	M00562830	EndoVive Low Profile Replacement Button	72
M00561480	Electrosurgical Bipolar Cable Adapter	53	M00562840	EndoVive Low Profile Replacement Button	72
M00561490	Electrosurgical Bipolar Cable Adapter	53	M00562850	EndoVive Low Profile Replacement Button	72
M00561500	Electrosurgical Bipolar Cable Adapter	53	M00563010	EndoVive One Step Button™ Low Profile PEG Kit, Pull	70
M00561600	Captivator EMR Device for Standard Gastrosopes	75	M00563020	EndoVive One Step Button Low Profile PEG Kit, Pull	70
M00561601	Captivator EMR Device for Standard Gastrosopes	75	M00563030	EndoVive One Step Button Low Profile PEG Kit, Pull	70
M00561610	Captivator EMR Device for Large Gastrosopes	75	M00563040	EndoVive One Step Button Low Profile PEG Kit, Pull	70
M00561611	Captivator EMR Device for Large Gastrosopes	75	M00563080	EndoVive One Step Button Low Profile PEG Kit, Push	70
M00561821	Single-Use Rotatable Snares – Small Oval	77	M00563090	EndoVive One Step Button Low Profile PEG Kit, Push	70
			M00563100	EndoVive One Step Button Low Profile PEG Kit, Push	70
			M00563120	EndoVive One Step Button Low Profile PEG Kit, Pull	70
			M00563130	EndoVive One Step Button Low Profile PEG Kit, Pull	70



Product Code Index

Product Code	Description	Page #	Product Code	Description	Page #
M00563140	EndoVive™ One Step Button Low Profile PEG Kit, Pull	70	M00566001	RX Accessories – Rapid Refill Continuous Injection System	12
M00563150	EndoVive One Step Button Low Profile PEG Kit, Pull	70	M00566320	EndoVive Through-the-PEG (TTP) Jejunal Feeding Tube – Bent Tip	73
M00563160	EndoVive One Step Button Low Profile PEG Kit, Pull	70	M00566330	EndoVive Through-the-PEG (TTP) Jejunal Feeding Tube – Pigtail Tip	73
M00563180	EndoVive One Step Button Low Profile PEG Kit, Push	70	M00566340	EndoVive Through-the-PEG (TTP) Jejunal Feeding Tube – Bent Tip	73
M00563190	EndoVive One Step Button Low Profile PEG Kit, Push	70	M00566350	EndoVive Through-the-PEG (TTP) Jejunal Feeding Tube – Pigtail Tip	73
M00564101	eXcelon™ Transbronchial Aspiration Needle	89	M00566461	EndoVive Safety PEG Kit – Pull	69
M00564111	eXcelon Transbronchial Aspiration Needle	89	M00566471	EndoVive Safety PEG Kit – Push	69
M00564121	eXcelon Transbronchial Aspiration Needle	89	M00566481	EndoVive Safety PEG Kit – Pull	69
M00565010	WallFlex™ Duodenal Stent System with Anchor Lock Delivery System	43	M00566491	EndoVive Safety PEG Kit – Push	69
M00565020	WallFlex Duodenal Stent System with Anchor Lock Delivery System	43	M00566670	Encore™ 26 Inflation Device	35
M00565030	WallFlex Duodenal Stent System with Anchor Lock Delivery System	43	M00567330	MaxForce™ Biliary Balloon Dilatation Catheter	35
M00565040	WallFlex Colonic Stent System with Anchor Lock Delivery System	43	M00567340	MaxForce Biliary Balloon Dilatation Catheter	35
M00565050	WallFlex Colonic Stent System with Anchor Lock Delivery System	43	M00567370	MaxForce Biliary Balloon Dilatation Catheter	35
M00565060	WallFlex Colonic Stent System with Anchor Lock Delivery System	43	M00567380	MaxForce Biliary Balloon Dilatation Catheter	35
M00565070	WallFlex Colonic Stent System with Anchor Lock Delivery System	43	M00567410	MaxForce Biliary Balloon Dilatation Catheter	35
M00565080	WallFlex Colonic Stent System with Anchor Lock Delivery System	43	M00567420	MaxForce Biliary Balloon Dilatation Catheter	35
M00565090	WallFlex Colonic Stent System with Anchor Lock Delivery System	43	M00568201	EndoVive Standard PEG Kit – Pull	69
M00565100	WallFlex Colonic Stent System with Anchor Lock Delivery System	43	M00568211	EndoVive Standard PEG Kit – Push	69
M00565110	WallFlex Colonic Stent System with Anchor Lock Delivery System	43	M00568241	EndoVive Standard PEG Kit – Pull	69
M00565120	WallFlex Colonic Stent System with Anchor Lock Delivery System	43	M00568251	EndoVive Standard PEG Kit – Push	69
M00565130	WallFlex Colonic Stent System with Anchor Lock Delivery System	43	M00569610	WALLSTENT™ RX Biliary Endoprosthesis Stent System	17
M00565140	WallFlex Colonic Stent System with Anchor Lock Delivery System	43	M00569620	WALLSTENT RX Biliary Endoprosthesis Stent System	17
M00565150	WallFlex Colonic Stent System with Anchor Lock Delivery System	43	M00569630	WALLSTENT RX Biliary Endoprosthesis Stent System	17
M00565560	WALLSTENT™ Enteral Endoprosthesis Colonic / Duodenal Stent with UNISTEP™ Plus Delivery System	44	M00569640	WALLSTENT RX Biliary Endoprosthesis Stent System	17
M00565570	WALLSTENT Enteral Endoprosthesis Colonic / Duodenal Stent with UNISTEP Plus Delivery System	44	M00569650	WALLSTENT RX Biliary Endoprosthesis Stent System	17
M00565580	WALLSTENT Enteral Endoprosthesis Colonic / Duodenal Stent with UNISTEP Plus Delivery System	44	M00569660	WALLSTENT RX Biliary Endoprosthesis Stent System	17
M00565590	WALLSTENT Enteral Endoprosthesis Colonic / Duodenal Stent with UNISTEP Plus Delivery System	44	M00569670	WALLSTENT RX Biliary Endoprosthesis Stent System	17
M00565620	WALLSTENT Enteral Endoprosthesis Colonic / Duodenal Stent with UNISTEP Plus Delivery System	44	M00569680	WALLSTENT RX Biliary Endoprosthesis Stent System	17
M00565630	WALLSTENT Enteral Endoprosthesis Colonic / Duodenal Stent with UNISTEP Plus Delivery System	44	M00569690	WALLSTENT RX Biliary Endoprosthesis Stent System	17
M00565640	WALLSTENT Enteral Endoprosthesis Colonic / Duodenal Stent with UNISTEP Plus Delivery System	44	M00569700	WALLSTENT RX Biliary Endoprosthesis Stent System	17
M00565650	WALLSTENT Enteral Endoprosthesis Colonic / Duodenal Stent with UNISTEP Plus Delivery System	44	M00569710	WALLSTENT RX Biliary Endoprosthesis Stent System	17
			M00569720	WALLSTENT RX Biliary Endoprosthesis Stent System	17
			M00569800	WALLSTENT RX Biliary Endoprosthesis Stent System	17
			M00569810	WALLSTENT RX Biliary Endoprosthesis Stent System	17
			M00570000	Polyflex™ Single-Use Self-Expanding Silicone Airway Stent	88
			M00570010	Polyflex Single-Use Self-Expanding Silicone Airway Stent	88
			M00570020	Polyflex Single-Use Self-Expanding Silicone Airway Stent	88
			M00570030	Polyflex Single-Use Self-Expanding Silicone Airway Stent	88
			M00570040	Polyflex Single-Use Self-Expanding Silicone Airway Stent	88
			M00570050	Polyflex Single-Use Self-Expanding Silicone Airway Stent	88
			M00570060	Polyflex Single-Use Self-Expanding Silicone Airway Stent	88
			M00570070	Polyflex Single-Use Self-Expanding Silicone Airway Stent	88
			M00570080	Polyflex Single-Use Self-Expanding Silicone Airway Stent	88



Product Code Index

Product Code	Description	Page #	Product Code	Description	Page #
M00570090	Polyflex™ Single-Use Self-Expanding Silicone Airway Stent....	88	M00572240	Epic Biliary Stent	16
M00570100	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	88	M00572250	Epic Biliary Stent	16
M00570110	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	88	M00572260	Epic Biliary Stent	16
M00570120	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	88	M00572270	Epic Biliary Stent	16
M00570130	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	88	M00572280	Epic Biliary Stent	16
M00570140	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	88	M00572290	Epic Biliary Stent	16
M00570150	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	88	M00572300	Epic Biliary Stent	16
M00570160	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	88	M00572310	Epic Biliary Stent	16
M00570170	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	88	M00572320	Epic Biliary Stent	16
M00570180	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	88	M00572330	Epic Biliary Stent	16
M00570190	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	88	M00572340	Epic Biliary Stent	16
M00570200	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	88	M00573000	Jagtome™ RX 49 Sphincterotome (Preloaded with Jagwire™ Guidewire)	5
M00570210	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	88	M00573010	Jagtome RX 49 Sphincterotome (Preloaded with Jagwire Guidewire).....	5
M00570220	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	88	M00573020	Jagtome RX 49 Sphincterotome (Preloaded with Jagwire Guidewire).....	5
M00570230	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	88	M00573030	Jagtome RX 49 Sphincterotome (Preloaded with Jagwire Guidewire).....	5
M00570240	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	88	M00573040	Jagtome RX 44 Sphincterotome (Preloaded with Jagwire Guidewire).....	5
M00570250	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	88	M00573050	Jagtome RX 44 Sphincterotome (Preloaded with Jagwire Guidewire).....	5
M00570260	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	88	M00573060	Jagtome RX 44 Sphincterotome (Preloaded with Jagwire Guidewire).....	5
M00570270	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	88	M00573070	Jagtome RX 44 Sphincterotome (Preloaded with Jagwire Guidewire).....	5
M00570280	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	88	M00573080	Jagtome RX 39 Sphincterotome (Preloaded with Jagwire Guidewire).....	5
M00570290	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	88	M00573090	Jagtome RX 39 Sphincterotome (Preloaded with Jagwire Guidewire).....	5
M00570300	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	88	M00573100	Jagtome RX 39 Sphincterotome (Preloaded with Jagwire Guidewire).....	5
M00570310	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	88	M00573110	Jagtome RX 39 Sphincterotome (Preloaded with Jagwire Guidewire).....	5
M00570320	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	88	M00576300	Ultraflex™ Single-Use Uncovered Tracheobronchial Stent System – Distal Release	86
M00570330	Polyflex Single-Use Self-Expanding Silicone Airway Stent....	88	M00576310	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Distal Release	86
M00570340	WallFlex™ Biliary RX Fully Covered Stent System RMV	15	M00576320	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Distal Release	86
M00570350	WallFlex Biliary RX Fully Covered Stent System RMV	15	M00576330	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Distal Release	86
M00570360	WallFlex Biliary RX Fully Covered Stent System RMV	15	M00576340	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Distal Release	86
M00570370	WallFlex Biliary RX Fully Covered Stent System RMV	15	M00576350	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Distal Release	86
M00570380	WallFlex Biliary RX Fully Covered Stent System RMV	15	M00576360	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Distal Release	86
M00570600	WallFlex Biliary RX Uncovered Stent System.....	15	M00576370	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Distal Release	86
M00570610	WallFlex Biliary RX Uncovered Stent System.....	15	M00576380	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Distal Release	86
M00570620	WallFlex Biliary RX Uncovered Stent System.....	15			
M00570630	WallFlex Biliary RX Uncovered Stent System.....	15			
M00570640	WallFlex Biliary RX Uncovered Stent System.....	15			
M00570650	WallFlex Biliary RX Uncovered Stent System.....	15			
M00570660	WallFlex Biliary RX Uncovered Stent System.....	15			
M00570670	Dynamic™ (Y) Stent	85			
M00570680	Dynamic (Y) Stent	85			
M00570690	Dynamic (Y) Stent	85			
M00570700	WallFlex Biliary RX Partially Covered Stent System	15			
M00570710	WallFlex Biliary RX Partially Covered Stent System	15			
M00570720	WallFlex Biliary RX Partially Covered Stent System	15			
M00570730	WallFlex Biliary RX Partially Covered Stent System	15			
M00570740	WallFlex Biliary RX Partially Covered Stent System	15			
M00570890	WallFlex Biliary RX Uncovered Stent System.....	15			
M00572230	Epic™ Biliary Stent	16			



Product Code Index

Product Code	Description	Page #	Product Code	Description	Page #
M00576410	Ultraflex™ Single-Use Uncovered Tracheobronchial Stent System – Proximal Release	86	M00580171	EndoVive™ Low Profile Replacement Button Accessories – 18Fr Bolus Feeding Set, Straight Adapter	70, 72
M00576420	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Proximal Release	86	M00580181	EndoVive Low Profile Replacement Button Accessories – 18Fr Right Angle Feeding Set, Continuous Adapter	70, 72
M00576430	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Proximal Release	86	M00580231	EndoVive Low Profile Replacement Button Accessories – 24Fr Bolus Feeding Set, Straight Adapter	70, 72
M00576440	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Proximal Release	86	M00580241	EndoVive Low Profile Replacement Button Accessories – 24Fr Right Angle Feeding Set, Continuous Adapter	70, 72
M00576450	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Proximal Release	86	M00580401	EndoVive Low Profile Replacement Button Accessories – 18Fr Decompression Tube	70, 72
M00576460	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Proximal Release	86	M00580411	EndoVive Low Profile Replacement Button Accessories – 18Fr Decompression Tube	70, 72
M00576470	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Proximal Release	86	M00580421	EndoVive Low Profile Replacement Button Accessories – 18Fr Decompression Tube	70, 72
M00576480	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Proximal Release	86	M00580431	EndoVive Low Profile Replacement Button Accessories – 18Fr Decompression Tube	70, 72
M00576490	Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Proximal Release	86	M00580441	EndoVive Low Profile Replacement Button Accessories – 18Fr Decompression Tube	70, 72
M00576520	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release	87	M00580461	EndoVive Low Profile Replacement Button Accessories – 24Fr Decompression Tube	70, 72
M00576530	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release	87	M00580471	EndoVive Low Profile Replacement Button Accessories – 24Fr Decompression Tube	70, 72
M00576540	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release	87	M00580481	EndoVive Low Profile Replacement Button Accessories – 24Fr Decompression Tube	70, 72
M00576550	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release	87	M00580491	EndoVive Low Profile Replacement Button Accessories – 24Fr Decompression Tube	70, 72
M00576560	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release	87	M00580641	EndoVive Initial Placement PEG Accessories – Y-Ports	69
M00576570	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release	87	M00580651	EndoVive Initial Placement PEG Accessories – Y-Ports	69
M00576580	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release	87	M00580731	EndoVive Initial Placement PEG Accessories – Round Bolsters	65, 69
M00576590	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release	87	M00580781	EndoVive Initial Placement PEG Accessories – C Clamps	65, 69
M00576600	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release	87	M00582040	EndoVive Replacement Gastrostomy Feeding Tubes, Straight Bolster	71
M00576610	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release	87	M00582050	EndoVive Replacement Gastrostomy Feeding Tubes, Straight Bolster	71
M00576620	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release	87	M00582060	EndoVive Replacement Gastrostomy Feeding Tubes, Straight Bolster	71
M00576630	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release	87	M00582070	EndoVive Replacement Gastrostomy Feeding Tubes, Straight Bolster	71
M00576640	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release	87	M00582080	EndoVive Replacement Gastrostomy Feeding Tubes, Straight Bolster	71
M00576650	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release	87	M00582090	EndoVive Replacement Gastrostomy Feeding Tubes, Straight Bolster	71
M00576660	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release	87	M00582100	EndoVive Replacement Gastrostomy Feeding Tubes, Straight Bolster	71
M00576670	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release	87	M00582110	EndoVive Replacement Gastrostomy Feeding Tubes, Right Angle Bolster	71
M00576680	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release	87	M00582120	EndoVive Replacement Gastrostomy Feeding Tubes, Right Angle Bolster	71
M00576690	Ultraflex Single-Use Partially Covered Tracheobronchial Stent System – Distal Release	87	M00582130	EndoVive Replacement Gastrostomy Feeding Tubes, Right Angle Bolster	71
			M00582140	EndoVive Replacement Gastrostomy Feeding Tubes, Right Angle Bolster	71



Product Code Index

Product Code	Description	Page #	Product Code	Description	Page #
M00582150	EndoVive™ Replacement Gastrostomy Feeding Tubes, Right Angle Bolster	71	M00584000	Dreamtome™ RX 49 Sphincterotome (Preloaded with Dreamwire High Performance Guidewire)	4
M00582160	EndoVive Replacement Gastrostomy Feeding Tubes, Right Angle Bolster	71	M00584010	Dreamtome RX 49 Sphincterotome (Preloaded with Dreamwire High Performance Guidewire)	4
M00582170	EndoVive Replacement Gastrostomy Feeding Tubes, Right Angle Bolster	71	M00584020	Dreamtome RX 49 Sphincterotome (Preloaded with Dreamwire High Performance Guidewire)	4
M00582860	EndoVive Replacement Gastrostomy Feeding Tubes, Straight Bolster	71	M00584030	Dreamtome RX 49 Sphincterotome (Preloaded with Dreamwire High Performance Guidewire)	4
M00582870	EndoVive Replacement Gastrostomy Feeding Tubes, Right Angle Bolster	71	M00584040	Dreamtome RX 44 Sphincterotome (Preloaded with Dreamwire High Performance Guidewire)	4
M00583000	Hydratome™ RX 49 Sphincterotome (Preloaded with Hydra Jagwire™ Guidewire)	4	M00584050	Dreamtome RX 44 Sphincterotome (Preloaded with Dreamwire High Performance Guidewire)	4
M00583010	Hydratome RX 49 Sphincterotome (Preloaded with Hydra Jagwire Guidewire)	4	M00584060	Dreamtome RX 44 Sphincterotome (Preloaded with Dreamwire High Performance Guidewire)	4
M00583020	Hydratome RX 49 Sphincterotome (Preloaded with Hydra Jagwire Guidewire)	4	M00584070	Dreamtome RX 44 Sphincterotome (Preloaded with Dreamwire High Performance Guidewire)	4
M00583030	Hydratome RX 49 Sphincterotome (Preloaded with Hydra Jagwire Guidewire)	4	M00584120	TRUEtome 39 Sphincterotome	8
M00583040	Hydratome RX 44 Sphincterotome (Preloaded with Hydra Jagwire Guidewire)	4	M00584150	TRUEtome 49 Sphincterotome	8
M00583050	Hydratome RX 44 Sphincterotome (Preloaded with Hydra Jagwire Guidewire)	4	M00584160	TRUEtome 49 Sphincterotome	8
M00583060	Hydratome RX 44 Sphincterotome (Preloaded with Hydra Jagwire Guidewire)	4	M00584170	TRUEtome 44 Sphincterotome	8
M00583070	Hydratome RX 44 Sphincterotome (Preloaded with Hydra Jagwire Guidewire)	4	M00584180	TRUEtome 44 Sphincterotome	8
M00583080	TRUEtome™ Jag 39 Cannulating Sphincterotome (Preloaded with High Performance Jagwire™ Guidewire)	7	M00584190	TRUEtome 39 Sphincterotome	8
M00583090	TRUEtome Jag 44 Cannulating Sphincterotome (Preloaded with High Performance Jagwire Guidewire)	7	M00589060	MaxForce TTS™ Esophageal High Performance Balloon Dilatation Catheter	59
M00583100	TRUEtome Jag 44 Cannulating Sphincterotome (Preloaded with High Performance Jagwire Guidewire)	7	M00589080	MaxForce TTS Esophageal High Performance Balloon Dilatation Catheter	59
M00583110	TRUEtome Jag 44 Cannulating Sphincterotome (Preloaded with High Performance Jagwire™ Guidewire)	7	M00589100	MaxForce TTS Esophageal High Performance Balloon Dilatation Catheter	59
M00583120	TRUEtome Jag 39 Cannulating Sphincterotome (Preloaded with High Performance Jagwire Guidewire)	7	M00589120	MaxForce TTS Esophageal High Performance Balloon Dilatation Catheter	59
M00583130	TRUEtome Jag 44 Cannulating Sphincterotome (Preloaded with High Performance Jagwire Guidewire)	7	M00589140	MaxForce TTS Esophageal High Performance Balloon Dilatation Catheter	59
M00583140	TRUEtome Jag 44 Cannulating Sphincterotome (Preloaded with High Performance Jagwire Guidewire)	7	M00589150	MaxForce TTS Esophageal High Performance Balloon Dilatation Catheter	59
M00583150	TRUEtome Jag 44 Cannulating Sphincterotome (Preloaded with High Performance Jagwire Guidewire)	7	M00589160	MaxForce TTS Esophageal High Performance Balloon Dilatation Catheter	59
M00583160	TRUEtome Hydra 44 Cannulating Sphincterotome (Preloaded with High Performance Hydra Guidewire)	7	M00589180	MaxForce TTS Esophageal High Performance Balloon Dilatation Catheter	59
M00583170	TRUEtome Hydra 44 Cannulating Sphincterotome (Preloaded with High Performance Hydra Guidewire)	7	M005ATS20000	Alair™ Bronchial Thermoplasty System Controller	85
M00583180	TRUEtome Dreamwire 44 Cannulating Sphincterotome (Preloaded with High Performance Dreamwire™ Guidewire)	7	M005ATS25010	Alair Bronchial Thermoplasty System Catheter	85
M00583190	TRUEtome Dreamwire 44 Cannulating Sphincterotome (Preloaded with High Performance Dreamwire Guidewire)	7	M005ATS20110	Alair RF Controller Accessory Kit	85
M00583200	TRUEtome Dreamwire 44 Cannulating Sphincterotome (Preloaded with High Performance Dreamwire Guidewire)	7	M005XXXX1	RX Cholangiogram Kit-Preloaded – Autotome™ Sphincterotome – RX Locking Device	6
M00583210	TRUEtome Dreamwire 44 Cannulating Sphincterotome (Preloaded with High Performance Dreamwire Guidewire)	7	M005XXXX1	RX Cholangiogram Kits – Autotome or Ultratome™ XL Sphincterotome	6
			SBB-119-50	HedgeHog™ Valve Brush	81
			SBB-382-50	HedgeHog Double Bundle (Double-End Brush + Valve Brush)	81
			SBD-228-50	HedgeHog Double-End Channel Brush	81
			SBD-289-50	HedgeHog Dual-End Channel / Valve Brush	81
			SBD-291-50	HedgeHog Sweeper Channel Brush	81
			SBP-236-40	Blox™ Endoscopic Bite Block with Elastic Cloth Strap (38 FR)	81



Product Code Index

Product Code	Description	Page #	Product Code	Description	Page #
SBS-227-50	HedgeHog™ Single-End Channel Brush.....	81	SIT-470-100	Hydra Disposable Auxiliary Water Jet Connector.....	82
SBT-114-100	Blox™ Endoscopic Bite Block with Elastic Cloth Strap (54 FR) ..	81	SIT-576-15	Hydra Irrigation Tubing (Compatible with Olympus Pumps) ...	82
SBT-321-100	Blox Endoscopic Bite Block with Elastic Cloth Strap (54 FR) (w/o Dental Retention).....	81	SLT-394-100	Endoglide™+ Sterile Bacteriostatic Lubricating Gel.....	82
SBT-546-100	Blox Endoscopic Bite Block with Elastic Cloth Strap (60 FR).....	81	SLT-395-50	Endoglide+ Sterile Bacteriostatic Lubricating Gel.....	82
SCT-350-10	Hydra™ Water Bottle Cap.....	82	SPX-017-25	CinchPad™ Endoscope Transport System.....	83
SCT-351-10	Hydra Water Bottle Cap.....	82	SSB-210-50	Koala™ Cleaning Sponge.....	83
SCT-466-10	Hydra Water Bottle Cap.....	82	SSB-240-50	Koala Cleaning Sponge Plus Pure Enzymatic	83
SCT-467-10	Hydra Water Bottle Cap.....	82	STE-297-50	TrapEase™ Single Chamber Polyp Trap	80
SEE-220-4	Pure™ Enzymatic Detergent - Spearmint Scent.....	81	STE-298-50	TrapEase Quad Chamber Polyp Trap	80
SEE-573-4	Pure Enzymatic Detergent	81	SUV-617-50	Orca™ Sterile Single-Use Air/Water & Suction Valves.....	83
SGE-368-100	EndoArmor™ Protective Gown.....	82			
SGE-389-100	EndoArmor+ Protective Gown.....	82			
SIK-242-24	SafeStart™ Cleaning Kit.....	83			
SIK-243-24	SafeStart CP Cleaning Kit.....	83			
SIT-355-15	Hydra Irrigation Tubing (Compatible with ERBE, Medivators Pumps).....	82			
SIT-361-100	Hydra Disposable Auxiliary Water Jet Connector.....	82			



Notes and Future Product Updates



Terms & Conditions



Terms & Conditions

🕒 Begins on Page..... 113

Certain Terms and Conditions

Ordering Information

All orders should be placed directly with Boston Scientific Endoscopy Customer Service
by mail, telephone, fax or Electronic Data Interchange (EDI).

Boston Scientific Endoscopy
500 Commander Shea Blvd.
Quincy, MA 02171
(888) 272-1001 Phone
(888) 272-9444 Fax
(877) 272-3340 EDI

Please contact EDI if interested in ordering your endoscopy products via the internet.

Please note that the above address is our mailing address only.

Customer Service is available between the hours of 8:30 am and 9:30 pm (Eastern Standard Time) weekdays excluding national holidays. Our toll-free fax number, (888) 272-9444, is available at all times. In addition, technical and emergency service is available 24 hours a day, 7 days a week, call (888) 272-1001 for further instructions.

Standing Orders

Boston Scientific Endoscopy accepts standing orders on all products. This convenience provides you with scheduled shipments using a single purchase order. Contact the Customer Service Database Department at (888) 325-1060 for details.

Terms

Terms are Net 30 Days. All shipments are F.O.B. Quincy, Massachusetts. Prices are subject to change without notice.

Returned Goods

Permission to return merchandise is required for all returns and, except as provided herein, may be withheld in the sole discretion of Boston Scientific Endoscopy. This authorization must be requested by calling or writing:

Customer Service

Boston Scientific Endoscopy

500 Commander Shea Blvd.

Quincy, MA 02171

(800) 862-1284

The request for authorization must include:

- Reason for the return – “Ordered in Error,” “Shipped in Error,” etc.
- Catalog number, size, quantity, lot number
- Purchase order number, Boston Scientific Endoscopy invoice number and date of invoice

All merchandise returned to Boston Scientific Endoscopy must be accompanied by a copy of the packing slip or invoice enclosed with the original shipment. The Returned Goods Authorization (RGA) number must be clearly marked on the shipping container. Merchandise returned without an RGA number will not be accepted or considered as returned product.

Merchandise that will not be accepted for return:

- Expired product or products with less than 6 months shelf life
- Merchandise manufactured as SPECIALS or CUSTOM
- Merchandise in unsaleable condition (opened, damaged, tampered, etc.)
- Merchandise not listed in this Price List or any Supplement
- Merchandise sterilized or used by customer
- Merchandise that is not packaged properly for shipping when returned
- Incomplete kits or boxes (unless defective)
- Unauthorized returns



Certain Terms and Conditions

Returned Goods (continued)

Allowance Schedule for Authorized Returns:

Allowance Schedule for Capital Equipment:

- Capital Equipment processed in error by Boston Scientific Endoscopy will receive full credit if returned in saleable condition
- Capital Equipment cannot be returned once the customer has taken title of the equipment

Allowance Schedule for Metal and Polyflex™ Stents and Alair™ Bronchial Thermoplasty Products:

- Merchandise processed in error by Boston Scientific Endoscopy will receive full credit if returned in saleable condition
- Merchandise returned within 30 days of invoice date in saleable condition will receive 100% credit
- Merchandise returned from 31 days to 90 days in saleable condition will receive full credit less a 30% restocking charge
- Merchandise returned from 91 days to 150 days in saleable condition will receive full credit less a 50% restocking fee
- Merchandise returned past 150 days will not receive any credit

Allowance Schedule for All Other Products:

- Merchandise processed in error by Boston Scientific Endoscopy will receive full credit if returned in saleable condition
- Other merchandise returned within 30 days of invoice date in saleable condition will receive full credit less a 15% restocking charge
- Merchandise returned past 31 days in saleable condition will receive full credit less a 40% restocking charge provided there is at least 6 months left on the expiration date

Pricing and Specification Changes

As we are continually introducing products and improving our existing products, this Price List does not purport to be complete, and the equipment supplied may not agree in all details with the description and specifications shown in the catalog and / or price list. All product offerings, prices and specifications are subject to change without notice.

Notice

Indications, contraindications and instructions for use, as well as product warranty information, if applicable, can be found in the product labeling supplied with each device. Caution: Federal law restricts this device to sale by or on the order of a physician.

Warning Re: Single-Use Devices

Products marked as “single-use” are for single use only. Do not reuse, reprocess or resterilize single-use products. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and / or lead to device failure which in turn may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and / or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. **BSC assumes no liability with respect to single-use devices reused, reprocessed or resterilized and makes no warranties, express or implied, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose, with respect to any such instrument.**

©2018 Boston Scientific Corporation or its Affiliates. All rights reserved.

INDICATIONS FOR USE

The WallFlex Biliary RX Fully Covered Stent System RMV is indicated for use in the palliative treatment of biliary strictures produced by malignant neoplasms, relief of malignant biliary obstruction prior to surgery and for indwell up to 12 months in the treatment of benign biliary strictures secondary to chronic pancreatitis.

LIMITATIONS

The sale, distribution, and use of the device are restricted to prescription use in accordance with 21 CFR 801.109.

Contraindications:

- The WallFlex Biliary RX Fully Covered Stent should not be placed in strictures that cannot be dilated enough to pass the delivery system, in a perforated duct, or in very small intrahepatic ducts.
- The WallFlex Biliary RX Fully Covered Stent System RMV should not be used in patients for whom endoscopic techniques are contraindicated.

Warnings:

- The safety and effectiveness of the stent has not been established for indwell periods exceeding 12 months, when used in the treatment of benign strictures secondary to chronic pancreatitis.
- The WallFlex Biliary RX Fully Covered Stent System RMV is for single-use only.
- The safety and effectiveness of the WallFlex Biliary RX Fully Covered Stent System RMV for use in the vascular system has not been established.
- The safety and effectiveness of the WallFlex Biliary RX Fully Covered Stent System RMV has not been established in the treatment of benign biliary anastomotic strictures in liver transplant patients and benign biliary post abdominal surgery strictures.
- Testing of overlapped stents has not been conducted.
- The stent contains nickel, which may cause an allergic reaction in individuals with nickel sensitivity.

PLEASE REFER TO THE LABELING FOR A MORE COMPLETE LIST OF WARNINGS, PRECAUTIONS AND CONTRAINDICATIONS.



Customer Ordering Information

Universal Product Number		
Radial Jaw™ 3 Single-Use Max Capacity Biopsy F		
Order Number	GTIN	Description
M005 15861	08714729096092	Radial Jaw 3 Max Capac
M005 15871	08714729085799	Radial Jaw 3 Max Capac
M005 15881	08714729085805	Radial Jaw 3 Max Capac
M005 15891	08714729096108	Radial Jaw 3 Max Capac
M005 15892	08714729718802	Radial Jaw 3 Max Capac

Bold type
is old order
number

For customer convenience, the Boston Scientific Endoscopy ordering system accepts the UPN (Universal Product Number) e.g., **M00515892** found on product labeling as well as Boston Scientific's old 4 digit order numbers e.g., **1589**.

Please use either one of these numbers when placing an order via phone, fax, internet, or Electronic Data Interchange (EDI).

For those customers who scan product today, Boston Scientific is moving to the GS1 GTIN barcode and labeling standards as part of our ongoing commitment to improve supply chain efficiency. For more information go to the website: www.bostonscientific.com/GS1.

Boston Scientific Endoscopy Customer Service Representatives will be happy to help you with any questions. We are eager to make the ordering procedure as simple and efficient as possible.

All trademarks are the property of their respective owners.



**Magnetic Resonance
Conditional¹**

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

Indications, Contraindications, Warnings and Instructions for Use can be found in the product labeling supplied with each device.
CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

**Boston
Scientific**
Advancing science for life™

Boston Scientific Corporation
300 Boston Scientific Way
Marlborough, MA 01752-1234
www.bostonscientific.com

Ordering Information
1.888.272.1001

©2018 Boston Scientific Corporation
or its affiliates. All rights reserved.

ENDO-519003-AA E February 2018