

Enteral & Esophageal Stents



Stents (Enteral and Esophageal)

Colonic and Duodenal Stents

➤ Ultraflex™ Precision Single-Use Colonic Stent System	44
➤ WallFlex™ Colonic Stent System with Anchor Lock Delivery System	43
➤ WallFlex Duodenal Stent System with Anchor Lock Delivery System	43
➤ WALLSTENT™ Enteral Endoprosthesis Colonic / Duodenal Stent with UNISTEP™ Plus Delivery System	44
➤ WALLSTENT Super Stiff Guidewire	44

Esophageal Stents

➤ Amplatz Super Stiff™ Guidewire	47
➤ Polyflex™ Single-Use Esophageal Stent System	46
➤ Ultraflex Single-Use Covered Esophageal NG Stent System – Distal Release	46
➤ Ultraflex Single-Use Covered Esophageal NG Stent System – Proximal Release	46
➤ Ultraflex Single-Use Covered Large Esophageal NG Stent System – Distal Release	47
➤ Ultraflex Single-Use Covered Large Esophageal NG Stent System – Proximal Release	47
➤ Ultraflex Single-Use Uncovered Esophageal NG Stent System – Distal Release	47
➤ Ultraflex Single-Use Uncovered Esophageal NG Stent System – Proximal Release	47
➤ WallFlex Fully Covered Esophageal Stent	45
➤ WallFlex Partially Covered Esophageal Stent	45



Enteral Stents

WallFlex™ Colonic and Duodenal Stents



 Magnetic Resonance Conditional¹

WallFlex Colonic Stent System with Anchor Lock Delivery System*



Order Number	GTIN	Stent Body Diameter (mm)	Stent Flare Diameter (mm)	Expanded Stent Length (mm)	Working Length (cm)	Delivery System Length (cm)	Delivery System Diameter (Fr)	Packaging
M00565040	08714729456513	25	30	60	230	270	10	Each
M00565050	08714729456520	25	30	90	230	270	10	Each
M00565060	08714729456537	25	30	120	230	270	10	Each
M00565070	08714729456544	25	30	60	135	175	10	Each
M00565080	08714729456551	25	30	90	135	175	10	Each
M00565090	08714729456568	25	30	120	135	175	10	Each
M00565100	08714729456575	22	27	60	230	270	10	Each
M00565110	08714729456582	22	27	90	230	270	10	Each
M00565120	08714729456599	22	27	120	230	270	10	Each
M00565130	08714729456605	22	27	60	135	175	10	Each
M00565140	08714729456612	22	27	90	135	175	10	Each
M00565150	08714729456629	22	27	120	135	175	10	Each

*Recommended Guidewire .035" WALLSTENT™ Super Stiff Guidewire, Order #18001,  (See page 44).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



 Magnetic Resonance Conditional¹

WallFlex Duodenal Stent System with Anchor Lock Delivery System*



Order Number	GTIN	Stent Body Diameter (mm)	Flare Diameter (mm)	Expanded Stent Length (mm)	Working Length (cm)	Delivery System Length (cm)	Delivery System Diameter (Fr)	Packaging
M00565010	08714729456483	22	27	60	230	270	10	Each
M00565020	08714729456490	22	27	90	230	270	10	Each
M00565030	08714729456506	22	27	120	230	270	10	Each

*Recommended Guidewire .035" Hydra Jagwire™ Guidewire Stiff Shaft 450cm Order #5606,  (See page 1).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



Section Table of Contents

Prev Page

43



Next Page

Main Table of Contents



Enteral Stents

WALLSTENT™ Colonic and Duodenal Stents (Continued)



 Magnetic Resonance Conditional¹

WALLSTENT Enteral Endoprosthesis Colonic / Duodenal Stent with UNISTEP™ Plus Delivery System*



Order Number	GTIN	Expanded Stent Diameter (mm)	Expanded Stent Length (mm)	Working Length (cm)	Delivery System Length (cm)	Delivery System Diameter (Fr)	Packaging
M00565560	08714729285052	20	60	230	255	10	Each
M00565570	08714729285069	20	90	230	255	10	Each
M00565580	08714729285076	22	60	230	255	10	Each
M00565590	08714729285083	22	90	230	255	10	Each
M00565620	08714729285113	20	60	135	160	10	Each
M00565630	08714729285120	20	90	135	160	10	Each
M00565640	08714729285137	22	60	135	160	10	Each
M00565650	08714729285144	22	90	135	160	10	Each


*Recommended Guidewire .035" WALLSTENT Super Stiff Guidewire, Order #18001, (see below).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

WALLSTENT Super Stiff Guidewire



Order Number	GTIN	Description	Distal Tip O.D. (in)	Length (cm)	Tip Style	Packaging
H965180011	08714729203070	WALLSTENT Super Stiff Guidewire	.035	500	Tapered	Box 5



 Magnetic Resonance Conditional¹

Ultraflex™ Precision Single-Use Colonic Stent System*



Order Number	GTIN	Stent Body Diameter (mm)	Stent Flare Diameter (mm)	Expanded Stent Length (mm)	Working Length (cm)	Delivery System Length (cm)	Delivery Tip Diameter (Fr)	Packaging
M00557330	08714729339175	25	30	57	100	105	22	Each
M00557360	08714729339199	25	30	87	100	105	22	Each
M00557380	08714729339205	25	30	117	100	105	22	Each

*Recommended Guidewire .038" Jagwire™ Guidewire, Order #5662,  (See page 2).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



Section Table of Contents


Prev Page  44  Next Page

Main Table of Contents



Esophageal Stents

WallFlex™ Esophageal Stents



 Magnetic Resonance Conditional¹


WallFlex Fully Covered Esophageal Stent*



Order Number	GTIN	Stent O.D. (mm)	Proximal / Distal Flares O.D. (mm)	Stent Length (cm)	Catheter Diameter (Fr / mm)	Working Length (cm)	Packaging
M00516700	08714729778035	18	25 / 23	10.3	18.5 / 6.2	78	Each
M00516710	08714729778042	18	25 / 23	12.3	18.5 / 6.2	78	Each
M00516720	08714729778059	18	25 / 23	15.3	18.5 / 6.2	78	Each
M00516730	08714729778066	23	28 / 28	10.5	18.5 / 6.2	78	Each
M00516740	08714729778073	23	28 / 28	12.5	18.5 / 6.2	78	Each
M00516750	08714729778080	23	28 / 28	15.5	18.5 / 6.2	78	Each

*Recommended Guidewire .038" Jagwire™ Guidewire, Order #5662,  (See page 2).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



 Magnetic Resonance Conditional¹

WallFlex Partially Covered Esophageal Stent*



Order Number	GTIN	Stent O.D. (mm)	Proximal / Distal Flares O.D. (mm)	Stent Length (cm)	Covered Length (cm)	Catheter Diameter (Fr / mm)	Working Length (cm)	Packaging
M00516900	08714729765240	18	23 / 23	10.3	7.1	18.5 / 6.2	78	Each
M00516910	08714729765257	18	23 / 23	12.3	9.1	18.5 / 6.2	78	Each
M00516920	08714729765264	18	23 / 23	15.3	12.1	18.5 / 6.2	78	Each
M00516930	08714729765271	23	28 / 28	10.5	7.3	18.5 / 6.2	78	Each
M00516940	08714729765288	23	28 / 28	12.5	9.3	18.5 / 6.2	78	Each
M00516950	08714729765295	23	28 / 28	15.5	12.3	18.5 / 6.2	78	Each

*Recommended Guidewire .038" Jagwire Guidewire, Order #5662,  (See page 2).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



Section Table of Contents

Prev Page

45

Next Page


Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Esophageal Stents

Polyflex™ Stents




Polyflex Single-Use Esophageal Stent System*



Order Number	GTIN	Stent I.D. (mm)	Proximal Flare I.D. (mm)	Stent Length (mm)	Delivery System Diameter (mm)	Packaging
M00514270	08714729457015	16	20	90	12	Each
M00514280	08714729457022	16	20	120	12	Each
M00514290	08714729457039	16	20	150	12	Each
M00514300	08714729457046	18	23	90	13	Each
M00514310	08714729457053	18	23	120	13	Each
M00514320	08714729457060	18	23	150	13	Each
M00514330	08714729457077	21	25	90	14	Each
M00514340	08714729457084	21	25	120	14	Each
M00514350	08714729457091	21	25	150	14	Each

*Recommended Guidewire .038" Jagwire™ Guidewire, Order #5662, (See page 2).

Ultraflex™ Stents



Magnetic Resonance Conditional¹

Ultraflex Single-Use Covered Esophageal NG Stent System – Distal Release*



Order Number	GTIN	Stent O.D. (mm)	Proximal Flare O.D. (mm)	Stent Length (cm)	Covered Length (cm)	Packaging
M00513730	08714729716150	18	23	10	7	Each
M00513740	08714729716167	18	23	12	9	Each
M00513750	08714729649113	18	23	15	12	Each

*Recommended Guidewire .038" Jagwire Guidewire, Order #5662, (See page 2).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

Magnetic Resonance Conditional¹

Ultraflex Single-Use Covered Esophageal NG Stent System – Proximal Release*



Order Number	GTIN	Stent O.D. (mm)	Proximal Flare O.D. (mm)	Stent Length (cm)	Covered Length (cm)	Packaging
M00513840	08714729716099	18	23	10	7	Each
M00513850	08714729716075	18	23	12	9	Each
M00513860	08714729716082	18	23	15	12	Each

*Recommended Guidewire .038" Jagwire Guidewire, Order #5662, (See page 2).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



Section Table of Contents

Prev Page 46 Next Page


Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Esophageal Stents

Ultraflex™ Stents (Continued)



 Magnetic Resonance Conditional¹

Ultraflex Single-Use Covered Large Esophageal NG Stent System – Distal Release*



Order Number	GTIN	Stent O.D. (mm)	Proximal Flare O.D. (mm)	Stent Length (cm)	Covered Length (cm)	Packaging
M00514200	08714729121190	23	28	10	7	Each
M00514210	08714729121206	23	28	12	9	Each

*Recommended Guidewire .038" Jagwire™ Guidewire, Order #5662,  (See page 2).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

 Magnetic Resonance Conditional¹


Ultraflex Single-Use Covered Large Esophageal NG Stent System – Proximal Release*



Order Number	GTIN	Stent O.D. (mm)	Proximal Flare O.D. (mm)	Stent Length (cm)	Covered Length (cm)	Packaging
M00514240	08714729716105	23	28	10	7	Each
M00514250	08714729716112	23	28	12	9	Each

*Recommended Guidewire .038" Jagwire Guidewire, Order #5662,  (See page 2).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.



 Magnetic Resonance Conditional¹

Ultraflex Single-Use Uncovered Esophageal NG Stent System – Distal Release*



Order Number	GTIN	Stent O.D. (mm)	Proximal Flare O.D. (mm)	Stent Length (cm)	Packaging
M00513700	08714729716129	18	23	7	Each
M00513710	08714729716136	18	23	10	Each
M00513720	08714729716143	18	23	15	Each

*Recommended Guidewire .038" Jagwire Guidewire, Order #5662,  (See page 2).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

 Magnetic Resonance Conditional¹

Ultraflex Single-Use Uncovered Esophageal NG Stent System – Proximal Release*



Order Number	GTIN	Stent O.D. (mm)	Proximal Flare O.D. (mm)	Stent Length (cm)	Packaging
M00513800	08714729716174	18	23	7	Each
M00513810	08714729716181	18	23	10	Each
M00513830	08714729716204	18	23	15	Each

*Recommended Guidewire .038" Jagwire Guidewire, Order #5662,  (See page 2).

¹Product with this symbol may be used in Magnetic Resonance environment according to the conditions described in Directions for Use.

 Magnetic Resonance Conditional¹

Amplatz Super Stiff™ Guidewire



Order Number	GTIN	Description	O.D. (in)	Length (cm)	Tip Style	Packaging
M00550090	08714729301929	Amplatz Super Stiff Guidewire	.038	260	Straight	Each