

Pulmonary Endoscopy

Alair™ Bronchial Thermoplasty System

- Alair Bronchial Thermoplasty System 85

Airway Stents

- Dynamic™ (Y) Stent..... 85
- Polyflex™ Single-Use Self-Expanding Silicone Airway Stent 88
- Ultraflex™ Single-Use Partially Covered Tracheobronchial Stent System – Distal Release 87
- Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Distal Release..... 86
- Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Proximal Release 86

Biopsy Forceps

- Radial Jaw™ 4 Single-Use Pulmonary Biopsy Forceps..... 89

Endobronchial Ultrasound Transbronchial Aspiration Needle (EBUS)

- Expect™ Pulmonary Endobronchial Ultrasound Transbronchial Aspiration Needle / Adaptor 89

Transbronchial Aspiration Needles

- eXcelon™ Transbronchial Aspiration Needle 89

Cytology Brushes

- Cellebri™ Single-Use Endoscopic Cytology Brush..... 90

Airway Retrieval Baskets

- Zero Tip™ Airway Retrieval Basket 90

Controlled Radial Expansion Pulmonary Balloon Dilatation Catheter

- CRE™ Pulmonary Balloon Dilatation Catheter..... 91

Pulmonary Guidewires


- Amplatz Super Stiff Guidewire 91
- Pulmonary Jagwire™ Pulmonary Guidewire 91

Inflation System

- Alliance™ II Integrated Inflation / Lithotripsy Device..... 91
- Alliance II Single-Use Syringe / Gauge Assembly 91

Pulmonary

Alair™ Bronchial Thermoplasty System



Alair Bronchial Thermoplasty System*


Order Number	GTIN	Model Number	Description	Packaging
M005ATS20000	08714729796695	Alair ATS 200	Alair Bronchial Thermoplasty Controller	Each
M005ATS25010	08714729802792	Alair ATS 2-5	⊗ Alair Bronchial Thermoplasty Catheter**	Each
M005ATS20110	08714729802778	Alair ATS 201	Alair RF Controller Accessory Kit	Each

*Service Agreement available with 1, 2, 3, and 4 year terms.

**Note: Initial stocking order requires a minimum order of 6 catheters (covering the complete treatment of 2 patients).

UPN	Description
SRVM005ATS1001	ENDO POS 1 Year Service Alair
SRVM005ATS1002	ENDO POS 2 Year Service Alair
SRVM005ATS1003	ENDO POS 3 Year Service Alair
SRVM005ATS1004	ENDO POS 4 Year Service Alair

Airway Stents



Dynamic™ (Y) Stent



Order Number	GTIN	Tracheal Width (ID) (mm)	Bronchial Width (ID) (mm)	Tracheal Length (mm)	Bronchial Lengths (R / L) (mm)	Packaging
M00570670	08714729764809	11	8	110	25 / 40	Each
M00570680	08714729764816	13	10	110	25 / 40	Each
M00570690	08714729764823	15	12	110	25 / 40	Each



Section Table of Contents

Prev Page

85


Next Page

Main Table of Contents



Pulmonary

Airway Stents (Continued)



 Magnetic Resonance Conditional¹

Ultraflex™ Single-Use Uncovered Tracheobronchial Stent System – Distal Release*



Order Number	GTIN	Description	Expanded Stent O.D. (mm)	Expanded Stent Length (mm)	Compressed Stent O.D. (Fr)**	Tip Max O.D. (mm)	Packaging
M00576300	08714729837688	Uncovered Distal Release	10	30	15	4.1	Each
M00576310	08714729837695	Uncovered Distal Release	12	30	17	4.1	Each
M00576320	08714729837701	Uncovered Distal Release	14	30	18	4.1	Each
M00576330	08714729837718	Uncovered Distal Release	16	40	20	5.3	Each
M00576340	08714729837725	Uncovered Distal Release	16	60	20	5.3	Each
M00576350	08714729837732	Uncovered Distal Release	18	40	20	5.3	Each
M00576360	08714729837749	Uncovered Distal Release	18	60	20	5.3	Each
M00576370	08714729837756	Uncovered Distal Release	20	40	20	5.3	Each
M00576380	08714729837763	Uncovered Distal Release	20	60	20	5.3	Each

*Recommended Guidewire .035" Pulmonary Jagwire™ Pulmonary Guidewire and Amplatz Super Stiff Guidewire, Order #1517,  (See page 91).
Delivery System working length is 95cm.

**Data on File at Boston Scientific.



 Magnetic Resonance Conditional¹

Ultraflex Single-Use Uncovered Tracheobronchial Stent System – Proximal Release*



Order Number	GTIN	Description	Expanded Stent O.D. (mm)	Expanded Stent Length (mm)	Compressed Stent O.D. (Fr)**	Tip Max O.D. (mm)	Packaging
M00576410	08714729837770	Uncovered Proximal Release	8	20	15	4.1	Each
M00576420	08714729837787	Uncovered Proximal Release	8	40	15	4.1	Each
M00576430	08714729837794	Uncovered Proximal Release	10	20	15	4.1	Each
M00576440	08714729837800	Uncovered Proximal Release	10	40	15	4.1	Each
M00576450	08714729837817	Uncovered Proximal Release	12	20	17	4.1	Each
M00576460	08714729837824	Uncovered Proximal Release	12	40	17	4.1	Each
M00576470	08714729837831	Uncovered Proximal Release	14	20	18	4.1	Each
M00576480	08714729837848	Uncovered Proximal Release	14	40	18	4.1	Each
M00576490	08714729837855	Uncovered Proximal Release	14	60	18	4.1	Each

*Recommended Guidewire .035" Pulmonary Jagwire Pulmonary Guidewire and Amplatz Super Stiff Guidewire, Order #1517,  (See page 91).
Delivery System working length is 95cm.

**Data on File at Boston Scientific.



Section Table of Contents


Prev Page  86  Next Page

Main Table of Contents



Pulmonary

Airway Stents (Continued)



 Magnetic Resonance
Conditional¹

Ultraflex™ Single-Use Partially Covered Tracheobronchial Stent System – Distal Release*




Order Number	GTIN	Description	Expanded Stent O.D. (mm)	Expanded Stent Length (mm)	Covered Compressed Stent O.D. (Fr)**	Cover Length (mm)	Tip Max O.D. (mm)	Packaging
M00576520	08714729842323	Partially Covered Distal Release	8	40	16	25	4.1	Each
M00576530	08714729842330	Partially Covered Distal Release	10	30	18	15	4.1	Each
M00576540	08714729842347	Partially Covered Distal Release	10	40	18	25	4.1	Each
M00576550	08714729842354	Partially Covered Distal Release	12	30	20	15	4.1	Each
M00576560	08714729842361	Partially Covered Distal Release	12	40	20	25	4.1	Each
M00576570	08714729842378	Partially Covered Distal Release	14	30	24	15	4.1	Each
M00576580	08714729842385	Partially Covered Distal Release	14	40	24	25	4.1	Each
M00576590	08714729842392	Partially Covered Distal Release	14	60	24	45	4.1	Each
M00576600	08714729842408	Partially Covered Distal Release	14	80	24	65	4.1	Each
M00576610	08714729842415	Partially Covered Distal Release	16	40	25	25	5.3	Each
M00576620	08714729842422	Partially Covered Distal Release	16	60	25	45	5.3	Each
M00576630	08714729842439	Partially Covered Distal Release	16	80	25	65	5.3	Each
M00576640	08714729842446	Partially Covered Distal Release	18	40	27	25	5.3	Each
M00576650	08714729842453	Partially Covered Distal Release	18	60	27	45	5.3	Each
M00576660	08714729842460	Partially Covered Distal Release	18	80	27	65	5.3	Each
M00576670	08714729842477	Partially Covered Distal Release	20	40	27	25	5.3	Each
M00576680	08714729842484	Partially Covered Distal Release	20	60	27	45	5.3	Each
M00576690	08714729842491	Partially Covered Distal Release	20	80	27	65	5.3	Each

*Recommended Guidewire .035" Pulmonary Jagwire™ Pulmonary Guidewire and Amplatz Super Stiff Guidewire, Order #1517,  (See page 91).
Delivery System working length is 95cm.

**Data on File at Boston Scientific.

Pulmonary

Airway Stents (Continued)




Polyflex™ Single-Use Self-Expanding Silicone Airway Stent					
Order Number	GTIN	Stent ID (mm)	Length (mm)	Delivery System Diameter (mm)	Packaging
M00570000	08714729457138	8	20	7	Each
M00570010	08714729457145	8	30	7	Each
M00570020	08714729457152	10	20	8	Each
M00570030	08714729457169	10	30	8	Each
M00570040	08714729457176	10	40	8	Each
M00570050	08714729457183	10	50	8	Each
M00570060	08714729457190	12	20	9	Each
M00570070	08714729457206	12	30	9	Each
M00570080	08714729457213	12	40	9	Each
M00570090	08714729457220	12	50	9	Each
M00570100	08714729457237	14	20	9	Each
M00570110	08714729457244	14	30	9	Each
M00570120	08714729457251	14	40	9	Each
M00570130	08714729457268	14	50	9	Each
M00570140	08714729457275	14	60	9	Each
M00570150	08714729457282	16	30	10	Each
M00570160	08714729457299	16	40	10	Each
M00570170	08714729457305	16	50	10	Each
M00570180	08714729457312	16	60	10	Each
M00570190	08714729457329	16	70	10	Each
M00570200	08714729457336	18	30	11	Each
M00570210	08714729457343	18	40	11	Each
M00570220	08714729457350	18	50	11	Each
M00570230	08714729457367	18	60	11	Each
M00570240	08714729457374	18	70	11	Each
M00570250	08714729457381	18	80	11	Each
M00570260	08714729457398	20	40	12	Each
M00570270	08714729457404	20	50	12	Each
M00570280	08714729457411	20	60	12	Each
M00570290	08714729457428	20	70	12	Each
M00570300	08714729457435	20	80	12	Each
M00570310	08714729457442	22	50	13	Each
M00570320	08714729457459	22	60	13	Each
M00570330	08714729457466	22	80	13	Each



Pulmonary

Biopsy Forceps



Radial Jaw™ 4 Single-Use Pulmonary Biopsy Forceps



Order Number	GTIN	Description	Differentiator	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Packaging
M00515181	08714729792864	Single-Use Biopsy Forceps	Pulmonary Standard Capacity	1.8	100	2.0	Box 5
M00515182	08714729799337	Single-Use Biopsy Forceps	Pulmonary Standard Capacity	1.8	100	2.0	Box 20
M00515191	08714729792871	Single-Use Biopsy Forceps	Pulmonary Standard Capacity w/ Needle	1.8	100	2.0	Box 5
M00515192	08714729818083	Single-Use Biopsy Forceps	Pulmonary Standard Capacity w/ Needle	1.8	100	2.0	Box 20
M00515201	08714729793281	Single-Use Biopsy Forceps	Pulmonary Large Capacity	2.2	100	2.8	Box 5
M00515202	08714729818076	Single-Use Biopsy Forceps	Pulmonary Large Capacity	2.2	100	2.8	Box 20

Expect™ Pulmonary Endobronchial Ultrasound Transbronchial Aspiration Needle




Expect Pulmonary Endobronchial Ultrasound Transbronchial Aspiration Needle



Order Number	GTIN	Needle Size	Minimum Working Channel (mm)	Sheath O.D. (mm)	Packaging
M00558220	08714729861409	Expect Pulmonary 22Ga – Olympus scope compatible	2.0	1.6	Please contact your Boston Scientific Sales Representative for ordering information
M00558250	08714729861416	Expect Pulmonary 25Ga – Olympus scope compatible	2.0	1.4	Please contact your Boston Scientific Sales Representative for ordering information
M00558731	08714729893189	Expect Pulmonary Needle Adaptor – Olympus scope compatible	N/A	N/A	Please contact your Boston Scientific Sales Representative for ordering information

Transbronchial Aspiration Needles



eXcelon™ Transbronchial Aspiration Needle*



Order Number	GTIN	Description	Needle Gauge	Needle Length (mm)	Catheter Length (cm)	Sheath O.D. (mm)	Packaging
M00564101	08714729456452	eXcelon Transbronchial Aspiration Needle	19	15	130	1.8	Box 5
M00564111	08714729456469	eXcelon Transbronchial Aspiration Needle	20	15	130	1.8	Box 5
M00564121	08714729456476	eXcelon Transbronchial Aspiration Needle	21	15	130	1.8	Box 5

*Packaging also includes a 20ml syringe.



Section Table of Contents

Prev Page ◀ 89 ▶ Next Page

Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at 1.888.272.1001

Pulmonary


Cytology Brushes

Cellebri™ Single-Use Endoscopic Cytology Brush



Order Number	GTIN	Required Working Channel (mm)	Bristle O.D. (mm)	Sheath Length (cm)	Packaging
M00516001	08714729138983	2.0	1.0	140	Box 10
M00516011	08714729138990	2.0	1.5	140	Box 10
M00516071	08714729139089	2.0	1.9	100	Box 10
M00516151	08714729014003	2.0	1.9	150	Box 10

Airway Retrieval Baskets



Zero Tip™ Airway Retrieval Basket



Order Number	GTIN	Description	Sheath O.D. (mm)	Working Length (cm)	Basket Opening (mm)	Packaging
M00513200	08714729414995	Zero Tip Airway Retrieval Basket	0.8	120	12	Each
M00513210	08714729415008	Zero Tip Airway Retrieval Basket	1.0	120	16	Each



Pulmonary

Controlled Radial Expansion Pulmonary Balloon Dilatation Catheter



CRE™ Pulmonary Balloon Dilatation Catheter*



Order Number	GTIN	Balloon Inflated O.D.		Balloon Length (cm)	Inflation Pressures		Catheter Length (cm)	Packaging
		(mm)	(Fr)		(atm)	(kPa)		
M00550300	08714729456186	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	110	Each
M00550310	08714729456193	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	110	Each
M00550320	08714729456209	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	110	Each
M00550330	08714729456216	8-9-10	24-27-30	3.0	3-5.5-9	304-557-912	110	Each
M00550340	08714729456223	10-11-12	30-33-36	3.0	3-5-8	304-507-811	110	Each
M00550350	08714729456230	12-13.5-15	36-40.5-45	3.0	3-4.5-8	304-456-811	110	Each

*Recommended Inflation Device: Alliance™ II Integrated Inflation / Lithotripsy Device Order #5062, Alliance II Single-Use Gauge / Syringe Assembly Order #5060 (see below).


Pulmonary Guidewires

Pulmonary Jagwire™ Pulmonary Guidewire



Order Number	GTIN	Description	O.D. (in)	Length (cm)	Packaging
M00515171	08714729455813	Pulmonary Jagwire Pulmonary Guidewire	.035	180	Box 2
M00550090	08714729301929	Amplatz Super Stiff Guidewire	.038	260	Each

Inflation System



Alliance™ II Integrated Inflation / Lithotripsy Device and Alliance II Single-Use Syringe / Gauge Assembly

Order Number	GTIN	Description	Packaging
M00550620	08714729283881	Alliance II Integrated Inflation / Lithotripsy Device	Each
M00550601	08714729129332	Alliance II Single-Use Syringe / Gauge Assembly	Box 5



Section Table of Contents

Prev Page **91** Next Page

Main Table of Contents



For pricing information contact your Boston Scientific Sales Representative or Customer Service at **1.888.272.1001**