

# Dreamwire™

High Performance Guidewire



RX Biliary System™  
**A difference**  
you can feel.™

**Boston  
Scientific**

*Delivering what's next.™*

## DREAM ACCESS

Facilitates deep cannulation of targeted ducts\*

- Forms an Alpha Loop to facilitate effective access into the targeted duct
- 10cm hydrophilic Dream Tip™ End is designed for soft, smooth access and excellent tactile feel
- Endoglide™ Coating promotes smooth tracking and enhances tactile sensation



## DREAM EXCHANGE

Designed to improve procedural efficiency and easy RX Biliary device exchange

- 260cm and 450cm guidewire lengths available to help meet physician preference for short-wire or traditional over-the-wire technique
- Kink-resistant Nitinol standard and stiff shaft, engineered for wire control and manipulation capabilities
- Striped movement markers designed to indicate wire positioning during wire advancement and device exchange
- Side loading Torque Device designed for effective guidewire tip manipulation (Angled-tip only)



## ORDERING INFORMATION

Product Number	Description	Length (cm)	Dream Tip Style	Shaft	Diameter (in)	Diameter (mm)
M00556101	Dreamwire	260	Straight	Standard	0.035	0.89
M00556111	Dreamwire	260	Angled	Standard	0.035	0.89
M00556121	Dreamwire	260	Straight	Stiff	0.035	0.89
M00556131	Dreamwire	260	Angled	Stiff	0.035	0.89
M00556141	Dreamwire	450	Straight	Standard	0.035	0.89
M00556151	Dreamwire	450	Angled	Standard	0.035	0.89
M00556161	Dreamwire	450	Straight	Stiff	0.035	0.89
M00556171	Dreamwire	450	Angled	Stiff	0.035	0.89

Packaged 2 per box

\*M.J. Bourke, G. Costamagna, M.L. Freeman: *Biliary cannulation during endoscopic retrograde cholangiopancreatography: core techniques and recent innovations.* Endoscopy 2009; 41: 612-617.

**European Distribution Centre – The Netherlands**  
T: +31 45 54 67 700 F: +31 45 54 67 800  
**European Headquarters – Paris**  
T: +33 1 57 66 80 00 F: +33 1 57 66 84 99  
**Argentina (Freephone)**  
T: 0800 555 2678 F: +5411 4896 8517  
**Australia/New Zealand**  
T: +61 2 8063 8100 F: +61 2 9330 1404  
**Austria**  
T: +43 1 60 810 F: +43 1 60 810 60  
**Belgium (Freephone)**  
T: 0800 94 494 F: 0800 93 343  
**Brazil**  
T: +55 11 5502 8500 F: +55 11 5103 2212  
**Canada**  
T: +1 888 359 9691 F: +1 888 575 7396  
**Chile**  
T: +562 445 4904 F: +562 445 4915  
**China – Beijing**  
T: +86 10 8525 1588 F: +86 10 8525 1566  
**China – Guangzhou**  
T: +86 20 8767 9791 F: +86 20 8767 9789  
**China – Shanghai**  
T: +86 21 6141 5959 F: +86 21 6141 5900

**Colombia**  
T: +57 1 629 5045 F: +57 1 612 4761  
**Czech Republic**  
T: +420 296 331 901 F: +420 296 331 935  
**Denmark (Freephone)**  
T: 80 30 80 02 F: 80 30 80 05  
**Eire (Freephone)**  
T: 1800 882 969 F: 1800 882 968  
**Finland**  
T: +358 20 762 8882 F: +358 20 762 8883  
**France**  
T: +33 1 39 30 49 00 F: +33 1 39 30 49 01  
**Germany (Freephone)**  
T: 08000 723300 F: 08000 723319  
**Greece**  
T: +30 210 9542 300 F: +30 210 9542 310  
**Hong Kong**  
T: +852 2960 7100 F: +852 2563 5276  
**Hungary**  
T: +36 1 456 30 40 F: +36 1 456 30 41  
**India – Bangalore**  
T: +91 80 2212 4928/9 F: +91 80 2207 5153  
**India – Chennai**  
T: +91 44 2220 1879 F: +91 44 2220 0716  
**India – Delhi**  
T: +91 11 4243 2222 F: +91 11 2610 0808

**India – Mumbai**  
T: +91 22 4030 9165 F: +91 22 4040 9199  
**Italy**  
T: +39 010 60 60 1 F: +39 010 60 60 200  
**Korea**  
T: +82 2 3476 2121 F: +82 2 3476 1776  
**Mexico**  
T: +52 55 5687 63 90 F: +52 55 5687 62 28  
**Middle East/Gulf/North Africa**  
T: +961 1 805 410 F: +961 1 805 445  
**The Netherlands**  
T: +31 30 602 55 44 F: +31 30 602 55 05  
**Norway (Freephone)**  
T: 800 104 04 F: 800 101 90  
**Poland**  
T: +48 22 435 14 14 F: +48 22 435 14 10  
**Portugal**  
T: +351 1 381 25 40 F: +351 21 381 25 58  
**South Africa**  
T: +27 11 840 8600 F: +27 11 463 6077  
**South East Asia – Malaysia**  
T: +60 3 2283 3813 F: +60 3 2284 3813  
**South East Asia – Philippines**  
T: +63 2 687 6994 F: +63 2 687 3047  
**South East Asia – Singapore**  
T: +65 6418 8888 F: +65 6418 8899

**South East Asia – Thailand**  
T: +66 2 6380 100 F: +66 2 6380 400  
**South East Europe**  
T: +30 210 9542 300 F: +30 210 9542 310  
**Spain – Barcelona**  
T: +34 93 444 72 00 F: +34 93 405 90 45  
**Spain – Madrid**  
T: +34 90 11 12 15 F: +34 91 319 50 03  
**Sweden**  
T: +46 42 25 69 00 F: +46 42 25 69 69  
**Switzerland**  
T: +41 32 626 57 00 F: +41 32 626 57 01  
**Taiwan**  
T: +886 2 2747 7278 F: +886 2 2747 7270  
**Turkey – Istanbul**  
T: +90 216 464 36 66 F: +90 216 464 36 67  
**UK (Freephone)**  
T: 0844 800 4512 F: 0844 800 4513  
**Uruguay**  
T: +59 82 900 6212 F: +59 82 900 6212  
**Venezuela**  
T: +58 212 959 6275 F: +58 212 959 5328

**Boston Scientific**

*Delivering what's next.™*

www.bostonscientific.com  
www.bostonscientific-international.com