



Dilation

CRE™ Balloon Dilators

CRE Single-Use Fixed Wire Esophageal Balloon Dilators	55
CRE Single-Use Wireguided Esophageal / Pyloric / Biliary Balloon Dilators	55
CRE Single-Use Wireguided Esophageal / Pyloric / Colonic / Biliary Balloon Dilators	56
CRE Single-Use Wireguided Esophageal / Pyloric Balloon Dilators	.59
CRE Single-Use Wireguided Esophageal / Pyloric / Colonic Balloon Dilators	59

MaxForce TTS™ Esophageal High Performance Balloon Dilatation Catheter

MaxForce TTS Esophageal High Performance Balloon Dilatation Catheter	57
---	----

Rigiflex™ II Single-Use Achalasia Balloon Dilators

Rigiflex II Single-Use Achalasia Balloon Dilators	57
---	----

Inflation Systems


Alliance™ II Inflation System	58
-------------------------------------	----

Dilation Accessories

Inflation / Deflation Devices and Accessories	58
---	----

Dilation

CRE™ Balloon Dilators





CRE Single-Use Fixed Wire Esophageal Balloon Dilators



Order Number	GTIN	Balloon Inflated O.D.		Balloon Length (cm)	Inflation Pressures		Minimum Working Channel (mm)	Catheter Size (Fr)	Working Length (cm)	Packaging
		(mm)	(Fr)		(ATM)	(kPa)				
M00558330	08714729201984	6-7-8	18-21-24	8	3-6-10	304-608-1013	2.8	6	180	Each
M00558340	08714729201991	8-9-10	24-27-30	8	3-5.5-9	304-557-912	2.8	6	180	Each
M00558350	08714729202004	10-11-12	30-33-36	8	3-5-8	304-507-811	2.8	6	180	Each
M00558360	08714729195986	12-13.5-15	36-40.5-45	8	3-4.5-8	304-456-811	2.8	6	180	Each
M00558370	08714729195993	15-16.5-18	45-49.5-54	8	3-4.5-7	304-456-709	2.8	6	180	Each
M00558380	08714729202011	18-19-20	54-57-60	8	3-4.5-6	304-456-608	2.8	6	180	Each
M00558361	08714729339298	12-13.5-15	36-40.5-45	8	3-4.5-8	304-456-811	2.8	6	180	Box 5
M00558371	08714729339304	15-16.5-18	45-49.5-54	8	3-4.5-7	304-456-709	2.8	6	180	Box 5
M00558381	08714729339311	18-19-20	54-57-60	8	3-4.5-6	304-456-608	2.8	6	180	Box 5

Recommended Inflation Device: Alliance™ II Integrated Inflation System Handle Order #5062, Sterile Syringe Order #5060, (See page 58).



CRE Single-Use Wireguided Esophageal / Pyloric / Biliary Balloon Dilators



Order Number	GTIN	Balloon Inflated O.D.		Balloon Length (cm)	Inflation Pressures		Minimum Working Channel (mm)	Catheter Size (Fr)	Working Length (cm)	Packaging
		(mm)	(Fr)		(ATM)	(kPa)				
M00558600	08714729797487	6-7-8	18-21-24	5.5	3-6-10	304-608-1013	2.8 / 3.2* / 3.2 (B1)	7.5	180	Each
M00558610	08714729797494	8-9-10	24-27-30	5.5	3-5.5-9	304-557-912	2.8 / 3.2* / 3.2 (B1)	7.5	180	Each
M00558620	08714729797500	10-11-12	30-33-36	5.5	3-5-8	304-507-811	2.8 / 3.2* / 3.2 (B1)	7.5	180	Each
M00558630	08714729797517	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	2.8 / 3.2* / 3.2 (B1)	7.5	180	Each
M00558640	08714729797524	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	2.8 / 3.2* / 3.7 / 4.2 (B2)	7.5	180	Each
M00558650**	08714729797531	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	2.8 / 3.2* / 3.7 / 4.2 (B2)	7.5	180	Each

Recommended Inflation Device: Alliance II Integrated Inflation System Handle Order #5062, Sterile Syringe Order #5060, (See page 58).

Guidewire included and preloaded in balloon catheter.

*Requires a 2.8mm Olympus™ Brand Endoscope channel or a 3.2mm Pentax™ or Fujinon™ Brand Endoscope channel.


**Esophageal and Biliary application only.

(B1) = Indication requires ≥ 3.2mm Duodenoscope channel.

(B2) = Indication requires ≥ 4.2mm Duodenoscope channel.

Dilation

CRE™ Balloon Dilators (Continued)



CRE Single-Use Wireguided Esophageal / Pyloric / Colonic / Biliary Balloon Dilators



Order Number	GTIN	Balloon Inflated O.D.		Balloon Length (cm)	Inflation Pressures		Minimum Working Channel (mm)	Catheter Size (Fr)	Working Length (cm)	Packaging
		(mm)	(Fr)		(ATM)	(kPa)				
M00558660	08714729797548	6-7-8	18-21-24	5.5	3-6-10	304-608-1013	2.8 / 3.2* / 3.2 (B1)	7.5	240	Each
M00558670	08714729797555	8-9-10	24-27-30	5.5	3-5.5-9	304-557-912	2.8 / 3.2* / 3.2 (B1)	7.5	240	Each
M00558680	08714729797562	10-11-12	30-33-36	5.5	3-5-8	304-507-811	2.8 / 3.2* / 3.2 (B1)	7.5	240	Each
M00558690	08714729797579	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	2.8 / 3.2* / 3.2 (B1)	7.5	240	Each
M00558700	08714729797586	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	2.8 / 3.2* / 3.7 / 4.2 (B2)	7.5	240	Each
M00558710**	08714729797593	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	2.8 / 3.2* / 3.7 / 4.2 (B2)	7.5	240	Each

Recommended Inflation Device: Alliance™ II Integrated Inflation System Handle Order #5062, Sterile Syringe Order #5060, (See page 58).

Guidewire included and preloaded in balloon catheter.

*Requires a 2.8mm Olympus™ Brand Endoscope channel or a 3.2mm Pentax™ or Fujinon™ Brand Endoscope channel.


**Esophageal, Colonic and Biliary application only.

(B1) = Indication requires ≥ 3.2mm Duodenoscope channel.

(B2) = Indication requires ≥ 4.2mm Duodenoscope channel.

Dilation

MaxForce TTS™ Esophageal High Performance Balloon Dilatation Catheter




MaxForce TTS Esophageal High Performance Balloon Dilatation Catheter



Order Number	GTIN	Balloon Inflated O.D.		Balloon Length (cm)	Inflation Pressures		Minimum Working Channel (mm)	Catheter Size		Working Length (cm)	Packaging
		(mm)	(Fr)		(ATM)	(kPa)		(mm)	(Fr)		
M00589060	08714729099581	6	18	6	10	1013	2.8	2.0	6	180	Each
M00589080	08714729099598	8	24	6	10	1013	2.8	2.0	6	180	Each
M00589100	08714729099895	10	30	6	7	709	2.8	2.0	6	180	Each
M00589120	08714729099901	12	36	6	7	709	2.8	2.0	6	180	Each
M00589140	08714729080602	14	42	6	7	709	2.8	2.0	6	180	Each
M00589150	08714729099918	15	45	6	5	507	2.8	2.0	6	180	Each
M00589160	08714729080619	16	48	6	5	507	2.8	2.0	6	180	Each
M00589180	08714729099925	18	54	6	5	507	2.8	2.0	6	180	Each

Recommended Inflation Device: Alliance™ II Integrated Inflation System Handle Order #5062, Sterile Syringe Order #5060, (See page 58).

Rigiflex™ II Single-Use Achalasia Balloon Dilators



Rigiflex II Single-Use Achalasia Balloon Dilators*




Order Number	GTIN	Balloon Inflated O.D.		Balloon Length (cm)	Inflation Pressures		Working Length (cm)	Catheter Size		Recommended Guidewire** (in)	Packaging
		(mm)	(Fr)		(psi)	(kPa)		(mm)	(Fr)		
M00554500	08714729719106	30	90	10	20	138	90	4.6	14	.038	Each
M00554510	08714729719113	35	105	10	20	138	90	4.6	14	.038	Each
M00554520	08714729719120	40	120	10	20	138	90	4.6	14	.038	Each

*Product Shipped Sterile. Recommended use with Achalasia Pneumatic Hand Pump and Monitor Order #5320, (See page 58).

**Guidewire included with Rigiflex II Balloon Dilatation Catheter.

Dilation

Inflation Systems



Alliance™ II Inflation System

Order Number	GTIN	Description	Packaging
M00550620	08714729283881	Inflation Handle	Each
M00550601	08714729129332	⊗ Single-Use Syringe / Gauge Assembly	Box 5


Dilation Accessories

Inflation / Deflation Devices and Accessories

Order Number	GTIN	Description	Packaging
M00553200	08714729195450	Achalasia Pneumatic Hand Pump and Monitor (not shown)	Each

Dilation

CRE™ Balloon Dilators



CRE Single-Use Wireguided Esophageal / Pyloric Balloon Dilators



Order Number	GTIN	Balloon Inflated O.D.		Balloon Length (cm)	Inflation Pressures		Minimum Working Channel (mm)	Catheter Size (Fr)	Working Length (cm)	Packaging
		(mm)	(Fr)		(ATM)	(kPa)				
M00558390	08714729339328	6-7-8	18-21-24	5.5	3-6-10	304-608-1013	2.8	7.5	180	Each
M00558400	08714729339335	8-9-10	24-27-30	5.5	3-5.5-9	304-557-912	2.8	7.5	180	Each
M00558410	08714729339342	10-11-12	30-33-36	5.5	3-5-8	304-507-811	2.8	7.5	180	Each
M00558420	08714729339359	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	2.8	7.5	180	Each
M00558430	08714729339366	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	2.8 / 3.2*	7.5	180	Each
M00558440**	08714729339373	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	2.8 / 3.2*	7.5	180	Each

Recommended Inflation Device: Alliance™ II Integrated Inflation System Handle Order #5062, Sterile Syringe Order #5060, (See page 58).

Guidewire included and preloaded in balloon catheter.

*Requires a 2.8mm Olympus™ Brand Endoscope channel or a 3.2mm Pentax™ or Fujinon™ Brand Endoscope channel.

**Esophageal application only.

CRE Single-Use Wireguided Esophageal / Pyloric / Colonic Balloon Dilators



Order Number	GTIN	Balloon Inflated O.D.		Balloon Length (cm)	Inflation Pressures		Minimum Working Channel (mm)	Catheter Size (Fr)	Working Length (cm)	Packaging
		(mm)	(Fr)		(ATM)	(kPa)				
M00558450	08714729339380	6-7-8	18-21-24	5.5	3-6-10	304-608-1013	2.8 / 3.2*	7.5	240	Each
M00558460	08714729339397	8-9-10	24-27-30	5.5	3-5.5-9	304-557-912	2.8 / 3.2*	7.5	240	Each
M00558470	08714729339403	10-11-12	30-33-36	5.5	3-5-8	304-507-811	2.8 / 3.2*	7.5	240	Each
M00558480	08714729339410	12-13.5-15	36-40.5-45	5.5	3-4.5-8	304-456-811	2.8 / 3.2*	7.5	240	Each
M00558490	08714729339427	15-16.5-18	45-49.5-54	5.5	3-4.5-7	304-456-709	2.8 / 3.2*	7.5	240	Each
M00558500**	08714729339434	18-19-20	54-57-60	5.5	3-4.5-6	304-456-608	2.8 / 3.2*	7.5	240	Each

Recommended Inflation Device: Alliance II Integrated Inflation System Handle Order #5062, Sterile Syringe Order #5060, (See page 58).

Guidewire included and preloaded in balloon catheter.

*Requires a 2.8mm Olympus Brand Endoscope channel or a 3.2mm Pentax or Fujinon Brand Endoscope channel.

**Esophageal or Colonic application only.

