




SpyGlass™ DS II Direct Visualization System featuring Sight Shield Technology

Accessory Information

Boston
Scientific

The **SpyGlass DS II Direct Visualization System** includes an expansive portfolio of devices that deliver innovative diagnostic and therapeutic approaches to biliary stone and stricture management to help deliver successful outcomes in complex procedures. These devices include:



Autolith™ Touch EHL System

The Autolith Touch 1.9Fr, 375cm EHL Probe and Autolith Touch EHL Generator are available to help physicians manage large biliary stones.

Slimline SIS GI™ Laser Fiber and VersaPulse® Holmium Lasers

Designed to help navigate difficult anatomies in the management of biliary stones and benign or malignant neoplasms.

SpyBite™ Max Biopsy Forceps

Helps enable greater tissue acquisition* which may result in faster more accurate diagnosis of malignancies.¹

SpyGlass™ Retrieval Basket

A durable 4-wire nitinol basket designed to grasp biliary and pancreatic stones while maintaining its shape in a difficult anatomy.

SpyGlass™ Retrieval Snare

Designed to efficiently capture and remove foreign bodies – such as migrated plastic stents – in the biliary and pancreatic ducts.

*The testing was performed by or on behalf of BSC. Data on file. Bench test results may not necessarily be indicative of clinical performance.

SpyGlass™ DS System

Order Number	Product Description	
M00546650	SpyGlass™ DS Digital Controller	Each
M00546610	SpyScope™ DS II Access & Delivery Catheter	Each

SpyGlass™ DS Accessory Devices

Order Number	Product Description	Cable Diameter (in/mm)	Jaw Outer Diameter (mm)	Working Length (cm)	Required Endoscope Working Channel (mm)
M00546470	SpyBite™ Max Biopsy Forceps	0.039 / 1.0	1	286	1.2
M00546270	SpyBite™ Biopsy Forceps	0.039 / 1.0	1	286	1.2
M00546550	SpyGlass™ Retrieval Basket	0.039 / 1.0	N/A	286	1.2
M00546560	SpyGlass™ Retrieval Snare	0.039 / 1.0	N/A	286	1.2

Biliary EHL Probe

Order Number	Product Description
M00546620	1.9 Fr 375cm Biliary EHL Probe

Autolith™ Touch System

Order Number	Product Description
M005466800	Autolith™ Touch EHL Generator

Autolith™ Touch System Accessories

Order Number	Product Description
M00546750	Autolith™ Touch Extender Cable
M00546760	Autolith™ Touch Foot Pedal

Slimline SIS GI™ Laser Fiber and VersaPulse® Holmium Lasers

Order Number	Product Description	Units
0644-007-01	1.8Fr SlimLine GITM Probe	Box of 3
GA-0006000	VersaPulse P20 Holmium Laser*	Each

*Includes generator, footswitch, fiber inspection scope, two laser safety glasses, one laser safety goggles, one remote interlock connector and one power cable.

1. Shah et al. Performance of a full disposable, digital, single-operator cholangiopancreatoinoscope. Endoscopy 2017; 49: 651-658

All photographs taken by Boston Scientific.

All trademarks are the property of their respective owners.

CAUTION: The law restricts these devices to sale by or on the order of a physician.

Indications, contraindications, warnings, and instructions for use can be found in the product labelling supplied with each device or at www.IFU-BSCI.com.

Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. This material not intended for use in France. RX Only.

2023 Copyright © Boston Scientific Corporation or its affiliates. All rights reserved.

ENDO-1519003-AA

**Boston
Scientific**
Advancing science for life™

Boston Scientific Corporation
300 Boston Scientific Way
Marlborough, MA 01752-1234
www.BostonScientific.com